Content

2 Welcome

5

9

Part 1: General

4 1	Ove	r 100	vears o	of ex	perience
T .	OVE		y Cui 3 C	, ,	Perience

- Past, Present and Future
- Mallinckrodt Baker's evolution over the last decade

8 2. Doing business with us

- Innovation: R&D, Six Sigma, S&OP
- Collaboration: Distribution Channel, Distributors
- Support: Service, Customer Service, Product and Application
 Specialists, Certificate of Analysis, MSDS, J.T.Baker Technical Library,
 Solv-It Center Knowledge Base, Product literature
- Quality: ISO 9001: 2000, ISO 14001, Responsible Care, CE mark, GMP, J.T.Baker Label, Packaging options

Part 2: The businesses

14 3. Analytical, Chemical reagents and Chromatography

- Chromatography products
- Reagents for trace element analysis
- Solvent applications
- CYCLE-TAINER High Purity Solvent Delivery System
- Reagents for pH Measurements, Titrimetry and Water Determination according to Karl Fischer
- Reagents for Molecular Biology and Biotechnology
- 28 4. Products for Hematology and Histopathology
- 32 5. Microelectronic materials
- 36 6. Pharmaceutical products

Part 3: A - Z

40 7. Product Specifications

Analytical reagents, Histopathology reagents,
 DNA synthesis reagents, Pharmaceutical products, Microelectronic materials,
 Biopharmaceutical products


• Chromatography products



Part 4: Miscellaneous

- 492 Trademarks
- 494 Risk and Safety phrases
- 496 General conditions
- 500 Catalogue index: Product number, CAS number and EINECS number
- cover Periodic table cover – A - Z Definitions





'Part of a pure process'

Dear Customer,

It is my pleasure to give you a brief introduction to Mallinckrodt Baker and our new J.T.Baker catalogue. The catalogue is designed to focus on your needs, allowing our quality products to be found quickly and easily for your applications.

Mallinckrodt Baker is a global organisation, with European headquarters in Deventer, The Netherlands, and employees working at five production facilities worldwide.

In Europe, a dedicated staff of people produces, researches and sells the highest technology chemistry on a daily basis. We offer high purity chemistry for analytical, biopharmaceutical, clinical and microelectronic applications.

At Mallinckrodt Baker we understand your needs. We employ specialist contacts that understand your business and offer technical and sales support on a daily basis. We are a customer focused organisation, where service, flexibility, innovation, research and development are fundamental to company growth.

We have selected the best dealers worldwide to strengthen our business with you. Providing an optimum service through the supply channel, is principal. We make every effort to have products available when you need them. No matter which communication or supply channel you choose, direct or via distribution, we do our utmost to assure the same high level of support.

Our mission is to deliver the best chemistry and be a leading global provider of high purity chemistry for research, analysis and demanding technology markets.

At Mallinckrodt Baker we continually demand the best for the customer. Should you have any questions or suggestions, please do not hesitate to contact us by e-mail at jtbaker.nl@emea.tycohealthcare.com or contact us through our Customer Service department.

Yours sincerely,



P.J.J. de Groot

Managing Director

Mallinckrodt Baker Europe



Past, Present and Future

In 1904, John Townsend Baker started a company to make the purest chemicals on earth. Along the way, he and his company achieved some major milestones.

A chemical manufacturer founded on one word: purity

Baker took every penny he had to invest in a new venture and founded the J.T.Baker Chemical Company.

Labelling that would influence an entire industry

To give his company a competitive difference, Baker focused on the label. He thought that if impurities were accurately listed on the label, it would benefit chemists. Accurate, informative labelling became an industry standard and is still considered a key J.T.Baker strength.

Innovative insight

Over the next thirty years, plant capacity steadily expanded. Baker had another innovative insight: pure laboratory chemicals could be adapted for manufacturing use. This was the genesis of our Beaker to Bulk packaging. In the early 1940s, the company was converted to wartime production. Penicillin, pesticides, chemicals for ammunition, film, x-ray plates and batteries were produced.

About Mallinckrodt Baker

Mallinckrodt was founded in 1867. In 1995,
Mallinckrodt Chemical acquired J.T.Baker to
create Mallinckrodt Baker, Inc. Tyco International
purchased Mallinckrodt, Inc., in October 2000.
We became one of the key businesses of
Tyco Healthcare. Today we represent both the
Mallinckrodt brand and the J.T.Baker brand in the
global market.

Both brands are now to be found in laboratories, pharmaceutical processes, clinical diagnostics laboratories and microelectronic processes around the world.

Global markets

Today, Mallinckrodt Baker employees work at five ISO 9000-certified manufacturing facilities located in the United States, Mexico, Europe and Malaysia, supplying more than 10,000 products worldwide.

Since its founding, the company has always sought markets that demand the highest purity. To that end, Mallinckrodt Baker, Inc has also established a strong presence as a supplier to the high-growth and high-tech markets of semiconductor manufacturing, biopharmaceutical production and clinical diagnostics.



At the time of Baker's retirement, the company's annual sales were a little over \$1 million and more than half of its business was in laboratory reagents and the rest in pharmaceutical and fine chemicals. About 100 employees worked at a single plant that supplied more than 1600 chemicals.

Mallinckrodt Baker's evolution over the last decade

1994

- The J.T.Baker CYCLE-TAINER solvent delivery system is introduced.
- Mallinckrodt Baker Europe achieves ISO 9001 certification.

1996

- BAKERBOND Speedisk solid phase extraction disk technology is introduced.
- OneShot SPE cartridge for individual samples is introduced.
- Histology and cytology reagents are developed and introduced.

1998

- 18-Litre CYCLE-TAINER solvent delivery system is now available.
- J.T.Baker guanidine hydrochloride in Flowmor containers is now available.
- Cyanide-free hematology reagents are developed for veterinary blood cell counting. Patented.
- Xylene/toluene free mounting medium for coverslipping in histopathology and cytopathology is developed. Our invention is patented.
- HYDRA-POINT pyridine free volumetric reagents are introduced.

1994

1995

1998

1999

1996

1997

1997

- J.T.Baker USP/NF test reagents are introduced.
- Low water THF introduced.
- New BAKERBOND Speedisk oil & grease extraction disks are introduced.

1999

- New Global Pharma line includes testing for international pharmacopoeia specifications.
- Multicompendial tested products help meet requirements in Europe and Asia.
- J.T.Baker BakerDRY Low Water Solvents are introduced.

1995

• J.T.Baker is acquired by Mallinckrodt. The company name is changed from J.T.Baker to Mallinckrodt Baker Inc. The J.T.Baker brand is maintained as a premium brand.

2000

- Mallinckrodt Inc. is acquired by Tyco International and becomes part of The Tyco Healthcare division.
- Mallinckrodt Baker Inc. is maintained as a company.
- Comprehensive amino acid line introduced.

2002

- Reagents for oligonucleotide synthesis are introduced.
- HYDRA-POINT pyridine- free coulometric reagents are developed and introduced Unique formulations are patent pending.

2004

- Versatile photoresist stripper and residue remover (CLK-288 cleaner) is introduced.
- Rapid Stat controls and quality assessment renewed.

2000 2001



 The new building of the European headquarters in Deventer, The Netherlands, in use.

2003

2002

- J.T.Baker special LC/MS solvents and reagents are developed and introduced.
- Mallinckrodt Baker Europe is ready with IVD CE marking.
- Mallinckrodt Baker Europe achieves ISO 14001 certification.

2005


2004

2003

Patent for xylene/toluene free mounting medium is filed and granted. Product is available as UltraKitt.

2005

- Mallinckrodt Baker Europe installed and validated a production tank for manufacturing of cGMP certified solutions. Begin multicompendial sodium hydroxide production.
- R&D group at Mallinckrodt Baker Europe developed a new SPE carbon column.



At Mallinckrodt Baker we understand your needs. We find it fundamental to introduce new products that embrace your existing research and production capabilities.

These competences translate into the four cornerstones of doing business with us:

Innovation, Collaboration, Support, Quality.

Innovation

At Mallinckrodt Baker innovation is fundamental to support our customers' increasing demands. For that reason, *R&D* is central in our business strategies. Our timeline, on the previous page, highlights Mallinckrodt Baker's innovation and vitality over the last decade.

Besides R&D, systematic improvements to our processes are made on a continual basis. Mallinckrodt Baker utilises the techniques of *Six Sigma project*

management and S&OP, Sales and Operational Planning.

Six Sigma projects have led to tighter specifications, more consistent results and innovative new products. Our goal is to provide you with products and services of the highest purity and value.

S&OP is used as our major business planning model to anticipate market changes and meet future customer demands.

Collaboration

We have selected the best *distribution channel* for you. On a daily basis we are in direct contact with you or we serve you through our official *distributors*. Our distribution partners are located throughout Europe, the Middle-East and Africa. An up-to-date list of

distributors can be found at www.jtbaker.com/
europe/distrib. For optimal service and more detailed information you can always contact our Customer
Service department. Address details on the rear cover of this cataloque.

Support

Mallinckrodt Baker's service includes on-time delivery and extensive customer support. Support incorporates dedicated centralised European Customer Service, application support by product and technical specialists and our up-to-date information programs on www.jtbaker.com

Customer Service

Our colleagues from *Customer Service* are always willing to help in any situation. A recent customer satisfaction survey describes: "Excellent performance, with a professional and responsible attitude, they are

customer, market and product oriented." A native speaking Customer Service department, with employees from throughout Europe, speak your native language for your enquiries.

Product and Application Specialists

Our *Product and Application Specialists* are professionals in the businesses we offer to you, answering specific inquiries and application issues, with the highest level of technical knowledge. Ask our Customer Service department to connect you to our specialists.

Website

On our *Website* <u>www.jtbaker.com</u> you have continual access to several of our services:


- Certificate of Analysis
 The Certificate of Analysis facilitates your use of incoming quality control.
- MSDS stands for Material Safety Data Sheets.
 We know the importance of the regulatory information regarding your raw materials. You can search on-line in our database of current MSDS information. Search by Product Number, by Product name or by CAS Number.
- J.T.Baker Technical Library
 Inside this digital library you find a large body of information about our products and applications using our products. You can either browse by Market Area and Document Type, or search documents by key words.
- Solv-IT Center Knowledge Base
 If you have a question about any of our chemicals,
 enter our digital Solv-IT Center Knowledge Base at
 www.jtbaker.com/solvit.

Browse or search the database for answers to frequently asked questions. If you cannot find the answer, a Mallinckrodt Baker representative will promptly supply one.



• Product literature

In addition to all the information you find in this catalogue and the website, we compile a large number of brochures each year specific to products from our different product groups. Visit the website for a summary of all available literature, which can be ordered on line.



Quality

Fundamental to Mallinckrodt Baker is quality. Offering the highest quality, we have established best practices:

ISO 9001: 2000

Certification guarantees continuous total quality.


ISO 14001

Certification highlights awareness and care for the environment









Responsible Care

Responsible Care is a worldwide initiative of the chemical industry for a continuous improvement of efforts on environment, safety and health.

At Mallinckrodt Baker these aspects are integrated in all our activities.

CE mark

The CE mark appears on all products from Mallinckrodt Baker that are registered in the IVD Directive 98/79/EC (In Vitro Diagnostic Medical Devices).

GMP Good Manufacturing Practice

GMP Good Manufacturing Practice encompasses the documented, validated procedures and quality systems applicable to the supply of our Bulk Pharmaceutical excipients, from one of our FDA-registered facilities. It is a guarantee to our customers that all of our Pharmaceutical grade materials are produced under consistent conditions according to validated, controlled processes, with complete traceability. Currently, our US manufacturing facilities are FDA-registered and we anticipate the achievement of FDA-registration for our European facility by the middle of 2006.

The J.T.Baker Label

Our approach to total quality ensures: That the product inside the package meets or exceeds the information stated on the package. On our label you will find all the required information.

- A Danger Class boxes.
- Information on storage, shipping, et cetera.
- С Product Name, multi-language.
- D Quality definition.
- Ε Lot Number. For traceability purposes.
- F Product number.
- G Pack size.
- H R/S Sentences. In 10 languages.
- Chemical formula.
- Molecular Mass.
- Certificate of Analysis/Specifications.
- Use Before Date.
- M UN Code.
- Production Batch Number. For traceability purposes.
- O Application

Packaging options

The packaging options of J.T.Baker products offer purity, safety, economical control of inventory and the assurance of consistent quality as you scale up from R&D through production.

Acids



Glass, narrow mouth 500 ml, 1, 2.5, 4 L.



<u>Plastic, narrow mouth</u> 500 ml, 1, 2.5, 5 L High Density Polyethylene (HDPE).



<u>Plastic, narrow mouth</u> 25, 200 L HDPE drums.



<u>Steel/HDPE, Combi-drum</u> 25, 200 L Combi Steel outside, HDPE inside drums.

Salts



Glass, narrow mouth

15, 30 and 60 ml containing 100 mg
up to 100 grams of dry salt.



Glass, wide mouth
100, 250, 500 ml, 2.5 L,
containing up to 5 kg of salt.



Plastic, wide mouth
100, 250, 500 ml, 1, 2.5 and 5 L
HDPE containing regular salts from 100 gram up to 5 kg.



Plastic, drum or pail
5 to 50 kg. The inside packaging
consists of a resealable Low Density
Polyethylene (LDPE) bag.



<u>Bag-in-Box</u>25 to 50 kg. Corrugated cardboard box contains an inner LDPE bag.

Solutions



<u>Ampoules</u>

DILUT-IT liquid concentrates of acids, bases and buffers are available in LDPE or glass ampoules.



Glass, narrow mouth
1 and 2.5 L.

Solutions



Plastic, narrow mouth
500 ml, 1, 2.5 and 5 L HDPE.



<u>Plastic, narrow mouth</u> 5, 10 and 20 L Jerrycan. Easy-to-carry HDPE cans.



Polycubes

5, 10 and 20 L polycube has an LDPE bag inside and a carton over pack.



Intermediate Bulk Container
1000 L for process solutions.

Solvents



<u>Glass, narrow mouth</u> 100, 250, 500 ml, 1 and 2.5 L.



<u>Plastic, narrow mouth</u> 1, 2.5, 5 L HDPE.



Aluminium, narrow mouth
1 and 5 L EcoTainer.



<u>Plastic, narrow mouth</u>
5, 10 and 20 L Jerrycan. Easy-to-carry
HDPE cans.



<u>Steel, narrow mouth</u>25 and 200 L Tri-Sure closureAnalytical reagent grade.




Returnable Container

10, 30 and 200 L for Analytical reagent grade and HPLC quality.



<u>CYCLE-TAINER System</u>
High Purity Solvent Delivery System
20, 50, 200 and 1400 L for high quality
HPLC and organic residue solvents.



Mallinckrodt Baker is recognised as one of the world's leading high purity brands of analytical reagents. 'Actual Lot Analysis' as proof of product purity is still valuable today. However in the past decade a lot has changed in analytical and research laboratories.

Highest technological standards

Traditional and classical methods are replaced by modern instrumentation, such as LC, GC and ICP with MS as a powerful detection principle. In other fields such as pH measurement, titrimetry and spectroscopy, equipment is modernised to the highest technological standards.

What we offer you

In this catalogue we offer you innovative, high quality and reliable products, comprising a broad chemical assortment. In the pharmaceutical, food, chemical, environmental and water industries our analytical and chemical competencies will meet your analytical and chemical needs. We have qualified products for monitoring and research laboratories in the previously mentioned industries. Besides reagents and chemicals for analytical and research applications, we can support you with chemicals in your preparative and process applications.



J.T.Baker analytical and chemical reagents with the J.T.Baker brand at a glance:

- Broad chemical assortment
- Innovative and state-of-the-art products
- High and consistent quality
- For monitoring and research laboratories and for preparative and process applications

Product assortment

The analytical and chemical reagent portfolio of Mallinckrodt Baker:

Salts

Acids

Solvents

Chromatography

Titrimetry

Calibration standards

Buffers

Molecular biology

Test strips

Safety products

Strategic contract manufacturing

General Quality definitions

In summary our product quality is described by:

Pure and Ultrapure Qualities

Chemicals and reagents, specially developed for specific analytical applications. Examples are

- ULTREX II
- BAKER INSTRA ANALYZED
- BAKER ANALYZED HPLC
- BAKER ANALYZED LC/MS
- ULTRA RESI ANALYZED
- BAKER BIO ANALYZED
- BakerDRY

Analytical Reagent Quality

Chemicals and reagents, for various analytical methods.

- BAKER ANALYZED may also include:
 - ACS specifications
 - ISO specifications
 - reagent specifications for testing pharmacopoeia monographs

Technical Quality

Technical Quality for production and general laboratory applications

• BAKER

Product qualities versus application

Our extensive product program is applicable to almost any analytical and process application.

See table 3a.

Table 3a Product qualities versus application

Product	Application
Salts	
ULTREX II	For trace element analysis Elements ppt level
BAKER ANALYZED	For preparation of solutions and buffers for critical analytical applications
BAKER	For preparation of solutions and buffers for production and general laboratory applications
Acids	
ULTREX II	For trace element analysis ICP-MS
	Trace impurities at ppt level
BAKER INSTRA ANALYZED	For trace element analysis ICP and AAS
	Trace impurities at ppb level
BAKER ANALYZED	For general analytical applications and AAS
	Trace impurities at ppm level
BAKER	For production and general laboratory applications
	Limited specifications
Solvents	
BAKER ANALYZED LC/MS	For analytical and preparative LC/MS
BAKER ANALYZED HPLC	For analytical and preparative (isocratic and gradient) LC
ULTRA RESI ANALYZED	For organic residue and environmental analysis
BAKER ANALYZED Spectroscopy quality	For UV and GC Spectroscopy
BakerDRY	For organometallic and other dry synthesis
	Low water
BAKER BIO-ANALYZED	For biotech applications, such as:
	DNA/RNA Synthesis (low water)
	Peptide synthesis (low amines and peroxides)
BAKER ANALYZED	For general analytical applications and extractions
BAKER	For production and general laboratory applications
	Limited specifications

Continuation of table 3a Product qualities versus application

Product	Application
Chromatography	
BAKERBOND Columns and Media	Liquid chromatography of small and large molecules
Solid Phase Extraction Products	Sample clean-up and concentration with BAKERBOND SPE columns and Speedisk columns
Biochromatography Products and Media	Designed for bioseparation and bioprocessing
Titrimetry	
HYDRA-POINT	For Karl Fischer water determination
Volumetric Solutions	For all kinds of titrimetric methods
DILUT-IT	For preparation of volumetric solutions
Primary Standards	Alkalimetric and reductometric standards for standardisation of volumetric solutions
Calibration standards	
Primary Standards	Alkalimetric and reductometric standards for standardisation of volumetric solutions
Standards for Trace Element Analysis	Standards for trace element analysis:
	Single and multi-element ICP standards
	Atomic spectral standards
Buffers	
Buffers	For pH calibration range from 1 to 14 (bottles), ready to use
DILUT-IT	For pH calibration range from 1 to 14 (ampoules), concentrates to be diluted
Molecular biology	
ULTRAPURE BIO REAGENTS	General laboratory reagents, biological buffers, amino acids, sugars, denaturants, detergents,
	biological stains, electrophoresis gels, chemicals and stains
Testrips	
Testrips	For quick testing, such as pH, semi-quantitative determinations of ions and water hardness
Safety products	
Spill Kits	For spill cleanup and removal
Packaging	Special packaging to eliminate risk to users of chemicals, such as EcoTainer, CYCLE-TAINER
5 5	System
Strategic contract manufacturing	
Special Solutions and mixtures	Does your company need to outsource a manufacturing process or produce customized
,	chemicals? Let J.T. Baker be Your Outsource Resource



More Information

For more information about our products, refer to our website www.jtbaker.com, contact our Customer Service or send us an e-mail at sales@www.jtbaker.nl

Chromatography products

In the field of separation science, Mallinckrodt Baker stands out with a prominent role understanding customer needs and offering products of unique separation capabilities for more than 25 years. Our reliable and highly efficient chromatography products for analysis and purification are designed to deliver optimum performance, reproducibility and easy scale-up, without changing the quality of established methods. Our chromatography product portfolio includes:

- BAKERBOND Columns and Media
- Solid Phase Extraction Products
- Biochromatography Products and Media
- Optimized Spherical Flash Silica

BAKERBOND Columns and Media

Liquid Chromatography of Small and Large Molecules

We are a basic manufacturer of BAKERBOND bonded phases, synthesised using trifunctional silane




Table 3b General product overview

Product	Silicia based SPE products		Polymer based SPE products
BAKERBOND spe column	BAKERBOND spe columns reverse phase normal phase ion exchange adsorption drug of abuse	TO MODEL SPEC	BAKERBOND spe columns Styrene Divinylbenzene copolymer (SDB)
BAKERBOND Speedisk column	BAKERBOND Speedisk columns reverse phase normal phase ion exchange adsorption drug of abuse	O sin Special Column	BAKERBOND Speedisk Divinylbenzene (DVB) columns H2O-Philic DVB H2O-Philic SC-DVB H2O-Phobic DVB H2O-Phobic SC-DVB H2O-Phobic SC-DVB
BAKERBOND Speedisk 96 - Well Plates	BAKERBOND <i>Speedisk</i> 96 Silica Columns (in holder)	SPEEDISK 96 ************************************	BAKERBOND <i>Speedisk</i> 96 Polymer Columns (in holder) BAKERBOND <i>Speedisk</i> 96- Well Plates Polymer Sorbents (compatible with automatic systems)
BAKERBOND Speedisk Extraction Disk	BAKERBOND Speedisk Extraction Disks for manual extraction stations and for automated extractors		BAKERBOND <i>Speedisk</i> Extraction Disks for manual extraction stations and for automated extractors

synthesis chemistry providing increased resistance to hydrolysis, greater stability at pH extremes and reduced secondary silanol interactions, enabled by consistent ligand density.

BAKERBOND chromatography media and columns demonstrate enhanced levels of performance meeting your most demanding needs in small and large molecule separations. A variety of functionalities from normal phase and reverse phase to ion exchange are available for separation of small organic molecules to large peptides and oligonucleotides.

Families of Silica Based Products

- Analytical HPLC columns and spherical media available with 3 µm and 5 µm particles and 120 Å pore diameter
- Chiral HPLC columns with 5 μm spherical or 10 μm irregular particles and 60 Å pore diameter
- C₁₈ preparative HPLC columns and media available with 10 µm irregular particles and 150 Å pore diameter
- Preparative LC media for low-pressure purification of grams to kilograms, available as 40 µm irregular particles with 60 Å pore diameters
- Biochromatography media incorporating our proprietary two-part bonding chemistry are also available with 5, 15, and 40 µm particle diameter and wide pore diameter (300 Å)

Solid Phase Extraction Products

For more than 25 years, Mallinckrodt Baker has been your essential companion in Solid Phase Extraction (SPE). Dr. Morris Zief, a R&D scientist at our company, first coined the term SPE in 1980. As a key to successful sample preparation, we offer a wide range of silica and polymer based SPE products. Following the demands in research and development for purity, as well as for improved detection and quantification limits in analytical techniques, we provide the market with a variety of J.T.Baker SPE products with guaranteed uniformity and reproducibility. Our wide choice of sorbents improves and simplifies sample clean-up and concentration.

The family of SPE products (Table 3b) includes:

- BAKERBOND spe Columns
- BAKERBOND Speedisk Columns
- BAKERBOND Speedisk 96 Columns
- BAKERBOND Speedisk 96-Well Plate
- BAKERBOND Speedisk Extraction Disk

NEW Polymeric sorbents

Solid Phase Extraction as a sample preparation technique often requires stationary phases stable over the whole pH range. Among the wide range of SPE sorbents, Mallinckrodt Baker offers you polymeric sorbents, which can improve the quality of your sample preparation. BAKERBOND polymeric resin particles have large surface areas, they are highly rigid and stable over pH range 1-14. These patented (US Patent 6,875,817 April 5, 2005) polymeric sorbents are available in hydrophilic and hydrophobic forms and as mixed mode ion exchangers. They are highly recommended when advanced detection methods have to be used, offering the unique properties of a polymer. Excellent structural, mechanical and chemical properties, packaging purity and loadability, as well as stability of the polymeric sorbents, offer improved sample clean up and preconcentration. Highly efficient and selective polymeric sorbents with high throughput performance will enable high recovery and reproducibility of the results, even in the SPE of compounds of different polarities.



Application note

J.T.Baker customers can obtain technical support quickly and easily by searching our web site www.jtbaker.com. In the J.T.Baker Technical Library, at www.jtbaker.com/chromatography/techlib You will find a large body of information about J.T.Baker chromatography products and SPE applications in your field of interest.



Processors

We offer a variety of J.T.Baker processors to meet your extraction needs in three ways:

- Standard vacuum processors for flexibility of processing different SPE devices in the same experiment
- Positive pressure processors for exceptional precision, control, and reliability and
- Speedisk extraction stations

Please refer to the Chromatography section for more detailled information about these products.

Biochromatography Products and media

BAKERBOND biochromatography products and media are specifically designed following the recent advances in the fields of bioseparation and bioprocessing.

Combining chromatographic properties of wide-pore silica with a proprietary surface chemistry has created excellent particle and surface properties. Offering confirmed stability, selectivity, packaging purity and loadability, BAKERBOND biochromatography products and media enable analysis and purification of proteins, peptides and polynucleotides with the possibilities of linear scale-up. The versatile family of biochromatography products and media offer stable, efficient and highly selective phases that meet your requirements in a variety of applications.

Optimised Spherical Flash Silica

Flash chromatography is a widely used preparative technique for separating a variety of primarily organic compounds quickly and cost effectively. In order to provide exceptional separation performance, Mallinckrodt Baker tailored spherical silica gel with stringent control of mechanical and structural properties. The separation is fully comparable to HPLC, higher flow rates with improved productivity are achieved as well as optimal physico-chemical properties.

Optimised Spherical Flash Silica gives:

1. better separation performance

A narrower particle size distribution as well as the spherical shape delivers separation performance closer to HPLC.

2. higher capacity

Better flow properties and pore distribution means more access to pores for the analyte of interest, more access delivers greater capacity.




◀ J.T.Baker Spherical Flash Silica Cartridges


3. higher flow rates – higher productivity


More consistent particle size as well as the spherical shape means a more even packing and distribution in the cartridge. These enable the solvent and analyte to move through faster at lower pressure!

4. higher mechanical strength

Spherical media means that there are no sharp edges to break off and create more small particles. Small particles create fines and increase backpressure.







▲ Spherical ▲ Irregular ▲ Surface of Spherical

Reagents for trace element analysis

Trace element analysis is a vital analytical parameter in several industries. In the table below components are listed in which trace elements are analysed.

Inorganic	Organic
High Purity Metals	• Liquids
Ferrous and Non-Ferrous Alloys	• Solids
Superalloys	Pharmaceuticals
Semiconductors	Nutrients
Superconductors	Agricultural Products
Ceramics	Biochemicals
• Glasses	
Carbon and Diamond	
• Salts	
Minerals and Geological Materials	
Synthetic Crystals	
Liquids	

Techniques available for analysing trace elements:

- Inductively Coupled Plasma Mass Spectrometry (ICPMS)
- 2. Glow Discharge Mass Spectrometry (GDMS)
- 3. Spark Source Mass Spectrometry (SSMS)
- 4. Atomic Absorption Spectrophotometry (AAS)
- 5. LECO Combustion Analysis

Analytical reagents with the J.T.Baker brand

For laboratory analysis we offer three qualities, for use in sample preparation:

- BAKER ANALYZED Reagent, including ACS, ISO and meets reagent specification for testing of USP/NF and/or Ph. Eur. Monographs.
- BAKER INSTRA-ANALYZED Reagent.
- ULTREX II Reagent.

Sample preparation

Analytical reagent acids with the trade name BAKER ANALYZED and high purity acids with trade names BAKER INSTRA-ANALYZED and ULTREX II used for dilution or solvation. Atomic spectroscopy, both AAS and ICP are methods that require liquefied samples. Therefore, all sample matrices must be dissolved or digested. For trace analysis, the reagents used for sample preparation must be of the highest available purity. As a basic rule, trace element impurities of the acids used for dissolving a sample must be a factor 10 lower than the analyte being analysed.

Application

See table 3c

Calibration

For calibration in atomic spectroscopy we offer a wide range of single and multi-element standards.

See table 3d

Mallinckrodt Baker offers the following with J.T.Baker products:

Built-in Excellence	Making life simpler
Ultra-low levels of metal impurities	Guarantee of trouble-free analysis
Advanced purification technology	Guarantee of consistent quality
Actual values for trace element impurities	Guarantee of your analytical results

Table 3c High purity acids

Quality definition	Trace element impurities, analysed at	Number of specifications	Application
ULTREX II	ppt level	60 elements tested to parts-per-trillion levels	ICP-MS
BAKER INSTRA-ANALYZED	ppb level	35 elements tested to 1 to 10 parts-per-billion levels	ICP-MS, ICP-OES and AAS
BAKER ANALYZED	ppm level	20 elements tested to parts per billion and parts-per-million levels	AAS



◀ One small step for a man, one giant leap for mankind (Neil Armstrong, July 20, 1969)

ULTREX II high purity acids were specifically developed for NASA to help complete trace element analysis of moon rock samples

Table 3d Calibration

Quality definition	Description	Purity	Application
BAKER INSTRA-ANALYZED	Single-element plasma standards	99.999%	ICP
BAKER INSTRA-ANALYZED	Multi-element plasma standards	N/A	ICP
BAKER ANALYZED	Atomic spectral standards	99.99%	ICP, AAS
BAKER	Atomic spectral standards	99.9%	AAS

Solvent applications

Solvents used for different analytical and process applications also differ in quality. In the schematic below, all our solvent qualities and applications offered are represented.

Besides the chemical purity, the solvent packaging is principal to guaranteeing the product integrity. Besides

glass bottles, we also offer other packaging options, such as 5 L EcoTainer and HDPE bottles. If you use bigger quantities of high purity solvents, we can offer you our CYCLE-TAINER High Purity Solvent Delivery System.

Mallinckrodt Baker offers

Built-in Excellence	Simplicity
Broad range of reagents and packaging	User-friendly and versatile
Multi-step purification technology	Trouble-free analysis
Specified and analysed to its application	Reproducible analytical results
Filtration through 0.2 micron filter	Ready-to-use

Analytical and preparative LC/MS	BAKER ANALYZED LC/MS	BAKER ANALYZED HPLC	Analytical and preparative isocratic and gradient LC
Organic residue and environmental analysis in gas chromatography	BAKER ULTRA RESI ANALYZED	BAKER ANALYZED Spectroscopy	UV and GC Spectroscopy
Organometallic and other dry synthesis	BakerDRY	BAKER BIO ANALYZED	Biotechnology applications Oligonucleotide and peptide synthesis
General analytical applications and extractions	BAKER ANALYZED	BAKER	Production and general laboratory applications

CYCLE-TAINER High Purity Solvent Delivery System

Preserve Purity, Protect People

Our unique solvent delivery system solves several major challenges related to solvents in laboratories, maintaining employee safety, assuring product purity, reducing packaging and disposal costs, reducing storage space, and complying with governmental regulations.

How does it work?

Using this system, solvents are delivered to the laboratory in stainless steel containers, empty drums collected, refilled and returned to the lab, as often as you require.

Safety stock at Containers in use at customer

Transport from Mallinckrodt Baker customer to to customer

Mallinckrodt Baker

Containers filled at Mallinckrodt Baker

▲ Flow chart CYCLE-TAINER System loop from Mallinckrodt Baker to customer and back

The CYCLE-TAINER System is available in a variety of sizes, ranging from 20 litres to 1400 litres. These containers can be placed in individual labs or accessed by several labs, through a system of stainless steel pipes that are configured and installed just like pipes for a lab's high purity gas system.

Because the CYCLE-TAINER System is completely contained, there is no risk of solvent contamination, or risks of solvent handling.

Why a High Purity Solvent Delivery System? Enhances employee safety:

The system reduces employee exposure hazards and eliminates concern for glass container breakage and storage.

Eliminates waste disposal costs

Eliminates waste packaging, no glass bottles, styrofoam or cardboard to dispose of. Without bulk packaging, storage space is maximised.

Enhances working convenience


Removes necessity to move heavy cases of solvent from storage room to the lab. Just press a button!!

Assures point-of-use solvent purity

The risk of solvent contamination is eliminated with a closed system. Also for each solvent, a different colour-coded quick connector is used. Shelf life is extended for at least two years, vastly reducing solvent degradation, even when solvents are dispensed from the container.



▲ Top of CYCLE-TAINER System with colourcoded quick connectors for the solvent and gas supply



▲ Example product integrity Acetonitrile BAKER BIO-ANALYZED

Easy connection

Delivery Hook up Piping Dispensing Solvent delivered in CYCLE-TAINER System Ouick connectors used Solvent delivered to point of use Solvent dispensed with tap or connected directly to equipment

200 L

201 L

201 L

▲ CYCLE-TAINER System sizes

Applications

- Solvent with a tap dispenser
- Direct connection to instrumentation HPLC equipment and DNA synthesizer

Reagents for pH Measurement, Titrimetry and Water Determination according to Karl Fischer

In several industries Titrimetry, pH measurement and Karl Fischer determination are required analytical techniques, for both research and quality control purposes.

Applications for all relevant areas:

pH measurement

- Cosmetics
- Detergents/surfactants
- Electroplating baths and dyes / pigments
- Water
- Soi
- · Surfaces, such as skin, paper, leather
- Biological materials, such as blood, urine
- Food / drinks / beverages / flavours / wine

Titrimetry, including Karl Fischer

- Water / environment
- Pharmaceutical / vitamins
- Petrochemistry
- Plastics / photo
- Food / drinks / beverages / flavours / wine
- Tobacco
- Biology / biochemistry
- Cosmetics
- Detergents / surfactants
- Metals / plating
- Fertilisers / base materials / explosives
- Textiles / paper / ceramics
- Paints / lacquers / solvents
- Mineral resources / cement
- Energy / power plants
- Miscellaneous: Kjehldal / hypochlorite / chlorides / chlorates / sulphates / et cetera



Mallinckrodt Baker offers

Analysis against NIST-SRM

Broad range of reagents and packaging Actual values for ready-to-use solutions

Analytical techniques available:

pH measurement

- pH testrips
- pH measurement with electrodes

Titrimetry, including Karl Fischer

- Acid/Base titrations with use of indicator
- Redox titrations, such as Potassium Permanganate
- Potentiometric titration
- Karl Fischer titration: Volumetric, Coulometric and Coulometric Oven

Analytical reagents with the J.T.Baker brand:

For pH measurement we offer

- BAKER ANALYZED Buffer Solutions, traceable to NIST
- DILUT-IT Buffer concentrates in ampoules from pH 1 and 14
- BAKER-pHIX pH Papers with colour scale

For Titrimetry we offer

- BAKER ANALYZED Volumetric solutions, traceable to NIST
- DILUT-IT Volumetric concentrates in ampoules
- BAKER ANALYZED Primary standards
- BAKER ANALYZED Reagent, including ACS, ISO and meets reagent specification for testing of USP/NF and/or Ph. Eur. Monographs

Karl Fischer titration products

 HYDRA-POINT Karl Fischer reagents are pyridine-free and especially formulated for existing application procedures.
 The HYDRA-POINT reagent line is comprised of the following products:

Volumetric reagents

- One component
- Two component
- Buffer, dry methanol, sodium tartrate standard

Coulometric reagents

- With diaphragm
- Without diaphragm

Simplicity
User-friendliness and versatility
Guarantee of your analytical results
Increased reliability of analytical results

Table 3e Applications pH Measurement, Titrimetry and Karl Fischer

Quality definition	Application	Comment
pH Measurement		
BAKER ANALYZED Buffer Solutions	Calibration pH meter	Traceable to NIST
DILUT-IT Buffer concentrates	Calibration pH meter	Traceable to NIST
BAKER-pHIX pH Papers	pH check/indication	Quick pH check
Titrimetry		
BAKER ANALYZED Reagents	Sample preparation prior to titration	Analytical reagent grade, ACS, ISO and meets reagent specification for testing of USP / NF and/or Ph. Eur. Monographs.
BAKER ANALYZED	Titrant	± 0.3%. When method is validated, 1 burette
Volumetric Solutions		is used and titration is done by 1 LAB technicar 0.1% is achievable.
DILUT-IT Volumetric concentrates	Titrant	± 1% is normal. Higher accuracy is possible, but depends on dilution and glassware used.
BAKER ANALYZED Primary standards	Primary Standard	99.95 – 100.05%
Karl Fischer		
HYDRA-POINT Volumetric reagents	Titrant for water concentrations > 100 ppm. Several reagents are available for different applications with 1 and 2 component reagents	N/A
HYDRA-POINT Coulometric reagents	Titrant for water concentrations < 100 ppm. Several reagents are available for different applications, for analysis with and without diaphragm	Excellent end point stability for water concentration below 10 ppm
HYDRA-POINT Buffer	For stabilising the pH	The buffer capacity is 5 mmol acid per ml.
HYDRA-POINT Dry Methanol	Sample solvation	Water maximum 0.01%
Sodium Tartrate Dihydrate	Standardisation	N/A

Reagents for Molecular Biology and Biotechnology

Biotechnology researchers expand the boundaries of science to benefit mankind by providing better healthcare, enhanced agriculture, and a cleaner and safer environment. BIO: The use of biological processes. TECHNOLOGY: To solve problems or make useful products. "NEW" BIOTECHNOLOGY: The use of cellular and biomolecular processes to solve problems or make useful products.

The Technologies and their Applications:

- Bioprocessing Technology
- Monoclonal Antibodies
- Cell Culture
- Recombinant DNA Technology
- Cloning
- Protein Engineering
- Biosensors
- Nanobiotechnology
- Microarrays

Mallinckrodt Baker offers

Built-in Excellence	Simplicity
Broad range of reagents and packaging	Reagents for your analytical and process application,
	including scale up from Beaker to Bulk
Biopharmaceutical reagents manufactured according to cGMP	Your ultimate guarantee of process consistency
and specified to all pharmacopoeias	

Reagents with the J.T.Baker brand

For applications in Molecular Biology and Biotechnology we offer the following reagents and solvents:

- Biopharmaceutical product line
- Bio reagents
- · Reagents for oligonucleotide synthesis
- BAKER BIO-ANALYZED Solvents
- BAKER ANALYZED HPLC and LC/MS Solvents and Reagents
- Chromatography products

Applications

Biopharmaceutical product line

From drug discovery to full scale manufacturing with cGMP-produced chemicals in Beaker to Bulk packaging options:

- USP, NF and FCC Products for Pharmaceutical and Biotech processes
- Global Pharma multicompendial-tested product line
 - Manufactured under cGMP requirements
 - Meets requirements of USP and NF
 - Meets the specifications of EP (PhEur),
 BP and JP
- Reagents for fermentation/cell culture
 - Sugars
 - Amino acids
 - Vitamins and minerals
- Reagents for product harvest/purification
 - Guanidine hydrochloride
 - Guanidine thiocyanate
 - ULTRAPURE BIO REAGENTS
- Biological buffers
 - ULTRAPURE BIO REAGENTS

- High purity solvents
 - BAKER BIO-ANALYZED solvents
- Chromatography products

Bio reagents

Our ULTRAPURE BIOREAGENT and BAKER ANALYZED Biochemical reagents are general laboratory reagents, biological buffers, amino acids, sugars, denaturants, detergents, biological stains, and electrophoresis gels, chemicals and stains.

Reagents for oligonucleotide synthesis

Reagents applicable to a wide range of synthesisers. Besides standard products, tailor made reagents are possible.

- Acetonitrile: see features under BAKER BIO-ANALYZED Solvents
- 2. Deblock
- Activator: substitutes for 1H-Tetrazole are available, such as 5-Ethylthio-1H-tetrazole (ETT) and 4,5-Dicyanoimidazole (DCI)
- 4. Capping A
- 5. Capping B
- 6. Oxidizer

BAKER BIO-ANALYZED Solvents

Features per product, see table 3f


BAKER ANALYZED HPLC and LC/MS Solvents and Reagents

Chromatography products

BAKERBOND spe and Speedisk technology

Table 3f BAKER BIO-ANALYZED Solvents

Solvent	Product code	Comment
Acetonitrile, for DNA/RNA synthesis	8134	Water < 10 ppm in CYCLE-TAINER System
Acetonitrile, for DNA/RNA synthesis	8144	Water < 30 ppm in Bottles
Dichloromethane	9316	Water < 30 ppm, Acidity <0.0001 meq/g
Dimethyl Sulfoxide	9234	Water < 250 ppm
Dimethylformamide	9344	Amines < 5 ppm, Water < 400 ppm
Ethyl Acetate	9276	Water < 300 ppm, Acidity <0.0008 meq/g
n-Heptane, 99% n-Heptane	9185	Water < 100 ppm
Methanol	9091	Water < 50 ppm
1-Methyl-2-pyrrolidone	9261	Amines < 100 ppm, Water <200 ppm
Pyridine HPLC 'BAKER ANALYZED'	9393	Water < 100 ppm
Tetrahydrofuran HPLC 'BAKER ANALYZED'	9439	Peroxide < 10 ppm, Water < 50 ppm



Mallinckrodt Baker's chemical expertise and over 30 years of experience are the foundation of our J.T.Baker products for **Hematology**. We understand not just cell lysing and chemical interference of compounds on cell membranes, but we also understand the characteristics of these compounds in aqueous solutions and their long-term behaviour in the presence of other chemicals.

Within our **Histopathology** product line we offer a comprehensive range of fixatives, reagents, solvents, stains, environmentally friendly cleaning and auxiliary products specially formulated and purified for use in the histopathology and cytopathology laboratories.

Our Hematology and Histopathology products are produced using state-of-the-art production techniques in modern, certified manufacturing facilities, to assure a minimum of impurities.

R&D and production is centralised at Mallinckrodt Baker's European headquarters in Deventer, The Netherlands.



Products for Hematology

Product range

Together with the broad range of reagents for cell analysers, Mallinckrodt Baker supplies a variety of hematology controls to monitor daily accuracy and precision.

Reagents for cell analysers

We offer a complete range of reagents. This range consists of Diluids (diluting fluids), CyMets (lysing reagents) and several types of rinsing and cleaning solutions. Diluid is a specially filtered non-sterile blood diluting fluid for use in cell counting and sizing.

J.T.Baker lysing reagents are available both as Cyanide (CN) free / SDS free solutions and as traditional



The hematology reagents manufacturing facility at Mallinckrodt Baker in Deventer, The Netherlands.

cyanide containing reagents. The CN/SDS free reagents are non-poisonous and more stable than the traditional reagents.

Hematology controls

Our range of controls covers 8 parameter control for most instruments and differential controls for various 3-part and 5-part differential analysis including reticulocyte counting. All controls are manufactured and delivered on a four or six times per year schedule (depending on the stability). Orders are guaranteed when received 9 weeks in advance. All deliveries are accompanied by assay values for the applicable instruments.

Rapid Stat, quality control and assessment program

Rapid Stat is a fast and convenient Quality Assessment program for hemocytometry. Users of all J.T.Baker Hematology Controls are allowed to participate in our external, between-lab and between-instrument

statistical program. This service is provided to inform users of our controls about their individual performance compared to other users.

Hematology stains

Besides May-Grünwald/Giemsa, Wright and Leishman staining solutions for traditional methods, we also offer an excellent fast alternative - Hemacolour for blood smear staining within 1 minute. Additionally, we offer a microscopic stain, Reticount, for staining of Reticulocytes.

All our products are distributed through our worldwide network of distributors. Some reagents are only available through specified instrument dealers.

Certification and legislation

Effective from June 2000, In Vitro Diagnostic Medical Devices (IVDs) have to comply with the IVD Directive 98/79/EC and have to show the CE mark when the registration is completed. Mallinckrodt Baker is in full compliance with this registration. All products subject to the directive are registered and CE marked.

Packaging to suit your needs

J.T.Baker hematology reagents are supplied in various packaging options. Diluids (diluting fluids) are packed in 5, 10 or 20 L polycubes. CyMets (lysing reagents) are provided in 0.5 or 1 L HDPE bottles or

5 L polycubes. Equipment dedicated bottles are also available. Hematology controls are supplied in 2.5 or 8 ml screw-cap vials or 4.5 ml pierceable vials.

Products for Histopathology

Product range

Mallinckrodt Baker supplies a variety of solvents, reagents and stains for use in histopathology and cytopathology laboratories.

Fixatives

Formaldehyde solutions

Effective and reproducible fixation requires buffered formaldehyde. Non-buffered formaldehyde solutions are acidic and therefore inconsistent for pH. When formaldehyde is acidic, formation of formaline-heme

pigments may occur. These pigments appear as black deposits in the tissue. To obtain optimised fixation the buffer concentration is also very important. We offer formaldehyde fixatives as a ready-to-use solution or as a concentrated solution.

All formaldehyde fixatives are phosphate buffered with an optimised buffer concentration.

Formalin Neutraliser

Mallinckrodt Baker developed a product that facilitates the safe disposal of used or spilled formalin. Formalin Neutraliser is a powerful chemical capable of neutralising formalin.

Paraffin

UltraPar, 54 - 56°C Paraffin

UltraPar is a special particle-free paraffin intended for use in tissue embedding and sectioning in histopathology.

When melted, UltraPar is a clear and bright solution.

The solubility of UltraPar in UltraClear is extremely good.



UltraClear - Environmentally and human friendly solvents to replace Xylene and Toluene UltraClear

UltraClear is an iso-paraffin based clearing reagent.

UltraClear may be used as a Xylene replacement during tissue embedding, deparaffination and staining processes. UltraClear is less-toxic, less-carcinogenic, less-flammable, odourless and easy to adapt with existing procedures.

UltraKitt

UltraKitt is a mounting medium which is miscible with UltraClear, xylene and toluene. UltraKitt is:

- Unique, patent for Europe and USA is granted and filed
- Safe, as it does not contain dangerous aromatic solvents

High quality staining solutions

Staining May-Grünwald and Giemsa, Wright and Leishman

These stains are used for tissue sections, cytology smears and bone marrow. Our production methods and formulations are modified to the highest standards, resulting in optimised colour intensity of stained cells. Staining Histology (H&E)

- Our hematoxyline is available according to Mayer or modified according to Gill II / Harris. Modern technology allows a time and temperature controlled manufacturing process. This process guarantees a final product with long stability and consistent, bright and clear staining of the cell nucleus. In addition we also offer a Scotch solution, which enables a stable and optimised bluing of the cell nucleus.
- J.T. Baker Eosine Y stains are offered as an alcoholic or as an aqueous solution. Our Eosine Y counter stains are optimised in combination with our hematoxyline stains. For advice on staining procedures, we offer product information bulletins on request or on our website.

Staining Cytology (Papanicolaou)

Experienced laboratory technicians claimed that
results obtained with the J.T.Baker PAP stains are
of superb colour intensity. All of our PAP staining
solutions are purified and free from contamination.
We offer Papanicolaou 1 (hematoxyline according to
Gill II / Harris), Papanicolaou 2A (Orange G6),
Papanicolaou 2B (Orange II) and Papanicolaou 3B
(EA 50). For more details and recommended
staining procedures, we offer product information
bulletins on our website.

Alcohols and Solvents

<u>Alcohols</u>

A broad range of ready-to-use alcohols and mixtures are available in convenient and easy-to-carry HDPE cans. Solvents

In addition to our widely used UltraClear, we offer both xylene and toluene in a Histo grade quality.



More Information

Our website, www.jtbaker.com (select Clinical) informs you in detail about our product program and suggested working protocols.



Our information society makes extensive use of semiconductor devices, in personal computers, cell phones, digital cameras and so on. Manufacturing of semiconductors that are the heart of these devices is very sensitive to contamination. As a result, manufacturing of semiconductor devices involves many cleaning steps and surface preparation steps.

For many of these steps, highly specialised, very pure chemicals are needed.

Mallinckrodt Baker, as a specialist in high purity chemicals, has been a natural partner for the semiconductor industry from the beginning. Cleaning and surface preparation are part of our core competencies. We have the technical knowledge and capabilities to support most lithographic processes. We offer three product lines:

- Wet Process Chemicals
- Photo Resist Strippers and Residue Removers
- Photo Resist Ancillaries and Specialty Products



Wet Process Chemicals

The finer the circuit geometry, the purer your chemicals need to be. To meet those needs, Mallinckrodt Baker has expanded its line of J.T.Baker fab process chemicals to include products with 1.0 and 0.1 ppb trace-metal impurity levels. In addition to J.T.Baker fab process solvents and acids, we offer a line of mixed acid, aluminium, and buffered oxide etchants. All of these products are available in a variety of container configurations including reusable containers.

Table 5a Wet process chemicals

J.T.Baker Grade	Trace Impurity Level	Device Geometry	SEMI Grade	SEMITier
CMOS	200-10-ppb	1.2 µm	1	-
Finyte/VLSI (Europe)	10-ppb	0.8 to 1.2 μm	2	А
Finyte 1	1.0-ppb	0.2 to 0.8 μm	3	В
ULTRYTE	0.1-ppb	0.09 to 0.2 μm	4	С

Photoresist Strippers and Residue Removers


As the industry drives to smaller geometries, device manufacturers are faced with the integration challenges associated with new materials such as copper interconnects and low-K dielectrics. To reduce time to market, device manufacturers are turning to their suppliers for technical solutions to these new cleaning challenges. The Mallinckrodt Baker Applications Engineering and Research staff has a thorough understanding of today's technical challenges and is

committed to developing the advanced chemistries required to meet the needs of the semiconductor, flat panel display, packaging (flip chip, bumps), disk drive, and MEMS industries.

Table 5b Overview by Product Family

luct lines provide versatile cleaning of bulk photoresist and ash/etch residue. ctive in traditional semiconductor (aluminium and silicon dioxide layers), flat panel
ctive in traditional semiconductor (aluminium and silicon dioxide layers), flat panel
hips/bumps), and disk drive applications.
e of residue removers are greater than 80% aqueous and are compatible with
on dioxide, and low-K dielectrics. They are effective in the removal of etch, ash residue,
om metal and via structures.
sist strippers and ash/etch residue removers has been specifically designed for use on
lus sensitive low-K and porous low-K dielectrics.
C .

New technology development: Copper and low-K



We are unique and especially strong in innovative products for the newest copper and low-K semiconductor devices, with incredibly small detail sizes below 100nm.

The drive for faster and smaller devices has presented a variety of challenges for device manufacturers.

One of the newest challenges is working with new materials of construction such as copper interconnects and low-K dielectrics. These new materials require strippers and residue removers that will not corrode the metal or damage the delicate dielectrics. Mallinckrodt Baker has developed new J.T.Baker strippers and residue removers for flat panel displays, compound semiconductor, and traditional aluminium/silicon dioxide technologies.

Photoresist Ancillaries and Specialty Products

In addition to the full line of J.T.Baker photoresist strippers and residue removers, Mallinckrodt Baker offers a line of photo resist ancillaries. Our photoresist ancillaries support every step of your process including priming, edge bead removal, development and rinse – you can depend on the J.T.Baker brand for all your chemical needs. All our photoresist ancillaries are available in reusable containers.














More Information

For the latest product developments or for more information about our products, you can visit our special Internet site www.jtbaker.com/micro or send an e-mail to sales@www.jtbaker.nl



Mallinckrodt Baker is the supplier of High Purity and Performance Excipients for the manufacture of traditional synthetic and biopharmaceutical therapeutics.

We manufacture and supply acids, solvents, salts and solutions in grades suitable for the pharmaceutical, nutritional, ophthalmics (contact lens manufacturers), cosmetics, food and beverage, dental health, special batteries and other industrial markets.



Our Pharmaceutical products

Our Pharmaceutical products are:

Global Pharma

GMP manufactured products that meet multi-compendial specifications: BP, Ph. Eur., USP, JP, NF, FCC.

ULTRAPURE BIOREAGENTS

General laboratory reagents, biological buffers, amino acids, sugars, denaturants, detergents, biological stains, and electrophoresis gels, chemicals and stains.

Synthesis (Process) chemicals

High purity chemicals in larger pack sizes.

Products Mallinckrodt Baker supplies for Biopharmaceutical processes

Upstream processing Fermentation/Cell Culture Amino Acids Minerals **Fermentation** Sugars Downstream processing Vitamins **Product Harvest** Ammonium Sulfate Biotech Reagent **Product Harvest** Urea USP ULTRAPURE BIOREAGENT Guanidine Hydrochloride Sodium Hydroxide Polysorbate 20, vegetable derived Polysorbate 80, vegetable derived **Initial Purification Buffer Component Salts** BAKERBOND Chromatography Media for Initial Purification Large Molecules Sodium Dodecyl Sulfate Ion Exchange Hydrophobic Interaction - Reverse Phase Affinity **Buffer Exchange** Tromethamine (Tris) **Buffer Exchange** Tris Hydrochloride **Buffer Component Salts** Biological Buffers ULTRAPURE BIOREAGENTS (UPR) Final Purification/Polishing BAKERBOND Chromatography Media for Final Purification/ Large Molecules **Polishing** Ion Exchange - Hydrophobic Interaction Reverse Phase Fill and Finish Ammonium Sulfate Biotech Reagent Fill and Finish Polysorbate 20, vegetable derived, low peroxide New! Polysorbate 80, vegetable derived, low peroxide Mannitol, low endotoxin Column Cleaning Tromethamine (Tris) Sodium Hydroxide Solutions Tris Hydrochloride Column Cleaning Phosphoric Acid Solutions Amino Acids Guanidine Hydrochloride Solutions Facility Cleaning Protocol C3 Sterile Cleansers **Facility Cleaning** Custom aqueous solutions can be made to your specifications Ask your sales person for details. **Quality Control** HPLC Solvents **Quality Control** HPLC Columns SPE Columns Products with a blue bullet (•) are Buffers Volumetric Solutions manufactured under PharmaTest Reagents and Solutions cGMP conditions ACS Grade Reagents

Beaker to Bulk packaging for easy scale-up



Beaker to Bulk packaging means we fill identical product in containers suitable for each stage of your development. Using chemicals that can be easy scaled up from R&D to production is fundamental to your business. Mallinckrodt Baker helps with scalable packaging optimised for the pharmaceutical environment. Whether you need 500 grams or 25,000 kilos, your raw-material comes from the same manufacturing process to ensure consistency.

Global Pharma

Global Pharma line is widely recognised as a "Best Practice". Developed to aid our customers as they prepare for production or clinical trials in Europe or Asia, the Global Pharma multicompendial product line will help you in meeting your regulatory requirements.





Quality Systems That You Can Depend On

A J.T.Baker and Mallinckrodt manufactured product that is labelled as USP, NF, FCC or cGMP-produced is

- Manufactured in a GMP compliant facility
- Made according to validated processes when required
- Documented via controlled procedures and records
- Traceable from raw materials to the customer
- We maintain extensive documentation, including



batch production record (BPR), label accountability, stability data, management of change (MOC) and QC test results

For BPEs (Bulk Pharmaceutical Excipients), the company adheres to Q7A and IPEC (International Pharmaceutical Excipient Council) guidelines.



More Information

Please request a copy of our complete pharmaceutical product list at one of our local offices. Visit us at www.jtbaker.com or www.mallchem.com



Analytical reagents
Histopathology reagents
DNA synthesis reagents
Pharmaceutical products
Microelectronic materials
Biopharmaceutical products

Acaci

Acacia

See for qualities meeting the requirements of cGMP, USP/NF, BP/Ph.Eur, JP or FCC, page 36

Assay

Density (g/ml) at 20°C

9020

Acetaldehyde

► CH₃CHO

M = 44.05 g/mol

1 I = 0.78 kg

FLASHPOINT −27 °C

CAS NO. 75-07-0

EINECS 200-836-8

NC CODE 2912 12 00

EC NO. 605 003 00 6

UN/ID NO. 1089

ADR/RID 3 F1

IMDG 3/I

R: 12-36/37-40

S: 16-33-36/37







Acetdimethylamide

See N,N-Dimethylacetamide

6903

Acetic Acid

99-100%, glacial / ULTREX II Ultrapure Reagent

► CH₃COOH

M = 60.05 g/mol

1 l = 1.05 kg

 $\textbf{FLASHPOINT} \ \ 40\ ^{\circ}\text{C}$

CAS NO. 64-19-7 **EINECS** 200-580-7

NC CODE 2915 21 00

EC NO. 607 002 00 6

UN/ID NO. 2789

ADR/RID 8 CF1

IMDG 8/II

R: 10-35

S: 23-26-45



Certificate Provided Reporting Actual Lot Analysis Assay (CH₃COOH) min. 99.0% Trace Impurities (in ppt) (pg/g): Aluminium (AI) max. 50

min. 99%

0.779-0.782

PRODUCT

9020.0500

PACKING

500 ml

CONT.

BOX

Aluminium (AI)	max. 50
Antimony (Sb)	max. 50
Arsenic (As)	max. 50
Barium (Ba)	max. 10
Beryllium (Be)	max. 10
Bismuth (Bi)	max. 10
Cadmium (Cd)	max. 10
Calcium (Ca)	max. 50
Cerium (Ce)	max. 10
Cesium (Cs)	max. 10
Chromium (Cr)	max. 10
Cobalt (Co)	max. 10
Copper (Cu)	max. 50
Dysprosium (Dy)	max. 1
Erbium (Er)	max. 1
Europium (Eu)	max. 1
Gadolinium (Gd)	max. 1
Gallium (Ga)	max. 10
Germanium (Ge)	max. 10
Hafnium (Hf)	max. 10
Holmium (Ho)	max. 1
Indium (In)	max. 1
Iron (Fe)	max. 50
Lanthanum (La)	max. 1
Lead (Pb)	max. 10
Lithium (Li)	max. 10
Lutetium (Lu)	max. 10
Magnesium (Mg)	max. 50
Manganese (Mn)	max. 10
Molybdenum (Mo)	max. 10
Neodymium (Nd)	max. 1
Nickel (Ni)	max. 50

Platinum (Pt)	max. 50
Potassium (K)	max. 50
Praseodymium (Pr)	max. 1
Rhenium (Re)	max. 10
Rhodium (Rh)	max. 50
Rubidium (Rb)	max. 10
Ruthenium (Ru)	max. 50
Samarium (Sm)	max. 1
Scandium (Sc)	max. 10
Selenium (Se)	act. value reported
Silver (Ag)	max. 50
Sodium (Na)	max. 100
Strontium (Sr)	max. 10
Tellurium (Te)	max. 1
Terbium (Tb)	max. 1
Γhallium (TI)	max. 10
Thorium (Th)	max. 1
Thulium (Tm)	max. 1
Γin (Sn)	max. 50
Fitanium (Ti)	max. 10
Γungsten (W)	max. 10
Jranium (U)	max. 1
Vanadium (V)	max. 10
Ytterbium (Yb)	max. 1
Yttrium (Y)	max. 1
Zinc (Zn)	max. 50
Zirconium (Zr)	max. 10

PRODUCT	PACKING	CONT.
NO.		вох
6903.0500	500 ml Teflon	

IMPORTANT: Material will freeze if stored below 17°C (63°F).

Acetic Acid

99-100%, glacial / 'BAKER HPLC ANALYZED' / for use in High Performance Liquid Chromatography

6152

•	CH ₃ COOH

M = 60.05 g/mol

1 l = 1.05 kg

 $\textbf{FLASHPOINT} \hspace{0.2cm} 40 \, ^{\circ}\text{C}$

CAS NO. 64-19-7

EINECS 200-580-7

NC CODE 2915 21 00

EC NO. 607 002 00 6

UN/ID NO. 2789

ADR/RID 8 CF1
IMDG 8/II

R: 10-35

S: 23-26-45

C corrosive

Assay (by thermometry)	min. 99.7%
Maximum Limits of Impurities:	
Residue after Evaporation (in ppm)	5
Water (H ₂ O)	0.1%
Physical Data (not specifications):	

Density (g/ml) at 20° C 1.049 Ultraviolet Absorbance (1.00-cm path vs water):

at 280 nm	max. 0.05
at 350 nm	max. 0.01
UV Cut-off, nm	max. 255

PRODUCT NO.	PACKING	CONT. BOX
6152.1000	11	6
6152.2500	2.5	4

Acetic Acid

99-100%, glacial / 'BAKER INSTRA-ANALYZED' / for Trace Metal Analysis / ACS

9524

► CH₃COOH

M = 60.05 g/mol

1 l = 1.05 kg

FLASHPOINT 40 °C

CAS NO. 64-19-7

EINECS 200-580-7

NC CODE 2915 21 00 **EC NO.** 607 002 00 6

UN/ID NO. 2789

ADR/RID 8 CF1

IMDG 8/II

R: 10-35

S: 23-26-45



Meets ACS Specifications	
Assay (by GC) (corrected for water)	min. 99.9%
Acetaldehyde	max. 0.005%
Acetic anhydride ((CH ₃ CO) ₂ O)	max. 0.01%
Color (APHA)	max. 10
Dilution Test	passes test
Residue after Ignition	max. 5 ppm
Solubility in Water	passes test
Specific Gravity at 20°C/20°C	min. 1.048
Substances Reducing Dichromate	passes test
Substances Reducing KMnO ₄	passes test
Titrable Base (meq/g)	max. 0.0004
Trace Impurities (in ppb):	
41 11 (41)	

Trace Impurities (in ppb):	
Aluminium (AI)	max. 50
Arsenic and Antimony (as As)	max. 5
Barium (Ba)	max. 20
Beryllium (Be)	max. 10
Bismuth (Bi)	max. 50
Boron (B)	max. 10
Cadmium (Cd)	max. 5
Calcium (Ca)	max. 100
Chromium (Cr)	max. 50
Cobalt (Co)	max. 10
Copper (Cu)	max. 10
Germanium (Ge)	max. 50
Heavy Metals (as Pb)	max. 300
Iron (Fe)	max. 50
Lead (Pb)	max. 10
Lithium (Li)	max. 50
Magnesium (Mg)	max. 50
Manganese (Mn)	max. 5
Mercury (Hg)	max. 5
Nickel (Ni)	max. 10
Potassium (K)	max. 200
Silicon (Si)	max. 50
Silver (Ag)	max. 5
Sodium (Na)	max. 500
Strontium (Sr)	max. 5
Tin (Sn)	max. 10
Titanium (Ti)	max. 300
Zinc (Zn)	max. 10

Trace Impurities (in ppm):

Chloride (CI)	max. 0.5
Phosphate (PO ₄)	max. 0.5
Sulfate (SO ₄)	max. 0.5

PRODUCT NO.	PACKING	CONT. BOX
9524.2500	2.5	4

(About 17.4N).

IMPORTANT: Material will freeze if stored below 17°C (63°F).

Aceti

Acetic Acid

6052 99-100%, glacial / 'BAKER ANALYZED' / ACS

► CH₃COOH

M = 60.05 g/mol **1** l = 1.05 kg

FLASHPOINT 40 °C

CAS NO. 64-19-7

EINECS 200-580-7

NC CODE 2915 21 00

EC NO. 607 002 00 6

UN/ID NO. 2789 **ADR/RID** 8 CF1

IMDG 8/II

R: 10-35

S: 23-26-45



Exceeds ACS Specifications	
Assay	min. 99.7%
Acetic anhydride ((CH ₃ CO) ₂ O)	max. 0.01%
Color (APHA)	max. 10
Dilution Test	passes tes
Residue after Evaporation	max. 0.001%
Substances Reducing Dichromate	passes tes
Substances Reducing KMnO ₄	passes tes
Titrable Base (meq/g)	max. 0.0004
Trace Impurities (in ppm):	
Aluminium (AI)	max. 0.
Arsenic (As)	max. 0.05
Chloride (CI)	max.
Copper (Cu)	max. 0.
Heavy Metals (as Pb)	max. 0.5
Iron (Fe)	max. 0.2
Nickel (Ni)	max. 0.
Sulfate (SO ₄)	max.

PRODUCT	PACKING	CONT.
NO.		вох
6052.1000	11	6
6052.2500	2.5	4
6052.9025	25 l	



Acetic Acid 99% MOS Grade

See for detailed information section Microelectronic materials page 32 or the website www.jtbaker.com/micro

Acetic Acid

4706

1 mol/l $\,/\,$ 1 equiv. = 60.05g, 1N Solution $\,/\,$ DILUT-IT

► CH₃COOH

M = 60.05 g/mol

1 l = 1.08 kg

CAS NO. 64-19-7 **EINECS** 200-580-7

NC CODE 2915 21 00

EC NO. 607 002 00 6

UN/ID NO. 2789

ADR/RID 8 C3

IMDG 8/II

R: 34

S: 23-26-36/37/39-45





Volumetric concentrate, for dilution to 500ml.



Acetic Acid, Glacial

See for qualities meeting the requirements of cGMP, USP/NF, BP/Ph.Eur, JP or FCC, page 36



Acetic Acid Solutions

See for qualities meeting the requirements of cGMP, USP/NF, BP/Ph.Eur, JP or FCC, page 36

6004

Acetic Anhydride

97% / 'BAKER ANALYZED' / ACS

► (CH₃CO)₂O

M = 102.09 g/mol

11= 1.08 kg

FLASHPOINT 49 °C **CAS NO.** 108-24-7

EINECS 203-564-8

NC CODE 2915 24 00

EC NO. 607 008 00 9

UN/ID NO. 1715

ADR/RID 8 CF1

IMDG 8/II

R: 10-20/22-34

S: 26-36/37/39-45



Meets ACS Specifications	
Assay	min. 97.0%
Phosphate (PO ₄)	max. 0.001%
Residue after Evaporation	max. 0.003%
Substances Reducing KMnO ₄	passes test
Trace Impurities (in ppm):	
Chloride (CI)	max. 5
Heavy Metals (as Pb)	max. 2
Iron (Fe)	max. 5
Sulfate (SO ₄)	max. 5

PRODUCT	PACKING	CONT.
NO.		вох
6004.1000	11	6
6004.2500	2.5	4
6004.9025	25 l	

Acetone

'BAKER ULTRA RESI-ANALYZED' / for Organic Residue Analysis

9254

► (CH₃)₂CO

M = 58.08 g/mol

1 l = 0.79 kg

FLASHPOINT −20 °C

CAS NO. 67-64-1

EINECS 200-662-2

NC CODE 2914 11 00

EC NO. 606 001 00 8

UN/ID NO. 1090

ADR/RID 3F1

IMDG 3/II

R: 11-36-66-67

S: 16-26-9





Assay (by GC) (corrected for water)	min. 99.4%	
Color (APHA)	max. 10	
Residue after Evaporation	max. 1 ppm	
Substances Reducing KMnO ₄	passes test	
Titrable Acid (µeq/g)	0.3	
Titrable base (µeq/g)	0.6	
Water (H ₂ O)	max. 0.5%	
ECD Sensitive Impurities (as Heptachlor Epoxide):		
Single Impurity Peak (pg/ml)	max. 10	

max. 5

FID-Sensitive Impurities (as 2-Octanol):

Single Impurity Peak (ng/ml)

PRODUCT **PACKING** CONT. NO. вох 9254.1000 11 6 9254.2500 2.5 l 4

Acetone

'BAKER HPLC ANALYZED' / for use in High Performance Liquid Chromatography

8142

► (CH₃)₂CO

M = 58.08 g/mol

1 l = 0.79 kg

FLASHPOINT −20 °C

CAS NO. 67-64-1

EINECS 200-662-2 NC CODE 2914 11 00

EC NO. 606 001 00 8

UN/ID NO. 1090

ADR/RID 3 F1

IMDG 3/II

R: 11-36-66-67

S: 16-26-9







Physical Data (not enseifications)	
Water (H ₂ O)	max. 0.2%
Titrable Base (meq/g)	max. 0.0006
Titrable Acid (meq/g)	max. 0.0003
Residue after Evaporation	max. 2 ppm
Assay (by GC) (corrected for water)	111111. 99.7%

Physical Data (not specifications):

0.791 Density (g/ml) at 20°C

Ultraviolet Absorbance (1.00-cm path vs water):

at 300 HH	max. 0.0 i
at 400 nm	max. 0.01
UV Cut-off, nm	max. 330

PRODUCT	PACKING	CONT
NO.		BOX
8142.1000	11	(
8142.2500	2.5	4

Filtered through a 0.2 micron filter. Packaged under Nitrogen.

Acetone

'BAKER ANALYZED' / Ultraviolet Spectrophotometry / ACS

8001

► (CH₃)₂CO

M = 58.08 g/mol

1 l = 0.79 kg

FLASHPOINT $-20 \, ^{\circ}\text{C}$

CAS NO. 67-64-1

EINECS 200-662-2

NC CODE 2914 11 00

EC NO. 606 001 00 8

UN/ID NO. 1090

ADR/RID 3 F1 IMDG 3/II

R: 11-36-66-67

S: 16-26-9





min. 99.5%
max. 0.002%
max. 10
max. 0.05%
max. 0.05%
max. 5 ppm
passes test
passes test
max. 0.0003
max. 0.0006
max. 0.5%

Ultraviolet Absorbance (1.00-cm path vs water; curve smooth throughout stated range with no extraneous impurity peaks):

at 330 nm	max. 1.00
at 340 nm	max. 0.10
at 350 nm	max. 0.02
at 400 nm	max. 0.01

PRODUCT	PACKING	CONT
NO.		BOX
8001.1000	11	
8001.2500	2.5	



Aceto

Acetone

8002 'BAKER ANALYZED' / ACS

► (CH₃)₂CO

M = 58.08 g/mol

1 l = 0.79 kg

FLASHPOINT $-20 \, ^{\circ}\text{C}$ **CAS NO.** 67-64-1

EINECS 200-662-2

NC CODE 2914 11 00 EC NO. 606 001 00 8

UN/ID NO. 1090

ADR/RID 3F1 IMDG 3/II

R: 11-36-66-67

S: 16-26-9





Exceeds ACS Specifications	
Assay (by GC)	min. 99.5%
Aldehyde (as HCHO)	max. 0.002%
Color (APHA)	max. 10
Diisopropylether((CH ₃) ₂ CHOCH(CH ₃) ₂)	max. 0.05%
Isopropyl Alcohol	max. 0.05%
Methanol (CH ₃ OH)	max. 0.05%
Residue after Evaporation	max. 0.001%
Solubility in Water	passes test
Substances Reducing KMnO ₄	passes test
Titrable Acid (meq/g)	max. 0.0003
Titrable Base (meq/g)	max. 0.0006
Water (H ₂ O)	max. 0.5%
Trace Impurities (in ppm):	
Aluminium (AI)	max. 0.5
Barium (Ba)	max. 0.1
Boron (B)	max. 0.02
Cadmium (Cd)	max. 0.05
Calcium (Ca)	max. 0.5
Chromium (Cr)	max. 0.02
Cobalt (Co)	max. 0.02
Copper (Cu)	max. 0.02
Iron (Fe)	max. 0.1
Lead (Pb)	max. 0.1
Magnesium (Mg)	max. 0.1

max. 0.02 max. 0.02

max. 0.1

max. 0.1

	CONT.
	вох
11	6
2.5	4
5 EcoTainer	4
25 l	
200 l	
e metal solvent can for more safe	ty in the lab.
	2.5 l 5 l EcoTainer 25 l

For safe handling of 25 l tin cans, see Self-closing tap.

3403

Acetone

99% / HISTO GRADE

► (CH₃)₂CO

M = 58.08 g/mol

11 = 0.79 kg

FLASHPOINT −20 °C **CAS NO.** 67-64-1

EINECS 200-662-2

NC CODE 2914 11 00

EC NO. 606 001 00 8

UN/ID NO. 1090

ADR/RID 3 F1

IMDG 3/II

R: 11-36-66-67

S: 16-26-9








Manganese (Mn)

Nickel (Ni) Tin (Sn)

Zinc (Zn)



Histo-Grade implicates that this reagent is specially tested and therefore solely intended for use in histo-pathology applications. This reagent is of an analytical quality.

Acetone

e for qualities meeting the requirements of cGMP, USP/NF, BP/Ph.Eur, JP or FCC, page 36



Acetone CMOS, Finyte Grade

www.jtbaker.com/europe

Acetonitrile

'BAKER ULTRA RESI-ANALYZED' / for Organic Residue Analysis

9255

► CH₃CN

M = 41.05 g/mol

1 l = 0.78 kg

FLASHPOINT $2\,^\circ\!\mathbb{C}$

CAS NO. 75-05-8

EINECS 200-835-2

NC CODE 2926 90 95

EC NO. 608 001 00 3

UN/ID NO. 1648

ADR/RID 3F1

IMDG 3/II

R: 11-20/21/22-36

S: 16-36/37





Assay (by GC) (corrected for water)	min. 99.8%
Color (APHA)	max. 10
Residue after Evaporation	max. 1 ppm
Titrable Acid (µeq/g)	max. 0.3
Titrable base (µeq/g)	max. 0.6
Water (H ₂ O) (by Coulometry)	max. 0.05%
ECD Sensitive Impurities (as Heptach	lor Epoxide):
Single Impurity Peak (pg/ml)	max. 10
FID-Sensitive Impurities (as 2-Octano	I):
Single Impurity Peak (ng/ml)	max. 5

PRODUCT	PACKING	CONT.
NO.		вох
9255.1000	11	6

Acetonitrile

BAKER ANALYZED LC-MS Reagent

9821

► CH₃CN **M** = 41.05 g/mol **1 l** = 0.78 kg FLASHPOINT 2°C **CAS NO.** 75-05-8 EINECS 200-835-2 NC CODE 2926 90 95 EC NO. 608 001 00 3 **UN/ID NO.** 1648

ADR/RID 3 F1 IMDG 3/II

R: 11-20/21/22-36 **S:** 16-36/37





Certificate Provided Reporting Actua	I Lot Analysis
Assay (by GC)	min. 99.8%
Residue after Evaporation	max. 1 ppm
Water (H ₂ O)	max. 0.01%
LC-Gradient-Diode Array Detection (a.u.), test
solution is modified with 0.1% (v/v) f	ormic acid:
at 220 nm	max. 0.002
at 254 nm	max. 0.001
LC-MS Gradient Suitability Test (TIC	, 100 to 2000
m/z), test solution is modified with 0	.1% (v/v)
formic acid:	
Positive ESI-MS Sensitive Impurities (as	
Reserpine)	max. 50 ng/ml
Product Information (not specification	ns):
Density (g/ml) at 20°C	0.78
Refractive Index $n_{\overline{D}}^{20}$	1.344
Trace Impurities (in ppb):	
Aluminium (AI)	max. 50
Calcium (Ca)	max. 50
Iron (Fe)	max. 50
Magnesium (Mg)	max. 50
Potassium (K)	max. 50
Sodium (Na)	max. 50
Ultraviolet Absorbance (1.00-cm path	

max. 0.05

max. 0.01

max. 0.01

PRODUCT	PACKING	CON
NO.		во
9821.1000GL	1 Glass	

Element concentrations are at time of lot release

The J.T.Baker CYCLE-TAINER High Purity Solvent Delivery System, preserves purity and protects people.

See chapter 3 of this catalogue for product details.

at 200 nm

at 220 nm

at 254 nm

Aceto

Acetonitrile

9017 'BAKER HPLC ANALYZED' / Ultra Gradient HPLC Grade

► CH₃CN

M = 41.05 g/mol

1 l = 0.78 kg

FLASHPOINT $2\,^\circ\!\mathbb{C}$

CAS NO. 75-05-8

EINECS 200-835-2 NC CODE 2926 90 95

EC NO. 608 001 00 3

UN/ID NO. 1648

ADR/RID 3 F1

IMDG 3/II

R: 11-20/21/22-36

S: 16-36/37





Meets Chromatography Reagent Specifications for testing Ph. Eur. monographs min. 99.8% Assay (by GC) Appearance passes test Distilling Range 80-82°C min. 95% Identity (by IR) passes test Transmittance 240 nm - 420 nm min. 98% Water (H₂O) max. 0.01%

Fluorescence Trace Impurities (as quinine base), ppb:

Measured at 450 nm max. 0.3 Measured at Emission Maximum for

Solvent Impurities max. 1.0

Gradient Elution Test (PAH Suitability test) Ultraviolet Absorbance (a.u.):

at 210 nm max. 0.002 at 254 nm max. 0.0005

Fluorescence (in ppb, measured as Benzo(a)Pyrene at EX/Em = 280nm/Total

Emission max. 0.5 Limits of Impurities:

Residue after Evaporation max. 1 ppm Titrable Acid (meq/g) max. 0.0008 Titrable Base (meq/g) max. 0.0006

Physical Data (not specifications):

0.786 Density (g/ml) at 20°C Refractive Index n²⁰ 1.344

Ultraviolet Absorbance (1.00-cm path vs water):

at 200 nm max. 0.05 at 220 nm max. 0.01 at 254-400 nm max. 0.01 UV Cut-off, nm max. 190

PRODUCT	PACKING	CONT.
NO.		вох
9017.1000	11	6
9017.2500	2.5	4
9017.5000	5 EcoTainer	4

EcoTainer, the metal solvent can for more safety in the lab.

Filtered through a 0.2 micron filter.

Packaged under Nitrogen.

9012

Acetonitrile

'BAKER HPLC ANALYZED' / HPLC Far UV / Gradient Grade

► CH₃CN

M = 41.05 g/mol

1 l = 0.78 kg

FLASHPOINT 2°C

CAS NO. 75-05-8

EINECS 200-835-2

NC CODE 2926 90 95

EC NO. 608 001 00 3

UN/ID NO. 1648

ADR/RID 3 F1

IMDG 3/II

R: 11-20/21/22-36

S: 16-36/37





Assay (by GC)	min. 99.8%	
Residue after Evaporation	max. 5 ppm	
Titrable Acid (meq/g)	max. 0.0008	
Titrable Base (meq/g)	max. 0.0006	
Water (H ₂ O)	max. 0.02%	
Fluorescence Trace Impurities (as quinine base),		
ppb:		
at Emission Maximum for Impurities	max. 1.0	

Gradient Elution Test (PAH Suitability test) Ultraviolet Absorbance (a.u.):

max. 0.001 at 254 nm

max. 0.3

Physical Data (not specifications):

Measured at 450 nm

Refractive Index n²⁰ 1.344 Ultraviolet Absorbance (1.00-cm path vs wate

Citaviolet Absorbance (1.00-cm path vs water).		
at 200 nm	max. 0.10	
at 210 nm	max. 0.05	
at 220 nm	max. 0.03	
at 254 nm	max. 0.01	
at 280 nm	max. 0.01	
at 350 nm	max. 0.01	
at 400 nm	max. 0.01	
UV Cut-off, nm	max. 190	

PRODUCI	PACKING CONT
NO.	B0
9012.1000GL	1 Glass
9012.2500GL	2.5 l Glass
9012.5000EC	5 EcoTainer
9012.9010RC	10 Returnable Container
9012.9030RC	30 Returnable Container
9012.9200RC	200 Returnable Container
EcoTainer, the	metal solvent can for more safety in the lab

Filtered through a 0.2 micron filter. Packaged under Nitrogen.

Acetonitrile

'BAKER HPLC ANALYZED' / HPLC Isocratic Grade

8257

CLLCN	
CH ₂ CN	

M = 41.05 g/mol

1 l = 0.78 kg

FLASHPOINT 2°C

CAS NO. 75-05-8

EINECS 200-835-2

NC CODE 2926 90 95

EC NO. 608 001 00 3

UN/ID NO. 1648

ADR/RID 3F1

IMDG 3/II

R: 11-20/21/22-36

S: 16-36/37





Assay (by GC)	min. 99.8%
Residue after Evaporation	max. 0.0005%
Water (H ₂ O)	max. 0.1%
Physical Data (not specification	ns):
Refractive Index n ²⁰ _D	1.344
Ultraviolet Absorbance (1.00-cr	n path vs water):
at 200 nm	max. 0.5
at 210 nm	max. 0.3
at 220 nm	max. 0.2
at 254 nm	max. 0.05
at 280 nm	max. 0.01

value on certificate

PRODUCT	PACKING	CONT.
NO.		вох
8257.1000	11	6
8257.2500	2.5	4
8257.5000	5 EcoTainer	4
8257.9010RC	10 l Returnable Container	
8257.9030RC	30 l Returnable Container	
8257.9200RC	200 l Returnable Container	
EcoTainer, the	metal solvent can for more safety	in the lab.

Acetonitrile

Ultra Low Water / BakerDRY / Low Water Solvent / ACS

UV Cut-off, nm

9035

► CH₃CN

M = 41.05 g/mol

1 l = 0.78 kg

FLASHPOINT 2 °C

CAS NO. 75-05-8

EINECS 200-835-2

NC CODE 2926 90 95

EC NO. 608 001 00 3

UN/ID NO. 1648

ADR/RID 3 F1

IMDG 3/II

R: 11-20/21/22-36

S: 16-36/37





Meets ACS Specifications		
Assay (by GC) (corrected for water)	min. 99.5%	
Appearance	passes test	
Color (APHA)	max. 10	
Residue after Evaporation	max. 1 ppm	
Titrable Acid (µeq/g)	max. 0.8	
Titrable base (µeq/g)	max. 0.6	
Water (by KF, coulometric)	max. 10 ppm	
Product Information (not specifications):		

Product Information (not specifications):		
Boiling Point (typical)	81.6°C	
Density (g/ml) at 25°C (typical)	0.778	

PRODUCT	PACKING	CONT.
NO.		вох
9035.1000	11	

Acetonitrile

'BAKER BIO-ANALYZED' / for DNA/RNA synthesis

8134



M = 41.05 g/mol

1 l = 0.78 kg

FLASHPOINT $2\,^{\circ}\mathrm{C}$

CAS NO. 75-05-8

EINECS 200-835-2

NC CODE 2926 90 95

EC NO. 608 001 00 3 **UN/ID NO.** 1648

ADR/RID 3 F1

IMDG 3/II

R: 11-20/21/22-36

S: 16-36/37





Elugraganas Traca Impuritias (as quining basa)		
Water (H ₂ O)	max. 10 ppm	
Titrable base (μeq/g)	max. 0.6	
Titrable Acid (µeq/g)	max. 0.5	
Residue after Evaporation	max. 1 ppm	
Assay (by GC)	min. 99.9%	

Fluorescence Trace Impurities (as quinine base),

ppb:

at Emission Maximum for Impurities max. 1.0 Measured at 450 nm max. 0.3 Ultraviolet Absorbance (1.00-cm path vs water):

at 200 nm	max. 0.05
at 220 nm	max. 0.01
at 254-400 nm	max. 0.01
UV Cut-off, nm	max. 190

PRODUCT	PACKING	CONT.
NO.		вох
8134	20 l Cycletainer	
8134	50 l Cycletainer	
8134	200 l Cycletainer	

D E F G

Aceto

Acetonitrile

8144 'BAKER BIO-ANALYZED' / for DNA/RNA synthesis

► CH₃CN

M = 41.05 g/mol

1 l = 0.78 kg

FLASHPOINT $2\,^{\circ}\mathrm{C}$

CAS NO. 75-05-8

EINECS 200-835-2

NC CODE 2926 90 95 EC NO. 608 001 00 3

UN/ID NO. 1648

ADR/RID 3F1

IMDG 3/II

R: 11-20/21/22-36

S: 16-36/37





Assay (corrected for H₂O) min. 99.9% Residue after Evaporation max. 1 ppm Titrable Acid (µeq/g) max. 0.5 Titrable base $(\mu eq/g)$ max. 0.6 Water (H₂O) max. 30 ppm

Fluorescence Trace Impurities (as quinine base),

Measured at 450 nm

Exceeds ACS Specifications

Zinc (Zn)

Measured at Emission Maximum for Solvent Impurities max. 1.0

max. 0.3

max. 0.1

Ultraviolet Absorbance (1.00-cm path vs water):

at 200 nm	max. 0.05
at 220 nm	max. 0.01
at 254-400 nm	max. 0.01
UV Cut-off, nm	max. 190

PRODUCT NO.	PACKING	CONT. BOX
8144.2500	2.5	
8144.4000	4 l Glass	

8004

Acetonitrile

'BAKER ANALYZED' / ACS

► CH₃CN **M** = 41.05 g/mol

1 l = 0.78 kg

FLASHPOINT 2°C

CAS NO. 75-05-8

EINECS 200-835-2

NC CODE 2926 90 95

EC NO. 608 001 00 3

UN/ID NO. 1648

ADR/RID 3 F1 IMDG 3/II

R: 11-20/21/22-36

S: 16-36/37





=noodad nood opcomodations	
Assay (by GC)	min. 99.5%
Appearance	passes test
Color (APHA)	max. 10
Residue after Evaporation	max. 0.005%
Titrable Acid (µeq/g)	max. 8
Titrable base (µeq/g)	max. 0.6
Water (H ₂ O)	max. 0.3%
Trace Impurities (in ppm):	
Aluminium (AI)	max. 0.5
Barium (Ba)	max. 0.1
Boron (B)	max. 0.02
Cadmium (Cd)	max. 0.05
Calcium (Ca)	max. 0.5
Chromium (Cr)	max. 0.02
Cobalt (Co)	max. 0.02
Copper (Cu)	max. 0.02
Iron (Fe)	max. 0.1
Lead (Pb)	max. 0.1
Magnesium (Mg)	max. 0.1
Manganese (Mn)	max. 0.02
Nickel (Ni)	max. 0.02
Tin (Sn)	max. 0.1

PRODUCT	PACKING	CONT.
NO.		вох
8004.1000	11	6
8004.2500	2.5	4
8004.9200	200 l	

Acetonitrile

2680

► CH₃CN

M = 41.05 g/mol

1 l = 0.78 kg

FLASHPOINT 2°C **CAS NO.** 75-05-8

EINECS 200-835-2

NC CODE 2926 90 95

EC NO. 608 001 00 3

UN/ID NO. 1648

ADR/RID 3 F1

IMDG 3/II

R: 11-20/21/22-36

S: 16-36/37





Specifications min. 99.8% Assay passes test Appearance Distillation Range 80-82°C min. 95% Transmittance 240 nm - 420 nm min. 98%

Meets Ph. Eur. Chromatography Reagent

PRODUCT	PACKING	CONT.
NO.		вох
2680.2500GL	2.5 Glass	

Acetylcysteine

e for qualities meeting the requirements of cGMP, USP/NF, BP/Ph.Eur, JP or FCC, page 36



Acetylene Tetrabromide

See 1,1,2,2-Tetrabromoethane



Acid Blue 90

See COOMASSIE Brilliant Blue R-250




1420

Acid Fuchsin

Identification (by IR) C₂₀H₁₇N₃Na₂O₉S₃ passes test **PRODUCT PACKING** CONT. NO.

M = 585.54 g/mol **CAS NO.** 3244-88-0 **EINECS** 221-816-5 NC CODE 3204 12 00



Acids for Trace Analysis

See for detailed information section Re agents for trace element analysis, page 21



Acrylamide



BOX

	BAKER ULTRAPURE BIOREAGEN	II' / Purified for Electrophoresis				4081
► CH ₂ :CHCONH	I_2	Assay	min. 99.9%	PRODUCT	PACKING	
M =	71.08 g/mol	Absorbance of a 1% Solution at 290 nm		NO.		
CAS NO.	79-06-1	(1-cm path vs water)(a.u.)	max. 0.1	4081.0100	100 g	
EINECS	201-173-7	Acrylic acid (CH ₂ :CHCOOH)	max. 0.001%	4081.0500	500 g	
NC CODE	2924 19 00	Conductivity of 35% Solution	max. 2.5 μmho	4081.2500	2.5 kg	
EC NO.	616 003 00 0	DNase Activity	none detected	4081.9012	12 kg	
UN/ID NO.	2074	Heavy Metals (as Pb)	max. 1 ppm			
ADR/RID	6.1 T2	Insoluble in Methanol	max. 0.005%	Store between	en 18 and 26°C.	
IMDG	6.1/III	Insoluble in Water	max. 0.005%			
R:	20/21-25-36/38-43-45-46-	Iron (Fe)	max. 1 ppm			
	48/23/24/25-62	pH of 10% Solution in 0.1 M NaCl	6.0-7.0			
S:	45-53	Protease Activity	none detected			

none detected

max. 0.5%

max. 0.01%

Acrylamide Gels for DNA Electrophoresis

RNase Activity



Activated Carbon

Acid-Washed, Steam-Activated / Powder / (DARCO G-60)

E343

▶ C	
M =	12.01 g/mol
CAS NO.	7440-44-0
EINECS	231-153-3
NC CODE	3802 10 00
UN/ID NO.	1362
ADR/RID	4.2 S2
IMDG	4.2/III

IMDG 4.2/III

▶ C

PRODUCI	PACKING	CONI.
NO.		вох
E343.0500	500 g	

Activated Carbon

max. 0.01% Chloride (CI) Heavy Metals (as Pb) **M** = 12.01 g/mol max. 0.003% **CAS NO.** 7440-44-0 Iron (Fe) max. 0.03% **EINECS** 231-153-3 Loss on Drying at 105°C max.12% NC CODE 3802 10 00 Residue after Ignition (600°C) max. 1% **UN/ID NO.** 1362 Substances soluble in Ethanol max. 0.2% **ADR/RID** 4.2 S2 Substances soluble in Hydrochloric Acid max. 1.0%

Sulfate (SO₄)

Substances soluble in Water

NO.	PACKING	CONT. BOX
1991.0500	500 g	
1991.0500	500 g	

Activ

Activator DCI

9479 BAKER ANALYZED' / for DNA/RNA synthesis

11 = 0.80 kg FLASHPOINT 5 °C NC CODE 2926 90 95 **UN/ID NO.** 1648 ADR/RID 3 F1

Appearance: clear, colorless solution free from visible particulates passes test Molarity (Dicyano-imidazole) 0.24-0.26 Water (H₂O) max. 30 ppm IMDG 3/II

Suitable for Oligonucleotide Synthesis

PRODUCT CONT. **PACKING** NO. вох 9479.0450 450 ml Glass 9479.2500 2.5 | Glass

S: 16-36/37





Activator ETT

R: 11-19-20/21/22-37/38-41

'BAKER ANALYZED' / for DNA/RNA synthesis

FLASHPOINT 5 °C **CAS NO.** 75-05-8 NC CODE 2926 90 95 **UN/ID NO.** 1648 ADR/RID 3 F1 IMDG 3/II

9478

R: 11-20/21/22-36 **S:** 16-26-33-36/37





Suitable for Oligonucleotide Synthesis Appearance: clear, colorless solution free from visible particulates passes test Molarity (5-Ethylthio-1H-Tetrazole) 0.28-0.32 Water (H₂O) max. 50 ppm

PRODUCT	PACKING	CONT.
NO.		вох
9478.0200	200 ml Glass	
9478.0450	450 ml Glass	
9478.2000GL	2 l Glass	

Activator reagents for use in DNA synthesis

e for detailed information section Reagents for DNA/RNA Synthesis, page 261

1892

Agar-Agar der / BAKER

CAS NO. 9002-18-0 EINECS 232-658-1 NC CODE 1302 31 00

Foreign Insoluble Matter max. 1% Foreign Organic Matter max. 1% Insoluble in HCI max. 0.5% Total Ash max. 5% Water (H₂O) max. 20%

PRODUCT	PACKING	CONT.
NO.		вох
1892.0250	250 g	
1892.1000	1 kg	

A426

Agarose, standard

smosis (EEO) / 'BAKER ULTRAPURE BIOREAGENT' / for electrophoresis

CAS NO. 9012-36-6 **EINECS** 232-731-8 NC CODE 3913 90 00 DNase Activity none detected Electroendosmosis (EEO)(-Mr) max. 0.13 min. 1200 Gel strength (g/cm 2)(MCI, 1.0%) Gelling Temperature of 1.5% (w/w) 34.5-37.5°C Loss on Drying at 70°C (in Vacuo) max. 7% Melting Temperature of 1.5% (w/w) solution 86.5-89.5°C Protease Activity none detected RNase Activity none detected Sulfate (SO₄) max. 0.2%

PRODUCT	PACKING	CONT
NO.		ВОХ
A426.0100	100 g	
A426.0500	500 g	
A426.0500	500 g	



L-Alanine

See for qualities meeting the requirements of cGMP, USP/NF, BP/Ph.Eur, JP or FCC, page 36

1002

Albumin, Egg

NC CODE 3502 90 90	PRODUCT	PACKING	CONT.
	NO.		вох
	1002.0500	500 g	

Albumin (Bovine) Fraction V

0604 'BAKER ANALYZED CAS NO. 9048-46-8 Total Protein (based on N) PRODUCT PACKING CONT. min. 95% EINECS 232-936-2 Water (H₂O) max. 5% NO. вох NC CODE 3502 90 70 0604.0010 10 g **EC NO.** 90 604 29 8 **Alcohol Absolute ALEG 310 Positive Resist Stripper**

ALEG 625 Positive Resist Stripper

See for detailed information section Microelectronic materials page 32 or the website www.jtbaker.com/micro

Nitrogen Compounds (as N)

See for detailed information section Microelectronic materials page 32 or the website www.jtbaker.com/micro



Alizarin Red S

1471 BAKER ANALYZED ► C₁₄H₇NaO₇S **Visual Transition Interval: PRODUCT PACKING** CONT. **M =** 342.26 g/mol pH 4.0 yellow BOX

CAS NO. 130-22-3 **EINECS** 204-981-8 NC CODE 2914 70 00

pH 6.0





Alum Iron

See Ammonium Iron(III) Sulfate Dodecahydrate




Aluminium

210 Sheets (6 x 6 x 0.001) / Foil 150x150x0.025 mm / 'BAKER'



► Al **M** = 26.98 g/mol **CAS NO.** 7429-90-5

EINECS 231-072-3 NC CODE 7606 11 90



9502.0450 450 g

1003.0500

PACKING

500 g

PACKING

25 g Glass

1471.0025

C.I. 58005.

red

CONT. BOX

Aluminium Shot / 'BAKER'

PACKING CONT. min. 95% **PRODUCT** Assay NO. BOX

M = 26.98 g/mol **CAS NO.** 7429-90-5 **EINECS** 231-072-3 NC CODE 7601 10 00

Aluminium 1000 μg/ml

(Matrix: 2% nitric acid) / 'BAKER INSTRA-ANALYZED' / Plasma Standard

5701

1003

► Al

► Al

M = 27.00 g/mol

NC CODE 3822 00 00

R: 36/38 **S**: 26



Certificate Provided Reporting Actual Lot Analysis Aluminium (AI) 998-1002 μg/ml

PRODUCT NO.

5701.0100 100 ml CONT.

Prepared from the highest purity raw material available, generally greater than 99.999% spectral purity. The content of the solution is confirmed to be accurate to within \pm 0.2%. Typically 1000 $\mu\text{g/ml.}$

The certificate of analysis provided reports actual lot analysis. The certificate also lists the trace impurities.

Alumi

Aluminium 1000 $\mu g/ml$ (Matrix: 2% nitric acid) / 'BAKER INSTRA-ANALYZED' / Atomic Absorption Standard 6917

Aluminium (AI)

► Al **M** = 27.00 g/mol

NC CODE 3822 00 00

R: 36/38

S: 26



PRODUCT 998-1002 μg/ml NO.

> 6917.0100 100 ml 6917.0500 500 ml

PACKING

Prepared by dissolution of high purity raw materials (min. 99.99% spectral purity). Assays are verified by ICP against standards traceable to NIST. Standard Reference Material numbers (SRM) are printed on each label.

CONT.

CONT

вох

вох

Aluminium 1000 μ g/ml 'BAKER ANALYZED' / Atomic Absorp

6801

► Al **M** = 27.00 g/mol

> NC CODE 3822 00 00 R: 36/38

> > S: 26



998-1002 μg/ml Aluminium (AI)

PRODUCT NO.

6801.0100 100 ml 6801.0500 500 ml

PACKING

Aluminium nitrate in nitric acid 0.5 mol/l.

Aluminium 10000 μ g/ml (Matrix: 2% nitric acid) / 'BAKER INSTRA-ANALYZED' / Plasma Standard

► Al **M** = 27.00 g/mol

> NC CODE 3822 00 00 R: 36/38

> > **S**: 26



Certificate Provided Reporting Actual Lot Analysis Aluminium (AI)

9980-10020 µg/ml

PRODUCT PACKING CONT. NO. BOX 5716.0100 100 ml

Prepared from the highest purity raw material available, generally greater than 99.999% spectral purity. The content of the solution is confirmed to be accurate to within 0.2 %. Typically 10000 μ g/ml. The certificate of analysis provided reports actual lot

analysis. The certificate also lists the trace impurities.

Aluminium Hydroxide

0005 BAKER ANALYZED

► Al(OH)₃ **M** = 78.00 g/mol

> **CAS NO.** 21645-51-2 **EINECS** 244-492-7

NC CODE 2818 30 00

Assay (by EDTA titrn.) min. 98.0% Chloride (CI) max. 0.002% Heavy Metals (as Pb) max. 0.005% max. 0.005% Iron (Fe) Sulfate (SO₄) max. 0.005%

PRODUCT PACKING CONT. вох 0005.0500 500 g

Aluminium Nitrate Nonahydrate 0006

'BAKER ANALYZED' / ACS

► AI(NO₃)₃.9H₂O

M = 375.13 g/mol

CAS NO. 7784-27-2

EINECS 236-751-8 NC CODE 2834 29 80

UN/ID NO. 1438

ADR/RID 5.1 O2

IMDG 5.1/III

R: 36/38-8 **S:** 17-26



Meets ACS Specifications	
Assay (by EDTA titrn.)	98.0-102.0%
Calcium (Ca)	max. 0.005%
Chloride (CI)	max. 0.001%
Heavy Metals (as Pb)	max. 0.001%
Insoluble Matter	max. 0.005%
Iron (Fe)	max. 0.002%
Magnesium (Mg)	max. 0.001%
pH of 5% Solution at 25°C	2.5-3.5
Potassium (K)	max. 0.002%
Sodium (Na)	max. 0.005%
Sulfate (SO ₄)	max. 0.005%

PRODUCT	PACKING	CONT.
NO.		вох
0006.0500	500 g	
0006.7300	300 lbs	

Alumi



Aluminium Oxide

 'BAKER ANALYZED'

 ► Al₂O₃
 Assay
 min. 98.5%
 PRODUCT
 PACKING
 CONTROL OF THE PACKING

 M = 101.96 g/mol
 Chloride (Cl)
 max. 0.005%
 NO.
 B

 M =
 101.96 g/mol
 Chloride (Cl)
 max. 0.005%

 CAS NO.
 1344-28-1
 Heavy Metals (as Pb)
 max. 0.001%

 EINECS
 215-691-6
 Iron (Fe)
 max. 0.02%

 NC CODE
 2818 20 00
 Loss on Ignition
 max. 1.0%

 Sulfate (SO₄)
 max. 0.005%

6	PRODUCT	PACKING	CONT.
6	NO.		вох
6	0007.0500	500 g	
6			
6			

Aluminium Oxide

'BAKER ANALYZED' / for Chromatography

8000

► Al ₂ O ₃		1
M =	101.96 g/mol	ı
CAS NO.	1344-28-1	ĺ
EINECS	215-691-6	,
NC CODE	2818 20 00	ı

Adsorption (mg o-nitroaniline/g)	min. 1.0
pH of 5% Slurry at 20°C	7.0-7.8
Physical Data (not specifications):	
Average Particle Diameter, µm (APD)	50-200
Bulk Density (g/cc)(typical)	0.9
Mean Pore Diameter, A	60

PRODUCT	PACKING	CONT.
NO.		вох
0008.0500	500 g	
Neutral, Bro	ckmann Activity Grade I.	

Aluminium Oxide

powder / 'BAKER ANALYZED' / for Chromatography

0537

Al_2U_3	
M =	101.96 g/mol
CAS NO.	1344-28-1
EINECS	215-691-6
NC CODE	2818 20 00

Adsorption (mg o-nitroaniline/g)	min. 0.05
Loss on Drying	max. 5%
pH of 5% Slurry at 20°C	7.0-7.8
Physical Data (not specifications):	
Average Particle Diameter, µm (APD)	50-200
Bulk Density (g/cc)(typical)	0.9
Mean Pore Diameter, A	60

PRODUCT	PACKING	CONT.
NO.		вох
0537.0500	500 g	

Aluminium Oxide

'BAKER ANALYZED' / for Chromatography

1848

► Al ₂ O ₃	
M =	101.96 g/mol
CAS NO.	1344-28-1
EINECS	215-691-6
NC CODE	2818 20 00

Adsorption (mg o-nitroaniline/g)	min 1.0
pH of 3% Slurry at 20°C	8.5-10.0
Physical Data (not specifications):	
Average Particle Diameter, µm (APD)	50-200
Bulk Density (g/cc)(typical)	0.9
Mean Pore Diameter, A	60

	вох
500 g	
ann Activity Grade I.	

Aluminium Oxide

'BAKER ANALYZED' / for Chromatography

1850

► Al ₂ O ₃	
M =	101.96 g/mol
CAS NO.	1344-28-1
EINECS	215-691-6
NC CODE	2818 20 00

PRODUCI	PACKING	CONT.
NO.		вох
1850.0500	500 g	
Acid, Brocki	mann Activity Grade I.	

Aluminium Oxide IB-F

Flexible TLC Sheets, 20 x 20 cm / 'BAKER-FLEX'

5006

PRODUCT	PACKING	CONT
NO.		ВОХ
5006	25 sheets	

A flexible sheet coated with high purity aluminium oxide containing an inert binder and a fluorescent indicator (activated at 2540 Å). Neither the binder nor the indicator chars on sulfuric acid heat-treatment.

Alumi

Aluminium Sulfate n-Hydrate

1889

Al₂(SO₄)₃.nH₂O **CAS NO.** 17927-65-0 **EINECS** 233-135-0

NC CODE 2833 22 00

R: 41

S: 26-39



Assay (Al ₂ (SO ₄) ₃)	51.0-59.0%
Alkali and Alkaline-earth Metals	max. 0.4%
Ammonium (NH ₄)	max. 500 ppm
Appearance of solution	passes test
Heavy Metals (as Pb)	max. 50 ppm
Identification	passes test
Iron (Fe)	max. 100 ppm
рН	2.5-4.0

PACKING	CONT.
	вох
1 kg	
50 kg	
	1 kg

Stored in an airtight container.



Aluminum Sulfate, n-Hydrate

ee for qualities meeting the requirements of cGMP, USP/NF, BP/Ph.Eur, JP or FCC, page 36

Aluminium Sulfate Octadecahydrate

'BAKER ANALYZED' / ACS

► Al₂(SO₄)₃. 18H₂O **M** = 666.42 g/mol **CAS NO.** 7784-31-8 EINECS 233-135-0 NC CODE 2833 22 00 **R:** 41 **S:** 26-39

100 - 110% Assav Ammonium (NH₄) max. 0.003% No. Appearance of solution max. 0.0005% 0010.7200 Arsenic (As) Chloride (CI) max. 0.005% Heavy Metals (as Pb) max. 0.002% Iron (Fe) max. 0.005% pH (2%, 20°C) 3 - 4

PRODUCT PACKING CONT. BOX passes test 0010.0500 500 g 200 lbs



Amidosulfonic Acid

See Sulfamic Acid



Aminoacetic Acid



p-Aminobenzenesulfonamide



p-Aminobenzenesulfonic Acid

e Sulfanilic Acid Anhydrous



2-Aminoethanesulfonic Acid



2-Aminoethanol



2-Amino-2-methyl-1-propanol

'BAKER ANALYZED'

► CH₃C(NH₂(CH₃)CH₂OH **M** = 89.14 g/mol

> **11=** 0.93 kg FLASHPOINT 67 °C

CAS NO. 124-68-5 **EINECS** 204-709-8

NC CODE 2922 1980 EC NO. 603 070 00 6

> R: 36/38-52/53 **S:** 61



Assay	93-97%
Appearance	passes test
Boiling Point	165°C
Freezing Point	24-28°C
pH of 0.1 M Solution at 25°C	11-12
Water (H ₂ O)	max. 1.0%

i koboci	1 ACKING	COIVI.
NO.		вох
3321.0500	500 g	
3321.1000	1 kg	
3321.2500	2.5 kg	



Ammonia Solutions

See for qualities meeting the requirements of cGMP, USP/NF, BP/Ph.Eur, JP or FCC, page 36



0390

Ammonium Acetate

'BAKER HPLC ANALYZED' / for use in High Performance Liquid Chromatography

► CH₃COONH₄

M = 77.08 g/mol CAS NO. 631-61-8 EINECS 211-162-9 NC CODE 2915 29 00 $\begin{tabular}{ll} Insoluble Matter & max. 0.005\% \\ Nitrate (NO_3) & max. 0.001\% \\ pH of 5\% Solution at 25°C & 6.7-7.3 \\ Residue after Ignition & max. 0.01\% \\ Sulfate (SO_4) & max. 0.001\% \\ \hline \textbf{Trace Impurities (in ppm):} \\ \end{tabular}$

Trace Impurities (in ppm):

Chloride (CI) max. 5

Heavy Metals (as Pb) max. 5

Iron (Fe) max. 5

Ultraviolet Absorbance (1M aqueous Solution; 1.00-cm path vs water):

at 254 nm max. 0.02 at 280 nm max. 0.01 at 350 nm max. 0.01

PRODUCT	PACKING	CONT.
NO.		вох
0390 1000	1 ka	

Ammonium Acetate

'BAKER ANALYZED' / ACS

0011

► CH₃COONH₄

M = 77.08 g/molCAS NO. 631-61-8EINECS 211-162-9NC CODE 2915 29 00

Meets ACS Specifications Assay min. 97% Insoluble Matter max. 0.005% Nitrate (NO₃) max. 0.001% pH of 5% Solution at 25°C 6.7-7.3 Residue after Ignition max. 0.01% Sulfate (SO₄) max. 0.001% Trace Impurities (in ppm): Chloride (CI) max. 5 Heavy Metals (as Pb) max. 5 Iron (Fe) max. 5

PRODUCT NO.	PACKING	CONT. BOX
0011.1000	1 kg	6
0011.9025	25 kg	
0011.9050	50 kg	

Ammonium Acetate

10 mM (pH 7) / BAKER ANALYZED LC-MS Reagent

9829

NC CODE 3822 00 00

9.5-10.5 mM Molarity pH at 20°C 7.00 ± 0.05 Product Information (not specifications): Density (g/ml) at 20°C 0.998 Trace Impurities (in ppb): Aluminium (AI) max. 50 Calcium (Ca) max. 50 Iron (Fe) max. 50 Magnesium (Mg) max. 50 Potassium (K) max. 50 Sodium (Na) max. 50

Ultraviolet Absorbance (1.00-cm path vs water):

max. 0.02

max. 0.01

at 254 nm

at 280 nm

Certificate Provided Reporting Actual Lot Analysis

PRODUCT	PACKING	CONT.
NO.		вох
9829.1000GL	1 Glass	6

Element concentrations are at time of lot release.

Ammonium Bicarbonate

See Ammonium Hydrogen Carbonate



Ammonium Bifluoride

See Ammonium Hydrogen Fluoride



Ammonium Biphosphate

See Ammonium Dihydrogen Phosphate



Ammonium Carbonate

0391 BAKER HPLC ANALYZED' / for use in High Performance Liquid Chromatography

CAS NO. 10361-29-2 **EINECS** 233-786-0 NC CODE 2836 10 00 **R:** 22-36/37 **S:** 26-36-47



min. 30.0% Assay (NH₃) Insoluble Matter max. 0.005% Non-volatile Matter max. 0.01% Sulfur Compounds (as SO₄) max. 0.002% Trace Impurities (in ppm):

Chloride (CI) max. 5 Heavy Metals (as Pb) max. 5 Iron (Fe) max. 3 Ultraviolet Absorbance (10% (w/v) aqueous

solution; 1.00-cm path vs distilled water; curve smooth throughout stated range with no extraneous impurity peaks):

at 254 nm	max. 0.02
at 280 nm	max. 0.01
at 350 nm	max. 0.01

PRODUCT NO.	PACKING	CONT BO
0391.1000	1 kg	
	_	

Ammonium Carbonate

0015

CAS NO. 10361-29-2 **EINECS** 233-786-0 NC CODE 2836 10 00 **R:** 22-36/37 **S:** 26-36-46 harmful

Meets ACS Specifications Assay (NH₃) min. 30.0% Additional Specification(s): Insoluble Matter max. 0.005% Non-volatile Matter max. 0.01% Sulfur Compounds (as SO₄) max. 0.002% Trace Impurities (in ppm): Chloride (CI) max. 5 Heavy Metals (as Pb) max. 5 Iron (Fe) max. 5

PRODUCT	PACKING	CONT
NO.		B0)
0015.1000	1 kg	
001311000	9	

0017

Ammonium Carbonate

CAS NO. 10361-29-2 EINECS 233-786-0 NC CODE 2836 10 00 R: 22-36/37 **S:** 26-36-46



Assay (NH₃) 30.0-34% Appearance passes test Chloride (CI) max. 0.002% Heavy Metals (as Pb) max. 0.001% Insoluble Matter max. 0.01% max. 0.1% Non-volatile Matter Potassium (K) max. 1 ppm Sodium (Na) max. 3 ppm Sulfur Compounds (as SO₄) max. 0.005%

PACKING	CONT.
	BOX
25 kg	



Ammonium Carbonate

e for qualities meeting the requirements of cGMP, USP/NF, BP/Ph.Eur, JP or FCC, page 36

Ammonium Cerium(IV) Nitrate 'BAKER ANALYZED' / ACS

► (NH₄)₂Ce(NO₃)₆

M = 548.23 g/mol **CAS NO.** 16774-21-3 EINECS 240-827-6 NC CODE 2846 10 00 **UN/ID NO.** 1477

ADR/RID 5.1 02 IMDG 5.1/II **R:** 41-8

S: 17-26-39



Meets ACS Specifications. Meets Reagent Specifications for testing USP/NF monographs Assay

min. 98.5% max. 0.01% Chloride (CI) Insoluble in Dilute H₂SO₄ max. 0.05% Iron (Fe) max. 0.005% Phosphate (PO₄) max. 0.02%

PRODUCT	PACKING	CONT.
NO.		вох
1050.0500GL	500 g Glass	

Ammonium Chloride

'BAKER ANALYZED' / ACS

0018

NH₄CI	
M =	53.49 g/mol
CAS NO.	12125-02-9
EINECS	235-186-4
NC CODE	2827 10 00
EC NO.	17 014 00 8
R:	22-36
S:	22
	Xn harmful

Meets ACS Specifications. Meets Reagent		
Specifications for testing USP/NF monographs		
Assay (argentometric titrn.)	min. 99.5%	
Calcium (Ca)	max. 0.001%	
Insoluble Matter	max. 0.005%	
pH of 5% Solution at 25°C	4.5-5.5	
Residue after Ignition	max. 0.01%	
Sulfate (SO ₄)	max. 0.002%	
Trace Impurities (in ppm):		

max. 5

max. 2

max. 5


max. 2

max. 150 ppm

PACKING	CONT.
	вох
500 g	6
1 kg	6
	500 g 1 kg

Ammonium Chloride

0019



Assay (dry basis)	99.5-100.5%
Acid or alkaline reacting substances	passes tes
Appearance of solution	passes tes
Bromide/lodide	passes tes
Calcium (Ca)	max. 200 ppm
Heavy Metals (as Pb)	max. 10 ppm
Identification	passes tes
Iron (Fe)	max. 20 ppm
Limit of thiocyanate (SCN)	passes tes
Loss on Drying	max. 0.5%
pH of 5% Solution at 25°C	4.6-6.0
Residue on Ignition	max 0.1%

PRODUCT	PACKING	CONT.
NO.		вох
0019.1000	1 kg	6
0019.5000	5 kg	4
0019.9050	50 kg	

Ammonium Chloride

See for qualities meeting the requirements of cGMP, USP/NF, BP/Ph.Eur, JP or FCC, page 36

Sulfate (SO₄)

Heavy Metals (as Pb)

Magnesium (Mg)

Phosphate (PO₄)

Iron (Fe)



Ammonium Citrate, Dibasic

See Ammonium Hydrogen Citrate



Ammonium Dichromate

'BAKER ANALYZED' / ACS

u	w	Z	ш	

(NH ₄) ₂ Cr ₂ O ₇	
M =	252.06 g/mol
CAS NO.	7789-09-5
EINECS	232-143-1
NC CODE	2841 50 00
EC NO.	24 003 00 1
UN/ID NO.	1439
ADR/RID	5.1 02
IMDG	5.1/II
R:	2-21-25-26-34-42/43-45-46
	48/23-50/53-60-61-8
ç.	45-53-60-61

Exceeds ACS Specifications	
Assay (dried)	min. 99.5%
Calcium (Ca)	max. 0.002%
Chloride (CI)	max. 0.005%
Insoluble Matter and NH ₄ OH Precipitate	max. 0.005%
Iron (Fe)	max. 0.002%
Loss on Drying at 105°C	max. 3.0%
pH of 5% Solution at 25°C	3.5-5.0
Sodium (Na)	max. 0.005%
Sulfate (SO ₄)	max. 0.005%

PRODUCI	PACKING	CONT.
NO.		вох
0021.0100	100 g	
0021.9050	50 kg	











Ammonium Dihydrogen Phosphate

4931 ULTREX Ultrapure Reagent

► NH₄H₂PO₄

M = 115.03 g/molCAS NO. 7722-76-1EINECS 231-764-5NC CODE 3105 40 00

Certificate Provided Reporting Actual Lot Analysis Actual Lot Analysis Lot. No. B08H23

Assay (NH₄H₂PO₄) 99.9% Particulate Matter < 0.001% pH of 5% Solution at 25°C 4.2 Metallic Impurities in parts per million ($\mu g/g$): Aluminium (AI) < 0.2 Barium (Ba) < 0.2 Bismuth (Bi) < 2.2 0.1 0.2 0.2 0.4

Cadmium (Cd) Calcium (Ca) Chromium (Cr) Cobalt (Co) Copper (Cu) 2 < 0.2 Iron (Fe) Lead (Pb) < 0.2 Magnesium (Mg) 0.1 Manganese (Mn) 0.1 Mercury (Hg) 7 Molybdenum (Mo) 0.1 Nickel (Ni) < 0.2 Niobium (Nb) 0.1 Potassium (K) 0.6 Silver (Ag) < 0.3 Sodium (Na) 0.3 Strontium (Sr) 0.2

 Vanadium (V)
 0.1

 Zinc (Zn)
 < 0.2</td>

 Zirconium (Zr)
 < 0.2</td>

 Non-Metallic Impurities in parts per million (μg/g):

 Arsenic (As)
 < 0.8</td>

 Fluoride (F)
 14.3

0.2

< 0.2

< 0.0005

min. 98.0%

max. 0.02

max. 0.01

PRODUCT

PACKING

 $\begin{tabular}{lll} Nitrate (NO_3) & < 0.0002 \\ Silicon (Si) & 0.3 \\ Sulfur Compounds (as <math>SO_4)$ & < 0.005 \\ \end{tabular}

Ultraviolet Absorbance (1.00-cm path vs water): at 215 nm at 220 nm

 at 220 nm
 0.01

 at 254 nm
 0.007

 at 260 nm
 0.007

 at 280 nm
 0.006

0.01

CONT.

 PRODUCT
 PACKING
 CONT.

 NO.
 BOX

 4931.0050
 50 g

 4931.0500
 500 g

For Laboratory, Research or Manufacturing Use.

Ammonium Dihydrogen Phosphate

crystal / 'BAKER HPLC ANALYZED' / for use in High Performance Liquid Chromatography

Assay

at 280 nm

at 350 nm

Halide (as CI)

Tin (Sn)

Titanium (Ti)

0392 ► NH₄H₂PO₄

M = 115.03 g/molCAS NO. 7722-76-1EINECS 231-764-5NC CODE 3105 40 00

Ammonium Hydroxide Precipitate max. 0.005% Insoluble Matter max. 0.005% max. 0.001% Iron (Fe) Nitrate (NO₃) max. 0.001% pH of 5% Solution at 25°C 3.8-4.4 Potassium (K) max. 0.005% max. 0.005% Sodium (Na) Sulfur Compounds (as SO₄) max. 0.005% Trace Impurities (in ppm): Arsenic (As) max. 0.5 Chloride (CI) max. 5 Heavy Metals (as Pb) max. 5 Ultraviolet Absorbance (1.00-cm path vs water): at 254 nm max. 0.03

NO.		BOX
0392.1000	1 kg	

Ammonium Dihydrogen Phosphate 'BAKER ANALYZED' / ACS

0029

$\sim NH_4H_2PO_4$	
M =	115.03 g/mo
CAS NO.	7722-76-1
EINECS	231-764-5
NC CODE	3105 40 00


Meets ACS Specifications. Meets Reagent Specifications for testing USP/NF monographs

Assay	min. 98.0%
Calcium (Ca)	max. 0.001%
Insoluble Matter	max. 0.005%
Iron (Fe)	max. 0.001%
Magnesium (Mg)	max. 0.0005%
Nitrate (NO ₃)	max. 0.001%
pH of 5% Solution at 25°C	3.8-4.4
Potassium (K)	max. 0.005%
Sodium (Na)	max. 0.005%
Sulfate (SO ₄)	max. 0.01%
Trace Impurities (in ppm):	
Chloride (CI)	max. 5
Heavy Metals (as Pb)	max. 5

PRODUCT	PACKING	CONT.
NO.		вох
0029.1000	1 kg	6
0029.9010	10 kg	

Ammonium Fluoride

0023



S: 1/2-26-45

Meets ACS Specifications	
Assay	min. 98.0%
Chloride (CI)	max. 0.001%
Heavy Metals (as Pb)	max. 5 ppm
Insoluble Matter	max. 0.005%
Iron (Fe)	max. 5 ppm
Residue after Ignition	max. 0.01%
Sulfate (SO ₄)	max. 0.005%

PACKING	CONT.
	вох
250 g	6

Ammonium Fluoride 40% MOS, VLSI Grade



0024

1890

Ammonium Heptamolybdate Tetrahydrate

'BAKER ANALYZED' / ACS

► (NH ₄) ₆ Mo ₇ O ₂₄	.4H ₂ O
M =	1235.86 g/mol
CAS NO.	12054-85-2
EINECS	234-722-4
NC CODE	2841 70 00

Meets ACS Specifications	
Assay (as MoO ₃)	81.0-83.0%
Arsenate, Phosphate and Silicate (as SiO ₂)	max. 0.001%
Chloride (CI)	max. 0.002%
Heavy Metals (as Pb)	max. 0.001%
Insoluble Matter	max. 0.005%
Magnesium (Mg)	max. 0.005%
Nitrate (NO ₃)	passes test
Potassium (K)	max. 0.01%
Sodium (Na)	max. 0.01%
Sulfate (SO ₄)	max. 0.02%
Trace Impurities (in ppm):	
Phosphate (PO ₄)	max. 5

PRODUCT	PACKING	CONT.
NO.		вох
0024.0250	250 g	6
0024.1000	1 kg	6
0024.5000	5 kg	
0024.9050	50 kg	

Ammonium Heptamolybdate Tetrahydrate

► (NH₄)₆Mo₇O₂₄.4H₂O

4/6.0.07024	
M =	1235.86 g/mol
CAS NO.	12054-85-2
EINECS	234-722-4
NC CODE	2841 70 00

Assay (as MoO ₃)	81.0-83.0%
Chloride (CI)	max. 0.005%
Heavy Metals (as Pb)	max. 0.003%
Insoluble Matter	max. 0.01%
Nitrate (NO ₃)	max. 0.01%
Phosphate (PO ₄)	max. 0.001%
Sulfate (SO ₄)	max. 0.05%

PRODUCT	PACKING	CONT.
NO.		вох
1890.9025	25 kg	
1890.9050	50 kg	

D E F G

Ammon

Ammonium Hydrogen Carbonate 'BAKER ANALYZED'

Iron (Fe)

0013

► NH₄HCO₃

M = 79.06 g/mol

CAS NO. 1066-33-7 **EINECS** 213-911-5

NC CODE 2836 10 00

R: 22-36/37

S: 26-36-46

Assay	99.0-101.0%
Heavy Metals (as Pb)	max. 0.001%
Insoluble Matter	max. 0.005%
pH of 5% Solution at 25°C	7.0-8.0
Residue after Ignition	max. 0.005%
Sulfate (SO ₄)	max. 0.001%
Trace Impurities (in ppm):	
Chloride (CI)	max. 5

max. 5

PRODUCT	PACKING	CONT.
NO.		вох
0013.1000	1 kg	
0013.9025	25 kg	
0013.9050	50 kg	

Ammonium Hydrogen Citrate

0020 'BAKER ANALYZED' / ACS

► (NH₄)₂HC₆H₅O₇

M = 226.19 g/mol **CAS NO.** 3012-65-5

EINECS 221-146-3

NC CODE 2918 15 00

R: 36



Meets ACS Specifications	
Assay (by formol method)	98.0-103.0%
Chloride (CI)	max. 0.001%
Insoluble Matter	max. 0.005%
Iron (Fe)	max. 0.001%
Oxalate (C ₂ O ₄)	passes test
Residue after Ignition	max. 0.01%
Sulfur Compounds (as SO ₄)	max. 0.005%
Trace Impurities (in ppm):	
Heavy Metals (as Pb)	max. 5
Phosphate (PO ₄)	max. 5

10.	ВО
020 1000 1 kg	
020.1000 1 kg	
020.9025 25 kg	

Ammonium Hydrogen Fluoride

0014 'BAKER ANALYZED'

NH₄F.HF

M = 57.04 g/mol

CAS NO. 1341-49-7

EINECS 215-676-4

NC CODE 2826 11 00

EC NO. 9 009 00 4 **UN/ID NO.** 1727

ADR/RID 8 C2

IMDG 8/II

R: 25-34

S: 22-26-37-45





Assay	min. 96.0%
Ammonium Fluoride (as NH ₄ F)	max. 3.0%
Heavy Metals (as Pb)	max. 0.002%
Hydrofluoric Acid	max. 1.0% (w/w)
Insoluble Matter	max. 0.02%
Moisture (as H ₂ O)	max. 0.4%
pH of 5% Solution at 25°C	4.0-6.5

0014.9045	100 lbs	
Product may	contain black specks that is indigenous to	

CONT.

BOX

product.

PACKING

PRODUCT

0014.5000

Ammonium Hydrogen Phosphate

0031 'BAKER ANALYZED' / ACS

► (NH₄)₂HPO₄

M = 132.06 g/mol

CAS NO. 7783-28-0

EINECS 231-987-8 NC CODE 3105 30 00

Exceeds ACS Specifications. Meets Reagents				
Specifications for testing USP/NF monographs				
Assay	min. 98.0%			
Calcium (Ca) max. 0.001				

7 100dy	111111. 00.070
Calcium (Ca)	max. 0.001%
Heavy Metals (as Pb)	max. 0.001%
Insoluble Matter	max. 0.005%
Iron (Fe)	max. 0.001%
Magnesium (Mg)	max. 0.0005%
Nitrate (NO ₃)	max. 0.003%
pH of 5% Solution at 25°C	7.7-8.1
Potassium (K)	max. 0.005%
Sodium (Na)	max. 0.005%
Sulfate (SO ₄)	max. 0.01%
Trace Impurities (in ppm):	
Chloride (CI)	max. 5

max. 5

PRODUCT	PACKING	CONT.
NO.		вох
0031.1000	1 kg	
0031.9050	50 kg	

Ammonium Hydroxide 33% / 'BAKER'

6125

► NH₄OH

M = 35.05 g/mol

1 l = 0.88 kg

CAS NO. 1336-21-6 **EINECS** 215-647-6

NC CODE 2814 20 00

EC NO. 7 001 01 2

UN/ID NO. 2672

ADR/RID 8 C5

IMDG 8/III

R: 34-50

S: 26-36/37/39-45-61





Trace Impurities (in ppm):	
Sulfate (SO ₄)	max. 0.003%
Substances Reducing KMnO ₄	max. 0.0008%
Residue after Ignition	max. 0.003%
Chloride (CI)	max. 0.002%
Assay (as NH ₃) (acidimetric)	30-33%

Residue after Ignition	max. 0.003%
Substances Reducing KMnO ₄	max. 0.0008%
Sulfate (SO ₄)	max. 0.003%
Trace Impurities (in ppm):	
Heavy Metals (as Pb)	max. 1

max. 5

max. 1

PRODUCT	PACKING	CONT.
NO.		вох
6125.1000	11	6
6125.2500	2.5	4
6125.9025	25 l	
For safe han	dling of 25 l tin cans, see S	elf-closing tap.

Ammonium Hydroxide
30% / 'BAKER INSTRA-ANALYZED' / for Trace Metal Analysis

Iron (Fe)

6162

► NH₄OH

M = 35.05 g/mol

1 l = 0.90 kg

CAS NO. 1336-21-6

EINECS 215-647-6

NC CODE 2814 20 00

EC NO. 7 001 01 2

UN/ID NO. 2672 ADR/RID 8 C5

IMDG 8/III

R: 34-50

S: 26-36/37/39-45-61





•	
Assay (NH ₃)	28.0-30.0%
Appearance	passes test
Color (APHA)	max. 10
Residue after Ignition	max. 3 ppm
Specific Gravity at 60°/60°F	0.896-0.902
Substances Reducing KMnO ₄	passes test
Trace Impurities (in ppb):	
Aluminium (AI)	max. 20
Arsenic and Antimony (as As)	max. 50
Barium (Ba)	max. 10
Boron (B)	max. 20
Cadmium (Cd)	max. 5
Calcium (Ca)	max. 200
Chromium (Cr)	max. 5
Cobalt (Co)	max. 1
Copper (Cu)	max. 20
Heavy Metals (as Pb)	max. 200
Iron (Fe)	max. 10
Lead (Pb)	max. 5
Lithium (Li)	max. 20
Magnesium (Mg)	max. 20
Manganese (Mn)	max. 5
Mercury (Hg)	max. 5
Nickel (Ni)	max. 5
Potassium (K)	max. 500
Silicon (Si)	max. 1000
Silver (Ag)	max. 5
Sodium (Na)	max. 300
Strontium (Sr)	max. 10
Tin (Sn)	max. 10
Zinc (Zn)	max. 5
Trace Impurities (in ppm):	
Chloride (CI)	max. 0.5
Phosphate (PO ₄)	max. 0.4

PRODUCT	PACKING	CONT.
NO.		вох
6162.0500	500 ml	6
6162.4000	4 I Glass	4

(Assay about 14.5N).

Assay value tends to be less than reported due to vapor loss, especially on opening container.

Ammonium Hydroxide 29% CMOS, Finyte, Finyte-1, Ultryte Grade

Total Sulfur (as SO₄)

See for detailed information section Microelectronic materials page 32 or the website www.jtbaker.com/micro





Liquid Chromatography (HPLC)

Fields of application for I.T.Baker HPLC solvents

- Pharmaceutical / vitamins
- Food / drinks / beverages / flavours / wine
- Water / environment
- Petrochemistry
- Biology / biochemistry
- Cosmetics
- Paints / lacquers / solvents
- Other chemical industries

We offer a complete product line:

- LC/MS solvents and reagents
- HPLC solvents
- HPLC acids
- HPLC buffer salts
- HPLC ion pair reagents

Quality definition of HPLC Acetonitrile, Methanol and Water

In analytical and preparative HPLC, acetonitrile, methanol and water are the most common and most critical mobile phases. Mallinckrodt Baker has an extensive HPLC solvent and reagent program. As an example, acetonitrile and water qualities versus the most common applications are represented on the next page.

Acetonitrile HPLC qualities versus applications	Preparative HPLC	Isocratic Grade HPLC	Far UV HPLC	Gradient elution 254 nm	Gradient elution 210 nm	Protein profiling	Pesticide Analysis	PAH analysis	Proteomics and LC/MS application
9821 LC/MS Grade									
9017 Ultra Gradient Grade									
9012 Far UV/Gradient Grade									
8275 Isocratic Grade									

Methanol HPLC qualities versus applications	Preparative HPLC	Isocratic Grade HPLC	Gradient elution 254 nm	Proteomics and LC/MS application
9822 LC/MS Grade				
8402 Gradient Grade				
8404 Isocratic Grade				

Water HPLC qualities versus applications	Preparative HPLC	Isocratic Grade HPLC	Gradient elution 254 nm	Gradient elution 210 nm	Proteomics and LC/MS application
9823 LC/MS Grade					
4218 Gradient Grade					

B C D E F G H

Ammon

Ammonium Hydroxide 28-30% / 'BAKER ANALYZED' / ACS

6005

► NH₄OH

M = 35.05 g/mol

1 l = 0.90 kg **CAS NO.** 1336-21-6

EINECS 215-647-6

NC CODE 2814 20 00

EC NO. 7 001 01 2

UN/ID NO. 2672 ADR/RID 8 C5

IMDG 8/III

R: 34-50

S: 26-36/37/39-45-61





Exceeds ACS Specifications	
Assay (as NH ₃) (acidimetric)	28.0-30.0%
Appearance	passes test
Carbon Dioxide (CO ₂)	max. 0.002%
Chloride (CI)	max. 0.5 ppm
Color (APHA)	max. 5
Residue after Ignition	max. 0.001%
Specific Gravity at 60°/60°F	0.896-0.902
Substances Reducing KMnO ₄	passes test
Trace Impurities (in ppm):	
Aluminium (AI)	max. 0.2
Arsenic and Antimony (as As)	max. 3
Barium (Ba)	max. 0.3
Boron (B)	max. 0.05
Chromium (Cr)	max. 0.1
Copper (Cu)	max. 0.1
Heavy Metals (as Pb)	max. 0.5
Iron (Fe)	max. 0.1
Lead (Pb)	max. 0.2
Magnesium (Mg)	max. 0.2
Manganese (Mn)	max. 0.1
Nickel (Ni)	max. 0.1
Nitrate (NO ₃)	max. 2
Phosphate (PO ₄)	max. 0.4
Sulfate (SO ₄)	max. 2
Tin (Sn)	max. 0.1
Titanium (Ti)	max. 0.1
Zinc (Zn)	max. 0.1

PRODUCT	PACKING	CONT.
NO.		вох
6005.1000	11	6
6005.2500	2.5	4

6051

Ammonium Hydroxide

25% / 'BAKER ANALYZED'

► NH₄OH

M = 35.05 g/mol **1 l** = 0.91 kg

CAS NO. 1336-21-6

EINECS 215-647-6

NC CODE 2814 20 00

EC NO. 7 001 01 2

UN/ID NO. 2672 ADR/RID 8 C5

IMDG 8/III

R: 34-50

S: 26-36/37/39-45-61





Assay (NH ₃)	25.0-26.0%
Carbon Dioxide (CO ₂)	max. 0.003%
Chloride (CI)	max. 0.5 ppm
Phosphate (PO ₄)	max. 1 ppm
Residue after Ignition	max. 0.003%
Substances Reducing KMnO ₄	passes test
Total Sulfur (as SO ₄)	max. 5 ppm
Trace Impurities (in ppm):	
Aluminium (AI)	max. 0.05
Arsenic (As)	max. 0.05
Barium (Ba)	max. 0.02
Beryllium (Be)	max. 0.01
Cadmium (Cd)	max. 0.02
Calcium (Ca)	max. 0.5
Chromium (Cr)	max. 0.1
Cobalt (Co)	max. 0.01
Copper (Cu)	max. 0.1
Iron (Fe)	max. 0.1
Lead (Pb)	max. 0.05
Lithium (Li)	max. 0.01
Magnesium (Mg)	max. 0.1
Manganese (Mn)	max. 0.01
Molybdenum (Mo)	max. 0.02
Nickel (Ni)	max. 0.05
Potassium (K)	max. 0.1
Silver (Ag)	max. 0.01
Sodium (Na)	max. 0.5
Strontium (Sr)	max. 0.01
Thallium (TI)	max. 0.05
Titanium (Ti)	max. 0.1
Vanadium (V)	max. 0.01

Zinc (Zn)

Zirconium (Zr)

6051.2500	2.5	
6051.9025	25 l	

PRODUCT

6051.1000

max. 0.05

max. 0.1

PACKING

-11

CONT.

BOX



Ammonium Hydroxide 25% MOS, VLSI Grade See for detailed information section Microelectronic materials page 32 or the website www.jtbaker.com/micro



Ammonium Hydroxide

20% / ULTREX II Ultrapure Reagent

4807

► NH₄OH	
M =	35.05 g/mol
CAS NO.	1336-21-6
EINECS	215-647-6
NC CODE	2814 20 00
EC NO.	7 001 01 2
UN/ID NO.	2672
ADR/RID	8 C5
IMDG	8/III
R:	34
ç.	20-26-36/37/

20-26-36/37/39-4
C C
corrosive

Certificate Provided Reporting Act	ual Lot Analysis
Assay (NH ₃)	20-22%
Trace Impurities (in ppt) (pg/g):	
Aluminium (AI)	max. 20
Antimony (Sb)	max. 10
Arsenic (As)	max. 10
Barium (Ba)	max. 10
Beryllium (Be)	max. 10
Bismuth (Bi)	max. 10
Cadmium (Cd)	max. 10
Calcium (Ca)	max. 20
Cerium (Ce)	max. 10
Cesium (Cs)	max. 10
Chromium (Cr)	max. 10
Cobalt (Co)	max. 10
Copper (Cu)	max. 20
Dysprosium (Dy)	max. 10
Erbium (Er)	max. 10
Europium (Eu)	max. 10
Gadolinium (Gd)	max. 10
Gallium (Ga)	max. 10
Germanium (Ge)	max. 10
Gold (Au)	max. 10
Hafnium (Hf)	act. value reported
Holmium (Ho)	max. 10
Indium (In)	max. 10
Iron (Fe)	max. 20
Lanthanum (La)	max. 10
Lead (Pb)	max. 10
Lithium (Li)	max. 10
Lutetium (Lu)	max. 10
Magnesium (Mg)	max. 20
Manganese (Mn)	max. 10
Molybdenum (Mo)	max. 10
Neodymium (Nd)	max. 10
Nickel (Ni)	max. 10
Niobium (Nb)	max. 10
Palladium (Pd)	act. value reported
Platinum (Pt)	act. value reported
Potassium (K)	max. 20
Praseodymium (Pr)	max. 10

Rhenium (Re)	act. value reported
Rhodium (Rh)	max. 10
Rubidium (Rb)	max. 10
Ruthenium (Ru)	act. value reported
Samarium (Sm)	max. 10
Scandium (Sc)	max. 10
Selenium (Se)	act. value reported
Silver (Ag)	max. 10
Sodium (Na)	max. 20
Strontium (Sr)	max. 10
Tellurium (Te)	max. 10
Terbium (Tb)	max. 10
Thallium (TI)	max. 10
Thorium (Th)	max. 10
Thulium (Tm)	max. 10
Tin (Sn)	max. 10
Titanium (Ti)	max. 10
Tungsten (W)	max. 10
Uranium (U)	max. 10
Vanadium (V)	max. 10
Ytterbium (Yb)	max. 10
Yttrium (Y)	act. value reported
Zinc (Zn)	max. 10
Zirconium (Zr)	max. 10

PRODUCT	PACKING	CONT.
NO.		вох
4807.0490	490 ml PE	

Ammonium Hydroxide Solution

See for qualities meeting the requirements of cGMP, USP/NF, BP/Ph.Eur, JP or FCC, page 36



1110

Ammonium Iron(III) Citrate

CAS NO. 1185-57-5 **EINECS** 214-686-6 NC CODE 2918 15 00

14.5-16.0% Assay (as Fe)

PRODUCT PACKING CONT. NO. BOX 1110.1000 1 kg 1110.9050 50 kg

Ammonium Iron (III) Citrate, brown
See for qualities meeting the requirements of cGMP, USP/NF, BP/Ph.Eur, JP or FCC, page 36



Ammonium Iron (III) Citrate, green

See for qualities meeting the requirements of cGMP, USP/NF, BP/Ph.Eur, JP or FCC, page 36





Ammonium Iron (II) Sulfate Hexahydrate

124 'BAKER ANALYZED' / ACS

(NH₄)₂Fe(SO₄)₂.6H₂O **M** = 392.14 g/mol **CAS NO.** 7783-85-9 **EINECS** 233-151-8 **NC CODE** 2842 90 90

Exceeds ACS Specifications. Meets Reagents Specifications for testing USP/NF monographs Assay (by KMnO₄ titrn.) Calcium (Ca) max. 0.005% Chloride (CI) max. 0.001% Copper (Cu) max. 0.003% Ferric Iron (Fe3+) max. 0.01% Insoluble Matter max. 0.01% max. 0.002% Magnesium (Mg) max. 0.01% Manganese (Mn) pH of 5% Solution at 25°C 3.0-5.0 Phosphate (PO₄) max. 0.003% Potassium (K) max. 0.002% Sodium (Na) max. 0.02% Zinc (Zn) max. 0.003%

PRODUCT	PACKING	CONT.
NO.		вох
0124.0250	250 g	
0124.1000	1 kg	6

Ai

Ammonium Iron(III) Sulfate Dodecahydrate

0118 'BAKER ANALYZED' / ACS

► NH₄Fe(SO₄)₂.12H₂O **M** = 482.19 g/mol **CAS NO.** 7783-83-7

EINECS 233-382-4 NC CODE 2833 30 00

Meets ACS Specifications	
Assay (by lodometry)	98.5-102.0%
Appearance	pale violet crystals
Chloride (CI)	max. 0.001%
Copper (Cu)	max. 0.003%
Ferrous Iron (Fe++)	passes test
Insoluble Matter	max. 0.01%
Nitrate (NO ₃)	max. 0.01%
Substances not Precipitated by NH ₄ OH	max. 0.05%
Zinc (Zn)	may 0.003%

	PACKING	CONT
NO.		ВОХ
0118.0100	100 g	
0118.1000	1 kg	

1010

Ammonium Metavanadate

Powder / 'BAKER ANALYZED' / ACS

► NH₄VO₃

M = 116.98 g/mol **CAS NO.** 7803-55-6

EINECS 232-261-3

NC CODE 2841 90 30

UN/ID NO. 2859

ADR/RID 6.1 T5

IMDG 6.1/II

R: 25-36/37/38

S: 37-45



Meets ACS Specifications. Meets Reagent Specifications for testing USP/NF monographs

Assay min. 99.0% Appearance passes test Carbonate (CO_3) passes test Chloride (CI) max. 0.2% Solubility in Ammonium Hydroxide passes test Sulfate (SO_4) max. 0.05%

 PRODUCT
 PACKING
 CONT.

 NO.
 BOX

 1010.0125
 125 g



Ammonium Molybdate

See Ammonium Heptamolybdate Tetrahydrate



Ammonium Monovanadate

See Ammonium Metavanadate

Certificates of Analysis are available at www.jtbaker.com/europe



Ammonium Nitrate

'BAKER ANALYZED' / ACS

0025

► NH₄NO₃	
M =	80.04 g/mol
CAS NO.	6484-52-2
EINECS	229-347-8
NC CODE	3102 30 90
UN/ID NO.	1942
ADR/RID	5.1 O2
IMDG	5.1/III
R:	8-9
S:	15-16-41
	_

Exceeds ACS Specifications. Meets Reagents Specifications for testing USP/NF monographs		
Insoluble Matter	max. 0.005%	
Nitrite (NO ₂)	passes test	
pH of 5% Solution at 25°C	4.5-6.0	
Residue after Ignition	max. 0.01%	
Sulfate (SO ₄)	max. 0.002%	
Trace Impurities (in ppm):		
Chloride (CI)	max. 5	
Heavy Metals (as Pb)	max. 5	
Iron (Fe)	max. 2	
Phosphate (PO ₄)	max. 5	

	PRODUCT	PACKING	CONT. BOX
6	0025.0100	100 g	
6	0025.1000	1 kg	
t			
0			
6			
6			
5			
5 5 2			
2			

Ammonium Peroxodisulfate

'BAKER ANALYZED' / ACS

0028

► (NH ₄) ₂ S ₂ O ₈	
M =	228.20 g/mol
CAS NO.	7727-54-0
EINECS	231-786-5
NC CODE	2833 40 00
EC NO.	16 060 00 6
UN/ID NO.	1444
ADR/RID	5.1 O2
IMDG	5.1/III
R:	22-36/37/38-42/43-8

weets ACS Specifications	
Assay	min. 98.0%
Chloride and Chlorate (as CI)	max. 0.001%
Heavy Metals (as Pb)	max. 0.005%
Insoluble Matter	max. 0.005%
Iron (Fe)	max. 0.001%
Residue after Ignition	max. 0.05%
Titrable Free Acid (meq/g)	max. 0.04
Trace Impurities (in ppm):	
Manganese (Mn)	max. 0.5

 PRODUCT
 PACKING
 CONT.

 NO.
 BOX

 0028.1000
 1 kg
 6


 0028.5000
 5 kg





Ammonium Peroxodisulfate

1899



EINECS 231-786-5
NC CODE 2833 40 00
EC NO. 16 060 00 6
UN/ID NO. 1444

ADR/RID 5.1 O2 IMDG 5.1/III

R: 22-36/37/38-42/43-8

S: 22-24-26-37





Assay min. 98.0% Chloride and Chlorate (as CI) max. 0.002% Heavy Metals (as Pb) max. 0.003% Iron (Fe) max. 0.003% Trace Impurities (in ppm): Manganese (Mn)

 PRODUCT NO.
 PACKING PACKING
 CONT. BOX

 1899.9050
 50 kg

Ammonium Persulfate

See Ammonium Peroxodisulfate



Ammonium Phosphate, Dibasic

See Ammonium Hydrogen Phosphat



Ammonium Phosphate, Monobasic

See Ammonium Dihydrogen Phosphate



Ammonium Purpurate

See Murexide



Ammonium Reineckate

See Reinecke Sal



Ammonium Rhodanide

e Ammonium Thiocvanate

1011

Ammonium Sulfamate

'BAKER ANALYZED' / ACS

► NH₂SO₃NH₄

M = 114.12 g/mol

CAS NO. 7773-06-0 **EINECS** 231-871-7

NC CODE 2842 90 90

Meets ACS Specifications. Meets Reagent Specifications for testing USP/NF monographs

Assay min. 98.0% Insoluble Matter max. 0.02% Melting Range max. 2.0°C Recorded Melting Point 133.0°C Residue after Ignition max. 0.10% Trace Impurities (in ppm):

max. 5

max. 0.5

PRODUCT PACKING CONT. BOX 500 g 1011.0500

4027

Ammonium Sulfate

► (NH₄)₂SO₄ **M** = 132.14 g/mol

CAS NO. 7783-20-2

EINECS 231-984-1

NC CODE 3102 21 00

min. 99.5% Assay DNase Activity none detected Insoluble Matter max. 0.005% pH of 5% Solution at 25°C max. 6.0 Protease Activity none detected Residue after Ignition max. 0.005% RNase Activity none detected Solution Test passes test

Absorbance of a 1 M Solution, Maximum (1-cm path vs water):

at 260 nm

Heavy Metals (as Pb)

at 280 nm 0.04

Trace Impurities (in ppm): Arsenic (As) max. 0.2 Cadmium (Cd) max. 2 Calcium (Ca) max. 2 Chloride (CI) max. 3 Copper (Cu) max. 0.2 Iron (Fe) max. 1 Lead (Pb) max. 0.1 Lithium (Li) max. 2 Magnesium (Mg) max. 2 Manganese (Mn) max. 0.6 Nitrate (NO₃) max. 10 Phosphate (PO₄) max. 3 Potassium (K) max. 2

PRODUCT	PACKING	CONT.
NO.		вох
4027.1000	1 kg	
4027.5000	5 kg	
4027.9012	12 kg	

For Protein Precipitation and Liquid chromatography.

Ammonium Sulfate BAKER ANALYZED' / ACS

► (NH₄)₂SO₄

0032

M = 132.14 g/mol

CAS NO. 7783-20-2 **EINECS** 231-984-1 NC CODE 3102 21 00

Meets ACS Specifications	
Assay (by formol method)	min. 99.0%
Insoluble Matter	max. 0.005%
Nitrate (NO ₃)	max. 0.001%
pH of 5% Solution at 25°C	5.0-6.0
Residue after Ignition	max. 0.005%
Trace Impurities (in ppm):	

Insoluble Matter	max. 0.005%
Nitrate (NO ₃)	max. 0.001%
pH of 5% Solution at 25°C	5.0-6.0
Residue after Ignition	max. 0.005%
Trace Impurities (in ppm):	
Arsenic (As)	max. 0.3
Chloride (CI)	max. 5
Heavy Metals (as Pb)	max. 5
Iron (Fe)	F
non (re)	max. 5

PRODUCT	PACKING	CONT.
NO.		вох
0032.1000	1 kg	6
0032.5000	5 kg	4
0032.9025BB	25 kg Bag-in-Box	
0032.9050	50 kg	



Ammonium Sulfate

See for qualities meeting the requirements of cGMP, USP/NF, BP/Ph.Eur, JP or FCC, page 36

Zinc (Zn)

Ammonium Sulfide

% / 'BAKER ANALYZED

► (NH₄)₂S

M = 68.14 g/mol

1 l = 1.00 kg

 $\textbf{FLASHPOINT} \quad 15~^{\circ}\text{C}$

CAS NO. 12135-76-1

EINECS 235-223-4

NC CODE 2830 90 80

UN/ID NO. 2683

ADR/RID 8 CFT

IMDG 8/II

R: 11-31-34

S: 16-23-26-36/37/39-45





Assay	min. 20.0%
Calcium (Ca)	max. 0.003%
Carbonate (CO ₃)	max. 0.005%
Chloride (CI)	max. 0.05%
Potassium (K)	max. 0.001%
Residue after Ignition	max. 0.05%
Sodium (Na)	max. 0.005%
Trace Impurities (in ppm):	
Copper (Cu)	max. 5
Iron (Fe)	max. 5
Magnesium (Mg)	max. 5

PRODUCT	PACKING	CONT
NO.		BOX
7003.1000	11	

Ammonium Thiocyanate

'BAKER ANALYZED' / ACS

0034

2662

7122

4650

CONT.

BOX

7003

NH₄SCN

M = 76.12 g/mol

CAS NO. 1762-95-4

EINECS 217-175-6

NC CODE 2838 00 00 **EC NO.** 615 004 00 3

R: 20/21/22-32

S: 13-2



Exceeds ACS Specifications, wieets neagents		
Specifications for testing USP/NF monographs		
Assay (argentometric titrn.)	min. 97.5%	
Appearance	passes test	

max. 0.005% Chloride (CI) Insoluble Matter max. 0.005% lodine-consuming Substances (as $Na_2S_2O_3$) (meq/g) max. 0.004

pH of 5% Solution at 25°C 4.5-6.0 Residue after Ignition max. 0.01% max. 0.005% Sulfate (SO₄)

Trace Impurities (in ppm):

Heavy Metals (as Pb) max. 5 Iron (Fe) max. 3

PRODUCT	PACKING	CONT.
NO.		вох
0034.0100	100 g	
0034.1000	1 kg	

Ammonium Thiocyanate

M = 76.12 g/mol

CAS NO. 1762-95-4 **EINECS** 217-175-6

► NH₄SCN

NC CODE 2838 00 00

EC NO. 615 004 00 3

R: 20/21/22-32

S: 13-2



Assay	min. 99.0%
Iron (Fe)	max. 3 ppm
pH of 5% Solution at 25°C	4.0-6.0

PRODUCT PACKING CONT. BOX 2662 9025 25 kg

Ammonium Thiocyanate

M = 76.12 g/mol **CAS NO.** 1762-95-4 **EINECS** 217-175-6 NC CODE 2838 00 00

Titer (mol/l) 0.0995-0.1005 **PRODUCT PACKING** NO. 7122.1000 11

Volumetric Solution, ready for use.

Ammonium Thiocyanate

0.1 mol/l / DILUT-IT

PRODUCT PACKING CONT. 4650 1 amp.

CAS NO. 1762-95-4 **EINECS** 217-175-6

M = 76.12 g/mol

NH₄SCN

NC CODE 2838 00 00

Volumetric Concentrate, for dilution to 1 l.

Ammoniummeta-Vanadate

e Ammonium Metavanadat

8009

n-Amyl Alcohol BAKER ANALYZED

► CH₃(CH₂)₄OH

M = 88.15 g/mol **11=** 0.82 kg

FLASHPOINT 33 °C

CAS NO. 71-41-0

EINECS 200-752-1

NC CODE 2905 15 00

EC NO. 603 006 00 7

UN/ID NO. 1105

ADR/RID 3 F1

IMDG 3/III

R: 10-20 **S:** 24/25

Assay (by GC) min. 98% Acids and Esters (as Amyl Acetate) max. 0.01% Aldehydes passes test Boiling Range 135-138°C Color (APHA) max. 5 max. 0.003% Residue after Evaporation Water (H₂O) max. 0.2%

PRODUCT PACKING CONT. BOX 8009.1000 11 8009.9025 25 I

8012

tert-Amyl Alcohol

BAKER ANALYZED

► CH₃CH₂C(CH₃)₂OH

M = 88.15 g/mol

11 = 0.81 kg

FLASHPOINT 19 °C

CAS NO. 75-85-4

EINECS 200-908-9

NC CODE 2905 15 00

EC NO. 603 007 00 2

UN/ID NO. 1105

ADR/RID 3F1 IMDG 3/II

R: 11-20

S: 16-24/25-9





Assay (by GC) min. 99.0% Acids and Esters (as Amyl Acetate) max. 0.06% Aldehydes passes test max. 0.004% Residue after Evaporation

PRODUCT PACKING CONT. вох 8012.0100 100 ml 8012.1000 11

Anhydrone

'BAKER ANALYZED' / magnesium perchlorate for drying / ACS

► Mg(ClO₄)₂

0035

M = 223.21 g/mol

CAS NO. 10034-81-8

EINECS 233-108-3

NC CODE 2829 90 10

UN/ID NO. 1475

ADR/RID 5.1 O2

IMDG 5.1/II

R: 36/37/38-9 **S:** 17-22-27





Meets ACS Specifications Loss on Drying at 190°C max. 8% Suitability for Moisture Absorption passes test max. 0.025 Titrable Base (meq/g) Titrable Free Acid (meq/g) max. 0.005

PRODUCT PACKING CONT. NO. BOX 0035.0500 500 g



Anhydrous Solvents

e for detailed information section Solvent applications, page 22

CONT.

Aniline

'BAKER ANALYZED' / ACS

► C₆H₅NH₂

M = 93.13 g/mol

1 l = 1.02 kg

FLASHPOINT 70 °C

CAS NO. 62-53-3

EINECS 200-539-3

NC CODE 2921 41 00

EC NO. 612 008 00 7

UN/ID NO. 1547

ADR/RID 6.1 T1

IMDG 6.1/II

R: 23/24/25-40-41-43-48/23/24/25-50-68

S: 26-27-36/37/39-45-46-61-63





7004

NO. вох

7004.0250 250 ml 7004.1000 11

PACKING

PRODUCT

min. 99.0%

max. 0.01%

max. 250

passes test

passes test

184.4°C

passes test

183-186°C

max. 0.005%

NOTE: This product darkens to a reddish-brown color on

Aniline

► C₆H₅NH₂ **M** = 93.13 g/mol

11 = 1.02 kg

FLASHPOINT 70 °C

CAS NO. 62-53-3 **EINECS** 200-539-3

NC CODE 29214100

EC NO. 612 008 00 7

UN/ID NO. 1547

ADR/RID 6.1 T1

IMDG 6.1/II

R: 23/24/25-40-41-43-48/23/24/25-

50-68

S: 26-27-36/37/39-45-46-61-63





dangerous for the environment

7005

8013

CONT.

NO. BOX 7005.1000 11

Stored protected from light.

PACKING

PRODUCT

Anisole

Meets ACS Specifications. Meets Reagent

Product Information (not specifications):

Assay (by GC) (corrected for water)

Chlorobenzene (C_6H_5CI)

Nitrobenzene (C₆H₅NO₂)

Residue after Ignition

Boiling Point (typical)

Color (APHA)

Hydrocarbons

Appearance

Boiling Point

Specifications for testing USP/NF monographs

► C₆H₅OCH₃

M = 108.14 g/mol

1 l = 0.99 kg

FLASHPOINT 43 °C

CAS NO. 100-66-3 **EINECS** 202-876-1

NC CODE 2909 30 90

UN/ID NO. 2222

ADR/RID 3 F1

IMDG 3.3/III

R: 10-36/37/38

S: 16-26-37/39



Assay (by GC) min. 99.0% **Boiling Point** 152-154°C Phenol (C₆H₅OH) max. 0.2% max. 0.5% Water (H₂O)

PRODUCT PACKING CONT. вох 8013.0500 500 ml

Anthr

Anthrone 1231

BAKER ANALYZED' / ACS

▶ Ç₆H₄COC₆H₄ÇH₂

M = 194.24 g/mol **CAS NO.** 90-44-8

EINECS 201-994-0 NC CODE 2914 39 00 Meets ACS Specifications Absorbance of reagent solution passes test Appearance passes test Melting range, within a 5°C range, including 156°C passes test Sensitivity to Carbohydrates passes test

passes test

PRODUCT CONT. **PACKING** вох 1231.0025 25 g Glass

Antifoam B

11 = 1.00 kg NC CODE 3822 00 00

PRODUCT	PACKING	CONT.
NO.		вох
7356.0100	100 ml	
7356.0500	500 ml	

Antimony 1000 µg/ml

5703 ice of tartaric acid) / 'BAKER INSTRA-ANALYZED' / Plasma Standard

Solubility in Ethyl Acetate

► Sb

M = 121.74 g/mol

NC CODE 3822 00 00

R: 36/38

S: 26-28-45



Certificate Provided Reporting Actual Lot Analysis Antimony (Sb) $998\text{-}1002~\mu\text{g/ml}$

PRODUCT PACKING CONT. BOX 5703.0100 100 ml

Prepared from the highest purity raw material available, generally greater than 99.999% spectral purity. The content of the solution is confirmed to be accurate to within \pm 0.2%. Typically 1000 $\mu\text{g/ml.}$ The certificate of analysis provided reports actual lot analysis. The certificate also lists the trace impurities.

Antimony 1000 µg/ml

(Matrix: 20% hydrochloric acid) / 'BAKER INSTRA-ANALYZED' / Atomic Absorption Standard

Antimony (Sb)

M = 121.74 g/mol **CAS NO.** 28300-74-5

NC CODE 3822 00 00

R: 36/38 S: 26



Antimony (Sb) $998\text{-}1002~\mu\text{g/ml}$

NO. 6968.0100 100 ml 6968.0500 500 ml

PRODUCT

 $998\text{-}1002~\mu\text{g/ml}$

PACKING

Prepared by dissolution of high purity raw materials (min. 99.99% spectral purity). Assays are verified by ICP against standards traceable to NIST. Standard Reference Material numbers (SRM) are printed on each label.

CONT

вох

6802

6968

Antimony 1000 µg/ml

'BAKER ANALYZED' / Atomic n Standard

► Sb **M =** 121.74 g/mol NC CODE 3822 00 00

UN/ID NO. 1789

ADR/RID 8 C1

IMDG 8/II

S: 26-28-45

R: 36/37/38

PRODUCT PACKING CONT. BOX 6802.0100 100 ml 6802.0500 500 ml

Antimony trichloride in hydrochloric acid 5 mol/l.



Antimony 10000 μg/ml

e of tartaric acid) / 'BAKER INSTRA-ANALYZED' / Plasma Standard

5717

► Sb

M = 121.74 g/mol

NC CODE 3822 00 00

R: 36/38 **S:** 26-28-45

Certificate Provided Reporting Actual Lot Analysis Antimony (Sb)

PRODUCT 9980-10020 μg/ml NO.

PACKING CONT. вох 5717.0100 100 ml

Prepared from the highest purity raw material available, generally greater than 99.999% spectral purity. The content of the solution is confirmed to be accurate to within 0.2 %. Typically 10000 μ g/ml.

The certificate of analysis provided reports actual lot analysis. The certificate also lists the trace impurities.

Antimony Potassium Tartrate, Trihydrate See for qualities meeting the requirements of GGMP, USP/NF, BP/Ph.Eur, JP or FCC, page 36



Antimony Trichloride 'BAKER ANALYZED' / ACS

0037

► SbCl₃

M = 228.11 g/mol

CAS NO. 10025-91-9

EINECS 233-047-2

NC CODE 2827 39 80

EC NO. 51001008

UN/ID NO. 1733

ADR/RID 8 C2 IMDG 8/II

R: 34-51/53

S: 26-45-61





Meets ACS Specifications	
Assay	min. 99.0%
Arsenic (As)	max. 0.02%
Calcium (Ca)	max. 0.005%
Copper (Cu)	max. 0.001%
Insoluble in Chloroform	max. 0.05%
Iron (Fe)	max. 0.002%
Lead (Pb)	max. 0.005%
Potassium (K)	max. 0.01%
Sodium (Na)	max. 0.02 %
Sulfate (SO ₄)	max. 0.005%

PRODUCT	PACKING	CONT.
NO.		вох
0037.0100	100 g	

Antimony(III) Chloride



4149

Apiezon Grease M

'BAKER' / for gaschromatography

NC CODE 2710 1999

PRODUCT	PACKING	CONT.
NO.		BOX
4149.0025	25 g Glass	

L-Arginine

See for qualities meeting the requirements of cGMP, USP/NF, BP/Ph.Eur, JP or FCC, page 36



L-Arginine Hydrochloride

See for qualities meeting the requirements of cGMP, USP/NF, BP/Ph.Eur, JP or FCC, page 36



Arsenic 1000 µg/ml

(Matrix: 2% nitric acid) / 'BAKER INSTRA-ANALYZED' / Plasma Standard



As

M = 74.92 g/mol

NC CODE 3822 00 00

EC NO. 33 001 00 0 R: 22-36/38-45

S: 20-26-36/37/39-53



Certificate Provided Reporting Actual Lot Analysis 998-1002 μg/ml Arsenic (As)

NO. 5704.0100 100 ml

PACKING

PRODUCT

CONT. BOX

Prepared from the highest purity raw material available, generally greater than 99.999% spectral purity. The content of the solution is confirmed to be accurate to within \pm 0.2%. Typically 1000 $\mu\text{g/ml}.$

The certificate of analysis provided reports actual lot analysis. The certificate also lists the trace impurities.

Arsen

Arsenic 1000 µg/ml

6919 (Matrix: 2% nitric acid) / 'BAKER INSTRA-ANALYZED' / Atomic Absorption Standard

Arsenic (As)

Arsenic (As)

► As **M** = 74.92 g/mol NC CODE 3822 00 00

> **EC NO.** 33 001 00 0 R: 22-36/38-45

> > **S:** 20-26-28-53



PRODUCT PACKING CONT. NO. вох 6919.0100 100 ml 6919.0500 500 ml

Prepared by dissolution of high purity raw materials (min. 99.99% spectral purity). Assays are verified by ICP against standards traceable to NIST. Standard Reference Material numbers (SRM) are printed on each label.

Arsenic 1000 µg/ml

6803 'BAKER ANALYZED' / Ā

► As **M** = 74.92 g/mol NC CODE 3822 00 00

R: 22-36/38-45

S: 20-26-36/37/39-45-53



CONT **PRODUCT PACKING** BOX NO. 6803.0100 100 ml 6803.0500 500 ml

Arsenic trioxide in nitric acid 0.5 mol/l.

Arsenic 10000 μ g/ml (Matrix: 2% nitric acid) / "BAKER INSTRA-ANALYZED" / Plasma Standard

► As **M** = 74.92 g/mol NC CODE 3822 00 00

> EC NO. 33 001 00 0 R: 25-36/38-45-52/53

> > **S:** 20-26-36/37/39-45-53



Certificate Provided Reporting Actual Lot Analysis Arsenic (As) 9980-10020 µg/ml

998-1002 μg/ml

 $998\text{-}1002~\mu\text{g/ml}$

PRODUCT **PACKING** CONT. NO. BOX 5718.0100 100 ml

Prepared from the highest purity raw material available, generally greater than 99.999% spectral purity. The content of the solution is confirmed to be accurate to within 0.2 %. Typically 10000 μ g/ml. The certificate of analysis provided reports actual lot analysis. The certificate also lists the trace impurities.

Ascorbic Acid

ee for qualities meeting the requirements of cGMP, USP/NF, BP/Ph.Eur, JP or FCC, page 36

1018

L(+)-Ascorbic Acid

BAKER ANALYZED' Bioche

► HOCH₂CH(OH)ÇHC(OH):C(OH)COQ

M = 176.13 g/mol **CAS NO.** 50-81-7 **EINECS** 200-066-2

NC CODE 2936 27 00

Assay (dried) min. 99.5% Heavy Metals (as Pb) max. 0.001% Iron (Fe) max. 5 Loss on Drying at 105°C max. 0.1% max. 0.05% Residue on Ignition Specific Rotation $\left[\alpha\right]^{\frac{20}{D}}$ (dried basis), c=10in water) $+20.5^{\circ}$ to $+21.5^{\circ}$ passes test Water Insoluble Material

PRODUCT	PACKING	CONT.
NO.		вох
1018.0100	100 g	6
1018.0500	500 g	6

1914

L(+)-Ascorbic Acid

► HOCH₂CH(OH)ÇHC(OH):C(OH)COQ

M = 176.13 g/mol CAS NO. 50-81-7 **EINECS** 200-066-2 NC CODE 2936 27 00

Appearance of solution passes test Assay 99.0 - 100.5% Copper (Cu) max. 5 ppm Heavy Metals (as Pb) max. 10 ppm Identification passes test Iron (Fe) max. 2 ppm Oxalic acid passes test +20.5 - +21.5° Specific Optical Rotation Sulfated Ash max. 0.1%

PRODUCT	PACKING	CONT.
NO.		вох
1914.0100	100 g	6
1914.0500	500 g	6

Stored in an airtight non-metallic container. Stored protected from light.



L-Asparagine, Monohydrate

nts of cGMP, USP/NF, BP/Ph.Eur, JP or FCC, page 36

DL-Aspartic Acid'BAKER ANALYZED' Bioch 1546 ► HOCOCH(NH₂)CH₂COOH PRODUCT PACKING CONT. min. 99% Assay **M** = 133.10 g/mol Heavy Metals (as Pb) max. 0.002% NO. вох **CAS NO.** 617-45-8 Homogenity (by TLC) 1546.0500 500 g passes test **EINECS** 210-513-3 Loss on Drying at 105°C max. 0.3% NC CODE 2922 49 95 Residue after Ignition max. 0.3% Trace Impurities (in ppm): Arsenic (As) max. 3 **Auric Chloride BAKER DUAL-TINT** 2867-01 pH 1-12 / pH Indicator Pape Product Information (not specifications): PRODUCT PACKING CONT. pH Gradation вох NO. 2867-01 Dispenser contains 5m For Refill see PN 2873-01. **BAKER DUAL-TINT** 2873-01 pH 1-12 / pH Indicator Pape Product Information (not specifications): **PRODUCT PACKING** CONT. BOX pH Gradation 2873-01 Refill contains 3x5m Refill for PN 2867-01. **BAKER DUAL-TINT** pH 1.0-4.3 / pH Indicator Pag 2878-01 Product Information (not specifications): **PRODUCT PACKING** CONT. pH Gradation 0.3 NO. BOX 2878-01 Dispenser contains 5m For Refill see PN 2879-01. **BAKER DUAL-TINT** 2879-01 pH 1.0-4.3 / pH Indicator Paper PRODUCT PACKING CONT. Product Information (not specifications): 0.3 pH Gradation BOX 2879-01 Refill contains 3x5m Refill for PN 2878-01. **BAKER DUAL-TINT** 2871-01 pH 3.5-6.8 / pH Indicator Paper PRODUCT PACKING CONT. Product Information (not specifications): pH Gradation 0.3 NO. вох 2871-01 Dispenser contains 5m For Refill see PN 2876-01. **BAKER DUAL-TINT** pH 3.5-6.8 / pH Indicator Paper 2876-01

Product Information (not specifications):

pH Gradation

PRODUCT

Refill for PN 2871-01.

NO. 2876-01

0.3

PACKING

Refill contains 3x5m

CONT. BOX



BAKER DUAL-TINT

2869-01 pH 5.0-8.0 / pH Indicator Pag

> PRODUCT CONT. Product Information (not specifications): **PACKING** pH Gradation 0.3 NO. вох 2869-01 Dispenser contains 5m

> > For Refill see PN 2875-01.

BAKER DUAL-TINT

2875-01 pH 5.0-8.0 / pH Indicator Pape

> Product Information (not specifications): **PRODUCT PACKING** CONT. pH Gradation 0.3 BOX NO. Refill contains 3x5m 2875-01

> > Refill for PN 2869-01.

BAKER DUAL-TINT

2868-01 pH 7.0-10.0 / pH Indicator Paper

> Product Information (not specifications): **PRODUCT PACKING** CONT. pH Gradation 0.3 BOX 2868-01 Dispenser contains 5m

> > For Refill see PN 2874-01.

BAKER DUAL-TINT

2874-01 pH 7.0-10.0 / pH Indicator Paper

> Product Information (not specifications): **PRODUCT PACKING** CONT. pH Gradation 0.3 BOX NO. 2874-01 Refill contains 3x5m

> > Refill for PN 2868-01.

BAKER DUAL-TINT 2880-01 pH 9.5-14.0 / pH Indicator Paper

> Product Information (not specifications): pH Gradation 0.5

PRODUCT PACKING CONT. вох NO. 2880-01 Dispenser contains 5m

For Refill see PN 2881-01.

BAKER DUAL-TINT 2881-01

4408

pH 9.5-14.0 / pH Indicator Pa

Product Information (not specifications): 0.5 pH Gradation

PRODUCT PACKING CONT. BOX 2881-01 Refill contains 3x5m

Refill for PN 2880-01.

BAKER TESTRIPS for Ammonium (NH₄*)

Range: 0-400 mg/l (ppm) / Test Strips for Semi-Qu

PRODUCT PACKING CONT. NO. BOX 4408 1 box

Each plastic box contains 100 strips, 6x95 mm. Gradations: 0-10-25-50-100-200-400 mg/l (ppm). Kits include the necessary reagents for preparation and test strips for testing. Test strips with color scale for quick and accurate

determinations.



BAKER TESTRIPS for Ascorbic Acid

Range: 0-2000 mg/l (ppm) / Test Strips for Semi-Quantitative Determinations of lons

4409

PRODUCT	PACKING	CONT.
NO.		вох
4409	1 box	

Each plastic box contains 100 strips, 6x95 mm.

Gradation: 0-50-100-200-300-500-1000-2000 mg/l (ppm).

Test strips with color scale for quick and accurate determinations.

BAKER TESTRIPS for Chromate (CrO₄²⁻)

Range: 0-100 mg/l (ppm) / Test Strips for Semi-Quantitative Determinations of lons

4410

PRODUCT NO.	PACKING	CONT.	
		вох	
4410	1 box		

Each plastic box contains 100 strips, 6x95 mm.
Gradations: 0-3-10-30-100 mg/L (ppm).
Kits include the necessary reagents for preparation and test strips for testing.

Test strips with color scale for quick and accurate determinations.

BAKER TESTRIPS for Cobalt (Co2+)

Range: 0-1000 mg/l (ppm) / Test Strips for Semi-Quantitative Determinations of lons

4411

PRODUCT	PACKING	CONT.
NO.		вох
4411	1 box	

Each plastic box contains 100 strips, 6x95 mm.

Gradations: 0-10-25-50-100-250-500-1000 mg/l (ppm).

Test strips with color scale for quick and accurate determinations.

BAKER TESTRIPS for Copper (Cu2+)

Range: 0-300 mg/l (ppm) / Test Strips for Semi-Quantitative Determinations of lons

4412

PRODUCT	PACKING	CONT.
NO.		вох
4412	1 box	

Each plastic box contains 100 strips, 6x95 mm. Gradation: 0-10-30-100-300 mg/l (ppm). Test strips with color scale for quick and accurate determinations.

BAKER TESTRIPS for Iron (Fe2+)

Range: 0-100 mg/l (ppm) / Test Strips for Semi-Quantitative Determinations of lons

4413

PACKING	CONT
	B0)
1 box	
c box contains 100 strips, 6x	95 mm.
	1 box

Gradations: 0-2-5-10-25-50-100 mg/l (ppm).

Kits include the necessary reagents for preparation and test strips for testing.

Test strips with color scale for quick and accurate

Test strips with color scale for quick and accurate determinations.



4414

4421

BAKER

BAKER TESTRIPS for Nickel (Ni $^{2+}$) Range: 0-1000 mg/l (ppm) / Test Strips for Semi-Qu

PRODUCT	PACKING	CONT.
NO.		вох
4414	1 hox	

Each plastic box contains 100 strips, 6x95 mm. Gradations: 0-10-25-50-100-250-500-1000 mg/l (ppm). Test strips with color scale for quick and accurate determinations.

BAKER TESTRIPS for Nitrate / Nitrite (N0 $_3$ -/N0 $_2$ -) Range: 0-500 mg/I NO $_3$ - (ppm) ; 0-80 mg/I NO $_2$ - (ppm) / Test Strips for Semi-

PRODUCT	PACKING	CONT.
NO.		вох
4421	1 box	

Each plastic box contains 100 strips, 6x95 mm. Gradations (NO_2^-) : 0-1-5-10-20-40-80 mg/l (ppm). Gradations (NO₃⁻): 0-10-25-50-100-250-500 mg/l (ppm). Test strips with color scale for quick and accurate determinations.

BAKER TESTRIPS for Nitrite (NO₂-)

Range: 0-80 mg/l (ppm) / Test Strips for S

PRODUCT	PACKING	CONT.
NO.		вох
4415	1 box	
5	1 00%	

Each plastic box contains 100 strips, 6x95 mm. Gradation: 0-1-5-10-20-40-80 mg/l (ppm). Test strips with color scale for quick and accurate determinations.

BAKER TESTRIPS for Peroxide (H₂O₂)

Range: 0-100 mg/I (ppm) / Test Strips for Semi-Qu

PRODUCT	PACKING	CONT.
NO.		вох
4416	1 box	

Each plastic box contains 100 strips, 6x95 mm. Gradation: 0-1-3-10-30-100 mg/l (ppm). Test strips with color scale for quick and accurate determinations.

PACKING

PRODUCT

determinations.

BAKER TESTRIPS for Potassium (K+)

Range: 0-1500 mg/l (ppm) / Test Strips for Semi-Qu

NO.		вох
4417	1 box	
Gradation	tic box contains 100 strips, 6x95 mm. ls: 0-200-400-700-1000-1500 mg/l (pp	,
Kits includ	de the necessary reagents for preparat	ion and
test strips	for testing.	

Test strips with color scale for quick and accurate

CONT.

80



BAKER TESTRIPS for Sulfite (\$0₃²⁻) Range: 0-1000 mg/I (ppm) / Test Strips for Semi-Quar

4418

PRODUCT	PACKING	CONT.
NO.		вох
4418	1 box	

Each plastic box contains 100 strips, 6x95 mm. Gradations: 0-10-25-50-100-250-500-1000 mg/l (ppm). Test strips with color scale for quick and accurate determinations.

BAKER TESTRIPS for Tin (Sn²⁺)

Test Strips for Semi-Quantitative Determinations of lons

PRODUCT	PACKING	CONT.	
NO.		вох	
4419	1 hox		

Each plastic box contains 100 strips, 6x95 mm. Gradations: 0-10-25-50-100-250-500 mg/l (ppm). Test strips with color scale for quick and accurate

BAKER TESTRIPS for Water Hardness

4420

PRODUCT	PACKING	CONT.
NO.		вох
4420	1 hox	

 $(1^{\circ}d = 17.8 \text{ mg/I CaCO}_3).$

Each plastic box contains 100 strips, 6x95 mm. Gradation (German Hardness): <3 - >5 - >10 - >15 ->20 - >25 °d.

Test strips with color scale for quick and accurate determinations.

BakerClean

See for detailed information section Microelectronic materials page 32 or the website www.jtbaker.com/micro



Baker-Clear

'BAKER' / Cleaning solution for general purpose

7341

1 =	0.95 kg
FLASHPOINT	79 °C
ç.	35

PACKING CONT. **PRODUCT** NO. 7341.1000 7341.2500 2.5 |

Baker-Clear 12

ning solution (pH 12) for general pupose

7340

NC CODE	3822 00 00
UN/ID NO.	1760
ADR/RID	8 C9
IMDG	8/111
R:	38-41

PACKING	CONT.
	вох
5 l Jerrycan	

S: 26-37/39-60

Use a 2% dilution for general use, 5% for intensive use or 20% for a strong cleaning effect (dilute with water).

4399

4400

4402

BAKER-pHIX pH Papers with color scale

pH 1-12 / pH Test Strips with color scale (also suitable for colored solutions)

Product Information (not specifications):

pH Gradation

 PRODUCT
 PACKING
 CONT.

 NO.
 BOX

 4399
 1 box

Each plastic box contains 200 strips, 11x100 mm. Indicator dyes and color scale fixed onto a single strip, for quick and accurate readings without referring to a seperate scale.

BAKER-pHIX pH Papers with color scale

pH 1.8-3.8 / pH Test Strips with color scale (also suitable for colored solutions

pH Gradation 0.2/0.3

 PRODUCT
 PACKING
 CONT.

 NO.
 BOX

 4400
 1 box

Each plastic box contains 200 strips, 11x100 mm. Indicator dyes and color scale fixed onto a single strip, for quick and accurate readings without referring to a seperate scale.

BAKER-pHIX pH Papers with color scale

pH 2.8-4.6 / pH Test Strips with color scale (also suitable for colored solutions)

Product Information (not specifications):
pH Gradation 0.2/0.3

 PRODUCT
 PACKING
 CONT.

 NO.
 BOX

 4401
 1 box

Each plastic box contains 200 strips, 11x100 mm. Indicator dyes and color scale fixed onto a single strip, for quick and accurate readings without referring to a seperate scale.

BAKER-pHIX pH Papers with color scale

pH 3.8-5.5 / pH Test Strips with color scale (also suitable for colored solutions)

Product Information (not specifications):

pH Gradation 0.2/0.3

 PRODUCT
 PACKING
 CONT.

 NO.
 BOX

 4402
 1 box

Each plastic box contains 200 strips, 11x100 mm. Indicator dyes and color scale fixed onto a single strip, for quick and accurate readings without referring to a seperate scale.

BAKER-pHIX pH Papers with color scale

pH 4.0-9.0 / pH Test Strips with color scale (also suitable for colored solutions)

Product Information (not specifications):

pH Gradation 0.5

 PRODUCT
 PACKING
 CONT.

 NO.
 BOX

 4403
 1 box

Each plastic box contains 200 strips, 11x100 mm. Indicator dyes and color scale fixed onto a single strip, for quick and accurate readings without referring to a seperate scale.

Innovation is principal to our business.



BAKER-pHIX pH Papers with color scale

pH 5.2-6.8 / pH Test Strips with color scale (also suitable for colored solutions)

4404

Product Information (not specifications):

pH Gradation

0.2/0.3

RODUCT	PACKING	CONT.
10.		вох
1404	1 hox	

Each plastic box contains 200 strips, 11x100 mm. Indicator dyes and color scale fixed onto a single strip, for quick and accurate readings without referring to a seperate scale.

BAKER-pHIX pH Papers with color scale

pH 6.0-8.1 / pH Test Strips with color scale (also suitable for colored solutions)

4405

Product Information (not specifications):

pH Gradation

0.3 **NO.**

0.2/0.3

0.2/0.3

1

PRODUCT	PACKING	CONT.	
NO.		вох	
4405	1 box		

Each plastic box contains 200 strips, 11x100 mm. Indicator dyes and color scale fixed onto a single strip, for quick and accurate readings without referring to a seperate scale. pH Gradation of 0.3 units.

BAKER-pHIX pH Papers with color scale

pH 7.2-8.8 / pH Test Strips with color scale (also suitable for colored solutions)

4406

Product Information (not specifications):

pH Gradation

PRODUCT	PACKING	CONT
NO.		ВОХ
4406	1 box	

Each plastic box contains 200 strips, 11x100 mm. Indicator dyes and color scale fixed onto a single strip, for quick and accurate readings without referring to a seperate scale.

BAKER-pHIX pH Papers with color scale

pH 8.0-9.7 / pH Test Strips with color scale (also suitable for colored solutions

4407

Product Information (not specifications):

pH Gradation

PRODUCT	PACKING	CONT
NO.		BOX
4407	1 box	

Each plastic box contains 200 strips, 11x100 mm. Indicator dyes and color scale fixed onto a single strip, for quick and accurate readings without referring to a seperate scale.

BAKER-pHIX pH Indicator Sticks

pH 0-14 / Non-bleeding Test Strips

4390

Product Information (not specifications):

pH Gradation

PRODUCT	PACKING	CONT.	
NO.		вох	
4390	1 box		

Each plastic box contains 100 sticks, 6x85 mm.

The pHIX indicator dyes are chemically bound to the cellulose fibers, avoiding bleeding even in strong alkaline samples.



BAKER-pHIX pH Indicator Sticks

4391

Product Information (not specifications): pH Gradation

0.5

0.3

PRODUCT	PACKING	CONT.
NO.		вох
/301	1 hov	

Each plastic box contains 100 sticks, 6x85 mm. The pHIX indicator dyes are chemically bound to the cellulose fibers, avoiding bleeding even in strong alkaline

BAKER-pHIX pH Indicator Sticks 4392

Product Information (not specifications):

pH Gradation

PACKING 4392 1 box

Each plastic box contains 100 sticks, 6x85 mm. The pHIX indicator dyes are chemically bound to the cellulose fibers, avoiding bleeding even in strong alkaline

BAKER-pHIX pH Indicator Sticks

Product Information (not specifications): pH Gradation

PRODUCT PACKING

1 box Each plastic box contains 100 sticks, 6x85 mm.

CONT.

The pHIX indicator dyes are chemically bound to the cellulose fibers, avoiding bleeding even in strong alkaline

BAKER-pHIX pH Indicator Sticks

4394

Product Information (not specifications):

pH Gradation 0.3/0.5

PRODUCT	PACKING	CONT	
NO.		вох	
4394	1 box		

Each plastic box contains 100 sticks, 6x85 mm. The pHIX indicator dyes are chemically bound to the cellulose fibers, avoiding bleeding even in strong alkaline

BAKER-pHIX pH Indicator Sticks

Product Information (not specifications):

pH Gradation 0.5

PACKING PRODUCT CONT. NO. BOX 4395 1 box

Each plastic box contains 100 sticks, 6x85 mm. The pHIX indicator dyes are chemically bound to the cellulose fibers, avoiding bleeding even in strong alkaline samples

www.jtbaker.com/europe

Bariu



BAKER-pHIX pH Indicator Sticks

4396

Product Information (not specifications):

pH Gradation

PRODUCT	PACKING	CONT	
NO.		ВОХ	
4396	1 hox		

Each plastic box contains 100 sticks, 6x85 mm. The pHIX indicator dyes are chemically bound to the cellulose fibers, avoiding bleeding even in strong alkaline samples

BAKER-pHIX pH Indicator Sticks

4397

Product Information (not specifications):

pH Gradation

0.5 4397

PRODUCT

0.3/0.4

PACKING 1 box

Each plastic box contains 100 sticks, 6x85 mm. The pHIX indicator dyes are chemically bound to the cellulose fibers, avoiding bleeding even in strong alkaline

BAKER-pHIX pH Indicator Sticks

4398

Product Information (not specifications):

pH Gradation

0.2/0.4

1 box

PACKING

PRODUCT

CONT.

Each plastic box contains 100 sticks, 6x85 mm. The pHIX indicator dyes are chemically bound to the cellulose fibers, avoiding bleeding even in strong alkaline

Barbituric Acid, Anhydrous

NHCONHCOCH₂ÇO

M = 128.09 g/mol

CAS NO. 67-52-7

EINECS 200-658-0

NC CODE 2933 54 00

Assay (acidimetric) min. 99% max. 0.001% Heavy Metals (as Pb) Iron (Fe) max. 0.002% Residue after Ignition max. 0.1%

PRODUCT PACKING CONT. NO. BOX 1547.0100 100 g

Barium 1000 μg/ml

'BAKER INSTRA-ANALYZED' / Plasma Standard

5705

1547

M = 137.33 g/mol

NC CODE 3822 00 00

EC NO. 7 004 00 1

R: 36/37/38 S: 26-28-45

Certificate Provided Reporting Actual Lot Analysis Barium (Ba)

998-1002 $\mu g/ml$

5705 0100 100 ml

PACKING

BOX

Prepared from the highest purity raw material available, generally greater than 99.999% spectral purity. The content of the solution is confirmed to be accurate to within \pm 0.2%. Typically 1000 μ g/ml.

The certificate of analysis provided reports actual lot analysis. The certificate also lists the trace impurities.

Innovation is principal to our business.



Bariu

Barium 1000 μg/ml

6920 (Matrix: 2% nitric acid) / 'BAKER INSTRA-ANALYZED' / Atomic Absorption Standard

Barium (Ba)

► Ba M = 137.33 g/mol

NC CODE 3822 00 00

R: 36/38 **S:** 26-37

Barium (Ba)

998-1002 μg/ml

PRODUCT PACKING CONT. NO. вох 6920.0100 100 ml 6920.0500 500 ml

Prepared by dissolution of high purity raw materials (min. 99.99% spectral purity). Assays are verified by ICP against standards traceable to NIST. Standard Reference Material numbers (SRM) are printed on each label.

Barium 1000 μg/ml 'BAKER ANALYZED' / Atomic A

6804

► Ba M = 137.33 g/mol NC CODE 3822 00 00

> R: 36/38 **S:** 26-37

998-1002 μg/ml

CONT **PRODUCT PACKING** NO. вох 6804.0100 100 ml 6804.0500 500 ml

Barium nitrate in nitric acid 0.5 mol/l.

Barium 10000 μ g/ml (Matrix: 2% nitric acid) / "BAKER INSTRA-ANALYZED" / Plasma Standard

► Ba M = 137.33 g/mol NC CODE 3822 00 00 **EC NO.** 7 004 00 1

R: 36/38 **S**: 26-37



Certificate Provided Reporting Actual Lot Analysis Barium (Ba) 9980-10020 µg/ml PRODUCT PACKING CONT. NO. BOX 5719.0100 100 ml

Prepared from the highest purity raw material available, generally greater than 99.999% spectral purity. The content of the solution is confirmed to be accurate to within 0.2 %. Typically 10000 μ g/ml. The certificate of analysis provided reports actual lot

analysis. The certificate also lists the trace impurities.

Barium Chloride

7248 0.5 mol/l / 'BAKER ANALYZED'

► BaCl₂ **M** = 208.25 g/mol

> **CAS NO.** 10326-27-9 NC CODE 2827 39 80

EC NO. 56 002 00 7 **UN/ID NO.** 3287

ADR/RID 6.1 T4

IMDG 6.1/III **R:** 20/22

S: 28

0.495-0.505 Titer (mol/l)

PRODUCT PACKING CONT. вох 7248.9010 10 l

Barium Chloride 7645

0.05 mol/l / 'BAKER ANALYZED'

► BaCl₂ M = 208.25 g/mol **CAS NO.** 10326-27-9

> NC CODE 2827 39 80 **UN/ID NO.** 3287

ADR/RID 6.1 T4 IMDG 6.1/III

R: 20/22

S: 36



0.0495-0.0505 Molarity (M) pH at 20°C 2.95-3.05

CONT. **PRODUCT** PACKING NO. вох 7645.1000 11



Barium Chloride Dihydrate 'BAKER ANALYZED' / ACS

0045

► BaCl₃.2H₃O	
2 2	244.28 g/mol
CAS NO.	10326-27-9
EINECS	233-788-1
NC CODE	2827 39 80
EC NO.	56 002 00 7
UN/ID NO.	1564
ADR/RID	6.1 T5
IMDG	6.1/III
R:	20-25
S:	45

Assay	99.0 - 101.0%
Calcium (Ca)	max. 0.05%
Heavy Metals (as Pb)	max. 5 ppm
Insoluble Matter	max. 0.005%
Iron (Fe)	max. 2 ppm
Loss on Drying at 150°C	14.0-16.0%
Oxidizing Substances (as NO ₃)	max. 0.005%
pH of 5% Solution at 25°C	5.2-8.2
Potassium (K)	max. 0.0025%
Sodium (Na)	max. 0.005%
Strontium (Sr)	max. 0.1%

PRODUCT	PACKING	CONT.
NO.		вох
0045.0100	100 g	
0045.1000	1 kg	6

Barium Nitrate

'BAKER ANALYZED' / ACS

0052

► Ba(NO ₃) ₂	
M =	261.35 g/mol
CAS NO.	10022-31-8
EINECS	233-020-5
NC CODE	2834 29 20
EC NO.	56 002 00 7
UN/ID NO.	1446
ADR/RID	5.1 OT2
IMDG	5.1/II
R:	20/22
S:	28

Assay	min. 99.0%
Calcium (Ca)	max. 0.05%
Insoluble Matter	max. 0.01%
pH of 5% Solution at 25°C	5.0-7.0
Potassium (K)	max. 0.005%
Sodium (Na)	max. 0.005%
Strontium (Sr)	max. 0.05%
Trace Impurities (in ppm):	
Chloride (CI)	max. 5
Heavy Metals (as Pb)	max. 5
Iron (Fe)	max. 2

	PRODUCT	PACKING	CONT.
6	NO.		вох
ó	0052.0500	500 g	
b	0052.9012PE	12 kg PE Pail	
0			
6			
b D			
6			
5			

Barium Sulfate

'BAKER ANALYZED'

M = 233.40 g/mol **CAS NO.** 7727-43-7 **EINECS** 231-784-4 NC CODE 2833 27 00

► BaSO₄

Assay	min. 97.0%
Acidity (as H ₂ SO ₄)	max. 0.01%
Chloride (CI)	max. 0.005%
Heavy Metals (as Pb)	max. 0.001%
Iron (Fe)	max. 0.001%
Phosphate (PO ₄)	max. 0.001%
Water Soluble Salts	max. 0.02%
Trace Impurities (in ppm):	
Arsenic (As)	max. 1

PACKING	CONT.
	вох
500 g	

Barium Sulfate

See for qualities meeting the requirements of cGMP, USP/NF, BP/Ph.Eur, JP or FCC, page 36 $\,$



0053

Basic Fuchsin



Bathophenanthroline

See 4,7-Diphenyl-1,10-phenanthroline



www.jtbaker.com/europe



Benze

Benzene

8014 'BAKER ANALYZED' / ACS

▶ C₆H₆ **M** = 78.11 g/mol **1 l** = 0.88 kg **FLASHPOINT** $-11\,^{\circ}\text{C}$ **CAS NO.** 71-43-2 **EINECS** 200-753-7 NC CODE 2902 20 00 **EC NO.** 601 020 00 8 **UN/ID NO.** 1114 ADR/RID 3F1 IMDG 3/II

R: 11-36/38-45-46-48/23/24/25-65

S: 45-53





Assay	min. 99.0%
Color (APHA)	max. 10
Residue after Evaporation	max. 0.001%
Substances Darkened by H ₂ SO ₄	passes tes
Sulfur Compounds (as S)	max. 0.005%
Thiophene	passes tes
Water (H ₂ O)	max. 0.05%
Trace Impurities (in ppm):	
Aluminium (AI)	max. 0.5
Barium (Ba)	max. 0.1
Boron (B)	max. 0.02
Cadmium (Cd)	max. 0.05
Calcium (Ca)	max. 0.5
Chromium (Cr)	max. 0.02
Cobalt (Co)	max. 0.02
Copper (Cu)	max. 0.02
Iron (Fe)	max. 0.1
Lead (Pb)	max. 0.1
Magnesium (Mg)	max. 0.1
Manganese (Mn)	max. 0.02
Nickel (Ni)	max. 0.02
Tin (Sn)	max. 0.1
Zinc (Zn)	max. 0.1

PRODUCT	PACKING	CONT.
NO.		вох
8014.1000	11	6
8014.2500	2.5	4
8014.9025	25 l	
8014.9200	200 l	

Thiophene Free.



1,2,3-Benzenetriol See Pyrogallol



Benzocaine

See for qualities meeting the requirements of cGMP, USP/NF, BP/Ph.Eur, JP or FCC, page 36

Benzoic Acid

C ₆ H ₅ COOH	
M =	122.12 g/mol
FLASHPOINT	121 °C
CAS NO.	65-85-0
EINECS	200-618-2
NC CODE	2916 31 00
R:	22-36
S:	24
	A A Xn

Meets ACS Specifications	
Assay	min. 99.5%
Chlorine Compounds (as CI)	max. 0.005%
Freezing Point	122.0-123.0°C
Heavy Metals (as Pb)	max. 5 ppm
Insoluble in Methanol	max. 0.005%
Residue after Ignition	max. 0.005%
Substances Reducing KMnO ₄	passes test
Sulfur Compounds (as S)	max. 0.002%

PRODUCT PACKING CONT. BOX 1270.0250 250 g



Benzoic Acid

See for qualities meeting the requirements of cGMP, USP/NF, BP/Ph.Eur, JP or FCC, page 36

Benzotriazole 'BAKER ANALYZED'

C ₆ H₄NHN:N	
M =	119.13 g/mol
CAS NO.	95-14-7
EINECS	202-394-1
NC CODE	2933 99 90
R:	22-36-52/53
S:	22-24
	Xn

Assay	min. 99%
Appearance	passes test
Color of 10% Solution in Ethanol (APHA)	max. 100
Heavy Metals (as Pb)	max. 0.001%
Iron (Fe)	max. 0.001%
Melting Point	min. 95°C
Water (H ₂ O)	max. 0.5%

PRODUCT	PACKING	CONT
NO.		ВОХ
2363.9040	40 kg	

7010

Benzyl Alcohol 'BAKER ANALYZED'

► C₆H₅CH₂OH

M = 108.14 g/mol

1 l = 1.04 kg

FLASHPOINT 101 °C

CAS NO. 100-51-6 **EINECS** 202-859-9

NC CODE 2906 21 00

EC NO. 603 057 00 5

R: 20/22

S: 26



Assay (by GC)	min. 99%
BDA/Diacetal	max. 10 ppm
Boiling Range	202.5-206.5°C
Chloride Compounds (as Cl)	max. 0.02%
Density (g/ml) at 25°C	1.038-1.042
Residue after Ignition	max. 0.005%
Solubility in Water	passes test

PRODUCT NO.	PACKING	CONT. BOX
7010.1000	11	
7010.2500	2.5	

Benzyl Alcohol

C₆H₅CH₂OH

M = 108.14 g/mol

1 l = 1.04 kg

FLASHPOINT 101 °C

CAS NO. 100-51-6

EINECS 202-859-9

NC CODE 2906 21 00

EC NO. 603 057 00 5

R: 20/22

S: 26



PRODUCT PACKING CONT. BOX NO. 3402.9025 25 l Jerrycan

Histo-Grade implicates that this reagent is specially tested and therefore solely intended for use in histo-pathology applications. This reagent is of an analytical quality.

Benzyl Alcohol

See for qualities meeting the requirements of cGMP, USP/NF, BP/Ph.Eur, JP or FCC, page 36

Assay (by GC)



3402

Beryllium 1000 μg/ml

(Matrix: 2% nitric acid) / 'BAKER INSTRA-ANALYZED' / Plasma Standard

5706

► Be

M = 9.01 g/mol

NC CODE 3822 00 00

EC NO. 4 002 00 2

R: 20-36/38-49

S: 26-45-53



Certificate Provided Reporting Actual Lot Analysis Beryllium (Be) 998-1002 µg/ml

CONT. **PRODUCT PACKING** NO. 5706.0100 100 ml

Prepared from the highest purity raw material available, generally greater than 99.999% spectral purity. The content of the solution is confirmed to be accurate to within \pm 0.2%. Typically 1000 $\mu\text{g/ml}.$

The certificate of analysis provided reports actual lot analysis. The certificate also lists the trace impurities.

Beryllium 1000 μg/ml

(Matrix: 2% nitric acid) / 'BAKER INSTRA-ANALYZED' / Atomic Absorption Standard

6921

CONT.

вох

► Be

M = 9.01 g/mol

NC CODE 3822 00 00

EC NO. 4 002 00 2

R: 20-36/38-49 **S:** 26-45-53

Beryllium (Be) 998-1002 μg/ml

min. 99%

6921.0100 100 ml 500 ml 6921 0500

PACKING

PRODUCT

Prepared by dissolution of high purity raw materials (min. 99.99% spectral purity). Assays are verified by ICP against standards traceable to NIST. Standard Reference Material numbers (SRM) are printed on each label.

Beryl

Beryllium 1000 μ g/ml

BAKER ANALYZED' / Ato 6805

► Be **M** = 9.01 g/mol

NC CODE 3822 00 00

R: 20-36/38-49 **S:** 26-45-53

Beryllium (Be) 995 - $1005~\mu g/ml$

PRODUCT CONT. **PACKING** NO. вох 6805.0100 100 ml 6805.0500 500 ml

Basic Beryllium Acetate in Nitric Acid 0.5 mol/l.

PACKING

Beryllium 10000 μg/ml

(Matrix: 5% nitric acid) / 'BAKER INSTRA-ANALYZED' / Plasma Standard

► Be

M = 9.01 g/mol

NC CODE 3822 00 00 **EC NO.** 4 002 00 2

UN/ID NO. 2031 ADR/RID 8 CO1

IMDG 8/II

R: 22-23-43-48/20-49-52/53

S: 36/37-45-53



Certificate Provided Reporting Actual Lot Analysis Beryllium (Be)

9980-10020 μg/ml

NO. 5720.0100 100 ml

PRODUCT

CONT. вох

Prepared from the highest purity raw material available, generally greater than 99.999% spectral purity. The content of the solution is confirmed to be accurate to within 0.2 %. Typically 10000 μ g/ml.

The certificate of analysis provided reports actual lot analysis. The certificate also lists the trace impurities.



Biopharmaceutical Products



Biotechnology Products

ation section Reagents for Molecular and Biotechnology, page 26

5707

Bismuth 1000 μg/ml

 $0.10\% \; \text{(w/v)} \; / \; \text{(Matrix: } 2\% \; \text{nitric acid)} \; / \; \text{'BAKER INSTRA-ANALYZED'} \; / \; \text{Plasma Standard}$

► Bi

M = 208.98 g/mol

NC CODE 3822 00 00

R: 36/38

S: 26



Certificate Provided Reporting Actual Lot Analysis Bismuth (Bi) 998-1002 $\mu g/ml$

NO. 5707.0100 100 ml

PRODUCT

PACKING CONT. BOX

Prepared from the highest purity raw material available, generally greater than 99.999% spectral purity. The content of the solution is confirmed to be accurate to within \pm 0.2%. Typically 1000 μ g/ml. The certificate of analysis provided reports actual lot

analysis. The certificate also lists the trace impurities.

Bismuth 1000 μg/ml

6922 (Matrix: 2% nitric acid) / 'BAKER INSTRA-ANALYZED' / Atomic Absorption Standard

► Bi

M = 208.98 g/mol

NC CODE 3822 00 00

R: 36/38



998-1002 μg/ml Bismuth (Bi)

PRODUCT PACKING CONT. NO. BOX 6922.0100 100 ml 6922.0500 500 ml

Prepared by dissolution of high purity raw materials (min. 99.99% spectral purity). Assays are verified by ICP against standards traceable to NIST. Standard Reference Material numbers (SRM) are printed on each label.

Bismuth 1000 μg/ml 'BAKER ANALYZED' / Atomic Abs

6806

CONT.

вох

► Bi

M = 208.98 g/mol

NC CODE 3822 00 00

R: 36/38 **S**: 26

998-1002 $\mu g/ml$ Bismuth (Bi)

PRODUCT PACKING NO.

6806.0100 100 ml 6806.0500 500 ml

Bismuth(III)nitrate in nitric acid 0.5 mol/l.

Bismuth 10000 μg/ml (Matrix: 10% nitric acid) / 'BAKER INSTRA-ANALYZED' Plasma Standard

5721

1029

CONT.

вох

► Bi

M = 208.98 g/mol

NC CODE 3822 00 00 **UN/ID NO.** 2031 ADR/RID 8 CO1

IMDG 8/II **R:** 34

S: 20-23-26-36/37/39-45



Certificate Provided Reporting Actual Lot Analysis Bismuth (Bi) 9980-10020 µg/ml

PRODUCT NO.

PACKING 5721.0100 100 ml

Prepared from the highest purity raw material available, generally greater than 99.999% spectral purity. The content of the solution is confirmed to be accurate to within 0.2 %. Typically 10000 μ g/ml.

The certificate of analysis provided reports actual lot analysis. The certificate also lists the trace impurities.

Bismuth(III) Nitrate Pentahydrate

► Bi(NO₃)₃.5H₂O

M = 485.07 g/mol **CAS NO.** 10035-06-0

EINECS 233-791-8 NC CODE 2834 29 80

UN/ID NO. 1477

ADR/RID 5.1 O2

IMDG 5.1/II

Assay (by EDTA titrn.)	min. 99.0%
Arsenic (As)	max. 0.001%
Chloride (CI)	max. 0.001%
Copper (Cu) (by AAS)	max. 0.002%
Insoluble in HNO ₃	max. 0.005%
Iron (Fe) (by AAS)	max. 0.001%
Lead (Pb) (by AAS)	max. 0.002%
Silver (Ag)	max. 0.001%
Sulfate (SO ₄)	max. 0.002%

PRODUCT	PACKING	CONT.
NO.		вох
1029.0500	500 g	

Bismuth(III) Oxide

► Bi₂O₃

M = 465.96 g/mol **CAS NO.** 1304-76-3 **EINECS** 215-134-7 NC CODE 2825 90 80

Assay (by EDTA titrn.) min. 99.0% Chloride (CI) max. 0.001% Iron (Fe) (by AAS) max. 0.005% Lead (Pb) (by AAS) max. 0.003% Nitrate (NO₃) max. 0.03% Sulfur Compounds (as SO₄) max. 0.002%

вох

Bismuth Subgallate

See for qualities meeting the requirements of cGMP, USP/NF, BP/Ph.Eur, JP or FCC, page 36



1032

Bismuth Subnitrate

See for qualities meeting the requirements of cGMP, USP/NF, BP/Ph.Eur, JP or FCC, page 36



Bismuth Subsalicylate

ents of cGMP, USP/NF, BP/Ph.Eur, JP or FCC, page 36



Вогах

um Tetraborate Decahydrate





Вогіс

Boric Acid

5168 **ULTREX Ultrapure Reagent**

► H₃BO₃

M = 61.83 g/mol **CAS NO.** 10043-35-3

EINECS 233-139-2 NC CODE 2810 00 90

R: 36/37/38-60

S: 45-53



Certificate Provided Reporting Actual Lot Analysis Actual Lot Analysis Lot No. Y36606

Assay (by acid-base titration) 99.9% 0.004% Nonvolatile with Methanol Particulate Matter < 0.002%

Metallic Impurities in parts per million ($\mu g/g$):

Aluminium (AI) < 20 Barium (Ba) < 20 Calcium (Ca) < 20 Chromium (Cr) < 20 Cobalt (Co) < 20 < 10 Copper (Cu) 120 Iron (Fe) Manganese (Mn) < 20 Nickel (Ni) < 10 Tin (Sn) < 100 Titanium (Ti) < 20 Vanadium (V) < 10 Zirconium (Zr) < 10 Non-Metallic Impurities in parts per million ($\mu g/g$):

< 0.04

Halide (as CI) 2 Phosphate (PO₄) < 0.4 Sulfur Compounds (as SO₄)

PRODUCT **PACKING** CONT. NO. вох 5168.0100 100 g PE

For application in Fiber Optics Manufacture.

Boric Acid

'BAKER ULTRAPURE BIOREAGENT'

► H₃BO₃

M = 61.83 g/mol **CAS NO.** 10043-35-3 **EINECS** 233-139-2

NC CODE 2810 00 90 **R:** 36/37/38-60

S: 45-53



For Molecular Biology Buffers	
Assay	min. 99.5%
Calcium (Ca)	max. 0.005%
Chloride (CI)	max. 0.001%
DNAase	none detected
Heavy Metals (as Pb)	max. 0.001%
Insoluble in Methanol	max. 0.005%
Iron (Fe)	max. 0.001%
Phosphate (PO ₄)	max. 0.001%
Protease	none detected
RNAase	none detected
Sulfate (SO ₄)	max. 0.005%
Trace Impurities (in ppm):	
Arsenic (As)	max. 1

Magnesium (Mg)

вох NO. 4035.0500 500 g 2.5 kg 4035.2500

CONT.

PACKING

PRODUCT

max. 5

Boric Acid 0055 'BAKER ANALYZED' / ACS

► H₃BO₃

M = 61.83 g/mol

CAS NO. 10043-35-3

EINECS 233-139-2 NC CODE 2810 00 90

R: 36/37/38-60

S: 45-53

Exceeds ACS Specifications	
Assay	min. 99.5%
Arsenic (As)	max. 1 ppm
Calcium (Ca)	max. 0.005%
Chloride (CI)	max. 0.001%
Heavy Metals (as Pb)	max. 0.001%
Insoluble in Methanol	max. 0.005%
Iron (Fe)	max. 0.001%
Nonvolatile with Methanol	max. 0.05%
Phosphate (PO ₄)	max. 0.001%
Sulfate (SO ₄)	max. 0.005%

PRODUCT	PACKING	CONT.
NO.		вох
0055.0500	500 g	6
0055.1000	1 kg	6
0055.9010	10 kg	
0055.9025	25 kg	
0055.9050	50 kg	

Boron



0501 'BAKER'

► H₃BO₃

M = 61.83 g/mol **CAS NO.** 10043-35-3

EINECS 233-139-2

NC CODE 2810 00 90 **R:** 36/37/38-60

S: 45-53



Appearance of solution	passes test
Assay	99.0-100.5%
Heavy Metals (as Pb)	max. 15 ppm
Identification	passes test
Loss on Drying	max. 0.5%
Organic Matter	passes test
рН	3.8-4.8
Solubility in Alcohol	passes test
Sulfates (as SO ₄)	max. 450 ppm

PRODUCT	PACKING	CONT.
NO.		вох
0501.1000	1 kg	6
0501.9050	50 kg	

Stored in a well closed container.

Boric Acid

See for qualities meeting the requirements of cGMP, USP/NF, BP/Ph.Eur, JP or FCC, page 36



Boric Anhydride



▶ B₂O₃

M = 69.62 g/mol CAS NO. 1303-86-2 **EINECS** 215-125-8

NC CODE 2810 00 10

Assay (B₂O₃) (ignited basis) min. 98.0% Loss on Ignition max. 10.0% max. 0.05% Sulfate (SO₄)

PRODUCT PACKING CONT. NO. BOX 1176 0500 500 g

Boron 1000 µg/ml

(Matrix 2% ammonium hydroxide) / 'BAKER INSTRA-ANALYZED' / Plasma Standard

5708

▶ B

M = 10.81 g/mol NC CODE 3822 00 00

Certificate Provided Reporting Actual Lot Analysis 998-1002 μg/ml **PRODUCT PACKING** CONT. вох 5708.0100 100 ml

Prepared from the highest purity raw material available, generally greater than 99.999% spectral purity. The content of the solution is confirmed to be accurate to within \pm 0.2%. Typically 1000 μ g/ml.

The certificate of analysis provided reports actual lot analysis. The certificate also lists the trace impurities.

Boron 1000 μg/ml

'BAKER INSTRA-ANALYZED' / Atomic Absorption Standard

6923

▶ B

M = 10.81 g/mol NC CODE 3822 00 00

Boron (B) 998-1002 µg/ml **PRODUCT PACKING** CONT. BOX 6923.0100 100 ml 6923.0500 500 ml

Prepared by dissolution of high purity raw materials (min. 99.99% spectral purity). Assays are verified by ICP against standards traceable to NIST. Standard Reference Material numbers (SRM) are printed on each label.

Boron 10000 μ g/ml

(Matrix 2% ammonium hydroxide) / 'BAKER INSTRA-ANALYZED' / Plasma Standard

5722

▶ B

M = 10.81 g/mol NC CODE 3822 00 00

Certificate Provided Reporting Actual Lot Analysis 9980-10020 $\mu g/ml$ **PRODUCT PACKING** CONT. 5722.0100 100 ml

Prepared from the highest purity raw material available, generally greater than 99.999% spectral purity. The content of the solution is confirmed to be accurate to within 0.2 %. Typically 10000 μ g/ml.

The certificate of analysis provided reports actual lot analysis. The certificate also lists the trace impurities.



Bouin

Bouin's Fixative

3880 HISTO GRADE

NC CODE 3822 00 00

R: 20/21/22-36/37/38-40-43

S: 26-36/37/39



Fixative according to Bouin

PRODUCT PACKING CONT. NO. вох 3880.1000 1 | Glass

Ready to use Fixative for use in Histo-Pathology



Bovine Albumin, Fraction V

See Albumin (Bovine) Fraction V

1240

Brij 35

CAS NO. 9002-92-0 NC CODE 3402 13 00

R: 22-36

S: 24

passes test

min. 99%

passes test

max. 0.002%

max. 0.02%

PRODUCT

9094.0100

9094.0500

PACKING

100 ml

500 ml

CONT.

вох

PRODUCT PACKING CONT. BOX NO. 1240.0250 250 g

9003

Bromine

'BAKER ANALYZED' / ACS

► Br

M = 79.90 g/mol

1 l = 3.11 kg **CAS NO.** 7726-95-6

EINECS 231-778-1 NC CODE 2801 30 90

EC NO. 35 001 00 5

UN/ID NO. 1744 ADR/RID 8 CT1 IMDG 8/I

> **R:** 26-35-50 **S:** 26-45-61-7/9







Meets ACS Specifications Assay min. 99.5% Chlorine (CI) max. 0.05% Organic Bromine Compounds passes test Sulfur Compounds (as S) max. 0.001% Additional Specification(s): lodine (I) max. 0.001% max. 0.005% Residue after Evaporation Trace Impurities (in ppm): Heavy Metals (as Pb) max. 2 Nickel (Ni) max. 5

Appearance of solution

PRODUCT PACKING CONT. NO. вох 9003.0100 100 ml 9003.0500 500 ml

Bromine

9094

▶ Br

M = 79.90 g/mol **1 l** = 3.11 kg

CAS NO. 7726-95-6

EINECS 231-778-1

NC CODE 2801 30 90 **EC NO.** 35 001 00 5

UN/ID NO. 1744 ADR/RID 8 CT1

IMDG 8/I

R: 26-35-50

S: 26-45-61-7/9





environment



Assay (Br)

lodine (I)

Appearance

Residue after Evaporation

Bromine (Bromide-Bromate) 0.05 mol Br₂/l

BAKER ANALYZED

7123

► Br₂

NC CODE 2801 30 90

R: 45

S: 37-45-53



PRODUCT PACKING CONT. NO. вох 7123.1000 11

Volumetric Solution, ready for use

PRODUCT

NO.

4708

passes test

min. 95.0%

max. 0.02%

min. 2.80

PRODUCT

7354 0500

PACKING

500 ml

Bromine (Bromide-Bromate) 0.05 mol Br $_2$ /l $_{^{1/}_{10}}$ equiv. = 7.99 g; 0.1N / DILUT-IT

4708

CONT

вох



NC CODE 2801 30 90

R: 45

S: 37-45-53



Volumetric Concentrate, for dilution to 11.

PACKING

1 amp.

Bromocresol Green

'BAKER ANALYZED' / ACS

1511

1330

7354

1369

CONT.

BOX

Ç₆H₄SO₂OÇ(C₆H-2-CH₃-3,5-Br₂-4-OH)₂

M = 698.02 g/mol **CAS NO.** 76-60-8

EINECS 200-972-8 NC CODE 2934 99 90 Clarity of Solution

Assay (by GC)

Density (g/ml) at 25°C

Residue after Evaporation

Visual Transition Interval: pH 3.8

Meets ACS Specifications

yellow pH 5.4 blue

PRODUCT PACKING CONT. NO. вох 1511.0005 5 g

Bromocresol Purple

 $C_6H_4SO_2OC(C_6H_2-3-CH_3-4-OH-5-Br)_2$

M = 540.23 g/mol **CAS NO.** 115-40-2 **EINECS** 204-087-8

Insoluble Matter passes test Visual Transition Interval: pH 5.2 vellow pH 6.8 purple **PRODUCT PACKING** CONT. NO. BOX 1330.0005 5 g

Bromoform

► CHBr₃

M = 252.73 g/mol

11 = 2.82 kg

NC CODE 2934 99 90

CAS NO. 75-25-2

EINECS 200-854-6

NC CODE 2903 30 36

EC NO. 602 007 00 0

UN/ID NO. 2515

ADR/RID 6.1 T1

IMDG 6.1/III

R: 23-36/38-51/53

S: 28-45-61





Bromophenol Blue

'BAKER ANALYZED' / ACS

Meets ACS Specifications	
Clarity of Solution	passes test
Visual Transition Interval:	
oH 3.0	yellow
oH 4.6	blue

PACKING BOX 1369 0005 5 g

Ç₆H₄SO₂ΟÇ(C₆H₂-3,5-Br₂-4-OH)₂

M = 669.97 g/mol

CAS NO. 115-39-9

EINECS 204-086-2

NC CODE 2934 99 90



Bromo

Bromophenol Blue Sodium Salt 'BAKER ANALYZED' / ACS

1472

► C₁9H₉Br₄O₅SNa **M** = 691.95 g/mol **CAS NO.** 34725-61-6 **EINECS** 252-170-2 NC CODE 2934 99 90

Meets ACS Specifications Clarity of Solution passes test Visual Transition Interval: pH 3.0 yellow pH 4.6 blue

NO.		BOX
1472.0010	10 g	

Bromothymol Blue

1513 'BAKER ANALYZED' / ACS

 $ightharpoonup C_6H_4SO_2OC[C_6H-2-CH_3-3-Br-4-OH-5-CH(CH_3)_2]_2$ **M =** 624.39 g/mol **CAS NO.** 76-59-5 **EINECS** 200-971-2 NC CODE 2934 99 90

Meets ACS Specifications Clarity of Solution passes test **Visual Transition Interval:** pH 6.0 yellow pH 7.6 blue

PRODUCT	PACKING	CONT.
NO.		вох
1513.0005	5 g	
1513.0025	25 g Glass	

Brucine

2459

C₂₃H₂₆N₂O₄

M = 394.48 g/mol **CAS NO.** 357-57-3 **EINECS** 206-614-7 NC CODE 2939 99 00 **EC NO.** 614 006 00 1 **UN/ID NO.** 1570 **ADR/RID** 6.1 T2

IMDG 6.1/I

R: 26/28-52/53 **S:** 13-45-61

Assay (by Perchloric Acid titrn.) min. 99.5% Melting Point 176-179°C Nitrate (NO₃) passes test Sensitivity to Nitrate passes test Specific Rotation [a] 20 -83 to -86° Sulfated Ash max. 0.1% Water (H₂O) max. 0.5%

BOX 2459.0010 10 g

CONT.

PACKING



Buffer

Buffer Concentrate

NC CODE 3822 00 00

pH at 20°C $1.00\,\pm\,0.02$ PRODUCT PACKING CONT. NO. вох 4850 1 amp.

Volumetric concentrate, for dilution to 500ml. Contains potassiumchloride / hydrochloric acid. This buffer is traceable to NIST Standard Reference

4851

Buffer Concentrate

pH 2 / DILUT-IT

NC CODE 3822 00 00

pH at 20°C $2.00\,\pm\,0.02$ **PRODUCT PACKING** CONT. NO. BOX 4851 1 amp.

Volumetric concentrate, for dilution to 500ml. Contains citric acid / hydrochloric acid. This buffer is traceable to NIST Standard Reference Buffer.



	Buffer Concentrate pH3 / DILUT-IT					4852
NC CODE	3822 00 00	pH at 20°C	3.00 ± 0.02	PRODUCT NO.	PACKING	CON BO
				4852	1 amp.	
				Contains cit	concentrate, for dilution to ric acid / sodium hydroxide is traceable to NIST Standa	e / sodium chloride
	Buffer Concentrate pH 4 / DILUT-IT					4795
NC CODE	3822 00 00	pH at 25°C	4.00 ± 0.01	PRODUCT	PACKING	CONT
	0022 00 00	Composition:	1.00 ± 0.01	NO.	T ACKING	ВО
		Buffer Capacity, B	0.016	4795	1 amp.	
		Dilution value, Δ pH $^1/_2$	+0.05			
		Ionic Strength, I	0.05	Volumetric	concentrate, for dilution to	500ml.
		NaOOC.C ₆ H ₄ COOH	0.05 M	Contains ph	ntalic acid / sodium hydroxi	ide.
		Temperature Coefficient, d (pH)/dT (25°C)	+0.001		is traceable to NIST Standa	rd Reference
		Useful Buffer Range, pH	3.0-5.0	Buffer.		
	Buffer Concentrate pH 4 RED / DILUT-IT					4860
NC CODE	3822 00 00	pH at 20°C	4.00 ± 0.02	PRODUCT	PACKING	CONT
				NO.		B0
				4860	1 amp.	
				Contains po	concentrate, for dilution to stassium hydrogen phtalate is traceable to NIST Standa	e / coloring agent.
	Buffer Concentrate ph 5.0 / 'BAKER ANALYZED'					7541
NC CODE	3822 00 00	pH 20°C (after dilution 7ml to 1L)	5.0 ± 0.05	PRODUCT	PACKING	CONT
		,		NO.		B0
				7541.9020	20 l Polycube	
					8 ml NaOH (50%) and 8.59 with water to 20 L.	02 kg NaH₂PO₄·2aq
	Buffer Concentrate ph 5 / DILUT-IT					4853
NC CODE	3822 00 00	pH at 20°C	5.00 ± 0.02	PRODUCT	PACKING	CONT
	0022 00 00	p11 4(20 0	0.00 ± 0.02	NO.		BOX
				4853	1 amp.	(
				Contains cit	concentrate, for dilution to ric acid / sodium hydroxidi is traceable to NIST Standa	е.
	Buffer Concentrate					4854
NC CODF	3822 00 00	pH at 20°C	6.00 ± 0.02	PRODUCT	PACKING	CONT
				NO.		B0
				4854	1 amp.	(
				Contains cit	concentrate, for dilution to ric acid / sodium hydroxid is traceable to NIST Standa	е.

This buffer is traceable to NIST Standard Reference

Buffer.



Buffe

Buffer Concentrate

nH 7 / DII UT-IT

4796

NC CODE 3822 00 00 $7.00\,\pm\,0.01$ pH at 25°C Composition: Buffer Capacity, B 0.028 Dilution value, $\Delta~\text{pH}^{\text{1}}\text{/}_{\text{2}}$ +0.08Ionic Strength, I 0.1 Na₂H PO₄ 0.029 M NaH₂PO₄ 0.021 M Temperature Coefficient, d (pH)/dT (25°C) -0.002 Useful Buffer Range, pH

PACKING	CONT.
	вох
1 amp.	6

Volumetric concentrate, for dilution to 500ml. Contains sodiumdihydrogenphosphate / disodium hydrogen phosphate.This buffer is traceable to NIST Standard Reference

Buffer Concentrate

4861 pH 7 GREEN / DILUT-IT

NC CODE 3822 00 00 pH at 20°C

PRODUCT	PACKING	CONT.
NO.		вох
4861	1 amp.	6

Volumetric concentrate, for dilution to 500ml. Contains potassium dihydrogen phosphate / disodiumhydrogen phosphate / coloring agent. This buffer is traceable to NIST Standard Reference Buffer.

Buffer Concentrate 4863

NC CODE 3822 00 00

pH 7.20, according to Weise / DILUT-IT

ţ	oH at 20°C	7.20 ± 0.05	PRODUCT	PACKING	CONT.
			NO.		BOX
			1042	1 2000	

6.5-7.5

 $7.00\,\pm\,0.02$

Volumetric Concentrate, for dilution to 1 l. Contains potassiumdihydrogenphosphate / disodiumhydrogenphosphate. This buffer is traceable to NIST Standard Reference Buffer.

4856

Buffer Concentrate

NC CODE 3822 00 00 pH at 20°C $8.00\,\pm\,0.02$

PRODUCT	PACKING	CONT.
NO.		вох
1856	1 amn	6

Volumetric concentrate, for dilution to 500ml. Contains sodium tetraborate / hydrochloric acid. This buffer is traceable to NIST Standard Reference

4878

Buffer Concentrate

pH9 / DILUT-IT

NC CODE 3822 00 00 pH at 20°C 9.00 ± 0.02

PRODUCT	PACKING	CONT.
NO.		вох
4878	1 amp.	6

Volumetric concentrate, for dilution to 500ml. Contains hydrochloric acid / disodium tetraborate. This buffer is traceable to NIST Standard Reference Buffer.



Buffer Concentrate pH 9 BLUE / DILUT-IT

4862

NC CODE 3822 00 00 pH at 20°C 9.00 ± 0.02 PRODUCT PACKING CONT. вох NO. 4862 1 amp. Volumetric concentrate, for dilution to 500ml. Contains sodium hydroxide / potassium chloride / boric acid / coloring agent.

Buffer Concentrate pH 10 / DILUT-IT

4797

NC CODE 3822 00 00	pH at 25°C	10.00 ± 0.01
	Composition:	
	Buffer Capacity, B	0.030
	Dilution value, Δ pH $^1/_2$	+0.08
	Ionic Strength, I	0.1
	Na ₂ CO ₃	0.025 M
	NaHCO ₃	0.025 M
	Temperature Coefficient, d (pH)/dT (25°C)	-0.010
	Useful Buffer Range, pH	9.5-10.5

PRODUCT	PACKING	CONT.
NO.		вох
4797	1 amp.	6

This buffer is traceable to NIST Standard Reference

Buffer.

Volumetric concentrate, for dilution to 500ml. Contains sodiumcarbonate / sodiumhydrogencarbonate. This buffer is traceable to NIST Standard Reference

Buffer Concentrate pH 11 / DILUT-IT

NC CODE 3822 00 00	pH at 20°C	11.00 ± 0.02	PRODUCT	PACKING	CONT.
			NO.		ВОХ
			4857	1 amp.	6

Volumetric concentrate, for dilution to 500ml. Contains sodium hydroxide / potassium chloride / boric

This buffer is traceable to NIST Standard Reference Buffer.

Buffer Concentrate

12.00 + 0.02 **PRODUCT** pH at 20°C

4858

4857

NC CODE 3822 00 00	pH at 20°C	12.00 ± 0.02	PRODUCT	PACKING	CONT.
			NO.		вох
			4858	1 amp.	6
			Contains sod phosphate.	oncentrate, for dilution to 500ml. ium hydroxide / disodium hydrogen traceable to NIST Standard Reference	ŧ

Buffered Oxide Etchants



Buffer Solution

NC CODE 3822 00 00	pH at 20°C	1.00 ± 0.02	PRODUCT	PACKING	CON
	Temperature (°C):		NO.		ВО
	0	pH 0.94	7261.1000	11	
	10	pH 0.99			
	20	pH 1.00	Contains po	tassiumchloride / hydrochl	oric acid.
	30	pH 1.00	This buffer is	s traceable to NIST Standar	rd Reference
	40	pH 1.01	Buffer.		
	50	pH 1.01			
	60	pH 1.01			
	70	pH 1.01			
	80	pH 1.02			
	90	pH 1.02			



Buffe

7144

Buffer Solution pH 2.00 [Citrate/HCI] / 'BAKER ANALYZED'

NC CODE 3822 00 00	pH at 20°C	2.00 ± 0.02
	Temperature (°C):	
	0	pH 1.99
	10	pH 1.99
	20	pH 2.00
	30	pH 2.00
	40	pH 2.00
	50	pH 2.00
	60	pH 2.00
	70	pH 2.00
	80	pH 2.00

90

PRODUCT	PACKING	CONT
NO.		B0)
7144.0500	500 ml	
7144.1000	11	6
Contains soc	lium citrate / hydrochloric ad traceable to NIST Standard	cid.

pH 2.00

Buffer Solution 7159

pH 3.00 [Citrate/HCI] / 'BAKER ANALYZED'

NC CODE 3822 00 00	pH at 20°C	3.00 ± 0.02
	Temperature (°C):	
	0	pH 3.03
	10	pH 3.02
	20	pH 3.00
	30	pH 3.00
	40	pH 2.99
	50	pH 2.98
	60	pH 2.98
	70	pH 2.98
	80	pH 2.98
	90	pH 2.97

l Polycube		ВО
,		
,		
adium banzas		
citrate / hydro	chloric ac	
		eable to NIST Standard

7262

Buffer Solution

pH 4.00 [Citrate/HCI] / 'BAKER ANALYZED'

NC CODE 3822 00 00	pH at 20°C	4.00 ± 0.02
	Temperature (°C):	
	0	pH 4.03
	10	pH 4.02
	20	pH 4.00
	30	pH 4.00
	40	pH 4.00
	50	pH 4.00
	60	pH 4.00
	70	pH 4.00
	80	pH 4.00

90

0 ± 0.02	PRODUCT	PACKING	CONT.
	NO.		BOX
pH 4.03	7262.1000	11	6
pH 4.02	7262.5000	5 I HDPE	
pH 4.00	7262.9010	10 l	
pH 4.00			
pH 4.00	Contains 0.0	1% sodium benzoate.	
pH 4.00	Contains soc	dium citrate / hydrochloric acid.	
pH 4.00	This buffer is	s traceable to NIST Standard Reference	e
pH 4.00	Buffer.		
pH 4.00			
pH 4.00			

5657

Buffer Solution pH 4.00 [Phthalate] / 'BAKER ANALYZED'

NC CODE 3822 00 00	pH at 25°C	4.00 ± 0.02
	Temperature (°C):	
	0	pH 4.00
	10	pH 4.00
	20	pH 4.00
	30	pH 4.01
	40	pH 4.03
	50	pH 4.05
	60	pH 4.08
	70	pH 4.12
	80	pH 4.16
	90	pH 4.21

PRODUCT	PACKING	CONT.
NO.		вох
5657.0500	500 ml	6
5657.1000	11	6
5657.2500PE	2.5 HDPE	4
5657.5000PE	5 I HDPE	4
5657.9010	10 l Polycube	

This buffer is traceable to NIST Standard Reference

Buffer.



Buffer Solution pH 4.00 [Phthalate] RED / 'BAKER ANALYZED'

5654

NC CODE 3822 00 00	pH at 20°C	4.00 ± 0.02	PRODUCT	PACKING	CONT.
	Temperature (°C):		NO.		вох
	0	pH 4.00	5654.0500	500 ml	6
	10	pH 4.00	5654.1000	1 HDPE	6
	20	pH 4.00	5654.9010	10 l	
	30	pH 4.01			
	40	pH 4.03	Contains 0.0	1% sodium benzoate.	
	50	pH 4.05	Contains pot	tassium hydrogen phtalate.	
	60	pH 4.08	This buffer is	s traceable to NIST Standard	d Reference
	70	pH 4.12	Buffer.		
	80	pH 4.16			
	90	pH 4.21			

Buffer Solution

pH 4.62 [Acetate] / 'BAKER ANALYZED'

7461

NC CODE 3822 00 00	pH at 20°C	4.62 ± 0.02	PRODUCT	PACKING	CONT.
	Temperature (°C):		NO.		вох
	20	pH 4.62	7461.0500	500 ml	
	30	pH 4.62	7461.1000	11	
	40	pH 4.62			
	50	pH 4.65	Contains 0.0	1% sodium benzoate.	
	60	pH 4.68	Contains so	dium acetate / acetic acid.	
	70	pH 4.69	This buffer is	traceable to NIST Standard Referen	се
	80	pH 4.71	Buffer.		
	90	pH 4.75			

Buffer Solution

NC CODE 3822 00 00	pH at 20°C	5.00 ± 0.05	PRODUCT	PACKING	CON
	Temperature (°C):		NO.		В
	0	pH 5.05	7263.1000	11	
	10	pH 5.02			
	20	pH 5.00	Contains citi	ric acid / sodium hydroxide.	
	30	pH 5.00		s traceable to NIST Standard Referen	ice
	40	pH 5.00			
	50	pH 5.02			
	60	pH 5.04			
	70	pH 5.07			
	80	pH 5.10			
	90	pH 5.13			

Buffer Solution

pH 6.00 [Citrate] / 'BAKER ANALYZED'

7264

NC CODE 3822 00 00	pH at 20°C	6.00 ± 0.05	PRODUCT	PACKING	(
	Temperature (°C):		NO.		
	0	pH 6.03	7264.1000	11	
	10	pH 6.01			
	20	pH 6.00	Contains cit	ric acid / sodium hydroxide	
	30	pH 6.02		d Reference	
	40	pH 6.04			
	50	pH 6.06			
	60	pH 6.09			
	70	pH 6.13			
	80	pH 6.18			
	90	pH 6.24			



Buffe

7462

Buffer Solution

pH 6.88 / 'BAKER ANALYZED'

NC CODE 3822 00 00 pH at 25°C 6.88 ± 0.02 Temperature (°C): 0 pH 6.98 10 pH 6.92 pH 6.88 20 30 pH 6.85 40 pH 6.83 50 pH 6.83 pH 6.84 60 70 pH 6.85 80 pH 6.86

90

PRODUCT	PACKING	CONT.
NO.		вох
7462.0500	500 ml	
7462.1000	11	6

Contains 0.01% sodium benzoate.

Contains potassiumdihydrogenphosphate / disodiumhydrogenphosphate.

This buffer is traceable to NIST Standard Reference

Buffer Solution
pH 7.00 [Phosphate] / 'BAKER ANALYZED'

NC CODE 3822 00 00 pH at 20°C 7.00 ± 0.02 Temperature (°C): 0 pH 7.13 10 pH 7.05 20 pH 7.00 pH 6.98 25 30 pH 6.98 40 pH 6.97 50 pH 6.96 pH 6.96 60 70 pH 6.97 80 pH 6.98

90

PRODUCT	PACKING	CONT.
NO.		вох
5656.0500	500 ml	6
5656.1000	11	6
5656.2500PE	2.5 I HDPE	4
5656.5000PE	5 I HDPE	4
5656.9010	10 l Polycube	

Contains 0.01% sodium benzoate.

Contains potassiumdihydrogenphosphate /
disodiumhydrogenphosphate.

This buffer is traceable to NIST Standard Reference

Buffer.

pH 7.00

pH 7.00

pH 6.88

Buffer Solution

5653

pH 7.00 [Phosphate] GREEN / 'BAKER ANALYZED'

NC CODE 3822 00 00 $7.00\,\pm\,0.02$ Temperature (°C): pH 7.13 0 10 pH 7.05 20 pH 7.00 25 pH 6.98 30 pH 6.98 40 pH 6.97 50 pH 6.96 60 pH 6.96 70 pH 6.97 80 pH 6.98

90

PRODUCI	PACKING	CON1.
NO.		вох
5653.0500	500 ml	6
5653.1000	1 HDPE	6
5653.9010	10 l	
5653.9020	20 l Polycube	

Contains potassiumdihydrogenphosphate /

disodium hydrogen phosphate.

This buffer is traceable to NIST Standard Reference Buffer.

Buffer Solution

7265 pH 8.00 [Phosphate] / 'BAKER ANALYZED'

NC CODE 3822 00 00	pH at 20°C	8.00 ± 0.05
	Temperature (°C):	
	0	pH 8.15
	10	pH 8.10
	20	pH 8.00
	30	pH 7.94
	40	pH 7.90
	50	pH 7.85
	60	pH 7.83
	70	pH 7.80
	80	pH 7.78
	90	pH 7.75

PRODUCT	PACKING	CONT.
NO.		вох
7265.1000	11	6

Contains potassium dihydrogen phosphate / sodium hydroxide.

This buffer is traceable to NIST Standard Reference Buffer.

7145

Buffer Solution
pH 9.00 [Borate] / 'BAKER ANALYZED'

NC CODE 3822 00 00	pH at 20°C	9.00 ± 0.02	PRODUCT	PACKING	CO
	Temperature (°C):		NO.		E
	0	pH 9.23	7145.0500	500 ml	
	10	pH 9.10	7145.1000	11	
	20	pH 9.00	7145.2500PE	2.5 I HDPE	
	30	pH 8.91	7145.5000PE	5 I HDPE	
	40	pH 8.84	7145.9020	20 l Polycube	
	50	pH 8.78			
	60	pH 8.73	Contains 0.01	% sodium benzoate.	
	70	pH 8.69	Contains bori	c acid / potassium chloride / sodi	ium
	80	pH 8.66	hydroxide.		
	90	pH 8.62	This buffer is	traceable to NIST Standard Refer	ence
			Buffer.		

Buffer Solution
pH 9.00 [Borate] BLUE / 'BAKER ANALYZED'

5652

NC CODE 3822 00 00	pH at 20°C	9.00 ± 0.02	PRODUCT	PACKING	CONT.	
	Temperature (°C):		NO.		вох	
		0	pH 9.23	5652.0500	500 ml	6
		10	pH 9.10	5652.1000	1 I HDPE	6
		20	pH 9.00	5652.9010	10 l	
		30	pH 8.91			
		40	pH 8.84	Contains 0.0	1% sodium benzoate.	
		50	pH 8.78	Contains bor	ric acid / potassium chlori	ide / sodium
		60	pH 8.73	hydroxide.		
		70	pH 8.69	This buffer is	s traceable to NIST Standa	ard Reference
		80	pH 8.66	Buffer.		
		90	pH 8.62			

Buffer Solution

pH 10.00 [Borate] / 'BAKER ANALYZED'

5655

NC CODE 3822 00 00	pH at 25°C	10.00 ± 0.05	PRODUCT	PACKING	
	Temperature (°C):		NO.		
	0	pH 10.30	5655.0500	500 ml	
	10	pH 10.17	5655.1000	11	
	20	pH 10.05	5655.9010	10 l Polycube	
	25	pH 10.00	5655.9020	20 l Polycube	
	30	pH 9.96			
	40	pH 9.89	Contains 0.0	1% sodium benzoate.	
	50	pH 9.82	Contains bo	ric acid / potassium chloride / sodium	1
	60	pH 9.76	hydroxide.		
	70	pH 9.72	This buffer is	s traceable to NIST Standard Reference	е
	80	pH 9.68	Buffer.		
	90	pH 9.64			

Buffer Solution

pH 10.00 [NH₄CI/EDTA] / 'BAKER ANALYZED'

7273

NC CODE 3822 00 00	pH at 20°C	10.0 - 10.6	PRODUCT	PACKING	CONT.
			NO.		вох
			7273.9020	20 l Polycube	
			Contains am	moniumchloride / magne	esium / potassium
			EDTA / sodi	umsulfide / ammonia.	
			This buffer is	traceable to NIST Stand	lard Reference
			Buffer.		



Buffe

Buffer Solution

pH 10.00 [Borate] BLUE / 'BAKER ANALYZED'

NC CODE 3822 00 00	pH at 25°C	10.00 ± 0.05
	Temperature (°C):	
	0	pH 10.30
	10	pH 10.17
	20	pH 10.05
	25	pH 10.00
	30	pH 9.96
	40	pH 9.89
	50	pH 9.82
	60	pH 9.76
	70	pH 9.72

80 90

pH at 20°C

PRODUCT	PACKING	CONT.
NO.		вох
5651.0500	500 ml	6
5651.1000	1 HDPE	6
5651.9010	10 l	

Contains 0.01% sodium benzoate.

Contains boric acid / potassium chloride / sodium

hydroxide

pH 9.68

pH 9.64

10-11

pH 10.68 pH 10.48 This buffer is traceable to NIST Standard Reference Buffer.

Buffer Solution

7268 pH 10-11 (NH₄CI/NH₄OH) / 'BAKER ANALYZED'

NC CODE 3822 00 00 **R:** 36/37/38

S: 26-36/37/39-45



 PRODUCT
 PACKING
 CONT.

 NO.
 BOX

 7268.1000
 1 l
 6

Contains ammoniumchloride / ammoniumhydroxide. This buffer is traceable to NIST Standard Reference Buffer.

Buffer Solution

pH 11.00 [Glycine] / 'BAKER ANALYZED'

NC CODE 3822 00 00 pH at 20°C 11.00 ± 0.05 Temperature (°C): 20 pH 11.00 30 pH 10.70 40 pH 10.42 50 pH 10.18 60 pH 9.92 70 pH 9.70 80 pH 9.50 90 pH 9.30

PRODUCT	PACKING	CONT.
NO.		вох
7267.1000	11	6

Contains glycine / sodium chloride / sodium hydroxide.
This buffer is traceable to NIST Standard Reference
Buffer.

Buffer Solution

pH 12.00 [Phosphate] / 'BAKER ANALYZED'

NC CODE 3822 00 00	pH at 20°C	12.00 ± 0.05
	Temperature (°C):	
	0	pH 12.58
	10	pH 12.26
	20	pH 12.00
	30	pH 11.75
	40	pH 11.53
	50	pH 11.31
	60	pH 11.09
	70	nH 10.88

80

	PACKING	CONT.
NO.		вох
7269.1000	11	
Contains 0.0)1% sodium benzoate.	
Contains so	dium hydroxide / disodium hydrog	en
phosphate.		
	. II . NICTO. I ID (
This buffer is	s traceable to NIST Standard Refer	ence

Innovation is principal to our business.



Buffer Solution

pH 13.00 [Glycine] / 'BAKER ANALYZED'

7270

NC CODE 3822 00 00	pH at 20°C	13.00 ± 0.05
	Temperature (°C):	
	0	pH 13.71
	10	pH 13.35
	20	pH 13.00
	30	pH 12.66
	40	pH 12.37
	50	pH 12.10
	60	pH 11.84
	70	pH 11.61
	80	pH 11.40
	90	pH 11.20

11	B02
11	
1% sodium benzoate.	
cine / sodium chloride / sod	dium hydroxide.
traceable to NIST Standard	d Reference

Buffer Solutions

See for detailed information section Reagents for pH Measurements, Titrimetry and Water Determination according to Karl Fischer, page 25



2-Butanone

See Methyl Ethyl Ketone



2-Butoxyethanol

8417

M = 118.18 g/mol

► CH₃(CH₂)₃OCH₂CH₂OH

1 l = 0.90 kg FLASHPOINT 61 °C

> **CAS NO.** 111-76-2 **EINECS** 203-905-0 NC CODE 2909 43 00

EC NO. 603 014 00 0

UN/ID NO. 2810 **ADR/RID** 6.1 T1

IMDG 6.1/III

R: 20/21/22-36/38

S: 36/37-46



Acidity (as CH ₃ COOH)	max. 0.01%
Boiling Range	169-173°C
Density (g/ml) at 20°C	0.900-0.904
Water (H ₂ O)	max. 0.2%

i koboci	1 ACKING	COIVI.
NO.		вох
8414.1000	11	
8414.2500	2.5	

2-(2-Butoxyethoxy)ethanol

► CH₃(CH₂)₃OCH₂CH₂OCH₂CH₂OH **M** = 162.23 g/mol

11 = 0.96 kg

FLASHPOINT 78 °C

CAS NO. 112-34-5

EINECS 203-961-6

NC CODE 2909 44 00

EC NO. 603 096 00 8

R: 36

S: 24-26



220-235°C	PRODUCT	PACKING	CONT.
0.950-0.956	NO.		вох
max. 0.5%	8417.1000	11	

www.jtbaker.com/europe

Boiling Range

Water (H₂O)

Density (g/ml) at 20°C



Butyl

n-Butyl Acetate

'BAKER

8018

M = 116.16 g/mol
1 = 0.88 kg

FLASHPOINT 24 °C

CAS NO. 123-86-4

EINECS 204-658-1

NC CODE 2915 33 00

EC NO. 607 025 00 1

UN/ID NO. 1123

ADR/RID 3 F1

IMDG 3/III

R: 10-66-67

Assay (by GC) min. 98% Acidity (as CH₃COOH) max. 0.005% Boiling Range 124-127°C Density (g/ml) at 20°C 0.878-0881 max. 0.5% n-Butanol Residue after Evaporation max. 0.005% Water (H₂O) max. 0.05% Trace Impurities (in ppm): Heavy Metals (as Pb) max. 0.5 Iron (Fe) max. 0.5

PRODUCT	PACKING	CONT.
NO.		вох
8018.1000	11	
8018.9025	25 l	
8018.9200	200	



n-Butyl Acetate MOS Grade

See for detailed information section Microelectronic materials page 32 or the website www.jtbaker.com/micro

8017

n-Butyl Alcohol

'BAKER ANALYZED' / ACS

► CH₃CH₂CH₂CH₂CH OH

M = 74.12 g/mol

1 I = 0.81 kg

FLASHPOINT 29 °C

CAS NO. 71-36-3

EINECS 200-751-6

NC CODE 2905 13 00

EC NO. 603 004 00 6

UN/ID NO. 1120

ADR/RID 3 F1 IMDG 3/III

R: 10-22-37/38-41-67 **S:** 13-26-37/39-46-7/9



Exceeds ACS Specifications	
Assay (by GC)	min. 99.4%
Aldehydes	passes test
Butyl Ether (by GC)	max. 0.2%
Carbonyl Compounds (as butyraldehyde)	max. 0.01%
Color (APHA)	max. 10
Residue after Evaporation	max. 0.001%
Titrable Acid (meq/g)	max. 0.0008
Water (H ₂ O)	max. 0.1%
Trace Impurities (in ppm):	
Aluminium (AI)	max. 0.5
Barium (Ba)	max. 0.1
Boron (B)	max. 0.02
Cadmium (Cd)	max. 0.05
Calcium (Ca)	max. 0.5
Chromium (Cr)	max. 0.02
Cobalt (Co)	max. 0.02
Copper (Cu)	max. 0.02
Iron (Fe)	max. 0.1
Lead (Pb)	max. 0.1
Magnesium (Mg)	max. 0.1
Manganese (Mn)	max. 0.02
Nickel (Ni)	max. 0.02
Tin (Sn)	max. 0.1
Zinc (Zn)	max. 0.1

PRODUCT	PACKING	CONT
NO.		ВОХ
8017.1000	11	ϵ
8017.2500	2.5	4
8017.5000	5 EcoTainer	
8017.9025	25 l	4
For safe han	dling of 25 l tin cans, see Se	If-closing tap.



n-Butylalcohol MOS Grade

See for detailed information section Microelectronic materials page 32 or the website www.jtbaker.com/micro

MSDS (Material Safety Data Sheets) are available in 16 languages at www.jtbaker.com/europe

8103

8019

sec-Butyl Alcohol 'BAKER ANALYZED'

► CH₃CH(OH)CH₂CH₃

M = 74.12 g/mol

1 l = 0.80 kg FLASHPOINT 24 °C

CAS NO. 78-92-2

EINECS 201-158-5

NC CODE 2905 14 90

EC NO. 603 004 01 3

UN/ID NO. 1120

ADR/RID 3F1

IMDG 3/III

R: 10-36/37-67

S: 13-24/25-26-46-7/9



Assay (by GC)	min. 99%
Boiling Range	98.5-100°C
Density (g/ml) at 25°C	0.801-0.803
Residue after Evaporation	max. 0.005%
Water (H ₂ O)	max. 0.2%

PRODUCT	PACKING	CONT.
NO.		вох
8103.1000	11	
8103.9025	25 l	
For safe han	dling of 25 l tin cans, see Self-closing tap	

tert-Butyl Alcohol

► (CH₃)₃COH

M = 74.12 g/mol

1 l = 0.77 kg

FLASHPOINT 11 °C

CAS NO. 75-65-0

EINECS 200-889-7 NC CODE 2905 14 10

EC NO. 603 005 00 1

UN/ID NO. 1120

ADR/RID 3 F1 IMDG 3/II

R: 11-20

S: 16-9





Assay (by GC)	min. 99%
Boiling Range	81.5-83°C.
Freezing Range	24.0-25.5°C.
Recorded Boiling Point	82.4°C
Residue after Evaporation	max. 0.006%
Trace Impurities (in ppm):	
Aluminium (AI)	max. 0.5
Barium (Ba)	max. 0.1
Boron (B)	max. 0.02
Cadmium (Cd)	max. 0.05
Calcium (Ca)	max. 0.5
Chromium (Cr)	max. 0.02
Cobalt (Co)	max. 0.02
Copper (Cu)	max. 0.02
Iron (Fe)	max. 0.1
Lead (Pb)	max. 0.1

max. 0.1

max. 0.02

max. 0.02

max. 0.1

max. 0.1

Magnesium (Mg)

Manganese (Mn)

Nickel (Ni)

Tin (Sn)

Zinc (Zn)

PRODUCT	PACKING	CONT.
NO.		вох
8019.1000	1 EcoTainer	6
8019.2500	2.5 EcoTainer	

tert-Butyl Alcohol 'BAKER'

M = 74.12 g/mol **1 l** = 0.77 kg

 $\textbf{FLASHPOINT} \quad 11~^{\circ}\text{C}$

► (CH₃)₃COH

CAS NO. 75-65-0

EINECS 200-889-7

NC CODE 2905 14 10

EC NO. 603 005 00 1

UN/ID NO. 1120

ADR/RID 3 F1

IMDG 3/II

R: 11-20

S: 16-9



81-83°C Boiling Range 23-25.5°C Freezing Point Residue after Evaporation max. 0.01%

PRODUCT PACKING CONT. вох NO. 8212.9025 25 l

For safe handling of 25 l tin cans, see Self-closing tap.

Butyl Chloride



8212

Butyl

tert-Butyl Methyl Ether

e Methyl tert-Butyl Eth



iso-Butyl Methyl Ketone



Butylparaben

eting the requirements of cGMP, USP/NF, BP/Ph.Eur, JP or FCC, page 36



Butyl Phosphate

See Tributyl Phosphate



Butyl Phthalate

See Dibutyl Phtalate

Cadmium 1000 µg/ml

5709 0.10% (w/v) / (Matrix: 2% nitric acid) / 'BAKER INSTRA-ANALYZED' / Plasma Standard

► Cd

M = 112.41 g/mol NC CODE 3822 00 00

R: 20/21/22-36/38

S: 26-36/37



Certificate Provided Reporting Actual Lot Analysis Cadmium (Cd) 998-1002 $\mu g/ml$

PRODUCT PACKING CONT. NO. BOX 5709.0100 100 ml

Prepared from the highest purity raw material available, generally greater than 99.999% spectral purity. The content of the solution is confirmed to be accurate to within \pm 0.2%. Typically 1000 μ g/ml. The certificate of analysis provided reports actual lot analysis. The certificate also lists the trace impurities.

Cadmium 1000 μg/ml

6924 (Matrix: 2% nitric acid) / 'BAKER INSTRA-ANALYZED' / Atomic Absorption Standard

► Cd

M = 112.41 g/mol

NC CODE 3822 00 00

R: 20/21/22-36/38 **S:** 26-36/37



Cadmium (Cd) 998-1002 μg/ml

NO. 6924.0100 100 ml 6924.0500 500 ml

PACKING

PRODUCT

Prepared by dissolution of high purity raw materials (min. 99.99% spectral purity). Assays are verified by ICP against standards traceable to NIST. Standard Reference Material numbers (SRM) are printed on each label.

CONT.

CONT.

BOX

BOX

Cadmium 1000 μg/ml

6807 'BAKER ANALYZED' / Atomic Absorption Standard

► Cd

M = 112.41 g/mol

NC CODE 3822 00 00

R: 20/21/22-36/38

S: 26-36/37



Cadmium (Cd)

998-1002 μg/ml

PRODUCT PACKING NO. 6807 0100 100 ml

6807.0500 500 ml Cadmium nitrate in nitric acid 0.5 mol/l.

Cadmium 10000 μg/ml

5723 (Matrix: 2% nitric acid) / 'BAKER INSTRA-ANALYZED' / Plasma Standard

► Cd

M = 112.41 g/mol NC CODE 3822 00 00

R: 20/21/22-36/38-52/53

S: 26-36/37



Certificate Provided Reporting Actual Lot Analysis

Cadmium (Cd) 9980-10020 $\mu g/ml$

PRODUCT PACKING CONT. NO. BOX 5723.0100 100 ml

Prepared from the highest purity raw material available, generally greater than 99.999% spectral purity. The content of the solution is confirmed to be accurate to within 0.2 %. Typically 10000 μ g/ml.

The certificate of analysis provided reports actual lot analysis. The certificate also lists the trace impurities.



Cadmium Chloride 2,5-Hydrate

'BAKER ANALYZED' / ACS 1208-01

► CdCl ₂ .2.5H ₂ O	
M =	228.34 g/mol
CAS NO. 7790-78-5	
FINECS	233-296-7

NC CODE 2827 39 80 **EC NO.** 48 008 00 3

UN/ID NO. 2570

ADR/RID 6.1 T5 IMDG 6.1/II

> R: 25-26-45-46-48/23/25-50/53-60-61

S: 45-53-60-61



environment





PRODUCT	PACKING	CONT.
NO.		вох
1208-01	500 g HDPE	

Cadmium Chloride Anhydrous

powder / 'BAKER ANALYZED' / AC

1039

1043

► CdCl₂ **M** = 183.31 g/mol **CAS NO.** 10108-64-2 **EINECS** 233-296-7 NC CODE 2827 39 80

EC NO. 48 008 00 3 **UN/ID NO.** 2570

ADR/RID 6.1 T5 IMDG 6.1/II

R: 25-26-45-46-48/23/25-50/53-60-61

S: 45-53-60-61



environment



Assay (argentometric titrn.)	min. 99.0%
Ammonium (NH ₄)	max. 0.01%
Calcium (Ca)	max. 0.01%
Copper (Cu)	max. 0.001%
Insoluble Matter	max. 0.01%
Iron (Fe)	max. 0.001%
Lead (Pb)	max. 0.005%
Nitrate and Nitrite (as NO ₃)	max. 0.003%
Potassium (K)	max. 0.02%
Sodium (Na)	max. 0.05%
Sulfate (SO ₄)	max. 0.01%
Zinc (Zn)	max. 0.05%

PRODUCT	PACKING	CONT.
NO.		вох
1039.0500	500 g	

Cadmium Sulfate Octahydrate

'BAKER ANALYZED' / ACS

M = 769.51 g/mol **CAS NO.** 7790-84-3 **EINECS** 233-331-6 NC CODE 2833 29 10 **EC NO.** 48 009 00 9 **UN/ID NO.** 2570 **ADR/RID** 6.1 T5

➤ 3CdSO₄.8H₂O

IMDG 6.1/III

R: 22-48/23/25-49-50/53

S: 45-53-60-61



dangerous for the environment



Meets ACS Specifications	
Assay	99.0-102.0%
Calcium (Ca)	max. 0.005%
Chloride (CI)	max. 0.001%
Copper (Cu)	max. 0.002%
Insoluble Matter	max. 0.005%
Iron (Fe)	max. 0.001%
Lead (Pb)	max. 0.001%
Nitrate and Nitrite (as NO ₃)	max. 0.003%
Potassium (K)	max. 0.01%
Sodium (Na)	max. 0.02%
Zinc (Zn)	max. 0.05%

PRODUCT	PACKING	CONT.
NO.		вох
1043.0100	100 g	
10.42 OOEO	EO ka	

Innovation is principal to our business.



Caffe

1514

Caffeine

► CH₃ŅCON(CH₃)COĆ:ÇN:CHNCH₃		
M =	194.19 g/mol	
CAS NO.	58-08-2	
EINECS	200-362-1	
NC CODE	2939 30 00	
EC NO.	613 086 00 5	
UN/ID NO.	1544	
ADR/RID	6.1 T2	
IMDG	6.1/III	

R: 22

Assay	98.5-101.5%
Acidity	passes test
Appearance of solution	passes test
Heavy Metals (as Pb)	max. 20 ppm
Identification	passes test
Loss on Drying	max. 0.5%
Sulfated Ash	max. 0.1%
Sulfates (as SO ₄)	max. 500 ppm
Related Substances:	
Chloroform (CHCl ₃)	max. 0.5%
Methanol (CH ₃ OH)	max. 0.5%

PRODUCT	PACKING	CONT.
NO.	ВОХ	
1514.0500	500 g	



Calamine

ee for qualities meeting the requirements of cGMP, USP/NF, BP/Ph.Eur, JP or FCC, page 36

Calcium 1000 μg/ml

(Matrix: 2% nitric acid) / 'BAKER INSTRA-ANALYZED' / Plasma Standard

► Ca **M** = 40.08 g/mol **NC CODE** 3822 00 00

R: 36/38



Certificate Provided Reporting Actual Lot AnalysisCalcium (Ca) 998-1002 μg/ml

 PRODUCT
 PACKING
 CONT.

 NO.
 BOX

 5710.0100
 100 ml

Prepared from the highest purity raw material available, generally greater than 99.999% spectral purity. The content of the solution is confirmed to be accurate to within \pm 0.2%. Typically 1000 $\mu g/ml$. The certificate of analysis provided reports actual lot analysis. The certificate also lists the trace impurities.

Calcium 1000 µg/ml

Matrix: 2% nitric acid) / 'BAKER INSTRA-ANALYZED' / Atomic Absorption Standard

Calcium (Ca)

M = 40.08 g/mol NC CODE 3822 00 00

6925

R: 36/38 **S:** 26

¥

Calcium (Ca) 998-1002 μg/ml

998-1002 μg/ml

 PRODUCT
 PACKING
 CONT.

 NO.
 BOX

 6925.0100
 100 ml

 6925.0500
 500 ml

Prepared by dissolution of high purity raw materials (min. 99.99% spectral purity). Assays are verified by ICP against standards traceable to NIST. Standard Reference Material numbers (SRM) are printed on each label.

Calcium 1000 µg/ml

'BAKER ANALYZED' / Atomic Absorption Standar

► Ca

M = 40.08 g/mol

NC CODE 3822 00 00

R: 36/38

S: 26.37

6808

Xi irritant
 PRODUCT NO.
 PACKING BOX

 6808.0100
 100 ml

 6808.0500
 500 ml

Calcium nitrate in nitric acid 0.5 mol/l.

www.jtbaker.com/europe



Calcium 10000 μ g/ml

(Matrix: 5% nitric acid) / 'BAKER INSTRA-ANALYZED' / Plasma Standard

5724

► Ca

M = 40.08 g/mol

NC CODE 3822 00 00 **UN/ID NO.** 2031

ADR/RID 8 CO2 IMDG 8/II

R: 34

S: 20-23-26-36/37/39-45



Certificate Provided Reporting Actual Lot Analysis Calcium (Ca)

9980-10020 μg/ml

PRODUCT	PACKING	CONT. BOX
NO.		
5724.0100	100 ml	

Prepared from the highest purity raw material available, generally greater than 99.999% spectral purity. The content of the solution is confirmed to be accurate to within 0.2 %. Typically 10000 μ g/ml.

The certificate of analysis provided reports actual lot analysis. The certificate also lists the trace impurities.

Calcium Acetate

See for qualities meeting the requirements of cGMP, USP/NF, BP/Ph.Eur, JP or FCC, page 36



0059

Calcium Acetate. 1/2 H₂0 'BAKER ANALYZED'

(CH₃COO)₂Ca.¹/₂ H₂O CAS NO. 62-54-4 **EINECS** 200-540-9 NC CODE 2915 29 00

Assay	98.0-101.0%
Alkalinity	passes test
Barium (Ba)	max. 0.005%
Chloride (CI)	max. 0.05%
Heavy Metals (as Pb)	max. 0.001%
Insoluble Matter	max. 0.005%
Iron (Fe)	max. 0.001%
Magnesium and Alkalies (as SO ₄)	max. 0.3%
Nitrate (NO ₃)	max. 0.003%
Sulfate (SO ₄)	max. 0.01%
Titrable Acid (meq/g)	max. 0.035

PRODUCT	PACKING	CONT
NO.		B0
0059.0500	500 g	

Calcium Acetate 1 mol/l

1N / 'BAKER ANALYZED

7117

CONT.

► (CH₂COO)₂Ca **M** = 158.17 g/mol **EINECS** 200-540-9 NC CODE 2915 29 00

Titer (mol/l) 0.995-1.005 **PRODUCT PACKING** 7117.9020 20 l Polycube

Volumetric Solution, ready for use.

Calibrate and standardise your analytical methods and equipment with J.T.Baker Volumetric and Buffer solutions.

Refer to the Analytical applications section of this catalogue for more details.



Calci

Calcium Carbonate ULTREX Ultrapure Reagent

4922

► CaCO₃ M = 100.09 g/mol CAS NO. 471-34-1 EINECS 207-439-9 NC CODE 2836 50 00

Certificate Provided Reporting Actual Lot Analysis Actual Lot Analysis of Lot. No. X02584

Actual Lot Analysis of Lot. No. No.	2001
Assay (dried basis)	100.00-100.00%
Alkalinity	passes test
Loss on Drying at 200°C	0.06%
Particulate Matter	0.2%
Metallic Impurities (in ppm):	
Aluminium (AI)	< 1
Barium (Ba)	< 10
Bismuth (Bi)	< 10
Cadmium (Cd)	< 5
Chromium (Cr)	< 1
Cobalt (Co)	< 1
Copper (Cu)	< 1
Iron (Fe)	< 1
Lead (Pb)	< 10
Lithium (Li)	< 100
Magnesium (Mg)	1
Manganese (Mn)	< 1
Mercury (Hg)	0.002
Molybdenum (Mo)	< 1
Nickel (Ni)	< 5
Niobium (Nb)	< 10
Potassium (K)	< 1
Silver (Ag)	< 1
Sodium (Na)	1
Strontium (Sr)	< 1
Tin (Sn)	< 10
Titanium (Ti)	< 1
Vanadium (V)	< 1
Zinc (Zn)	< 10
Zirconium (Zr)	< 1
Non Metallic Impurities (in ppm):	
Arsenic (As)	< 0.2
Fluoride (F)	4
Halide (as CI)	3
Nitrogen Compounds (as N)	1300
Di (DO)	

PRODUCT	PACKING	CONT.
NO.		вох
4922.0025	25 g Glass	
4922.0100	100 g Glass	

1464 P

Calcium Carbonate

Powder / 'BAKER INSTRA-ANALYZED'

► CaCO₃

M = 100.09 g/mol CAS NO. 471-34-1 EINECS 207-439-9 NC CODE 2836 50 00

Assay	min. 99.0%
Chloride (CI)	max. 0.001%
Fluoride (F)	max. 0.001%
Loss on Drying at 285°C	max. 1.0%
Sulfate (SO ₄)	max. 0.01%

Phosphate (PO₄)

Sulfur Compounds (as SO₄)

Silicon (Si)

Sulfate (SO ₄)	max. 0.01%
Trace Impurities (in ppm):	
Aluminium (AI)	max. 10
Cadmium (Cd)	max. 50
Chromium (Cr)	max. 5
Cobalt (Co)	max. 10
Copper (Cu)	max. 10
Iron (Fe)	max. 2
Lead (Pb)	max. 10
Magnesium (Mg)	max. 50
Manganese (Mn)	max. 10
Nickel (Ni)	max. 20
Potassium (K) (by FES)	max. 50
Silicon (Si)	max. 5
Sodium (Na) (by FES)	max. 80
Tin (Sn)	max. 5
Titanium (Ti)	max. 20
Vanadium (V)	max. 5

PRODUCT	PACKING	CONT.
NO.		вох
1464.0500	500 g	

Flux Grade.

0.5

< 1

2



Calcium Carbonate

Powder / 'BAKER ANALYZED' / low in Alkalies / ACS

0061

CaCO ₃	
	M =

M = 100.09 g/mol CAS NO. 471-34-1 EINECS 207-439-9 NC CODE 2836 50 00

Exceeds ACS Specifications	
Assay (dried basis)	min. 99.0%
Ammonium (NH ₄)	max. 0.003%
Average Particle Diameter, μ m (APD)	act. value reported
Barium (Ba)	max. 0.01%
Bulk Density (g/cc)(typical)	act. value reported
Chloride (CI)	max. 0.001%
Fluoride (F)	max. 0.0015%
Heavy Metals (as Pb)	max. 0.001%
Insoluble in Dilute HCl	max. 0.01%
Iron (Fe)	max. 0.002%
Magnesium (Mg)	max. 0.01%
Potassium (K)	max. 0.01%
Sodium (Na)	max. 0.01%
Specific Surface Area, m ² /g (typical)	act. value reported
Strontium (Sr)	max. 0.1%
Sulfate (SO ₄)	max. 0.005%
Mesh (Wet Screen Analysis):	
On U.S. No. 325 Sieve	act. value reported

PRODUCT	PACKING	CONT
NO.		BO)
0061.0500	500 g	

Calcium Carbonate

Precipitated / 'BAKER ANALYZED

0503

► CaCO₃

M = 100.09 g/mol **CAS NO.** 471-34-1

EINECS 207-439-9 **NC CODE** 2836 50 00

Assay	min. 99.0%
Barium (Ba)	max. 0.005%
Chloride (CI)	max. 0.005%
Insoluble in HCl	max. 0.01%
Iron (Fe)	max. 0.002%
Lead (Pb)	max. 0.0005%
Magnesium (Mg)	max. 0.05%
Nitrogen Compounds (as N)	max. 0.001%
Sulfate (SO ₄)	max. 0.01%

 PRODUCT
 PACKING
 CONT.

 NO.
 BOX

 0503.1000
 1 kg

Calcium Carbonate

See for qualities meeting the requirements of cGMP, USP/NF, BP/Ph.Eur, JP or FCC, page 36



Calcium Chloride Anhydrous

Granular < 6mm / 'BAKER'

M = 110.99 g/mol
CAS NO. 10043-52-4
EINECS 233-140-8
NC CODE 2827 20 00

R: 36

S: 22-24

ssay	min. 95%	PRODUCT	PACKING	CONT.
		NO.		BOX
		0070.1000	1 kg	6
		0070.5000	5 kg	

Calcium Chloride, Anhydrous

See for qualities meeting the requirements of cGMP, USP/NF, BP/Ph.Eur, JP or FCC, page 36



Find more Chromatography information at www.jtbaker.com/chromatography

Calci

Calcium Chloride Dihydrate 'BAKER ANALYZED'

0504

► CaCl₂.2H₂O

M = 147.02 g/mol **CAS NO.** 10035-04-8

EINECS 233-140-8

NC CODE 2827 20 00

EC NO. 17 013 00 2

R: 36

S: 22-24



Assay	min. 99.0%
Ammonium (NH ₄)	max. 0.005%
Barium (Ba)	max. 0.003%
Insoluble in Water	max. 0.005%
Magnesium (Mg)	max. 0.05%
Oxidizing Substances (as NO ₃)	max. 0.003%
pH of 5% Solution at 20°C	4.5-6.5
Potassium (K)	max. 0.01%
Sodium (Na)	max. 0.01%
Strontium (Sr)	max. 0.05%
Sulfate (SO ₄)	max. 0.005%
Trace Impurities (in ppm):	
Heavy Metals (as Pb)	max. 5

PRODUCT	PACKING	CONT.
NO.		вох
0504.0100	100 g	
0504.1000	1 kg	6
0504.9050	50 kg	

Appearance white to light cream colored

PACKING

Calcium Chloride Dihydrate

0064

► CaCl₂.2H₂O

M = 147.02 g/mol **CAS NO.** 10035-04-8

EINECS 233-140-8

NC CODE 2827 20 00

EC NO. 17 013 00 2

R: 36 **S:** 22-24

Assay	99.0-103.0%
Acid or alkaline reacting substances	passes test
Aluminium (AI)	passes test
Appearance of solution	passes test
Barium (Ba)	passes test
Heavy Metals (as Pb)	max. 0.001%
Identification	passes test
Iron (Fe)	max. 0.001%
Iron, Aluminium and Phosphate	passes test
Magnesium and Alkalies (as SO ₄)	max. 0.3%
Organic Volatile Impurities	passes test
pH of 5% Solution at 25°C	4.5-9.2
Sulfate (SO ₄)	max. 200 ppm

BOX 0 1 kg 6
) 1 kg 6
,
) 5 kg
) 50 kg

CONT.

CONT.

вох



Calcium Chloride, Dihydrate
See for qualities meeting the requirements of cGMP, USP/NF, BP/Ph.Eur, JP or FCC, page 36

Iron (Fe)

Calcium Chloride

7232 0.0125 mol/l / 'BAKER ANALYZED'

► CaCl₂

M = 110.99 g/mol **CAS NO.** 10043-52-4

EINECS 233-140-8 NC CODE 2827 20 00

0.0120-0.0130 Titer (mol/l)

PACKING PRODUCT

7232.9020 20 l Polycube

Volumetric Solution, ready for use.

Calcium Chloride

0.01 mol/l / 'BAKER ANALYZED'

► CaCl₂

7116

M = 110.99 g/mol

CAS NO. 10043-52-4

EINECS 233-140-8 NC CODE 2827 20 00 Titer (mol/l)

0.0095-0.0105

max. 5

PRODUCT

PRODUCT NO.

PACKING CONT. BOX

7116.1000 11

Volumetric Solution, ready for use.

Each lot of this product is standardized potentiometrically against Calcium Carbonate (NIST tracable reference standard).

The concentration of the solution is confirmed by titration with EDTA.



Calcium Citrate

See for qualities meeting the requirements of cGMP, USP/NF, BP/Ph.Eur, JP or FCC, page 36





► CaF₂

Assay (by EDTA titrn.) min. 95.0% **M** = 78.08 g/mol Chloride (CI) max. 0.05% **CAS NO.** 7789-75-5 Heavy Metals (as Pb) max. 0.005%

PRODUCT PACKING 1046.0500

max. 0.005%

max. 0.05%

Calcium Gluconate, Anhydrous

See for qualities meeting the requirements of cGMP, USP/NF, BP/Ph.Eur, JP or FCC, page 36

Iron (Fe)

Sulfate (SO₄)



1046

CONT.

вох

Calcium Gluconate, Monohydrate

See for qualities meeting the requirements of cGMP, USP/NF, BP/Ph.Eur, JP or FCC, page 36



Calcium Hydroxide 'BAKER ANALYZED'



EINECS 232-188-7

NC CODE 2826 19 00

► Ca(OH)₂ **M** = 74.09 g/mol **CAS NO.** 1305-62-0 **EINECS** 215-137-3 NC CODE 2825 90 19 **S:** 22-24-26-39

Assay	min. 95.0%
Chloride (CI)	max. 0.03%
Heavy Metals (as Pb)	max. 0.003%
Insoluble in HCI	max. 0.03%
Iron (Fe)	max. 0.05%
Magnesium and Alkali Salts (as SO ₄)	max. 1.0%
Sulfur Compounds (as SO ₄)	max. 0.1%

PRODUCT	PACKING	CONT.
NO.		вох
0073.1000	1 kg	

Calcium Hydroxide

0074

► Ca(OH) ₂	
M =	74.09 g/mol
CAS NO.	1305-62-0
EINECS	215-137-3
NC CODE	2825 90 19
R:	41
S:	22-24-26-39

Assay	95.0-100.5%
Acid Insoluble Substances	max. 0.5%
Carbonate (CO ₃)	passes test
Heavy Metals (as Pb)	max. 20 ppm
Identification	passes test
Magnesium and Alkali salts	max. 4.8%
Organic Volatile Impurities	passes test

PRODUCT	PACKING	CONT.
NO.		вох
0074.1000	1 kg	
0074.9025	25 kg	

Calcium Hydroxide

See for qualities meeting the requirements of cGMP, USP/NF, BP/Ph.Eur, JP or FCC, page 36



Calcium Lactate, 5-Hydrate

See for qualities meeting the requirements of cGMP, USP/NF, BP/Ph.Eur, JP or FCC, page 36



Calcium Lactate-Calcium Acetate solution



CONT.

NC CODE 3822 00 00

PRODUCT NO.	PACKING	
7233.9020	20 l Polycube	

Volumetric Solution, ready for use.

Innovation is principal to our business.

Calci

Calcium Nitrate Tetrahydrate 'BAKER ANALYZED' 0076

Ca(NO₃)₂.4H₂O **M** = 236.15 g/mol **CAS NO.** 13477-34-4 **EINECS** 233-332-1 NC CODE 2834 29 80 **UN/ID NO.** 1454 ADR/RID 5.1 02 IMDG 5.1/III **R:** 36-8

Assay (by EDTA titrn.) 99.0-103.0% Barium (Ba) max. 0.005% Chloride (CI) max. 0.005% Insoluble Matter and NH₄OH Precipitate max. 0.005% max. 0.2% Magnesium and Alkali Salts (as SO₄) pH of 5% Solution at 25°C 5.0-7.0 Sulfate (SO₄) max. 0.002% Trace Impurities (in ppm): Heavy Metals (as Pb) max. 5 Iron (Fe) max. 5

PRODUCT	PACKING	CONT.
NO.		вох
0076.0100	100 g	
0076.1000	1 kg	6

S: 15-17-24-7

Calcium Oxide

0078 **'BAKER ANALYZED**

► CaO **M** = 56.08 g/mol **CAS NO.** 1305-78-8 **EINECS** 215-138-9 NC CODE 2825 90 19 **R:** 41 **S:** 22-24-26-39

Chloride (CI) max. 0.005% Heavy Metals (as Pb) max. 0.005% Insoluble in CH₃COOH and NH₄OH max. 1.0% Precipitate Iron (Fe) max. 0.1% max. 5.0% Loss on Ignition Nitrate (NO₃) max. 0.05% Sulfate (SO₄) max. 0.1%

PRODUCT	PACKING	CONT.
NO.		вох
0078.1000	1 kg	6
0078.7100	100 lbs	
0078.9012	12 kg	

Calcium Oxide

0154

► CaO

M = 56.08 g/mol **CAS NO.** 1305-78-8 **EINECS** 215-138-9 NC CODE 2825 90 19

> **R:** 41 **S:** 22-24-26-29

Xi

Assay	min. 90%
Carbonate	passes test
Insoluble Matter in HCI	max. 1.0%
Loss on Ignition	max. 10%
Magnesium and Alkalies	passes test

PACKING	CONT.
	вох
1 kg	6
50 kg	
	1 kg

Calcium Oxide

eting the requirements of cGMP, USP/NF, BP/Ph.Eur, JP or FCC, page 36



Calcium Pantothenate

See for qualities meeting the requirements of cGMP, USP/NF, BP/Ph.Eur, JP or FCC, page 36



Calcium Phosphate, Dibasic, Anhydrous

See for qualities meeting the requirements of cGMP, USP/NF, BP/Ph.Eur, JP or FCC, page 36



Calcium Phosphate, Tribasic

e for qualities meeting the requirements of cGMP, USP/NF, BP/Ph.Eur, JP or FCC, page 36

1701

Calcium Sulfate Anhydrous 1-3mm / 'BAKER' / for drying

► CaSO₄ **M** = 136.14 g/mol **CAS NO.** 7778-18-9 **EINECS** 231-900-3 NC CODE 2833 29 90

Appearance passes test

PACKING CONT. **PRODUCT** BOX 1701.1000 1 kg

Calcium Sulfate Anhydrous 3-5mm / 'BAKER' / for drying

1702

► CaSO ₄	Appearance	passes test	PRODUCT	PACKING	CONT.
M = 136.14 g/mol			NO.		вох
CAS NO. 7778-18-9			1702.1000	1 kg	
EINECS 231-900-3			1702.7100	100 lbs	
NC CODE 3833 30 00					

Calcium Sulfate Anhydrous

	with indicator, 3-15mm / 'E					1703
► CaSO ₄		Appearance	passes test	PRODUCT	PACKING	CONT.
M =	136.14 g/mol			NO.		вох
CAS NO.	7778-18-9			1703.1000	1 kg	
EINECS	231-900-3			1703.9050	50 kg	
NC CODE	2833 29 90					
R:	42/43-49-51/53			Contains co	baltchloride.	
S:	20-22-24-53					
	N Quantity T toxic					
	for the					

Calcium Sulfate Dihydrate

CaSO ₄ .2H ₂ O	Meets ACS Specifications		PRODUCT	PACKING	CON
M = 172.17 g/mol	Assay	98.0-102.0%	NO.	1 Acking	BC
CAS NO. 10101-41-4	Carbonate (CO ₃)	passes test	0081.1000	1 kg	
EINECS 231-900-3	Chloride (CI)	max. 0.005%	0081.9050	50 kg	
NC CODE 2833 29 90	Heavy Metals (as Pb)	max. 0.002%			
	Insoluble in Dilute HCI	max. 0.02%			
	Iron (Fe)	max. 0.001%			
	Magnesium (Mg)	max. 0.02%			
	Nitrate (NO ₃)	passes test			
	Potassium (K)	max. 0.005%			
	Sodium (Na)	max. 0.02%			
	Strontium (Sr)	max. 0.05%			

Calcium Sulfate Hemihydrate Powder / 'BAKER TLC' / for Thin Layer Chrom

0536

► CaSO ₄ .0,5H ₂ O M = 145.15 g/mol	Chloride (CI) Iron (Fe)	max. 0.01% max. 0.01%	PRODUCT	PACKING	CONT. BOX
CAS NO. 7778-18-9 EINECS 231-900-3 NC CODE 2833 29 90	Suitability for TLC	passes test	0536.0500	500 g	

Capping A (ABI)
'BAKER ANALYZED' / for DNA/RNA synthesis

9497

CONT.

11=	0.91 kg
FLASHPOINT	−21 °C
NC CODE	3822 00 00
UN/ID NO.	2924
ADR/RID	3 FC
IMDG	3/II
R:	11-19-20/21/22-37/38-41
S:	26-28-36/37/39-43A-45

	8	F
h	ighly	
fl	amma	ble

Suitable for Oligonucleotide Synthes	is
Assay Acetic Anhydride	9.0-11.0% (v/v)
Assay Pyridine	9.0-11.0% (v/v)
Appearance: clear, colorless solution free	
from visible particulates	passes test

NO.		BOX
9497.0180	180 ml Glass	
9497.0450	450 ml Glass	
9497.2000	2 l Glass	
9497.2500	2.5 l Glass	
9497.4000GL	4 Glass	

PACKING

PRODUCT



Cappi

Capping A (ABI)

9510 'BAKER ANALYZED' / for DNA/RNA synthesis

11 = 0.91 kg

FLASHPOINT −21 °C

NC CODE 3822 00 00

UN/ID NO. 2924

ADR/RID 3 FC

IMDG 3/II

R: 11-19-37/38-41 **S:** 26-36/37/39-43A-45



 Suitable for Oligonucleotide Synthesis

 Assay Acetic Anhydride
 8.2-10.0% (v/v)

 Appearance: clear, colorless solution free from visible particulates
 passes test

 PRODUCT
 PACKING
 CONT.

 NO.
 BOX

 9510.0200
 200 ml Glass

 9510.0450
 450 ml Glass

 9510.0900
 900 ml Glass

 9510.2500
 2.5 l Glass

 9510.4000GL
 4 l Glass

(

Capping B (ABI)

'BAKER ANALYZED' / for DNA/RNA synthesis

1 I = 0.91 kg

FLASHPOINT −21 °C

NC CODE 3822 00 00

UN/ID NO. 2924

ADR/RID 3 FC

IMDG 3/II

R: 11-19-34

S: 16-20-26-9





 Suitable for Oligonucleotide Synthesis

 Assay (1-Methylimidazole)
 15.2-16.8% (v/v)

 Appearance: clear, colorless solution free from visible particulates
 passes test

 Water (H₂O)
 max. 120 ppm

PRODUCT	PACKING	CONT.
NO.		вох
9485.0180	180 ml Glass	
9485.0450	450 ml Glass	
9485.2000	2 l Glass	
9485.2500	2.5 Glass	
9485.4000GL	4 Glass	

Capping B (ABI)

'BAKER ANALYZED' / for DNA/RNA synthesis

1 I = 0.91 kg
FLASHPOINT -21 °C
NC CODE 3822 00 00
UN/ID NO. 2924
ADR/RID 3 FC
IMDG 3/II
R: 11-19-20/21/22-34

9512

S: 26-36/37/39-43A-45





PRODUCT NO.	PACKING	CONT. BOX
9512.0200	200 ml Glass	
9512.0450	450 ml Glass	
9512.0900	900 ml Glass	
9512.2500	2.5 Glass	
9512.4000GL	4 Glass	



Capping A/B reagents for use in DNA synthesis

See for detailed information section Reagents for DNA/RNA Synthesis, page 261



Capryl Alcohol

See 2-Octanol



CAPS

'BAKER ULTRAPURE BIOREAGENT'

► C₉H₁₉NO₃S **M** = 221.32 g/mol **CAS NO.** 1135-40-6 **EINECS** 214-492-1 **NC CODE** 2922 50 00

min. 98.0% Assav Appearance passes test DNase Activity none detected Heavy Metals (as Pb) max. 5 ppm Insoluble Matter max. 0.5% Protease Activity none detected RNase Activity none detected Product Information (not specifications): pKa at 20°C 10.40

PRODUCT	PACKING	CONT.
NO.		вох
4118.0025	25 g Glass	
4118.0100	100 g	
4118.1000	1 kg	



Carbamide

See Urea

Carbon Tetrachloride

See Tetrachloromethane

Castor Oil

See for qualities meeting the requirements of cGMP, USP/NF, BP/Ph.Eur, JP or FCC, page 36



Caustic Potash

See Potassium Hydroxide



CCV Solution I & II

'BAKER INSTRA-ANALYZED' / Plasma Standard

6106-01

CONT.

NC CODE	3822 00 00
UN/ID NO.	3264
ADR/RID	8 C1
IMDG	8/III
R:	36/38
S:	26-37
	A A Xi

Kit contains one bottle of each	solution
Element Concentrations of Solu	ıtion l (μg/ml):
Barium (Ba)	500
Beryllium (Be)	200
Cadmium (Cd)	250
Cobalt (Co)	500
Copper (Cu)	500
Iron (Fe)	500
Lead (Pb)	500
Manganese (Mn)	500
Nickel (Ni)	500
Silver (Ag)	100
Thallium (TI)	500
Zinc (Zn)	500
Element Concentrations of Solu	ıtion II (μg/ml):
Aluminium (AI)	500
Antimony (Sb)	500
Arsenic (As)	500
Calcium (Ca)	5000
Chromium (Cr)	500
Magnesium (Mg)	5000
Potassium (K)	5000
Selenium (Se)	500
Sodium (Na)	5000
Vanadium (V)	500

NO.		BO
6106-01	100 ml x 2	
For use in F	PA Contract Laboratory Progra	am (CLP)
		iiii (CLI).
Traceable to	NIST.	

PRODUCT

PACKING

CDTA, Sodium Salt

e-N,N,N,N-tetra- acetic acid, Sod



Cedarwood Oil

7397

1 = 0.98 kg	Refractive Index n ²⁰	1.501-1.511	PRODUCT	PACKING	CONT.
NC CODE 3302 90 90			NO.		вох
			7397.0025	25 ml	6
			For Immersion).	

	Celite 503 'BAKER' / analytical filter aid					0902
EINECS	272-489-0	Appearance	passes test	PRODUCT	PACKING	CONT.
NC CODE	3822 00 00			NO.		вох

0902.1000

1 kg

R: 68/20 **S**: 22



119

Celit

Celite 545 3376 'BAKER ANALYZED

CAS NO. 68855-54-9 **EINECS** 272-489-0

R: 68/20

S: 22



Suitable for use in pesticide residue analysis after extraction with petroleum ether

Loss on Drying at 110°C max. 0.2% Petroleum Ether Extractables value on label pH of 5% Slurry at 20°C 8.0-10.0

PRODUCT	PACKING	CONT.
NO.		вох
3376.2500	2.5 kg	

Cellulose

/ Flexible TLC Sheets, 20 X 20 cm 'BAKER-FLEX'

NC CODE 3912 11 00

5007

PRODUCT	PACKING	CONT.
NO.		вох
5007	25 sheets	

A flexible sheet coated with high purity cellulose powder.

5036

Cellulose

'BAKER TLC' / for Thin Layer Chromatography

NC CODE 3912 11 00 Suitability for TLC Physical Data (not specifications): Average Particle Diameter, μ m (APD) 2-20 Bulk Density (g/cc) (untapped)

PRODUCT PACKING CONT. BOX NO. 5036.0500 500 g

Cellulose DEAE

'BAKER-FLEX' / Flexible TLC Sheets, 20 X 20 cm

NC CODE 3912 39 10

PRODUCT	PACKING	CONT.
NO.		вох
5016	25 sheets	

A flexible sheet coated with (2-(diethylamino)ethyl) cellulose powder.



Cellulose 2-(Diethylamino)ethyl Ether See Cellulose DEAE

5012

Cellulose PEI

'BAKER-FLEX' / Flexible TLC Sheets, 20 X 20 cm

NC CODE 3912 39 80

PRODUCT	PACKING	CONT.
NO.		вох
5012	25 sheets	

A flexible sheet coated with polyetylenimine cellulose



Cerium(IV) Ammonium Nitrate

See Ammonium Cerium(IV) Nitrate

2008

Cerium(IV) Sulfate 'BAKER ANALYZED'

► Ce(SO ₄) ₂ .4H ₂ O			
M =	404.30 g/mol		
CAS NO.	10294-42-5		
EINECS	237-029-5		
NC CODE	2846 10 00		

Assay	min. 98%
Chloride (CI)	max. 0.001%
Heavy Metals (as Pb)	max. 0.005%
Insoluble in H ₂ SO ₄	max. 0.005%
Iron (Fe)	max. 0.005%
Phosphate (PO ₄)	max. 0.01%

PRODUCT	PACKING	CONT.
NO.		вох
2008.1000PE	1 kg HDPE	6

7114

3869

CONT.

CONT.

вох



Cerium(IV) Sulfate 0.1 mol/l

0.1 mol/l / 'BAKER ANALYZED'

► Ce(SO₄)₂.4H₂O **M** = 404.30 g/mol **CAS NO.** 10294-42-5 **EINECS** 237-029-5

NC CODE 2846 10 00

Titer (mol/l)

0.0997-0.1003

PRODUCT PACKING NO.

7114.1000 7114.9020 20 l Polycube

PRODUCT

3869.1200

NO.

Volumetric Solution, ready for use.

PACKING

125 ml HDPE x 12

Cervix Fixative

NC CODE 3822 00 00

UN/ID NO. 1993 ADR/RID 3 F1

IMDG 3/II

R: 11-36

S: 16-24/25-26-7/9





Cervix Spray for Cytology

r.com and select Clinical



Cesium Chloride

BAKER ULTRAPURE BIOREAGENT

000.	M =	168.37 g/mol

CAS NO. 7647-17-8 EINECS 231-600-2 NC CODE 2827 39 80

For Density Gradient Centrifugation	
Assay	min. 99.9%
Absorbance of 50% (w/v) Aqueous	
Solution at 260 nm (1-cm path)	max. 0.02
DNase Activity	none detected
Protease Activity	none detected
RNase Activity	none detected
Trace Impurities (in ppm):	

FIDLEASE ACTIVITY	none detected
RNase Activity	none detected
Trace Impurities (in ppm):	
Aluminium (AI)	max. 10
Barium (Ba)	max. 5
Calcium (Ca)	max. 10
Chromium (Cr)	max. 5
Copper (Cu)	max. 3
Iron (Fe)	max. 2
Lead (Pb)	max. 5
Lithium (Li)	max. 5
Magnesium (Mg)	max. 3
Manganese (Mn)	max. 5
Potassium (K)	performance only
Rubidium (Rb)	max. 10

performance only

max. 5

max. 5

PRODUCT	PACKING	CONT
NO.		ВОХ
4042.0100	100 g	
4042.1000	1 kg	

Analytical applications are available in our Technical Library at www.jtbaker.com/europe

Sodium (Na)

Strontium (Sr)

Sulfate (SO₄)



► CsCl

Cesiu

Cesium Chloride 1955 'BAKER INSTRA-ANALYZED'

M = 168.37 g/mol **CAS NO.** 7647-17-8 **EINECS** 231-600-2 NC CODE 2827 39 80

min. 99.5% Assay Absorbance of 50% (w/v) Aqueous Solution at 260 nm (1-cm path) max. 0.02 Trace Impurities (in ppm): Aluminium (AI) max. 0.1 Barium (Ba) max. 5 Calcium (Ca) max. 0.3 Chromium (Cr) max. 0.1 Iron (Fe) max. 0.2 Lithium (Li) max. 0.5 Magnesium (Mg) max. 0.1 Manganese (Mn) max. 0.1 Phosphorus (as P₂O₅) max. 0.5 Potassium (K) max. 1 Rubidium (Rb) max. 5 Silicon Dioxide (SiO₂) max. 1 Sodium (Na) max. 0.2 Strontium (Sr) max. 0.2 Sulfate (SO₄) max. 5

PRODUCT	PACKING	CONT.
NO.		вох
1955.0025	25 g Glass	
1955.0250	250 g	

2012

Cesium Chloride


BAKER ANALYZED

► CsCl

M = 168.37 g/mol **CAS NO.** 7647-17-8 **EINECS** 231-600-2

NC CODE 2827 39 80

Assay	min. 99.5%
Nitrogen Compounds (as N)	max. 0.001%
Trace Impurities (in ppm):	
Aluminium (AI)	max. 5
Heavy Metals (as Pb)	max. 5
Iron (Fe)	max. 3
Magnesium (Mg)	max. 5
Potassium (K)	max. 20
Sodium (Na)	max. 20
Sulfate (SO ₄)	max. 50





Cetane



cGMP Manufactured Products

e for detailed information section Biopharmaceutical Products, page 36

4145

CHAPS

'BAKER ULTRAPURE BIOREAGENT'

C₃₂H₅₈N₂O₇S

M = 614.89 g/mol **CAS NO.** 75621-03-3 NC CODE 3822 00 00

Assay (by TLC) min. 98% Conductivity, µmhos max. 100 Insoluble Matter max. 0.2% pH of 10% Slurry at 25°C 5.0-7.0 Solubility passes test

PRODUCT PACKING CONT. NO. BOX 4145.0005 5 g 4145.0025 25 g Glass 100 g 4145.0100 4145.1000 1 kg



Charcoal, Activated

See for qualities meeting the requirements of cGMP, USP/NF, BP/Ph.Eur, JP or FCC, page 36

4147

CHES, Sodium salt

'BAKER ULTRAPURE BIOREAGENT'

► C₈H₁₆NNaO₃S **M** = 229.27 g/mol **CAS NO.** 3076-05-9

> **EINECS** 221-352-3 NC CODE 2921 30 99

Assay (dried basis) min. 99.0% Appearance passes test DNase Activity none detected Heavy Metals (as Pb) max. 5 ppm Insoluble Matter max. 0.005% Protease Activity none detected RNase Activity none detected Water (by KF, volumetric) max. 8.0%

9.30

Product Information (not specifications): pKa at 20°C

PRODUCI	PACKING	CONT.
NO.		вох
4147.0025	25 g Glass	
4147.0100	100 g	
4147.1000	1 kg	

Chloral Hydrate

'BAKER'

min. 24.5%

max. 2.0%

PRODUCT

1255.0250

NO.



M = 165.40 g/mol

CAS NO. 302-17-0 **EINECS** 206-117-5

NC CODE 2905 59 99

EC NO. 605 014 00 6

UN/ID NO. 2811

ADR/RID 6.1 T2

IMDG 6.1/II

R: 25-36/38

S: 25-45

Tayin T

Assay	98.5 - 101.0%
Appearance of solution	passes test
Chloral-alcoholate	passes test
Chlorides (as CI)	max. 100 ppm
Heavy Metals (as Pb)	max. 20 ppm
Identification	passes test
Nonvolatile Residue	max. 0.1%
pH (10%; water)	3.5 - 5.5

PRODUCT	PACKING	CONT.	
NO.		вох	
1052.0500	500 g		

Stored in an airtight container.

PACKING

250 g

Chloramine T

'BAKER ANALYZED

► CH₃C₆H₄SO₂NCINa.3H₂O

M = 281.69 g/mol

CAS NO. 7080-50-4

EINECS 204-854-7

NC CODE 2935 00 90

EC NO. 616 010 00 9

UN/ID NO. 3263 **ADR/RID** 8 C8

IMDG 8/III

R: 22-31-34-42

S: 22-26-36/37/39-45-7



$\hbox{N-Chloro-} p-toluene sulfonamide, N-Sodium\ Derivate$

Assay (as Chlorine)

Insoluble in Ethanol

See Chloramine T



1255

CONT.

вох

Chloroacetic Acid 'BAKER ANALYZED'

M = 94.50 g/mol

FLASHPOINT 126 °C

► CICH₂COOH

CAS NO. 79-11-8

EINECS 201-178-4

NC CODE 2915 40 00

EC NO. 607 003 00 1

UN/ID NO. 1751

ADR/RID 6.1 TC2

IMDG 6.1/II

R: 25-34-50

S: 23-37-45-61



dangerous for the



Assay min. 99.0% max. 0.001% Chloride (CI) Heavy Metals (as Pb) max. 0.001% Insoluble Matter max. 0.01% Iron (Fe) max. 0.002% 62-64°C Melting Point Residue after Ignition max. 0.02% Substances Darkened by H₂SO₄ passes test Sulfate (SO_4) max. 0.02%

PRODUCT	PACKING	CON
NO.		BO
0177.1000	1 kg	

Chloroauric Acid

See Tetrachloroauric(III) Acid



Chlor

Chlorobenzene 8023 'BAKER ANALYZED

M = 112.56 g/mol

1 l = 1.10 kg

 $\textbf{FLASHPOINT} \ \ 28\,^{\circ}\text{C}$

CAS NO. 108-90-7

EINECS 203-628-5

NC CODE 2903 61 00

EC NO. 602 033 00 1

UN/ID NO. 1134

ADR/RID 3F1

IMDG 3/III

R: 10-20-51/53

S: 24/25-61





Assay (by GC) min. 99% Boiling Range max. 1.0°C Density (g/ml) at 25°C 1.101-1.104 Recorded Boiling Point 132°C

PRODUCT	PACKING	CONT.
NO.		вох
8023.1000	11	6
8023.2500	2.5	4
8023.9200	200 l	

Chloroform

'BAKER ULTRA RESI-ANALYZED' / for Organic Residue Analysis 9257

► CHCl₃

M = 119.38 g/mol

1 l = 1.47 kg

CAS NO. 67-66-3

EINECS 200-663-8

NC CODE 2903 13 00 EC NO. 602 006 00 4

UN/ID NO. 1888

ADR/RID 6.1 T1

IMDG 6.1/III

R: 22-38-40-48/20/22

S: 36/37



Assay (by GC) (corrected for water) (exclusive of preservative) min. 99.8% Chloride (CI) (in ppm) max. 10 Color (APHA) max. 10 Residue after Evaporation (in ppm) max. 2 Substances Darkened by H₂SO₄ passes test Titrable Acid (meq/g) max. 0.0005 Water (H₂O) max. 0.05% ECD Sensitive Impurities (as Heptachlor Epoxide):

Single Impurities (pg/ml) max. 10

FID-Sensitive Impurities (as 2-Octanol):

Assay (by GC) (exclusive of preservative)

Single Impurities (ng/ml) max. 10

min. 99.8%

PRODUCT NO.	PACKING	CONT. BOX
9257.1000	11	6
9257.4000	4 Glass	4

Stabilized with about 0.75% ethanol.

Chloroform

'BAKER HPLC ANALYZED' / for use in High Performance Liquid Chromatography

► CHCl₃

9174

M = 119.38 g/mol

1 l = 1.47 kg

CAS NO. 67-66-3

EINECS 200-663-8

NC CODE 2903 13 00

EC NO. 602 006 00 4

UN/ID NO. 1888

ADR/RID 6.1 T1 IMDG 6.1/III

R: 22-38-40-48/20/22

S: 36/37



Preservative (as Amylene) max. 50 ppm Water (H₂O) max. 0.01% Fluorescence Trace Impurities (as quinine base), ppb: Measured at 450 nm max. 1.0 Measured at Emission Maximum for max. 2.5 Solvent Impurities Physical Data (not specifications): Density (g/ml) at 20°C 1.483 R

Trace Impurities:	
Chloride (CI) (in ppm)	max. 10
Residue after Evaporation (in ppm)	max. 5
itrable Acid (meq/g)	max. 0.0005
Iltraviolet Absorbance (1.00-cm path vs water):	
at 254 nm	max. 0.15
at 280 nm	max. 0.05
at 350 nm	max. 0.01
JV Cut-off, nm	max. 245

PRODUCT I	PACKING	CONT.
NO.		вох
9174.1000	11	6
9174.2500	2.5	4

Stabilized with about 50 ppm Amylene. Filtered through a 0.2 micron filter. Packaged under Nitrogen.



Chloroform

'BAKER HPLC ANALYZED' / for use in High Performance Liquid Chromatography

9175

\triangleright	CHCI ₃
	Of FOIg

M = 119.38 g/mol **1 l** = 1.47 kg **CAS NO.** 67-66-3

EINECS 200-663-8 NC CODE 2903 13 00 EC NO. 602 006 00 4

UN/ID NO. 1888

ADR/RID 6.1 T1 IMDG 6.1/III

R: 22-38-40-48/20/22

S: 36/37



Assay (by GC) (exclusive of preservative)	min. 99.8%
Preservative (C ₂ H ₅ OH) (by GC)	0.5-1.0% (w/w)
Residue after Evaporation	max. 2 ppm
Substances Darkened by H ₂ SO ₄	passes test
Titrable Acid (meq/g)	max. 0.0005
Water (H ₂ O)	max. 0.01%

Fluorescence Trace Impurities (as quinine base), ppb:

Measured at 450 nm max. 0.25 Measured at Emission Maximum for Solvent Impurities max. 1.0 Physical Data (not specifications):

1.478

Density (g/ml) at 20°C Trace Impurities (in ppm):

Chloride (CI) max. 10

Ultraviolet Absorbance (1.00-cm path vs water):		
at 254 nm	max. 0.15	
at 280 nm	max. 0.0	
at 350 nm	max. 0.0	
LIV Cut-off nm	max 24 ^p	

PRODUCT NO.	PACKING	CONT. BOX
9175.2500	2.5	4

Stabilized with about 0.75% ethanol. Filtered through a 0.2 micron filter. Packaged under Nitrogen.

Chloroform

'BAKER ANALYZED' / for Dithizon / ACS



M = 119.38 g/mol

1 l = 1.47 kg

CAS NO. 67-66-3 EINECS 200-663-8

NC CODE 2903 13 00

EC NO. 602 006 00 4

UN/ID NO. 1888

ADR/RID 6.1 T1

IMDG 6.1/III

R: 22-38-40-48/20/22

S: 36/37



Exceeds ACS Specifications	
Assay (by GC)	min. 99.8%
Acetone and Aldehyde (as (CH ₃) ₂ CO)	max. 0.005%
Acid and Chloride	passes test
Color (APHA)	max. 10
Free Chlorine	passes test
Preservative (C ₂ H ₅ OH)	about 0.75%
Residue after Evaporation	max. 0.001%
Substances Darkened by H ₂ SO ₄	passes test
Suitability for Use in Dithizone Test	passes test
Trace Impurities (in ppm):	
Aluminium (AI)	max. 0.5
Barium (Ba)	max. 0.1
Boron (B)	max. 0.02
Cadmium (Cd)	max. 0.05

Aluminium (AI)	max. 0.5
Barium (Ba)	max. 0.1
Boron (B)	max. 0.02
Cadmium (Cd)	max. 0.05
Calcium (Ca)	max. 0.5
Chromium (Cr)	max. 0.02
Cobalt (Co)	max. 0.02
Copper (Cu)	max. 0.02
Iron (Fe)	max. 0.1
Lead (Pb)	max. 0.05
Magnesium (Mg)	max. 0.1
Manganese (Mn)	max. 0.02
Nickel (Ni)	max. 0.02
Tin (Sn)	max. 0.1
Zinc (Zn)	max. 0.1

PRODUCT	PACKING	CONT.
NO.		вох
7018.1000	11	6
7018.1000S	1 EcoTainer	
7018.2500	2.5	4
7018.5000	5 EcoTainer	
7018.9200	200	
EcoTainer, the	e metal solvent can for more s	afety in the lab.

Stabilized with about 0.75% ethanol.

Solv-IT center - A new solution to chemistry related technical questions.

Visit us at www.jtbaker.com/solvit



Chlor

Chloroform

7019 'BAKER ANALYZED' / Ultraviolet Spectrophotometry / ACS

at 260 nm

at 270 nm

at 290-400 nm

► CHCl₃

M = 119.38 g/mol **1 l** = 1.47 kg

CAS NO. 67-66-3 **EINECS** 200-663-8

NC CODE 2903 13 00

EC NO. 602 006 00 4

UN/ID NO. 1888 **ADR/RID** 6.1 T1

IMDG 6.1/III

R: 22-38-40-48/20/22

S: 36/37



Exceeds ACS Specifications	
Assay (by GC)	min. 99.8%
Acetone and Aldehyde (as (CH ₃) ₂ CO)	max. 0.001%
Acid and Chloride	passes test
Color (APHA)	max. 10
Free Chlorine	passes test
Lead (Pb)	max. 0.05 ppm
Preservative (C ₂ H ₅ OH)	about 0.75%
Residue after Evaporation	max. 5 ppm
Substances Darkened by H ₂ SO ₄	passes test
Suitability for Use in Dithizone Test	passes test
Ultraviolet Absorbance (1.00-cm pati	n vs water):
at 245 nm	max. 1.00
at 255 nm	max. 0.25

max. 0.15

max. 0.05

max. 0.01

PRODUCT	PACKING	CONT.
NO.		вох
7019.0500	500 ml	
7019.1000	11	6
7019.2500	2.5	4

Stabilized with about 0.75% ethanol.

Chloroform

7386 'BAKER ANALYZED' / for Dithizon

► CHCl₃

M = 119.38 g/mol

1 l = 1.47 kg

CAS NO. 67-66-3

EINECS 200-663-8 NC CODE 2903 13 00

EC NO. 602 006 00 4

UN/ID NO. 1888

ADR/RID 6.1 T1

IMDG 6.1/III **R:** 22-38-40-48/20/22

S: 36/37



Assay (by GC)	min. 99%
Acetone and Aldehyde (as (CH ₃) ₂ CO)	max. 0.005%
Acid and Chloride	passes test
Color (APHA)	max. 10
Density (g/ml) at 25°C	1.467-1.474
Free Chlorine	passes test
Preservative (C ₂ H ₅ OH)	about 0.75%
Residue after Evaporation	max. 0.001%
Substances Darkened by H ₂ SO ₄	passes test
Trace Impurities (in ppm):	
Aluminium (AI)	max. 0.5
Barium (Ba)	max. 0.1
Boron (B)	max. 0.02
Cadmium (Cd)	max. 0.05
Calcium (Ca)	max. 0.5
Chromium (Cr)	max. 0.02
Cobalt (Co)	max. 0.02
Copper (Cu)	max. 0.02
Iron (Fe)	max. 0.1
Lead (Pb)	max. 0.1
Magnesium (Mg)	max. 0.1
Manganese (Mn)	max. 0.02
Nickel (Ni)	max. 0.02
Tin (Sn)	max. 0.1
Zinc (Zn)	max. 0.1

PRODUCT	PACKING	CONT.
NO.		вох
7386.1000	11	6
7386.2500	2.5	4
7386.5000	5 EcoTainer	
7386.9025	25 l	
7386.9200	200 l	
EcoTainer, th	e metal solvent can for mo	ore safety in the lab.
For safe han	dling of 25 l tin cans, see S	self-closing tap.

Stabilized with about 0.75% ethanol.

Chloroform

7152

► CHCl₃ **M =** 119.38 g/mol

1 l = 1.47 kg

CAS NO. 67-66-3

EINECS 200-663-8

NC CODE 2903 13 00

EC NO. 602 006 00 4

UN/ID NO. 1888

ADR/RID 6.1 T1

IMDG 6.1/III

R: 22-38-40-48/20/22

S: 36/37



Assay (by GC)	min. 99%
Appearance	passes test
Color (APHA)	max. 10
Ethanol (C ₂ H ₅ OH)	0.4-1.0% (m/m)
Residue after Evaporation	max. 0.002%

PRODUCI	PACKING	CON1.
NO.		вох
7152.1000	11	6
7152.2500	2.5	4
7152.9025	25 l	
7152.9200	200 l	
For safe han	dling of 25 l tin cans, see S	elf-closing tap.

Stabilized with about 0.75% ethanol.

7413

R: 22-38-40-48/20/22

S: 36/37

Chloroform-d₁ min. 99.8 Atom% D / 'BAKER'

► CDCl₃ PRODUCT PACKING CONT. min. 99.8% Assay **M** = 120.38 g/mol NO. вох **CAS NO.** 865-49-6 7413.0100 100 g **EINECS** 212-742-4 NC CODE 2903 13 00 **EC NO.** 602 006 00 4

Chloroplatinic Acid Hexahydrate

2890-03

► H ₂ PtCl ₆ .6H ₂ O	
M =	517.92 g/mol
CAS NO.	18497-13-7
EINECS	241-010-7
NC CODE	2843 90 90
EC NO.	78 009 00 4
UN/ID NO.	2507
ADR/RID	8 C2
IMDG	8/111
R:	25-34-42/43

UN/ID NO. 1888 **ADR/RID** 6.1 T1 IMDG 6.1/III

Meets ACS Specifications. Meets Reagent Specifications for testing USP/NF monographs Assay (Pt) min. 37.50% Alkalies and other Salts (as SO₄) max. 0.05% Solubility in Alcohol passes test Suitability for Potassium Determination passes test

PRODUCT	PACKING	CONT.
NO.		вох
2890-03	1 g	

Keep in a cool place and tightly closed.

3-(3-cholamidopropyl)dimethylammonio-1-propanesulfonate

S: 22-26-36/37/39-45



Chromic Acid Anhydride



Chromium 1000 µg/ml

0.10% (w/v) / (Matrix: 2% nitric acid) / 'BAKER INSTRA-ANALYZED' / Plasma Standard



Certificate Provided Reportin	g Actual Lot Analysis
Chromium (Cr)	998-1002 μg/ml

PRODUCT	PACKING	CONT
NO.		B0
5711.0100	100 ml	

Prepared from the highest purity raw material available, generally greater than 99.999% spectral purity. The content of the solution is confirmed to be accurate to within \pm 0.2%. Typically 1000 μ g/ml. The certificate of analysis provided reports actual lot analysis. The certificate also lists the trace impurities.

Chromium 1000 µg/ml

(Matrix: 2% nitric acid) / 'BAKER INSTRA-ANALYZED' / Atomic Absorption Standard

6926

BOX

M =	52.00 g/mol
NC CODE	3822 00 00
R:	36/38
S:	26-37
	Xi

► Cr

Chromium (Cr)	998-1002 μg/ml	PRODUCT	PACKING
		NO.	
		6926.0100	100 ml
		6926.0500	500 ml

Prepared by dissolution of high purity raw materials (min. 99.99% spectral purity). Assays are verified by ICP against standards traceable to NIST. Standard Reference Material numbers (SRM) are printed on each label.

Chrom

Chromium 1000 μ g/ml 'BAKER ANALYZED' / Atomic Absol

6809

► Cr **M** = 52.00 g/mol

NC CODE 3822 00 00

R: 36/38 **S:** 26-37

Chromium (Cr)

998-1002 $\mu g/ml$

PRODUCT CONT. **PACKING** NO. вох 6809.0100 100 ml 6809.0500 500 ml

Chromium nitrate in nitric acid 0.5 mol/l.

Chromium 10000 µg/ml (Matrix: 2% nitric acid) / 'BAKER INSTRA-ANALYZED' / Plasma Standard

► Cr

M = 52.00 g/mol

NC CODE 3822 00 00

R: 36/38

Certificate Provided Reporting Actual Lot Analysis 9980-10020 μg/ml Chromium (Cr)

PRODUCT PACKING NO. 5727.0100 100 ml

Prepared from the highest purity raw material available, generally greater than 99.999% spectral purity. The content of the solution is confirmed to be accurate to within 0.2 %. Typically 10000 μ g/ml.

CONT.

вох

The certificate of analysis provided reports actual lot analysis. The certificate also lists the trace impurities.

Chromium(III) Chloride Hexahydrate

1055

► CrCl₃.6H₂O

M = 266.45 g/mol **CAS NO.** 10060-12-5

EINECS 233-038-3

NC CODE 2827 39 80 R: 21/22

S: 24/25-28

Assay (by lodometry) 99.0-102.0% Insoluble Matter max. 0.01% max. 0.01% Iron (Fe) pH of 5% Solution at 25°C 2.0-3.0 Sulfate (SO₄) max. 0.01%

PRODUCT PACKING CONT. BOX 1055.0250 250 g

Chromium(III) Oxide

▶ Cr₂O₃

0083

M = 151.99 g/mol **CAS NO.** 1308-38-9 EINECS 215-160-9

NC CODE 2819 90 90

max. 0.01% Ammonium (NH₄) Chloride (CI) max. 0.02% Iron (Fe) max. 0.05% Sulfate (SO₄) max. 0.7%

PRODUCT PACKING CONT. BOX 0083.1000 1 kg

Chromium(VI) Oxide

0086 'BAKER ANALYZED' / ACS

► CrO₃

M = 99.99 g/mol**CAS NO.** 1333-82-0

EINECS 215-607-8

NC CODE 2819 10 00

EC NO. 24 001 00 0

UN/ID NO. 1463

ADR/RID 5.1 OC2

IMDG 5.1/II

R: 25-35-43-49-50/53-8

S: 45-53-60-61









Meets ACS Specifications Assay (by lodometry) min. 98.0% Chloride (CI) max. 0.005% max. 0.01% Insoluble Matter Iron, Aluminium and Barium max. 0.03% max. 0.05% Nitrate (NO₃) Sodium (Na) max. 0.2% max. 0.005% Sulfate (SO₄)

PRODUCT CONT. **PACKING** вох 0086.0500



Chromium Sesquioxide

See Chromium(III) Oxide

Chromium Trioxide

See Chromium(VI) Oxide



Chromotropic Acid, Disodium Salt 'BAKER ANALYZED'

► (HO)₂C₁₀H₄(SO₃Na)₂.2H₂O **M** = 400.30 g/mol

CAS NO. 5808-22-0 EINECS 204-972-9 NC CODE 2908 20 00

Insoluble Matter	max. 0.01%
Suitability for Formaldehyde Determination	passes test
Water (H ₂ O)	about 9%

PRODUCT PACKING CONT. вох NO. 1261.0025 25 g Glass

Citric Acid Anhydrous

'BAKER ANALYZED' / ACS

0090

HOC(COOH)(CH ₂ COOH) ₂
M =	192.13 g/mol
CAS NO.	77-92-9
EINECS	201-069-1
NC CODE	2918 14 00

min. 99.5%
max. 0.001%
max. 0.005%
passes tes
max. 0.001%
max. 0.02%
passes tes
max. 0.002%
max. 3

max. 2

Mosts ACE Enskitions

PRODUCT	PACKING	CONT.
NO.		вох
0090.1000	1 kg	6
0090.9025	25 kg	
0090.9050	50 kg	

Citric Acid Anhydrous

1995

► HOC(COOH)(CH ₂ COOH) ₂		
M =	192.13 g/mol	
CAS NO.	77-92-9	
EINECS	201-069-1	
NC CODE	2918 14 00	

Assay	99.5-101.0%
Appearance of solution	passes test
Heavy Metals (as Pb)	max. 10 ppm
Identification	passes test
Oxalic acid	max. 350 ppm
Readily Carbonisable Substances	passes test
Sulfated Ash	max. 0.1%
Sulfates (as SO ₄)	max. 150 ppm
Water (H ₂ O)	max. 1.0%

PRODUCT	PACKING	CONT.
NO.		вох
1995.1000	1 kg	6
1995.9050	50 kg	

Citric Acid, Anhydrous

See for qualities meeting the requirements of cGMP, USP/NF, BP/Ph.Eur, JP or FCC, page 36

Lead (Pb)



Citric Acid Monohydrate crystal / 'BAKER ANALYZED' / ACS

0088

► HOC(COOH)(CH ₂ COOH) ₂ .H ₂ O			
M = 210.14 g/mol			
CAS NO. 5949-29-1			
EINECS 201-069-1			
NC CODE 2918 14 00			

Meets ACS Specifications	
Assay	99.0-102.0%
Chloride (CI)	max. 0.001%
Insoluble Matter	max. 0.005%
Oxalate (C ₂ O ₄)	passes test
Phosphate (PO ₄)	max. 0.001%
Residue after Ignition	max. 0.02%
Substances Carbonizable by Hot H ₂ SO ₄	
(Tartrates, etc.)	passes test
Sulfate (SO ₄)	max. 0.002%
Trace Impurities (in ppm):	
Iron (Fe)	max. 3
Lead (Pb)	max. 2

PRODUCT NO.	PACKING	CONT. BOX
0088.1000	1 kg	6
0088.5000	5 kg	4
0088.9050	50 kg	



Citri

Citric Acid Monohydrate 2064

stal / 'BAKER

► HOC(COOH)(CH₂COOH)₂.H₂O **M** = 210.14 g/mol

> **CAS NO.** 5949-29-1 **EINECS** 201-069-1 NC CODE 2918 14 00

Assay	99.5-100.5%
Arsenic (As)	max. 3 ppm
Heavy Metals	max. 0.001%
Identification	passes test
Limit of oxalate	passes test
Organic Volatile Impurities	passes test
Readily Carbonisable Substances	passes test
Residue after Ignition	max. 0.05%
Sulfate (SO ₄)	passes test
Water (H ₂ O)	max. 8.8%

PRODUCT	PACKING	CONT.
NO.		вох
2064.1000	1 kg	6
2064.5000	5 kg	
2064.9050	50 kg	

Citric Acid, Monohydrate

See for qualities meeting the requirements of cGMP, USP/NF, BP/Ph.Eur, JP or FCC, page 36

Citric Acid Solutions

See for qualities meeting the requirements of cGMP, USP/NF, BP/Ph.Eur, JP or FCC, page 36



Cleaner for Histology



Cleaners for Use in Electronic Industries

See for detailed information section Microelectronic materials page 32 or the website www.jtbaker.com/micro



Cleaning

7272

Cleaning Solution

NC CODE 3822 00 00

PRODUCT PACKING CONT. NO. BOX 7272,1000 11 7272.9020 20 I

Solution, ready for use.

F780

Cleland's Reagent

'BAKER ULTRAPURE BIOREAGENT'

► HSCH₂(CHOH)₂CH₂SH **M** = 154.25 g/mol **CAS NO.** 27565-41-9 NC CODE 2930 90 70 **R:** 22-36/38

Assay min. 99.5% Appearance passes test 40-43°C Melting Point Absorbance of a 0.1M Solution (1-cm path vs water): at 260 nm max. 0.40 at 280 nm max. 0.06

PRODUCT PACKING CONT. BOX F780.0005 5 g F780.0025 25 g Glass F780.0100 100 g F780.1000 1 kg

Keep material refrigerated between 2-8°C (36-46°F).

CLk-222 Positive Resist Stripper

See for detailed information section Microelectronic materials page 32 or the website www.jtbaker.com/micro

CLk-288 Positive Resist Stripper

See for detailed information section Microele ctronic materials page 32 or the website www.jtbaker.com/micro

CLk-820 Positive Resist Stripper

See for detailed information section Microelectronic materials page 32 or the website www.jtbaker.com/micro



CLk-888 Positive Resist Stripper

See for detailed information section Micros tronic materials page 32 or the website www.jtbaker.com/micro



CMOS chemicals



Gas Chromatography

ULTRA RESI ANALYZED solvents and salts for the industries:

- Pharmaceutical
- Food and Agriculture
- Environmental

We offer ultra high purity solvents and salts for:

- Organic residue extraction/concentration procedures
- Purge and trap analysis

Check our extensive line of ULTRA RESI ANALYZED solvent and reagents at www.jtbaker.com

Application fields of Ultra Resi Analyzed	Pesticide residues	PCB and Dioxine	PAH (Polyaromatic Hydrocarbons)	EOX	Various pharma applications	Capillairy GC	GC/MS
Pharmaceutical							
Food and Agriculture							
Environmental (water and soil)							

pCMP

pCMP 800, pCMP 850 cleaners

rmation section Microelectronic materials page 32 or the website www.jtbaker.com/micro

Cobalt 1000 µg/ml 5712

'BAKER INSTRA-ANALYZED' / Plasma Standard

M = 58.93 g/mol

NC CODE 3822 00 00

EC NO. 7 004 00 1

R: 36/38

S: 26



Certificate Provided Reporting Actual Lot Analysis Cobalt (Co) 998-1002 µg/ml **PRODUCT PACKING** CONT. BOX 5712.0100 100 ml

Prepared from the highest purity raw material available, generally greater than 99.999% spectral purity. The content of the solution is confirmed to be accurate to within \pm 0.2%. Typically 1000 μ g/ml.

The certificate of analysis provided reports actual lot analysis. The certificate also lists the trace impurities.

PACKING

6927

Cobalt 1000 µg/ml

'BAKER INSTRA-ANALYZED' / Atomic Absorption Standard

► Co

M = 58.93 g/mol

NC CODE 3822 00 00

R: 36/38

S: 26



Cobalt (Co) $998\text{-}1002~\mu\text{g/ml}$

NO.		вох
6927.0100	100 ml	
6927.0500	500 ml	

CONT

CONT.

CONT.

вох

BOX

Prepared by dissolution of high purity raw materials (min. 99.99% spectral purity). Assays are verified by ICP against standards traceable to NIST. Standard Reference Material numbers (SRM) are printed on each label.

6810

Cobalt 1000 µg/ml

'BAKER ANALYZED' /

► Co

M = 58.93 g/mol

NC CODE 3822 00 00

R: 36/38 **S:** 26-37

Cobalt (Co) 998-1002 μg/ml

6810.0100 100 ml 6810.0500 500 ml

PRODUCT

PRODUCT

Cobalt nitrate in nitric acid 0.5 mol/l.

PACKING

5728

Cobalt 10000 µg/ml

(Matrix: 2% nitric acid) / 'BAKER INSTRA-ANALYZED' / Plasma Standard

► Co

M = 58.93 g/mol

NC CODE 3822 00 00

EC NO. 7 004 00 1 **R:** 36/38-43

S: 24-36/37



Certificate Provided Reporting Actual Lot Analysis 9980-10020 μg/ml Cobalt (Co)

PRODUCT PACKING NO. 100 ml 5728.0100

Prepared from the highest purity raw material available, generally greater than 99.999% spectral purity. The content of the solution is confirmed to be accurate to within 0.2 %. Typically 10000 $\mu g/ml$.

The certificate of analysis provided reports actual lot analysis. The certificate also lists the trace impurities.

Cobalt Acetate Tetrahydrate'BAKER ANALYZED'

1057

► (CH₃COO)₂Co.4H₂O

M = 249.08 g/mol

CAS NO. 6147-53-1 **EINECS** 200-755-8

NC CODE 2915 23 00

R: 22-42/43 **S:** 22-24-37

Assay (by EDTA titrn.)	min. 99.0%
Chloride (CI)	max. 0.001%
Copper (Cu) (by AAS)	max. 0.002%
Insoluble Matter	max. 0.01%
Iron (Fe) (by AAS)	max. 0.001%
Lead (Pb) (by AAS)	max. 0.001%
Nickel (Ni) (by AAS)	max. 0.15%
Nitrate (NO ₃)	max. 0.015%
pH of 5% Solution at 25°C	6.0-7.5
Sulfate (SO ₄)	max. 0.005%
Zinc (Zn) (by AAS)	max. 0.01%

PRODUCT	PACKING	CONT.
NO.		вох
1057.0100	100 g	

1059



Cobalt Chloride Hexahydrate 'BAKER ANALYZED' / ACS

► CoCl₂.6H₂O **M** = 237.93 g/mol

> **CAS NO.** 7791-13-1 **EINECS** 231-589-4

NC CODE 2827 34 00 **EC NO.** 27 004 00 5

UN/ID NO. 3288

ADR/RID 6.1 T5

IMDG 6.1/III

R: 22-42/43-49-50/53

S: 22-45-53-60-61



environment





Exceeds ACS Specifications. Meets Reagents Specifications for testing USP/NF monographs	
Calcium (Ca)	max. 0.005%
Copper (Cu)	max. 0.002%
Insoluble Matter	max. 0.01%
Iron (Fe)	max. 0.005%
Lead (Pb)	max. 0.002%
Magnesium (Mg)	max. 0.005%
Nickel (Ni)	max. 0.1%
Nitrate (NO ₃)	max. 0.01%
Potassium (K)	max. 0.01%
Sodium (Na)	max. 0.05%
Sulfate (SO ₄)	max. 0.01%
Zinc (Zn)	max. 0.02%

PRODUCT NO.	PACKING	CONT. BOX
1059.0100	100 g	
1059.0250	250 g	6

Cobalt Nitrate Hexahydrate

0091

Co(NO₃)₂.6H₂O

M = 291.04 g/mol

CAS NO. 10026-22-9

EINECS 233-402-1

NC CODE 2834 29 80

UN/ID NO. 1477 **ADR/RID** 5.1 O2

IMDG 5.1/II

R: 22-40-43

S: 36-37



Specifications for testing USP/NF monographs		
98.0-102.0%		
max. 0.005%		
max. 0.002%		
max. 0.002%		
max. 0.01%		
max. 0.001%		
max. 0.002%		

Exceeds ACS Specifications. Meets Reagents

Insoluble Matter	max. 0.01%
Iron (Fe)	max. 0.001%
Lead (Pb)	max. 0.002%
Magnesium (Mg)	max. 0.005%
Nickel (Ni)	max. 0.15%
pH of 5% Solution at 25°C	3.0-6.0
Potassium (K)	max. 0.01%
Sodium (Na)	max. 0.05%
Sulfate (SO ₄)	max. 0.005%
Zinc (Zn)	max. 0.01%

PRODUCT NO.	PACKING	CONT
		ВОХ
0091.0500	500 g	

Cobalt Sulfate Heptahydrate

0093

CoSO₄.7H₂O

M = 281.10 g/mol

CAS NO. 10026-24-1

EINECS 233-334-2

NC CODE 2833 29 30

EC NO. 27 005 00 0

R: 22-42/43-49-50/53

S: 22-45-53-60-61



dangerous for the environment



Assay (by EDTA titrn.)	97.0-103.0%
Chloride (CI)	max. 0.001%
Copper (Cu)	max. 0.002%
Insoluble Matter	max. 0.01%
Iron (Fe)	max. 0.005%
Lead (Pb)	max. 0.002%
Nickel (Ni)	max. 500 ppm
Nitrate (NO ₃)	max. 0.005%
pH of 5% Solution at 25°C	3.0-7.0
Zinc (Zn)	max. 0.02%

PACKING	CONT.
	вох
500 g	
50 kg	
	500 g

Sym-Colidine

ee 2,4,6-Trimethylpyridine



See for qualities meeting the requirements of cGMP, USP/NF, BP/Ph.Eur, JP or FCC, page 36





Congo

Congo Red 1262

► [(NaOSO₂)(NH₂)C₁₀H₅N:NC₆H₄-]₂

M = 696.68 g/mol

CAS NO. 573-58-0 **EINECS** 209-358-4

NC CODE 3204 14 00

EC NO. 611 027 00 8

R: 45-63

S: 45-53



COOMASSIE Brilliant Blue G-250

1473

C₄₇H₅₀N₃NaO₇S₂ Identification (by IR) **PRODUCT PACKING** CONT. passes test **M** = 854.03 g/mol NO. BOX **CAS NO.** 6104-58-1 1473.0025 25 g Glass **EINECS** 228-058-4 NC CODE 3204 12 00 C.I.42655.

COOMASSIE Brilliant Blue R-250

M = 825.98 g/mol

1474

PRODUCT C₄₅H₄₄N₃NaO₇S₂ **PACKING**

CAS NO. 6104-59-2 **EINECS** 228-060-5 NC CODE 3204 12 00

C.I. 42660.

1474.0005

5 g

NO.

max. 5

min. 99.90%

max. 0.005%

min. 99.0%

max. 0.05%

PRODUCT

1262.0100

NO.

PACKING

100 g

CONT.

CONT.

CONT.

вох

вох

PRODUCT

9503 Heavy Foil 0.125 mm (0.005") / 'BAKER ANALYZED'

► Cu Assay min. 99.90% **M** = 63.55 g/mol Antimony and Tin (as Sn) max. 0.005% **CAS NO.** 7440-50-8 Insoluble in HNO₃ max. 0.02%

EINECS 231-159-6 Iron (Fe) max. 0.01% NC CODE 7408 19 90 Lead (Pb) max. 0.005% Manganese (Mn) max. 0.001% Phosphorus (P) max. 0.001%

Trace Impurities (in ppm): Arsenic (As)

PRODUCT PACKING CONT. NO. BOX 9503.0500 500 g

Copper 9504

Wire 0.50 mm Ø / 'BAKER ANALYZED' / ACS

Exceeds ACS Specifications

M = 63.55 g/mol Assay **CAS NO.** 7440-50-8 Antimony and Tin (as Sn) EINECS 231-159-6 Arsenic (As) NC CODE 7408 19 90 Insoluble in diluted HNO₃

max. 5 ppm max. 0.02% max. 0.005% Iron (Fe) Lead (Pb) max. 0.005% Manganese (Mn) max. 0.001% Phosphorus (P) max. 0.001% max. 0.002% Silver (Ag)

NO. BOX 9504.0500 500 g

PACKING

1728

Copper er / Purified

EINECS 231-159-6

NC CODE 7408 19 90

► Cu Assay (Cu) **M** = 63.55 g/mol Insoluble Matter in HNO₃ **CAS NO.** 7440-50-8

PRODUCT PACKING CONT. NO. BOX 1728.0500 500 a 1728.9012 12 kg

Copper 1000 µg/ml

rix: 2% nitric acid) / 'BAKER INSTRA-ANALYZED' / Plasma Standard

5713

► Cu

M = 63.55 g/mol

NC CODE 3822 00 00

R: 36/38

S: 26



Certificate Provided Reporting Actual Lot Analysis

Copper (Cu) 998-1002 μg/ml

PRODUCT PACKING CONT. NO. вох 5713.0100 100 ml

Prepared from the highest purity raw material available, generally greater than 99.999% spectral purity. The content of the solution is confirmed to be accurate to within \pm 0.2%. Typically 1000 μ g/ml.

The certificate of analysis provided reports actual lot analysis. The certificate also lists the trace impurities.

Copper 1000 μ g/ml (Matrix: 2% nitric acid) / 'BAKER INSTRA-ANALYZED' / Atomic Absorption Standard Copper (Cu)

6928



M = 63.55 g/mol

NC CODE 3822 00 00

R: 36/38 **S**: 26



PRODUCT	PACKING	CONT.
NO.		вох
6928.0100	100 ml	
6928.0500	500 ml	

Prepared by dissolution of high purity raw materials (min. 99.99% spectral purity). Assays are verified by ICP against standards traceable to NIST. Standard Reference Material numbers (SRM) are printed on each label.

Copper 1000 µg/ml 'BAKER ANALYZED' / Atomic Absorption Standard

6811

CONT.

вох

► Cu

M = 63.55 g/mol

NC CODE 3822 00 00

R: 36/38

S: 26



998-1002 μg/ml Copper (Cu)

998-1002 μg/ml

NO. 6811.0100 100 ml

PRODUCT

6811.0500 500 ml

Copper(II)nitrate in nitric acid 0.5 mol/l.

PACKING

Copper 10000 µg/ml

(Matrix: 2% nitric acid) / 'BAKER INSTRA-ANALYZED' / Plasma Standard

5729

► Cu

M = 63.55 g/mol

NC CODE 3822 00 00

R: 36/38



Certificate Provided Reporting Actual Lot Analysis Copper (Cu)

9980-10020 μg/ml

PRODUCT PACKING CONT. NO. 5729.0100 100 ml

Prepared from the highest purity raw material available, generally greater than 99.999% spectral purity. The content of the solution is confirmed to be accurate to within 0.2 %. Typically 10000 µg/ml.

The certificate of analysis provided reports actual lot analysis. The certificate also lists the trace impurities.

Copper(II) Acetate Monohydrate

0095

► (CH ₃ COO) ₂ Cu	ı.H ₂ O
M =	199.65 g/m
CAS NO.	6046-93-1

35 g/mol

EINECS 205-553-3

NC CODE 2915 29 00

R: 22



Exceeds ACS Specifications	
Assay (by lodometry)	98.0-102.0%
Calcium (Ca)	max. 0.005%
Chloride (CI)	max. 0.003%
Insoluble Matter	max. 0.01%
Iron (Fe)	max. 0.002%
Nickel (Ni)	max. 0.01%
pH of 5% Solution at 25°C	5.0-6.0
Potassium (K)	max. 0.01%
Sodium (Na)	max. 0.05%

max. 0.005%

Sulfate (SO₄)

PRODUCT NO.	PACKING	CONT. BOX
0095.9050	50 kg	

Coppe

Copper(II) Bromide

➤ CuBr ₂	
M =	223.35 g/mol
CAS NO.	7789-45-9
EINECS	232-167-2
NC CODE	2827 59 00
UN/ID NO.	1759
AND/DID	8 C10

IMDG 8/III R: 22-34 **S:** 26-36/39-45

Assay	111111 3070
Chloride (CI)	max. 0.3%
Insoluble Matter	max. 0.01%
Iron (Fe)	max. 0.005%
pH of 5% Solution at 25°C	2.5-4.5
Sulfate (SO ₄)	max. 0.01%

	PRODUCT	PACKING	CONT.
-	NO.		вох
-	1067.0500	500 g	

Copper(I) Chloride

0108 / 'BAKER ANALYZED' / ACS

► CuCl **M** = 99.00 g/mol **CAS NO.** 7758-89-6 **EINECS** 231-842-9 NC CODE 2827 39 90 **EC NO.** 29 001 00 4 **UN/ID NO.** 2802 ADR/RID 8 C2 IMDG 8/III

R: 22-50/53 **S:** 22-60-61

dangerous for the



Meets ACS Specifications min. 90.0% Assay (CuCI) Calcium (Ca) max. 0.01% Insoluble in Acid max. 0.02% Iron (Fe) max. 0.005% max. 0.02% Potassium (K) Sodium (Na) max. 0.05% Sulfate (SO₄) max. 0.1%

PRODUCI	PACKING	CONT.
NO.		вох
0108.0500	500 g	

Copper(II) Chloride Dihydrate

► CuCl₂.2H₂O **M** = 170.48 g/mol CAS NO. 10125-13-0 **EINECS** 231-210-2 NC CODE 2827 39 80 **UN/ID NO.** 2802 ADR/RID 8 C2 IMDG 8/III **R:** 22-36/38-50/53 **S:** 22-26-61

0097

dangerous for the environment



99.0-101.0% Calcium (Ca) max. 0.005% Insoluble Matter max. 0.01% Iron (Fe) max. 0.005% Nickel (Ni) max. 0.01% Nitrate (NO₃) max. 0.005% pH of 5% Solution at 25°C 2.0-4.0 Potassium (K) max. 0.01% Sodium (Na) max. 0.02% Sulfate (SO₄) max. 0.005%

PACKING вох 0097.0250 250 g

Up-to-date product specifications are available online at www.jtbaker.com/europe

Copper(II) Chloride Dihydrate

► CuCl₂.2H₂O Assay min. 95% PRODUCT PACKING CONT. **M** = 170.48 g/mol max. 0.05% вох Insoluble Matter NO. **CAS NO.** 10125-13-0 0515.1000 1 kg **EINECS** 231-210-2 NC CODE 2827 39 80 **UN/ID NO.** 2802

Copper(II) Nitrate 2.5-Hydrate

0099

0515

Cu(110/3/2.2,31120		
M =	232.59 g/mol	
CAS NO.	3251-23-8	
EINECS 221-838-5		
NC CODE	2834 29 30	
UN/ID NO.	1477	
ADR/RID 5.1 O2		
IMDG	5.1/II	
R:	22-36/38	
S:	24/25	

► Cu/NO \ 25H O

ADR/RID 8 C2 IMDG 8/III

> **R:** 22-36/38-50/53 **S:** 22-26-61

> > dangerous for the environment

Exceeds ACS Specifications	
Assay	98.0-102.0%
Calcium (Ca)	max. 0.005%
Chloride (CI)	max. 0.002%
Insoluble Matter	max. 0.01%
Iron (Fe)	max. 0.005%
Lead (Pb)	max. 0.001%
Nickel (Ni)	max. 0.01%
pH of 5% Solution at 25°C	3.0-4.0
Potassium (K)	max. 0.005%
Sodium (Na)	max. 0.01%
Sulfate (SO ₄)	max. 0.005%

PRODUCT	PACKING	CONT.
NO.		вох
0099.0250	250 g	
0099.9050	50 kg	

Copper(I) Oxide

► Cu₂O **M** = 143.09 g/mol **CAS NO.** 1317-39-1 **EINECS** 215-270-7 NC CODE 2825 50 00 **EC NO.** 29 002 00 0

R: 22 **S:** 22-36-60

Assay	min. 96.0%
Chloride (CI)	max. 0.5%
Free Copper (Cu)	act. value reported
Insoluble in HNO ₃	max. 0.3%
Iron (Fe)	max. 0.05%
Preservative (Zinc Stearate)	0.2-0.5%
Sulfate (SO ₄)	max. 0.05%

PACKING	CONT.
	вох
500 g HDPE	

Copper(II) Oxide

Wire / 'BAKER ANALYZED' / ACS

0102

1878-01

M =	79.55 g/mol
CAS NO.	1317-38-0
EINECS	215-269-1
NC CODE	2825 50 00
R:	22

► CuO

S: 36-60

M =	79.55 g/mol
S NO.	1317-38-0
INECS	215-269-1
CODE	2825 50 00

Exceeds ACS Specifications Carbon Compounds (as C) max. 0.002% max. 0.001% Nitrogen Compounds (as N) max. 0.012% Sulfate (SO₄)

PRODUCI	PACKING	CONT.
NO.		вох
0102.9050	50 kg	



Coppe

Copper(II) Oxide

0101 r / 'BAKER ANALYZED' / ACS

► CuO

M = 79.55 g/mol **CAS NO.** 1317-38-0

EINECS 215-269-1 NC CODE 2825 50 00

R: 22

S: 36-60



Exceeds ACS Specifications	
Assay	min. 99.0%
Average Particle Diameter, µm (APD) (by	
Sedigraph)(typical)	VO
Bulk Density (g/cc)(typical)	VC
Calcium (Ca)	max. 0.01%
Carbon Compounds (as C)	max. 0.01%
Chloride (CI)	max. 0.005%
Free Alkali	passes tes
Insoluble in Dilute HCl	max. 0.02%
Iron (Fe)	max. 0.05%
Nitrogen Compounds (as N)	max. 0.002%
Potassium (K)	max. 0.02%
Silicon (Si)	VO
Sodium (Na)	max. 0.05%
Specific Surface Area, m²/g (typical)	VO
Sulfate (SO ₄)	max. 0.02%
Mesh (Wet Screen Analysis):	
On U.S. No. 325 Sieve	VO

PRODUCT	PACKING	CONT.
NO.		вох
0101.0500	500 g	
0101.9050	50 kg	

vol = value on label.

0107

Copper(II) Sulfate Anhydrous

► CuSO₄

M = 159.61 g/mol

CAS NO. 7758-98-7 **EINECS** 231-847-6

NC CODE 2833 25 00 **EC NO.** 29 004 00 0

UN/ID NO. 3077

ADR/RID 9 M7

IMDG 9/III

R: 22-36/38-50/53

S: 22-60-61





Assay (by lodometry) min. 99.0% Ammonium Sulfide Metals, other than Iron max. 0.01% (as Ni) Chloride (CI) max. 0.002% Insoluble Matter max. 0.01% Iron (Fe) max. 0.005% pH of 5% Solution at 25°C 3.0-5.0

PRODUCT PACKING CONT. NO. BOX 0107.1000 1 kg 50 kg

0104

Copper(II) Sulfate Pentahydrate

Fine Crystals / 'BAKER ANALYZED' / A

► CuSO₄.5H₂O

M = 249.68 g/mol

CAS NO. 7758-99-8 **EINECS** 231-847-6

NC CODE 2833 25 00

EC NO. 29 004 00 0

UN/ID NO. 3077

ADR/RID 9 M7

IMDG 9/III

R: 22-36/38-50/53

S: 22-60-61





Exceeds ACS Specifications		
Assay	98.0-102.0%	
Calcium (Ca)	max. 0.005%	
Chloride (CI)	max. 0.001%	
Insoluble Matter	max. 0.005%	
Iron (Fe)	max. 0.003%	
Nickel (Ni)	max. 0.005%	
Nitrogen Compounds (as N)	max. 0.001%	
Potassium (K)	max. 0.01%	
Sodium (Na)	max. 0.02%	

PRODUCT	PACKING	CONT.
NO.		вох
0104.0250	250 g	6
0104.1000	1 kg	6
0104 9050	50 kg	



J.T.Baker

Copper(II) Sulfate Pentahydrate

0105 ► CuSO₄.5H₂O 98.5-100.5% PRODUCT PACKING CONT. Assay M = 249.68 g/mol NO. вох **CAS NO.** 7758-99-8 0105.1000 1 kg 6 **EINECS** 231-847-6 0105.5000 5 kg NC CODE 2833 25 00 50 kg 0105.9050 **EC NO.** 29 004 00 0 **UN/ID NO.** 3077 ADR/RID 9 M7 IMDG 9/III **R:** 22-36/38-50/53 **S:** 22-60-61

Copper(II) Sulfate

harmful

0.1 mol/l / 'BAKER ANALYZED'

7251

0.0995-0.1005 CONT. ► CuSO₄ Titer (mol/l) **PRODUCT PACKING** вох **M** = 159.61 g/mol NC CODE 2833 25 00 7251.1000 11 6 **R:** 52/53

o-Cresol

S: 61

dangerous for the

environment

Assay (by GC) min. 98% **PRODUCT PACKING**

► CH₃C₆H₄OH CONT. 189-191°C **M =** 108.14 g/mol **Boiling Point** NO. BOX FLASHPOINT 81 °C 29-31°C Freezing Point 7085,1000 11 **CAS NO.** 95-48-7 **EINECS** 202-423-8 NC CODE 2907 12 00 EC NO. 604 004 00 9 **UN/ID NO.** 2076 ADR/RID 6.1 TC2 IMDG 6.1/II R: 24/25-34 **S:** 36/37/39-45

Crystal Violet

1065 'BAKER ANALYZED'

[(CH₃)₂NC₆H₄]₂C:C₆H₄:N(CH₃)₂CI **M** = 408.00 g/mol

CAS NO. 548-62-9

EINECS 208-953-6

NC CODE 3204 13 00

R: 22-40-41-50/53

S: 22-26-36/37/39-61





Assay (dried basis) min. 90% Adsorbance passes test Insoluble Matter passes test Loss on Drying max. 4% Sensitivity passes test **Visual Transition Interval:** 0.8 pH (green)

pH (violet) 2.6

NO. 1065.0025 25 g Glass C.I. 42555

PACKING

PRODUCT

Cupric Sulfate, 5-Hydrate

See for qualities meeting the requirements of cGMP, USP/NF, BP/Ph.Eur, JP or FCC, page 36



CONT.

BOX

7085

Cyclo

Cyclohexane

AKER ULTRA RESI-ANALYZED' / for Organic Residue Analysis

► ÇH₂(CH₂)₄CH₂

9258

M = 84.16 g/mol

1 l = 0.78 kg

FLASHPOINT $-18\,^{\circ}\text{C}$

CAS NO. 110-82-7

EINECS 203-806-2

NC CODE 2902 11 00

EC NO. 601 017 00 1

UN/ID NO. 1145

ADR/RID 3F1

IMDG 3/II

R: 11-38-50/53-65-67

S: 16-33-60-61-62-9







Assay (by GC) (corrected for water)	min. 99.0%
Color (APHA)	max. 10
Residue after Evaporation	max. 2 ppm
Substances Darkened by H ₂ SO ₄	passes test
Water (by KF, coulometric)	max. 100 ppm
ECD Sensitive Impurities (as Heptac	hlor Epoxide):
Single Impurity Peak (pg/ml)	max. 10
FID-Sensitive Impurities (as 2-Octan	iol):

max. 10

min. 99.5%

min. 99.0%

Single Impurity Peak (ng/ml)

Assay (by GC) (corrected for water)

Exceeds ACS Specifications

Assay (by GC)

PRODUCT	PACKING	CONT.
NO.		вох
9258.1000	11	6
9258.2500	2.5	4

Cyclohexane

'BAKER HPLC ANALYZED' / for Use in Liquid Chromatography and Spectrophotometry

► ÇH₂(CH₂)₄CH₂

M = 84.16 g/mol

1 l = 0.78 kg

FLASHPOINT −18 °C

CAS NO. 110-82-7

EINECS 203-806-2

NC CODE 2902 11 00

EC NO. 601 017 00 1

UN/ID NO. 1145

ADR/RID 3F1

IMDG 3/II

R: 11-38-50/53-65-67

S: 16-33-60-61-62-9









highly flammable

Residue after Evaporation	max. 2 ppm
Substances Darkened by H ₂ SO ₄	passes test
Water (by KF, coulometric)	max. 0.01%
Fluorescence Trace Impurities (as quinine base)	
ppb:	
Measured at 450 nm	max. 0.5
Measured at Emission Maximum	max. 1.0
Ultraviolet Absorbance (1.00-cm pa	th vs water):
at 220 nm	max. 0.30
at 230 nm	max. 0.20
at 240 nm	max. 0.08
at 254-400 nm	max. 0.01
UV Cut-off, nm	max. 205

PRODUCT	PACKING	CONT.
NO.		вох
9292.4000	4 l Glass	4

Filtered through a 0.2 micron filter. Packaged under Nitrogen.

Cyclohexane 8026

► ÇH₂(CH₂)₄ÇH₂

M = 84.16 g/mol

1 l = 0.78 kg **FLASHPOINT** -18 °C

CAS NO. 110-82-7

EINECS 203-806-2

NC CODE 2902 11 00

EC NO. 601 017 00 1

UN/ID NO. 1145

ADR/RID 3 F1

IMDG 3/II

R: 11-38-50/53-65-67

S: 16-33-60-61-62-9



environment





Appearance	clear
Color (APHA)	max. 10
Residue after Evaporation	max. 0.002%
Substances Darkened by H ₂ SO ₄	passes test
Water (H ₂ O)	max. 0.02%
Trace Impurities (in ppm):	
Aluminium (AI)	max. 0.5
Barium (Ba)	max. 0.1
Boron (B)	max. 0.02
Cadmium (Cd)	max. 0.05
Calcium (Ca)	max. 0.5
Chromium (Cr)	max. 0.02
Cobalt (Co)	max. 0.02
Copper (Cu)	max. 0.02
Iron (Fe)	max. 0.1
Lead (Pb)	max. 0.1
Magnesium (Mg)	max. 0.1
Manganese (Mn)	max. 0.02
Nickel (Ni)	max. 0.02
Tin (Sn)	max. 0.1
Zinc (Zn)	max. 0.1

PRODUCT	PACKING	CONT.
NO.		вох
8026.1000	11	6
8026.2500	2.5	4
8026.5000	5 EcoTainer	4
8026.9200	200 l	
EcoTainer, th	e metal solvent can for more	e safety in the lab.
		,

Cyclohexane

'BAKER ANALYZED' / Ultraviolet Spectrophotometry / ACS

8706

8104

8028

CONT.

BOX

► ÇH₂(CH₂)₄ÇH₂

M =	84.16 g/mol	
W -	04.10 9/11101	

1 l = 0.78 kg

FLASHPOINT $-18\,^{\circ}\text{C}$

CAS NO. 110-82-7

EINECS 203-806-2

NC CODE 2902 11 00

EC NO. 601 017 00 1

UN/ID NO. 1145

ADR/RID 3F1

IMDG 3/II

R: 11-38-50/53-65-67

S: 16-33-60-61-62-9







Assay (by GC)	min. 99.0%
Appearance	clear
Color (APHA)	max. 10
Residue after Evaporation	max. 0.002%
Substances Darkened by H ₂ SO ₄	passes test
Water (H ₂ O)	max. 0.02%
Ultraviolet Absorbance (1.00-cm pa	th vs water):
at 210 nm	max. 1.00
at 220 nm	max. 0.50
at 230 nm	max. 0.20
at 240 nm	max. 0.08
at 250 nm	max. 0.03
at 260 nm	max. 0.02
at 300-400 nm	max. 0.01

Meets ACS Specifications

PRODUCT NO.	PACKING	CONT. BOX
8706.1000	11	
8706.2500	2.5	

Cyclohexane

M = 84.16 g/mol

1 l = 0.78 kg

FLASHPOINT − 18 °C

CAS NO. 110-82-7

EINECS 203-806-2

NC CODE 2902 11 00

EC NO. 601 017 00 1

UN/ID NO. 1145

ADR/RID 3F1

IMDG 3/II

R: 11-38-50/53-65-67

S: 16-33-60-61-62-9







Assay (by GC) min. 99% Residue after Evaporation max. 0.005% Water (H₂O) max. 0.02%

153-157°C

max. 0.5%

PRODUCT

8028.2500

NO.

PACKING

2.5 |

PRODUCT	PACKING	CONT
NO.		BOX
8104.1000	11	(
8104.2500	2.5	4
8104.9025	25 l	
8104.9200	200 l	
For safe han	dling of 25 l tin cans, see S	elf-closing tap.

Cyclohexanone

► ÇH₂(CH₂)₄ÇO

M = 98.15 g/mol

1 I = 0.94 kg

FLASHPOINT 43 °C

CAS NO. 108-94-1

EINECS 203-806-2

NC CODE 2914 22 00

EC NO. 606 010 00 7

UN/ID NO. 1915 ADR/RID 3 F1

IMDG 3/III

R: 10-20

S: 25

Cyclohexanone MOS Grade

for detailed information section Microelectronic materials page 32 or the website www.jtbaker.com/micro

Boiling Range

Water (H₂O)





Cyclo

Cyclohexene

G045-01

► ÇH:CH(CH₂)₃ÇH₂ **M** = 82.15 g/mol

1 l = 0.81 kg

 $\textbf{FLASHPOINT} \quad < \text{-20 °C}$

CAS NO. 110-83-8

EINECS 203-807-8

NC CODE 2902 1980

UN/ID NO. 2256

ADR/RID 3 F1

IMDG 3/II

R: 11-22-36/37

S: 16-23-33-9





2-(cyclohexylamino)ethanesulfonic acid, Sodium salt

Assay (by GC)

3-(Cyclohexylamino)-1-propanesulfonic acid



L-Cysteine

L-Cysteine Hydrochloride, Monohydrate

See for qualities meeting the requirements of cGMP, USP/NF, BP/Ph.Eur, JP or FCC, page 36

L(+)-Cysteinium Chloride Monohydrate

► HSCH₂CH(NH₂)COOH.HCI.H₂O **M** = 175.64 g/mol

CAS NO. 52-89-1 **EINECS** 200-157-7

NC CODE 2930 90 16

Meets NAS/NRC Specificat	tions and Criteria for	
Biochemical Compounds		
Assay	98.5-100.5%	
Heavy Metals (as Pb)	max. 0.001%	
Homogenity (by TLC)	passes test	
Iron (Fe)	max. 0.003%	
Residue after Ignition	max. 0.1%	
Specific Rotation $[\alpha]^{\frac{26}{5}}(c=2)$ in	5M HCI,	
calcd. as cysteine)	$+6.03^{\circ}$ to $+7.03^{\circ}$	
Sulfate (SO ₄)	max. 0.03%	
Water Insoluble Matter	passes test	

min. 99%

PRODUCT

G045-01

PRODUCT

1576.0100

1576.1000

NO.

max. 1

PACKING

100 g

1 kg

NO.

PACKING

1 | Glass

CONT.

вох

CONT.

вох

Trace Impurities (in ppm): Arsenic (As)



DCYTA ee trans-1,2-Diaminocyclohexane-N,N,N',N'-tetraacetic Acid

DEAE Cellulose



Deblock reagents for use in DNA synthesis

www.jtbaker.com/europe



Deblock-DCA (ABI)

'BAKER ANALYZED' / for DNA/RNA synthesis

9518

1 I = 1.34 kg NC CODE 3822 00 00 UN/ID NO. 2927 ADR/RID 6.1 TC1

R/RID 6.1 IC1 IMDG 6.1/II

R: 36/38-40 **S:** 23-25-26-3



 Suitable for Oligonucleotide Synthesis

 Assay (Dichloroacetic Acid)
 2.90-3.10% (w/v)

 Appearance: clear, colorless solution free from visible particulates
 passes test

 Water (H₂O)
 max. 150 ppm

 PRODUCT
 PACKING
 CONT.

 NO.
 BOX

 9518.0900
 900 ml Glass

 9518.2500
 2.5 l Glass

 9518.4000GL
 4 l Glass

Deblock-TCA (ABI)

'BAKER ANALYZED' / for DNA/RNA synthesis

9448

9451

1 l = 1.34 kg NC CODE 3822 00 00 UN/ID NO. 2927 ADR/RID 6.1 TC1 IMDG 6.1/II R: 36/37/38-40-52/53

S: 23-25-26-3-36/37-61



 Suitable for Oligonucleotide Synthesis

 Assay (Trichloroacetic Acid)
 2.90-3.10% (w/v)

 Appearance: clear, colorless solution free from visible particulates
 passes test

 Water (H₂O)
 max. 150 ppm

 PRODUCT NO.
 PACKING PACKING
 CONT. ROX

 9448.0450
 450 ml Glass

 9448.2000
 2 l Glass

 9448.2500
 2.5 l Glass

Deblock-TCA (ABI)

'BAKER ANALYZED' / for DNA/RNA synthesis

1 l = 1.34 kg NC CODE 3822 00 00 UN/ID NO. 2927

ADR/RID 6.1 TC1 IMDG 6.1/II

R: 36/37/38-40-52/53 **S:** 23-25-26-3-36/37-61

Xn harmful
 Suitable for Oligonucleotide Synthesis

 Assay (Trichloroacetic Acid)
 2.90-3.10% (w/v)

 Appearance: clear, colorless solution free from visible particulates
 passes test

 Water (H₂O)
 max. 150 ppm

 PRODUCT NO.
 PACKING PACKING
 CONT.

 9451.0180
 180 ml Glass

 9451.0900
 900 ml Glass

Deblock TCA-DCE (ABI)

'BAKER ANALYZED' / for DNA/RNA synthesis

9449

11 = 1.25 kg
FLASHPOINT 13 °C
NC CODE 3822 00 00
UN/ID NO. 2924
ADR/RID 3 FC
IMDG 3/II

R: 11-22-36/37/38-45-52/53

S: 16-20-26-33-53-9





 Suitable for Oligonucleotide Synthesis

 Assay (Trichloroacetic Acid)
 2.90-3.10%

 Appearance: clear, colorless solution free from visible particulates
 passes test

 Water (H₂O)
 max. 100 ppm

 PRODUCT NO.
 PACKING SOME PACKING PACKI

Find our up-to-date Product Literature at www.jtbaker.com/europe



7080

Decah

Decahydronaphthalene cis-trans mixture / 'BAKER'

PRODUCT PACKING CONT. ► C₁₀H₁₈ Assay (by GC) min. 98% **M** = 138.25 g/mol NO. вох **1 l** = 0.88 kg 7080.2500 2.5 l FLASHPOINT 58 °C **CAS NO.** 91-17-8 **EINECS** 202-046-9 NC CODE 2902 1980 **UN/ID NO.** 1147 ADR/RID 3 F1 IMDG 3/III **R:** 10-65 **S:** 16-36-60-9

Decalcifier, rapid HISTO GRADE

PRODUCT PACKING CONT. NC CODE 3822 00 00 **UN/ID NO.** 3264 NO. ADR/RID 8 C1 3930.2500 2.5 | HDPE IMDG 8/III Ready to use for Histo-Pathology applications. Rapid Decalcifier contains an acid mixture.

Decalin

See Decahydronaphthalene

8168

n-Decane 'BAKER'

► CH₃(CH₂)₈CH₃ Assay (by GC) **M** = 142.29 g/mol **Boiling Point 1** I = 0.73 kg FLASHPOINT 46 °C **CAS NO.** 124-18-5 EINECS 204-686-4

NC CODE 2901 10 00 **UN/ID NO.** 2247 ADR/RID 3F1 IMDG 3/III R: 10-65



S: 23A-24-62

min. 99% 173-174°C

CONT. PRODUCT PACKING NO. вох 500 ml

вох

Decanedicyl Chloride

See Sebacoyl Chloride

Deuterochloroform

See Chloroform-d₁

Devarda's Alloy 0112 'BAKER ANALYZED'

NC CODE 7403 29 00 Aluminium (AI) 44 - 46% **PRODUCT PACKING** CONT. Copper (Cu) 49 - 51% NO. BOX Nitrogen (N) max. 0.005% 0112.0500 500 g Zinc (Zn) 4-6%

PRODUCT

4893.0100

NO.

PACKING

100 g Glass



Dextrose, Anhydrous

ULTREX Ultrapure Reagent

4893

CONT.

вох

► HOCH₂ÇH(CHOH)₄Q			
180.16 g/mol			
50-99-7			
200-075-1			
1702 30 51			

Actual Analysis Lot. No. N23H06	
Ash (sulfated)	< 0.01%
Loss on Drying at 105°C	0.04%
Particulate Matter	< 0.001%
Specific Rotation [α] $\frac{25}{D}$ ($c = 10$ in water)	52.7 '
Starch	passes tes
Titrable Acid (µeq/g)	0.0008
Metallic Impurities in parts per million	(μ g/g) :
Aluminium (AI)	0.1
Cadmium (Cd)	< 0.1
Calcium (Ca)	0.3
Chromium (Cr)	0.1
Cobalt (Co)	< 0.0
Copper (Cu)	0.3
Iron (Fe)	0.3
Lead (Pb)	0.1
Magnesium (Mg)	0.3
Manganese (Mn)	< 0.008
Molybdenum (Mo)	< 0.0
Nickel (Ni)	< 0.0
Silver (Ag)	0.3
Sodium (Na)	3
Tin (Sn)	
Titanium (Ti)	0.1
Vanadium (V)	< 0.0
Zinc (Zn)	< 0.
Non-Metallic Impurities in parts per m	illion (μg/g):
Arsenic (As)	< '
Boron (B)	0.0
Chloride (CI)	50
Nitrogen Compounds (as N)	< 2
Silicon (Si)	0.3
Sulfate and Sulfite (as SO ₄)	< 50



See for qualities meeting the requirements of cGMP, USP/NF, BP/Ph.Eur, JP or FCC, page 36



Dextrose, Monohydrate

See for qualities meeting the requirements of cGMP, USP/NF, BP/Ph.Eur, JP or FCC, page 36



2,7-Diamino-10-ethyl-9-phenyl- phenanthridinium bromide

See Ethidium Bromide



1269

trans-1,2-Diaminocyclohexane-N,N,N',N'-tetraacetic Acid 'BAKER ANALYZED'

► (HOOCCH₂)₂N<u>C</u>H(CH₂)₄CHN(CH₂COOH)₂.H₂O

M = 364.36 g/mol CAS NO. 13291-61-7 EINECS 236-308-9 NC CODE 2922 49 95
 Assay
 min. 99%

 Heavy Metals (as Pb)
 max. 0.001%

 Iron (Fe)
 max. 0.001%

 Residue after Ignition
 max. 0.1%

 Sulfate (SO₄)
 max. 0.01%

 Water (H₂O)
 4.8-5.6%



trans-1,2-Diaminocyclohexane-N,N,N',N'-tetraacetic acid, Sodium Salt

0.1 mol/l / 'BAKER ANALYZEI



NC CODE 2922 49 70	Titer	0.095-0.105 mol/l	PRODUCT	PACKING	CONT.
			NO.		вох
			7244.1000	11	

1,2-Diaminoethane

See Ethylenediamine



Diamm

Diammonium Hydrogen Citrate

See Ammonium Hydrogen Citrat



Diammonium Hydrogen Phosphate

See Ammonium Hydrogen Phosphate

0610

Diamond Fuchsin

'BAKER ANALYZED

► (CH₃)(NH₂)C₆H₃C(C₆H₄NH₂):C₆H₄:NH.HCl **M** = 337.86 g/mol

CAS NO. 632-99-5 EINECS 211-189-6 NC CODE 2925 20 00 Insoluble Matterpasses testVisual Transition Interval:purplepH 1.0purplepH 3.1red

min. 96%

clear

max. 30

min. 99%

PRODUCT

8578.1000

NO.

PACKING

11

PRODUCT

7368.1000

PACKING

11

CONT.

CONT.

вох

BOX

 PRODUCT NO.
 PACKING BOX

 0619.0025
 25 g Glass

 6
 C.I.42510.



1,3-Diazole

ee Imidazole



5',5-Dibromo-o-cresolsulfonphthalein

Assay (by GC)

Appearance

Color (APHA)

See Bromocresol Purple

7368

1,2-Dibromoethane

'BAKER

7500

1 I = 2.15 kg

CAS NO. 106-93-4

--- ...

EINECS 203-444-5

NC CODE 2903 30 36 EC NO. 602 010 00 6

UN/ID NO. 1605

ADR/RID 6.1 T1

, IMDG 6.1/I

R: 23/24/25-36/37/38-45-51/53

S: 45-53-61







3',3-Dibromothymolsulfonphthalein

Assay

See Bromothymol Blue

8578

Dibutylamine 'BAKER'

ightharpoonup (CH₃CH₂CH₂CH₂)₂NH $m M = 129.25 \ g/mol$

1 I = 0.76 kg FLASHPOINT 39 °C

FLASHPOINT 39 °C

CAS NO. 111-92-2

EINECS 203-921-8

NC CODE 2921 19 80

EC NO. 612 049 00 0

UN/ID NO. 2248

ADR/RID 8 CF1 IMDG 8/II

R: 10-20/21/22



146

Dibutyl Phtalate

7024

► C₆H₄(COOC₄H₉)₂

M = 278.35 g/mol

1 l = 1.04 kg

FLASHPOINT 157 °C

CAS NO. 84-74-2

EINECS 201-557-4

NC CODE 2917 31 00

EC NO. 607 318 00 4

UN/ID NO. 3082

ADR/RID 9 M6

IMDG 9/III

R: 50-61-62

S: 45-53-61





Assay min. 99.0% Density (g/ml) at 25°C 1.041-1.043 PRODUCT **PACKING** CONT. NO. вох 7024.1000 11

o-Dichlorobenzene

'BAKER HPLC ANALYZED' / for use in High Performance Liquid Chromatography

9233

7025

► C₆H₄Cl₂

M = 147.00 g/mol

1 l = 1.30 kg

FLASHPOINT 66 °C

CAS NO. 95-50-1

EINECS 202-425-9

NC CODE 2903 61 00

EC NO. 602 034 00 7

UN/ID NO. 1591 **ADR/RID** 6.1 T1

IMDG 6.1/III

R: 22-36/37/38-50/53

S: 23-60-61





Assay (by GC) (corrected for water) min. 98.5% Residue after Evaporation max. 8 ppm Water (H₂O) max. 0.02% Ultraviolet Absorbance (1.00-cm path vs water): at 300 nm max. 0.30 at 305 nm max. 0.10 max. 0.05 at 335 nm at 375 nm max. 0.01 UV Cut-off, nm max. 296

PRODUCT	PACKING	CONT.
NO.		вох
9233.1000	11	
9233.4000	4 Glass	

Filtered through a 0.2 micron filter. Packaged under Nitrogen.

o-Dichlorobenzene

► C₆H₄Cl₂

M = 147.00 g/mol

1 I = 1.30 kg

FLASHPOINT 66 °C

CAS NO. 95-50-1

EINECS 202-425-9

NC CODE 2903 61 00

EC NO. 602 034 00 7

UN/ID NO. 1591

ADR/RID 6.1 T1

IMDG 6.1/III

R: 22-36/37/38-50/53

S: 23-60-61



environment



Assay (by GC)	min. 99%
Acidity (as HCI)	max. 0.005%
Boiling Range	179.5-181.5°C
Density (g/ml) at 25°C	1.299-1.302
Residue after Evaporation	max. 0.005%
Water (H ₂ O)	max. 0.05%

PRODUCT	PACKING	CONT
NO.		ВОХ
7025.1000	11	



Dichl

o-Dichlorobenzene-Phenol Mix

1:1 / 'BAKER ANALYZED' / for the viscosity dete 7227 rmination of plastics, DIN 53728

Color (APHA)

Water (H₂O)

Meets ACS Specifications

1 l = 1.18 kg FLASHPOINT 79 °C NC CODE 3822 00 00 **UN/ID NO.** 2810 ADR/RID 6.1 T1

IMDG 6.1/II

R: 23/24/25-34-37-48/20/21/22-50/53-68

S: 26-28A-36/37/39-45-57

dangerous for the environment



PRODUCT PACKING CONT. NO. вох 7227.2500 2.5 l

Volumetric Solution, ready for use.

PACKING

500 ml

CONT.

вох

PRODUCT

8707.0500

max. 200

max. 300 ppm

1,2-Dichloroethane

8707 PHOTREX Reagent / For Spectrophotometry / ACS

► CICH₂CH₂CI **M** = 98.96 g/mol **1 l** = 1.25 kg FLASHPOINT 13 °C

CAS NO. 107-06-2 **EINECS** 203-458-1 NC CODE 2903 15 00 **EC NO.** 602 012 00 7

UN/ID NO. 1184 ADR/RID 3 FT1 IMDG 3/II

R: 11-22-36/37/38-45

S: 45-53





Assay (by GC) (corrected for water)	min. 99.0%
Appearance	passes test
Color (APHA)	max. 10
Residue after Evaporation	max. 0.0005%
Titrable Acid (meq/g)	max. 0.0003
Water (H ₂ O)	max. 0.03%
Physical Data (not specifications):	
Boiling Point (typical)	83.5°C
Density (g/ml) at 25°C (typical)	1.246
Ultraviolet Absorbance (1.00-cm pat	h vs water):
at 228 nm	max. 1.00
at 230 nm	max. 0.50
at 235 nm	max. 0.20
at 240 nm	max. 0.10
at 245 nm	max. 0.05
at 250 nm	max. 0.02
at 255 nm	max. 0.01
at 400 nm	max. 0.01
Windows of Infrared Transmittance	(0.1-mm path,
50-100% T), μm:	
2.5-3.3	passes test
3.4-6.7	passes test
7.2-7.4	passes test
8.3-9.6	passes test
9.8-10.4	passes test
11.6-12.7	passes test

1,2-Dichloroethane

'BAKER ANALYZED' / ACS

8042 ► CICH₂CH₂CI

M = 98.96 g/mol

1 l = 1.25 kg FLASHPOINT 13 °C **CAS NO.** 107-06-2

> **EINECS** 203-458-1 NC CODE 2903 15 00

EC NO. 602 012 00 7

UN/ID NO. 1184 ADR/RID 3 FT1 IMDG 3/II

R: 11-22-36/37/38-45

S: 45-53





Meets ACS Specifications	
Assay (by GC)	min. 99.0%
Appearance	clear
Color (APHA)	max. 10
Residue after Evaporation	max. 0.002%
Titrable Acid (meq/g)	max. 0.0003
Water (H ₂ O)	max. 0.03%

PRODUCT	PACKING	CONT.
NO.		вох
8042.1000	11	6
8042.2500	2.5	4



1,2-Dichloroethane

See for qualities meeting the requirements of cGMP, USP/NF, BP/Ph.Eur, JP or FCC, page 36

Dichloromethane

'BAKER ULTRA RESI-ANALYZED' / for Organic Residue Analysis

926/

M = 84.93 g/mol

1 l = 1.32 kg

CAS NO. 75-09-2

EINECS 200-838-9

NC CODE 2903 12 00

EC NO. 602 004 00 3

UN/ID NO. 1593

ADR/RID 6.1 T1

IMDG 6.1/III

R: 40

S: 23-24/25-36/37



A (b CC) (b	:- 00.00/	
Assay (by GC) (exclusive of preservative)	min. 99.8%	
Chloride (CI)	max. 10 ppm	
Color (APHA)	max. 10	
Residue after Evaporation	max. 1 ppm	
Titrable Acid (meq/g)	max. 0.0003	
Water (H ₂ O) (by Coulometry)	max. 0.02%	
ECD Sensitive Impurities (as Heptachlor Epoxide):		
01 1 1 1 1 1 1 1 1 1 1 1	4.0	

Single Impurity Peak (pg/ml) max. 10

FID-Sensitive Impurities (as 2-Octanol):

Single Impurity Peak (ng/ml) max. 5

PRODUCT NO.	PACKING	CONT
		ВОХ
9264.1000	11	ϵ
9264.2500	2.5	4

Dichloromethane

'BAKER HPLC ANALYZED' / for use in High Performance Liquid Chromatography

9315



M = 84.93 g/mol

1 l = 1.32 kg

CAS NO. 75-09-2

EINECS 200-838-9 **NC CODE** 2903 12 00

EC NO. 602 004 00 3

UN/ID NO. 1593

ADR/RID 6.1 T1

IMDG 6.1/III

R: 40

S: 23-24/25-36/37



Assay (by GC)	min. 99.8%
Chloride (CI)	max. 10 ppm
Density (g/ml) at 25°C	1.315-1.321
Preservative (Amylene) (ppm)	about 20
Residue after Evaporation (in ppm)	max. 2
Titrable Acid (meq/g)	max. 0.0003
Water (H ₂ O)	max. 0.02%

Ultraviolet Absorbance (1.00-cm path vs water):	
at 254 nm	max. 0.01
at 280 nm	max. 0.01
at 350 nm	max. 0.01
UV Cut-off, nm	max. 233

PRODUCT NO.	PACKING	CONT. BOX
9315.1000	11	6
9315.2500	2.5	4

Filtered through a 0.2 micron filter. Packaged under Nitrogen.

Dichloromethane

BakerDRY / Low Water Solvent / ACS

9295

► CH₂Cl₂

M = 84.93 g/mol

1 l = 1.32 kg

CAS NO. 75-09-2

EINECS 200-838-9

NC CODE 2903 12 00

EC NO. 602 004 00 3 **UN/ID NO.** 1593

ADR/RID 6.1 T1

IMDG 6.1/III R: 40

S: 23-24/25-36/37



Meets ACS Specifications	
Assay (CH ₂ CL ₂)(by GC, corrected for	
water)	min. 99.5%
Appearance	passes test
Color (APHA)	max. 10
Free Halogens	passes test
Preservative (cyclohexane)	information only
Residue after Evaporation	max. 2 ppm
Titrable Acid (µeq/g)	max. 0.2

max. 30 ppm

PRODUCT	PACKING	CONT.
NO.		вох
9295.1000	11	

Innovation is principal to our business.

Water (by KF, coulometric)



Dichl

Dichloromethane

9316 'BAKER BIO-ANALYZED

► CH₂Cl₂

M = 84.93 g/mol

1 l = 1.32 kg

CAS NO. 75-09-2 **EINECS** 200-838-9

NC CODE 2903 12 00

EC NO. 602 004 00 3

UN/ID NO. 1593

ADR/RID 6.1 T1

IMDG 6.1/III

R: 40 **S:** 23-24/25-36/37

Xn harmful

Assay (CH ₂ Cl ₂)(by GC)	min. 99.5%
Color (APHA)	max. 10
Preservative (as Amylene)	about 20 ppm
Residue after Evaporation	max. 5 ppm
Titrable Acid (meq/g)	max. 0.0001
Water (by KF, coulometric)	max. 30 ppm

PRODUCT	PACKING	CONT.
NO.		вох
9316.2500GL	2.5 Glass	
9316.4000GL	4 l Glass	

7053 T

Dichloromethane

'BAKER ANALYZED' / ACS

► CH₂Cl₂

M = 84.93 g/mol

1 l = 1.32 kg

CAS NO. 75-09-2

EINECS 200-838-9

NC CODE 2903 12 00

EC NO. 602 004 00 3

UN/ID NO. 1593 **ADR/RID** 6.1 T1

IMDG 6.1/III

R: 40

S: 23-24/25-36/37



Exceeds ACS Specifications	
Assay (by GC)	min. 99.5%
Appearance	clear
Color (APHA)	max. 10
Free Halogens	passes test
Preservative (as Amylene)	about 20 ppm
Residue after Evaporation	max. 0.002%
Titrable Acid (meq/g)	max. 0.0002
Water (H ₂ O)	max. 0.02%
Trace Impurities (in ppm):	
Aluminium (AI)	max. 0.5
Barium (Ba)	max. 0.1
Boron (B)	max. 0.02
Cadmium (Cd)	max. 0.05
Calcium (Ca)	max. 0.5
Chromium (Cr)	max. 0.02
Cobalt (Co)	max. 0.02
Copper (Cu)	max. 0.02
Iron (Fe)	max. 0.1
Lead (Pb)	max. 0.1
Magnesium (Mg)	max. 0.1
Manganese (Mn)	max. 0.02
Nickel (Ni)	max. 0.02
Tin (Sn)	max. 0.1
Zinc (Zn)	max. 0.1

NO.		B0)
7053.1000	11	(
7053.2500	2.5	4
7053.5000	5 EcoTainer	2
7053.9025	25	
7053.9200	200	
EcoTainer, the	e metal solvent can for more safety	in the lab
For safe hand	lling of 25 l tin cans, see Self-closing	n tap.
rui saie Ilaliu	illing of 23 i till calls, see self-closing	ј тар.

PRODUCT PACKING

Dichloromethane

'BAKER ANALYZED' / Ultraviolet Spectrophotometry / ACS

7305
► CH₂Cl₂

M = 84.93 g/mol **1** l = 1.32 kg

CAS NO. 75-09-2

EINECS 200-838-9

NC CODE 2903 12 00

EC NO. 602 004 00 3

UN/ID NO. 1593

ADR/RID 6.1 T1 **IMDG** 6.1/III

R: 40

S: 23-24/25-36/37



Exceeds ACS Specifications	
Assay (by GC)	min. 99.5%
Appearance	clear
Boiling Range (initial to dry point)	max. 2.0°C
Color (APHA)	max. 10
Density (g/ml) at 25°C	1.315-1.321
Free Halogens	passes test
Preservative (as Amylene)	about 20 ppm
Recorded Boiling Point	39.8°C
Residue after Evaporation	max. 5 ppm
Titrable Acid (meq/g)	max. 0.0002
Water (H ₂ O)	max. 0.02%
Ultraviolet Absorbance (1.00-cm pa	th vs water):
at 235 nm	max. 1.00
at 240 nm	max. 0.35
at 250 nm	max. 0.10
at 260 nm	max. 0.04
at 340-400 nm	max. 0.01

PRODUCI	PACKING	CONT
NO.		B0)
7305.1000	11	(

Dichloromethane

'BAKER' 7153

► CH₂Cl₂

M = 84.93 g/mol

1 I = 1.32 kg

CAS NO. 75-09-2

EINECS 200-838-9

NC CODE 2903 12 00

EC NO. 602 004 00 3

UN/ID NO. 1593

ADR/RID 6.1 T1

IMDG 6.1/III R: 40

S: 23-24/25-36/37



Assay (by GC) min. 99%
Boiling Range 38-41°C
Preservative (as Amylene) about 20 ppm
Residue after Evaporation max. 0.005%
Water (H₂O) max. 0.1%

PRODUCT NO.	PACKING CONT. BOX
7153.2500	2.5 4
7153.5000	5 EcoTainer
7153.9025	25
7153.9200	200
EcoTainer, th	e metal solvent can for more safety in the lab.

For safe handling of 25 l tin cans, see Self-closing tap.

PACKING

100 ml

Stabilized with tert-Butylpyrocatechol.

166-173°C

PRODUCT

8599.0100

Dichromium Trioxide

See Chromium(VI) Oxide



8599

7026

CONT.

вох

Dicyclopentadiene

► ÇH:CHCHCH₂ÇHÇHCHCH₂CH:ÇH

M = 132.21 g/mol

1 l = 0.98 kg

FLASHPOINT 32 $^{\circ}\mathrm{C}$

CAS NO. 77-73-6

EINECS 201-052-9

NC CODE 2902 19 80

EC NO. 601 044 00 9

UN/ID NO. 2048

ADR/RID 3 F1 IMDG 3/III

R: 11-20/22-36/37/38-51/53

S: 36/37-61







highly

Boiling Point

Diethanolamine

BAKER ANALYZED

► (HOCH₂CH₂)₂NH

M = 105.14 g/mol

1 l = 1.09 kg

FLASHPOINT 152 °C

CAS NO. 111-42-2

EINECS 203-868-0 NC CODE 2922 12 00

EC NO. 603 071 00 1

R: 22-38-41-48/22

S: 26-36/37/39-46

X

Assay	min. 98.5%
Density (g/ml) at 30°C	1.088-1.093
Heavy Metals (as Pb)	max. 0.002%
Iron (Fe)	max. 0.001%
Residue after Ignition	max. 0.005%
Trace Impurities (in ppm):	
Chloride (CI)	max. 5

PRODUCT PACKING	CONT.	
NO.		вох
7026.2500	2.5 l	4

www.jtbaker.com/europe



Dieth

Diethylamine 8030 BAKER ANALYZED

► (C₂H₅)₂NH

M = 73.14 g/mol

1 l = 0.70 kg

FLASHPOINT -39 °C

CAS NO. 109-89-7

EINECS 203-716-3

NC CODE 2921 12 00

EC NO. 612 003 00 0

UN/ID NO. 1154

ADR/RID 3 FC IMDG 3/II

R: 11-20/21/22-35

S: 16-26-29-3-36/37/39-45





Assay (by GC)	min. 99.0%
Density (g/ml) at 25°C	0.698-0.702
Insoluble Matter	passes test
Residue after Evaporation	max. 0.01%

PRODUCT	PACKING	CONT.
NO.		вох
8030.1000	11	6

Diethyl Ether

'BAKER HPLC ANALYZED' / for use in Liquid Chromatography and Spectrophotometry

9237 ► (C₂H₅)₂O

M = 74.12 g/mol

1 I = 0.71 kg

FLASHPOINT −45 °C

CAS NO. 60-29-7

EINECS 200-467-2

NC CODE 2909 11 00

EC NO. 603 022 00 4

UN/ID NO. 1155

ADR/RID 3 F1

IMDG 3/I

R: 12-19-22-66-67

S: 16-29-33-9





Assay (by GC) (corrected for water)	
(exclusive of preservative)	min. 99.5%
Peroxide (as H ₂ O ₂)	max. 5 ppm
Preservative (C ₂ H ₅ OH)	1.5-2.5%
Residue after Evaporation	max. 5 ppm
Titrable Acid (µeq/g)	max. 0.2

Fluorescence Trace Impurities (as quinine base), ppb:

Measured at 450 nm max. 0.5 Measured at Emission Maximum for Solvent Impurities max. 1.0

Ultraviolet Absorbance (1.00-cm path vs water):

at 231 nm	max. 0.4
at 254 nm	max. 0.07
at 280 nm	max. 0.02
at 350-400 nm	max. 0.01
UV Cut-off, nm	max. 220

PRODUCT	PACKING	CONT.
NO.		вох
9237.1000	11	6
9237.2500	2.5	

Contains Alcohol as a preservative. Filtered through a 0.2 micron filter. Packaged under Nitrogen.

Diethyl Ether 8033 'BAKER ANALYZED' / ACS

► (C₂H₅)₂O

M = 74.12 g/mol

11 = 0.71 kg

FLASHPOINT −45 °C **CAS NO.** 60-29-7

EINECS 200-467-2

NC CODE 2909 11 00

EC NO. 603 022 00 4

UN/ID NO. 1155

ADR/RID 3 F1

IMDG 3/I

R: 12-19-22-66-67

S: 16-29-33-9





Meets ACS Specifications. Meets Reagent Specifications for testing USP/NF monographs

min. 98.0% Assay (by GC) (exclusive of preservative) Carbonyl Compounds (as HCHO) (by polarography) max. 0.001% Color (APHA) max. 10 Peroxide (as H₂O₂) max. 1 ppm Preservative (BHT) min. 7 ppm Preservative (C₂H₅OH) 1.5-3.5% max. 0.5% (w/w) Preservative (H₂O) Residue after Evaporation max. 0.001% Titrable Acid (µeq/g) max. 0.2

Product Information (not specifications):

Density (g/ml) at 25°C (typical) 0.713

PRODUCT PACKING CONT. NO. BOX 8033.1000 1 | EcoTainer 8033.5000 5 | EcoTainer EcoTainer, the metal solvent can for more safety in the lab.

Contains Alcohol, Water and BHT as Preservatives.

Diethyl Ether, Anhydrous'BAKER ULTRA RESI-ANALYZED' / for Organic Residue Analysi

9259

► (C ₂ H ₅) ₂ O	
M =	74.12 g/mol
1 i =	0.71 kg
FLASHPOINT	<-20 °C
CAS NO.	60-29-7
EINECS	200-467-2
NC CODE	2909 11 00
EC NO.	603 022 00 4
UN/ID NO.	1155
ADR/RID	3 F1
IMDG	3/I
R:	12-19-66-67
S:	16-29-33-9

Assay (by GC) (corrected for water)	
(exclusive of preservative)	min. 99.0%
Diethyl Ether Interferences:	
Color (APHA)	max. 10
Peroxide (as H ₂ O ₂)	max. 5 ppm
Residue after Evaporation	max. 1 ppm
Substances Darkened by H ₂ SO ₄	passes test
Suitability for U.S. EPA Method 8151A	passes test
Titrable Acid (µeq/g)	max. 0.2
Water (by KF, coulometric)	max. 0.08%

PRODUCT	PACKING	CONT.
NO.		вох
9259.1000	11	6
9259.4000	4 Glass	

Diethyl Ether, Anhydrous

9250

C ₂ H ₅) ₂ O	
M =	74.12 g/mol
1 =	0.71 kg
FLASHPOINT	−20 °C
CAS NO.	60-29-7
EINECS	200-467-2
NC CODE	2909 11 00
EC NO.	603 022 00 4
UN/ID NO.	1155
ADR/RID	3 F1
IMDG	3/I
R:	12-19-22-66-67
S:	16-29-33-9

Meets ACS Specifications		
Assay $(C_2H_5)_2O$) (by GC)	min. 99.0%	
Alcohol (C ₂ H ₅ OH)	max. 0.05%	
Carbonyl (HCHO) (by polarography)	max. 0.001%	
Color (APHA)	max. 10	
Peroxide (as H ₂ O ₂)	max. 1 ppm	
Residue after Evaporation	max. 5 ppm	
Titrable Acid (µeq/g)	max. 0.2	
Water (by KF, coulometric)	max. 10 ppm	
Product Information (not specifications):		
Density (g/ml) at 25°C (typical)	0.7079	

PRODUCT	PACKING	CO
NO.		BO
9250.1000	11	

Mallinckrodt Baker's cGMP Manufactured Chemicals for the Biopharmaceutical industry are a necessity for uncomplicated scale-up.

See chapter 6 of this catalogue.



Dieth

Diethyl Ether, Anhydrous 'BAKER ANALYZED' / ACS

8254

► (C₂H₅)₂O **M** = 74.12 g/mol **1 l** = 0.71 kg **FLASHPOINT** −20 °C **CAS NO.** 60-29-7 **EINECS** 200-467-2 NC CODE 2909 11 00 **EC NO.** 603 022 00 4 **UN/ID NO.** 1155 ADR/RID 3F1 **IMDG** 3/I

R: 12-19-22-66-67 **S:** 16-29-33-9



Exceeds ACS Specifications	
Assay (by GC)	min. 99.5%
Alcohol (C ₂ H ₅ OH)	max. 0.01%
Carbonyl Compounds (as HCHO)	max. 0.001%
Color (APHA)	max. 10
Density (g/ml) at 20°C	0.714-0.716
Foreign Odor	passes test
Peroxide (as H ₂ O ₂)	max. 3 ppm
Residue after Evaporation	max. 0.001%
Substances Darkened by H ₂ SO ₄	passes test
Titrable Acid (meq/g)	max. 0.0002
Water (H ₂ O)	max. 0.01%
Trace Impurities (in ppm):	
Aluminium (AI)	max. 0.5
Barium (Ba)	max. 0.1
Boron (B)	max. 0.02
Cadmium (Cd)	max. 0.05
Calcium (Ca)	max. 0.5
Chromium (Cr)	max. 0.02
Cobalt (Co)	max. 0.02
Copper (Cu)	max. 0.02
Iron (Fe)	max. 0.1
Lead (Pb)	max. 0.1
Magnesium (Mg)	max. 0.1
Manganese (Mn)	max. 0.02
Nickel (Ni)	max. 0.02
Tin (Sn)	max. 0.1
Zinc (Zn)	max. 0.1

PRODUCT	PACKING	CONT.
NO.		вох
8254.1000	11	6
8254.1000S	1 EcoTainer	6
8254.5000	5 EcoTainer	
8254.9025	25	
FcoTainer the	e metal solvent can for more safety	in the lab

For safe handling of 25 l tin cans, see Self-closing tap.

Stabilized with 7 ppm BHT.

8401

► (C₂H₅)₂O

Diethyl Ether, Anhydrous

BAKER ANALYZED

M = 74.12 g/mol **1 l =** 0.71 kg FLASHPOINT −20 °C **CAS NO.** 60-29-7 **EINECS** 200-467-2 NC CODE 2909 11 00

EC NO. 603 022 00 4 **UN/ID NO.** 1155

ADR/RID 3 F1 **IMDG** 3/I

R: 12-19-22-66-67

S: 16-29-33-9





Assay (by GC)	min. 99.5%
Alcohol (C ₂ H ₅ OH)	max. 0.01%
Carbonyl Compounds (as HCHO)	max. 0.001%
Color (APHA)	max. 10
Density (g/ml) at 20°C	0.714-0.716
Foreign Odor	passes test
Peroxide (as H ₂ O ₂)	max. 3 ppm
Residue after Evaporation	max. 0.001%
Substances Darkened by H ₂ SO ₄	passes test
Water (H ₂ O)	max. 0.01%
Trace Impurities (in ppm):	
Aluminium (AI)	max. 0.5
Barium (Ba)	max. 0.1
Boron (B)	max. 0.02
Cadmium (Cd)	max. 0.05
Calcium (Ca)	max. 0.5
Chromium (Cr)	max. 0.02
Cobalt (Co)	max. 0.02
Copper (Cu)	max. 0.02
Iron (Fe)	max. 0.1
Lead (Pb)	max. 0.1
Magnesium (Mg)	max. 0.1
Manganese (Mn)	max. 0.02
Nickel (Ni)	max. 0.02
Tin (Sn)	max. 0.1
Zinc (Zn)	max. 0.1

PRODUCT	PACKING	CONT
NO.		B0
8401.1000	11	
8401.5000	5 EcoTainer	
EcoTainer, the	e metal solvent can for more safety in	the lab



1,4-Diethylene Dioxide



Diethylene Glycol



Diethylene Glycol Monobutyl Ether



min. 98.0%

1.054-1.055

passes test

PRODUCT

7028.1000

PACKING

11

Assay (by GC)

Density (g/ml) at 20°C

Identification (by IR)

7028

CONT.

BOX

4



M = 160.17 g/mol

1 l = 1.07 kg

FLASHPOINT 93 °C

CAS NO. 105-53-3

EINECS 203-305-9

NC CODE 2917 19 10 **ADR/RID** 3,32c

S: 24/25

Diglycol



9,10-Dihydro-9-oxoanthracene



3,4-Dihydroxy-2-anthraquinonesulfonic acid, sodium salt



2,2'-Dihydroxydiethylamine



2,2'-Dihydroxyethyl Ether

See 2,2'-Oxydiethanol



3',6'-Dihydroxyfluoran

e Fluorescein Disodium Salt



4,5-Dihydroxy-2,7-naphthalenedisulfonic acid, disodium salt

Magnesium (Mg)

Manganese (Mn)

Nickel (Ni)

Tin (Sn)

Zinc (Zn)



Diiron Trioxide

e Iron(III) Oxid



Diisopropyl Ether



CONT.



M = 102.18 g/mol

1 l = 0.72 kg

FLASHPOINT −20 °C

CAS NO. 108-20-3

EINECS 203-560-6

NC CODE 2909 19 00

EC NO. 603 045 00 0

UN/ID NO. 1159

ADR/RID 3 F1

IMDG 3/II

R: 11-19-66-67

S: 16-33-9



min. 98%
max. 0.005%
66-69°C
0.716-0.720
max. 0.05%
max. 0.005%
max. 0.5
max. 0.1
max. 0.02
max. 0.05
max. 0.5
max. 0.02
max. 0.02
max. 0.02
max. 0.1
max. 0.1

NO.		ВОХ
8072.1000	11	ϵ

PACKING

PRODUCT

max. 0.1

max. 0.02

max. 0.02

max. 0.1

max. 0.1



Dimet



Dimethoxystrychnine

7030

N,N-Dimethylacetamide

► CH₃CON(CH₃)₂

M = 87.12 g/mol

1 l = 0.93 kg FLASHPOINT 70 °C

CAS NO. 127-19-5

EINECS 204-826-4

NC CODE 2924 10 00 **EC NO.** 616 011 00 4

R: 20/21-61

S: 45-53

Boiling Range	165-167°C
Density (g/ml) at 25°C	0.933-0.937
pH of 20% Solution at 25°C	4.0-7.0
Residue after Ignition	max. 0.005%
Water (H ₂ O)	max. 0.05%

PRODUCT	PACKING	CONT.
NO.		вох
7030.1000	11	6



p-(Dimethylamino)azobenzene-o-carboxylic Acid

1279

p-Dimethylaminobenzaldehyde

► (CH₃)₂NC₆H₄CHO **M** = 149.19 g/mol

> **CAS NO.** 100-10-7 EINECS 202-819-0 NC CODE 2922 30 00

Assay (by Perchloric Acid titrn.) min 99% Heavy Metals (as Pb) max. 0.001% Iron (Fe) max. 0.001% Residue after Ignition max. 0.1%

PRODUCT PACKING CONT. вох NO.



o-[(p-(Dimethylamino(phenyl]azo] benzoic Acid



Dimethylbenzene

Dimethylformamide

BakerDRY / For Use in Organic Synthesis / ACS

► HCON(CH₃)₂

M = 73.10 g/mol **1 l** = 0.95 kg

FLASHPOINT 58 °C

CAS NO. 68-12-2

EINECS 200-679-5

NC CODE 2924 19 00

EC NO. 616 001 00 0

UN/ID NO. 2265

ADR/RID 3 F1

IMDG 3/III R: 20/21-36-61

S: 45-53



Meets ACS Specifications		
Assay (HCON(CH ₃) ₂) (by GC,corrected for		
water)	min. 99.8%	
Appearance	passes test	
Color (APHA)	max. 15	
Residue after Evaporation	max. 0.005%	
Titrable Acid (meq/g)	max. 0.0005	
Titrable Base (meq/g)	max. 0.003	
Water (by KF, coulometric)	max. 20 ppm	
Product Information (not specifications):		
Boiling Point (typical)	153.0°C	

Dimethylformamide can have a teratogenic effect. Women in childbearing age are warned.

CONT.

PACKING

PRODUCT

9213.1000

0.944

Innovation is principal to our business.

Density (g/ml) at 25°C (typical)

9344

Dimethylformamide'BAKER BIO-ANALYZED' / For Biotech Applications

HCON(CH ₃)
M -

M = 73.10 g/mol

1 l = 0.95 kg

FLASHPOINT 58 °C

CAS NO. 68-12-2

EINECS 200-679-5

NC CODE 2924 19 00 EC NO. 616 001 00 0

UN/ID NO. 2265

ADR/RID 3F1

IMDG 3/III

R: 20/21-36-61

S: 45-53

Assay (by GC) (corrected for water)	min. 99.8%	
Amines (as dimethylamines)	max. 5 ppm	
Appearance	passes test	
Color (APHA)	max. 15	
Residue after Evaporation	max. 5 ppm	
Titrable Acid (meq/g)	max. 0.0005	
Titrable Base (meq/g)	max. 0.003	
Water (H ₂ O)	max. 400 ppm	
Ultraviolet Absorbance (1.00-cm path vs water):		

at 270 nm max. 1.00

at 275 nm	max. 0.30
at 295 nm	max. 0.10
at 310 nm	max. 0.05
at 340 nm	max. 0.01
at 400 nm	max. 0.01

PRODUCT	PACKING	CONT
NO.		BO)
9344.4000	4 Glass	

Dimethylformamide can have a teratogenic effect. Women in childbearing age are warned.

Dimethylformamide

'BAKER ANALYZED' / ACS

7032

► HCON(CH₃)₂

M = 73.10 g/mol

1 l = 0.95 kg

FLASHPOINT 58 °C

CAS NO. 68-12-2

EINECS 200-679-5

NC CODE 2924 19 00

EC NO. 616 001 00 0

UN/ID NO. 2265

ADR/RID 3 F1

IMDG 3/III

R: 20/21-36-61

S: 45-53



Assay (by GC)	min. 99.8%
Appearance	clear
Color (APHA)	max. 15
Residue after Evaporation	max. 0.005%
Titrable Acid (meq/g)	max. 0.0005
Titrable Base (meq/g)	max. 0.003
Water (H ₂ O)	max. 0.1%
Trace Impurities (in ppm):	
Aluminium (AI)	max. 0.5
Barium (Ba)	max. 0.1
Boron (B)	max. 0.02
Cadmium (Cd)	max. 0.05
Calcium (Ca)	max. 0.5
Chromium (Cr)	max. 0.02
Cobalt (Co)	max. 0.02
Copper (Cu)	max. 0.02
Iron (Fe)	max. 0.1
Lead (Pb)	max. 0.1
Magnesium (Mg)	max. 0.1
Manganese (Mn)	max. 0.02
Nickel (Ni)	max. 0.02

max. 0.1

max. 0.1

PRODUCT	PACKING	CONT
NO.		B0
7032.1000	11	
7032.2500	2.5	
7032.9025	25 l	
For safe han	dling of 25 l tin cans, see Self-c	losing tap.

Dimethylformamide can have a teratogenic effect. Women in childbearing age are warned.

Dimethylformamide

'BAKER ANALYZED' / Ultraviolet Spectrophotometry / ACS

Tin (Sn)

Zinc (Zn)

Exceeds ACS Specifications

7400

M = 73.10 g/mol **1 l** = 0.95 kg FLASHPOINT 58 °C **CAS NO.** 68-12-2 **EINECS** 200-679-5 NC CODE 2924 19 00 **EC NO.** 616 001 00 0

► HCON(CH₃)₂

UN/ID NO. 2265 ADR/RID 3 F1

IMDG 3/III

R: 20/21-36-61

S: 45-53



Assay (by GC)	mın. 99.8%
Appearance	clear
Color (APHA)	max. 15
Residue after Evaporation	max. 5 ppm
Titrable Acid (meq/g)	max. 0.0005
Titrable Base (meq/g)	max. 0.003
Water (H ₂ O)	max. 0.05%
Trace Impurities (in ppm):	
Iron (Fe)	max. 0.05
Ultraviolet Absorbance (1.00-cm	path vs water):
at 270 nm	max. 1.00
at 275 nm	max. 0.30
at 295 nm	max. 0.10
at 310 nm	max. 0.05
at 340-400 nm	max. 0.01

PRODUCT	PACKING	CONT
NO.		ВОХ
7400.0500	500 ml	ϵ
7400.2500	2.5	

Dimethylformamide can have a teratogenic effect. Women in childbearing age are warned.

Dimet

2,9-Dimethyl-1,10 phenanthroline Hydrochloride

2690 BAKER ANALYZED

► (CH₃)₂C₁₂H₆N₂.HCl.H₂O

M = 262.75 g/mol

CAS NO. 7296-20-1

EINECS 230-732-8 **NC CODE** 2933 99 90

 Sensitivity to Copper
 passes test
 PRODUCT PACKING
 CONT.

 NO.
 BOX

 2690.0001
 1 g



2',2-Dimethyl-5',5-di-iso-propylphenolphthalein

See Thymolphthalein

9234

Dimethyl Sulfoxide

'BAKER BIO-ANALYZED'

M = 78.13 g/mol
1 = 1.10 kg
FLASHPOINT 95 °C
CAS NO. 67-68-5

CAS NO. 67-68-5 EINECS 200-664-3 NC CODE 2930 90 70 R: 36/38

S: 26

Assay ((CH₃)₂SO)(by GC, corrected for water) min. 99.9% Appearance passes test Residue after Evaporation max. 5 ppm Titrable Acid (meq/g) max. 0.001 Water (by KF, coulometric) max. 250 ppm Ultraviolet Absorbance (1.00-cm path vs water):

 Ultraviolet Absorbance (1.00-cm path vs water):

 at 270 nm
 max. 0.40

 at 290 nm
 max. 0.18

 at 310 nm
 max. 0.06

 at 330 nm
 max. 0.02

 at 350-400 nm
 max. 0.01

 PRODUCT NO.
 PACKING PACKING
 CONT.

 9234.4000GL
 4 l Glass

D

7093

Dimethyl Sulfoxide

'BAKER ANALYZED' / Ultraviolet Spectrophotometry

M = 78.13 g/mol 1 I = 1.10 kg FLASHPOINT 95 °C CAS NO. 67-68-5 EINECS 200-664-3 NC CODE 2930 90 70 R: 36/38

S: 26

 Assay
 min. 99.5%

 Color (APHA)
 max. 10

 Residue after Evaporation
 max. 0.005%

 Water (H₂O)
 max. 0.05%

Ultraviolet Absorbance (1.00-cm path vs water; curve smooth throughout stated range with no extraneous impurity peaks):

at 270 nm max. 1.0
at 275 nm max. 0.40
at 280 nm max. 0.30
at 300 nm max. 0.10
at 335 nm max. 0.03
at 400 nm max. 0.03

 PRODUCT
 PACKING
 CONT.

 NO.
 BOX

 7093.1000
 1 l
 6

Dimethyl Sulfoxide

7033 'BAKER ANALYZED

► (CH₃)₂SO

M = 78.13 g/mol **1** l = 1.10 kg

FLASHPOINT 95 °C **CAS NO.** 67-68-5

EINECS 200-664-3 NC CODE 2930 90 70

R: 36/38 **S:** 26

X

Assay	min. 99.0%
Color (APHA)	max. 10
Density (g/ml) at 25°C	min. 1.095
Freezing Point	18.0-19.0°C.
Residue after Evaporation	max. 0.01%
Water (H ₂ O)	max. 0.2%

PRODUCT	PACKING	CONT.
NO.		вох
7033.1000	11	6
7033.5000JC	5 l Jerrycan	
7033.9200	200	

Din See f

Dimethyl Sulfoxide

See for qualities meeting the requirements of cGMP, USP/NF, BP/Ph.Eur, JP or FCC, page 36

CONT.

вох

2,4-Dinitrophenylhydrazine

1293

198-200°C

max. 0.1%

PRODUCT

1293.0100

PRODUCT

8434.1000

NO.

NO.

PACKING

PACKING

11

► (NO₂)₂C₆H₃NHNH₂ **M** = 198.14 g/mol **CAS NO.** 119-26-6 **EINECS** 204-309-3

NC CODE 2928 00 90 **UN/ID NO.** 1325

ADR/RID 4.1 D IMDG 4.1/II

R: 1-22-36/38

S: 35



1,4-Dioxane PHOTREX Reagent / For Spectrophotmetry and Liquid Scintillation Counting / ACS

Melting Point

Residue after Ignition

8434

CONT. вох

QCH₂CH₂OCH₂CH₂ **M** = 88.11 g/mol 11 = 1.03 kg FLASHPOINT 11 °C **CAS NO.** 123-91-1 **EINECS** 204-661-8 NC CODE 2932 99 85 **EC NO.** 603 024 00 5 **UN/ID NO.** 1165 ADR/RID 3 F1 IMDG 3/II **R:** 11-19-36/37-40-66

S: 16-36/37-46





Meets ACS Specifications		
Assay (by GC) (corrected for water)	min. 99.0%	
at 420 nm	max. 0.005	
Carbonyl Compounds (as HCHO) (by		
polarography)	max. 0.01%	
Color (APHA)	max. 10	
Counting Efficiency for ³ H in Prepared		
'Cocktail'	act value reported	
Freezing Point	min. 11.0°C.	
Peroxide (as H ₂ O ₂)	max. 0.005%	
Residue after Evaporation	max. 0.001%	
Titrable Acid (meq/g)	max. 0.0016	
Water (by KF, coulometric)	max. 0.02%	
Ultraviolet Absorbance (1.00-cm path vs water):		
at 215 nm	max. 1.00	
at 225 nm	max. 0.60	
at 250 nm	max. 0.25	
at 270 nm	max. 0.10	
at 280 nm	max. 0.05	
at 295-400 nm	max. 0.01	
Windows of Infrared Transmittane	ce (0.1 mm path,	
60-100%T) , μ m :		
2.0-3.3	info only (not spec)	
3.8-5.0	info only (not spec)	
5.1-6.7	info only (not spec)	
9.9-10.7	info only (not spec)	
12.1-15.4	info only (not spec)	

Use J.T.Baker Ultrex II and BAKER INSTRA-ANALYZED acids for low level trace element analysis.

See chapter 3 of this catalogue for more details.



Dioxa

1,4-Dioxane 8031 BAKER ANALYZED

► QCH₂CH₂OCH₂CH₂

M = 88.11 g/mol

1 l = 1.03 kg

FLASHPOINT 11 $^{\circ}\mathrm{C}$ **CAS NO.** 123-91-1

EINECS 204-661-8

NC CODE 2932 99 85

EC NO. 603 024 00 5

UN/ID NO. 1165

ADR/RID 3F1 IMDG 3/II

R: 11-19-36/37-40-66

S: 16-36/37-46





Assay (by GC)	min. 99.0%
Carbonyl (as HCHO)	max. 0.05%
Color (APHA)	max. 10
Freezing Point	min. 11.0°C
Peroxide (as H ₂ O ₂)	max. 0.01%
Residue after Evaporation	max. 0.005%
Titrable Base (meq/g)	max. 0.0002
Water (H ₂ O)	max. 0.1%
Trace Impurities (in ppm):	
Aluminium (AI)	max. 0.5
Barium (Ba)	max. 0.1
Boron (B)	max. 0.02
Cadmium (Cd)	max. 0.05
Calcium (Ca)	max. 0.5
Chromium (Cr)	max. 0.02
Cobalt (Co)	max. 0.02
Copper (Cu)	max. 0.02
Iron (Fe)	max. 0.1
Lead (Pb)	max. 0.1
Magnesium (Mg)	max. 0.1
Manganese (Mn)	max. 0.02
Nickel (Ni)	max. 0.02
Tin (Sn)	max. 0.1
Zinc (Zn)	max. 0.1

PRODUCT	PACKING	CONT.
NO.		вох
8031.1000	11	6
8031.1000S	1 EcoTainer	
8031.2500	2.5	4
8031.5000	5 EcoTainer	
8031.9200	200 l	

Stabilized.

Diphenylamine 'BAKER ANALYZED' / ACS 1071

► (C₆H₅)₂NH

M = 169.23 g/mol

CAS NO. 122-39-4

EINECS 204-539-4 NC CODE 2921 44 00

EC NO. 612 026 00 5

UN/ID NO. 2811

ADR/RID 6.1 T2

IMDG 6.1/II

R: 23/24/25-33-50/53

S: 28-36/37-45-60-61





dangerous for the environment

Meets ACS Specifications	
Melting Point	52.5-54.0°C
Nitrate (NO ₃)	passes test
Residue after Ignition	max. 0.03%
Sensitivity to Nitrate	passes test
Solubility in Alcohol	passes test

PRODUCT PACKING CONT. вох NO. 1071.0100 100 g

Diphenylcarbazide 'BAKER ANALYZED' / ACS

► C₆H₅NHNHCONHNHC₆H₅ **M** = 242.28 g/mol

CAS NO. 140-22-7

1295

EINECS 205-403-7 NC CODE 2928 00 90

Meets ACS Specifications Melting Point

173-176°C Residue after Ignition max. 0.05% Sensitivity for Chromate passes test Solubility in Aqueous Acetone passes test

PRODUCT PACKING CONT. BOX NO. 1295.0025 25 g Glass

Diphenylcarbazone

1296

► C₆H₅NHNHCON:NC₆H₅ **M** = 240.27 g/mol

CAS NO. 538-62-5

EINECS 208-698-0

NC CODE 2928 00 90

Residue after Ignition	max. 0.05%
Sensitivity for Mercury Detection	passes test
Solubility	passes test

PRODUCT PACKING CONT. NO. BOX 1296.0005 5 g



1,5-Diphenylcarbohydrazide

4,7-Diphenyl-1,10-phenanthroline

1691

 $ightharpoonup (C_6H_5)_2C_{12}H_6N_2$ **M** = 332.41 g/mol **CAS NO.** 1662-01-7 **EINECS** 216-767-1

NC CODE 2933 90 95

PRODUCT PACKING CONT. Sensitivity to Iron passes test NO. вох 1691.0002 2 g

min 00 00/

PRODUCT

0367.5000

PRODUCT

0268.1000

0268.9010

NO.

PACKING

5 kg

Dipotassium Hydrogen Phosphate



Disodium (Ethylenedinitrilo)tetraacetate



Disodium Tetraborate Anhydrous'BAKER ANALYZED' / for X-ray Fluorescence Analysis

0367

CONT.

BOX

► Na₂B₄O₇ **M** = 201.22 g/mol

CAS NO. 1330-43-4 **EINECS** 215-540-4 NC CODE 2840 11 00 **R:** 36/37/38-62-63

S: 22-26-36/37/39-45

Assay (acidimetric)	min. 99.0%
Barium (Ba)	max. 0.002%
Calcium (Ca)	max. 0.005%
Chloride (CI)	max. 0.005%
Chromium (Cr)	max. 0.0002%
Cobalt (Co)	max. 0.0002%
Copper (Cu)	max. 0.0002%
Iron (Fe)	max. 0.001%
Lead (Pb)	max. 0.0002%
Lithium (Li)	max. 0.005%
Loss on ignition (800°C)	max. 1.0%
Magnesium (Mg)	max. 0.001%
Manganese (Mn)	max. 0.0005%
Nickel (Ni)	max. 0.0002%
Phosphate (as PO ₄)	max. 0.002%
Potassium (K)	max. 0.005%
Silicate (as SiO ₂)	max. 0.05%
Strontium (Sr)	max. 0.001%
Sulfate (SO ₄)	max. 0.005%

Disodium Tetraborate Decahydrate

crystal / 'BAKER ANALYZED' / ACS

0268

CONT. BOX

Na₂B₄O₇. 10H₂O **M =** 381.37 g/mol **CAS NO.** 1303-96-4

EINECS 215-540-4 NC CODE 2840 19 90

R: 36/37/38-62-63 **S:** 22-26-36/37/39-45



Meets ACS Specifications	
Assay	99.5-105.0%
Calcium (Ca)	max. 0.005%
Chloride (CI)	max. 0.001%
Heavy Metals (as Pb)	max. 0.001%
Insoluble Matter	max. 0.005%
pH of 0.01 M Solution at 25°C	9.15-9.20
Phosphate (PO ₄)	max. 0.001%
Sulfate (SO ₄)	max. 0.005%
Trace Impurities (in ppm):	
Iron (Fe)	max. 5

Disodium Tetraborate Decahydrate

PACKING

1 kg

10 kg

► Na₂B₄O₇.10H₂O

M = 381.37 g/mol **CAS NO.** 1303-96-4 **EINECS** 215-540-4 NC CODE 2840 19 90 **R:** 22

Assay	99.0-103.0%
Ammonium (NH ₄)	max. 10 ppm
Appearance of solution	passes test
Arsenic (As)	max. 5 ppm
Calcium (Ca)	max. 100 ppm
Heavy Metals (as Pb)	max. 25 ppm
Identification	passes test
рН	9.0 - 9.6
Sulfates (as SO ₄)	max. 50 ppm

PRODUCT	PACKING	CONT.
NO.		вох
0270.1000	1 kg	

Dithiothreitol



0270

Dithi



1,4-Dithiothreitol (DTT, Cleland's Reagent)

e for qualities meeting the requirements of cGMP, USP/NF, BP/Ph.Eur, JP or FCC, page 36



DMSO

See Dimethyl Sulfoxide



DNA/RNA Synthesis reagents

mation section Reagents for DNA/RNA Synthesis, page 261

8172

n-Dodecane

CH ₃ (CH ₂) ₁₀	CH	3
M	1 =	170.

.34 g/mol **1 l** = 0.75 kg FLASHPOINT 74 °C **CAS NO.** 112-40-3 **EINECS** 203-967-9

Assay (by GC) 215-216°C **Boiling Point** Density (g/ml) at 20°C 0.748-0.749 Freezing Point -9 to -10°C

PRODUCT PACKING CONT. NO. BOX 8172.0500 500 ml

L034-07

1-Dodecanol

'BAKER'



M = 186.34 g/mol **1 l =** 0.83 kg FLASHPOINT 107 °C

NC CODE 2901 10 00

CAS NO. 112-53-8 **EINECS** 203-982-0 NC CODE 2905 17 00

min. 98% **PACKING** CONT. Assay (by GC) **PRODUCT** Appearance passes test NO. BOX Water (H₂O) max. 0.1% L034-07 500 ml Glass

min. 99%

Dodecyl Alcohol

4533

Drum opener for 10, 25 and 200 litre solvents cans

PRODUCT PACKING CONT. BOX 4533 1 unit



Edetate Disodium

e for qualities meeting the requirements of cGMP, USP/NF, BP/Ph.Eur, JP or FCC, page 36



Edge bead removers

See for detailed information section Microelectronic materials page 32 or the website www.jtbaker.com/micro

1109

EDTA

Powder / 'BAKER ANALYZED'

► (HOCOCH₂)₂NCH₂CH₂N(CH₂COOH)₂ **M** = 292.25 g/mol

CAS NO. 60-00-4 **EINECS** 200-448-9 NC CODE 2922 49 95

> **R:** 36-52/53 **S:** 61

X

Assay	99.4-100.6%
Heavy Metals (as Pb)	max. 0.001%
Insoluble in Dilute NH ₄ OH	max. 0.005%
Iron (Fe)	max. 0.005%
Residue after Ignition	max. 0.2%

PRODUCT PACKING CONT. BOX NO. 1109.0100 100 g 1109.1000 1 kg



EDTA

See Ethylenediamine Tetraacetic Acid

4040

1073

2060

CONT.

вох



EDTA, Disodium Salt

crystal / 'BAKER ULTRAPURE BIOREAGENT'

► (CH₂)₂N₂(CH₂COOH)₂(CH₂COONa)₂.2H₂O

M = 372.24 g/mol **CAS NO.** 6381-92-6

EINECS 205-358-3

NC CODE 2922 49 95

R: 20/21/22 **S:** 36/37-9

99.5-100.5% Assay Calcium (Ca) max. 0.005% DNase Activity none detected Heavy Metals (as Pb) max. 0.001% Insoluble Matter max. 0.005% Iron (Fe) max. 0.005% Nitrilotriacetic Acid max. 0.1% pH of 5% Solution at 25°C 40-60 Protease Activity none detected

none detected

min 98%

PRODUCT

2060.1000

2060.5000

NO.

PRODUCT	PACKING	CONT.
NO.		вох
4040.0100	100 g	
4040.0500	500 g	
4040.1000	1 kg	
4040.5000	5 kg	

EDTA, Disodium Salt

BAKER ANALYZED

► (CH₂)₂N₂(CH₂COOH)₂(CH₂COONa)₂.2H₂O

M = 372.24 g/mol

CAS NO. 6381-92-6

EINECS 205-358-3 NC CODE 2922 49 95

R: 20/21/22

S: 36/37-9



99.0-101.0% Assav Heavy Metals (as Pb) max.0.001% Insoluble Matter max.0.005% pH of 5% Solution at 25°C 4.0-6.0 Sulfate (SO₄) max. 0.01%

Trace Impurities (in ppm):		
Copper (Cu)	max. 3	
Iron (Fe)	max. 5	

PRODUCT	PACKING	CONT
NO.		ВОХ
1073.0250	250 g	ϵ
1073.0500	500 g	ϵ
1073.1000	1 kg	ϵ

PACKING

1 kg

5 kg

EDTA, Disodium Salt

(CH₂)₂N₂(CH₂COOH)₂(CH₂COONa)₂.2H₂O

M = 372.24 g/mol

CAS NO. 6381-92-6

FINECS 205-358-3

NC CODE 2922 49 95

R: 20/21/22

S: 36/37-9



EDTA, Disodium Salt

See for qualities meeting the requirements of cGMP, USP/NF, BP/Ph.Eur, JP or FCC, page 36

Assay

RNase Activity



EDTA, Disodium Salt

(1ml = 1mg CaCO₃) / standard solution / 'BAKER ANALYZED'

7092

(CH₂)₂N₂(CH₂COOH)₂(CH₂COONa)₂.2H₂O

M = 372.24 g/mol

EINECS 231-791-2

NC CODE 2922 49 95

Heavy Metals (as Pb) max. 5 ppm Iron (Fe) max. 5 ppm 0.98-1.02mg CaCO3/ml Titer

PRODUCT PACKING BOX 11 7092.1000 7092.9020 20 l Polycube

Volumetric Solution, ready for use. The titer found falls within the range 0.98-1.02 mg CaCO₃/ml.

EDTA, Disodium Salt 0.1 mol/l / 'BAKER ANALYZED'

7124

CONT.

вох

EINECS 205-358-3 NC CODE 2922 49 95 Titer (mol/l) 0.0997-0.1003

NO. 7124.1000

PRODUCT

11 5 I HDPE 7124.5000 7124.9020 20 l Polycube

Volumetric Solution, ready for use.

PACKING

Each lot of this product is standardized potentiometrically against NIST tracable reference standard.



EDTA

EDTA, Disodium Salt

olution A, 0.1 mol/l / 'BAKER ANALYZED' / 0.1M solution with Mg- and Zn-complex for the determination of waterhardness. 1 ml represents 5.6°dH/100 ml water.

NC CODE 2922 49 95

7126

PRODUCT PACKING CONT. NO. вох 7126.1000 11

Volumetric Solution, ready for use.

EDTA, Disodium Salt

7481 0.05 mol/l / 'BAKER ANALYZED'

EINECS 231-791-2 NC CODE 2922 49 95 Titer (mol/l)

0.0495-0.0505

PRODUCT PACKING NO. 20 l Polycube

Volumetric Solution, ready for use.

7481.9020

Each lot of this product is standardized potentiometrically against NIST tracable reference standard.

EDTA, Disodium Salt 0.02 mol/I / 'BAKER ANALYZED' 7125

CAS NO. 139-33-3 **EINECS** 205-358-3 NC CODE 2922 49 95

0.0195-0.0205 Titer (mol/l)

PRODUCT PACKING

7125.1000

CONT. вох

6

CONT.

BOX

Volumetric Solution, ready for use.

11

EDTA, Disodium Salt

Solution B, 0.01783 mol/I / 'BAKER ANALYZED' / 0.01783M solution with Mg- and Zn-complex for the determination of waterhardness. 1 ml represents 1°dH/100 ml water.

NC CODE 2922 49 95

7127

PRODUCT PACKING CONT. BOX NO. 7127,1000 11 7127.5000 5 I HDPE 7127.9010 10 l Polycube 7127.9020 20 l Polycube

Volumetric Solution, ready for use

EDTA, Disodium Salt

0.01 mol/l / 'BAKER ANALYZED'

CAS NO. 139-33-3 **EINECS** 205-358-3 NC CODE 2922 49 95

0.0099-0.0101 Molarity (M)

PRODUCT PACKING CONT. NO. BOX 7613.1000 11

Volumetric Solution, ready for use

4871

EDTA, Disodium Salt

Solution A, 0.1 mol/l / DILUT-IT

NC CODE 2922 49 95

PRODUCT PACKING CONT. BOX NO. 1 amp.

Volumetric Concentrate, for dilution to 1 l.

4653

EDTA, Disodium Salt

 $0.1 \, \text{mol/l}$; $^{1}/_{10} \, M = 33.62 g$, 0.1 M / DILUT-IT

(CH₂)₂N₂(CH₂COOH)₂(CH₂COONa)₂

M = 336.21 g/mol **CAS NO.** 139-33-3 **EINECS** 205-358-3 NC CODE 2922 49 95

PRODUCT PACKING CONT. вох NO. 4653 500 ml Sealed Bottle 6

Volumetric Concentrate, for dilution to 1 l.

4652

1883

1303

3874

3871



EDTA, Disodium Salt $0.02 \text{ mol/I}; ^{1}/_{50} M = 6.724g; 0.02M / DILUT-IT$

(CH₂)₂N₂(CH₂COOH)₂(CH₂COONa)₂ PRODUCT PACKING CONT. **M** = 336.21 g/mol NO. вох **CAS NO.** 139-33-3 4652 1 amp. **EINECS** 205-358-3

NC CODE 2922 49 95 Volumetric Concentrate, for dilution to 1 l.

EDTA, Disodium Salt

4870 0.01 mol/l / DILUT-IT EINECS 205-358-3 **PRODUCT PACKING** CONT. NC CODE 2922 49 95 BOX NO. 4870 1 amp.

Volumetric Concentrate, for dilution to 1 l.

EDTA, Iron(III) Sodium Salt

C₁₀H₁₂FeN₂NaO₈ Appearance passes test **PRODUCT PACKING** CONT. **M** = 367.05 g/mol Iron (Fe) 12.5-13.5% BOX **CAS NO.** 15708-41-5 Total EDTA min. 67.5% 1883.0500 500 g **EINECS** 239-802-2

Ehrlich's Reagent

NC CODE 2922 49 90

See p-Dimethylaminobenzaldehyde

Eosin Y (Yellowish)

Ç₆H₄COOÇC₆H-2,4-Br₂-ONaOC₆H-5,7-Br₂-6-ONa Assay min. 88% PRODUCT PACKING CONT. M = 691.88 g/molNO. вох **CAS NO.** 15086-94-9 1303.0025 25 g Glass

EINECS 239-138-3 NC CODE 3204 12 00 C.I. 45380

Eosine-Y (aquaous)

NC CODE 3822 00 00 **PRODUCT PACKING** CONT. NO. BOX 3874.1000 1 | Glass 3874.2500 2.5 | Glass

Eosine-Y-solution

11 = 0.91 kg Eosine Y Staining Solution for use in combination PRODUCT PACKING CONT. FLASHPOINT 29 °C with Hematoxylin NO. вох NC CODE 3822 00 00 3871.1000 1 | Glass **UN/ID NO.** 1170 3871 2500 2.5 | Glass ADR/RID 3 F1 IMDG 3.2/III **R:** 10

Eosin-Methyleneblue dye according Leishman

S: 16-7/9

Eosin-Methyleneblue solution according Leishman

Eosin-Methyleneblue solution according Wright

Eosin

Eosine Solution for HE staining



1,2-Epoxypropane

Eriochrome Black T

► HOC₁₀H₆N:NC₁₀H₄(OH)(NO₂)SO₃Na **M =** 461.40 g/mol

CAS NO. 1787-61-7

EINECS 217-250-3 NC CODE 3204 19 00

UN/ID NO. 3077

ADR/RID 9 M7

IMDG 9/III

R: 36-51/53

S: 26-61







Etching Chemicals

See for detailed information section Microelectronic materials page 32 or the website www.jtbaker.com/micro

Sensitivity as Indicator



Ethanal



1,2-Ethanediol



N,N'-1,2-Ethanediylbis[N-(carboxymethyl)glycine], Disodium salt

Ethanol

Absolute / 'BAKER HPLC ANALYZED' / for use in High Performance Liquid Chromatography

► C₂H₅OH

M = 46.07 g/mol

11 = 0.79 kg

FLASHPOINT 12 °C

CAS NO. 64-17-5

EINECS 200-578-6

NC CODE 2207 10 00

 $\pmb{\mathsf{EC N0.}} \ \ 603\ 002\ 00\ 5$

UN/ID NO. 1170

ADR/RID 3 F1

IMDG 3/II

R: 11 **S:** 16-7



Assay (by GC) min. 99.5% Density (g/ml) at 20°C 0.789 - 0.790 Free Acid (as CH₃COOH) max. 0.001% 8462.1000 Residue after Evaporation max. 0.001% 8462.2500 max. 0.2% Excluding excise.

Ultraviolet Absorbance (1.00-cm path vs water): at 210 nm max. 0.7 at 240 nm at 260 nm max. 0.01

PACKING CONT. BOX 11

PRODUCT

1107.0025

passes test

PACKING

25 g Glass

CONT.

BOX

www.jtbaker.com/europe



Ethanol

Absolute / 'BAKER ANALYZED' / ACS

► C_2H_5OH **M** = 46.07 g/mol **1** I = 0.79 kg

FLASHPOINT 12 °C

CAS NO. 64-17-5 EINECS 200-578-6

NC CODE 2207 10 00

EC NO. 603 002 00 5

UN/ID NO. 1170 **ADR/RID** 3 F1

IMDG 3/II

R: 11 **S:** 16-7

highly

min.99.9%(by volume)
passes test
max. 10
max. 0.1%
max. 0.001%
passes test
passes test
passes test
max. 0.0005
max. 0.0002
max. 0.2%
max. 0.5
max. 0.1
max. 0.02
max. 0.05
max. 0.5
max. 0.02
max. 0.02
max. 0.2
max. 0.1
max. 0.1

PRODUCT	PACKING	CONT
NO.		BOX
8006.0250GL	250 ml Glass	
8006.1000	11	ϵ
8006.1000S	1 EcoTainer	
8006.1000PE	1 HDPE	ϵ
8006.2500	2.5	4
8006.2500PE	2.5 HDPE	4
8006.5000	5 EcoTainer	4
8006.5000PE	5 HDPE	
8006.9025	25	
EcoTainer, the	metal solvent can for more safety i	n the lab.

For safe handling of 25 l tin cans, see Self-closing tap.

Excluding excise.

max. 0.1

max. 0.02 max. 0.02

max. 0.1

max. 0.1

Fthanol

Absolute / 'BAKER ANALYZED' / Ultraviolet Spectrophotometry / ACS

Magnesium (Mg)

Manganese (Mn)

Nickel (Ni) Tin (Sn)

Zinc (Zn)

8098

M = 46.07 g/mol 1 I = 0.79 kg FLASHPOINT 12 °C CAS NO. 64-17-5 EINECS 200-578-6

► C₂H₅OH

EINECS 200-578-6 NC CODE 2207 10 00 EC NO. 603 002 00 5

EC NO. 603 002 00 5

UN/ID NO. 1170 **ADR/RID** 3 F1 **IMDG** 3/II

R: 11

S: 16-7



Meets ACS Specifications		
Assay	min. 99.5%	
Residue after Evaporation	max. 0.001%	
Water (H ₂ O)	max. 0.2%	
Ultraviolet Absorbance (1.00-cm path vs water;		
curve smooth throughout stated range with no		

extraneous impurity peaks):	
at 210 nm	max. 1.00
at 220 nm	max. 0.50
at 230 nm	max. 0.20
at 250 nm	max. 0.05
at 270-400 nm	max 0.01

PRODUCT	PACKING	CONT.
NO.		вох
8098.0500	500 ml	6
8098.2500	2.5	
Excluding ex	cise.	

J.T.Baker HYDRA-POINT Karl Fischer Reagents are fast and accurate.

Refer to the Hydra-Point section of this catalogue .



Ethan

Ethanol

Absolute / 'BAKER

► C ₂ H ₅ OH		Assay	min. 99.5%	PRODUCT	PACKING	CONT.
M =	46.07 g/mol	Color (APHA)	max. 10	NO.		вох
1 l =	0.79 kg			8228.1000	11	6
FLASHPOINT	12 °C			8228.2500	2.5	4
CAS NO.	64-17-5			8228.9025	25 l	
EINECS	200-578-6			8228.9200	200 l	
NC CODE	2207 10 00			Excluding exc	cise.	
EC NO.	603 002 00 5			For safe hand	dling of 25 l tin cans, see Self-closing t	ap.
UN/ID NO.	1170					
ADR/RID	3 F1					
IMDG	3/II					
R:	11					
S:	16-7					
	F highly					

3406	Ethanol 99.8% / HISTO GRADE					
► C ₂ H ₅ OH		Assay	min. 99.8%	PRODUCT	PACKING	CONT.
M =	46.07 g/mol	Calcium (Ca)	max. 5 ppm	NO.		вох
11=	0.79 kg	Iron (Fe)	max. 1 ppm	3406.2500PE	2.5 HDPE	
FLASHPOINT	12 °C	Residue after Evaporation	max. 0.0015%	3406.5000	5 l Jerrycan	
CAS NO.	64-17-5	Substances Darkened by H ₂ SO ₄	passes test	3406.9010	10 l Jerrycan	
EINECS	200-578-6	Titrable Acid (meq/g)	max. 0.001	3406.9025	25 l Jerrycan	
NC CODE	2207 10 00	Water (H ₂ O)	max. 0.2%	Excluding exci	se.	
EC NO.	603 002 00 5					
UN/ID NO.	1170			Histo-Grade i	mplicates that this reagent is spe	cially tested
ADR/RID	3 F1			and therefore	solely intended for use in histo-p	pathology
IMDG	3/II			applications.	This reagent is of an analytical qu	uality.
R:	11					
S:	16-7					
	highly flammable					

0

P Q R S T U

Ethanol, absolute MOS, VLSI Grade

Calibrate and standardise your analytical methods and equipment with J.T.Baker Volumetric and Buffer solutions.

Refer to the Analytical applications section of this catalogue for more details.



Ethanol

95% / 'BAKER ANALYZED' / ACS

	95% / 'BAKER ANALYZED' / ACS					8007
► C ₂ H ₅ OH		Exceeds ACS Specifications		PRODUCT	PACKING	CONT.
M =	46.07 g/mol	Assay	min. 95.0% by volume	NO.		вох
1 i =	0.81 kg	Acetone, isopropylalcohol	passes test	8007.1000	11	6
FLASHPOINT	12 °C	Color (APHA)	max. 10	8007.1000S	1 EcoTainer	
CAS NO.	64-17-5	Methanol (CH ₃ OH)	max. 0.1%	8007.2500	2.5	4
EINECS	200-578-6	Residue after Evaporation	max. 0.001%	8007.2500PE	2.5 HDPE	4
NC CODE	2207 10 00	Solubility in Water	passes test	8007.5000	5 EcoTainer	4
EC NO.	603 002 00 5	Substances Darkened by H ₂ SO ₄	passes test	8007.9025	25 l	
UN/ID NO.	1170	Substances Reducing KMnO ₄	passes test	8007.9200	200 l	
ADR/RID	3 F1	Titrable Acid (meq/g)	max. 0.0005	EcoTainer, the	metal solvent can for r	more safety in the lab.
IMDG	3/II	Titrable Base (meq/g)	max. 0.0002	Excluding exc	ise.	
R:	11	Trace Impurities (in ppm):		For safe hand	ling of 25 l tin cans, see	Self-closing tap.
S:	16-7	Aluminium (AI)	max. 0.5			
	 F	Barium (Ba)	max. 0.1			
	10	Boron (B)	max. 0.02			
	highly flammable	Cadmium (Cd)	max. 0.05			
	nammable	Calcium (Ca)	max. 0.5			
		Chromium (Cr)	max. 0.02			
		Cobalt (Co)	max. 0.02			
		Copper (Cu)	max. 0.02			
		Iron (Fe)	max. 0.1			
		Lead (Pb)	max. 0.1			
		Magnesium (Mg)	max. 0.1			
		Manganese (Mn)	max. 0.02			
		Nickel (Ni)	max. 0.02			
		Tin (Sn)	max. 0.1			

Ethanol
96% / 'BAKER'
8229

max. 0.1

	30 % / BAKEN
► C ₂ H ₅ OH	
M =	46.07 g/mol
1 =	0.81 kg
FLASHPOINT	12 °C
CAS NO.	64-17-5
EINECS	200-578-6
NC CODE	2207 10 00
EC NO.	603 002 00 5
UN/ID NO.	1170
ADR/RID	3 F1
IMDG	3/11
R:	11
S:	16-7
	F highly
	flammable

Assay (96.0-97.2%)	passes test
Appearance	passes test
Boiling Point	78.0-79.0°C
Density (g/ml) at 20°C	0.804-0.809
Furfural	passes test
Fuseloil, aldehyde	passes test
Heavy Metals, Zinc (as Pb)	max. 2 ppm
Identification (by IR)	passes test
Iron (Fe)	max. 1 ppm
Methanol (CH ₃ OH)	max. 500 ppm
Residue after Evaporation	max. 0.0015%
Substances Reducing KMnO ₄	passes test
Titrable acid	passes test
Titrable base	passes test
UV-absorbance (1cm vs hexane):	
220-270 nm (no maxima)	passes test
at 220 nm	max. 0.30
at 230 nm	max. 0.18
at 240 nm	max. 0.08
at 270 nm	max. 0.02

PRODUCT PA NO. 8229.1000 1 I 8229.2500 2.5	ACKING CONT BO
8229.1000 1 l	
8229.2500 2.5	
	51
8229.5000 5 l	EcoTainer
8229.9025 25	I .
EcoTainer, the me	etal solvent can for more safety in the lab
Excluding excise.	
For safe handling	of 25 l tin cans, see Self-closing tap.

Mallinckrodt Baker's chemistry is Part of a pure process TM .

Zinc (Zn)

Ethan

Ethanol

3405 96% / HISTO GRADE

► C₂H₅OH

M = 46.07 g/mol

1 l = 0.81 kg

 $\textbf{FLASHPOINT} \quad 12~^{\circ}\text{C}$

CAS NO. 64-17-5

EINECS 200-578-6

NC CODE 2207 10 00

EC NO. 603 002 00 5

UN/ID NO. 1170

ADR/RID 3 F1

IMDG 3/II

R: 11

S: 16-7

highly

Assay	min. 96%
Calcium (Ca)	max. 5 ppm
Iron (Fe)	max. 1 ppm
Residue after Evaporation	max. 0.0015%
Substances Darkened by H ₂ SO ₄	passes test
Titrable Acid (meq/g)	max. 0.001

PRODUCT NO.	PACKING	CONT. BOX
3405.2500PE	2.5 I HDPE	
3405.5000	5 Jerrycan	
3405.9010	10 l Jerrycan	
3405.9025	25 l Jerrycan	
Excluding exci	se.	

Histo-Grade implicates that this reagent is specially tested and therefore solely intended for use in histo-pathology applications. This reagent is of an analytical quality.

3408

Ethanol 99%, denatured

408 HISTO GRAD

► C₂H₅OH

M = 46.07 g/mol

1 l = 0.79 kg

FLASHPOINT 12 °C

LASHPUINI 12°C

CAS NO. 64-17-5 **EINECS** 200-578-6

NC CODE 2207 10 00

EC NO. 603 002 00 5

UN/ID NO. 1170

ADR/RID 3 F1

IMDG 3/II

R: 11-20/21/22-68/20/21/22

S: 16-23-36/37-45-7/9





 Assay (exclusive about 5% methanol)
 min. 99%

 Methanol (CH₃OH)
 4-6%

 Water (H₂O)
 max. 0.5%

Denatured with 5% methanol.

Histo-Grade implicates that this reagent is specially tested and therefore solely intended for use in histo-pathology applications. This reagent is of an analytical quality.

Ethanol 96%, denatured

3407 HISTO GRADE

► C₂H₅OH

M = 46.07 g/mol

1 l = 0.79 kg

FLASHPOINT 12 °C **CAS NO.** 64-17-5

EINECS 200-578-6

NC CODE 2207 10 00

EC NO. 603 002 00 5

UN/ID NO. 1170

ADR/RID 3 F1

IMDG 3/II

R: 11-20/21/22-68/20/21/22

S: 16-23-36/37-45-7/9





 $\begin{array}{lll} \mbox{Assay (exclusive about 5\% methanol)} & \mbox{min. 96\%} \\ \mbox{Methanol (CH}_3\mbox{OH)} & \mbox{4-6\%} \\ \mbox{Water (H}_2\mbox{O)} & \mbox{4-6\%} \end{array}$

Denatured with 5% methanol.

Histo-Grade implicates that this reagent is specially tested and therefore solely intended for use in histo-pathology applications. This reagent is of an analytical quality.

Ethanol 70% Solution Denatured, Sterile

See for qualities meeting the requirements of cGMP, USP/NF, BP/Ph.Eur, JP or FCC, page 36

вох

Ethanolamine

7055 BAKER ANALYZED

► HOCH₂CH₂NH₂

M = 61.08 g/mol

1 l = 1.01 kg

 $\textbf{FLASHPOINT} \;\; 85 \, ^{\circ}\text{C}$

CAS NO. 141-43-5 **EINECS** 205-483-3

NC CODE 2922 11 00

EC NO. 603 030 00 8


UN/ID NO. 2491

ADR/RID 8 C7

IMDG 8/III **R:** 20-36/37/38



min. 99.0% Assay 1.011-1.015 Density (g/ml) at 25°C Heavy Metals (as Pb) max. 0.002% max. 0.001% Iron (Fe) max. 0.005% Residue after Ignition



Ether

See Diethyl Ether



Ether

See for qualities meeting the requirements of cGMP, USP/NF, BP/Ph.Eur, JP or FCC, page 36



Ethidium Bromide



'BAKER ULTRAPURE BIOREAGENT'



.31 g/mol
9-45-8
-984-6
3 99 90

UN/ID NO. 2811 **ADR/RID** 6.1 T2 IMDG 6.1/I

R: 21/22-26-36/37/38-40 **S:** 22-26-36/37-43A-45



For Nucleic Acid Detection	
Assay (dried basis)	min. 98%
DNase Activity	none detected
Protease Activity	none detected
RNase Activity	none detected
UV Absorbance at Maximum	min. 1.25
UV Absorbance, nm	283-287
Visible Absorbance at Maximum	min. 0.65
Visible Absorbance, nm	475-481

PRODUCT	PACKING	CONT.
NO.		вох
4007.0001	1 g	
4007.0005	5 g	

2-Ethoxyethanol

8083 'BAKER ANALYZED'



M = 90.12 g/mol

1 I = 0.93 kg FLASHPOINT 40 °C

CAS NO. 110-80-5

EINECS 203-804-1

NC CODE 2909 44 00

EC NO. 603 012 00 0

UN/ID NO. 1171

ADR/RID 3 F1

IMDG 3/III

R: 10-20/21/22-60-61

S: 45-53



Assay (by GC)	min. 98%	PRODUCT	PACKING
Boiling Range (95%)	133-137°C	NO.	
Density (g/ml) at 20°C	0.929-0.930	8083.1000	11
Residue after Evaporation	max. 0.01%		
Water (H ₂ O)	max. 0.1%		



Ethyl

Ethyl Acetate

'BAKER ULTRA RESI-ANALYZED' / for Organic Residue Analysis 9260

► CH₃COOC₂H₅

M = 88.11 g/mol

1 l = 0.90 kg

 $\textbf{FLASHPOINT} \quad -4~^{\circ}\text{C}$

CAS NO. 141-78-6 **EINECS** 205-500-4

NC CODE 2915 31 00

EC NO. 607 022 00 5

UN/ID NO. 1173

ADR/RID 3F1

IMDG 3/II

R: 11-36-66-67

S: 16-26-33





Assay (by GC) (corrected for water)	min. 99.6%
Color (APHA)	max. 10
Residue after Evaporation	max. 1 ppm
Substances Darkened by H ₂ SO ₄	passes test
Titrable Acid (meq/g)	max. 0.0008
Water (H ₂ O)	max. 0.05%

ECD Sensitive Impurities (as Heptachlor Epoxide): max. 10

Single Impurity Peak (pg/ml) FID-Sensitive Impurities (as 2-Octanol):

Single Impurity Peak (ng/ml) max. 5

PRODUCT	PACKING	CONT.
NO.		вох
9260.1000	11	6
9260.2500	2.5	4

9831

Ethyl Acetate

BAKER ANALYZED LC-MS Reagent

► CH₃COOC₂H₅ **M** = 88.11 g/mol **1 l** = 0.90 kg FLASHPOINT $-4\,^{\circ}\text{C}$ **CAS NO.** 141-78-6

EINECS 205-500-4 NC CODE 2915 31 00 EC NO. 607 022 00 5

UN/ID NO. 1173 ADR/RID 3F1

> IMDG 3/II **R:** 11-36-66-67





Certificate Provided Reporting Assay (by GC)	min. 99.6%
Residue after Evaporation	max. 1 ppm
Water (H ₂ O)	max. 0.05%
Product Information (not speci	fications):
Density (g/ml) at 20°C	0.90
Trace Impurities (in ppb):	
Aluminium (AI)	max. 50
Calcium (Ca)	max. 50
Iron (Fe)	max. 50
Magnesium (Mg)	max. 50
Potassium (K)	max. 50
Sodium (Na)	max. 50
Ultraviolet Absorbance (1.00-ci	n path vs water):
at 265 nm	max. 0.05
at 280 nm	max. 0.02

PRODUCT PACKING CONT. BOX 9831.1000GL 1 | Glass 6

Element concentrations are at time of lot release.

Ethyl Acetate

9282 'BAKER HPLC ANALYZED' / for Use in Liquid Chromatography and Spectrophotometry

UV Cut-off, nm

► CH₃COOC₂H₅ **M** = 88.11 g/mol

> 11 = 0.90 kg FLASHPOINT $-4\,^{\circ}\mathrm{C}$ **CAS NO.** 141-78-6

EINECS 205-500-4 NC CODE 2915 31 00

EC NO. 607 022 00 5 **UN/ID NO.** 1173

ADR/RID 3 F1 IMDG 3/II

R: 11-36-66-67

S: 16-26-33





Residue after Evaporation	max. 2 ppm
Titrable Acid (µeq/g)	max. 0.8
Water (by KF, volumetric)	max. 0.04%
Fluorescence Trace Impurities (as quinine base),	
ppb:	
Measured at 450 nm	max. 0.3
Measured at Emission Maximum	max. 1.0
Ultraviolet Absorbance (1.00-cm pat	h vs water):
at 265 nm	max. 0.05
at 280 nm	max. 0.02
at 330-400 nm	max. 0.01

min. 99.6%

max. 255

Assay (by GC) (corrected for water)

PRODUCT NO.	PACKING	CONT. BOX
9282.1000	11	6
9282.4000	4 l Glass	4

Filtered through a 0.2 micron filter. Packaged under Nitrogen.

Ethyl Acetate

'BAKER INSTRA-ANALYZED' / GC-Spectrophotometric quality / ACS

8204

► CH₃COOC₂H₅

M = 88.11 g/mol

1 l = 0.90 kg

 $\textbf{FLASHPOINT} \quad -4~^{\circ}\text{C}$

CAS NO. 141-78-6

EINECS 205-500-4

NC CODE 2915 31 00

EC NO. 607 022 00 5

UN/ID NO. 1173

ADR/RID 3 F1

IMDG 3/II R: 11-36-66-67

S: 16-26-33





Meets ACS Specifications	
Assay (by GC)	min. 99.5%
Color (APHA)	max. 10
Residue after Evaporation	max. 0.003%
Substances Darkened by H ₂ SO ₄	passes test
Titrable Acid (meq/g)	max. 0.0009
Water (H ₂ O)	max. 0.2%

Ultraviolet Absorbance (1.00-cm path vs water; curve smooth throughout stated range with no extraneous impurity peaks):

at 255 nm	max. 1.00
at 257 nm	max. 0.50
at 263 nm	max. 0.10
at 275 nm	max. 0.05
at 330-400 nm	max. 0.01

PRODUCT	PACKING	CONT.
NO.		вох
8204.1000	1	

Ethyl Acetate

BAKER BIO-ANALYZED

9276

► CH₃COOC₂H₅

M = 88.11 g/mol **1 l** = 0.90 kg

FLASHPOINT $-4\,^{\circ}\text{C}$

CAS NO. 141-78-6

EINECS 205-500-4

NC CODE 2915 31 00

EC NO. 607 022 00 5

UN/ID NO. 1173

ADR/RID 3 F1 IMDG 3/II

R: 11-36-66-67

S: 16-26-33





min. 99.6% Assay (by GC) (corrected for water) Residue after Evaporation max. 1 ppm Titrable Acid (meq/g) max. 0.0008 Water (H₂O) max. 300 ppm

Fluorescence Trace Impurities (as quinine base), ppb:

at 450 nm Emission max. 0.3 at Emission Maximum for Impurities max. 1.0

Ultraviolet Absorbance (1.00-cm path vs water):

at 265 nm	max. 0.05
at 280 nm	max. 0.02
at 330-400 nm	max. 0.01
UV Cut-off, nm	max. 255

PRODUCT	PACKING	CONT.
NO.		вох
9276 1000GL	1 Glass	

Ethyl Acetate

'BAKER ANALYZED' / ACS

8037

► CH₃COOC₂H₅

M = 88.11 g/mol

11 = 0.90 kg

FLASHPOINT $-4\,^{\circ}\mathrm{C}$

CAS NO. 141-78-6

EINECS 205-500-4 NC CODE 2915 31 00

EC NO. 607 022 00 5

UN/ID NO. 1173

ADR/RID 3 F1

IMDG 3/II

R: 11-36-66-67 **S:** 16-26-33





Assay (by GC)

Color (APHA)

Residue after Evaporation max. 0.003% Substances Darkened by H₂SO₄ passes test k. 0.0009

min. 99.5%

max. 10

Exceeds ACS Specifications

Titrable Acid (meq/g)	max. 0.0009
Water (H ₂ O)	max. 0.2%
Trace Impurities (in ppm):	
Aluminium (AI)	max. 0.5
Barium (Ba)	max. 0.1
Boron (B)	max. 0.02
Cadmium (Cd)	max. 0.05
Calcium (Ca)	max. 0.5
Chromium (Cr)	max. 0.02
Cobalt (Co)	max. 0.02
Copper (Cu)	max. 0.02
Iron (Fe)	max. 0.1
Lead (Pb)	max. 0.1
Magnesium (Mg)	max. 0.1
Manganese (Mn)	max. 0.02
Nickel (Ni)	max. 0.02
Tin (Sn)	max. 0.1
Zinc (Zn)	max. 0.1

PRODUCT	PACKING	CONT.
NO.		вох
8037.1000	11	6
8037.2500	2.5	4
8037.2500PE	2.5 I HDPE	4
8037.5000	5 EcoTainer	
8037.9025	25 l	4
8037.9200	200	

EcoTainer, the metal solvent can for more safety in the lab. For safe handling of 25 l tin cans, see Self-closing tap.

Ethyl

Ethyl Acetate

8038 'BAI

► CH₃COOC₂H₅

M = 88.11 g/mol

1 l = 0.90 kg

 $\textbf{FLASHPOINT} \quad -4~^{\circ}\text{C}$

CAS NO. 141-78-6

EINECS 205-500-4

NC CODE 2915 31 00

EC NO. 607 022 00 5

UN/ID NO. 1173

ADR/RID 3 F1

IMDG 3/II

R: 11-36-66-67 **S:** 16-26-33







Ethyl Acetate

ee for qualities meeting the requirements of cGMP, USP/NF, BP/Ph.Eur, JP or FCC, page 36

Assay (by GC)

Residue after Evaporation

Appearance

Water (H₂O)



Ethylene Chloride

See 1,2-Dichloroethane

8409

Ethylenediamine

BAKER'

► NH₂CH₂CH₂NH₂

M = 60.10 g/mol

1 l = 0.90 kg

FLASHPOINT 34 $^{\circ}\mathrm{C}$

CAS NO. 107-15-3

EINECS 203-468-6

NC CODE 2921 21 00

EC NO. 612 006 00 6

UN/ID NO. 1604

ADR/RID 8 CF1

IMDG 8/II

R: 10-21/22-34-42/43

S: 23-26-36/37/39-45





min. 99%

passes test

max. 0.003%

max. 0.1%

PRODUCT

8038.9010

8038.9025

8038.9200

NO.

PACKING

10 l

25 I

200 l

For safe handling of 25 l tin cans, see Self-closing tap.

CONT.

вох

PRODUCT	PACKING	CONT.
NO.		вох
8409.1000	11	
8409.9025	25 l	
For safe han	dling of 25 l tin cans, see Self-clos	sing tap.

Fth

Ethylenediamine Tetraacetic Acid

See EDTA

Ethylenediamine Tetraacetic Acid, Disodium Salt

Ethylenediamine Tetraacetic Acid, Iron(III) Sodium Salt

See EDTA, Iron(III) Sodium Salt



Ethylene Dibromide

see 1,2-Dibromoethane



Ethylene Dichloride

See for qualities meeting the requirements of cGMP, USP/NF, BP/Ph.Eur, JP or FCC, page 36

7037

CONT.

вох

6

4



Ethylene Glycol 'BAKER ANALYZED'

► HOCH₂CH₂OH

M = 62.07 g/mol **1 l** = 1.11 kg

FLASHPOINT 111 °C

CAS NO. 107-21-1

EINECS 203-473-3

NC CODE 2905 31 00 EC NO. 603 027 00 1

R: 22



Assay (by GC) min. 99.0% Acidity (as CH₃COOH) max. 0.01% **Boiling Point** 195-198°C 194.0-200.0°C Boiling Range Color (APHA) max. 10 Density (g/ml) at 20°C 1.110-1.114 Residue after Ignition max. 0.005% Water (H₂O) max. 0.2% Trace Impurities (in ppm):

PRODUCT PACKING NO. 7037.1000 11 7037.2500 2.5 l 7037.9010 10 l 7037.9025 25 I 7037.9200 200 l For safe handling of 25 l tin cans, see Self-closing tap.

max. 5

max. 0.2

130-131°C

value on label

value on label

value on label

passes test

PRODUCT

8782.0250

NO.

Ethylene Glycol MOS Grade

Chloride (CI)

Iron (Fe)



Ethylene Glycol Monoethyl Ether



Ethylene Glycol Monomethyl Ether

See 2-Methoxyethanol



Ethylene Tetrachloride

See Tetrachloroethylene



Ethylene Trichloride



Ethyl Malonate

See Diethyl Malonate



Ethyl Methyl Ketone



Ethylparaben

See for qualities meeting the requirements of cGMP, USP/NF, BP/Ph.Eur, JP or FCC, page 36

Boiling Point



1-Ethylpiperidine

'BAKER



CONT.

вох

C₂H₅N(CH₂)₄CH₂ **M** = 113.20 g/mol

11 = 0.82 kg

FLASHPOINT 19 °C

CAS NO. 766-09-6

EINECS 212-161-6

NC CODE 2933 39 99

UN/ID NO. 2386

ADR/RID 3 FC

IMDG 3/II

R: 11-20/22-34

S: 16-26-36/37/39-45









Fast Green FCF 'BAKER ANALYZED

1476

C₃₇H₃₄N₂Na₂O₁₀S₃

M = 808.86 g/mol

Absorbance at maximum (0.6 mg/200 ml in 50% C₂H₅OH, 1-cm path) **CAS NO.** 2353-45-9 Absorbance Maximum, nm NC CODE 3204 19 00 Biological Test Total Dye Content

PRODUCI	PACKING	CONT
NO.		ВОХ
1476.0025	25 g Glass	
	10.1	
For Histolog	y and Cytology.	

PACKING

250 ml



Fast



Fast staining in Hematology

Fehling's Reagent I

7478

1 l = 1.05 kg CAS NO. 7758-98-7

EINECS 231-847-6 NC CODE 2833 25 00

R: 51/53

S: 57

N N dangerous for the Suitability for determination of reducing

passes test

7478.0100 100 ml

PRODUCT

7478.0500 500 ml 7478.2500 2.5 l

PACKING

CONT.

BOX

6

Copper(II)Sulfate solution.

7479

Fehling's Reagent II

environment

11 = 1.24 kg **CAS NO.** 6381-59-5 **EINECS** 205-698-2


NC CODE 2918 13 00 UN/ID NO. 3266

ADR/RID 8 C5 IMDG 8/II

R: 35

S: 20-26-36/37/39-45





passes test

PRODUCT PACKING CONT. вох 7479.0100 100 ml 7479.0500 500 ml 7479.2500 251

Potassium sodium tartrate solution, alkaline.



Ferric Ammonium Citrate, Brown

peting the requirements of cGMP, USP/NF, BP/Ph.Eur, JP or FCC, page 36



Ferric Ammonium Citrate, Green

See for qualities meeting the requirements of cGMP, USP/NF, BP/Ph.Eur, JP or FCC, page 36



Ferric Ammonium Sulfate

See Ammonium Iron(III) Sulfate Dodecahydrate



Ferric Subsulfate Solutions

ee for qualities meeting the requirements of cGMP, USP/NF, BP/Ph.Eur, JP or FCC, page 36

7333

Ferroin Solution 1/40 M

'BAKER ANALYZED'

C₃₆H₂₄FeN₆.O₄S **1 l** = 1.00 kg **CAS NO.** 14634-91-4 **EINECS** 238-676-6 NC CODE 3822 00 00

Sensitivity as a redox indicator

passes test

PRODUCT PACKING 7333.0100 100 ml

500 ml

CONT. BOX

Redox Indicator.

7333.0500



Ferrous Ammonium Sulfate



Ferrous Fumarate

See for qualities meeting the requirements of cGMP, USP/NF, BP/Ph.Eur, JP or FCC, page 36



Ferrous Sulfate

See for qualities meeting the requirements of cGMP, USP/NF, BP/Ph.Eur, JP or FCC, page 36



Ferrous Sulfate, 7-Hydrate

See for qualities meeting the requirements of cGMP, USP/NF, BP/Ph.Eur, JP or FCC, page 36

Finyte, Finyte-1 chemicals



Fixatives

See Bouin's Fixative or Formalin solution



Florisil

Activated at 675°C / 60-100 mesh / 'BAKER ANALYZED'

3369

NC CODE 2839 90 00	Suitable for use in chromatogrape pesticide residues	hic cleanup of
	Mesh:	
	On U.S. No. 60 Sieve	max. 10%
	Thru U.S. No. 100 Sieve	max. 10%
	Physical Data (not specifications)	:

Bulk Density (g/cc)

PRODUCT	PACKING	CONT.
NO.		вох
3369.0500	500 g	
3369.2000	2 kg	

Florisil

60-100 mesh / 'BAKER ANALYZED' / For Hydrocarbon Oil Index determination according to ISO 9377-2 and NEN 5733

Average Particle Diameter, µm (APD)

7061-00

	1343-88-0 215-681-1	Florisil recovery test with stearyl-stearate/mineral oil standard acc. to		PRODUCT	PACKING	CONT. BOX
NC CODE	2839 90 00	ISO 9377-2	passes test	7061-00	100 g	
		Physical Data (not specifications):				
		Average Particle Diameter, µm (APD)	140-250	Activated ma	agnesium silicate.	
		Bulk Density (g/cc)(typical)	0.5	After opening	g: Store in desiccator.	

140-250

0.5

Florisil

60-100 mesh / 'BAKER' / Gas Chromatographic Adsorbent

0567

CAS NO. 1343-88-0	Physical Data (not specifications):		PRODUCT	PACKING	CONT.
EINECS 215-681-1	Average Particle Diameter, μm (APD)	140-250	NO.		вох
NC CODE 2839 90 00	Bulk Density (g/cc)(typical)	0.5	0567.0250	250 g	6

Activated magnesium silicate.

Fluorescein Disodium Salt

1314

► 2-NaOCOC ₆ H ₄ C:C ₆ H ₃ -3(:O)OC ₆ H ₃ -6-ONa	Appearance	passes test	PRODUCT	PACKING	CONT.
M = 376.28 g/mol	Appearance of solution	passes test	NO.		вох
CAS NO. 518-47-8			1314.0025	25 g Glass	6
EINECS 208-253-0					
NC CODE 3204 12 00			C.L. 45350		

Formaldehyde solution 38% (w/v) / HISTO GRADE / Fixative for histology

3859

► CH ₂ O	
M =	30.03 g/mol
CAS NO.	50-00-0
EINECS	200-001-8
NC CODE	2912 11 00
EC NO.	605 001 00 5
UN/ID NO.	2209
ADR/RID	8 C9
IMDG	8/III
R:	23/24/25-34-40-43
S:	26-36/37-45-51

Formaldhyde (HCHO)(after dilution)	36 - 40 g/l
pH 25°C (after dilution)	6.8 - 7.2

PRODUCT	PACKING	CONT.
NO.		вох
3859.5000	5 I HDPE	
3859.9010	10 l Jerrycan	

Store at room temperature (18-30°C). Cooling down formaldehyde solution below 18°C may cause irreversible polymerisation.



Forma

Formaldehyde solution

7040 / 'BAKER ANALYZED' / stabilized with 10-15% methanol

Assay

Color (APHA)

Sulfate (SO₄)

Preservative (CH₃OH)

Residue after Ignition

► HCHO **M** = 30.03 g/mol **1 l** = 1.09 kg $\textbf{FLASHPOINT} \ \ 60\ ^{\circ}\text{C}$ **CAS NO.** 50-00-0 **EINECS** 200-001-8 NC CODE 2912 11 00 EC NO. 605 001 00 5 **UN/ID NO.** 2209

Titrable Acid (meq/g) Trace Impurities (in ppm): Chloride (CI) Heavy Metals (as Pb) ADR/RID 8 C9 Iron (Fe) IMDG 8/III R: 23/24/25-34-40-43 **S:** 26-36/37-45-51

36.5-38.0%	PRODUCT	PACKING	CONT.
max. 10	NO.		вох
10-15%	7040.1000	11	6
max. 0.01%	7040.2500	2.5	4
max. 0.002%	7040.5000	5 I HDPE	
max. 0.006			

max. 5

max. 5

max. 5

Formaldehyde solution

7041 37% / 'BAKER' / stabilized with 9-15% methanol

► HCHO **M** = 30.03 g/mol 11 = 1.09 kg FLASHPOINT 60 °C **CAS NO.** 50-00-0 EINECS 200-001-8 NC CODE 2912 11 00 **EC NO.** 605 001 00 5 **UN/ID NO.** 2209 ADR/RID 8 C9

IMDG 8/III R: 23/24/25-34-40-43 **S:** 26-36/37-45-51

37.0-38.0% Assay Acidity passes test Appearance of solution passes test Identification passes test Methanol (CH₃OH) 9.0-15.0% (v/v) Sulfated Ash max. 0.1%

PRODUCT PACKING CONT. NO. BOX 7041.1000 11 7041.2500 2.5 l 7041.2500PE 2.5 | HDPE 7041.9010 10 I 7041.9025 25 I 200 l 7041.9200 For safe handling of 25 l tin cans, see Self-closing tap.

Stored at a temperature of 15°C to 25°C. Stored in an airtight container. Stored protected from light.

PACKING

PRODUCT

Formaldehyde solution 3858

19% (w/v) / HISTO GRADE / Fixative for histology

F

► HCHO M = 30.03 g/mol**1 l** = 1.04 kg **CAS NO.** 50-00-0

EINECS 200-001-8 NC CODE 2912 11 00 **EC NO.** 605 001 00 5

R: 20/21/22-36/37/38-40-43

S: 26-36/37/39-51



Formaldhyde (HCHO)(after dilution)	38 - 42 g/I
pH 25°C (after dilution)	6.8 - 7.2
Sodium (Na)	61 - 71 mmol/l

NO. 3858.9020 20 L jerrycan

CONT.

BOX

Store at room temperature (18-30°C). Cooling down formaldehyde solution below 18°C may cause irreversible polymerisation.

Formaldehyde Solutions

eting the requirements of cGMP, USP/NF, BP/Ph.Eur, JP or FCC, page 36

Formalin Neutralizer HISTO GRADE

EINECS 231-673-0 NC CODE 2832 10 00

R: 22-31-37-41 **S:** 26-39



Contains: Sodium Disulfite min. 80% PRODUCT PACKING CONT. NO. BOX 3934.5000 5 kg Plastic Pail

Tested by an independent Sterlab Qualified environmental laboratory.

Formalin solution

13% (v/v) / 'BAKER' / For Histology

► HCHO

M = 30.03 g/mol **CAS NO.** 50-00-0

EINECS 200-001-8 NC CODE 2912 11 00

EC NO. 605 001 00 5

R: 20/21/22-36/37/38-40-43

S: 26-36/37/39-51



Formaldehyde (HCHO) pH at 25°C

49-53 g/l 6.8-7.2

9.7-10.3 %

PRODUCT PACKING NO. 7390.9010 10 l Polycube

(4.0 g/l sodium dihydrogen phosphate, 6.5 g/l sodium monohydrogen phosphate and 130 ml/l formaldehyde sol. (37-38%) in water).

Formalin Solution

10% (w/v) neutralized / 'BAKER ANALYZED' / stabilized with methanol

Assay

7424

7390

CONT.

вох

► HCHO

M = 30.03 g/mol

1 l = 1.04 kg **CAS NO.** 50-00-0

EINECS 200-001-8 NC CODE 2912 11 00

EC NO. 605 001 00 5 R: 20/21/22-36/37/38-40-43

S: 26-36/37/39-51



PRODUCT	PACKING	CONT.
NO.		вох
7424.1000	11	6

Store at room temperature (18-30°C).

Cooling down formaldehyde solution below 18°C may cause irreversible polymerisation.

Formalin Solution

10% (v/v) [=4% (w/v)] neutralized / HISTO GRADE / Fixative for histology

3933

CONT.

BOX

► HCHO

M = 30.03 g/mol

CAS NO. 50-00-0

EINECS 200-001-8

NC CODE 2912 11 00 EC NO. 605 001 00 5

R: 40-43

S: 36/37



38-42 g/l Formaldehyde (HCHO) pH at 25°C 6.8-7.2 Sodium (Na) 61-71 mmol/l

NO. 3933.1000

PRODUCT

3933.9020

PRODUCT

1 | HDPE 3933.5000PC 5 l Polycube 3933.9010 10 l Polycube 3933.9010PE 10 l Jerrycan

PACKING

3933.9020PE 20 L jerrycan This solution is phosphate buffered at 40 mmol/l.

20 l Polycube

Formalin solution

7385

CONT.

BOX

EINECS 200-001-8 NC CODE 2912 11 00

R: 40-43

S: 26-36/37/39-45



Formaldehyde (HCHO) pH at 25°C

38-42 g/l 6.8-7.2

NO. 7385.1000 11 7385.5000 5 l Polycube 7385.9020 20 l Polycube 200 l 7385.9200

PACKING

Phosphate Buffered with 40 mmol/l.

Formalin Solution See Formaldehyde solution

Formaline Phosphate buffered 38% (w/v) for Histology



See for detailed information www.jtbaker.com and select Cli



Formaline Phosphate buffered 19% (w/v) for Histology



Formaline Phosphate buffered 10% (v/v) [=4% (w/v)] for Histology

See for detailed information www.jtbaker.com and select Clini





Forma

Formamide 4028

'BAKER ULTRAPURE BIOREAGENT'

► HCONH₂

M = 45.04 g/mol

1 l = 1.13 kg

FLASHPOINT 155 °C

CAS NO. 75-12-7

EINECS 200-842-0

NC CODE 2924 19 00

R: 61

S: 24/25-37-45-53



For denaturing nucleic acids prior to analysis Assay (by GC) min. 99.5% Color (APHA) max. 10 Freezing Point 2.0-3.0°C

0.1

Absorbance of a 1 M Solution, Maximum (1-cm path vs water):

Zinc (Zn)

at 280 nm Trace Impurities (in ppm): Copper (Cu) max. 0.1 Iron (Fe) max. 0.5 Lead (Pb) max. 0.5 max. 0.5 PRODUCT **PACKING** CONT. NO. вох 4028.0100 100 ml 4028.0500 500 ml

Packaged under Nitrogen. Deionization of material is not

Formamide 7042 'BAKER ANALYZED'

► HCONH₂

M = 45.04 g/mol

1 l = 1.13 kg

FLASHPOINT 155 °C

CAS NO. 75-12-7 **EINECS** 200-842-0

NC CODE 2924 19 00

EC NO. 616 052 00 8

R: 61 **S:** 45-53

Assay	min. 99.5 %
Color (APHA)	max. 25
Density (g/ml) at 20°C	1.130-1.135
Freezing Point	2 - 3°C
Heavy Metals (as Pb)	max. 0.001%
Iron (Fe)	max. 0.001%
Residue after Ignition	max. 0.05%
Water (H ₂ O)	max. 0.5%

PRODUCT	PACKING	CONT.	
NO.		вох	
7042.1000	11	6	
7042.1000PE	1 HDPE		
7042.9025	25 l		
7042.9200	200		

Formic Acid

9820 99% / BAKER ANALYZED LC-MS Reagent

▶ НСООН

M = 46.03 g/mol

1 l = 1.22 kg

FLASHPOINT 69 °C

CAS NO. 64-18-6

EINECS 200-579-1

NC CODE 2915 11 00

EC NO. 607 001 00 0

UN/ID NO. 1779

ADR/RID 8 C3

IMDG 8/II

R: 35 **S:** 23-26-45



Certificate Provided Reporting Actual Lot Analysis min. 99% Residue after Evaporation max. 20 ppm LC-Gradient with Diode Array Detection (a.u.), 0.1% (v/v) aqueous solution: at 220 nm max. 0.002

at 254 nm max. 0.002 LC-MS Gradient Suitability Test (TIC, 100 to 2000

m/z), 0.1% (v/v) aqueous solution:

Positive ESI-MS Sensitive Impurities (as

Reserpine) max, 50 ppl Product Information (not specifications):

Density (g/ml) at 20°C 1.22

Trace Impurities (in ppb), 0.1% (v/v) aqueous solution:

Aluminium (AI) Calcium (Ca) max. 50 Iron (Fe) max. 50 Magnesium (Mg) max. 50 Potassium (K) max. 50 Sodium (Na) max 50

Ultraviolet Absorbance, 0.1% (v/v) aqueous solution (1.00-path vs water):

at 240 nm max. 0. at 260 nm max. 0.01

PRODUCT	PACKING	CONT.
NO.		вох
9820.0010AM	1 ml Amp. x 10	

6037

6166

8699

CONT.

BOX

Formic Acid

3% / 'BAKER ANALYZED'

M = 46.03 g/mol

1 l = 1.22 kg

FLASHPOINT 69 $^{\circ}\mathrm{C}$

▶ НСООН

CAS NO. 64-18-6

EINECS 200-579-1 NC CODE 2915 11 00

EC NO. 607 001 00 0

UN/ID NO. 1779

ADR/RID 8 C3

IMDG 8/II **R:** 35

S: 23-26-45

min. 98% Acetic Acid (CH₃COOH) max. 0.4% Ammonium (NH₄) max. 0.005% Chloride (CI) max. 0.001% Color (APHA) max. 15 Heavy Metals (as Pb) max. 0.001% Iron (Fe) max. 0.001% Residue after Evaporation max. 0.003% Sulfate (SO₄) max. 0.003%

passes test

89.5-90.5%

PRODUCT

6166.9025

NO.

Sulfite (SO₃)

Assay

PRODUCT	PACKING	CONT.
NO.		вох
6037.1000	11	6
6037.2500	2.5	
6037.9060	60 l	

Stored at a temperature of 15°C to 25°C.

PACKING

25 I

Formic Acid

90% / 'BAKER ANALYZED'

▶ нсоон

M = 46.03 g/mol

1 l = 1.19 kg

FLASHPOINT 69 °C

CAS NO. 64-18-6 EINECS 200-579-1

NC CODE 2915 11 00

EC NO. 607 001 00 0

UN/ID NO. 1779

ADR/RID 8 C3

IMDG 8/II **R:** 35

S: 23-26-36/37/39-45

Formic Acid

90% (PCS 7) / 'BAKER PCS'

▶ нсоон

M = 46.03 g/mol

1 l = 1.19 kg

FLASHPOINT 69 °C

CAS NO. 64-18-6 **EINECS** 200-579-1

NC CODE 2915 11 00

EC NO. 607 001 00 0

UN/ID NO. 1779

ADR/RID 8 C3

IMDG 8/II **R:** 34

S: 23-26-45

Assay 89.8-91.5% Clarity passes test Color (APHA)

PACKING CONT. **PRODUCT** вох 500 ml

Polymer Characterization Solvent.

Innovation is principal to our business.



Formi

Formic Acid

6008 88% / 'BAKER ANALYZED' / ACS

Meets ACS Specifications. Meets Reagent Specifications for testing USP/NF monographs Acetic Acid (CH₃COOH) max. 0.4 Ammonium (NH₄) max. 0.005 Color (APHA) max. Dilution Test passes te Residue after Evaporation max. 0.002 Sulfate (SO₄) max. 0.002 Sulfite (SO₃) passes te Trace Impurities (in ppm): Chloride (CI) max Heavy Metals (as Pb) max

PRODUCT	PACKING	CONT
NO.		BO
6008.0500	500 ml	

Formic Acid

9832 10

10% / BAKER ANALYZED LC-MS Reagent

Iron (Fe)

■ HCOOH

1 I = 1.02 kg

CAS NO. 64-18-6

NC CODE 2915 11 00

UN/ID NO. 1779

ADR/RID 8 C3

IMDG 8/II

R: 34

S: 26-36/37/39-45

Certificate Provided Reporting Actual Lot Analysis		
Assay	$10 \pm 0.5\% (v/v)$	
Product Information (not s	pecifications):	
Density (g/ml) at 20°C	1.028	
Trace Impurities (in ppb):		
Aluminium (AI)	max. 50	
Calcium (Ca)	max. 50	
Iron (Fe)	max. 50	
Magnesium (Mg)	max. 50	
Potassium (K)	max. 50	
Sodium (Na)	max. 50	
Ultraviolet Absorbance (1.00-cm path vs water):		
at 260 nm	max. 0.1	
at 280 nm	max. 0.05	

9832.1000PE 1 | HDPE 6

Element concentrations are at time of lot release.

CONT.

PACKING

PRODUCT

Formic Acid

corrosive

9826 0.1% / BAKER ANALYZED LC-MS Reagent

► HCOOH

M = 46.03 g/mol **1 l** = 1.0 kg **NC CODE** 2915 11 00

Certificate Provided Reporting Actual Lot Analysis $0.1 \pm 0.005\% (v/v)$ Assay Residue after Evaporation max. 5 ppm Gradient Elution Test (a.u.): at 220 nm max. 0.002 at 254 nm max. 0.002 LC-Gradient-MS Suitability Test (TIC, 100 to 2000 Positive ESI-MS Sensitive Impurities (as Reserpine) max. 50 ppb Product Information (not specifications): Density (g/ml) at 20°C 1.01 Trace Impurities (in ppb): Aluminium (AI) max. 50 Calcium (Ca) max. 50 Iron (Fe) max 50 Magnesium (Mg) max. 50 Potassium (K)

PRODUCT	PACKING	CONT.
NO.		вох
9826.1000GL	1 Glass	6
Element cond	entrations are at time of lot release.	

Sodium (Na) max. 50

Ultraviolet Absorbance (1.00-cm path vs water):
at 240 nm max. 0.1
at 260 nm max. 0.01



Formic Acid in Acetonitrile

0.1% / BAKER ANALYZED I C-MS Reagent

9824

11=	0.78 kg
FLASHPOINT	5 ℃
NC CODE	2926 90 95
UN/ID NO.	1648
ADR/RID	3 F1
IMDG	3/II
_	44 00 104 100

R: 11-20/21/22 **S:** 16-23-26-9





Certificate Provided Reporting Actual Lot Analysis		
Assay	$0.1 \pm 0.005\% (v/v)$	
Residue after Evaporation	max. 5 ppm	
Water (H ₂ O)	max. 0.02%	
Gradient Elution Test (a.u.):		
at 220 nm	max. 0.002	
at 254 nm	max. 0.002	
LC-Gradient-MS Suitability Test (TIC, 100 to 2000		
/ \		

m/zj:	
Positive ESI-MS Sensitive Impurities (as	
Reserpine)	max. 50 ppb
Product Information (not specifications	s):
Density (g/ml) at 20°C	0.78
Trace Impurities (in ppb):	
Aluminium (AI)	max. 50
Calcium (Ca)	max. 50

Calcium (Ca)	max. 50	
Iron (Fe)	max. 50	
Magnesium (Mg)	max. 50	
Potassium (K)	max. 50	
Sodium (Na)	max. 50	
Ultraviolet Absorbance (1.00-cm path vs water):		
at 260 nm	max. 0.1	
at 280 nm	max. 0.05	

PRODUCT	PACKING	CONT
NO.		BOX
9824.1000GL	1 Glass	(

Element concentrations are at time of lot release.

p-Fructopyranose

See D(-)Fructo



1563

D(-)Fructose

▶ QCH₂(CHOH)₃ÇOHCH₂OH

M = 180.16 g/mol **CAS NO.** 57-48-7 **EINECS** 200-333-3 NC CODE 1702 50 00

Lead (Pb) max. 0.001% Loss on Drying max. 0.5% Residue after Ignition max. 0.05% Specific Rotation [α] $\frac{20}{D}$ (c = 10 in H_2O) -91 - -93° Trace Impurities (in ppm): Arsenic (As) max. 1

PRODUCT PACKING 1563.0500 500 g

Fuchsin or Fuchsin S

See Acid Fuchsin



Gallium 1000 μ g/ml (Matrix: 1% nitric acid) / 'BAKER INSTRA-ANALYZED' / Plasma Standard

CONT.

BOX

M = 69.72 g/mol NC CODE 3822 00 00 R: 36/38 **S:** 26-37

► Ga



Certificate Provided Reporting Actual Lot Analysis Gallium (Ga) 998-1002 μg/ml PRODUCT PACKING 5714.0100 100 ml

Prepared from the highest purity raw material available, generally greater than 99.999% spectral purity. The content of the solution is confirmed to be accurate to within \pm 0.2%. Typically 1000 μ g/ml. The certificate of analysis provided reports actual lot

analysis. The certificate also lists the trace impurities.

www.jtbaker.com/europe



Galli

Gallium 10000 μ g/ml (Matrix: 1% nitric acid) / 'BAKER INSTRA-ANALYZED' / Plasma Standard 5758

► Ga

M = 69.72 g/mol

NC CODE 3822 00 00

R: 36/38

S: 26-37



Certificate Provided Reporting Actual Lot Analysis Gallium (Ga)

9980-10020 μg/ml

PRODUCT CONT. **PACKING** NO. вох 5758.0100 100 ml

Prepared from the highest purity raw material available, generally greater than 99.999% spectral purity. The content of the solution is confirmed to be accurate to within 0.2 %. Typically 10000 μ g/ml.

The certificate of analysis provided reports actual lot analysis. The certificate also lists the trace impurities.



Gallotannic Acid



Gelatin

ee for qualities meeting the requirements of cGMP, USP/NF, BP/Ph.Eur, JP or FCC, page 36



Gentian Violet

Germanium 1000 μ g/ml

oric acid) / 'BAKER INSTRA-ANALYZED' / Plasma Standard (Matrix: H₂O plus a trace of hyd

► Ge

M = 72.61 g/mol

NC CODE 3822 00 00

R: 20/21/22-36

S: 26-36/37



Certificate Provided Reporting Actual Lot Analysis

Germanium (Ge) 998-1002 μg/ml NO. 5762.0100 100 ml

PACKING

PRODUCT

Prepared from the highest purity raw material available, generally greater than 99.999% spectral purity. The

CONT.

CONT.

BOX

content of the solution is confirmed to be accurate to within \pm 0.2%. Typically 1000 μ g/ml. The certificate of analysis provided reports actual lot analysis. The certificate also lists the trace impurities.

Germanium 10000 µg/ml

rofluoric acid) / 'BAKER INSTRA-ANALYZED' / Plas

► Ge

5759

M = 72.61 g/mol

NC CODE 3822 00 00 R: 20/21/22-36

S: 26-36/37



Certificate Provided Reporting Actual Lot Analysis 9980-10020 µg/ml **PRODUCT PACKING** 5759.0100 100 ml

Prepared from the highest purity raw material available, generally greater than 99.999% spectral purity. The content of the solution is confirmed to be accurate to within 0.2 %. Typically 10000 μ g/ml.

The certificate of analysis provided reports actual lot analysis. The certificate also lists the trace impurities.

Giemsa

3856

HEMATOLOGY/CYTOLOGY/HISTOLOGY

1 I = 0.79 kg FLASHPOINT 11 °C

NC CODE 3822 00 00

UN/ID NO. 1992

ADR/RID 3 FT1 IMDG 3/II

R: 11-23/24/25-39-39/23/24/25

S: 16-36/37/39-45





PRODUCT PACKING CONT. BOX 3856 0100 100 ml Glass 3856.0500 500 ml Glass 3856.1000 1 | Glass 2.5 | Glass 3856.2500

Giemsa

Powder / HEMATOLOGY/CYTOLOGY/HISTOLOGY / Giemsa's Azur Eosin Methylene Blue Stain, suitable for staining blood smean

3815

CAS NO. 51811-82-6 **EINECS** 257-438-2

NC CODE 3204 19 00

R: 20/21/22-40-41 **S:** 22-26-36





Gill II Hematoxyline

See Hematoxyline (Papanicolaou 1)



Glucitol

See D(-)-Sorbitol



D(+)-Glucose Anhydrous

KÉR ANALYZED' / ACS




► HOCH₂ÇH(CHOH)₄Q **M** = 180.16 g/mol

CAS NO. 50-99-7 **EINECS** 200-075-1 NC CODE 1702 30 51

Chloride (CI) max. 0.01% Insoluble Matter max. 0.005% Loss on Drying at 105°C max. 0.2% Residue after Ignition max. 0.015% Specific Rotation $[\alpha]^{\frac{26}{D}}$ +52.5-+53.0° Starch passes test Sulfate and Sulfite (as SO₄) max. 0.005% Titrable Acid (meq/g) max. 0.002 Trace Impurities (in ppm):

Exceeds ACS Specifications

Heavy Metals (as Pb) max. 5 Iron (Fe) max. 5



D(+)-Glucose Anhydrous Ultrex



0113

D(+)-Glucose Monohydrate

'BAKER ANALYZED

► HOCH₂CH(CHOH)₄O.H₂O **M** = 198.17 g/mol **CAS NO.** 50-99-7 **EINECS** 200-075-1 NC CODE 1702 30 51

Chloride (CI) max. 0.005% Insoluble Matter max. 0.005% S Т

Residue after Ignition	max. 0.015%
Starch	passes test
Sulfate and Sulfite (as SO ₄)	max. 0.005%
Trace Impurities (in ppm):	
Arsenic (As)	max. 0.2
Heavy Metals (as Pb)	max. 5
Iron (Fe)	max. 5

PACKING	CONT.
	вох
1 kg	6
50 kg	
	1 kg

L-Glutamic Acid

See for qualities meeting the requirements of cGMP, USP/NF, BP/Ph.Eur, JP or FCC, page 36



L-Glutamine

See for qualities meeting the requirements of cGMP, USP/NF, BP/Ph.Eur, JP or FCC, page 36



See for qualities meeting the requirements of cGMP, USP/NF, BP/Ph.Eur, JP or FCC, page 36





Glyce

Glycerol

4043

/ 'BAKER ULTRAPURE BIOREAGENT'

► HOCH₂CHOHCH₂OH

M = 92.10 g/mol **1 l** = 1.23 kg

FLASHPOINT 160 $^{\circ}\mathrm{C}$

CAS NO. 56-81-5

EINECS 200-289-5

NC CODE 1520 00 00

For Enzyme stabilization, Freezing and Electrophoresis Buffers

min. 99.5% Assay Color (APHA) max. 10 DNAase none detected Fatty Acid Esters (as Butyric Acid) max.0.05% Protease none detected Residue after Ignition max.0.005% RNAase none detected

Absorbance of a 1 M Solution, Maximum (1-cm path vs water):

Zinc (Zn)

Heavy Metals (as Pb)

at 280 nm 0.05 Trace Impurities (in ppm): Iron (Fe) max.5 Lead (Pb) max. 2 Mercury (Hg) max. 10

max. 2

max. 2

84-89%

passes test

0.1

PRODUCT

PRODUCT	PACKING	CONT.
NO.		вох
4043.0500	500 ml	
4043.4000	4 I Glass	

7044

Glycerol

BAKER ANALYZED

► HOCH₂CHOHCH₂OH

M = 92.10 g/mol

1 l = 1.26 kg

FLASHPOINT 160 °C **CAS NO.** 56-81-5

EINECS 200-289-5

NC CODE 1520 00 00

Assay	min. 99.5%
Color	passes test
Density (g/ml) at 25°C	min. 1.255
Fatty Acid Esters	max. 0.05%
Neutrality	passes test
Residue after Ignition	max. 0.005%
Silver Reducing Substances	passes test
Substances Darkened by H ₂ SO ₄	passes test
Sulfate (SO ₄)	max. 0.001%
Trace Impurities (in ppm):	
Chloride (CI)	max. 5

NO. BOX 7044.2500 2.5 l 7044.9010 10 I 7044.9025 25 I 7044.9200 For safe handling of 25 l tin cans, see Self-closing tap.

CONT.

CONT.

вох

PACKING

7158

Glycerol

► HOCH₂CHOHCH₂OH

M = 92.10 g/mol

11 = 1.23 kg

FLASHPOINT 160 °C CAS NO. 56-81-5

EINECS 200-289-5

NC CODE 1520 00 00

Assay

PRODUCT PACKING NO. 7158.1000 11 7158.5000 5 I HDPE

Glycine 4059

'BAKER ULTRAPURE BIOREAGENT'

► NH₂CH₂COOH

M = 75.07 g/mol

CAS NO. 56-40-6 **EINECS** 200-272-2

NC CODE 2922 49 10

For Electrophoresis, liquid Chromatography and Molecular Biology

Assay (dried basis, by non-aqueous

min. 99.5% acid-base titration) DNAase none detected

Homogenity (by TLC) no extraneous spots Iron (Fe) max. 0.003%

Loss on Drying at 105°C max. 0.2% Protease none detected

Readily Carbonisable Substances passes test Residue after Ignition max. 0.1% RNAase none detected

Absorbance of a 1 M Solution, Maximum (1-cm path vs water):

at 280 nm Trace Impurities (in ppm):

Solution Test

Arsenic (As) max. 3 Heavy Metals (as Pb) max. 5

PRODUCT	PACKING	CONT
NO.		ВОХ
4059.0250	250 g	
4059.1000	1 kg	
4059.5000	5 kg	

1504

Glycine
'BAKER ANALYZED' Biochemical

► NH₂CH₂COOH Assay (by Perchloric Acid titrn.) min. 99% **M** = 75.07 g/mol Heavy Metals (as Pb) max. 0.002% **CAS NO.** 56-40-6 Homogenity (by TLC) passes test **EINECS** 200-272-2 Loss on Drying at 105°C max. 0.2% NC CODE 2922 49 10 Readily Carbonisable Substances passes test Residue after Ignition max. 0.1% Trace Impurities (in ppm):

Arsenic (As)

Assay (HOCH₂COOH)

Assay

Iron (Fe)

PRODUCT	PACKING	CONT.	
NO.		вох	
1504.0100	100 g		
1504.0500	500 g	6	

Glycine 'BAKER'

1933

max. 3

max. 5

68-72%

35-45%

PRODUCT

1930.0500

PRODUCT

M834-09

NO.

PACKING

500 ml


PACKING

4 | Glass

 NH_2CH_2COOH M = 75.07 g/mol

CAS NO. 56-40-6 **EINECS** 200-272-2 **NC CODE** 2922 49 10

98.5-101.5% Assay passes test Appearance Chloride (CI) max. 0.007% Heavy Metals (as Pb) max. 0.002% Hydrolyzable substances passes test Identification passes test Loss on Drying at 105°C max. 0.2% Organic Volatile Impurities passes test pH of 5% Solution at 20°C 5.9-6.4 Residue after Ignition max. 0.1% Sulfate (SO₄) max. 0.0065%



Glycocoll



Glycol

See Ethylene Glycol



1930

CONT.

вох

BOX

Glycollic Acid

68-72% / 'BAKER'

► HOCH₂COOH

M = 76.05 g/mol

CAS NO. 79-14-1

EINECS 201-180-5

NC CODE 2918 19 80 UN/ID NO. 3265 ADR/RID 8 C3

IMDG 8/II

R: 34

S: 26-36/37/39-45



Glyoxal solution, 40% in H₂O / Technical M834-09

► OHCCHO **M** = 58.04 g/mol

1 l = 1.27 kg

CAS NO. 107-22-2 EINECS 203-474-9

NC CODE 2912 19 00 **EC NO.** 605 016 00 7

R: 20-36/38-40-43 **S:** 36/37



_.

Glyoxaline



Gold

Gold 1000 μg/ml (Matrix: 5% hydrochloric acid

5763

ric acid) / 'BAKER INSTRA-ANALYZED' / Plasma Standard

► Au M = 196.97 g/mol NC CODE 3822 00 00

Certificate Provided Reporting Actual Lot Analysis Gold (Au) 998-1002 μg/ml)

998-1002 μg/ml

PRODUCT CONT. **PACKING** NO. вох 5763.0100 100 ml

Prepared from the highest purity raw material available, generally greater than 99.999% spectral purity. The content of the solution is confirmed to be accurate to within \pm 0.2%. Typically 1000 $\mu g/ml$ The certificate of analysis provided reports actual lot analysis. The certificate also lists the trace impurities.

Gold 1000 μg/ml (Matrix: 5% hydrochloric aci

ric acid) / 'BAKER INSTRA-ANALYZED' / Atomic Absorption Standard

Gold (Au)

► Au M = 196.97 g/mol NC CODE 3822 00 00

6952

PRODUCT PACKING CONT. NO. BOX 6952.0100 100 ml 6952.0500 500 ml

Prepared by dissolution of high purity raw materials (min. 99.99% spectral purity). Assays are verified by ICP against standards traceable to NIST. Standard Reference Material numbers (SRM) are printed on each label.

Gold 10000 μg/ml

5730 'BAKER INSTRA-ANALYZED' / Plasma Standard

► Au M = 196.97 g/mol NC CODE 3822 00 00

Certificate Provided Reporting Actual Lot Analysis 9980-10020 μg/ml Gold (Au)

PRODUCT PACKING CONT. NO. вох 5730.0100 100 ml

Prepared from the highest purity raw material available, generally greater than 99.999% spectral purity. The content of the solution is confirmed to be accurate to within 0.2 %. Typically 10000 $\mu g/ml$. The certificate of analysis provided reports actual lot analysis. The certificate also lists the trace impurities.

Gold Chloride

e Tetrachloroauric(III) Acid

Granusic A

0506 anular / 'BAKER ANALYZED' / ACS

► P₂O₅

M = 141.94 g/mol **CAS NO.** 1314-56-3 **EINECS** 215-236-1 NC CODE 2809 10 00 **EC NO.** 15 010 00 0

UN/ID NO. 1807

ADR/RID 8 C2 IMDG 8/II

> **R:** 35 **S:** 22-26-45

Exceeds ACS Specifications		
Assay	min. 99.0%	
Ammonium (NH ₄)	max. 0.005%	
Heavy Metals (as Pb)	max. 0.005%	
Insoluble Matter	max. 0.01%	
Phosphorus Trioxide (P _o O _o)	max. 0.02%	

PRODUCT	PACKING	CONT.
NO.		вох
0506.0500	500 g	
0506.1000	1 kg	6

Guanidine Hydrochloride

e Guanidinium Chloride



Guanidine Hydrochloride

ements of cGMP, USP/NF, BP/Ph.Eur, JP or FCC, page 36

PRODUCT

4168.9012

PRODUCT

4075.0500

4075.9012

PRODUCT

4045.4000

PACKING

500 g

12 kg

PACKING

4 | Glass

NO.

PACKING

12 kg

J.T.Baker

Guanidine Hydrochloride SolutionsSee for qualities meeting the requirements of cGMP, USP/NF, BP/Ph.Eur, JP or FCC, page 36

Guanidinium Chloride

'BAKER ULTRAPURE BIOREAGENT' / For Protein Denaturation and Nucleic Acid Purification

GMP Manufactured Product

4168

CONT.

вох

•	NH ₂ (:NH)N	H ₂	.HC	:I
	M	_	QE.	53

M = 95.53 g/mol **CAS NO.** 50-01-1 **EINECS** 200-002-3

NC CODE 2925 20 00 **EC NO.** 607 148 00 0

R: 22-36/38

S: 22



GIVIP IVIanutacturea Product	
Assay	99.5-101.0%
Acidity (as HCI)	max. 0.01%
DNase Activity	none detected
Loss on Drying at 105°C	max. 0.5%
Melamine and Related	
Compounds (g/l)	act. value reported
Nitrate (NO ₃)	max. 0.005%
Particulate Matter	passes test
Protease Activity	none detected
Residue after Ignition	max. 0.05%
RNase Activity	none detected
Sulfate (SO ₄)	max. 0.005%
Absorbance of a 6M Aqueous S	Solution:
at 230 nm	act. value reported
at 260 nm	max. 0.03
at 275 nm	max. 0.03
Trace Impurities (in ppm):	
Arsenic and Antimony (as As)	max. 0.5
Copper (Cu)	max. 0.2
Cyanide (CN)	max. 0.1
Iron (Fe)	max. 2
Lead (Pb)	max. 0.5

Guanidinium Chloride

4075

CONT.

BOX



M = 95.53 g/mol **CAS NO.** 50-01-1

EINECS 200-002-3 NC CODE 2925 20 00

EC NO. 607 148 00 0 R: 22-36/38

S: 22



GMP Manufactured Product	
Assay (NH ₂ C(:NH)NH ₂ .HCI)	99.5-101.0%
Appearance	passes test
Iron (Fe)	max. 2 ppm
Loss on Drying at 105°C	max. 0.5%
Melamine and Related	
Compounds (g/l)	act. value reported
Melting Range	information only
Water-Insoluble Matter	max. 0.3%
Ultraviolet Absorbance (1.00-cm path vs water):	

at 260 nm act. value reported at 275 nm act. value reported

Guanidinium Chloride

8 mol/l / 'BAKER ULTRAPURE BIOREAGENT'

4045

CONT.

► NH₂C(:NH)NH₂.HCl

M = 95.53 g/mol

CAS NO. 50-01-1 EINECS 200-002-3

NC CODE 2925 20 00

EC NO. 607 148 00 0

R: 22-36/38

S: 26-36/37



GMP Manufactured Product	
Cyanide (CN)	passes test
Density (g/ml) at 25°C	1.1-1.3
DNase Activity	none detected
Molarity (by titrimetry)	7.8-8.2
Protease Activity	none detected
RNase Activity	none detected
Solution Test	passes test
Sulfate (SO ₄)	max. 0.005%
Trace Impurities (in ppm):	
Copper (Cu)	max. 0.1
Iron (Fe)	max. 1
Ultraviolet Absorbance (1.00-cr	n path vs water;
curve emooth throughout state	d range with no

curve smooth throughout stated range with no extraneous impurity peaks):

at 260 nm max. 0.04

Harris Hematoxyline

See Hematoxyline (Papanicolaou 1)



Hemac

Hemacolour

3719

1 l = 0.79 kg FLASHPOINT 11 °C

NC CODE 3822 00 00

UN/ID NO. 1992

ADR/RID 3 FT1

IMDG 3/II

R: 11-23/24/25-36/38-39/23/24/25

S: 16-36/37-45-7





Staining solution for blood smears

PRODUCT **PACKING** CONT. вох 3719.0300 300 ml Glass 3719.1000 1 | Glass

Hematoxyline

3870

NC CODE 3822 00 00 **UN/ID NO.** 2810 ADR/RID 6.1 T1

IMDG 6.1/III

R: 22

S: 36



Hematoxyline Solution according to Mayer for use in HE staining

PRODUCT **PACKING** CONT. вох NO. 3870.1000 1 | Glass 3870.2500 2.5 | Glass

Hematoxyline (Papanicolaou 1)

NC CODE 3822 00 00

R: 22

S: 36



Hematoxyline (Papanicolaou I), modified to match Gill II and Harris

PRODUCT PACKING BOX 3873.1000 1 | Glass 3873.2500 2.5 | Glass



Hematoxyline Solution for HE staining

d information www.jtbaker.com and select Cl

HEPES

'BAKER ULTRAPURE BIOREAGENT' / For Liquid Chromatography and Molecular Biology Buffers

► CH2CH2N(CH2CH2OH)CH2CH2NCH2CH2SO3H **M** = 238.31 g/mol **CAS NO.** 7365-45-9

EINECS 230-907-9

NC CODE 2933 59 95

Absorbance of a 0.1M Solution /1 cm path vs	
RNAase	none detected
Protease	none detected
pH of 5% Solution at 25°C	5.0-6.5
DNAase	none detected
Ash (sulfated)	max. 0.2%
Appearance	passes test
Assay	min. 99%

water): at 260 nm max. 0.

Trace Impurities (in ppm):

Iron (Fe)

Arsenic and Antimony (as As) Heavy Metals (as Pb)

NO.	
4018.0025	
4018.0100	
4018.0500	
4018.5000	
Ref. Good, N., e	9
	4018.0025 4018.0100 4018.0500 4018.5000

PRODUCT PACKING CONT вох 25 g Glass 100 g 500 g 5 kg

et al. (1966) Biochemistry, 5:467.

Innovation is principal to our business.



HEPES sodium salt

4153 'BAKER UI TRAPURE BIOREAGENT

min. 99.0%

min. 99.0%

0.684

C₈H₁₇N₂NaO₄S

M = 260.29 g/mol **CAS NO.** 75277-39-3 NC CODE 2933 59 95

min. 99.0%
passes test
none detected
max. 1 ppm
max. 1.0%
max. 3.0%
act. value reported
none detected
none detected
passes test

CONT.
вох

n-Heptane

6 n-Heptane / 'BAKER ULTRA RESI-ANALYZED' / for Organic Residue Analysis

9338

► CH₃(CH₂)₅CH₃

M = 100.21 g/mol

1 l = 0.68 kg

FLASHPOINT −4 °C

CAS NO. 142-82-5

EINECS 205-563-8

NC CODE 2901 10 00

EC NO. 601 008 00 2

UN/ID NO. 1206

ADR/RID 3 F1 IMDG 3/II

R: 11-38-50/53-65-67

S: 16-29-33-60-61-62-9







Color (APHA) max. 10 Residue after Evaporation max. 1 ppm Substances Darkened by H2SO4 passes test Water (H₂O) max. 0.01% ECD Sensitive Impurities (as Heptachlor Epoxide): Single Impurity Peak (pg/ml) max. 10 FID-Sensitive Impurities (as 2-Octanol): Single Impurity Peak (ng/ml) max. 5 Neat solvent front characterization: ECD-Sensitive Impurities (as Ethylene Dibromide): Single Impurities (ng/ml) max. 5

Assay (by GC) (corrected for water)

PRODUCT	PACKING	CONT.
NO.		ВОХ
9338.1000	11	
9338 4000	4 l Glass	1

99% n-Heptane / 'BAKER HPLC ANALYZED' / for Use in Liquid Chromatography and Spectrophotometry

Density (g/ml) at 20°C

Assay (by GC) (corrected for water)

9177

CH₃(CH₂)₅CH₃

M = 100.21 g/mol

11 = 0.68 kg

FLASHPOINT −4 °C

CAS NO. 142-82-5

EINECS 205-563-8

NC CODE 2901 10 00

EC NO. 601 008 00 2

UN/ID NO. 1206

ADR/RID 3 F1 IMDG 3/II

R: 11-38-50/53-65-67

S: 16-29-33-60-61-62-9









PRODUCT	PACKING	CONT.
NO.		вох
9177.1000	11	
9177.2500	2.5	4
9177.5000	5 EcoTainer	

Filtered through a 0.2 micron filter Packaged under Nitrogen.

Questions or suggestions, please contact us at jtbaker.nl@emea.tycohealthcare.com



Hepta

n-Heptane

9365 99% n-Heptane / BakerDRY / Low Water Solvent

► CH₃(CH₂)₅CH₃

M = 100.21 g/mol

1 l = 0.68 kg

 $\textbf{FLASHPOINT} \quad -4~^{\circ}\text{C}$

CAS NO. 142-82-5

EINECS 205-563-8

NC CODE 2901 10 00

EC NO. 601 008 00 2

UN/ID NO. 1206

ADR/RID 3F1

IMDG 3/II

R: 11-38-50/53-65-67

S: 16-29-33-60-61-62-9







Assay (CH₃(CH₂)₅CH₃)(by GC)

Residue after Evaporation

Water (by KF, coulometric)

Assay (CH₃(CH₂)₅CH₃)(by GC, corrected for

Color (APHA)

water)

Color (APHA)

Assay (by GC)

Boiling Range

Color (APHA)

Recorded Boiling Point

Residue after Evaporation

Water (by KF, coulometric)

n-Heptane

99% n-Heptane / 'BAKER BIO-ANALYZED'

► CH₃(CH₂)₅CH₃

M = 100.21 g/mol

1 l = 0.68 kg

FLASHPOINT -4 °C

CAS NO. 142-82-5

EINECS 205-563-8

NC CODE 2901 10 00

EC NO. 601 008 00 2

UN/ID NO. 1206

ADR/RID 3 F1

IMDG 3/II

R: 11-38-50/53-65-67

S: 16-29-33-60-61-62-9







n-Heptane

99% n-Heptane / 'BAKER ANALYZED'

8110 ► CH₃(CH₂)₅CH₃

M = 100.21 g/mol

1 I = 0.68 kg

 $\textbf{FLASHPOINT} \quad -4~^{\circ}\text{C}$

CAS NO. 142-82-5 **EINECS** 205-563-8

NC CODE 2901 10 00

EC NO. 601 008 00 2

UN/ID NO. 1206

ADR/RID 3 F1

IMDG 3/II

R: 11-38-50/53-65-67

S: 16-29-33-60-61-62-9







Residue after Evaporation	max. 0.001%
Total Sulfur Compounds (as S)	max. 0.005%
Water (H ₂ O)	max. 0.01%
Physical Data (not specifications):	
Density (g/ml) at 20°C	0.683-0.684
Trace Impurities (in ppm):	
Aluminium (AI)	max. 0.5
Barium (Ba)	max. 0.1
Boron (B)	max. 0.02
Cadmium (Cd)	max. 0.05
Calcium (Ca)	max. 0.5
Chromium (Cr)	max. 0.02
Cobalt (Co)	max. 0.02
Copper (Cu)	max. 0.02
Iron (Fe)	max. 0.1
Lead (Pb)	max. 0.1
Magnesium (Mg)	max. 0.1
Manganese (Mn)	max. 0.02
Nickel (Ni)	max. 0.02
Tin (Sn)	max. 0.1
Zinc (Zn)	max. 0.1

min. 98.5%

max. 0.001%

max. 20 ppm

min. 99.0%

max. 1 ppm

max. 100 ppm

max. 10

min. 99%

97.5-99°C

max.10

98.4°C

PRODUCT

PACKING

max. 10

PRODUCT

9365.1000

PRODUCT

9185.1000GL 1 | Glass

PACKING

NO.

PACKING

11

CONT.

вох

CONT.

CONT.

вох

NO.		BOX
8110.1000	11	
8110.2500	2.5	4
8110.9025	25 l	
3110.9200	200	
For safe han	dling of 25 l tin cans, see Self-clo	sing tap.
	<u> </u>	J - 1

PRODUCT

8662.1000

8662.2500

8662.9025

8662.9200

PRODUCT

8111.1000

8111.5000

8111.9200

NO.

min. 95%

98-99°C.

94-99°C

max. 0.05

max. 0.05

PACKING

11

2.5 l

25 I

PACKING

5 | EcoTainer

200 l

For safe handling of 25 l tin cans, see Self-closing tap.

8662

8111

CONT.

вох

CONT.

вох



n-Heptane

Assay (by GC)

Boiling Point

Boiling Range

► CH₃(CH₂)₅CH₃

M = 100.21 g/mol

1 l = 0.68 kg

 $\textbf{FLASHPOINT} \quad -4~^{\circ}\text{C}$

CAS NO. 142-82-5

EINECS 205-563-8

NC CODE 2901 10 00

EC NO. 601 008 00 2

UN/ID NO. 1206

ADR/RID 3 F1

IMDG 3/II

R: 11-38-50/53-65-67

S: 16-29-33-60-61-62-9







n-Heptane

► CH₃(CH₂)₅CH₃ **M** = 100.21 g/mol

1 l = 0.68 kg

FLASHPOINT -4 °C

CAS NO. 142-82-5

EINECS 205-563-8

NC CODE 2901 10 00

EC NO. 601 008 00 2

UN/ID NO. 1206

ADR/RID 3F1

IMDG 3/II

R: 11-38-50/53-65-67

S: 16-29-33-60-61-62-9







1-Heptanesulfonic Acid Sodium Salt

► CH₃(CH₂)₆SO₃Na **M =** 220.26 g/mol CAS NO. 22767-50-6

EINECS 245-210-5

NC CODE 2904 10 00

For Ion-Pair Chromatography of Basic Compounds		
Assay (acidimetric)	min. 98.0%	
UV Absorbance of 0.25M Solution:		
at 200 nm	max. 0.2	
at 210 nm	max. 0.08	
at 220 nm	max. 0.06	
at 230 nm	max. 0.05	

PRODUCT NO.	PACKING	CONT BOX
2173.0025	25 g Glass	

Hexachlorplatinum(IV) Acid Hexahydrate

at 240 nm

at 250 nm

See Chloroplatinic Acid Hexahydrate



2173

www.jtbaker.com/europe



Hexad

n-Hexadecane

7194

► CH₃(CH₂)₁₄CH₃

M = 226.45 g/mol

1 l = 0.77 kg

 $\textbf{FLASHPOINT} \quad 135~^{\circ}\text{C}$ **CAS NO.** 544-76-3

EINECS 208-878-9

NC CODE 2901 10 00

R: 65

S: 23A-24-62

1,1,1,3,3,3-Hexafluoro-2-propanol

7108

► CF₃CHOHCF₃

M = 168.04 g/mol

1 l = 1.62 kg

CAS NO. 920-66-1

EINECS 213-059-4

NC CODE 2905 59 10

UN/ID NO. 3265

ADR/RID 8 C3

IMDG 8/II

R: 20-34-37

S: 26-36/37/39-45



Hexamethyldisilazane (HMDS) CMOS, Finyte Grade

Assay (by GC)

Freezing Point

Assay (by GC) (corrected for water)

Product Information (not specifications):

Residue after Evaporation

Boiling Point (typical)

Color (APHA)

Water (H₂O)

Hexamethylenetetramine

'BAKER ANALYZED' / ACS

► (CH₂)₆N₄

3371

M = 140.19 g/mol

CAS NO. 100-97-0 EINECS 202-905-8

NC CODE 2933 69 20

EC NO. 612 101 00 2

UN/ID NO. 1328

ADR/RID 4.1 F1 IMDG 4.1/III

R: 11-42/43

S: 16-22-24-37





Exceeds ACS Specifications Assay (dried) min. 99.0% Heavy Metals (as Pb) max. 0.001% Insoluble Matter max. 0.005% max. 2.0% Loss on Drying Residue after Ignition max. 0.1%

PRODUCT PACKING CONT. вох

CONT.

вох

CONT.

вох

PRODUCT

7194.0250

PRODUCT

7108.0100

NO.

min. 99%

min. 99%

max. 0.01%

max. 0.05%

max. 10

58.2°C

18-19°C

PACKING

250 ml

PACKING

100 ml

J.T.Baker: over 100 years of experience.

See chapter 1 of this catalogue.



In several industries water determination with Karl Fischer reagents is an important analytical technique, for both research and quality control purposes.

Major fields of application for Karl Fischer reagents:

- Pharmaceutical/vitamins
- Food / drinks / beverages / flavours / wine
- Petrochemistry
- Other chemical industries

We offer a complete product line:

- Volumetric reagents (for > 100 ppm water)
 - One component reagents
 - Two component reagents
- Coulometric reagents (for < 100 ppm water)
 - With diaphragm
 - Without diaphragm

Program highlights:

Volumetric

- Safer: pyridine free
- Productive: end point reached quickly, accurately and with high reproducibility
- Accurate: stable end point and less drift.

Coulometric

- Accurate: less background noise, especially below 10 ppm
- Trouble free: a lower amount of reduction products
- Fast: high imidazole buffer concentration
- Safer: HYDRA-POINT coulometric gen and oven are halogen free reagents
- Easy: only five products for a broad range of applications
- Innovative: coulometric reagent line is patent pending (reference Prof. A. Cedergren, Department of Chemistry, Umea University, Sweden)

Our trade name for the Karl Fischer reagents line is HYDRA-POINT. For our complete reagent line, see the alphabetical listings.



Hexam

Hexamethylenetetramine

1113

► (CH₂)₆N₄

M = 140.19 g/mol **CAS NO.** 100-97-0 **EINECS** 202-905-8

NC CODE 2933 69 20 **EC NO.** 612 101 00 2


UN/ID NO. 1328

ADR/RID 4.1 F1

IMDG 4.1/III

R: 11-42/43 **S:** 16-22-24-37





PRODUCT	PACKING	CONT.
NO.		вох
1113.9050	50 kg	



Hexamethyl p-rosanilinium Chloride



Hexamine

Hexane

99% n-Hexane / 'BAKER ANALYZED'

► CH₃(CH₂)₄CH₃ **M** = 86.18 g/mol

1 l = 0.66 kg FLASHPOINT < -23 °C **CAS NO.** 110-54-3

EINECS 203-777-6

NC CODE 2901 10 00 EC NO. 601 037 00 0

UN/ID NO. 1208 ADR/RID 3 F1 IMDG 3/II

R: 11-38-48/20-51/53-62-65-67

S: 16-29-33-36/37-61-62-9









Assay (by GC)	min. 99%
Aromatic Compounds (as C ₆ H ₆)	max. 0.02%
Boiling Range	max. 2°C
Density (g/ml) at 25°C	max. 0.657
Recorded Boiling Point	68.7°C
Residue after Evaporation	max. 5 ppm
Sulfur Compounds (as S)	max. 0.005%
Thiophene	passes test
Water (H ₂ O)	max. 0.02%
Water Soluble Titrable Acid (meq/g)	max. 0.0003
Trace Impurities (in ppm):	
Aluminium (AI)	max. 0.5
Barium (Ba)	max. 0.1
Boron (B)	max. 0.02
Cadmium (Cd)	max. 0.05
Calcium (Ca)	max. 0.5
Chromium (Cr)	max. 0.02
Cobalt (Co)	max. 0.02
Copper (Cu)	max. 0.02
Iron (Fe)	max. 0.1
Lead (Pb)	max. 0.1
Magnesium (Mg)	max. 0.1
Manganese (Mn)	max. 0.02
Nickel (Ni)	max. 0.02
Tin (Sn)	max. 0.1
Zinc (Zn)	max. 0.1

PRODUCT	PACKING	CONT
NO.		ВОХ
8668.1000	11	ϵ
8668.2500	2.5	4
8668.5000	5 EcoTainer	
8668.9025	25	
EcoTainer, th	e metal solvent can for mor	re safety in the lab
For safe han	dling of 25 l tin cans, see Se	elf-closing tap.
TOT SUIC HUIT	aning of 25 i till calls, see se	in closing tap.

Certificates of Analysis are available at www.jtbaker.com/europe

8205

► CH₃(CH₂)₄CH₃

		~~				
М	=	-86.	18	0/	/mol	ı

1 l = 0.66 kg

FLASHPOINT $-23 \, ^{\circ}\text{C}$

CAS NO. 110-54-3

EINECS 203-777-6

NC CODE 2901 10 00

EC NO. 601 037 00 0

UN/ID NO. 1208

ADR/RID 3F1

IMDG 3/II

R: 11-38-48/20-51/53-62-65-67

S: 16-29-33-36/37-61-62-9







Exceeds ACS Specifications	
Assay (n-hexane)	min. 99%
Assay (sum of 5 isomers, total hexanes	
plus methylcyclopentane)	min. 99%
Aromatic Compounds (as C ₆ H ₆)	max. 0.02%
Color (APHA)	max. 10
Residue after Evaporation	max. 5 ppm
Sulfur Compounds (as S)	max. 0.005%
Thiophene	passes test
Water Soluble Titrable Acid (meq/g)	max. 0.0002
Illtraviolet Absorbance (1.00 em noth	io seintori

Ultraviolet Absorbance (1.00-cm path vs water; curve smooth throughout stated range with no extraneous impurity peaks):

at 210 nm	max. 1.00
at 220 nm	max. 0.20
at 230 nm	max. 0.10
at 240 nm	max. 0.04
at 250 nm	max. 0.02
at 280-400 nm	max. 0.01

Meets ACS Specifications

Residue after Evaporation

Sulfur Compounds (as S)

Water (by KF, coulometric)

Color (APHA)

Thiophene

Assay (C₆H₁₄)(by GC)(total isomers)

Water Soluble Titrable Acid (meq/g)

PRODUCT NO.	PACKING	CONT. BOX
8205.0500	500 ml	
8205.1000	11	

PRODUCT

9277.1000

NO.

min. 98.5%

max. 2 ppm

max. 0.005%

passes test

max. 20 ppm

max. 0.0003

max. 10

PACKING

11

Hexane

BakerDRY / Low Water Solvent / ACS

9277

CONT.

BOX

► CH₃(CH₂)₄CH₃

M = 86.18 g/mol

1 I = 0.66 kg

FLASHPOINT <-23 °C

CAS NO. 110-54-3

EINECS 203-777-6

NC CODE 2901 10 00

EC NO. 601 037 00 0

UN/ID NO. 1208 ADR/RID 3 F1

IMDG 3/II

R: 11-38-48/20-51/53-62-65-67

S: 16-29-33-36/37-61-62-9



environment





Hexane

95% n-Hexane / 'BAKER ULTRA RESI-ANALYZED' / for Organic Residue Analysis

► CH₃(CH₂)₄CH₃ **M** = 86.18 g/mol

1 l = 0.66 kg

 $\textbf{FLASHPOINT} \quad < \text{-}20~^{\circ}\text{C}$

CAS NO. 110-54-3

EINECS 203-777-6

NC CODE 2901 10 00

EC NO. 601 037 00 0 **UN/ID NO.** 1208

ADR/RID 3 F1

IMDG 3/II

R: 11-38-48/20-51/53-62-65-67

S: 16-29-33-36/37-61-62-9









Assay (total saturated C_6 isomers)(by GC,		
corrected for water)	min. 99.5%	
Color (APHA)	max. 10	
n-Hexane (by GC) (corrected for water)	min. 95%	
Residue after Evaporation	max. 1 ppm	
Substances Darkened by H ₂ SO ₄	passes test	
Water (H ₂ O) (by Coulometry)	max. 0.05%	
FCD Sensitive Impurities (as Hentachlor Epoxide):		

Single Impurity Peak (pg/ml) max. 10

FID-Sensitive Impurities (as 2-Octanol):

Single Impurity Peak (ng/ml) max. 5

Neat solvent front characterization: ECD-Sensitive Impurities (as Ethylene Dibromide):

Single Impurities (ng/ml)

PRODUCT	PACKING	CONT BOX
9262.1000	11	(
9262.2500	2.5	
7202.2300	2.31	



Hexan

Hexane

% n-Hexane / 'BAKER HPLC ANALYZED' / for use in High Performance Liquid Chromatograph

► CH₃(CH₂)₄CH₃ **M** = 86.18 g/mol **1 l** = 0.66 kg **FLASHPOINT** $-23 \, ^{\circ}\text{C}$ CAS NO. 110-54-3 **EINECS** 203-777-6

9304

NC CODE 2901 10 00 EC NO. 601 037 00 0 **UN/ID NO.** 1208 ADR/RID 3F1

IMDG 3/II

R: 11-38-48/20-51/53-62-65-67 **S:** 16-29-33-36/37-61-62-9



Assay (by GC)	111111. 3070		
Residue after Evaporation	max. 2 ppm		
Substances Darkened by H ₂ SO ₄	passes test		
Water (H ₂ O)	max. 0.01%		
Ultraviolet Absorbance (1.00-cm path vs water;			
curve smooth throughout stated range with no			
extraneous impurity peaks):			
at 210 nm	max. 0.2		

Oltraviolet Absorbance (1.00-cm path vs water;		
curve smooth throughout stated range with no extraneous impurity peaks):		
at 220 nm	max. 0.08	
at 254 nm	max. 0.008	
at 280 nm	max. 0.005	
at 350 nm	max. 0.005	
UV Cut-off, nm	max. 192	

PRODUCT	PACKING	CONT.
NO.		вох
9304.1000	11	6
9304.2500	2.5	4
9304.5000	5 EcoTainer	

EcoTainer, the metal solvent can for more safety in the lab.

Filtered through a 0.2 micron filter. Packaged under Nitrogen.

PACKING

PRODUCT

Hexane

95% n-Hexane / 'BAKER ANALYZED' / Ultraviolet Spectrophotometry

► CH₃(CH₂)₄CH₃ **M** = 86.18 g/mol **1 l** = 0.66 kg FLASHPOINT <-22 °C **CAS NO.** 110-54-3 **EINECS** 205-563-8 NC CODE 2901 10 00 EC NO. 601 037 00 0 **UN/ID NO.** 1208 ADR/RID 3 F1 IMDG 3/II

R: 11-38-48/20-51/53-62-65-67





Assay (by GC)	min. 95%	
Boiling Range (initial to dry point)	68.0-70.0°C	
Color (APHA)	max. 10	
Recorded Boiling Point	68.7°C	
Residue after Evaporation	max. 5 ppm	
Sulfur Compounds (as S)	max. 0.005%	
Thiophene	passes test	
Water (H ₂ O)	max. 0.05%	
Water Soluble Titrable Acid (meq/g)	max. 0.0003	
Physical Data (not specifications):		
Density (g/ml) at 25°C	0.687	
Ultraviolet Absorbance (1.00-cm path vs water:		

curve smooth throughout stated range with no extraneous impurity peaks):

at 210 nm	max. 1.00
at 220 nm	max. 0.20
at 230 nm	max. 0.10
at 240 nm	max. 0.04
at 250 nm	max. 0.02
at 280-400 nm	max. 0.01

BOX 8710.2500 2.5 |

Mallinckrodt Baker's cGMP Manufactured Chemicals for the Biopharmaceutical industry are a necessity for uncomplicated scale-up.

See chapter 6 of this catalogue.

Hexane

95% n-Hexane / 'BAKER ANALYZED'

8044

min. 9



M = 86.18 g/mol

1 l = 0.66 kg

FLASHPOINT $-23 \, ^{\circ}\text{C}$

CAS NO. 110-54-3

EINECS 203-777-6

NC CODE 2901 10 00

EC NO. 601 037 00 0

UN/ID NO. 1208

ADR/RID 3F1 IMDG 3/II

R: 11-38-48/20-51/53-62-65-67

S: 16-29-33-36/37-61-62-9







Assay (by GC)	min. 95%
Aromatic Compounds (as C ₆ H ₆)	max. 0.02%
Color (APHA)	max. 10
Density (g/ml) at 25°C	max. 0.657
Residue after Evaporation	max. 0.001%
Sulfur Compounds (as S)	max. 0.005%
Thiophene	passes test
Water Soluble Titrable Acid (meq/g)	max. 0.0003
Trace Impurities (in ppm):	
Aluminium (AI)	max. 0.5
Barium (Ba)	max. 0.1
Boron (B)	max. 0.02
Cadmium (Cd)	max. 0.05
Calcium (Ca)	max. 0.5
Chromium (Cr)	max. 0.02
Cobalt (Co)	max. 0.02
Copper (Cu)	max. 0.02
Iron (Fe)	max. 0.1
Lead (Pb)	max. 0.1
Magnesium (Mg)	max. 0.1
Manganese (Mn)	max. 0.02
Nickel (Ni)	max. 0.02
Tin (Sn)	max. 0.1
Zinc (Zn)	max. 0.1

PRODUCT	PACKING	CONT.
NO.		вох
8044.1000	11	6
8044.2500	2.5	4
8044.5000	5 EcoTainer	4
8044.9025	25 l	
8044.9200	200	
market and Ale	a secretarily and the second secretarily	and the state of the

EcoTainer, the metal solvent can for more safety in the lab. For safe handling of 25 l tin cans, see Self-closing tap.

Hexane

95% n-Hexane / 'BAKER'

8669

► CH₃(CH₂)₄CH₃

M = 86.18 g/mol

1 l = 0.66 kg

FLASHPOINT $-23 \, ^{\circ}\text{C}$

CAS NO. 110-54-3 **EINECS** 203-777-6

NC CODE 2901 10 00

EC NO. 601 037 00 0

UN/ID NO. 1208

ADR/RID 3F1

IMDG 3/II

R: 11-38-48/20-51/53-62-65-67

S: 16-29-33-36/37-61-62-9









95%	PRODUCT	PACKING	CONT
	NO.		B0
	8669.1000	11	(
	8669.2500	2.5	
	8669.5000	5 EcoTainer	
	8669.9025	25 l	
	8669.9200	200 l	
	EcoTainor th	a motal solvent can for more safet	ty in the lab

For safe handling of 25 I tin cans, see Self-closing tap.

iso-Hexane

'BAKER HPLC ANALYZED' / for use in High Performance Liquid Chromatography

Assay (by GC)

9305

► C₆H₁₄

M = 86.18 g/mol

1 l = 0.66 kg

FLASHPOINT −23 °C

CAS NO. 64742-49-0 NC CODE 2901 10 00

UN/ID NO. 1208

ADR/RID 3 F1

IMDG 3/II

R: 11-38-48/20-51/53-62-65-67

S: 16-29-33-36/37-61-62-9



environment





Assay (as C ₆ H ₁₄ isomers)	min. 95%	
n-Hexane	max. 3%	
Residue after Evaporation	max. 0.0005%	
Substances Darkened by H ₂ SO ₄	passes test	
Water (H ₂ O)	max. 0.01%	
Ultraviolet Absorbance (1.00-cm path vs water):		
at 210 pm	may 0.2	

at 210 nm	max. 0.2
at 220 nm	max. 0.08
at 254 nm	max. 0.008
at 280 nm	max. 0.005
at 350 nm	max. 0.005
UV Cut-off, nm	max. 192

PRODUCT	PACKING	CONT.
NO.		вох
9305.2500	2.5	4
9305.5000	5 EcoTainer	
FeeTelman 4h	a santal california con for social	andaturia tha lab

EcoTainer, the metal solvent can for more safety in the lab.

Filtered through a 0.2 micron filter. Packaged under Nitrogen.

Hexan

iso-Hexane

9267 'BAKER ANALYZED' / For Hydrocarbon Oil Index determination

C₆H₁₄

M = 86.18 g/mol

1 l = 0.66 kg

FLASHPOINT $-23 \, ^{\circ}\text{C}$

CAS NO. 64742-49-0 NC CODE 2901 10 00

UN/ID NO. 1208

ADR/RID 3 F1

IMDG 3/II

R: 11-38-48/20-51/53-62-65-67

S: 16-29-33-36/37-61-62-9








highly flammable

Color (APHA) max. 10 Residue after Evaporation max. 5 ppm Water (H₂O) max. 0.01%

Hydrocarbon oil index concentration (as RIVM oil reference standard):

Total peaks between n-decane (C₁₀H₂₂)

and n-tetracontane (C₄₀H₈₂) max. 0.5 mg/l



Suitable for determination of Hydrocarbon Oil Index according to ISO 9377-2.

iso-Hexane

environment

'BAKER ANALYZED

8043

C₆H₁₄

M = 86.18 g/mol

11 = 0.66 kg

FLASHPOINT −23 °C

CAS NO. 64742-49-0

NC CODE 2901 10 00

UN/ID NO. 1208

ADR/RID 3F1

IMDG 3/II

R: 11-38-48/20-51/53-62-65-67

S: 16-29-33-36/37-61-62-9







Assay (as C ₆ H ₁₄ isomers)	min. 95%
Aromatic Compounds (as C ₆ H ₆)	max. 0.01%
Boiling Range (initial to dry point)	56 - 63°C
Color (APHA)	max. 10
n-Hexane	max. 3%
Residue after Evaporation	max. 0.001%
Sulfur Compounds (as S)	max. 0.005%
Thiophene	passes test
Water (H ₂ O)	max. 0.01%
Water Soluble Titrable Acid (meq/g)	max. 0.0003
Trace Impurities (in ppm):	
Aluminium (AI)	max. 0.5
Barium (Ba)	max. 0.1
Boron (B)	max. 0.02
Cadmium (Cd)	max. 0.05
Calcium (Ca)	max. 0.5
Chromium (Cr)	max. 0.02
Cobalt (Co)	max. 0.02
Copper (Cu)	max. 0.02
Iron (Fe)	max. 0.1
Lead (Pb)	max. 0.1
Magnesium (Mg)	max. 0.1
Manganese (Mn)	max. 0.02

max. 0.02

max. 0.1

max. 0.1

PRODUCT	PACKING	CONT
NO.		B0)
8043.2500	2.5	4
8043.9200	200 l	

1-Hexanesulfonic Acid Sodium Salt

BAKER HPLC ANALYZED

► CH₃(CH₂)₅SO₃Na

M = 188.22 g/mol

CAS NO. 2832-45-3 EINECS 220-601-3

NC CODE 2904 10 00

For Ion-Pair Chromatography of Basic Compounds Assay (acidimetric) min. 98.0% UV Absorbance of 0.25M Solution:

at 200 nm max. 0.2 at 210 nm max. 0.08 at 220 nm max. 0.06 at 230 nm max 0.05 at 240 nm max. 0.05 at 250 nm max. 0.05

PACKING	CONT.
	вох
25 g Glass	



L-Histidine

See for qualities meeting the requirements of cGMP, USP/NF, BP/Ph.Eur, JP or FCC, page 36

Nickel (Ni)

Tin (Sn)

Zinc (Zn)



CONT.

вох

CONT.

CONT.

BOX

6

BOX

6

HYDRA

L(+)-Histidine
'BAKER ANALYZED' Biochemical 1603 NHCH:NCH:CCH2CH(NH2)COOH 98.5-101.5% PRODUCT PACKING Assay (dried basis) **M** = 155.16 g/mol Arsenic (As) max. 0.0003% **CAS NO.** 71-00-1 Ash (sulfated) max. 0.1% 1603.0025 25 g Glass **EINECS** 200-745-3 Heavy Metals (as Pb) max. 0.002% NC CODE 2933 29 90 Iron (Fe) max. 0.003% Loss on Drying at 105°C max. 0.2% Specific Rotation [α] 2 5 (dried basis, c=1in 6N HCI) $+ 12.6^{\circ} \text{ to } + 14.0^{\circ}$ L-Histidine Monohydrochloride nents of cGMP, USP/NF, BP/Ph.Eur, JP or FCC, page 36 Homidium bromide **HPLC Columns** See for detailed information section Chromatography products, page 448 **HPLC Solvents and Reagents** See for detailed information section HPLC Solvents and Reagents, page 64 Hyamine 1622 0.004 mol/I / 'BAKER ANALYZED' / for determination of anionic tensides 7255 **CAS NO.** 121-54-0 0.0039-0.0041 **PRODUCT PACKING** Titer (mol/l) **EINECS** 204-479-9 NO. NC CODE 2923 90 00 7255.1000 11 7255.2500 2.5 I Volumetric Solution, ready for use. **HYDRA-POINT Buffer** 8899 'BAKER ANALYZED' / Buffer for stabilizing pH during Karl Fischer titration **1 l** = 1.0 kg Buffer capacity min. 5 mmol acid/ml **PRODUCT PACKING** NC CODE 3822 00 00 8899.0500 UN/ID NO. 3286 500 ml ADR/RID 3 FTC 8899.1000 11 IMDG 3/I R: 11-23/24/25-34-39/23/24/25-63 Contains: imidazole, sulfur dioxide, methanol. **S:** 16-26-36/37/39-45-7/9

HYDRA-POINT Composite	2
-----------------------	---

BAKER AWALTZED / Fyndine nee one-component reagent for volumetric kan Fischer thration			0071		
1 l = 1.1 kg		vith HYDRA-POINT Methanol dry	PRODUCT	PACKING	CONT.
NC CODE 3822 00	00 product no. 8	8898	NO.		вох
UN/ID NO. 2922	Titer (mg H ₂ O/	(ml)1) min. 2.00	8891.1000	11	6
ADR/RID 8 CT1	1) At time of lot re	elease.	8891.2500	2.5	4
IMDG 8/III					
R: 20-34-6	3		Contains: in	nidazole, iodine, sulfur d	ioxide.
S: 26-36/3	7/39-45				
corrosive					

HYDRA

HYDRA-POINT Composite 5

8890 'BAKER ANALYZED' / Pyridine free one-component reagent for volumetric Karl Fischer 1

1 l = 1.2 kg NC CODE 3822 00 00 UN/ID NO. 2922 ADR/RID 8 CT1 IMDG 8/III R: 20-34-63

To be used with HYDRA-POINT Methanol dry product no. 8898 $\label{eq:product} \mbox{Titer (mg H}_2\mbox{O/ml})^1) \qquad \mbox{min. 5.00}$ $^1) At time of lot release.$

 PRODUCT NO.
 PACKING
 CONT.

 8890.1000
 1 l
 6

 8890.2500
 2.5 l
 4

Contains: imidazole, iodine, sulfur dioxide.



S: 26-36/37/39-45

HYDRA-POINT Composite 5K

8892 (BAKER ANALYZED) / Pyridine free one-component reagent for volumetric Karl Fischer titration of ketones

NC CODE 3822 00 00 UN/ID NO. 1760 ADR/RID 8 C9 IMDG 8/III R: 20-34-63 Suitability for coulometric Karl Fischer

Suitability for coulometric Karl Fischer

 PRODUCT
 PACKING
 CONT.

 NO.
 BOX

 8892.1000
 1 l
 6

 8892.2500
 2.5 l
 4

Contains: imidazole, iodine, sulfur dioxide

PACKING

corrosive

S: 26-36/37/39-45

HYDRA-POINT Coulometric An

BAKER ANALYZED' / Anolyte for coulometric Karl Fischer determinations with diafragr

1 I = 1.20 kg

FLASHPOINT <21 °C
NC CODE 3822 00 00

UN/ID NO. 3286

ADR/RID 3/II

R: 11-23/24/25-3

R: 11-23/24/25-34-39-40 **S:** 26-36/37/39-45

F F

highly flammable



passes test

8862.0500 500 ml

BOX

Contains: methanol, imidazole, sulphur dioxide, iodine, chloroform.

8863

HYDRA-POINT Coulometric Cat

'BAKER ANALYZED' / Catholyte for coulometric Karl Fischer determinations with diafragm

R: 11-23/24/25-34-39/23/24/25-40-48/23-59-63

S: 26-36/37/39-45

dangerou

highly flammable



 PRODUCT NO.
 PACKING PACKING
 CONT. PACKING PACKING

 8863.0025
 25 ml
 6

Contains: methanol, imidazole, sulphur dioxide, diethanolamine, tetrachloromethane.

Innovation is principal to our business.

HYDRA-POINT Coulometric Gen

'BAKER ANALYZED' / Halogen free reagent for coulometric Karl Fischer determinations without diafragm

titration

titration

Suitability for coulometric Karl Fischer

passes test

passes test

max. 0.01%

passes test

PRODUCT

8898.1000

8898.2500

NO.

PACKING

11

2.5 |

8860

1 l =	0.98 kg
FLASHPOINT	<21°0

NC CODE 3822 00 00

UN/ID NO. 3286 ADR/RID 3 FTC

IMDG 3/II

R: 11-23/24/25-34-39/23/24/25-63

S: 26-36/37/39-45





PRODUCT	PACKING	CONT
NO.		ВОХ
8860.0500	500 ml	ϵ

Contains: methanol, 1-hexanol, imidazole, sulphur dioxide, iodide.

HYDRA-POINT Coulometric Oven

'BAKER ANALYZED' / Halogen free reagent for coulometric Karl Fischer determinations without diafragm

Suitability for coulometric Karl Fischer

To be used with HYDRA-POINT Composite 2 product no. 8891 or HYDRA-POINT Composite 5

To be used with HYDRA-POINT Titrant 2 product

no. 8845 or HYDRA-POINT Titrant 5 product no.

Suitability for Karl Fischer titration

8861

11 = 1.04 kg FLASHPOINT <21 °C NC CODE 3822 00 00

UN/ID NO. 3286 ADR/RID 3 FTC IMDG 3/II

R: 11-23/24/25-34-39/23/24/25-63

S: 26-36/37/39-45





PRODUCT PACKING CONT. вох NO. 500 ml 8861.0500

Contains: methanol, imidazole, sulphur dioxide, iodine, propylene glycol.

HYDRA-POINT Methanol Dry

orking medium for Karl Fischer titration 'BAKER ANALYZED' / Dry methanol is a

product no. 8890

Water (H₂O)

8898

CONT

вох

4

► CH₃OH

M = 32.04 g/mol

1 I = 0.79 kg

FLASHPOINT 11 °C **CAS NO.** 67-56-1

EINECS 200-659-6

NC CODE 3822 00 00

EC NO. 603 001 00 0

UN/ID NO. 1230

ADR/RID 3 FT1

IMDG 3/II

R: 11-23/24/25-39/23/24/25

S: 16-36/37-45-7





HYDRA-POINT Solvent G

'BAKER ANALYZED' / Pyridine free solvent for volumetric Karl Fischer titration

8844

8855

11 = 0.93 kg $\textbf{FLASHPOINT} \quad 11~^{\circ}\text{C}$ NC CODE 3822 00 00 **UN/ID NO.** 1992 ADR/RID 3 FT1

IMDG 3/II

R: 11-23/24/25-36/38-39/23/24/25-63

S: 16-26-36/37/39-45-7/9





PRODUCT NO.	PACKING	CONT. BOX
8855.1000	11	6
8855.2500	2.5	4

Contains: imidazole, sulfur dioxide, methanol.



HYDRA

HYDRA-POINT Titrant 2 mg H₂0/ml

BAKER ANALYZED' / Titrant for volumetric Karl Fischer titration

FLASHPOINT 11 °C **UN/ID NO.** 1992 ADR/RID 3 FT1 IMDG 3/II

To be used with HYDRA-POINT Solvent K product no. 8840 or HYDRA-POINT Solvent G product no. 8855

PRODUCT PACKING CONT. вох 8845.1000 11 8845.2500 2.5 l

R: 11-23/24/25-39/23/24/25

S: 16-36/37-45-7/9



Titer (mg H₂O/ml)¹) min. 2.00 1) At time of lot release

Contains: methanol, iodine.

HYDRA-POINT Titrant 5 mg H₂0/ml

'BAKER ANALYZED' / Titrant for volumetric Karl Fischer titration

1) At time of lot release

FLASHPOINT 11 °C NC CODE 3822 00 00 **UN/ID NO.** 1992 ADR/RID 3 FT1 IMDG 3/II

To be used with HYDRA-POINT Solvent K product no. 8840 or HYDRA-POINT Solvent G product no. 8855 Titer (ma H₂O/ml)¹) min. 5.00

PRODUCT PACKING CONT. NO. BOX 8844 1000 11 6 8844.2500 2.5 |

R: 11-23/24/25-39/23/24/25

S: 16-36/37-45-7/9





Contains: methanol, iodine.

Hydrobromic Acid

6010 47-49% / 'BAKER ANALYZED' / ACS

► HBr M = 80.92 g/mol **1 I =** 1.49 kg **CAS NO.** 10035-10-6 **EINECS** 233-113-0 NC CODE 2811 19 10 EC NO. 35 002 01 8 UN/ID NO. 1788 ADR/RID 8 C1 IMDG 8/II R: 34-37 **S:** 26-45-7/9

Exceeds ACS Specifications	
Assay	47.0-49.0%
Chloride (CI)	max. 0.05%
lodide (I)	max. 0.003%
Phosphate (PO ₄)	max. 0.001%
Residue after Ignition	max. 0.002%
Sulfate and Sulfite (as SO ₄)	max. 0.003%
Trace Impurities (in ppm):	
Arsenic (As)	max. 0.5
Heavy Metals (as Pb)	max. 5
Iron (Fe)	max. 1
Selenium (Se)	max. 0.01

PRODUCI	PACKING	CONI
NO.		ВОХ
6010.1000	11	
6010.9200	200 l	

Use J.T.Baker Ultrex II and BAKER INSTRA-ANALYZED acids for low level trace element analysis.

See chapter 3 of this catalogue for more details.



Hydrochloric Acid

37-38% (max. 5 ppb Hg) / 'BAKER ANALYZED'

6081

► HCI	
M =	36.46 g/mol
1 =	1.19 kg
CAS NO.	7647-01-0
EINECS	231-595-7
NC CODE	2806 10 00
EC NO.	17 002 01 0
UN/ID NO.	1789
ADR/RID	8 C1
IMDG	8/11
R:	34-37
S:	26-45
	C corrosive

Assay	37-38%
Ammonium (NH ₄)	max. 3 ppm
Free Chlorine	max. 0.5 ppm
Mercury (Hg)	max. 5 ppb
Residue after Ignition	max. 5 ppm
Sulfate (SO ₄)	max. 1 ppm
Sulfite (SO ₃)	max. 2 ppm
Trace Impurities (in ppm):	
Aluminium (AI)	max. 0.05
Arsenic (As)	max. 0.05
Barium (Ba)	max. 0.02
Beryllium (Be)	max. 0.01
Cadmium (Cd)	max. 0.02
Calcium (Ca)	max. 0.1
Chromium (Cr)	max. 0.05
Cobalt (Co)	max. 0.01
Copper (Cu)	max. 0.01
Germanium (Ge)	max. 0.1
Iron (Fe)	max. 0.05
Lead (Pb)	max. 0.05
Lithium (Li)	max. 0.01
Magnesium (Mg)	max. 0.05
Manganese (Mn)	max. 0.01
Molybdenum (Mo)	max. 0.02
Nickel (Ni)	max. 0.05
Potassium (K)	max. 0.05
Silver (Ag)	max. 0.01
Sodium (Na)	max. 0.05
Strontium (Sr)	max. 0.01
Thallium (TI)	max. 0.05
Titanium (Ti)	max. 0.05
Vanadium (V)	max. 0.01
Zinc (Zn)	max. 0.5

max. 0.05

Zirconium (Zr)

PRODUCT	PACKING	CONT.
NO.		вох
6081.1000	11	6
6081.2500	2.5	4
6081.2500PE	2.5 I HDPE	4
6081.9025	25 l	
6081.9200	200	

Calibrate and standardise your analytical methods and equipment with J.T.Baker Volumetric and Buffer solutions.

Refer to the Analytical applications section of this catalogue for more details.



Hydro

Hydrochloric Acid

11 37% / 'BAKER ANALYZED' / ACS

M = 36.46 g/mol 1 I = 1.19 kg CAS NO. 7647-01-0 EINECS 231-595-7 NC CODE 2806 10 00 EC NO. 17 002 01 0 UN/ID NO. 1789 ADR/RID 8 C1 IMDG 8/II R: 34-37 S: 26-45

Specifications for testing USP/NF monographs Appearance passes test Bromide (Br) max. 0.005% Color (APHA) max. 10 Extractable Organic Substances max. 5 ppm Free Chlorine max. 1 ppm Residue after Ignition max. 3 ppm 1.185-1.192 Specific Gravity at 60°/60°F Trace Impurities (in ppb): Aluminium (AI) max. 100 Arsenic and Antimony (as As) max. 5 Boron (B) max. 50 Calcium (Ca) max. 200 Chromium (Cr) max. 100 Copper (Cu) max. 100

max. 100

max. 100

max. 100

max. 50

max. 300

max. 300

max. 100

max. 300

max. 300

max. 300

max. 300

max. 100

max. 3

max. 1

max. 0.5

max. 0.8

max. 5

Meets ACS Specifications. Meets Reagent

PRODUCT	PACKING	CONT.
NO.		BOX
6011.0500	500 ml	6



Hydrochloric Acid 37% MOS, VLSI Grade

ee for detailed information section Microelectronic materials page 32 or the website www.jtbaker.com/micro

Trace Impurities (in ppm): Ammonium (NH₄)

Phosphate (as PO₄)

Sulfate (SO₄)

Sulfite (SO₃)

Gold (Au)

Iron (Fe)

Lead (Pb)

Heavy Metals (as Pb)

Magnesium (Mg)

Manganese (Mn)

Mercury (Hg)

Potassium (K)

Sodium (Na)

Titanium (Ti)

Tin (Sn)

Zinc (Zn)

Nickel (Ni)

The J.T.Baker CYCLE-TAINER
High Purity Solvent Delivery System,
preserves purity and protects people.

See chapter 3 of this catalogue for product details.

Hydrochloric Acid36.5-38% / 'BAKER INSTRA-ANALYZED' / for Trace Metal Analysis / ACS

► HCI

M = 36.46 g/mol

1 l = 1.19 kg

CAS NO. 7647-01-0 **EINECS** 231-595-7

NC CODE 2806 10 00

EC NO. 17 002 01 0

UN/ID NO. 1789

ADR/RID 8 C1 IMDG 8/II

R: 34-37

S: 26-45



Meets ACS Specifications	
Assay (by acid-base titration)	36.5-38.0%
Appearance	passes test
Bromide (Br)	max. 0.005%
Color (APHA)	max. 10
Extractable Organic Substances	max. 5 ppm
Free Chlorine	max. 0.5 ppm
Residue after Ignition	max. 3 ppm
Specific Gravity at 60°/60°F	1.185-1.192
Trace Impurities (in ppb):	
Aluminium (AI)	max. 10
Arsenic and Antimony (as As)	max. 5
Barium (Ba)	max. 1
Beryllium (Be)	max. 1
Bismuth (Bi)	max. 10
Boron (B)	max. 20
Cadmium (Cd)	max. 1
Calcium (Ca)	max. 50
Chromium (Cr)	max. 1
Cobalt (Co)	max. 1
Copper (Cu)	max. 1
Gallium (Ga)	max. 1
Germanium (Ge)	max. 3
Gold (Au)	max. 4
Heavy Metals (as Pb)	max. 100
Iron (Fe)	max. 15
Lead (Pb)	max. 1
Lithium (Li)	max. 1
Magnesium (Mg)	max. 10
Manganese (Mn)	max. 1
Mercury (Hg)	max. 0.5
Molybdenum (Mo)	max. 10
Nickel (Ni)	max. 4
Niobium (Nb)	max. 1
Potassium (K)	max. 9
Selenium (Se)	act. value reported
Silicon (Si)	max. 100
Silver (Ag)	max. 1
Sodium (Na)	max. 100
Strontium (Sr)	max. 1
Tantalum (Ta)	max. 1

Thallium (TI)	max. 5
Tin (Sn)	max. 5
Titanium (Ti)	max. 1
Vanadium (V)	max. 1
Zinc (Zn)	max. 5
Zirconium (Zr)	max. 1
Trace Impurities (in ppm):	
Ammonium (NH ₄)	max. 3
Phosphate (PO ₄)	max. 0.05
Sulfate (SO ₄)	max. 0.5
Sulfite (SO ₃)	max. 0.8

PRODUCT	PACKING	CONT.
NO.		вох
9530.0500	500 ml	6
9530.2500	2.5	4

Hydrochloric Acid

35-39% / 'BAKER'

► HCI

M = 36.46 g/mol

1 l = 1.19 kg

CAS NO. 7647-01-0

EINECS 231-595-7

NC CODE 2806 10 00

EC NO. 17 002 01 0

UN/ID NO. 1789

ADR/RID 8 C1

IMDG 8/II

R: 34-37

S: 26-45



Assay	35.0-39.0%
Appearance of solution	passes test
Free Chlorine	max. 4 ppm
Heavy Metals (as Pb)	max. 2 ppm
Identification	passes test
Residue after Evaporation	max. 0.01%
Sulfates (as SO ₄)	max. 20 ppm

PRODUCT	PACKING	CONT.
NO.		вох
6012.1000	11	
6012.2500	2.5	4
6012.5000	5 l	4
6012 9025	25	

6012



Hydro

Hydrochloric Acid

33-36% / ULTREX II Ultrapure Reagent

M = 36.46 g/mol
1 I = 1.19 kg
CAS NO. 7647-01-0
EINECS 231-595-7
NC CODE 2806 10 00
EC NO. 17 002 01 0
UN/ID NO. 1789
ADR/RID 8 C1
IMDG 8/II
R: 34-37



S: 26-45

Certificate Provided Reporting Actual Lot Analysis		
Assay	33-36% (w/w)	
Trace Impurities (in ppt) (pg/g):		
Aluminium (AI)	max. 20	
Antimony (Sb)	max. 20	
Arsenic (As)	max. 50	
Barium (Ba)	max. 10	
Beryllium (Be)	max. 10	
Bismuth (Bi)	max. 10	
Boron (B)	max. 100	
Cadmium (Cd)	max. 10	
Calcium (Ca)	max. 20	
Cerium (Ce)	max. 10	
Cesium (Cs)	max. 10	
Chromium (Cr)	max. 20	
Cobalt (Co)	max. 10	
Copper (Cu)	max. 20	
Dysprosium (Dy)	max. 1	
Erbium (Er)	max. 1	
Europium (Eu)	max. 1	
Gadolinium (Gd)	max. 1	
Gallium (Ga)	max. 10	
Gold (Au)	max. 100	
Hafnium (Hf)	max. 10	
Holmium (Ho)	max. 1	
Indium (In)	max. 1	
Iron (Fe)	max. 20	
Lanthanum (La)	max. 1	
Lead (Pb)	max. 10	
Lithium (Li)	max. 10	
Lutetium (Lu)	max. 10	
Magnesium (Mg)	max. 10	
Manganese (Mn)	max. 10	
Mercury (Hg)	max. 100	
Molybdenum (Mo)	max. 10	
Neodymium (Nd)	max. 1	
Nickel (Ni)	max. 50	
Niobium (Nb)	max. 1	
Palladium (Pd)	act. value reported	
Platinum (Pt)	act. value reported	
Potassium (K)	max. 10	

Praseodymium (Pr)	max. 1
Rhenium (Re)	max. 10
Rhodium (Rh)	max. 10
Rubidium (Rb)	max. 10
Ruthenium (Ru)	max. 10
Samarium (Sm)	max. 1
Scandium (Sc)	max. 10
Selenium (Se)	act. value reported
Silver (Ag)	max. 10
Sodium (Na)	max. 10
Strontium (Sr)	max. 10
Tantalum (Ta)	act. value reported
Tellurium (Te)	max. 1
Terbium (Tb)	max. 1
Thallium (TI)	max. 10
Thorium (Th)	max. 1
Thulium (Tm)	max. 1
Tin (Sn)	max. 20
Titanium (Ti)	max. 20
Tungsten (W)	max. 10
Jranium (U)	max. 1
Vanadium (V)	max. 10
Ytterbium (Yb)	max. 1
Yttrium (Y)	max. 1
Zinc (Zn)	max. 20
Zirconium (Zr)	max. 10

PACKING	CONT.
	вох
500 ml Fluoropolymer, pre-leached	
1 l Fluoropolymer, pre-leached	
2 l Fluoropolymer, pre-leached	
	500 ml Fluoropolymer, pre-leached 1 l Fluoropolymer, pre-leached

Mallinckrodt Baker's cGMP Manufactured Chemicals for the Biopharmaceutical industry are a necessity for uncomplicated scale-up.

See chapter 6 of this catalogue.



Hydrochloric Acid 32% / 'BAKER ANALYZED'

6070

► HCI	
M =	36.46 g/mol
1 =	1.16 kg
CAS NO.	7647-01-0
EINECS	231-595-7
NC CODE	2806 10 00
EC NO.	17 002 01 0
UN/ID NO.	1789
ADR/RID	8 C1
IMDG	8/11
R:	34-37
S:	26-45
	C C

Assay	min. 32%
Ammonium (NH ₄)	max. 1 ppm
Free Chlorine	max. 0.5 ppm
Residue after Ignition	max. 5 ppm
Sulfate (SO ₄)	max. 1 ppm
Sulfite (SO ₃)	max. 2 ppm
Trace Impurities (in ppm):	
Aluminium (AI)	max. 0.5
Arsenic (As)	max. 0.1
Barium (Ba)	max. 0.1
Beryllium (Be)	max. 0.1
Cadmium (Cd)	max. 0.1
Calcium (Ca)	max. 0.5
Chromium (Cr)	max. 0.1
Cobalt (Co)	max. 0.1
Copper (Cu)	max. 0.1
Germanium (Ge)	max. 0.1
Iron (Fe)	max. 0.5
Lead (Pb)	max. 0.1
Lithium (Li)	max. 0.1
Magnesium (Mg)	max. 0.5
Manganese (Mn)	max. 0.1
Molybdenum (Mo)	max. 0.1
Nickel (Ni)	max. 0.1
Potassium (K)	max. 0.1
Silver (Ag)	max. 0.1
Sodium (Na)	max. 0.5
Strontium (Sr)	max. 0.1
Thallium (TI)	max. 0.1
Titanium (Ti)	max. 0.1
Vanadium (V)	max. 0.1
Zinc (Zn)	max. 0.5
Zirconium (Zr)	max. 0.1

PRODUCT	PACKING	CONT.
NO.		вох
6070.1000	11	
6070.2500	2.5	4
6070.9025	25 l	
6070.9200	200 l	

Use J.T.Baker Ultrex II and BAKER INSTRA-ANALYZED acids for low level trace element analysis.

See chapter 3 of this catalogue for more details.



Hydro

Hydrochloric Acid

9600 30% / 'BAKER INSTRA-ANALYZED' / for Trace Metal Analysis

M = 36.46 g/mol
1 I = 1.19 kg
CAS NO. 7647-01-0
EINECS 231-595-7
NC CODE 2806 10 00
EC NO. 17 002 01 0
UN/ID NO. 1789
ADR/RID 8 C1
IMDG 8/II
R: 34-37
S: 26-45

C

Actual lot analysis at time of lot	release
Assay	29-31%
Color	max. 10 Hazen
Trace Impurities:	
Aluminium (AI)	max. 1 ppb
Ammonium (NH ₄)	max. 0.5 ppm
Antimony (Sb)	max. 1 ppb
Arsenic (As)	max. 1 ppb
Barium (Ba)	max. 1 ppb
Beryllium (Be)	max. 1 ppb
Bismuth (Bi)	max. 1 ppb
Boron (B)	max. 1 ppb
Bromide (Br)	max. 10 ppm
Cadmium (Cd)	max. 1 ppb
Calcium (Ca)	max. 1 ppb
Chromium (Cr)	max. 1 ppb
Cobalt (Co)	max. 1 ppb
Copper (Cu)	max. 1 ppb
Free Chlorine	max. 0.5 ppm
Gallium (Ga)	max. 1 ppb
Germanium (Ge)	max. 1 ppb
Gold (Au)	max. 1 ppb
Indium (In)	max. 1 ppb
Iron (Fe)	max. 1 ppb
Lead (Pb)	max. 1 ppb
Lithium (Li)	max. 1 ppb
Magnesium (Mg)	max. 1 ppb
Manganese (Mn)	max. 1 ppb
Mercury (Hg)	max. 1 ppb
Molybdenum (Mo)	max. 1 ppb
Nickel (Ni)	max. 1 ppb
Potassium (K)	max. 1 ppb
Selenium (Se)	max. 1 ppb
Silver (Ag)	max. 1 ppb
Sodium (Na)	max. 1 ppb

max. 1 ppb

max. 1 ppb

max. 1 ppb max. 1 ppb

max. 1 ppb

max. 0.01 ppm

Total Sulfur		max. 0.3 ppm
Uranium (U)		max. 1 ppb
Vanadium (V)		max. 1 ppb
Zinc (Zn)		max. 1 ppb
Zirconium (Zr	٠١	max. 1 ppb
Ziicomum (Zi)	тих. тррь
PRODUCT	PACKING	CONT. BOX
PRODUCT NO.	,	CONT.

Calibrate and standardise your analytical methods and equipment with J.T.Baker Volumetric and Buffer solutions.

Strontium (Sr)

Thallium (TI)

Thorium (Th)

Titanium (Ti)

Total Phosphorus

Refer to the Analytical applications section of this catalogue for more details.



6167

Hydrochloric Acid

/ 'BAKER ANALYZED

► HCI **M** = 36.46 g/mol **1 l** = 1.12 kg **CAS NO.** 7647-01-0 **EINECS** 231-595-7 NC CODE 2806 10 00 **EC NO.** 17 002 01 0 **UN/ID NO.** 1789

> ADR/RID 8 C1 IMDG 8/II R: 34-37 **S:** 26-45

min. 25% Assay Ammonium (NH₄) max. 3 ppm Free Chlorine max. 0.5 ppm Residue after Ignition max. 5 ppm Sulfate (SO₄) max. 1 ppm Sulfite (SO₃) max. 2 ppm Trace Impurities (in ppm): Aluminium (AI) max. 0.5 Arsenic (As) max. 0.1 max. 0.1 Barium (Ba) Beryllium (Be) max. 0.1 Cadmium (Cd) max. 0.1 Calcium (Ca) max. 0.5 max. 0.1 Chromium (Cr) Cobalt (Co) max. 0.1 Copper (Cu) max. 0.1

Germanium (Ge) max. 0.1 Iron (Fe) max. 0.5 Lead (Pb) max. 0.1 Lithium (Li) max. 0.1 Magnesium (Mg) max. 0.5 Manganese (Mn) max. 0.1 Molybdenum (Mo) max. 0.1 Nickel (Ni) max. 0.1 Potassium (K) max. 0.1 Silver (Ag) max. 0.1 Sodium (Na) max. 0.5 Strontium (Sr) max. 0.1 Thallium (TI) max. 0.1 Titanium (Ti) max. 0.1 Vanadium (V) max. 0.1

max. 0.5

max. 0.1

18-19%

PRODUCT

7655.1000PE 1 | HDPE

PACKING

PRODUCT	PACKING	CONT. BOX
6167.2500	2.5	
6167.9025	25 l	

Hydrochloric Acid 18.5% VLSI

for detailed information section Microelectronic materials page 32 or the website www.jtbaker.com/micro

Zinc (Zn)

Assav

Zirconium (Zr)



7655

CONT.

вох

6

Hydrochloric Acid 18% / 'BAKER'

M = 36.46 g/mol**CAS NO.** 7647-01-0 EINECS 231-595-7 NC CODE 2806 10 00 **UN/ID NO.** 1789 ADR/RID 8 C1 IMDG 8/II

R: 36/37/38 **S:** 26-45

Find more Chromatography information at www.jtbaker.com/chromatography



Hydro

Hydrochloric Acid

7219 10% / 'BAKER

M = 36.46 g/mol
CAS NO. 7647-01-0
EINECS 231-595-7
NC CODE 2806 10 00
EC NO. 17 002 01 0

UN/ID NO. 1789 **ADR/RID** 8 C1 **IMDG** 8/III

R: 36/37/38 **S:** 26-45

 Assay
 9.5 - 10.5%

 Appearance
 passes test

 Free Chlorine
 max. 1 ppm

 Heavy Metals (as Pb)
 max. 2 ppm

 Identification
 passes test

 Residue after Evaporation
 max. 0.01%

 Sulfates (as SO₄)
 max. 5 ppm

 PRODUCT NO.
 PACKING PACKING
 CONT. BOX

 7219.9010
 10 | Polycube

 7219.9060
 60 |

Hydrochloric Acid

7212 3.57 mol/l, 1/0.28 mol/l / 'BAKER ANALYZED'

Titer (mol/l)

Titer (mol/l)

► HCI

M = 36.46 g/mol

CAS NO. 7647-01-0

EINECS 231-595-7 NC CODE 2806 10 00 UN/ID NO. 1789

ADR/RID 8 C1

R: 36/37/38 **S:** 26-45

*

Volumetric Solution, ready for use.

Hydrochloric Acid

7230 2 mol/l / 'BAKER ANALYZED

▶ HCI

M = 36.46 g/molCAS NO. 7647-01-0EINECS 231-595-7NC CODE 2806 10 00

1.99-2.01

 PRODUCT
 PACKING
 CONT.

 NO.
 BOX

 7230.1000
 1 l
 6

Volumetric Solution, ready for use.

Hydrochloric Acid

7088 1 mol/l / 'BAKER ANALYZED'

M = 36.46 g/mol
CAS NO. 7647-01-0
EINECS 231-595-7
NC CODE 2806 10 00

 PRODUCT
 PACKING
 CONT.

 NO.
 BOX

 7088.1000
 1 l
 6

 7088.9010
 10 l Polycube

 7088.9020
 20 l Polycube

Volumetric Solution, ready for use.

Each lot of this product is standardized potentiometrically against NIST tracable reference standard.

Hydrochloric Acid

7111 0.5 mol/l / 'BAKER ANALYZED'

► HCI

M = 36.46 g/mol

CAS NO. 7647-01-0

EINECS 231-595-7

NC CODE 2806 10 00

 Titer (mol/I)
 0.4975-0.5025

 Trace Impurities (in ppm):

 Ammonium (NH4)
 max. 3

 Heavy Metals (as Pb)
 max. 0.5

 Iron (Fe)
 max. 0.1

 Sulfate (SO4)
 max. 1

PRODUCT	PACKING	CONT.
NO.		вох
7111.1000	11	6
7111.9010	10 l Polycube	
7111.9020	20 l Polycube	

Volumetric Solution, ready for use.

Each lot of this product is standardized potentiometrically against NIST tracable reference standard.



Hydrochloric Acid 0.357 mol/l, 13.021 g HCl/l, 1/2.8 N / 'BAKER ANALYZED'

7477

► HCI	Titer (mol/l)	0.355-0.359	PRODUCT	PACKING	CONT.
M = 36.46 g/mol			NO.		вох
CAS NO. 7647-01-0			7477.9010	10 l Polycube	
EINECS 231-595-7					
NC CODE 2806 10 00			Volumetric S	Solution, ready for use.	

Hydrochloric Acid 0.2 mol/I / 'BAKER ANALYZED'

7651

► HCI	Titer (mol/l)	0.1990-0.2010	PRODUCT	PACKING	CONT.
M = 36.46 g/mol			NO.		вох
CAS NO. 7647-01-0			7651.9010	10 l	
EINECS 231-595-7					
NC CODE 2806 10 00			Volumetric S	olution, ready for use.	

Hydrochloric Acid

0.1 mol/l / 'BAKER ANALYZED'

7038

► HCI	Titer (mol/I)	0.0997-0.1003	PRODUCT	PACKING	CONT.
M = 36.46 g/mol	Trace Impurities (in ppm):		NO.		вох
CAS NO. 7647-01-0	Ammonium (NH ₄)	max. 3	7038.1000	11	6
EINECS 231-595-7	Heavy Metals (as Pb)	max. 0.5	7038.2500	2.5	4
NC CODE 2806 10 00	Iron (Fe)	max. 0.1	7038.9010	10 l Polycube	
	Sulfate (SO ₄)	max. 1	7038.9020P	20 l Polycube	
			Volumetric S	olution ready for use	

Hydrochloric Acid 0.0357 mol/l, 1/28 mol/l / 'BAKER ANALYZED'

► HCI	Titer (mol/l)	0.0350 - 0.0370	PRODUCT	PACKING	CONT.
M = 36.46 g/mol			NO.		вох
CAS NO. 7647-01-0			7211.9010	10 l Polycube	
EINECS 231-595-7			7211.9020	20 l Polycube	
NC CODE 2806 10 00					
			Volumetric S	Solution ready for use	

Hydrochloric Acid 0.01 mol/l / 'BAKER ANALYZED'

7493

► HCI	Titer (mol/l)	0.0095-0.0105	PRODUCT	PACKING	CONT.
CAS NO. 7647-01-0			NO.		вох
EINECS 231-595-7			7493.9020	20 l Polycube	
NC CODE 2806 10 00					
			Volumetric S	Solution, ready for use.	
			Each lot of the	his product is standardized	l potentiometrically
			against NIST	tracable reference standa	ırd.

Hydrochloric Acid 1 mol/I / 1 equiv. = 36.46g; 1N / DILUT-IT

	:, , ,,			
HCI		PRODUCT	PACKING	CONT
M =	: 36.46 g/mol	NO.		вох
CAS NO.	. 7647-01-0	4657	1 amp.	6
EINECS	231-595-7			
NC CODE	: 2806 10 00	Volumetric	Concentrate, for dilution to 1 l.	
UN/ID NO.	. 1789			
ADR/RID	8 C1			
IMDG	i 8/III			
R:	: 36/37/38			
S:	: 26-45			
	xi irritant			



Hydro

Hydrochloric Acid

4 0.5 mol/l / 1/2 equiv. = 18.23g; 0.5 N / DILUT-IT

M = 36.46 g/mol

CAS NO. 7647-01-0

EINECS 231-595-7

NC CODE 2806 10 00

UN/ID NO. 1789

UN/ID NO. 1789 ADR/RID 8 C1 IMDG 8/II

R: 36/37/38 **S:** 26-45

irritant

 PRODUCT NO.
 PACKING BOX

 4654
 1 amp.
 6

Volumetric Concentrate, for dilution to 11.

Hydrochloric Acid

4655 0.1 mol/l / DILUT-I

M = 36.46 g/mol
CAS NO. 7647-01-0
EINECS 231-595-7
NC CODE 2806 10 00
UN/ID NO. 1789
ADR/RID 8 C1

IMDG 8/III

 PRODUCT NO.
 PACKING PACKING
 CONT. BOX

 4655
 1 amp.
 6

Volumetric Concentrate, for dilution to 1 l.

Hydrochloric Acid

4867 0.01 mol/l / DILUT-IT

► HCI

M = 36.46 g/mol

EINECS 231-595-7

NC CODE 2806 10 00

 PRODUCT NO.
 PACKING PACKING
 CONT. BOX PACKING

 4867
 1 amp.
 6

Volumetric Concentrate, for dilution to 11.

Hydrochloric Acid

See for qualities meeting the requirements of cGMP, USP/NF, BP/Ph.Eur, JP or FCC, page 36

Hydrochloric Acid Solutions See for qualities meeting the requirements of

See for qualities meeting the requirements of cGMP, USP/NF, BP/Ph.Eur, JP or FCC, page 36

Hydrocortisone, Micronized

See for qualities meeting the requirements of cGMP, USP/NF, BP/Ph.Eur, JP or FCC, page 36

The J.T.Baker CYCLE-TAINER
High Purity Solvent Delivery System,
preserves purity and protects people.

See chapter 3 of this catalogue for product details.

Hydrofluoric Acid

/ ULTREX II Ultrapure Reagent

6904

>	HF	
	M =	20.01 g/mol
	11=	1.16 kg
	CAS NO.	7664-39-3
	EINECS	231-634-8
	NC CODE	2811 11 00
	EC NO.	9 003 00 1
	UN/ID NO.	1790
	ADR/RID	8 CT1
	IMDG	8/11
	D.	26/27/28 35





Certificate Provided Reportin	g Actual Lot Analysis
Assay	47-51% (w/w)
Trace Impurities (in ppt) (pg/	g):
Aluminium (AI)	max. 20
Antimony (Sb)	max. 20
Arsenic (As)	max. 50
Barium (Ba)	max. 10
Beryllium (Be)	max. 10
Bismuth (Bi)	max. 10
Boron (B)	max. 100
Cadmium (Cd)	max. 10
Calcium (Ca)	max. 20
Cerium (Ce)	max. 10
Cesium (Cs)	max. 10
Chromium (Cr)	max. 20
Cobalt (Co)	max. 10
Copper (Cu)	max. 20
Dysprosium (Dy)	max. 1
Erbium (Er)	max. 1
Europium (Eu)	max. 1
Gadolinium (Gd)	max. 1
Gallium (Ga)	max. 10
Germanium (Ge)	max. 10
Gold (Au)	max. 20
Hafnium (Hf)	max. 10
Holmium (Ho)	max. 1
Indium (In)	max. 1
Iron (Fe)	max. 20
Lanthanum (La)	max. 10
Lead (Pb)	max. 10
Lithium (Li)	max. 10
Lutetium (Lu)	max. 1
Magnesium (Mg)	max. 10
Manganese (Mn)	max. 10
Mercury (Hg)	max. 100
Molybdenum (Mo)	max. 10
Neodymium (Nd)	max. 1
Nickel (Ni)	max. 50
Niobium (Nb)	max. 10
Palladium (Pd)	max. 20
Platinum (Pt)	max. 20

PRODUCT PACKING	CONT.
Zirconium (Zr)	max. 10
Zinc (Zn)	max. 20
Yttrium (Y)	max. 1
Ytterbium (Yb)	max. 1
Vanadium (V)	max. 10
Uranium (U)	max. 1
Tungsten (W)	max. 20
Titanium (Ti)	max. 50
Tin (Sn)	max. 20
Thulium (Tm)	max. 1
Thorium (Th)	max. 1
Thallium (TI)	max. 10
Terbium (Tb)	max. 1
Tellurium (Te)	max. 1
Tantalum (Ta)	act. value reported
Strontium (Sr)	max. 10
Sodium (Na)	max. 10
Silver (Ag)	max. 10
Selenium (Se)	act. value reported
Scandium (Sc)	max. 10
Samarium (Sm)	max. 1
Ruthenium (Ru)	max. 20
Rubidium (Rb)	max. 20
Rhodium (Rh)	max. 20
Rhenium (Re)	max. 10
Praseodymium (Pr)	max. 1
Potassium (K)	max. 10

. Koboci	Tricking	COIL
NO.		B0
6904.0500	500 ml Fluoropolymer, pre-leached	
6904.1000	1 Fluoropolymer, pre-leached	

Mallinckrodt Baker's cGMP Manufactured Chemicals for the Biopharmaceutical industry are a necessity for uncomplicated scale-up.

See chapter 6 of this catalogue.



Hydro

Hydrofluoric Acid

48-51% / 'BAKER INSTRA-ANALYZED' / for Trace Metal Analysis / ACS

M = 20.01 g/mol
1 l = 1.16 kg
CAS NO. 7664-39-3
EINECS 231-634-8
NC CODE 2811 11 00
EC NO. 9 003 00 1
UN/ID NO. 1790
ADR/RID 8 CT1
IMDG 8/II
R: 26/27/28-35

9563

C C

S: 26-36/37-45-7/9



Meets ACS Specifications. Meets Reagent		
Specifications for testing USP/NF monographs		
Assay (acidimetric)	48.0-51.0%	
Color (APHA)	max. 10	
Fluorsilicic Acid (H ₂ SiF ₆)	max. 0.010%	
Residue after Ignition	max. 5 ppm	
Trace Impurities (in ppb):		
Aluminium (AI)	max. 50	
Arsenic and Antimony (as As)	max. 30	
Barium (Ba)	max. 50	
Beryllium (Be)	max. 20	
Bismuth (Bi)	max. 100	
Boron (B)	max. 50	
Cadmium (Cd)	max. 10	
Calcium (Ca)	max. 50	
Chromium (Cr)	max. 10	
Cobalt (Co)	max. 5	
Copper (Cu)	max. 5	
Gallium (Ga)	max. 20	
Germanium (Ge)	max. 50	
Gold (Au)	max. 10	
Heavy Metals (as Pb)	max. 100	
Iron (Fe)	max. 100	
Lead (Pb)	max. 20	
Lithium (Li)	max. 20	
Magnesium (Mg)	max. 50	
Manganese (Mn)	max. 50	
Mercury (Hg)	max. 10	
Molybdenum (Mo)	max. 50	
Nickel (Ni)	max. 10	
Potassium (K)	max. 100	
Silver (Ag)	max. 20	
Sodium (Na)	max. 200	
Strontium (Sr)	max. 20	
Tantalum (Ta)	max. 100	
Tin (Sn)	max. 50	
Titanium (Ti)	max. 50	
Vanadium (V)	max. 50	

max. 20

max. 20

Chloride (CI)		max.	
Nitrate (NO ₃)	max.	
Phosphate (PO ₄) Sulfate and Sulfite (as SO ₄)		max.	
		max. !	
PRODUCT	PACKING	CONT	
NO.		BO	
9563.0500	500 ml		



Hydrofluoric Acid 49% CMOS, Finyte, Finyte-1 Grade

Zinc (Zn)

Zirconium (Zr)

ee for detailed information section Microelectronic materials page 32 or the website www.jtbaker.com/micr

Use J.T.Baker Ultrex II and BAKER INSTRA-ANALYZED acids for low level trace element analysis.

See chapter 3 of this catalogue for more details.



Hydrofluoric Acid

48% / 'BAKER ANALYZED' / ACS

6013

HF	
M =	20.01 g/mol
11=	1.16 kg
CAS NO.	7664-39-3
EINECS	231-634-8
NC CODE	2811 11 00
EC NO.	9 003 00 1
UN/ID NO.	1790
ADR/RID	8 CT1
IMDG	8/11
R:	26/27/28-35
S:	26-36/37-45-7

Exceeds ACS Specifications	
Assay	min. 48.0-51.0%
Chloride (CI)	max. 5 ppm
Fluorsilicic Acid (H ₂ SiF ₆)	max. 0.01%
Heavy Metals (as Pb)	max. 0.5 ppm
Phosphate (PO ₄)	max. 1 ppm
Residue after Ignition	max. 5 ppm
Sulfate and Sulfite (as SO ₄)	max. 5 ppm
Trace Impurities (in ppm):	
Aluminium (AI)	max. 0.05
Arsenic and Antimony (as As)	max. 0.05
Barium (Ba)	max 0.1
Beryllium (Be)	max 0.01
Cadmium (Cd)	max 0.02
Calcium (Ca)	max 0.5
Chromium (Cr)	max 0.1
Cobalt (Co)	max 0.01
Copper (Cu)	max 0.01
Iron (Fe)	max 0.1
Lead (Pb)	max 0.05
Lithium (Li)	max 0.02
Magnesium (Mg)	max 0.1
Manganese (Mn)	max 0.01
Molybdenum (Mo)	max 0.02
Nickel (Ni)	max 0.05
Potassium (K)	max 0.1
Silver (Ag)	max 0.02
Sodium (Na)	max 0.5
Strontium (Sr)	max 0.02
Thallium (TI)	max 0.05
Titanium (Ti)	max 0.1
Vanadium (V)	max 0.01
Zinc (Zn)	max. 0.05
Zirconium (Zr)	max. 0.1

	BOX
11	(
	11

PACKING

PRODUCT

Hydrofluoric Acid 40% VLSI

See for detailed information section Microelectronic materials page 32 or the website www.jtbaker.com/micro



Calibrate and standardise your analytical methods and equipment with J.T.Baker Volumetric and Buffer solutions.

Refer to the Analytical applications section of this catalogue for more details.



Hydro

Hydrofluoric Acid 38% / 'BAKER ANALYZED'

► HF

M = 20.01 g/mol

1 I = 1.13 kg

CAS NO. 7664-39-3 **EINECS** 231-634-8

NC CODE 2811 11 00

EC NO. 9 003 00 1

UN/ID NO. 1790 **ADR/RID** 8 CT1

IMDG 8/II

R: 26/27/28-35

S: 26-36/37-45-7/9





A	00.400/
Assay	38-40%
Chloride (CI)	max. 5 ppm
Fluorsilicic Acid (H ₂ SiF ₆)	max. 0.005%
Phosphate (PO ₄)	max. 2 ppm
Residue after Ignition	max. 5 ppm
Sulfate and Sulfite (as SO ₄)	max. 5 ppm
Trace Impurities (in ppm):	
Aluminium (AI)	max. 0.05
Arsenic and Antimony (as As)	max. 0.05
Barium (Ba)	max. 0.1
Beryllium (Be)	max. 0.01
Cadmium (Cd)	max. 0.02
Calcium (Ca)	max. 0.5
Chromium (Cr)	max. 0.1
Cobalt (Co)	max. 0.01
Copper (Cu)	max. 0.01
Iron (Fe)	max. 0.1
Lead (Pb)	max. 0.05
Lithium (Li)	max. 0.02
Magnesium (Mg)	max. 0.1
Manganese (Mn)	max. 0.01
Molybdenum (Mo)	max. 0.02
Nickel (Ni)	max. 0.05
Potassium (K)	max. 0.1
Silver (Ag)	max. 0.02
Sodium (Na)	max. 0.5
Strontium (Sr)	max. 0.02
Thallium (TI)	max. 0.05
Titanium (Ti)	max. 0.1
Vanadium (V)	max. 0.01
Zinc (Zn)	max. 0.05
Zirconium (Zr)	max. 0.1

PRODUCT	PACKING	CONT.
NO.		вох
6079.1000	11	6
6079.2500	2.5	4



Hydrofluoric Acid 25% VLSI

See for detailed information section Microelectronic materials page 32 or the website www.jtbaker.com/micro

7048

Hydrofluoric Acid 5% VLSI

See for detailed information section Microelectronic materials page 32 or the website www.jtbaker.com/micro

Assay1)

1) At time of lot release.

Hydrogen Peroxide 35% / 'BAKER'

► H₂O₂

M = 34.01 g/mol

1 I = 1.13 kg CAS NO. 7722-84-1

EINECS 231-765-0

NC CODE 2847 00 00

EC NO. 8 003 00 9

UN/ID NO. 2014

ADR/RID 5.1 OC1

IMDG 5.1/II

R: 22-37/38-41

S: 26-36/37/39



min. 35%	PRODUCT	PACKING	CONT.
	NO.		BOX
	7048.1000	11	6
	7048.2500PE	2.5 HDPE	4
	7048.9025	25	
	For safe hand	ling of 25 Ltin cans see Self-closing tan	



5155

Hydrogen Peroxide

30% / ULTREX II Ultrapure Reagent

M = 34.01 g/mol 1 I = 1.11 kg CAS NO. 7722-84-1 EINECS 231-765-0 NC CODE 2847 00 00 EC NO. 8 003 00 9 UN/ID NO. 2014 ADR/RID 5.1 OC1 IMDG 5.1/II R: 22-41 S: 26-36/39

Meets ACS Ultratrace Specifications	
Assay	25.0-35.0%
Appearance	passes test
Free acid (µeq/g)	max. 0.5
Trace Impurities (in ppm):	
Ammonium (NH ₄)	max. 5
Chloride (CI)	max. 1
Nitrate (NO ₃)	max. 2
Phosphate (PO ₄)	max. 2
Sulfate (SO ₄)	max. 5
Trace Impurities (in ppt):	
Aluminium (AI)	max. 100
Antimony (Sb)	max. 1000
Arsenic (As)	max. 1000
Barium (Ba)	max. 50
Beryllium (Be)	max. 1000
Bismuth (Bi)	max. 1000
Boron (B)	max. 1000
Cadmium (Cd)	max. 50
Calcium (Ca)	max. 100
Chromium (Cr)	max. 50
Cobalt (Co)	max. 50
Copper (Cu)	max. 50
Gallium (Ga)	max. 1000
Iron (Fe)	max. 100
Lead (Pb)	max. 100
Lithium (Li)	max. 50
Magnesium (Mg)	max. 50
Manganese (Mn)	max. 50
Mercury (Hg)	max. 1000
Molybdenum (Mo)	max. 1000
Nickel (Ni)	max. 50
Potassium (K)	max. 100
Silicon (Si)	max. 10000
Silver (Ag)	max. 50
Sodium (Na)	max. 100
Strontium (Sr)	max. 50
Tin (Sn)	max. 50

max. 1000 max. 1000

Zinc	: (Zn)		max. 50
Zirco	onium (Zr		max. 1000
PRO	DUCT	PACKING	CONT.
5155	5.0450	450 ml PE	

The J.T.Baker CYCLE-TAINER
High Purity Solvent Delivery System,
preserves purity and protects people.

Titanium (Ti)

Vanadium (V)

See chapter 3 of this catalogue for product details.



Hydro

Hydrogen Peroxide 30% / 'BAKER ANALYZED' / ACS

► H₂O₂ **M** = 34.01 g/mol **1 l** = 1.11 kg **CAS NO.** 7722-84-1 **EINECS** 231-765-0 NC CODE 2847 00 00 **EC NO.** 8 003 00 9 **UN/ID NO.** 2014 **ADR/RID** 5.1 O C1 IMDG 5.1/II R: 22-41 **S:** 26-36/39

7047

Exceeds ACS Specifications	
Assay1)	30.0-32.0%
Ammonium (NH ₄)	max. 5 ppm
Chloride (CI)	max. 1 ppm
Color (APHA)	max. 10
Heavy Metals (as Pb)	max. 1 ppm
Nitrate (NO ₃)	max. 2 ppm
Phosphate (PO ₄)	max. 2 ppm
Residue after Evaporation	max. 0.002%
Sulfate (SO ₄)	max. 2 ppm
Titrable Acid (meq/g)	max. 0.0006
Trace Impurities (in ppm):	
Aluminium (AI)	max. 0.05
Arsenic (As)	max. 0.05
Barium (Ba)	max. 0.02
Beryllium (Be)	max. 0.01
Cadmium (Cd)	max. 0.02
Calcium (Ca)	max. 0.5
Chromium (Cr)	max. 0.1
Cobalt (Co)	max. 0.01
Copper (Cu)	max. 0.01
Germanium (Ge)	max. 0.05
Iron (Fe)	max. 0.1
Lead (Pb)	max. 0.05
Lithium (Li)	max. 0.01
Magnesium (Mg)	max. 0.1
Manganese (Mn)	max. 0.01
Molybdenum (Mo)	max. 0.02
Nickel (Ni)	max. 0.05
Potassium (K)	max. 0.1
Silver (Ag)	max. 0.01
Sodium (Na)	max. 0.5
Strontium (Sr)	max. 0.01

max. 0.05

max. 0.01 max. 0.05

max. 0.1

PRODUCT	PACKING	CONT.
NO.		вох
7047.0250	250 ml	
7047.1000	11	6
7047.2500PE	2.5 HDPE	



Hydrogen Peroxide 30% CMOS, Finyte, Finyte-1, Ultryte Grade

1) At time of lot release.

Thallium (TI)

Titanium (Ti) Vanadium (V)

Zinc (Zn) Zirconium (Zr)

iled information section Microelectronic materials page 32 or the website www.jtbaker.com/micro



Hydrogen Peroxide 30% MOS, VLSI, SLSI Grade



Hydrogen Peroxide Solution

lities meeting the requirements of cGMP, USP/NF, BP/Ph.Eur, JP or FCC, page 36

Mallinckrodt Baker's Microelectronic Materials are Engineered for the future™.

See chapter 5 of this catalogue or visit us at www.jtbaker.com/micro

CONT.

вох

1116

2195

CONT.

BOX

Hydroquinone

1115

max. 0.001%

172-174°C

max. 0.05%

PRODUCT

1115.0250

PRODUCT

1116.9050

NO.

min. 99.0%

max. 0.1%

passes test

max. 0.01%

max. 0.005%

max 5

max. 5

max. 0.05

max. 0.25 meq/g

PACKING

50 kg

NO.

PACKING

250 g

Heavy Metals (as Pb)

Residue after Ignition

Melting Point



M = 110.11 g/mol

FLASHPOINT $165\,^{\circ}\text{C}$ **CAS NO.** 123-31-9

> **EINECS** 204-617-8 NC CODE 2907 22 00

EC NO. 604 005 00 4

UN/ID NO. 2662

ADR/RID 6.1 T2

IMDG 6.1/III

R: 22-40-41-43-50

S: 26-36/37/39-61







Hydroxyacetic Acid



(2-Hydroxyethl)mercaptan

Hydroxylamine Hydrochloride

e Hydroxylammonium Chloride

Hydroxylammonium Chloride

(max. 0.05 ppm Hg) / 'BAKER ANALYZED' / ACS

► HONH₃CI Exceeds ACS Specifications

M = 69.49 g/mol **CAS NO.** 5470-11-1 **EINECS** 226-798-2

NC CODE 2825 10 00 **EC NO.** 612 123 00 2

UN/ID NO. 3260 ADR/RID 8 C2 IMDG 8/III

R: 22-36/38-43-48/22-50

S: 22-24-37-61





Hydroxylammonium Chloride

M = 69.49 g/mol **CAS NO.** 5470-11-1

Assay

Iron (Fe)

Mercury (Hg)

Ammonium (NH₄)

Titrable Free Acid

Clarity of alcohol solution

Sulfur Compounds (as SO₄)

Trace Impurities (in ppm):

Residue after Ignition

Heavy Metals (as Pb)

EINECS 226-798-2 NC CODE 2825 10 00 **EC NO.** 612 123 00 2

UN/ID NO. 3260

ADR/RID 8 C2 IMDG 8/III

► HONH₃CI

R: 22-36/38-43-48/22-50

S: 22-24-37-61






Exceeds ACS Specifications	
Assay	min. 99.0%
Ammonium (NH ₄)	max. 0.1%
Clarity of alcohol solution	passes test
Residue after Ignition	max. 0.01%
Sulfur Compounds (as SO ₄)	max. 0.005%
Titrable Free Acid	max. 0.25 meq/g
Trace Impurities (in ppm):	
Heavy Metals (as Pb)	max. 5
Iron (Fe)	max. 5

PRODUCT	PACKING	CONT.
NO.		вох
2195.0100	100 g	
2195.0250	250 g	6
2195.1000	1 kg	
2195.9050	50 kg	



Hydro

Hydroxylammonium Chloride





4-Hydroxy-4-Methyl-2-Pentanone MOS Grade

6016

Hypophosphorous Acid

► H ₃ PO ₂		Assay	49-52%	PRODUCT	PACKING	CONT.
M =	66.00 g/mol	Heavy Metals (as Pb)	max. 0.002%	NO.		вох
11=	1.25 kg	Identification	passes test	6016.0500	500 ml	
CAS NO.	6303-21-5	Phosphate (PO ₄)	max. 0.01%			
EINECS	228-601-5					
NC CODE	2811 19 80					
UN/ID NO.	3264					
ADR/RID	8 C1					
IMDG	8/III					
R:	34					
S:	26-27-36/37/39-45					
	C C					



ICP CCV Solution I & II

See CCV Solution I & II



ICP ICS Solution AB-1&AB-2

See ICS Solution AB-1 & AB-2



ICP ICV Solution I & II

See ICV Solution I & II



ICP Interference Check Standard I

ee Interference Check Standard I



ICP Interference Check Standard II

See Interference Check Standard II



ICP Interference Check Standard IV

See Interference Check Standard IV



ICP Mixed Calibration Standard III



ICP Mixed Calibration Standard IV

See Mixed Calibration Standard IV



ICP Mixed Calibration Standard V

See Mixed Calibration Standard V

ICP Primary Drinking Water Std I See Primary Drinking Water Standard I

ICP Trace Metal Standard I

See Trace Metal Standard I

ICP Trace Metal Standard III



ICP Water Standard I & II

See Water Standard I & II



ICS Solution AB-1 & AB-2

Matrix I : 5% nitric acid - Matrix II: 20% HCI / 'BAKER INSTRA-ANALYZED' / Plasma Standard

6119-01

NC CODE 3822 00 00 UN/ID NO. 3264 ADR/RID 8 C1 IMDG 8/III **R:** 34 **S:** 23-26-36/37/39-45

KIT CONTAINS ONE DOTTIE OF EACH SOLUTION	
Element Concentrations of Solution I (µg/ml):	
Beryllium (Be)	50
Cadmium (Cd)	100
Cobalt (Co)	50
Copper (Cu)	50
Lead (Pb)	100
Manganese (Mn)	50
Nickel (Ni)	100
Silver (Ag)	100
Zinc (Zn)	100
Element Concentrations of Solution II ($\mu g/mI$):	
Barium (Ba)	50
Chromium (Cr)	5
Vanadium (V)	50

PRODUCT	PACKING	CONT.
NO.		вох
5119-01	100 ml x 2	

For use in EPA Contract Laboratory Program (CLP). Traceable to NIST.

ICV Solution I & II

Matrix I : 5% nitric acid - Matrix II: 20% HCI / 'BAKER INSTRA-ANALYZED' / Plasma Standard

6103-01

NC CODE 3822 00 00 **UN/ID NO.** 3264 ADR/RID 8 C1 IMDG 8/III **R:** 34 **S:** 23-26-36/37/39-45

Kit contains one bottle of each solution	
Element Concentrations of Solution I (µg	g/ml):
Barium (Ba)	100
Beryllium (Be)	40
Cadmium (Cd)	50
Cobalt (Co)	100
Copper (Cu)	100
Iron (Fe)	100
Lead (Pb)	100
Manganese (Mn)	100
Nickel (Ni)	100
Silver (Ag)	20
Thallium (TI)	100
Zinc (Zn)	100
Element Concentrations of Solution II (μ	g/ml):
Aluminium (AI)	100
Antimony (Sb)	100
Arsenic (As)	100
Calcium (Ca)	1000
Chromium (Cr)	100
Magnesium (Mg)	1000
Potassium (K)	1000
Selenium (Se)	100
Sodium (Na)	1000
Vanadium (V)	100

PRODUCT	PACKING	CONT.
NO.		вох
6103-01	100 ml x 2	

For use in EPA Contract Laboratory Program (CLP). Traceable to NIST.

Imida

Imidazole

1747

'BAKER ANALYZED

► NHCH:NCH:ÇH **M** = 68.08 g/mol

CAS NO. 288-32-4

EINECS 206-019-2 NC CODE 2933 29 90

UN/ID NO. 3263

ADR/RID 8 C8

IMDG 8/III

R: 22-34-63

S: 22-26-36/37/39-45



min. 99.5% Assay Heavy Metals (as Pb) max. 0.001% Iron (Fe) max. 0.001% Melting Point 88-90°C Residue after Ignition max. 0.1%

PRODUCT	PACKING	CONT.
NO.		вох
1747.0100	100 g	
1747.9010	10 kg	

Immuno PBS

3058

NC CODE 3822 00 00

PRODUCT PACKING BOX 3058.9010 10 l Polycube

Dilute 20 times before use.

PRODUCT

Interference Check Standard I

(Matrix: 5% nitric acid) / 'BAKER INSTRA-ANALYZED' / Plasma Standard

NC CODE 3822 00 00 **UN/ID NO.** 3264 ADR/RID 8 C1 IMDG 8/III

R: 36/38-45

S: 26-37-45-53



Elemental Concent (µg/ml):	
Arsenic (As)	1000
Barium (Ba)	300
Beryllium (Be)	100
Cadmium (Cd)	300
Chromium (Cr)	300
Cobalt (Co)	300
Copper (Cu)	300
Lead (Pb)	1000
Manganese (Mn)	200
Mercury (Hg)	50
Nickel (Ni)	300
Potassium (K)	20000
Selenium (Se)	500
Silver (Ag)	300
Thallium (TI)	1000
Vanadium (V)	300
Zinc (Zn)	300

NO.		вох
6011-01	100 ml	

For use in EPA SW-846 Methods 6010 and 200.7. Traceable to NIST.

PACKING

Interference Check Standard II

6012-01 (Matrix: 5~% nitric acid plus a trace of hydrofluoric acid) / 'BAKER INSTRA-ANALYZED' / Plasma Standard

UN/ID NO. 3264 ADR/RID 8 C1 IMDG 8/III **R:** 34

S: 23-26-36/37/39-45



Elemental Concent (µg/ml):	
Boron (B)	500
Molybdenum (Mo)	300
Silicon (Si)	230
Titanium (Ti)	1000

PRODUCT	PACKING	CONT.
NO.		вох
6012-01	100 ml	

For use in EPA SW-846 Methods 6010 and 200.7. Traceable to NIST.

www.jtbaker.com/europe



J.T.Baker

Interference Check Standard IV

(Matrix: 5% nitric acid) / 'BAKER INSTRA-ANALYZED' / Plasma Standard

6014-01

NC CODE	3822 00 00
UN/ID NO.	3264
ADR/RID	8 C1
IMDG	8/III
R:	34
S:	23-26-36/37/39-45

Elemental Concent (µg/ml):	
Aluminium (AI)	1200
Calcium (Ca)	6000
Iron (Fe)	5000
Magnesium (Mg)	3000
Sodium (Na)	1000

PRODUCT NO.	PACKING	CONT BO)
6014-01	100 ml	

For use in EPA SW-846 Methods 6010 and 200.7. Traceable to NIST.

Iodide/Iodate

4876

EINECS 231-659-4 NC CODE 3822 00 00

PRODUCT	PACKING	CONT.
NO.		вох
4876	1 amp.	6

Volumetric Concentrate, for dilution to 1 l.

Iodine

sublimed / 'BAKER ANALYZED' / ACS

1118

I_2	

M = 253.81 g/mol

CAS NO. 7553-56-2

EINECS 231-442-4

NC CODE 2801 20 00 **EC NO.** 53 001 00 3

UN/ID NO. 1759

ADR/RID 8 C10

IMDG 8/III

R: 20/21-50 **S:** 23-25-61

¥2

dangerous for the environment

Exceeds ACS Specifications	
Assay	min. 99.8%
Chloride and Bromine (as CI)	max. 0.003%
Jon volatila Matter	may 0.01%

NO.	PACKING	CONI	
		BO)	
1118.0100	100 g	(
1118.0250	250 g	(
1118.1000	1 kg		
1118.5000	5 kg		

Iodine

0132



M = 253.81 g/mol

CAS NO. 7553-56-2

EINECS 231-442-4

NC CODE 2801 20 00

EC NO. 53 001 00 3

UN/ID NO. 1759

ADR/RID 8 C10

IMDG 8/III **R:** 20/21-50

S: 23-25-61



dangerous for the environment



Assay 99.5-100.5% Bromides and Chlorides max. 250 ppm Identification passes test Non-volatile substances max. 0.1%

PRODUCT PACKING CONT. NO. вох 0132.9050 50 kg

Iodine

See for qualities meeting the requirements of cGMP, USP/NF, BP/Ph.Eur, JP or FCC, page 36





Iodin

Iodine Solution

0.5 mol/l 1N / 'BAKER ANALYZED' 7577

M = 253.81 g/mol **1 l** = 1.26 kg **CAS NO.** 7553-56-2 **EINECS** 231-441-9 NC CODE 2801 20 00

R: 20/21/22

S: 26-36/37

 $0.5M~\pm~0.5\%$

Molarity (titration with Na₂S₂O₃)

PRODUCT **PACKING** CONT. NO. 7577.1000 11

вох

CONT.

BOX

6

Volumetric Solution, ready for use.

Iodine Solution

7463 0.1 mol IBr/l; 0.2N / 'BAKER ANALYZED' / According to Hanus

Titer (N)

Titer (N)

FLASHPOINT 40 °C NC CODE 2801 20 00 **UN/ID NO.** 2789 ADR/RID 8 CF1 IMDG 8/II **R:** 10-35

S: 26-36/37/39-45

0.19-0.21 PRODUCT **PACKING** CONT. вох NO. 7463.1000 11

Volumetric Solution, ready for use.

Each lot of this product is standardized potentiometrically against NIST tracable reference standard.

Iodine Solution

6124 0.1 mol ICI/I, 0.2N solution / 'BAKER' / According to Wijs

1 l = 1.21 kg FLASHPOINT 40 °C UN/ID NO. 2920 ADR/RID 8 CF1 IMDG 8/II

R: 23/24/25-34-40-48/23-52/53-59

S: 26-36/37/39-45





Keep in a cool place and tightly closed.

11

PACKING

0.2

PRODUCT

6124.1000

Iodine Solution

0.05 mol/I 0.1N / 'BAKER ANALYZED'

EINECS 231-441-9 NC CODE 2801 20 00

Titer (N) 0.0997-0.1003 **PRODUCT PACKING** CONT. NO. BOX 7210.1000 11 7210.2500 2.5 |

Volumetric Solution, ready for use. Each lot of this product is standardized potentiometrically against Sodium Thiosulfate.

Iodine Solution

7128 0.05 mol/l 0.1N / 'BAKER ANALYZED'

EINECS 231-441-9 NC CODE 2801 20 00

0.0997-0.1003 Titer (N)

PRODUCT PACKING CONT. NO. вох 7128.1000 11 7128.2500 2.5 l

Volumetric Solution, ready for use.

Each lot of this product is standardized potentiometrically against Arsenic Trioxide (primary standard).

$\begin{array}{c} \textbf{Iodine (Iodine-Iodide)} \\ \textbf{0.05 mol I}_2/I \ / \ \textbf{DILUT-IT} \ / \ \textbf{1}/_{10} \ \textbf{equiv.} \end{array}$

v. = 12.6905g; 0.1N

4662

Þ	
	-

M = 253.81 g/mol **CAS NO.** 7553-56-2 **EINECS** 231-441-9

NC CODE 3822 00 00

PRODUCT **PACKING** CONT. NO. вох 4662 1 amp.

Volumetric Concentrate, for dilution to 1 l.

Volumetric Concentrate, for dilution to 1 l.

PACKING

PRODUCT

Iodine (Iodine-Iodide)

 $0.005 \text{ mol } I_2/I / DILUT-IT / \frac{1}{100} \text{ equiv.} = 1.269g; 0.01N$

4660



M = 253.81 g/mol **CAS NO.** 7553-56-2 **EINECS** 231-441-9 NC CODE 3822 00 00

PRODUCT PACKING CONT. BOX NO. 4660 1 amp.

Ion Exchange Resin, Anion, IONAC NA-38

4601

1099

1100

CONT.

NC CODE	3914 00 00	

(strong base, alkyl quaternary ammonius	m
polystyrene copolymer)	
Moisture	55-65%
Total Exchange Capacity (wet volume),	
(meq/ml)	min. 0.9
Mesh (Wet Screen Analysis):	
On U.S. No. 16 Sieve	max. 5%
Thru U.S. No. 50 Sieve	max. 5%

NO.		ROX
4601.0500	500 g	
OH- Form,	Type I, Beads (16-50 Mesh).	

Ion exchange Resin, Cation, Dowex C-211 'BAKER ANALYZED'

NC CODE 3914 00 00

NC CODE 3914 00 00

Identification	passes test
Total Exchange Capacity (wet volume),	
(meq/ml)	min. 1.8
Mesh (Wet Screen Analysis):	
On U.S. No. 16 Sieve	max. 5%
Thru U.S. No. 50 Sieve	max. 1%

PRODUCT	PACKING	CONT.
NO.		вох
1099.0500	500 g	
H+ Form Si	oherical Beads (16-50 Mesh).	

lon exchange Resin, Cation, Dowex HCR-W2 'BAKER ANALYZED'

(strong acid; styrene-DVB copolymer; nuclear sulfonic acid active group)		
Identification	passes test	
Total Exchange Capacity (dry basis),		
meq/g	act. value reported	
Total Exchange Capacity (wet volume)),	
(meq/ml)	min. 1.9	
Mesh (Wet Screen Analysis):		
On U.S. No. 16 Sieve	max. 5%	

max. 5%

PRODUCT	PACKING	CONT
NO.		B0)
1100.0500	500 g	

Na+ Form, Spherical Beads (16-40 Mesh).

Ion exchange Resin, Mixed bed, IONAC NM-60 BAKER ANALYZED

Thru U.S. No. 40 Sieve

4631

(strong acid/strong base; sulfonated/alkyl	
quaternary ammonium polystyrenes)	
Column Capacity, kgr/cf	min. 12.0
Column Capacity, meq/ml	min. 0.55
Effluent Greater than 15 Megohms	min. 50%
Moisture	50-60%
Mesh (Wet Screen Analysis):	
On U.S. No. 16 Sieve	max. 5%
Thru U.S. No. 50 Sieve	max. 5%

PRODUCI	PACKING	CONT
NO.		ВОХ
4631.0500	500 g	

H+/OH- Form, Type I, Beads (16-50 Mesh).

Iron

Iron

0134

M = 55.85 g/mol **CAS NO.** 7439-89-6

EINECS 231-096-4 NC CODE 7205 29 00

min. 96.0% Insoluble in H₂SO₄ max. 0.2% Lead (Pb) max. 0.005% Sulfide (S) max. 0.03% Water Soluble Substances max. 0.03% Trace Impurities (in ppm):

CONT. **PRODUCT PACKING** вох 0134.0500GL 500 g Glass

/ 'BAKER' / reduced by hydrogen

/ 'BAKER ANALYZED'

M = 55.85 g/mol CAS NO. 7439-89-6 **EINECS** 231-096-4 NC CODE 7205 29 00

min. 96% Assav Acid Insoluble Substances max. 1.25% Lead (Pb) max. 0.003% Trace Impurities (in ppm): Arsenic (As) max. 8 Mercury (Hg) max. 5

PRODUCT PACKING CONT. 0135.2500 2.5 kg

Iron 1000 µg/ml

(Matrix: 2% nitric acid) / 'BAKER INSTRA-ANALYZED' / Plasma Standard

Arsenic (As)

► Fe

M = 55.85 g/mol NC CODE 3822 00 00

R: 36/38



Certificate Provided Reporting Actual Lot Analysis 998-1002 μg/ml Iron (Fe)

PRODUCT PACKING CONT. 5764.0100 100 ml

Prepared from the highest purity raw material available, generally greater than 99.999% spectral purity. The content of the solution is confirmed to be accurate to within \pm 0.2%. Typically 1000 $\mu\text{g/ml.}$ The certificate of analysis provided reports actual lot analysis. The certificate also lists the trace impurities.

Iron 1000 μ g/ml

(Matrix: 2% nitric acid) / 'BAKER INSTRA-ANALYZED' / Atomic Absorption Standard

Iron (Fe)

M = 55.85 g/mol

NC CODE 3822 00 00

6929

R: 36/38

S: 26-37



998-1002 µg/ml **PRODUCT** NO.

6929.0100 100 ml 6929.0500 500 ml

PACKING

Prepared by dissolution of high purity raw materials (min. 99.99% spectral purity). Assays are verified by ICP against standards traceable to NIST. Standard Reference Material numbers (SRM) are printed on each label.

BOX

CONT.

BOX

Iron 1000 μ g/ml

'BAKER ANALYZED' / Atomic Absorption Standard

M = 55.85 g/mol

NC CODE 3822 00 00 R: 36/38

S: 26-37



Iron (Fe) 998-1002 µg/ml

NO. 6812.0100 100 ml 6812.0500 500 ml

PRODUCT

Iron(III)nitrate in nitric acid 0.5 mol/l.

PACKING

Innovation is principal to our business.

5731

▶ Fe

M = 55.85 g/mol

NC CODE 3822 00 00

R: 36/38

S: 26



Certificate Provided Reporting Actual Lot Analysis

9980-10020 μg/ml Iron (Fe)

PRODUCT	PACKING	CONT.
NO.		вох
5731 0100	100 ml	

Prepared from the highest purity raw material available, generally greater than 99.999% spectral purity. The content of the solution is confirmed to be accurate to within 0.2 %. Typically 10000 $\mu g/ml$.

The certificate of analysis provided reports actual lot analysis. The certificate also lists the trace impurities.

Iron Ammonium Sulfate

See Ammonium Iron (II) Sulfate Hexahydrate



Iron(II) Chloride Tetrahydrate

0125



M = 198.81 g/mol

CAS NO. 13478-10-9

EINECS 213-843-4 NC CODE 2827 33 00

UN/ID NO. 3260

ADR/RID 8 C2

IMDG 8/III

R: 22-38-41

S: 26-39

Assay	min. 99.0%
Appearance of solution	passes test
Arsenic (As)	max. 0.001%
Copper (Cu)	max. 0.005%
Insoluble Matter	max. 0.01%
Phosphate (PO ₄)	max. 0.001%

Sulfate (SO₄)

Zinc (Zn)

Assay

max. 0.01% max. 0.005%

CONT. **PRODUCT PACKING** NO. вох 0125.0250 250 g 0125.9050 50 kg

Iron(III) Chloride Anhydrous

M = 162.22 g/mol CAS NO. 7705-08-0

EINECS 231-729-4

NC CODE 2827 33 00

UN/ID NO. 1773

► FeCl₃

ADR/RID 8 C2 IMDG 8/II

R: 22-34-41

S: 26-36/39-45

min. 98%	PRODUCT	PACKING
	NO.	
	9051.1000	1 kg

Iron(III) Chloride Hexahydrate

Lump / 'BAKER ANALYZED' / ACS

0119

9051

CONT.

BOX

► FeCl₃.6H₂O

M = 270.30 g/mol

CAS NO. 10025-77-1

EINECS 231-729-4 NC CODE 2827 33 00

R: 22-38-41

S: 26-39

Meets ACS Specifications	
Assay (by lodometry)	97.0-102.0%
Calcium (Ca)	max. 0.01%
Copper (Cu) (by AAS)	max. 0.003%
Ferrous Iron (as Fe)	max. 0.002%
Insoluble Matter	max. 0.01%
Magnesium (Mg)	max. 0.005%
Nitrate (NO ₃)	max. 0.01%
Phosphorus Compounds (as PO ₄)	max. 0.01%
Potassium (K)	max. 0.005%
Sodium (Na)	max. 0.05%
Sulfate (SO ₄)	max. 0.01%
Zinc (Zn) (by AAS)	max. 0.003%

PRODUCT	PACKING	CONT
NO.		ВОХ
0119.0100	100 g	
0119.1000	1 kg	ϵ



Ironi

Iron(III) Nitrate Nonahydrate 0121

BAKER ANALYZED' / ACS

► Fe(NO₃)₃.9H₂O **M** = 404.00 g/mol **CAS NO.** 7782-61-8 **EINECS** 233-899-5 NC CODE 2834 29 80 **UN/ID NO.** 1466 **ADR/RID** 5.1 O2 IMDG 5.1/III **R:** 36/38-8

S: 26

Exceeds ACS Specifications Assay (by lodometry) 98.0-101.0 Insoluble Matter max. 0.005^o Manganese (Mn) max. 0.010 pH of 5% Solution at 25°C 1.5-2. Substances not Precipitated by NH₄OH (as max. 0.19 Sulfate (SO₄) max. 0.019 Trace Impurities (in ppm): Chloride (CI) max.

PRODUCT NO.	PACKING	CONT. BOX
0121.1000	1 kg	

Iron(III) Oxide

0122

► Fe₂O₃ **M =** 159.69 g/mol **CAS NO.** 1309-37-1 **EINECS** 215-168-2 NC CODE 2821 10 00

Assay (by lodometry) min. 98.0% max. 0.01% Copper (Cu) (by AAS) Insoluble in HCI max. 0.2% Manganese (Mn) (by AAS) max. 0.05% Phosphate (PO₄) max. 0.02% Sulfate (SO₄) max. 0.2% Zinc (Zn) (by AAS) max. 0.01%

PRODUCT	PACKING	CONT.
NO.		вох
0122.0500	500 g	
0122.9012	12 kg PE Pail	



Iron Reduced

See for qualities meeting the requirements of cGMP, USP/NF, BP/Ph.Eur, JP or FCC, page 36

0123

Iron(III) Sulfate n-Hydrate

► Fe₂(SO₄)₃.nH₂O **CAS NO.** 15244-10-7 **EINECS** 233-072-9 NC CODE 2833 29 50 **R:** 22-41 **S**: 26

Assay (Fe₂(SO₄)₃) min. 73.0% Chloride (CI) max. 0.002% Copper (Cu) max. 0.005% Ferrous Iron (as Fe) max. 0.02% Insoluble Matter max. 0.02% Nitrate (NO₃) max. 0.01% Substances not Precipitated by NH₄OH (as SO₄) max. 0.1% max. 0.005% Zinc (Zn)

PRODUCT PACKING CONT. NO. BOX 1 kg 0123,1000

0126

Iron(II) Sulfate Heptahydrate

► FeSO₄.7H₂O **M =** 278.02 g/mol **CAS NO.** 7782-63-0 **EINECS** 231-753-5 NC CODE 2833 29 50 R: 22-41 **S**: 26



Meets ACS Specifications	
Assay (by KMnO ₄ titrn.)	min. 99.0%
Calcium (Ca)	max. 0.005%
Chloride (CI)	max. 0.001%
Copper (Cu)	max. 0.005%
Ferric ion (Fe+)	max. 0.1%
Insoluble Matter	max. 0.01%
Magnesium (Mg)	max. 0.002%
Manganese (Mn)	max. 0.05%
Phosphate (PO ₄)	max. 0.001%
Potassium (K)	max. 0.002%
Sodium (Na)	max. 0.02%
Zinc (Zn)	max. 0.005%

PRODUCI	PACKING	CON1.
NO.		вох
0126.0500	500 g	6

Iron(II) Sulfate Heptahydrate

0127

► FeSO₄.7H₂O

M = 278.02 g/mol **CAS NO.** 7782-63-0

EINECS 231-753-5

NC CODE 2833 29 50

R: 22-41 **S:** 26

Assay	98.0-105.0%
Appearance of solution	passes test
Chlorides (as CI)	max. 300 ppm
Ferric Ions	max. 0.5%
Heavy Metals (as Pb)	max. 50 ppm
Identification	passes test
Manganese (Mn)	max. 0.1%
рН	3.0-4.0
Zinc (Zn)	max 500 nnm

PRODUCT NO.	PACKING	CONT. BOX
0127.1000	1 kg	
0127.9050	50 kg	

Iron Sesquinoxide

See Iron(III) Oxide



8010

Isoamyl Alcohol

► (CH₃)₂CHCH₂CH₂OH **M** = 88.15 g/mol

11= 0.81 kg

FLASHPOINT 43 °C

CAS NO. 123-51-3

EINECS 204-633-5

NC CODE 2905 15 00

UN/ID NO. 1105

ADR/RID 3 F1

IMDG 3/III

R: 10-20

S: 24/25



Assay (by GC)	min. 98%
Acids and Esters (as Amyl Acetate)	max. 0.06%
Aldehydes	passes test
Boiling Range	128-132°C
Color (APHA)	max. 5
Residue after Evaporation	max. 0.003%
Substances Darkened by H ₂ SO ₄	passes test
Titrable Acid (meq/g)	max. 0.002
Water (H ₂ O)	max. 0.2%
Trace Impurities (in ppm):	
Aluminium (AI)	max. 0.5
Barium (Ba)	max. 0.1
Boron (B)	max. 0.02
Cadmium (Cd)	max. 0.05
Calcium (Ca)	max. 0.5
Chromium (Cr)	max. 0.02
Cobalt (Co)	max. 0.02
Copper (Cu)	max. 0.02
Iron (Fe)	max. 0.1
Lead (Pb)	max. 0.1
Magnesium (Mg)	max. 0.1
Manganese (Mn)	max. 0.02
Nickel (Ni)	max. 0.02
Tin (Sn)	max. 0.1
Zinc (Zn)	max. 0.1

PRODUCT	PACKING	CON
NO.		BC
8010.1000	11	
8010.9010	10 l Jerrycan	
8010.9025	25 l	
8010.9200	200 l	
For safe hand	dling of 25 l tin cans, see Self-	closing tap.

Isobutyl Alcohol

'BAKER HPLC ANALYZED' / For use in High Performance Liquid Chromatography

9048-03

► (CH₃)₂CHCH₂OH

M = 74.12 g/mol

1 l = 0.80 kg

 $\textbf{FLASHPOINT} \ \ 27\ ^{\circ}\text{C}$ **CAS NO.** 78-83-1

EINECS 201-148-0

NC CODE 2905 14 90

EC NO. 603 004 00 6 **UN/ID NO.** 1212

ADR/RID 3 F1

IMDG 3/III

R: 10-22-37/38-41-67

S: 13-26-37/39-46-7/9



Water (by KF, coulometric)	max. 0.05%
Ultraviolet Absorbance (1.00-cm	path vs water):
at 254 nm	max. 0.05
at 280 nm	max. 0.02
at 350 nm	max. 0.01

min. 99.5%

max. 5 ppm

max. 220

PF

N

Assay (by GC) (corrected for water)

Residue after Evaporation

UV Cut-off, nm

RODUCT	PACKING	CONT.
0.		вох
048-03	4 Glass	

Isobu

Isobutyl Alcohol 'BAKER ANALYZED'

8056 ► (CH₃)₂CHCH₂OH

M = 74.12 g/mol

1 l = 0.80 kg

 $\textbf{FLASHPOINT} \ \ 27\ ^{\circ}\text{C}$ **CAS NO.** 78-83-1

EINECS 201-148-0

NC CODE 2905 14 90

EC NO. 603 004 00 6

UN/ID NO. 1212

ADR/RID 3F1 IMDG 3/III

R: 10-22-37/38-41-67

S: 13-26-37/39-46-7/9

Assay (by GC)	min. 99.0 %
Carbonyl Compounds	max. 0.1%
Color (APHA)	max. 10
Residue after Evaporation	max. 0.001%
Solubility in Water	passes test
Titrable Acid (meq/g)	max. 0.0005
Water (H ₂ O)	max. 0.1%
Trace Impurities (in ppm):	
Aluminium (AI)	max. 0.5
Barium (Ba)	max. 0.1
Boron (B)	max. 0.02
Cadmium (Cd)	max. 0.05
Calcium (Ca)	max. 0.5
Chromium (Cr)	max. 0.02
Cobalt (Co)	max. 0.02
Copper (Cu)	max. 0.02
Iron (Fe)	max. 0.1
Lead (Pb)	max. 0.1
Magnesium (Mg)	max. 0.1
Manganese (Mn)	max. 0.02
Nickel (Ni)	max. 0.02
Tin (Sn)	max. 0.1
Zinc (Zn)	max. 0.1

PRODUCT	PACKING	CONT.
NO.		вох
8056.1000	11	6

8210

Isobutyl Alcohol

► (CH₃)₂CHCH₂OH

M = 74.12 g/mol

11= 0.80 kg

 $\textbf{FLASHPOINT} \ \ 27\ ^{\circ}\text{C}$

CAS NO. 78-83-1 **EINECS** 201-148-0

NC CODE 2905 14 90

EC NO. 603 004 00 6

UN/ID NO. 1212

ADR/RID 3 F1

IMDG 3/III

R: 10-22-37/38-41-67

S: 13-26-37/39-46-7/9



Boiling Range 106-109°C Residue after Evaporation max. 0.005% Water (H₂O) max. 0.5%

PRODUCT PACKING CONT. BOX NO. 2.5 l 8210.2500 8210.9025 25 l For safe handling of 25 l tin cans, see Self-closing tap.



See for qualities meeting the requirements of cGMP, USP/NF, BP/Ph.Eur, JP or FCC, page 36



Iso-Octane

See 2,2,4-Trimethylpentane



Isopropyl Alcohol



Isopropyl Alcohol

ting the requirements of cGMP, USP/NF, BP/Ph.Eur, JP or FCC, page 36



Isopropyl Alcohol 70% Solution, Sterile

See for qualities meeting the requirements of cGMP, USP/NF, BP/Ph.Eur, JP or FCC, page 36



Isopropylamine

8821 ► (CH₃)₂CHNH₂ 31-32°C PRODUCT PACKING CONT. **Boiling Point M** = 59.11 g/mol NO. вох **1 l** = 0.69 kg 8821.0100 100 ml **FLASHPOINT** -37 °C **CAS NO.** 75-31-0 **EINECS** 200-860-9 NC CODE 2921 19 30 EC NO. 612 007 00 1 **UN/ID NO.** 1221 ADR/RID 3 FC IMDG 3/I **R:** 12-36/37/38 **S:** 16-26-29

Kaolin

See for qualities meeting the requirements of cGMP, USP/NF, BP/Ph.Eur, JP or FCC, page 36



Karl Fischer Reagent

5.0 mg H₂O/ml / 'BAKER ANALYZED' / Titrant for moisture determination

		,				
1 =	1.18 kg	Titer (mg H ₂ O/ml)	min. 5.0	PRODUCT	PACKING	CONT.
FLASHPOINT	>21 °C			NO.		вох
NC CODE	3822 00 00			8792.1000	11	6
UN/ID NO.	1993			8792.2500	2.5	4
ADR/RID	3 F1					
IMDG	3/III			Solution Stab	ilized.	
R:	10-20/21/22-60-61			Contains Pyrio	dine, 2-Methoxyethanol and S	Sulfur Dioxide.
S:	23-36/37-45-53			One compon	ent reagent.	
	T					

Karl Fischer Reagent

2.0 mg H₂O/ml / 'BAKER ANALYZED' / Titrant for moisture determination

9 Q		

FLASHPOINT 21 °C	Titer (mg H ₂ O/ml)	min. 2.0	PRODUCT	PACKING	CONT.
NC CODE 3822 00 00			NO.		ВОХ
UN/ID NO. 1993			8834.1000	11	
ADR/RID 3 F1			8834.2500	2.5	
IMDG 3/III					
R: 10-20/21/22-60-61			Solution stal	bilized.	
S: 23-36/37-45-53			Contains Pyr	ridine, 2-Methoxyethano	l and Sulfur Dioxide.
T TOXIC			One compo	nent reagent.	

Karl Fischer Reagents

See for detailed information section Karl Fischer reagents, page 195



Find our up-to-date Product Literature at www.jtbaker.com/europe



Keros

Kerosene

8387 'BAKER'

1 l = 0.79 kg FLASHPOINT 72 °C **CAS NO.** 64742-96-7 **EINECS** 265-200-4

> NC CODE 2710 11 25 **EC NO.** 649 406 00 5

UN/ID NO. 1223 ADR/RID 3F1

IMDG 3/III

R: 65

S: 23-24-62



PRODUCT PACKING CONT. NO. вох 8387.5000 5 | EcoTainer 8387.9200 200 l

194-250°C

passes test

Kieselguhr 0116

'BAKER' / For filtration

NC CODE 3802 90 00

R: 40/20

S: 22



PACKING CONT. **PRODUCT** NO. вох 0116.1000 1 kg

DL-Lactic Acid 6069 'BAKER ANALYZED

► CH₃CHOHCOOH **1 l** = 1.25 kg

CAS NO. 598-82-3

EINECS 201-196-2

NC CODE 2918 11 00

R: 36/38 **S:** 24/25



Assay (C₃H₆O₃) 85.0-90.0% Chloride (CI) passes test Citric, Oxalic, Phosphoric or Tartaric Acids passes test Identification passes test Readily Carbonisable Substances passes test Residue after Ignition max. 0.05% Specific Rotation [α]²⁵ -0.05 to +0.05° Sugars passes test Sulfate (SO₄) passes test Trace Impurities (in ppm): Heavy Metals (as Pb) max. 10 Iron (Fe) max. 10

PRODUCT PACKING CONT. BOX 6069.4000 4 | Glass 6069.9060 60 l

DL-Lactic Acid

See for qualities meeting the requirements of cGMP, USP/NF, BP/Ph.Eur, JP or FCC, page 36

Boiling Range

Appearance

6017

Lactic Acid

'BAKER ANALYZED' / ACS

► CH₃CHOHCOOH

M = 90.08 g/mol **1 l** = 1.25 kg

FLASHPOINT 110 $^{\circ}\mathrm{C}$

CAS NO. 79-33-4

EINECS 201-196-2 NC CODE 2918 11 00

R: 38-41

S: 24-26-37/39



Meets ACS Specifications	
Assay	85.0-90.0%
Chloride (CI)	max. 0.001%
Residue after Ignition	max. 0.02%
Substances Darkened by H ₂ SO ₄	passes test
Sulfate (SO ₄)	max. 0.002%
Trace Impurities (in ppm):	
Heavy Metals (as Pb)	max. 5
Iron (Fe)	max. 5

PRODUCT	PACKING	CONT	
NO.		вох	
6017.1000	11		

Lanth



Lactic Acid 'BAKER

► CH₃CHOHCOOH **M** = 90.08 g/mol

1 l = 1.25 kg

FLASHPOINT 110 °C

CAS NO. 79-33-4

EINECS 201-196-2

NC CODE 2918 11 00

R: 38-41

S: 24-26-37/39



Assay (C ₃ H ₆ O ₃)	88.0-90.5%
Appearance	passes test
Calcium (Ca)	max. 200 ppm
Citric, oxalic or phosphoric acid	passes test
Ether-insoluble Substances	passes tes
Identification	passes test
Sugars and other Reducing Substances	passes test
Sulfated Ash	max. 0.05%
Sulfates (as SO ₄)	max. 200 ppm
Trace Impurities (in ppm):	
Heavy Metals (as Pb)	max. 10 ppm

PRODUCT	PACKING	CONT.
NO.		вох
6060.1000	11	

Lactic Acid

g the requirements of cGMP, USP/NF, BP/Ph.Eur, JP or FCC, page 36

Exceeds ACS Specifications



6060

Lactic Acid, Sodium Salt

See Sodium Lactate



D(+)-Lactose Monohydrate

'BAKER ANALYZED' Biochemical / ACS



C₁₂H₂₂O₁₁.H₂O

M = 360.32 g/mol **CAS NO.** 10039-26-6

EINECS 200-559-2 NC CODE 1702 11 00

Dextrose passes test Insoluble Matter max. 0.005% Residue after Ignition max. 0.03% Sucrose passes test Water (H₂O) 4.0-6.0% Trace Impurities (in ppm):

Dextrin, Starch max. 5 Heavy Metals (as Pb) max. 5 Iron (Fe) max. 5



Lactose, Monohydrate

See for qualities meeting the requirements of cGMP, USP/NF, BP/Ph.Eur, JP or FCC, page 36



Lanolin

See for qualities meeting the requirements of cGMP, USP/NF, BP/Ph.Eur, JP or FCC, page 36



Lanthanum 1% (w/v)
lanthanum chloride (99.999%) in 0.3M hydrochloric acid / 'BAKER INSTRA-ANALYZED' / Atomic Absorption Standard



NC CODE 3822 00 00

Chloride (CI) passes test Lanthanum Content (La) 0.95-1.05% **PRODUCT PACKING** CONT. NO. BOX 6947.0500 500 ml

Prepared by dissolution of high purity raw materials (min. 99.99% spectral purity). Assays are verified by ICP against standards traceable to NIST. Standard Reference Material numbers (SRM) are printed on each label.

Lanthanum 1% (w/v)

lanthanum nitrate (99.999%) in 0.3M nitric acid / 'BAKER INSTRA-ANALYZED' / Atomic Absorption Standard



NC CODE 3822 00 00 R: 36/38



Lanthanum (La) 0.95-1.05% Nitrate (NO₃) passes test

PACKING CONT. **PRODUCT** вох NO. 500 ml 6948.0500

Prepared by dissolution of high purity raw materials (min. 99.99% spectral purity). Assays are verified by ICP against standards traceable to NIST. Standard Reference Material numbers (SRM) are printed on each label.



Lanth

Lanthanum 1000 μg/ml

6949 0.10% (w/v) / (Matrix: 1% nitric acid) / 'BAKER INSTRA-ANALYZED' / Atomic Absorption Standard Lanthanum (La)

► La **M** = 138.90 g/mol

NC CODE 3822 00 00

R: 36/38 **S**: 26

PRODUCT PACKING CONT. NO. вох 6949.0100 100 ml 6949.0500 500 ml

Prepared by dissolution of high purity raw materials (min. 99.99% spectral purity). Assays are verified by ICP against standards traceable to NIST. Standard Reference Material numbers (SRM) are printed on each label.

Lanthanum Chloride 7-Hydrate

0383 'BAKER ANALYZED' / for Atomic Absorption / Flame Photometric Methods / ACS

► LaCl₃.7H₂O Exceeds ACS Specifications **M** = 371.38 g/mol 64.5-70.0% Assay (LaCl₂) **CAS NO.** 10099-58-8 Calcium (Ca) max. 0.001% **EINECS** 233-237-5 Insoluble Matter max. 0.01% NC CODE 2846 90 00

Trace Impurities (in ppm): Magnesium (Mg) max. 1

998-1002 μg/ml

PRODUCT	PACKING	CONT.
NO.		вох
0383.0025	25 g Glass	6
0383.0100	100 g	

Lanthanum Nitrate Hexahydrate

0597 'BAKER' / for Atomic Absorption / Flame Photometric Methods

► La(NO₃)₃.6H₂O **M** = 433.02 g/mol

> **CAS NO.** 10277-43-7 **EINECS** 233-238-0 NC CODE 2846 90 00

UN/ID NO. 1477 ADR/RID 5.1 02 IMDG 5.1/II

> R: 41-8 **S:** 26-39



Calcium (Ca) (by AAS) max. 0.001%

NO.

PRODUCT PACKING CONT. вох 0597.0100 100 g

Lanthanum Oxide

► La₂O₃ **M =** 325.84 g/mol **CAS NO.** 1312-81-8 **EINECS** 215-200-5

NC CODE 2846 90 00

4991

Certificate Provided Reporting Actual Lot Analysis Actual Lot Analysis Lot. No. B32431 Assay (by EDTA titrn.) 99.3% Loss on Ignition at 1150°C 0.6% Particulate Matter < 0.0001%

Metallic Impurities in parts per million ($\mu g/g$): < 1 Aluminium (AI) Calcium (Ca) < 0.5 Copper (Cu) < 1 Iron (Fe) < 1 Lead (Pb) < 10 Magnesium (Mg) < 0.3 Potassium (K) < 0.5 Rare Earths < 10 < 1 Silver (Ag)

< 10 Tin (Sn) Non-Metallic Impurities in parts per million ($\mu g/g$): Halide (as CI) Silicon (Si) < 1

< 0.3

5

3

Sodium (Na)

Sulfur Compounds (as SO₄)

PRODUCT PACKING CONT. NO. BOX 4991.0100 100 g PE

For Laboratory, Research or Manufacturing Use.

Lauryl Alcohol

(Matrix: 2% nitric acid) / 'BAKER INSTRA-ANALYZED' / Plasma Standard

5765

► Pb

M = 207.21 g/mol

NC CODE 3822 00 00

R: 36/38

S: 26



Certificate Provided Reporting Actual Lot Analysis

Lead (Pb) 998-1002 μg/ml

998-1002 μg/ml

998-1002 μg/ml

PRODUCT PACKING CONT. NO. вох 5765.0100 100 ml

Prepared from the highest purity raw material available, generally greater than 99.999% spectral purity. The content of the solution is confirmed to be accurate to within \pm 0.2%. Typically 1000 μ g/ml.

The certificate of analysis provided reports actual lot analysis. The certificate also lists the trace impurities.

Lead (Pb)

Lead (Pb)

6930

▶ Pb

M = 207.20 g/mol

NC CODE 3822 00 00

R: 36/38

S: 26-37



PRODUCT PACKING CONT. NO. BOX 6930.0100 100 ml 6930.0500 500 ml

Prepared by dissolution of high purity raw materials (min. 99.99% spectral purity). Assays are verified by ICP against standards traceable to NIST. Standard Reference Material numbers (SRM) are printed on each label.

Lead 1000 μg/ml'BAKER ANALYZED' / Atomic Absorption Standard

6813

▶ Ph

M = 207.21 g/mol

NC CODE 3822 00 00

R: 36/38

S: 26-37



CONT. PRODUCT PACKING NO. вох 6813.0100 100 ml 6813.0500 500 ml

Lead(II)nitrate in nitric acid 0.5 mol/l.

Lead 10000 μg/ml

(Matrix: 2% nitric acid) / 'BAKER INSTRA-ANALYZED' / Plasma Standard

5732

► Pb

M = 207.21 g/mol

NC CODE 3822 00 00

R: 22-33-36/38-52/53-61

S: 26-45-61



Certificate Provided Reporting Actual Lot Analysis 9980-10020 μg/ml Lead (Pb)

PRODUCT PACKING CONT. NO. 5732.0100 100 ml

Prepared from the highest purity raw material available, generally greater than 99.999% spectral purity. The content of the solution is confirmed to be accurate to within 0.2 %. Typically 10000 µg/ml.

The certificate of analysis provided reports actual lot analysis. The certificate also lists the trace impurities.

Lead Acetate Basic Anhydrous

'BAKER ANALYZED' / For sugar analysis / ACS

0142

CAS NO. 1335-32-6 **EINECS** 215-630-3 NC CODE 2915 29 00 EC NO. 82 007 00 9 **UN/ID NO.** 1616 **ADR/RID** 6.1 T5 IMDG 6.1/III **R:** 33-40-48/22-50/53-61-62 **S:** 45-53-60-61

dangerous for the



Exceeds ACS Specifications	
Basic Lead (PbO)	min. 33.0%
Calcium (Ca)	max. 0.01%
Chloride (CI)	max. 0.002%
Copper (Cu)	max. 0.002%
Insoluble in Dilute Acetic Acid	max. 0.02%
Insoluble in Water	max. 1.0%
Iron (Fe)	max. 0.002%
Loss on Drying at 105°C	max. 1.5%
Nitrate and Nitrite (as NO ₃)	max. 0.003%
Potassium (K)	max. 0.02%
Sodium (Na)	max. 0.05%

PRODUCT NO.	PACKING	CON1 BO
0142.1000	1 kg	

LeadI

Lead(II) Acetate Trihydrate

'BAKER ANALYZED' / ACS

► Pb(CH₃COO)₂.3H₂O **M** = 379.33 g/mol **CAS NO.** 6080-56-4 **EINECS** 206-104-4 NC CODE 2915 29 00 **EC NO.** 82 005 00 8 **UN/ID NO.** 1616

> **ADR/RID** 6.1 T5 IMDG 6.1/III

> > **R:** 33-48/22-50/53-61-62

S: 45-53-60-61

dangerous for the

environment



Exceeds ACS Specifications	
Assay	99.0-103.0%
Calcium (Ca)	max. 0.005%
Copper (Cu)	max. 0.002%
Insoluble Matter	max. 0.01%
Iron (Fe)	max. 0.001%
Nitrate and Nitrite (as NO ₃)	max. 0.005%
pH of 5% Solution at 25°C	5.5-6.5
Potassium (K)	max. 0.005%
Sodium (Na)	max. 0.01%
Trace Impurities (in ppm):	
Chloride (CI)	max 5

min. 99.0%

max. 0.005%

max. 0.01%

max. 0.001%

max. 0.005%

max. 0.01%

Assay (by EDTA titrn.)

Copper (Cu)

Nitrate (NO₃)

Sulfate (SO₄)

Iron (Fe)

Insoluble Matter

PRODUCT	RODUCT PACKING	CONT.		
NO.		вох		
0141.0250	250 g			

Lead(II) Chloride

0145

► PbCl₂

M = 278.10 g/mol **CAS NO.** 7758-95-4

EINECS 231-845-5

NC CODE 2827 39 80 **EC NO.** 82 001 00 6

UN/ID NO. 2291 **ADR/RID** 6.1 T5

IMDG 6.1/III

R: 20/22-33-50/53-61-62

S: 45-53-60-61





PRODUCT PACKING CONT. BOX 0145.0500 500 g

Lead(II) Nitrate

► Pb(NO₃)₂

M = 331.20 g/mol **CAS NO.** 10099-74-8 **EINECS** 233-245-9 NC CODE 2834 29 80 **EC NO.** 82 001 00 6 **UN/ID NO.** 1469

ADR/RID 5.1 OT2

IMDG 5.1/II

R: 20/22-33-50/53-61-62

S: 45-53-60-61





min. 99.0%
111111111111111111111111111111111111111
max. 0.005%
max. 0.001%
max. 0.002%
max. 0.005%
max. 0.001%
3.5-4.2
max. 0.005%
max. 0.02%

Exceeds ACS Specifications

PRODUCT	PACKING	CONT.
NO.		вох
0148.0100	100 g	
0148.9050	50 kg	

Mallinckrodt Baker's chemistry is Part of a pure process™.

Lead(II) Oxide (Litharge) powder / 'BAKER ANALYZED' / ACS

2338

CONT. вох

M = 223.19 g/mol Assay		PRODUCT	PACKING
220.10 9,11101	min. 99.0%	NO.	
CAS NO. 1317-36-8 Average Particle Diameter, μm (APD) a	act. value reported	2338.0500PE	500 ml HDPE
EINECS 215-267-0 Bulk Density (g/cc)(typical) a	act. value reported		
NC CODE 2824 10 00 Calcium (Ca)	max. 0.005%		
EC NO. 82 001 00 6 Chloride (CI)	max. 0.002%		
UN/ID NO. 2291 Copper (Cu)	max. 0.005%		
ADR/RID 6.1 T5 Insoluble in CH ₃ COOH	max. 0.02%		
IMDG 6.1/III Nitrate (NO ₃)	max. 0.01%		
R: 20/22-33-50/53-61-62 Potassium (K)	max. 0.005%		
S: 45-53-60-61 Silicon (Si)	act. value reported		
N Sodium (Na)	max. 0.02%		
Specific Surface Area, m²/g (typical) a	act. value reported		
dangerous toxic Mesh (Wet Screen Analysis):			
	act. value reported		
Trace Impurities (in ppm):			
Iron (Fe)	max. 5		

Lead Subacetate



3879

Leishman

11 = 0.79 kg Leihsmans stain for differential blood picture

staining

Silver (Ag)

PRODUCT PACKING CONT. NO. вох 3879.1000 1 | Glass 3879.2500 2.5 | Glass

FLASHPOINT 11C °C **UN/ID NO.** 1230 ADR/RID 3 FT1 IMDG 3/II

R: 11-23/24/25-39/23/24/25

S: 16-36/37-45-7





Leishman Stain

CAS NO. 12627-53-1 **PRODUCT PACKING** Leihsmans stain for differential blood picture

EINECS 235-732-1 **R:** 36 **S:** 26



staining

NO. 3813.0025 25 g Glass

max. 5

CONT. вох

3813

L-Leucine

See for qualities meeting the requirements of cGMP, USP/NF, BP/Ph.Eur, JP or FCC, page 36



Levulose

See D(-)Fructose



Lithium 1000 μ g/ml (Matrix: 2% nitric acid) / 'BAKER INSTRA-ANALYZED' / Plasma Standard

5766

▶ Li

M = 6.94 g/mol

NC CODE 3822 00 00

R: 36/38



S: 26

Certificate Provided Reporting Actual Lot Analysis 998-1002 μg/ml Lithium (Li)

NO. 5766.0100 100 ml

PACKING

PRODUCT

CONT. вох

Prepared from the highest purity raw material available, generally greater than 99.999% spectral purity. The content of the solution is confirmed to be accurate to within \pm 0.2%. Typically 1000 μ g/ml.

The certificate of analysis provided reports actual lot analysis. The certificate also lists the trace impurities.

Lithi

Lithium 1000 μg/ml

6931 (Matrix: 2% nitric acid) / 'BAKER INSTRA-ANALYZED' / Atomic Absorption Standard

▶ Li

M = 6.94 g/mol NC CODE 3822 00 00

R: 36/38 **S:** 26-37

Lithium (Li) 998-1002 μg/ml

PRODUCT CONT. **PACKING** вох 6931.0100 100 ml 6931.0500 500 ml

Prepared by dissolution of high purity raw materials (min. 99.99% spectral purity). Assays are verified by ICP against standards traceable to NIST. Standard Reference Material numbers (SRM) are printed on each label.

Lithium 1000 μg/ml 'BAKER ANALYZED' / Atomic A

▶ Li

M = 6.94 g/molNC CODE 3822 00 00

998-1002 μg/ml Lithium (Li)

CONT **PRODUCT PACKING** BOX 6814.0100 100 ml 6814.0500 500 ml

Lithium chloride in water.

Lithium 10000 μg/ml

5733 (Matrix: 7% nitric acid) / 'BAKER INSTRA-ANALYZED' / Plasma Standard

▶ Li

M = 6.94 g/mol

NC CODE 3822 00 00

UN/ID NO. 2031

ADR/RID 8 CO1 IMDG 8/II

R: 34

S: 20-23-26-36/37/39-45



Certificate Provided Reporting Actual Lot Analysis Lithium (Li)

9980-10020 $\mu g/ml$

PRODUCT PACKING CONT. вох NO. 5733.0100 100 ml

Prepared from the highest purity raw material available, generally greater than 99.999% spectral purity. The content of the solution is confirmed to be accurate to within 0.2 %. Typically 10000 $\mu g/ml$.

The certificate of analysis provided reports actual lot analysis. The certificate also lists the trace impurities.

Lithium meta-Borate

BAKER ANALYZED' / ACS

► LiBO₂

0382

M = 49.75 g/mol **CAS NO.** 13453-69-5 **EINECS** 236-631-5 NC CODE 2840 20 90

Exceeds ACS Specifications. Meets Reagents Specifications for testing USP/NF monographs

Assay (LiBO ₂)	98.0-102.0%
Aluminium (AI)	max. 0.001%
Bulk Density	min. 0.25 g/ml
Calcium (Ca)	max. 0.01%
Heavy Metals (as Pb)	max. 0.001%
Insoluble Matter	max. 0.01%
Iron (Fe)	max. 0.001%
Loss on Fusion at 950°C	max. 2.0%
Phosphorus Compounds (as PO ₄)	max. 0.004%
Potassium (K)	max. 0.005%
Silicon (Si)	max. 0.005%
Sodium (Na)	max. 0.005%
Trace Impurities (in ppm):	
Magnesium (Mg)	max. 5

PRODUCT PACKING CONT. вох NO. 0382.0100 100 g

The formula weight of this reagent is likely to deviate from the value cited above, since the natural distribution of 6Li and 7Li isotopes is often altered in current sources of lithium compounds.

www.jtbaker.com/europe

Lithium Carbonate

'BAKER ANALYZED' 0155

► Li₂CO₃

M = 73.89 g/mol

CAS NO. 554-13-2

EINECS 209-062-5 NC CODE 2836 91 00

R: 22-38

S: 22



Assay	min. 99.0%
Calcium (Ca)	max. 0.01%
Chloride (CI)	max. 0.005%
Heavy Metals (as Pb)	max. 0.002%
Insoluble in HCI	max. 0.01%
Iron (Fe)	max. 0.002%
Potassium (K)	max. 0.005%
Sodium (Na)	max. 0.01%
Total Sulfur (as SO ₄)	max. 0.2%
Trace Impurities (in ppm):	
Ammonium (NH ₄)	max. 5
Nitrate (NO ₃)	max. 5

PRODUCT NO.	PACKING	CONT BO
0155.0100	100 g	
0155.9050	50 kg	

Lithium Chloride

'BAKER ULTRAPURE BIOREAGENT'

4002

► LiCl

M = 42.39 g/mol

CAS NO. 7447-41-8

EINECS 231-212-3

NC CODE 2827 39 80

R: 22-36/38

S: 24



Assay (argentometric titrn.)	min. 99.0%
Calcium (Ca) (by FES)	max. 0.01%
DNase Activity	none detected
Heavy Metals (as Pb)	max. 0.002%
Insoluble Matter	max. 0.01%
Iron (Fe)	max. 0.001%
Loss on Drying at 105°C	max. 1.0%
Potassium (K) (by FES)	max. 0.01%
Protease Activity	none detected
RNase Activity	none detected
Sodium (Na) (by FES)	max. 0.2%
Titrable Base (meq/g)	max. 0.008

PRODUCT	PACKING	CONT.
NO.		вох
4002.0500	500 g	
4002.2500	2.5 kg	

Lithium Chloride

'BAKER ANALYZED' / ACS

0156

► LiCl

M = 42.39 g/mol

CAS NO. 7447-41-8 **EINECS** 231-212-3

NC CODE 2827 39 80

R: 22-36/38

S: 24



Assay	min. 99%
Barium (Ba)	max. 0.003%
Calcium (Ca)	max. 0.01%
Heavy Metals (as Pb)	max. 0.002%
Insoluble Matter	max. 0.01%
Iron (Fe)	max. 0.001%
Loss on Drying at 105°C	max. 1.0%
Nitrate (NO ₃)	max. 0.001%
Potassium (K)	max. 0.01%
Sodium (Na)	max. 0.20%
Sulfate (SO ₄)	max. 0.01%
Titrable Base (meq/g)	max. 0.008

PRODUCI	PACKING	CONT.
NO.		вох
0156.0100	100 g	6
0156.1000	1 kg	
0156.5000	5 kg	
0156.9050	50 kg	

Lithium Chloride

0157

► LiCl

M = 42.39 g/mol

CAS NO. 7447-41-8

EINECS 231-212-3

NC CODE 2827 39 80

R: 22-36/38

S: 24



Assay (argentometric titrn.)	min. 98.0%
Heavy Metals (as Pb)	max. 0.01%
Insoluble Matter	max. 0.05%
pH of 5% Solution at 25°C	6.0-9.0
Sodium (Na)	max. 0.2%

PRODUCI	PACKING	CONT.
NO.		вох
0157.0250	250 g	
0157.9050	50 kg	



Lithi

Lithium Nitrate 'BAKER ANALYZED

► LiNO₃ **M** = 68.94 g/mol **CAS NO.** 7790-69-4

> **EINECS** 232-218-9 NC CODE 2834 29 80

UN/ID NO. 2722

ADR/RID 5.1 02

IMDG 5.1/III

R: 8 **S:** 17-24/25

Assay (acidimetric) min. 97.0% Barium (Ba) max. 0.002% Calcium (Ca) (by FES) max. 0.02% Chloride (CI) max. 0.002% Heavy Metals (as Pb) max. 0.001% Insoluble Matter max. 0.01% Iron (Fe) max. 0.001% Potassium (K) (by FES) max. 0.05% Sodium (Na) (by FES) max. 0.05% Sulfate (SO₄) max. 0.1%

PRODUCT	PACKING	CONT.
NO.		вох
0158.0500	500 g	6

Lithium Tetraborate

1465

► Li₂B₄O₇

M = 169.12 g/mol

CAS NO. 12007-60-2 **EINECS** 234-514-3

NC CODE 2840 20 90

Assay (dried basis) min. 99.5% Loss on Drying at 285°C max. 1% Metallic and Non-metallic Impurities (by DC-arc

Spectrography) (in ppm): Aluminium (AI)

Calcium (Ca) max. 100 Chromium (Cr) max. 20 Cobalt (Co) max. 50 Copper (Cu) max. 20 Iron (Fe) max. 30 Lead (Pb) max. 20 Magnesium (Mg) max. 10 Manganese (Mn) max. 50 Nickel (Ni) max. 20 **PRODUCT PACKING** BOX 1465.0500 500 g

Flux Grade.

max. 30

max. 80

max. 150

min. 99%

L-Lysine Hydrochloride

r qualities meeting the requirements of cGMP, USP/NF, BP/Ph.Eur, JP or FCC, page 36

Assay

Potassium (K) Silicon (Si)

Sodium (Na)

Magnesium

M = 24.30 g/mol

CAS NO. 7439-95-4

EINECS 231-104-6

NC CODE 8104 30 00

EC NO. 12 002 00 9 **UN/ID NO.** 1869

ADR/RID 4.1 F3

IMDG 4.1/III

R: 11-15 **S:** 43A-7/8

PRODUCT PACKING BOX 1120.0500 500 g

Innovation is principal to our business.

Magnesium 1000 μg/ml (Matrix: 2% nitric acid) / 'BAKER INSTRA-ANALYZED' / Plasma Standard

5767

► Mg

M = 24.31 g/mol

NC CODE 3822 00 00

R: 36/38 **S**: 26

Certificate Provided Reporting Actual Lot Analysis

Magnesium (Mg) 998-1002 μg/ml

PRODUCT NO.	PACKING	CONT.	
		вох	
5767 0100	100 ml		

Prepared from the highest purity raw material available, generally greater than 99.999% spectral purity. The content of the solution is confirmed to be accurate to within \pm 0.2%. Typically 1000 μ g/ml.

The certificate of analysis provided reports actual lot analysis. The certificate also lists the trace impurities.

Magnesium 1000 µg/ml (Matrix: 2% nitric acid) / 'BAKER INSTRA-ANALYZED' / Atomic Absorption Standard

Magnesium (Mg)

6932



M = 24.31 g/mol

NC CODE 3822 00 00

R: 36/38

S: 26-37



PRODUCT NO.	PACKING	CONT.
		вох
6932.0100	100 ml	
6932.0500	500 ml	

Prepared by dissolution of high purity raw materials (min. 99.99% spectral purity). Assays are verified by ICP against standards traceable to NIST. Standard Reference Material numbers (SRM) are printed on each label.

Magnesium 1000 μg/ml 'BAKER ANALYZED' / Atomic Absorpt

6815

► Mg

M = 24.31 g/mol NC CODE 3822 00 00

R: 36/38

S: 26-37



Magnesium (Mg)

998-1002 $\mu g/ml$

998-1002 μg/ml

PRODUCT	PACKING	CONT
NO.		BOX
6815.0100	100 ml	
6815.0500	500 ml	

Magnesium nitrate in nitric acid 0.5 mol/l.

Magnesium 10000 µg/ml

(Matrix: 2% nitric acid) / 'BAKER INSTRA-ANALYZED' / Plasma Standard

5734

► Mg

M = 24.31 g/mol

NC CODE 3822 00 00

R: 36/38

S: 24/25-26-37



Certificate Provided Reporting Actual Lot Analysis

9980-10020 μg/ml Magnesium (Mg)

PRODUCT PACKING CONT. NO. 5734.0100 100 ml

Prepared from the highest purity raw material available, generally greater than 99.999% spectral purity. The content of the solution is confirmed to be accurate to within 0.2 %. Typically 10000 µg/ml.

The certificate of analysis provided reports actual lot analysis. The certificate also lists the trace impurities.

PACKING

PRODUCT

max. 5

max. 5

Magnesium Acetate Tetrahydrate

'BAKER ANALYZED' / ACS

0160

CONT.

1vig(Cr13COO)2.4r12O		
M =	214.46 g/mol	
CAS NO.	16674-78-5	
EINECS	205-554-9	

NC CODE 2915 29 00

Ma(CH COO) 4H O

Meets ACS Specifications	
Assay	98.0-102.0%
Barium (Ba)	max. 0.001%
Calcium (Ca)	max. 0.01%
Chloride (CI)	max. 0.001%
Insoluble Matter	max. 0.005%
Manganese (Mn)	max. 0.001%
Potassium (K)	max. 0.005%
Sodium (Na)	max. 0.005%
Strontium (Sr)	max. 0.005%
Sulfate (SO ₄)	max. 0.005%
Trace Impurities (in ppm):	

Heavy Metals (as Pb)

Iron (Fe)

NO.		ВО
0160.0250	250 g	

Magne



Magnesium Carbonate

or qualities meeting the require nents of cGMP, USP/NF, BP/Ph.Eur, JP or FCC, page 36

0161

Magnesium Carbonate n-Hydrate

► 4MgCO₃.Mg(OH)₂.nH₂O **CAS NO.** 12125-28-9 **EINECS** 235-192-7 NC CODE 2836 99 11

Assay (as MgO) 40.0-45.0% Appearance of solution passes test Arsenic (As) max. 2 ppm Calcium (Ca) max. 0.75% Chlorides (as CI) max. 0.07% Heavy Metals (as Pb) max. 20 ppm Iron (Fe) max. 400 ppm Soluble Substances max. 1.0% max. 0.05% Substances Insoluble in Acetic Acid Sulfates (as SO₄) max. 0.3%

B0)

4003

Magnesium Chloride Hexahydrate

BAKER ULTRAPURE BIOREAGENT

▶ MgCl₂.6H₂O

M = 203.30 g/mol **CAS NO.** 7791-18-6

EINECS 232-094-6 NC CODE 2827 31 00

99.0 - 102.0% Assay Ammonium (NH₄) max. 0.002% Calcium (Ca) max. 0.01% **DNase Activity** none detected max. 0.005% Insoluble Matter max. 0.005% Potassium (K) Protease Activity none detected RNase Activity none detected

Sodium (Na) Trace Impurities (in ppm): Heavy Metals (as Pb)

PRODUCT PACKING CONT BOX 4003.0500 500 g

Iron (Fe)

max. 5 max. 5 Manganese (Mn) max. 5

max. 0.005%

0162

Magnesium Chloride Hexahydrate

'BAKER ANALYZED' / ACS

▶ MgCl₂.6H₂O M = 203.30 g/mol **CAS NO.** 7791-18-6 **EINECS** 232-094-6

NC CODE 2827 31 00

Meets ACS Specifications 99.0-102.0% Assay Ammonium (NH₄) max. 0.002% Barium (Ba) max. 0.005% Calcium (Ca) max. 0.01% Insoluble Matter max. 0.005% Nitrate (NO₃) max. 0.001% Potassium (K) max. 0.005% max. 0.005% Sodium (Na)

Strontium (Sr) max. 0.005% Sulfate (SO₄) max. 0.002%

Trace Impurities (in ppm):

Heavy Metals (as Pb) max. 5 Iron (Fe) max. 5 Manganese (Mn) max. 5 Phosphate (PO₄) max. 5

PRODUCT	PACKING	CONT.
NO.		вох
0162.0250	250 g	6
0162.1000	1 kg	6

Magnesium Chloride Hexahydrate

0163

▶ MgCl₂.6H₂O

M = 203.30 g/mol **CAS NO.** 7791-18-6

EINECS 232-094-6 NC CODE 2827 31 00

Assay min. 98.5% Heavy Metals (as Pb) max. 0.001% Insoluble Substances max. 0.01% pH of 5% Solution at 25°C 4.5-7.0

PRODUCT PACKING CONT. NO. BOX 0163.1000 1 ka

Magnesium Chloride, 6-Hydrate

nts of cGMP, USP/NF, BP/Ph.Eur, JP or FCC, page 36



Magnesium Hydroxide

See for qualities meeting the requirements of cGMP, USP/NF, BP/Ph.Eur, JP or FCC, page 36

Magne



Magnesium Nitrate Hexahydrate 'BAKER ANALYZED' / ACS

0164

Mg(NO ₃) ₂ .6H ₂ O		
M =	256.41 g/mol	
CAS NO.	13446-18-9	
EINECS	233-826-7	
NC CODE	2834 29 80	
UN/ID NO.	1474	
ADR/RID	5.1 O2	
IMDG	5.1/III	
R:	8	
S:	24/25	

Exceeds ACS Specifications	
Assay	98.0-102.0%
Ammonium (NH ₄)	max. 0.002%
Barium (Ba)	max. 0.002%
Calcium (Ca)	max. 0.01%
Chloride (CI)	max. 0.001%
Insoluble Matter	max. 0.005%
pH of 5% Solution at 25°C	5.0-7.0
Potassium (K)	max. 0.005%
Sodium (Na)	max. 0.005%
Strontium (Sr)	max. 0.005%
Sulfate (SO ₄)	max. 0.005%
Trace Impurities (in ppm):	
Heavy Metals (as Pb)	max. 5
Iron (Fe)	max. 5
Manganese (Mn)	max. 5
Phosphate (PO ₄)	max. 5

PRODUCT	PACKING	CONT
NO.		BO)
0164.0500	500 g	Ć

Magnesium Oxide 'BAKER ANALYZED' / ACS

MgO	

M = 40.30 g/mol **CAS NO.** 1309-48-4 **EINECS** 215-171-9 NC CODE 2519 90 10 **S:** 22

Exceeds ACS Specifications. Meets Re	eagents
Specifications for testing USP/NF mor	nographs
Barium (Ba)	max. 0.005%
Calcium (Ca)	max. 0.05%
Chloride (CI)	max. 0.01%
Heavy Metals (as Pb)	max. 0.003%
Insoluble in Dilute HCl	max. 0.02%
Iron (Fe)	max. 0.01%
Loss on Ignition	max. 2.0%
Manganese (Mn)	max. 5 ppm
Nitrate (NO ₃)	max. 0.005%
Potassium (K)	max. 0.005%
Sodium (Na)	max. 0.5%
Strontium (Sr)	max. 0.005%
Sulfate and Sulfite (as SO ₄)	max. 0.02%
Water Soluble Substances	max. 0.4%
Product Information (not specifications	s):
Average Particle Diameter, µm (APD) (by	
Sedigraph)(typical)	4

PRODUCI	PACKING	CONT
NO.		BO
1122.1000	1 kg	
1122.9050	50 kg	

Magnesium Oxide

quirements of cGMP, USP/NF, BP/Ph.Eur, JP or FCC, page 36

Bulk Density (g/cc)(typical) Specific Surface Area, m²/g (typical)



Magnesium Perchlorate



Magnesium Silicate, Activated



Magnesium Stearate

See for qualities meeting the requirements of cGMP, USP/NF, BP/Ph.Eur, JP or FCC, page 36



www.jtbaker.com/europe



Magne

Magnesium Sulfate Anhydrous

0366 'BAKER ANALYZED

► MgSO₄

M = 120.37 g/molCAS NO. 7487-88-9EINECS 231-298-2NC CODE 2833 21 00

Assay (after ignition) min. 99.0% Ammonium (NH₄) max. 0.004% Calcium (Ca) max. 0.06% Chloride (CI) max. 0.001% Heavy Metals (as Pb) max. 0.001% Iron (Fe) max. 0.001% Loss on Ignition max. 2.0% max. 0.001% Manganese (Mn) Nitrate (NO₃) max. 0.004% Trace Impurities (in ppm): Arsenic (As) max. 1

PRODUCT	PACKING	CONT.
NO.		вох
0366.0500	500 g	6



Magnesium Sulfate Anhydrous

See for qualities meeting the requirements of cGMP, USP/NF, BP/Ph.Eur, JP or FCC, page 36

0168

Magnesium Sulfate Heptahydrate

'BAKER ANALYZED' / ACS

MgSO₄.7H₂O

M = 246.47 g/mol

CAS NO. 10034-99-8

EINECS 231-298-2

NC CODE 2833 21 00

Exceeds ACS Specifications 99.0-102.0% Assay (by EDTA titrn.) Ammonium (NH₄) max. 0.002% max. 0.02% Calcium (Ca) Insoluble Matter max. 0.005% Nitrate (NO₃) max. 0.002% 5.0-7.0 pH of 5% Solution at 25°C Potassium (K) max. 0.005% Sodium (Na) max. 0.005% Strontium (Sr) max. 0.005% Trace Impurities (in ppm): Chloride (CI) max. 5 Heavy Metals (as Pb) max. 5 Iron (Fe) max. 5

30X
6

CONT.

PACKING

PRODUCT

max. 5

0169

Magnesium Sulfate Heptahydrate

BAKER'

► MgSO₄.7H₂O M = 246.47 g/mol CAS NO. 10034-99-8 EINECS 231-298-2 NC CODE 2833 21 00

Assay (dried basis) 99.0-100.5% Acidity or Alkalinity passes test Appearance of solution passes test Arsenic (As) max. 2 ppm Chlorides (as CI) max. 300 ppm Heavy Metals (as Pb) max. 10 ppm Identification passes test Iron (Fe) max. 20 ppm 48.0-52.0% Loss on Drying

PRODUCT	PACKING	CONT.
NO.		вох
0169.5000	5 kg	4
0169.9050	50 kg	



Magnesium Sulfate, 7-Hydrate

See for qualities meeting the requirements of cGMP, USP/NF, BP/Ph.Eur, JP or FCC, page 36

Appearance

Manganese (Mn)

1325

Malachite Green Oxalate

'BAKER ANALYZED'

 $\begin{array}{c} \textbf{C}_{62}\textbf{H}_{54}\textbf{N}_{4}\textbf{O}_{12} \\ \textbf{M} = 927.00 \text{ g/mol} \\ \textbf{CAS NO.} \ \ \, 2437-29-8 \\ \textbf{EINECS} \ \ \, 219-441-7 \\ \textbf{NC CODE} \ \ \, 3204 \ \, 13 \ \, 00 \\ \textbf{UN/ID NO.} \ \ \, 2811 \\ \end{array}$

ADR/RID 6.1 T2 **IMDG** 6.1/III

R: 21/22 **S:** 24/25

X

passes test PRODUCT
NO.

 PRODUCT
 PACKING
 CONT.

 NO.
 BOX

 1325.0100
 100 g

C.I. 42000

Certified Stain for Use in Histology and Bacteriology.

2790

Maleic Acid

'BAKER ANALYZED



M = 116.07 g/mol **CAS NO.** 110-16-7

EINECS 203-742-5

NC CODE 2917 1990

EC NO. 607 095 00 3

UN/ID NO. 3261

ADR/RID 8 C4

IMDG 8/III

R: 22-36/37/38

S: 26-28-37



Assay (acidimetric)	min. 99.0%
Chloride (CI)	max. 0.002%
Heavy Metals (as Pb)	max. 0.001%
Iron (Fe)	max. 0.001%
Melting Point	132-135°C
Residue after Ignition	max. 0.1%
Sulfate (SO ₄)	max. 0.01%

PRODUCT	PACKING	CONT.
NO.		вох
2790.0500	500 g	

Malonic Acid

'BAKER ANALYZED

► CH₂(COOH)₂

M = 104.06 g/mol

CAS NO. 141-82-2

EINECS 205-503-0

NC CODE 2917 19 10

R: 22-36 **S:** 22-24



min. 99.0% Assay Heavy Metals (as Pb) max. 0.001% Iron (Fe) max. 0.005% Melting Range 133-134°C Residue after Ignition max. 0.3%

PRODUCT	PACKING	CONT.
NO.		вох
1124.0100	100 g	

Manganese 1000 μg/ml

(Matrix: 2% nitric acid) / 'BAKER INSTRA-ANALYZED' / Plasma Standard

5793

1124

► Mn

M = 54.94 g/mol

NC CODE 3822 00 00

R: 36/37/38

S: 26-28-45



Certificate Provided Reporting Actual Lot Analysis

998-1002 μg/ml Manganese (Mn)

PRODUCT PACKING CONT. вох NO. 5793.0100 100 ml

Prepared from the highest purity raw material available, generally greater than 99.999% spectral purity. The content of the solution is confirmed to be accurate to within \pm 0.2%. Typically 1000 $\mu g/ml$

The certificate of analysis provided reports actual lot analysis. The certificate also lists the trace impurities

Manganese 1000 μg/ml

cid) / 'BAKER INSTRA-ANALYZED' / Atomic Absorption Standard 0.10% (w/v) / (Matrix: 2% nit

6933

CONT.

вох

► Mn

M = 200.59 g/mol

NC CODE 3822 00 00

R: 36/38 **S:** 26-37



Manganese (Mn) $998\text{-}1002~\mu\text{g/ml}$

PRODUCT NO.

6933.0100 100 ml 6933.0500 500 ml

PACKING

Prepared by dissolution of high purity raw materials (min. 99.99% spectral purity). Assays are verified by ICP against standards traceable to NIST. Standard Reference Material numbers (SRM) are printed on each label.

Manganese 1000 μ g/ml

'BAKER ANALYZED' / Atomic A bsorption Standard 6816

► Mn

M = 54.94 g/mol

NC CODE 3822 00 00

R: 36/38

S: 26-37



Manganese (Mn) 998-1002 $\mu g/ml$

PRODUCT NO.	PACKING	CONT.		
		вох		
6816.0100	100 ml			
6816 0500	500 ml			

Manganese(II)nitrate in nitric acid 0.5 mol/l.



Manga

Manganese 10000 $\mu g/ml$ (Matrix: 2% nitric acid) / 'BAKER INSTRA-ANALYZED' / Plasma Standard 5735

► Mn

M = 54.94 g/mol

NC CODE 3822 00 00

R: 36/38 **S**: 26



Certificate Provided Reporting Actual Lot Analysis

9980-10020 μg/ml Manganese (Mn)

PRODUCT CONT. **PACKING** NO. вох 5735.0100 100 ml

Prepared from the highest purity raw material available, generally greater than 99.999% spectral purity. The content of the solution is confirmed to be accurate to within 0.2 %. Typically 10000 $\mu\text{g/ml.}$

The certificate of analysis provided reports actual lot analysis. The certificate also lists the trace impurities.

Manganese(II) Carbonate 'BAKER ANALYZED'

0172

► MnCO₃

M = 114.95 g/mol **CAS NO.** 598-62-9 EINECS 209-942-9 NC CODE 2836 99 18

S: 22-24/25

Assay (as Mn)	min. 43.0%
Average Particle Diameter, μm (APD)	max. 3
Calcium (Ca)	max. 0.01%
Chloride (CI)	max. 0.02%
Insoluble in HCl	max. 0.01%
Iron (Fe)	max. 0.02%
Magnesium (Mg)	max. 0.01%
Other Heavy Metals (as Pb)	max. 0.005%
Potassium (K)	max. 0.01%
Sodium (Na)	max. 0.02%
Sulfate (SO ₄)	max. 0.005%
Zinc (Zn)	max. 0.05%

PACKING	CONT.
	вох
500 g	
12 kg	
	500 g

0173

Manganese(II) Chloride Tetrahydrate

'BAKER ANALYZED' / ACS

MnCl₂.4H₂O **M** = 197.91 g/mol **CAS NO.** 13446-34-9 **EINECS** 231-869-6

NC CODE 2827 39 90

R: 22



Meets ACS Specifications	
Assay (by EDTA titrn.)	98.0-101.0%
Calcium (Ca)	max. 0.005%
Insoluble Matter	max. 0.005%
Magnesium (Mg)	max. 0.005%
pH of 5% Solution at 25°C	3.5-6.0
Potassium (K)	max. 0.01%
Sodium (Na)	max. 0.05%
Sulfate (SO ₄)	max. 0.005%
Zinc (Zn)	max. 0.005%
Trace Impurities (in ppm):	
Heavy Metals (as Pb)	max. 5
Iron (Fe)	max. 5

CONT.

вох

6

PACKING

500 g

PRODUCT

0173.0500

NO.



Manganese Dioxide

ee Manganese(IV) Oxide

Manganese(IV) Oxide

0170 BAKER ANALYZED

► MnO₂

M = 86.94 g/mol **CAS NO.** 1313-13-9

EINECS 215-202-6 NC CODE 2820 10 00

EC NO. 25 001 00 3 **R:** 20/22

S: 25



Assay (by EDTA titrn.)	min. 99.0%
Chloride (CI)	max. 0.01%
Insoluble in HCI	max. 0.03%
Iron (Fe)	max. 0.05%
Nitrate (NO ₃)	max. 0.05%
Sulfate (SO ₄)	max. 0.08%

PRODUCT	PACKING	CONT.
NO.		ВОХ
0170.1000	1 kg	

0174

Manganese(II) Sulfate Monohydrate 'BAKER ANALYZED' / ACS

98.0-101.0%

max. 0.005%

max. 0.001%

max. 0.002%

max. 0.01%

max. 0.001%

10.0-12.0%

max. 0.005%


max. 0.005%

max. 0.005% max. 0.005%

max. 5

max. 5

max. 5



PACKING	CONT
	B0
1 kg	

Mannitol

dangerous for the

environment

'BAKER ANALYZED' / ACS

harmful

1125

1998

► HOCH₂(CHO	H) ₄ CH ₂ OH
M =	182.17 g/mo
CACNO	CO CE O

CAS NO. 69-65-8 **EINECS** 200-711-8 NC CODE 2905 43 00

Meets ACS Specifications	
Insoluble Matter	max. 0.01%
Loss on Drying at 105°C	max. 0.05%
Reducing Sugars	passes test
Residue after Ignition	max. 0.01%
Specific Rotation [α] 25/10	$+23.3^{\circ}$ to $+24.3^{\circ}$
Titrable Acid (meq/g)	max. 0.0008
Trace Impurities (in ppm):	
Heavy Metals (as Pb)	max. 5

Trace Impurities (in ppm):

Substances Reducing KMnO₄

Copper (Cu)

Nickel (Ni)

PRODUCT	PACKING	CONT
NO.		ВОХ
1125.0250	250 g	
1125.1000	1 kg	6

Mannitol 'BAKER'

► HOCH₂(CHOH)₄CH₂OH

M = 182.17 g/mol **CAS NO.** 69-65-8

EINECS 200-711-8 NC CODE 2905 43 00

Assay (dried basis)	96.0-101.5%
Acidity	passes test
Arsenic (As)	max. 1 ppm
Chloride (CI)	max. 0.007%
Identification	passes test
Loss on Drying at 105°C	max. 0.3%
Melting Range	164-169°C
Reducing Sugars	passes test
Specific rotation	+ 137- + 145°
Sulfate (SO ₄)	max. 0.01%

PRODUCT	PACKING	CONT.
NO.		вох
1998.1000	1 kg	6
1998.9050	50 kg	

Mannitol

See for qualities meeting the requirements of cGMP, USP/NF, BP/Ph.Eur, JP or FCC, page 36



Marshall's Reagent

See N-(1-Naphthyl) Ethy



May Grünwald

Powder / HEMATOLOGY/CYTOLOGY/HISTOLOGY



NC CODE 3204 12 00 **R:** 22-36

S: 22-26



Non oxidized Methylene Blue/Eosine Y Stain,
suitable for staining blood smears

PRODUCT	PACKING	CONT.
١0.		вох
8814.0025	25 g Glass	

May

May Grünwald

3855 HEMATOLOGY/CYTOLOGY/HISTOLOGY

1 l = 0.79 kg FLASHPOINT 11 °C NC CODE 3204 12 00 **UN/ID NO.** 1230 ADR/RID 3 FT1

IMDG 3/II

modified for Microscopy (contains Methanol)

PRODUCT CONT. **PACKING** NO. вох 3855.0100 100 ml Glass 3855.0500 500 ml Glass 3855.1000 1 | Glass 3855.2500 2.5 | Glass 3855.9020PE 20 L jerrycan

R: 11-23/24/25-39-39/23/24/25 **S:** 16-36/37/39-45











Mayer's Hematoxyline



Menthol

See for qualities meeting the requirements of cGMP, USP/NF, BP/Ph.Eur, JP or FCC, page 36



levo-Menthol

See for qualities meeting the requirements of cGMP, USP/NF, BP/Ph.Eur, JP or FCC, page 36

4049

2-Mercaptoethanol

'BAKER ULTRAPURE BIOREAGENT'

► HSCH₂CH₂OH

M = 78.13 g/mol

1 | = 1.11 kg FLASHPOINT 74 °C

CAS NO. 60-24-2 **EINECS** 200-464-6

NC CODE 2930 90 70

UN/ID NO. 2966 **ADR/RID** 6.1 T1

IMDG 6.1/II

R: 22-24-34-51/53

S: 26-28A-36/37/39-45-61







For Electrophoresis and other Molecular Biology Applications Assay (by GC) min. 99.0% Appearance passes test DNase Activity none detected Protease Activity none detected RNase Activity none detected Trace Impurities (in ppm):

May-Grünwalds Eosine-Methylene Blue Solution

	,
Copper (Cu)	max. 10
Iron (Fe)	max. 5
Lead (Pb)	max. 5
Zinc (Zn)	max. 5

CONT. **PRODUCT PACKING** NO. вох 4049.0100 100 g 4049.0500 500 g

environment

Mercury
'BAKER ANALYZED' / ACS 9006

► Hg

M = 200.59 g/mol

CAS NO. 7439-97-6

EINECS 231-106-7

NC CODE 2805 40 90 EC NO. 80 001 00 0

UN/ID NO. 2809

ADR/RID 8 C9

IMDG 8/III

R: 23-33-50/53

S: 45-60-61-7





dangerous for the environment

Meets ACS Specifications	
Appearance	passes test
Non-volatile Matter	max. 5 ppm

PRODUCT	PACKING	CONT.
NO.		ВОХ
9006.0500	500 g	
9006.1000	1 kg	



Mercury 1000 μg/ml

nitric acid) / 'BAKER INSTRA-ANALYZED' / Plasma Standard

5768

► Hg

M = 200.59 g/mol

NC CODE 3822 00 00 **EC NO.** 7 004 00 1

R: 20/21/22-33-36/38

S: 26-36/37



Certificate Provided Reporting Actual Lot Analysis

Mercury (Hg) 998-1002 μg/ml

PRODUCT CONT. **PACKING** NO. вох 5768.0100 100 ml

Prepared from the highest purity raw material available, generally greater than 99.999% spectral purity. The content of the solution is confirmed to be accurate to within \pm 0.2%. Typically 1000 μ g/ml.

The certificate of analysis provided reports actual lot analysis. The certificate also lists the trace impurities.

Mercury 1000 μg/ml 0.10% (w/v) / (Matrix: 2% nitric a

/ 'BAKER INSTRA-ANALYZED' / Atomic Absorption Standard

Mercury (Hg)

6934

► Hg

M = 200.59 g/mol

NC CODE 3822 00 00

R: 20/21/22-33-36/38

S: 26-36/37



PRODUCT PACKING CONT. 6934.0100 100 ml 6934.0500 500 ml

Prepared by dissolution of high purity raw materials (min. 99.99% spectral purity). Assays are verified by ICP against standards traceable to NIST. Standard Reference Material numbers (SRM) are printed on each label.

Mercury 1000 μg/ml

BAKER ANALYZED' / At

6817

CONT.

вох

CONT.

► Hg

NC CODE 3822 00 00

R: 20/21/22-33-36/38

M = 200.59 g/mol

S: 26-36/37



Mercury (Hg) 998-1002 $\mu g/ml$

PRODUCT NO.

6817.0100 100 ml 6817.0500 500 ml

Mercury(II)nitrate in nitric acid 0.5 mol/l.

PACKING

PACKING

Mercury 10000 μg/ml

(Matrix: 2% nitric acid) / 'BAKER INSTRA-ANALYZED' / Plasma Standard

5736

► Hq

M = 200.59 g/mol

NC CODE 3822 00 00

EC NO. 7 004 00 1

R: 23/24/25-33-36/38-52/53

S: 26-36/37-45



Certificate Provided Reporting Actual Lot Analysis Mercury (Hg) 9980-10020 $\mu g/ml$

PRODUCT

998-1002 µg/ml

5736.0100 100 ml

Prepared from the highest purity raw material available, generally greater than 99.999% spectral purity. The content of the solution is confirmed to be accurate to within 0.2 %. Typically 10000 µg/ml.

The certificate of analysis provided reports actual lot analysis. The certificate also lists the trace impurities.

Make a clear diagnosis with J.T.Baker Histopathology reagents and stains.

See chapter 4 of this catalogue for more details.

Mercu

Mercury(II) Acetate

1128 'BAKER ANALYZED' / ACS

► Hg(CH₃COO)₂ **M** = 318.68 g/mol

CAS NO. 1600-27-7 **EINECS** 216-491-1

NC CODE 2915 29 00 **EC NO.** 80 004 00 7

UN/ID NO. 1629

ADR/RID 6.1 T5 IMDG 6.1/II

R: 26/27/28-33-50/53

S: 13-28-36-45-60-61





dangerous

very toxic

Meets ACS Specifications	
Assay	min. 98.0%
Chloride (CI)	max. 0.005%
Insoluble Matter	max. 0.01%
Iron (Fe)	max. 0.001%
Mercurous Mercury (as Hg)	max. 0.4%
Nitrate (NO ₃)	max. 0.005%
Other Heavy Metals (as Pb)	max. 0.002%
Residue after Reduction	max. 0.02%
Sulfate (SO ₄)	max. 0.005%

Meets ACS Specifications

Residue after Reduction

Solution in Ethyl Ether

Assay

Iron (Fe)

PRODUCT	PACKING	CONT
NO.		BO)
1128.0500	500 g	

PRODUCT

1131.0100

1131.0500

min. 99.5%

max. 0.002%

max. 0.02%

passes test

PACKING

100 g

500 g

CONT.

BOX

1131

Mercury(II) Chloride

► HgCl₂

M = 271.50 g/mol

CAS NO. 7487-94-7

EINECS 231-299-8

NC CODE 2827 39 80 **EC NO.** 80 010 00 0

UN/ID NO. 1624

ADR/RID 6.1 T5

IMDG 6 1/II

R: 28-34-48/24/25-50/53





dangerous for the

Mercury(II) Iodide, Red

M = 454.40 g/mol CAS NO. 7774-29-0

EINECS 231-873-8

NC CODE 2827 60 00

EC NO. 80 002 00 6

UN/ID NO. 1638

ADR/RID 6.1 T5

IMDG 6.1/II

R: 26/27/28-33-50/53

S: 13-28-45-60-61





Meets ACS Specifications		
Assay	min. 99.0%	
Mercurous Mercury (as Hg)	max. 0.1%	
Solubility in Potassium Iodide Solution	passes test	
Soluble Mercury Salts (as Ha)	max. 0.05%	

PRODUCT NO.	PACKING	CONT.
		вох
1133.0100	100 a	

Questions or suggestions, please contact us at jtbaker.nl@emea.tycohealthcare.com

1134

7252

1136

CONT.

вох

Mercury(II) Nitrate Monohydrate

Titer (mol/l)

► Hg(NO₃)₂.H₂O M = 342.62 g/mol

CAS NO. 10045-94-0 **EINECS** 233-152-3

NC CODE 2834 29 30

EC NO. 80 002 00 6

UN/ID NO. 1625

ADR/RID 6.1 T5

IMDG 6.1/II

R: 26/27/28-33-50/53

S: 13-28-45-60-61





dangerous

Assay	min. 99.0%
Chloride (CI)	max. 0.002%
Iron (Fe)	max. 0.001%
Mercurous Mercury (as Hg)	max. 0.1%
Residue after Ignition	max. 0.01%
Sulfate (SO ₄)	max. 0.002%

0.0495-0.0505

PRODUCT

7252.1000

PACKING

11

PRODUCT PACKING CONT. вох 1134.0100

Mercury (II) Nitrate

► Hg(NO₃)₂

M = 342.62 g/mol

1 l = 1.00 kg

NC CODE 2834 29 30

UN/ID NO. 2024 **ADR/RID** 6.1 T4

IMDG 6.1/II

R: 23/24/25-33-52/53

S: 26-36/37-45



Mercury(II) Oxide, Red

► HgO **M** = 216.59 g/mol **CAS NO.** 21908-53-2

EINECS 244-654-7

NC CODE 2825 90 50

EC NO. 80 002 00 6

UN/ID NO. 1641 **ADR/RID** 6.1 T5

IMDG 6.1/II

R: 26/27/28-33-50/53

S: 13-28-45-60-61





Meets ACS Specifications	
Assay	min. 99.0%
Chloride (CI)	max. 0.025%
Insoluble in Dilute HCI	max. 0.03%
Iron (Fe)	max. 0.005%
Nitrogen Compounds (as N)	max. 0.005%
Residue after Reduction	max. 0.025%
Sulfate (SO ₄)	max. 0.015%

PACKING CONT. **PRODUCT** 1136.0100 100 g

Only the peaks that matter are analysed with J.T.Baker LC/MS solvents and blends!

Refer to LC/MS section of this catalogue for more details.



► H

Mercu

Mercury(II) Sulfate 'BAKER ANALYZED' / ACS

1140

lgSO₄	
M =	296.65 g/mol
CAS NO.	7783-35-9
EINECS	231-992-5
NC CODE	2833 29 70
EC NO.	80 002 00 6
UN/ID NO.	1645
ADR/RID	6.1 T5

IMDG 6.1/II **R:** 26/27/28-33-50/53

environment

S: 13-28-45-60-61 dangerous for the very toxic

Meets ACS Specifications	
Assay	min. 99.0%
Chloride (CI)	max. 0.003%
Iron (Fe)	max. 0.003%
Mercurous Mercury (as Hg)	max. 0.15%
Nitrate (NO ₃)	max. 0.005%
Residue after Reduction	max. 0.02%
Suitable for COD det. acc. NEN (max. 0.8 g	
Hg(II)SO ₄ /20 mI)	passes test

PRODUCT	PACKING	CONT. BOX
1140.0100	100 g	6
1140.0500	500 g	6

Mercury(II) Sulfate

7464 te sulfuric acid / 'BAKER ANALYZED' / for the COD determination according DIN 38409-H43

EINECS 231-992-5 NC CODE 2833 29 70 UN/ID NO. 3289 ADR/RID 6.1 TC3 IMDG 6.1/II

R: 26/27/28-33-36/38-51/53

S: 28-29-36/37-45



PRODUCT PACKING CONT. NO. BOX 7464.1000 11 6 7464.2500 2.5 l

Volumetric Solution, ready for use.

PACKING

11

2.5 |

Volumetric Solution, ready for use

CONT.

вох

6

PRODUCT

7135.1000

7135.2500

NO.

passes test

0.0198-0.0202

Mercury(II) Sulfate

80 g/l in 0.020 mol/l $K_2Cr_2O_7$ solution in sulfuric acid. / 'BAKER ANALYZED' / For COD determination according DIN 38409-H41.

Identity (Hg+)

Molarity (K₂Cr₂O₇) (M)

11 = 1.19 kg **CAS NO.** 7783-35-9 **EINECS** 231-992-5 NC CODE 2833 29 70 UN/ID NO. 3289 ADR/RID 6.1 TC3 IMDG 6.1/II

R: 26/27/28-33-35-43-45-46-51/53

S: 26-28A-29-36/37/39-45-53





dangerous for the environment

very toxic

Molarity (M)

Mercury(II) Sulfate 80 g/l in 0.04 mol/l K₂Cr₂O₂ solution in sulfuric acid / 'BAKER ANALYZED' / for the COD determiniation according AFNOR NF T 90-101

1 l = 1.19 kg **CAS NO.** 7783-35-9 NC CODE 2833 29 70 UN/ID NO. 3289 ADR/RID 6.1 TC3 IMDG 6.1/II

7642

R: 26/27/28-33-35-43-46-49-51/53

S: 26-28-29-36/37/39-45-53-57







0.0396 - 0.0404

PRODUCT PACKING BOX 7642 1000 11 6 7642.2500 2.5 |

 $Mercury (II) Sulfate \ 80 \ g/I \ in \ Sulfuric \ acid \ / \ Potassium$ dichromate 0.04 mol/l.



MES Monohydrate, Free Acid

See 2-(N-Morpholino)-ethanesulfonic Acid Monohydrate

PACKING

PRODUCT

9077.1000

Metaphosphoric Acid





Methanoic Acid

See Formaldehyde solution

See Formic Acid



CONT.

Methanol

'BAKER ULTRA RESI-ANALYZED' / Absolute, suitable for Purge and Trap analysis

CH₂OH	

M = 32.04 g/mol

11 = 0.79 kg

FLASHPOINT 11 °C

CAS NO. 67-56-1

EINECS 200-659-6 NC CODE 2905 11 00

EC NO. 603 001 00 0

UN/ID NO. 1230

ADR/RID 3 FT1

IMDG 3/II

R: 11-23/25

S: 16-24-45-7





Assay (by GC) (corrected for water)	min. 99.9%	
Residue after Evaporation	max. 1 ppm	
Titrable Acid (meq/g)	max. 0.0003	
Titrable Base (meq/g)	max. 0.0001	
Water (H ₂ O)	max. 0.08%	
Volatile Organic Trace Analysis. Gas		
Chromatography with Purge and Trap		
Concentration. EPA Contract Required Quantitation		
Limit-CRQL.:		
Electroconductivity Detection, Below CRQL	passes test	
Photoionization Detection, Below CROI	passes test	

Methanol

AKER ULTRA RESI-ANALYZED' / for Organic Residue Analys

9263

BOX

CHI3OTT	
M =	32.04 g/mol
11=	0.79 kg
FLASHPOINT	11 °C
CAS NO.	67-56-1
EINECS	200-659-6
NC CODE	2905 11 00
EC NO.	603 001 00 0
UN/ID NO.	1230

► CH OH

ADR/RID 3 FT1 IMDG 3/II **R:** 11-23/25

> highly flammable



Assay (by GC) (corrected for water)	mm. 99.8%
Acetone	max. 0.001%
Color (APHA)	max. 10
Residue after Evaporation	max. 1 ppm
Substances Darkened by H ₂ SO ₄	passes test
Titrable Acid (meq/g)	max.0.0003
Titrable Base (meq/g)	max.0.0001
Water (H ₂ O)	max. 0.08%
ECD Sensitive Impurities (as Heptach	lor Epoxide):
Single Impurity Peak (pg/ml)	max. 10
FID-Sensitive Impurities (as 2-Octano	I):
Single Impurity Peak (ng/ml)	max. 5

Assay (by GC) (corrected for water)

Mallinckrodt Baker's cGMP Manufactured
Chemicals for the Biopharmaceutical industry
are a necessity for uncomplicated scale-up.

min 99.8%

PRODUCT

9263.1000

9263.2500

PACKING

11

2.5 l

See chapter 6 of this catalogue.



Metha

Methanol

9822 BAKER ANALYZED LC-MS Reagent

M = 32.04 g/mol 1 I = 0.79 kg FLASHPOINT 11 °C CAS NO. 67-56-1 EINECS 200-659-6 NC CODE 2905 11 00 EC NO. 603 001 00 0 UN/ID NO. 1230 ADR/RID 3 FT1

IMDG 3/II

R: 11-23/24/25-39/23/24/25

S: 16-36/37-45-7





Certificate Provided Reporting	g Actual Lot Analysis
Assay (by GC)	min. 99.8%
Residue after Evaporation	max. 1 ppn
Water (H ₂ O)	max. 0.02%
LC-Gradient-Diode Array Det	ection (a.u.), test
solution is modified with 0.1%	(v/v) formic acid:
at 235 nm	max. 0.005
at 254 nm	max. 0.00
LC-MS Gradient Suitability Te	est (TIC, 100 to 2000
m/z), test solution is modified	l with 0.1% (v/v)
formic acid:	
Positive ESI-MS Sensitive Impuri	ties (as
Reserpine)	max. 50 ng/m
Product Information (not spec	cifications):
Density (g/ml) at 20°C	0.79
Trace Impurities (in ppb):	
Aluminium (AI)	max. 50
Calcium (Ca)	max. 50
Iron (Fe)	max. 50
Magnesium (Mg)	max. 50
Sodium (Na)	max. 50
Ultraviolet Absorbance (1.00-	cm path vs water):
at 225 nm	max. 0.20
at 254 nm	max. 0.02
at 280 nm	max. 0.0

PACKING	CONT
	BO
1 Glass	(

Element concentrations are at time of lot release

Methanol

'BAKER HPLC ANALYZED' / HPLC Gradient Grade

M = 32.04 g/mol

1 I = 0.79 kg

FLASHPOINT 11 °C

CAS NO. 67-56-1

EINECS 200-659-6

NC CODE 2905 11 00

EC NO. 603 001 00 0

UN/ID NO. 1230 ADR/RID 3 FT1 IMDG 3/II

R: 11-23/24/25-39/23/24/25 **S:** 16-36/37-45-7

highly



Acetone	max. 0.001%	
Residue after Evaporation (in ppm)	max. 2	
Titrable Acid (meq/g)	max. 0.0003	
Titrable Base (meq/g)	max. 0.0001	
Water (H ₂ O)	max. 0.02%	
Fluorescence Trace Impurities (as quinine base),		
ppb:		
Measured at 450 nm	max. 0.3	
Measured at Emission Maximum for		
Solvent Impurities	max. 1.0	
Gradient Elution Test Ultraviolet Absorbance (a.u.):		
at 235 nm	max. 0.005	
at 254 nm	max. 0.001	
Ultraviolet Absorbance (1.00-cm path vs water):		
at 225 nm	max. 0.17	

min. 99.8%

max. 0.01

max. 0.01

PRODUCT

8402.1000

8402.2500

Assay (by GC) (corrected for water)

8402.5000 5 | EcoTainer

Filtered through a 0.2 micron filter.

PACKING

11

Packaged under Nitrogen.

2.5 |

J.T.Baker SPE and Speedisk products for separation specialists.

at 280 nm

at 350 nm

UV Cut-off, nm

Refer to the special Chromatography section in the rear of this catalogue.

Methanol

'BAKER HPLC ANALYZED' / HPLC Isocratic Grade

8404

► CH₃OH

M = 32.04 g/mol

1 l = 0.79 kg

FLASHPOINT 11 °C

CAS NO. 67-56-1 **EINECS** 200-659-6

NC CODE 2905 11 00

EC NO. 603 001 00 0

UN/ID NO. 1230

ADR/RID 3 FT1

IMDG 3/II

R: 11-23/24/25-39/23/24/25

S: 16-36/37-45-7





Assay (by GC)	min. 99.8%	
Color (APHA)	max. 10	
Residue after Evaporation	max. 5 ppm	
Titrable Acid (meq/g)	max. 0.0003	
Titrable Base (meq/g)	max. 0.0002	
Water (H ₂ O)	max. 0.05%	
Ultraviolet Absorbance (1.00-cm path vs water):		
at 225 nm	max. 0.30	
at 240 nm	max. 0.10	

max. 0.01

min. 99.8%

passes test

max. 10

PRODUCT	PACKING	CONT.
NO.		вох
8404.1000	11	6
8404.2500	2.5	4
8404.5000	5 EcoTainer	4
8404.9010RC	10 Returnable Container	
8404.9030RC	30 l Returnable Container	
8404.9200RC	200 l Returnable Container	

Methanol

'BAKER ANALYZED' / Ultraviolet Spectrophotometry / ACS

at 265 nm

8046

9097

> (CH₃OH
---------------	-------

M = 32.04 g/mol

1 I = 0.79 kg

FLASHPOINT 11 °C

CAS NO. 67-56-1

EINECS 200-659-6

NC CODE 2905 11 00

EC NO. 603 001 00 0

UN/ID NO. 1230

ADR/RID 3 FT1

IMDG 3/II

R: 11-23/24/25-39/23/24/25

S: 16-36/37-45-7





Exceeds ACS Specifications		
Assay (by GC)	min. 99.8%	
Appearance	clear	
Carbonyl Compounds	passes test	
Color (APHA)	max. 10	
Residue after Evaporation	max. 5 ppm	
Solubility in Water	passes test	
Substances Darkened by H ₂ SO ₄	passes test	
Substances Reducing KMnO ₄	passes test	
Titrable Acid (meq/g)	max. 0.0003	
Titrable Base (meq/g)	max. 0.0002	
Water (H ₂ O)	max. 0.10%	
Ultraviolet Absorbance (1.00-cm path vs water):		

Oftraviolet Absorbance (1.00-cm path vs water):		
at 205 nm	max. 1.00	
at 210 nm	max. 0.80	
at 220 nm	max. 0.40	
at 230 nm	max. 0.20	
at 240 nm	max. 0.10	
at 260 nm	max. 0.04	
at 280-400 nm	max. 0.01	

PRODUCT NO.	PACKING	CONT. BOX
8046.1000	11	6
8046.2500	2.5	4

Methanol

BakerDRY / Low Water Solvent / ACS

► CH₃OH

M = 32.04 g/mol

1 l = 0.79 kg

FLASHPOINT 11 °C **CAS NO.** 67-56-1

EINECS 200-659-6

NC CODE 2905 11 00

EC NO. 603 001 00 0

UN/ID NO. 1230

ADR/RID 3 FT1

IMDG 3/II

R: 11-23/24/25-39/23/24/25

S: 16-36/37-45-7





Appearance
Color (APHA)
Residue after Evaporation
Solubility in Water

Assay (CH₃OH)(by GC,corrected for water)

Meets ACS Specifications

max. 1 ppm passes test Substances Darkened by H₂SO₄ passes test Substances Reducing Permanganate passes test Titrable Acid (µeq/g) max. 0.3 max. 0.2 Titrable base (μ eq/g) Water (by KF, coulometric) max. 30 ppm

Carbonyl Compounds:	
Acetaldehyde	max. 0.001%
Acetone	max. 0.001%
Formaldehyde (HCHO)	max. 0.001%

PRODUCT	PACKING	CONT
NO.		BO)
9097.1000	11	



Metha

Methanol

'BAKER BIO-ANALYZED' / For Biotechnology Applications

M = 32.04 g/mol 11 = 0.79 kg FLASHPOINT 11 °C CAS NO. 67-56-1 EINECS 200-659-6 NC CODE 2905 11 00

EC NO. 603 001 00 0
UN/ID NO. 1230
ADR/RID 3 FT1
IMDG 3/II

R: 11-23/24/25-39/23/24/25

S: 16-36/37-45-7





Assay (CH ₃ OH)(by GC,corrected for water)	min. 99.8%
Acetone	max. 0.001%
Residue after Evaporation	max. 2.0 ppm
Titrable Acid (meq/g)	max. 0.0003
Titrable Base (meq/g)	max. 0.0001
Water (by KF, coulometric)	max. 0.005%
Gradient Elution Test (a.u.):	
at 254 nm	max. 0.005
Physical Data (not specifications):	
Density (g/ml) at 20°C	0.791

 Ultraviolet Absorbance (1.00-cm path vs water):

 at 225 nm
 max. 0.15

 at 254-400 nm
 max. 0.01

 UV Cut-off, nm
 max. 205

Exceeds ACS Specifications

PRODUCT NO.	PACKING	CON B(
9091.4000GL	4 I Glass	

Methanol

max. 0.01% H₂O / 'BAKER ANALYZED' / Dried / ACS

M = 32.04 g/mol 1 I = 0.79 kg FLASHPOINT 11 °C CAS NO. 67-56-1 EINECS 200-659-6 NC CODE 2905 11.00

9049

NC CODE 2905 11 00 EC NO. 603 001 00 0 UN/ID NO. 1230

ADR/RID 3 FT1 IMDG 3/II

R: 11-23/24/25-39/23/24/25

S: 16-36/37-45-7





Exceeds ACS Specifications	
Assay	min. 99.8%
Appearance	clear
Carbonyl Compounds	passes test
Color (APHA)	max. 10
Residue after Evaporation	max. 0.001%
Solubility in Water	passes test
Substances Darkened by H ₂ SO ₄	passes test
Substances Reducing KMnO ₄	passes test
Titrable Acid (meq/g)	max. 0.0003
Titrable Base (meq/g)	max. 0.0002
Water (H ₂ O)	max. 0.01%
Trace Impurities (in ppm):	
Aluminium (AI)	max. 0.5
Barium (Ba)	max. 0.1
Boron (B)	max. 0.02
Cadmium (Cd)	max. 0.05
Calcium (Ca)	max. 0.5
Chromium (Cr)	max. 0.02
Cobalt (Co)	max. 0.02
Copper (Cu)	max. 0.02
Iron (Fe)	max. 0.1
Lead (Pb)	max. 0.1
Magnesium (Mg)	max. 0.1
Manganese (Mn)	max. 0.02
Nickel (Ni)	max. 0.02
Tin (Sn)	max. 0.1
. ,	

RODUCI	FACKING	CONT.
NO.		вох
9049.1000	11	6
9049.1000S	1 EcoTainer	
9049.2500	2.5	4
coTainer, the	metal solvent can for mo	re safety in the lab.

PRODUCT

PACKING

Prepare yourself for the future with J.T. Baker molecular biology and biotechnology products.

max. 0.1

See chapter 3 of this catalogue for more details.

Zinc (Zn)

2685

Methanol

8045 'BAKER ANALYZED' / ACS

► CH₃OH

M = 32.04 g/mol

1 l = 0.79 kg

FLASHPOINT 11 $^{\circ}\mathrm{C}$

CAS NO. 67-56-1

EINECS 200-659-6 NC CODE 2905 11 00

EC NO. 603 001 00 0

UN/ID NO. 1230

ADR/RID 3 FT1 IMDG 3/II

R: 11-23/24/25-39/23/24/25

S: 16-36/37-45-7





Assay	min. 99.8%
Appearance	clear
Carbonyl Compounds	passes test
Color (APHA)	max. 10
Residue after Evaporation	max. 0.001%
Resistivity (megaohm-cm)	min. 0.1
Solubility in Water	passes test
Substances Darkened by H ₂ SO ₄	passes test
Substances Reducing KMnO ₄	passes test
Titrable Acid (meq/g)	max. 0.0003
Titrable Base (meq/g)	max. 0.0002
Water (H ₂ O)	max. 0.05%
Trace Impurities (in ppm):	
Aluminium (AI)	max. 0.5
Barium (Ba)	max. 0.1
Boron (B)	max. 0.02
Cadmium (Cd)	max. 0.05
Calcium (Ca)	max. 0.5
Chromium (Cr)	max. 0.02
Cobalt (Co)	max. 0.02
Copper (Cu)	max. 0.02
Iron (Fe)	max. 0.1
Lead (Pb)	max. 0.1
Magnesium (Mg)	max. 0.1
Manganese (Mn)	max. 0.02
Nickel (Ni)	max. 0.02
Tin (Sn)	max. 0.1
Zinc (Zn)	max. 0.1

PRODUCT NO.	PACKING CON	
8045.1000	11	6
8045.2500	2.5	4
8045.2500PE	2.5 I HDPE	4
8045.5000	5 EcoTainer	4
8045.9010RC	10 l Returnable Container	
8045.9025	25	4
8045.9030RC	30 l Returnable Container	
8045.9200	200	
EcoTainer, the	metal solvent can for more safety in the la	b.
For safe handl	ing of 25 l tin cans, see Self-closing tap.	

Methanol

'BAKER ANALYZED' / ACS, PhEur

► CH₃OH

M = 32.04 g/mol

11= 0.79 kg

FLASHPOINT 11 °C

CAS NO. 67-56-1 **EINECS** 200-659-6

NC CODE 2905 11 00

EC NO. 603 001 00 0

UN/ID NO. 1230

ADR/RID 3 FT1

IMDG 3/II

R: 11-23/24/25-39/23/24/25

S: 16-36/37-45-7





Exceeds ACS Specifications. Meets	Ph.Eur.
Reagent Specifications	
A.C.S. Specification(s):	
Appearance	clear
Assay	min. 99.8%
Carbonyl Compounds	passes test
Color (APHA)	max. 10
Residue after Evaporation	max. 0.001%
Resistivity (megaohm-cm)	min. 0.1
Solubility in Water	passes test
Substances Darkened by H ₂ SO ₄	passes test
Substances Reducing KMnO ₄	passes test
Titrable Acid (meq/g)	max. 0.0003
Titrable Base (meq/g)	max. 0.0002
Water (H ₂ O)	max. 0.05%
Ph. Eur. Reagent Specification(s):	
Acetone, Aldehydes (as Acetone)	max. 0.001%
Appearance	passes test

PRODUCT NO.	PACKING	CONT. BOX
2685.9200ST	200 l Steel Drum	

J.T.Baker: over 100 years of experience.

See chapter 1 of this catalogue.



Methanol 'BAKER'

8047

► CH₃OH

M = 32.04 g/mol **1 l** = 0.79 kg

FLASHPOINT 11 °C

CAS NO. 67-56-1

EINECS 200-659-6

NC CODE 2905 11 00

EC NO. 603 001 00 0

UN/ID NO. 1230 ADR/RID 3 FT1

IMDG 3/II

R: 11-23/24/25-39/23/24/25

S: 16-36/37-45-7





Assay	min. 99.5%
Acetone, Aldehydes (as Acetone)	max. 0.002%
Boiling Range (initial to dry point)	max. 2°C
Residue after Evaporation	max. 0.002%
Titrable Acid (meq/g)	max. 0.0003
Titrable Base (meq/g)	max. 0.0002
Water (H ₂ O)	max. 0.1%

PRODUCT	PACKING	CONT.
NO.		вох
8047.1000	11	
8047.2500	2.5	4
8047.5000	5 EcoTainer	
8047.9025	25 l	4
8047.9200	200	
FeeTainer th	a motal coluent can for mor	o cafety in the lab

For safe handling of 25 l tin cans, see Self-closing tap.

3400

Methanol

99.5% / HISTO GRADE

► CH₃OH

M = 32.04 g/mol

1 l = 0.79 kg

FLASHPOINT 11 °C

CAS NO. 67-56-1

EINECS 200-659-6

NC CODE 2905 11 00 EC NO. 603 001 00 0

UN/ID NO. 1230

ADR/RID 3 FT1

IMDG 3/II

R: 11-23/24/25-39/23/24/25

S: 16-36/37-45-7





PRODUCT PACKING CONT. NO. BOX 3400.0500PE 500 ml HDPE 3400.5000 5 | Jerrycan 3400.9010 10 l Jerrycan 3400.9025 25 l Jerrycan

Histo-Grade implicates that this reagent is specially tested and therefore solely intended for use in histo-pathology applications. This reagent is of an analytical quality.

Methanol

See for qualities meeting the requirements of cGMP, USP/NF, BP/Ph.Eur, JP or FCC, page 36



Methanol MOS, VLSI, ULSI Grade

ee for detailed information section Microelectronic materials page 32 or the website www.jtbaker.com/micro



Methenamine

See for qualities meeting the requirements of cGMP, USP/NF, BP/Ph.Eur, JP or FCC, page 36



L-Methionine

See for qualities meeting the requirements of cGMP, USP/NF, BP/Ph.Eur, JP or FCC, page 36



Methotrexate

e for qualities meeting the requirements of cGMP, USP/NF, BP/Ph.Eur, JP or FCC, page 36



Methoxybenzene



Major fields of application for reagents:

Biotechnology

We offer a complete range:

- Synthesis reagents
- DNA/RNA purification and analysis

Synthesis reagents available for the following DNA/RNA Synthesisers:

- ABI series 390Z, 392, 394, 3900 and 3948
- Milligen Biosearch System 8800
- Expedite 8905 and 8909

We offer a complete range of synthesis solvents and reagents:

- Acetonitrile
 (<30 ppm bottles; <10 ppm CYCLE-TAINER System)
- Dichloromethane (< 30 ppm)
- Deblock
- Activator: besides 1H-Tetrazole we have available 5-Ethylthio-1H-tetrazole (ETT) and 4,5-Dicyanoimidazole (DCI)
- Capping A
- Capping B
- Oxidising

Reagents for DNA/RNA purification and analysis

- HPLC Solvents and reagents
- BAKERBOND spe
- BAKERBOND Speedisk

The critical flushing solvent: Acetonitrile

Product Name	Description	Product Code	Content
Acetonitrile	Acetonitrile low water in CYCLE-TAINER (< 10 ppm)	8134	20, 50, 200, 1400 L
Acetonitrile	Acetonitrile low water in bottles (< 30 ppm)	8144	2,5 and 4 L



Metho

2-Methoxyethanol 'BAKER ANALYZED'

8084

► CH₃OCH₂CH₂OH

M = 76.10 g/mol

1 l = 0.96 kg

FLASHPOINT 39 °C

CAS NO. 109-86-4 **EINECS** 203-713-7

NC CODE 2909 42 00

EC NO. 603 011 00 4

UN/ID NO. 1188

ADR/RID 3F1 IMDG 3/III

R: 10-20/21/22-60-61

S: 45-53



Assay (by GC)	min. 98%
Boiling Range	123-125°C
Color (APHA)	max. 10
Density (g/ml) at 20°C	0.964-0.966
Peroxide	passes test
Residue after Evaporation	max. 0.002%
Water (H ₂ O)	max. 0.2%
Trace Impurities (in ppm):	
Aluminium (AI)	max. 0.5
Barium (Ba)	max. 0.1
Boron (B)	max. 0.02
Cadmium (Cd)	max. 0.05
Calcium (Ca)	max. 0.5
Chromium (Cr)	max. 0.02
Cobalt (Co)	max. 0.02
Copper (Cu)	max. 0.02
Iron (Fe)	max. 0.1
Lead (Pb)	max. 0.2
Magnesium (Mg)	max. 0.1
Manganese (Mn)	max. 0.02
Nickel (Ni)	max. 0.02
Tin (Sn)	max. 0.1
Zinc (Zn)	max. 1.0

PRODUCT	PACKING	CONT
NO.		B0)
8084.1000	11	6
8084.9025	25 l	
For safe han	dling of 25 l tin cans, see Self-closing	tap.



Methyl Alcohol

See Methanol



Methyl Alcohol

e for qualities meeting the requirements of cGMP, USP/NF, BP/Ph.Eur, JP or FCC, page 36

Assay

Meets ACS Specifications

Suitability for Phopshate Determination

Residue after Ignition

2794

p-(Methylamino)phenol Sulfate

'BAKER ANALYZED' / ACS

► (CH₃NHC₆H₄OH)₂.H₂SO₄ **M** = 344.39 g/mol

> **CAS NO.** 55-55-0 **EINECS** 200-237-1

NC CODE 2922 29 00

EC NO. 650 031 00 4

UN/ID NO. 3077

ADR/RID 9 M7

IMDG 9/III

R: 22-43-48/22-50/53

S: 36/37-46-60-61



environment



harmful

Methyl Benzoate

8123 ► C₆H₅COOCH₃

M = 136.15 g/mol

1 l = 1.09 kg

FLASHPOINT 83 °C **CAS NO.** 93-58-3

EINECS 202-259-7

NC CODE 2916 31 00

R: 22

S: 23-24/25



min. 99%
passes test
197-200°C

PRODUCT	PACKING	CONT.
NO.		вох
8123.1000	11	

PRODUCT

2794.0500

99.0-101.5%

max. 0.1%

passes test

PACKING

500 g

CONT.

CONT.

вох

2-Methylbutane

8113

PRODUCT

8113.1000

8113.9025

NO.

min. 99%

27.5-28.5°C

max. 0.005%

PACKING

For safe handling of 25 l tin cans, see Self-closing tap.

11

25 l

► CH₃CH₂CH(CH₃)₂

M = 72.15 g/mol

1 l = 0.62 kg

FLASHPOINT $-55\,^{\circ}\text{C}$

CAS NO. 78-78-4

EINECS 201-142-8

NC CODE 2901 10 00

FC NO. 601 006 00 1

UN/ID NO. 1265

ADR/RID 3 F1

IMDG 3/I

R: 12-51/53-65-66-67

S: 16-29-33-61-62-9







Assay (by GC)

Residue after Evaporation

Boiling Point

3-Methyl-1-butanol

See Isoamyl Alcohol



Methyl tert-Butyl Ether

'BAKER ULTRA RESI-ANALYZED' / for Organic Residue Analysis

9043

► (CH₃)COCH₃ **M** = 88.15 g/mol **11 =** 0.74 kg

FLASHPOINT $-28\,^{\circ}\text{C}$

CAS NO. 1634-04-4

EINECS 216-653-1

NC CODE 2909 19 00

UN/ID NO. 2398

ADR/RID 3 F1

IMDG 3/II

R: 12

S: 16-29-33-9



Assay (corrected for H₂O) min. 99.0% max. 1 ppm Peroxide (as H₂O₂) (in ppm) Residue after Evaporation (in ppm) max. 2 ppm Water (H₂O) max. 0.05% ECD Sensitive Impurities (as Heptachlor Epoxide):

Single Impurities (pg/ml)

FID-Sensitive Impurities (as 2-Octanol):

Single Impurity Peak (ng/ml) max. 10 CONT. вох

NO. 9043.1000 9043.4000 4 | Glass

PACKING

PRODUCT

Methyl tert-Butyl Ether

'BAKER HPLC ANALYZED' / for use in High Performance Liquid Chromatography

Assav

ppb:

9042

CONT.

вох

► (CH₃)₃COCH₃

M = 88.15 g/mol

1 I = 0.74 kg

FLASHPOINT −28 °C

CAS NO. 1634-04-4

EINECS 216-653-1

NC CODE 2909 19 00

UN/ID NO. 2398

ADR/RID 3 F1

IMDG 3/II

R: 12

S: 16-29-33-43A-9



Water (H₂O) Product Information (not specifications):

Density (g/ml) at 20°C

UV Cut-off, nm

Residue after Evaporation

Measured at 450 nm

Peroxide (as H₂O₂)

Measured at Emission Maximum

Maximum Limits of Impurities:

Ultraviolet Absorbance (1.00-cm path vs water): at 254 nm max. 0.1 at 280 nm max. 0.02 at 350 nm max. 0.01

Fluorescence Trace Impurities (as quinine base), NO. 9042.1000 value on certificate 9042,2500 9042.5000 value on certificate

min. 99.0%

PRODUCT

EcoTainer, the metal solvent can for more safety in the lab. 1 ppm 3 ppm 0.05%

0.740

max. 215

Packaged under Nitrogen.

PACKING

11

2.5 [

5 | EcoTainer

Methyl Cyanide



Methy

Methylcyclohexane

050 'BAKER ANALYZED

► CH₃CH(CH₂)₄CH₂

M = 98.19 g/mol

1 l = 0.77 kg

FLASHPOINT -4 °C CAS NO. 108-87-2

EINECS 203-624-3

NC CODE 2902 19 80

EC NO. 601 018 00 7

UN/ID NO. 2296

ADR/RID 3 F1

IMDG 3/II

R: 11-38-51/53-65-67

S: 16-33-61-62-9







7

Assay (by GC)

Identification (by IR)

Boiling Point (typical)

Residue after Evaporation

Density (g/ml) at 25°C (typical)

Product Information (not specifications):

4-Methyl-1,3-dioxolan-2-one

See Propylene Carbonate

4031

N,N'-Methylenbisacrylamide

BAKER ULTRAPURE BIOREAGENT

► CH₂(NHCOCH:CH₂)₂ **M** = 154.20 g/mol **CAS NO.** 110-26-9

EINECS 203-750-9

NC CODE 2924 19 00

R: 22



Crosslinking Agent for Polyacrylamide Gels min. 99.5% Assay Acrylic acid (CH₂:CHCOOH) max. 0.01% Conductivity of 2% Solution (µmho) max. 10 DNAase none detected pH of 1% Solution at 25°C min. 5.3 Protease none detected RNAase none detected Absorbance of a 1% Solution (1-cm path vs water):

nax. 10 4031.0100 etected nin. 5.3 etected ete

PRODUCT

4031.0025

NO.

PACKING

25 g Glass

100 g

PRODUCT

8050.0500

min. 99%

101.0°C

0.765

passes test

max. 0.001%

PACKING

379

Methylene Blue

1329 'BAKER ANALYZED

► C₁₆H₁₈CIN₃S.nH₂O

M = 319.86 g/mol **CAS NO.** 61-73-4

EINECS 200-515-2

NC CODE 2934 30 90

R: 22

S: 22-24/25





PRODUCT	PACKING	CONT.
NO.		вох
1329.0025	25 g Glass	6
1329.0100GL	100 g Glass	

CONT.

CONT.

вох

вох



Methylenechloride

See Dichloromethane

Certificates of Analysis are available at www.jtbaker.com/europe

Methyl Ethyl Ketone

'BAKER HPLC ANALYZED' / for use in High Performance Liquid Chromatography

9214

► CH₃CH₂COCH₃

M = 72.11 g/mol

1 l = 0.80 kg

 $\textbf{FLASHPOINT} \quad -\, 1\,\,^{\circ}\text{C}$

CAS NO. 78-93-3

EINECS 201-159-0

NC CODE 2914 12 00

EC NO. 606 002 00 3

UN/ID NO. 1193

ADR/RID 3 F1

IMDG 3/II

R: 11-36-66-67

S: 16-9





Assay (by GC)	min. 99.5%
Residue after Evaporation (in ppm)	max. 3
Water (H ₂ O)	max. 0.03%

Fluorescence Trace Impurities (as quinine base), ppb:

Measured at 450 nm information only Measured at Emission Maximum for

Solvent Impurities information only

Physical Data (not specifications):

Density (g/ml) at 20°C 0.805

Ultraviolet Absorbance (1.00-cm path vs water):		
at 340 nm	max. 0.07	
at 350 nm	max. 0.01	
at 362 nm	max. 0.01	
UV Cut-off, nm	max. 330	

PRODUCT	PACKING	CONT.
NO.		вох
9214.2500	2.5	

Filtered through a 0.2 micron filter. Packaged under Nitrogen.

Methyl Ethyl Ketone

8052

► CH₃CH₂COCH₃

M = 72.11 g/mol

1 l = 0.80 kg

FLASHPOINT −1°C

CAS NO. 78-93-3

EINECS 201-159-0 NC CODE 2914 12 00

EC NO. 606 002 00 3

UN/ID NO. 1193

ADR/RID 3 F1 IMDG 3/II

R: 11-36-66-67

S: 16-9



Assay (by GC)	min. 99%
Acidity (as CH ₃ COOH)	max. 0.003%
Boiling Range	max. 1.5°C
Density (g/ml) at 25°C	0.797-0801
Recorded Boiling Point	79.6°C
Residue after Evaporation	max. 0.001%
Water (H ₂ O)	max. 0.2%
Trace Impurities (in ppm):	

Trace	Impurities	(in ppm):
Alumin	ium (AI)	

Zinc (Zn)

Trace Impurities (in ppm):		
Aluminium (AI)	max. 0.5	
Barium (Ba)	max. 0.1	
Boron (B)	max. 0.02	
Cadmium (Cd)	max. 0.05	
Calcium (Ca)	max. 0.5	
Chromium (Cr)	max. 0.02	
Cobalt (Co)	max. 0.02	
Copper (Cu)	max. 0.02	
Iron (Fe)	max. 0.1	
Lead (Pb)	max. 0.1	
Magnesium (Mg)	max. 0.1	
Manganese (Mn)	max. 0.02	
Nickel (Ni)	max. 0.02	
Tin (Sn)	max. 0.1	

PRODUCT	PACKING	CONT.
NO.		вох
8052.1000	11	6
8052.2500	2.5	4
8052.9025	25	
8052.9200	200	
For cafe han	dline of 25 Ltip cape, con 5	alf clasing tan

Methyl Ethyl Ketone 8053

max. 0.1

► CH₃CH₂COCH₃

M = 72.11 g/mol

11 = 0.80 kg

FLASHPOINT −1°C

CAS NO. 78-93-3

EINECS 201-159-0

NC CODE 2914 12 00

EC NO. 606 002 00 3

UN/ID NO. 1193

ADR/RID 3 F1

IMDG 3/II **R:** 11-36-66-67

S: 16-9





Acidity (as CH₃COOH) max. 0.005% Boiling Range max. 2°C Density (g/ml) at 25°C 0.797-0.801 Residue after Evaporation max. 0.005%

PRODUCT	PACKING	CONT.
NO.		вох
8053.1000	11	6
8053.2500	2.5	
8053.9025	25 l	
8053.9200	200 l	

For safe handling of 25 l tin cans, see Self-closing tap.

Methyl Ethyl Ketone MOS, VLSI Grade

See for detailed information section Microelectronic materials page 32 or the website www.jtbaker.com/micro





Methy

Methyl Isobutyl Ketone

ULTREX Ultrapure Reagent / for Trace Metal Analysis

► CH₃COCH₂CH(CH₃)₂

M = 100.16 g/mol

1 I = 0.80 kg

FLASHPOINT 17 °C

CAS NO. 108-10-1

EINECS 203-550-1

NC CODE 2914 13 00

EC NO. 606 004 00 4

UN/ID NO. 1245

ADR/RID 3 F1

IMDG 3/II

R: 11-20-36/37-66

S: 16-29-9





PRODUCT Certificate Provided Reports Actual Lot Analysis Assay (by GC) (corrected for water) 99.0-100.0% Density (g/ml) at 25°C act. value reported 4855.1000 Residue after Evaporation act. value reported Titrable Acid (µeq/g) act, value reported Water (H₂O) act. value reported Metallic Impurities (in ppb)(ng/g): Bismuth (Bi) Cadmium (Cd) max. 1 Chromium (Cr) Cobalt (Co) Copper (Cu) max. 1 Gallium (Ga) max. 1 Iron (Fe) Lead (Pb) Manganese (Mn) max. 1 Mercury (Hg) Molybdenum (Mo) max. 1 Nickel (Ni) max. 1 Silver (Ag) max. 1 Thallium (TI) Tin (Sn) Vanadium (V) max. 1 Zinc (Zn) max. 10 Zirconium (Zr) Ultraviolet Absorbance (1.00-cm path vs water; curve smooth throughout stated range with no extraneous impurity peaks): at 335 nm act. value reported at 340 nm act, value reported at 350 nm act. value reported at 360 nm act. value reported at 380 nm act, value reported at 400 nm act, value reported Windows of Infrared Transmittance (0.1-mm path, **50-100% T)**, μm: 2.5-2.9 passes test 3.6-5.5 passes test 6.2-6.4 passes test 92-99 passes test 10.9-11.7 passes test

CONT.

вох

PACKING

Calibrate and standardise your analytical methods and equipment with J.T.Baker Volumetric and Buffer solutions.

Refer to the Analytical applications section of this catalogue for more details.

12.0-15.0

Methyl Isobutyl Ketone PHOTREX Reagent / For Spectrophoto

8224

► CH ₃ COCH ₂ CH(CH ₃) ₂			
M =	100.16 g/mol		
1 l =	0.80 kg		
FLASHPOINT	17 °C		
CAS NO.	108-10-1		
EINECS	203-550-1		
NC CODE	2914 13 00		
EC NO.	606 004 00 4		
UN/ID NO.	1245		
ADR/RID	3 F1		

IMDG 3/II **R:** 11-20-36/37-66 **S:** 16-29-9





Assay (by GC) (corrected for water)	min. 98.5%
Appearance	passes test
Color (APHA)	max. 10
Residue after Evaporation	max. 0.001%
Titrable Acid (meq/g)	max. 0.002
Water (H ₂ O)	max. 0.05%
Product Information (not specificati	ons):
Boiling Point (typical)	115.5°C
Density (g/ml) at 25°C (typical)	0.796
Ultraviolet Absorbance (1.00-cm pat	th vs water):
at 335 nm	max. 1.00
at 340 nm	max. 0.50
at 350 nm	max. 0.25
at 360 nm	max. 0.15
at 380 nm	max. 0.02
at 400 nm	max. 0.01
Windows of Infrared Transmittance	(0.1-mm path,
50-100% T), μm:	
2.5-3.2	passes test
3.6-5.5	passes test
6.2-6.4	passes test
9.2-9.9	passes test
10.9-11.7	passes test
12.0-15.0	passes test

PRODUCI	PACKING	CONT.
NO.		вох
8224.0500	500 ml	

Methyl Isobutyl Ketone

'BAKER ANALYZED' / ACS

8049

► CH₃COCH₂CH(CH₃)₂

M = 100.16 g/mol

1 l = 0.80 kg

FLASHPOINT 17 °C

CAS NO. 108-10-1

EINECS 203-550-1 NC CODE 2914 13 00

EC NO. 606 004 00 4

UN/ID NO. 1245

ADR/RID 3 F1

IMDG 3/II

R: 11-20-36/37-66

S: 16-29-9





Exceeds ACS Specifications	
Assay (by GC)	min. 99%
Appearance	clear
Color (APHA)	max. 15
Residue after Evaporation	max. 0.002%
Titrable Acid (meq/g)	max. 0.002
Water (H ₂ O)	max. 0.1%
Trace Impurities (in ppm):	
Aluminium (AI)	max. 0.5
Barium (Ba)	max. 0.1
Boron (B)	max. 0.02
Cadmium (Cd)	max. 0.05
Calcium (Ca)	max. 0.5
Chromium (Cr)	max. 0.02
Cobalt (Co)	max. 0.02
Copper (Cu)	max. 0.02
Iron (Fe)	max. 0.1
Lead (Pb)	max. 0.1
Magnesium (Mg)	max. 0.1
Manganese (Mn)	max. 0.02
Nickel (Ni)	max. 0.02
Tin (Sn)	max. 0.1
Zinc (Zn)	max. 0.1

PRODUCT	PACKING	CONT.
NO.		вох
8049.1000	11	6
8049.2500	2.5	4
8049.5000	5 EcoTainer	
8049.9025	25 l	
8049.9200	200	
a		6

EcoTainer, the metal solvent can for more safety in the lab. For safe handling of 25 l tin cans, see Self-closing tap.

Methyl Isobutyl Ketone MOS, VLSI Grade
See for detailed information section Microelectronic materials page 32 or the website www.jtbaker.com/micro



Methyl Orange

0.1% in H₂O / BAKER A	NALYZED				/ 136
NC CODE 3204 13 00	Visual Transition Interval:		PRODUCT	PACKING	CONT.
	pH 3.2	pink or red	NO.		вох
	pH 4.4	yellow	7136.0250	250 ml	
			7136.1000	11	
			Solution, rea	ady for use.	



Methy

Methyl Orange Sodium Salt

1145 BAKER ANALYZED' / ACS

► 4-(CH₃)₂NC₆H₄N:NC₆H₄-4-SO₃Na **M** = 327.34 g/mol

> **CAS NO.** 547-58-0 **EINECS** 208-925-3

NC CODE 2927 00 00

UN/ID NO. 3143 **ADR/RID** 6.1 T2

IMDG 6.1/III

R: 25

S: 22-24/25-45



Meets ACS Specifications Clarity of Solution passes test Visual Transition Interval: pH 3.2 pink or red yellow pH 4.4

PRODUCT CONT. **PACKING** NO. вох 1145.0100 100 g

C.I. 13025

PRODUCT

PRODUCT

9345.4000GL 4 | Glass

min. 99.5%

max. 0.01%

max. 200 ppm

max. 20

9261.4000GL 4 | Glass

PACKING

PACKING

CONT.

CONT.

BOX

BOX



Methyloxirane



Methylparaben

quirements of cGMP, USP/NF, BP/Ph.Eur, JP or FCC, page 36



4-Methyl-2-pentanone

See Methyl Isobutyl Ketone



Methyl Phenyl Ether



9261

2-Methyl-1-propanol

1-Methyl-2-pyrrolidone

'BAKER BIO-ANALYZED

► C₅H₉NO

M = 99.13 g/mol

1 l = 1.03 kg

FLASHPOINT 95 °C

CAS NO. 872-50-4

EINECS 212-828-1

NC CODE 2933 79 00

EC NO. 606 021 00 7

R: 36/38 **S:** 41



Assay (C₅H₅NO)	min. 99.5%
Color (APHA)	max. 20
Free Amines (as CH ₃ NH ₂)	max. 0.01%
Heavy Metals (as Pb)	max. 0.1 ppm
Residue after Ignition	max. 10 ppm
Water (H₂O)	max. 200 ppm
Ultraviolet Absorbance (1.00-cr	n path vs water):
at 300 nm	max. 0.5
at 325 nm	max. 0.1

at 350 nm max. 0.03 at 400 nm max. 0.01 max. 285 UV Cut-off, nm

1-Methyl-2-pyrrolidone 9345 'BAKER ANALYZED' / For peptide synthe

► C₅H₉NO

M = 99.13 g/mol

11 = 1.03 kg

 $\textbf{FLASHPOINT} \ 95\,^{\circ}\text{C}$

S:

CAS NO. 872-50-4

EINECS 212-828-1 NC CODE 2933 79 00

EC NO. 606 021 00 7

R: 36/38

41	
Xi	
inside and	

1-Methyl-2-Pyrrolidone (NMP) CMOS, Finyte Grade

Assay

Color (APHA)

Water (H₂O)

Free Amines (as CH₃NH₂)

1146

7137

J.T.Baker

Methyl Red

 $(CH_3)_2NC_6H_4N:NC_6H_4COOH$ **M** = 269.31 g/mol

CAS NO. 493-52-7 **EINECS** 207-776-1 **NC CODE** 2927 00 00

11 = 0.80 kg

NC CODE 2927 00 00

R: 11 S: 16-7/9

FLASHPOINT 11 °C

UN/ID NO. 1170 ADR/RID 3 F1 IMDG 3/II

Insoluble Matter in Alcohol	passes test
Visual Transition Interval:	
pH 4.2	pink
pH 6.2	yellow

PRODUCT	PACKING	CONT
NO.		вох
1146.0025	25 g Glass	

C.I. 13020.

pink

yellow

Methyl Red

0.1% in ethanol / 'BAKER ANALYZED'

Visual Transition Interval:

pH 4.2 pH 6.2

PRODUCT	PACKING	CONT
NO.		BOX
7137.0250	250 ml	
Excluding ex	cise	

Solution, ready for use.

Methyl Salicylate

See for qualities meeting the requirements of cGMP, USP/NF, BP/Ph.Eur, JP or FCC, page 36

Methylsulfoxide See Dimethyl Sulfoxide

•

Metol

See p-(Methylamino)phenol Sulfate



Mineral Oil

See for qualities meeting the requirements of cGMP, USP/NF, BP/Ph.Eur, JP or FCC, page 36



Mixed Calibration Standard III

(Matrix: 2% nitric acid) / 'BAKER INSTRA-ANALYZED' / Plasma Standard



R: 36/38 S: 26-37

NC CODE 3822 00 00

Elemental Concent (µg/ml):		
Arsenic (As)	500	
Molybdenum (Mo)	100	
Silicon (Si)	100	

PRODUCT	PACKING	CONT.
NO.		ВОХ
6003-01	100 ml	

For use in EPA SW-846 Methods 6010 and 200.7. Traceable to NIST.

Mixed Calibration Standard IV

(Matrix: 5% nitric acid) / 'BAKER INSTRA-ANALYZED' / Plasma Standard

6004-01

NC CODE 3822 00 00 UN/ID NO. 3264 ADR/RID 8 C1 IMDG 8/III R: 34 S: 23-26-36/37/39-45

corrosive

Elemental Concent (μg/ml):		
Aluminium (AI)	200	
Calcium (Ca)	1000	
Chromium (Cr)	20	
Nickel (Ni)	20	
Potassium (K)	400	
Sodium (Na)	200	

PRODUCT	PACKING	CONT.
NO.		вох
6004-01	100 ml	
	TDA CIA/ 040 A4 // / 0010	1000 7
For use in E	PA SW-846 Methods 6010 a	nd 200.7.
Tracoable to	TOIM	

Innovation is principal to our business.



Mixed

Mixed Calibration Standard V

Matrix: 5% HNO₃ and a trace of tartaric acid / 'BAKER INSTRA-ANALYZED' / Plasma Standard 6005-01

Elemental Concent (µg/ml):

NC CODE 3822 00 00 UN/ID NO. 3264 ADR/RID 8 C1 IMDG 8/III **R:** 34

Antimony (Sb) Boron (B) Magnesium (Mg) Silver (Ag) **S:** 23-26-36/37/39-45 Thallium (TI)

PRODUCT	PACKING	CONT
NO.		BO)
6005-01	100 ml	

For use in EPA SW-846 Methods 6010 and 200.7. Traceable to NIST.

200

100

1000

50

200



Mohr's Salt

e Ammonium Iron (II) Sulfate Hexahydrate

Molecular Sieve 3Å

8-12 mesh / 'BAKER ANALYZED'

NC CODE 3824 90 15 Loss on Drying at 300°C max. 5% **PRODUCT PACKING**

CONT. вох NO. 1569.1000 1 kg

Molecular Sieve 4Å

NC CODE 3824 90 15 Loss on Drying at 300°C max. 5% **PRODUCT PACKING** CONT. NO. BOX 0810.1000 1 kg Activated

5769

Molybdenum 1000 μg/ml

'BAKER INSTRA-ANALYZED' / Plasma Standard (Matrix 2% ammonium hydroxide)

► Mo M = 95.94 g/molNC CODE 3822 00 00

Certificate Provided Reporting Actual Lot Analysis Molybdenum (Mo) 998-1002 µg/ml

PRODUCT PACKING CONT. BOX NO. 5769.0100 100 ml

Prepared from the highest purity raw material available, generally greater than 99.999% spectral purity. The content of the solution is confirmed to be accurate to within \pm 0.2%. Typically 1000 $\mu\text{g/ml}.$ The certificate of analysis provided reports actual lot analysis. The certificate also lists the trace impurities.

6935

Molybdenum 1000 μg/ml

(Matrix 2% ammonium hydroxide) / 'BAKER INSTRA-ANALYZED' / Atomic Absorption Standard

► Mo **M** = 95.94 g/mol NC CODE 3822 00 00

Molybdenum (Mo) 998-1002 $\mu g/ml$ **PRODUCT PACKING** CONT. вох 6935.0100 100 ml 6935.0500 500 ml

Prepared by dissolution of high purity raw materials (min. 99.99% spectral purity). Assays are verified by ICP against standards traceable to NIST. Standard Reference Material numbers (SRM) are printed on each label.

6818

Molybdenum 1000 μg/ml

'BAKER ANALYZED' / Atomic Absorption Standard

► Mo **M** = 95.94 g/mol NC CODE 3822 00 00

998-1002 μg/ml Molybdenum (Mo)

PACKING CONT. PRODUCT NO. BOX 6818.0100 100 ml 6818.0500 500 ml

Ammonium heptamolybdate in ammoniumhydroxide 0.5

Morph



Molybdenum 10000 μ g/ml

(Matrix: 4% Ammonium Hydroxide) / 'BAKER INSTRA-ANALYZED' / Plasma Standard

5737

► Mo

M = 95.94 g/mol **NC CODE** 3822 00 00

Certificate Provided Reporting Actual Lot Analysis

Molybdenum (Mo) 9980-10020 μg/ml

 PRODUCT
 PACKING
 CONT.

 NO.
 BOX

 5737.0100
 100 ml

Prepared from the highest purity raw material available, generally greater than 99.999% spectral purity. The content of the solution is confirmed to be accurate to within 0.2 %. Typically 10000 $\mu g/ml$.

The certificate of analysis provided reports actual lot analysis. The certificate also lists the trace impurities.

Molybdic Acid

powder / 'BAKER ANALYZED' / ACS

0176

CONT.

CAS NO. 7782-91-4
EINECS 231-970-5
NC CODE 2825 70 00
UN/ID NO. 3288
ADR/RID 6.1 T5
IMDG 6.1/III

R: 36/37-48/20/22 S: 22-25

> Xn harmful

Meets ACS Specifications	
Assay (as MoO ₃)	min. 85.0%
Arsenate, Phosphate and Silicate (as SiO ₂)	max. 0.001%
Chloride (CI)	max. 0.002%
Heavy Metals (as Pb)	max. 0.003%
Insoluble in Dilute NH ₄ OH	max. 0.01%
Sulfate (SO ₄)	max. 0.2%
Trace Impurities (in ppm):	
Phosphate (PO ₄)	max. 5

0176.0250 250 g

Molybdic Acid Reagent is commonly known as Molybdic Acid 85% and consists largely of ammonium molybdate.

PACKING

PRODUCT

Monochloroacetic Acid

See Chloroacetic Acid



Monochlorobenzene

See Chlorobenzen



Monoethanolamine

See Ethanolamine



3-(4-Morpholine)propanesulfonic acid, sodium salt

See MOPS, Sodium salt



2-(N-Morpholino)-ethanesulfonic Acid Monohydrate

crystalline / 'BAKER ULTRAPURE BIOREAGENT

4014

► ÇH₂CH₂OCH₂CH₂NCH₂CH₂SO₃H.H₂O

M = 213.25 g/mol
CAS NO. 4432-31-9
EINECS 224-632-3
NC CODE 2934 90 96

Buffer for Liquid Chromatography and other Molecular Biology Applications min. 98% Appearance passes test Colour of a 1M Alkaline Solution passes test DNase Activity none detected Identification (by IR) passes test 7-10% Loss on Drying at 130°C Protease Activity none detected none detected

Product Information (not specifications):
pKa at 20°C 6.15
Trace Impurities (in ppm):

Heavy Metals (as Pb)

 PRODUCT NO.
 PACKING
 CONT.

 4014.0020
 20 g

 4014.0200
 200 g

 4014.1000
 1 kg

www.jtbaker.com/europe



Morph

Morpholinopropane Sulfonicacid

BAKER ULTRAPURE BIOREAGENT

► C₇H₁₅NO₄S

4004

M = 209.26 g/mol **CAS NO.** 1132-61-2 **EINECS** 214-478-5 NC CODE 2934 99 90

Assay (dried basis) min. 99% Appearance passes test DNase Activity none detected Loss on Drying max. 1% pH of 1% Solution at 25°C 2.5-4.5 Protease Activity none detected RNase Activity none detected Solution (10% in water) passes test Absorbance of a 0.1M Solution (1-cm path vs

CONT. **PRODUCT PACKING** NO. вох 4004.0100 100 g 4004.0500 500 g 4004.2500 2.5 kg

PACKING

25 g Glass

100 g

2.5 kg

PRODUCT

4163.0025

4163.0100

4163.2500

CONT.

BOX

water):

at 260 nm max. 0.02 Trace Impurities (in ppm): Heavy Metals (as Pb) max. 5 Iron (Fe) max. 5

MOPS, Sodium salt

4163

'BAKER ULTRAPURE BIOREAGENT'

► C₇H₁₄NNaO₄S

M = 231.25 g/mol **CAS NO.** 7119-22-7 NC CODE 2934 99 90

min. 99.0% Assay (dried basis) Appearance passes test **DNase Activity** none detected Heavy Metals (as Pb) max. 5 ppm Identification (by IR) passes test Insoluble Matter max. 1.0% Loss on Drying max. 1.0% pH of 0.1 M Solution at 25°C 10.0-11.0 Protease Activity none detected RNase Activity none detected Solubility passes test Absorbance of a 0.1 M Solution, Maximum (1-cm path vs water): max. 0.03 at 260 nm

Product Information (not specifications): pKa at 20°C 7.20



MOS chemicals

See for detailed information section Microelectronic materials page 32 or the website www.jtbaker.com/micro



Mounting medium

MountingClear

3904 HISTO GRADE

11 = 0.69 kg **FLASHPOINT** -12 °C **CAS NO.** 540-84-1 NC CODE 2901 10 00 **EC NO.** 601 009 00 8 **UN/ID NO.** 1262

ADR/RID 3 F1 IMDG 3/II

R: 11-38-50/53-65-67

S: 16-29-33-9



dangerous for the environment



Cytology Aromates (by UV) not detected Assay (by GC)(2,2,4 Trimethylpentane) min. 90% (w/w) Color (APHA) max. 10

intermediate during coverslipping in Histology and

Intended for use as a Alcohol-Mountant

PRODUCT PACKING CONT. NO. BOX 3904 2500 2.5 I Glass

Multicompendial Products



Multi-Element CCV Solution I & II



Multi-Element ICS Solution AB-1&AB-2 ee ICS Solution AB-1 & AB-2 Multi-Element ICV Solution I & II See ICV Solution I & II Multi-Element Interference Check Standard I **Multi Element Interference Check Standard II** See Interference Check Standard II **Multi-Element Interference Check Standard IV** See Interference Check Standard IV **Multi-Element Mixed Calibration Standard III Multi-Element Mixed Calibration Standard IV** See Mixed Calibration Standard IV **Multi-Element Mixed Calibration Standard V** See Mixed Calibration Standard V Multi-Element Primary Drinking Water Std I See Primary Drinking Water Standard I Multi-Element Trace Metal Standard I See Trace Metal Standard I **Multi-Element Trace metal Standard III** See Trace Metal Standard III Multi-Element Water Standard I & II See Water Standard I & II Murexide 1149 'BAKER ANALYZED' C₈H₈N₆O₆.H₂O Sensitivity as Indicator **PRODUCT PACKING** CONT. passes test **M** = 302.22 g/mol **CAS NO.** 3051-09-0 1149.0025 25 g Glass **EINECS** 221-266-6 NC CODE 2933 54 00 1-Naphthol 1150 ► C₁₀H₇OH min 99% PACKING CONT. Assav PRODUCT **M** = 144.17 g/mol Melting Point 93-96°C вох NO. **CAS NO.** 90-15-3 Residue after Ignition max. 1% 1150.0100 **EINECS** 201-969-4 NC CODE 2907 15 10 **EC NO.** 604 029 00 5 **R:** 21/22-37/38-41 **S:** 22-26-37/39

Napht

N-(1-Naphthyl) Ethylenediamine Dihydrochloride

C₁₀H₇NHCH₂CH₂NH₂.2HCl M = 259.18 g/mol

1341

CAS NO. 1465-25-4 **EINECS** 215-981-2

NC CODE 2921 59 90

R: 36/38 **S:** 24/25-26



Meets ACS Specifications. Meets Reagent Specifications for testing USP/NF monographs Sensitivity for Sulfanilamide passes test Solubility passes test Suitability for Nitrogen Dioxide Determination(ASTM D1607) passes test Water (H₂O) max. 5% PRODUCT CONT. **PACKING** вох 1341.0025 25 g Glass

Marshall's Reagent



Neocuproine Hydrochloride



ee for qualities meeting the requirements of cGMP, USP/NF, BP/Ph.Eur, JP or FCC, page 36

Nickel 1000 μg/ml

(Matrix: 2% nitric acid) / 'BAKER INSTRA-ANALYZED' / Plasma Standard

► Ni

M = 58.69 g/mol

NC CODE 3822 00 00 EC NO. 7 004 00 1

R: 36/37/38

S: 26-37

Certificate Provided Reporting Actual Lot Analysis

998-1002 $\mu g/ml$

5770.0100 100 ml

PRODUCT PACKING CONT.

Prepared from the highest purity raw material available, generally greater than 99.999% spectral purity. The content of the solution is confirmed to be accurate to within \pm 0.2%. Typically 1000 μ g/ml. The certificate of analysis provided reports actual lot analysis. The certificate also lists the trace impurities.

Nickel 1000 μg/ml

(Matrix: 2% nitric acid) / 'BAKER INSTRA-ANALYZED' / Atomic Absorption Standard

Nickel (Ni)

► Ni

M = 58.69 g/mol

NC CODE 3822 00 00

EC NO. 7 004 00 1

R: 36/38



998-1002 μg/ml Nickel (Ni)

NO.

998-1002 μg/ml

PRODUCT

6936.0100 100 ml 6936.0500 500 ml

PACKING

Prepared by dissolution of high purity raw materials (min. 99.99% spectral purity). Assays are verified by ICP against standards traceable to NIST. Standard Reference Material numbers (SRM) are printed on each label.

CONT.

BOX

Nickel 1000 μg/ml

'BAKER ANALYZED' / Atomic Absorption Standard

► Ni

M = 58.69 g/mol

NC CODE 3822 00 00

EC NO. 7 004 00 1 **R:** 36/38

S: 26-37



PRODUCT PACKING CONT. BOX 6819.0100 100 ml 6819.0500 500 ml

Nickel nitrate in nitric acid 0.5 mol/l.

Innovation is principal to our business.

Nickel 10000 μ g/ml 1.00% (w/v) / (Matrix: 2% nitric acid) / 'BAKER INSTRA-ANALYZED' / Plasma Standard

5738

► Ni

M = 58.69 g/mol

NC CODE 3822 00 00

EC NO. 7 004 00 1

R: 36/38 **S:** 26

Certificate Provided Reporting Actual Lot Analysis

9980-10020 μg/ml Nickel (Ni)

PRODUCT **PACKING** CONT. NO. вох 5738.0100 100 ml

Prepared from the highest purity raw material available, generally greater than 99.999% spectral purity. The content of the solution is confirmed to be accurate to within 0.2 %. Typically 10000 $\mu\text{g/ml.}$

The certificate of analysis provided reports actual lot analysis. The certificate also lists the trace impurities.

Nickel(II) Acetate 4-hydrate

0179

► (CH₃COO)₂Ni.4H₂O

M = 248.86 g/mol

CAS NO. 373-02-4

EINECS 206-761-7 NC CODE 2915 29 00

R: 22

S: 24



Assay (by EDTA titrn.)	min. 97.0%
Chloride (CI)	max. 0.001%
Cobalt (Co)	max. 0.2%
Copper (Cu)	max. 0.01%
Insoluble Matter	max. 0.005%
Iron (Fe)	max. 0.002%
Lead (Pb)	max. 0.005%
Nitrogen Compounds (as N)	max. 0.005%
pH of 5% Solution at 25°C	6.0-7.0
Sulfate (SO ₄)	max. 0.005%
Zinc (Zn)	max. 0.02%

PRODUCT	PACKING	CONT.
NO.		вох
0179.9050	50 kg	

Nickel(II) Carbonate

CAS NO. 3333-67-3 **EINECS** 222-068-2

NC CODE 2836 99 18

EC NO. 28 010 00 0

R: 22-40-43-50/53

S: 22-36/37-60-61





Assay (as Ni) (by EDTA titrn.)	min. 45.0%
Chloride (CI)	max. 0.002%
Cobalt (Co) (by AAS)	max. 0.1%
Copper (Cu) (by AAS)	max. 0.01%
Insoluble in HCl	max. 0.01%
Iron (Fe) (by AAS)	max. 0.004%
Lead (Pb) (by AAS)	max. 0.002%
Nitrogen Compounds (as N)	max. 0.02%
Sulfate (SO ₄)	max. 0.01%
Zinc (Zn) (by AAS)	may 0.01%

PRODUCT PACKING	CONT.	
NO.		вох
0181.0500	500 g	
0181.0300	300 g	

Nickel(II) Chloride Hexahydrate

0182

CONT.

0181

► NiCl₂.6H₂O

M = 237.71 g/mol

CAS NO. 7791-20-0

EINECS 231-743-0

NC CODE 2827 35 00

UN/ID NO. 3288

ADR/RID 6.1 T5

IMDG 6.1/III

R: 23/24/25E-45

S: 22-36/37-44



Assay (by EDTA titrn.)	min. 97.0%
Barium (Ba)	max. 0.005%
Calcium (Ca)	max. 0.005%
Cobalt (Co)	max. 0.002%
Insoluble Matter	max. 0.005%
Iron (Fe)	max. 0.002%
Lead (Pb)	max. 0.001%
Magnesium (Mg)	max. 0.005%
Nitrogen Compounds (as N)	max. 0.005%
pH of 5% Solution at 25°C	4.0-7.0
Potassium (K)	max. 0.001%
Sodium (Na)	max. 0.01%
Sulfate (SO ₄)	max. 0.005%
Trace Impurities (in ppm):	
Copper (Cu)	max. §

Lithium (Li)

Zinc (Zn)

NO.		B0)
0182.9050	50 kg	

PACKING

PRODUCT

max. 1

max. 50



Nicke

Nickel(II) Chloride Hexahydrate

0183 'BAKER

NiCl₂6H₂0

M = 237.71 g/mol

CAS NO. 7791-20-0

EINECS 231-743-0

NC CODE 2827 35 00

UN/ID NO. 3288

UN/ID NO. 3288 ADR/RID 6.1 T5 IMDG 6.1/III

R: 23/24/25E-45 **S:** 22-36/37-44

toxic

Assay (by EDTA titrn.)	min. 96%
Insoluble Substances	max. 0.01%
Lead (Pb)	max. 0.01%
pH of 5% Solution at 25°C	3.0-7.0

PRODUCT	PACKING	CONT.
NO.		вох
0183.2500	2.5 kg	

Nickel(II) Nitrate Hexahydrate

0184 'BAKER ANALYZED

M = 290.81 g/mol

CAS NO. 13478-00-7

EINECS 236-068-5

NC CODE 2834 29 20

UN/ID NO. 2725

ADR/RID 5.1 O2

IMDG 5.1/III

R: 22-43-8

S: 24-36/37

Xn

Assay (by EDTA titrn.) min. 99.0% max. 0.05% Ammonium (NH₄) max. 0.001% Chloride (CI) Cobalt (Co) (by AAS) max. 0.05% Insoluble Matter max. 0.005% Iron (Fe) (by AAS) max. 0.001% Lead (Pb) (by AAS) max. 0.002% pH of 5% Solution at 25°C 3.5-5.5 max. 0.003% Sulfate (SO₄) Zinc (Zn) (by AAS) max. 0.01% Trace Impurities (in ppm): Copper (Cu) (by AAS) max. 5

 PRODUCT NO.
 PACKING BOX

 0184.0500
 500 g

Nickel(II) Oxide Green

Powder / 'BAKER ANALYZED'

M = 74.71 g/mol

CAS NO. 1313-99-1

EINECS 215-215-7

NC CODE 2825 40 00

EC NO. 28 003 00 2

R: 43-49-53

S: 45-53-61

min 99.0% Assay (by EDTA titrn.) Aluminium (AI) max. 0.005% Average Particle Diameter, μm (APD) (by Sedigraph)(typical) max. 3 max. 0.005% Calcium (Ca) Chloride (CI) max. 0.001% Cobalt (Co) max. 0.2% max. 0.005% Copper (Cu) max. 0.01% Lead (Pb) max. 0.01% max. 0.005% Nitrogen Compounds (as N) max. 0.005% Potassium (K) Silicon (Si) max. 0.01% Sodium (Na) max 0.02% Total Sulfur (S) max. 0.01% max. 0.02% Zinc (Zn)

PRODUCT PACKING CONT.

NO. BOX

1153.2500 2.5 kg

Would you like to receive an extra copy of this catalogue?

Call our Customer Service department. Contact details on the rear cover of this catalogue.

Nickel(II) Sulfate Hexahydrate

'BAKER ANALYZED' / ACS 0185

► NiSO₄.6H₂O

M = 262.86 g/molCAS NO. 10101-97-0EINECS 232-104-9

NC CODE 2833 24 00 EC NO. 28 009 00 5

R: 22-40-42/43-50/53

S: 22-36/37-60-61



environment



Exceeds ACS Specifications	
Assay	98.0-102.0%
Calcium (Ca)	max. 0.005%
Chloride (CI)	max. 0.001%
Cobalt (Co)	max. 0.002%
Copper (Cu)	max. 0.005%
Insoluble Matter	max. 0.005%
Iron (Fe)	max. 0.001%
Lithium (Li)	max. 1 ppm
Magnesium (Mg)	max. 0.005%
Manganese (Mn)	max. 0.002%
Nitrogen Compounds (as N)	max. 0.002%
pH of 5% Solution at 25°C	3.0-5.0
Potassium (K)	max. 0.001%
Sodium (Na)	max. 0.01%
Zinc (Zn)	max. 50 ppm

PRODUCT NO.	PACKING	CONT. BOX
0185.0250	250 g	
0185.1000	1 kg	

Niobium 1000 μ g/ml

(Matrix: H₂O plus a trace of hydrofluoric acid) / 'BAKER INSTRA-ANALYZED' / Plasma Standard

5771

► Nb

M = 92.91 g/mol **NC CODE** 3822 00 00

R: 20/21/22-36

S: 26-36/37



Certificate Provided Reporting Actual Lot AnalysisNiobium (Nb) 998-1002 μg/ml

 PRODUCT
 PACKING
 CONT.

 NO.
 BOX

 5771.0100
 100 ml

Prepared from the highest purity raw material available, generally greater than 99.999% spectral purity. The content of the solution is confirmed to be accurate to within $\pm~0.2\%$. Typically 1000 $\mu\text{g/ml}$.

The certificate of analysis provided reports actual lot analysis. The certificate also lists the trace impurities.

Niobium 10000 μg/ml

1.00% (w/v) / (Matrix: H_2O plus a trace of hydrofluoric acid) / 'BAKER INSTRA-ANALYZED' / Plasma Standard

5760

CONT

вох

► Nb

M = 92.91 g/mol **NC CODE** 3822 00 00

R: 20/21/22-36

S: 26-36/37



Certificate Provided Reporting Actual Lot AnalysisNiobium (Nb)9980-10020 μg/ml

 PRODUCT
 PACKING

 NO.
 5760.0100

 100 ml

Prepared from the highest purity raw material available, generally greater than 99.999% spectral purity. The content of the solution is confirmed to be accurate to within 0.2 %. Typically 10000 μ g/ml.

The certificate of analysis provided reports actual lot analysis. The certificate also lists the trace impurities.

The J.T.Baker CYCLE-TAINER High Purity Solvent Delivery System, preserves purity and protects people.

See chapter 3 of this catalogue for product details.



Nitri

Nitric Acid

6901 70% / ULTREX II Ultrapure Reagent

NC CODE 2808 00 00 EC NO. 7 004 00 1 UN/ID NO. 2031

R: 35

ADR/RID 8 CO1

S: 23-26-36-45



Assay	67-70% (w/w)
Trace Impurities (in ppt) (pg/g):	
Aluminium (AI)	max. 20
Antimony (Sb)	max. 10
Arsenic (As)	max. 20
Barium (Ba)	max. 10
Beryllium (Be)	max. 10
Bismuth (Bi)	max. 10
Boron (B)	max. 20
Cadmium (Cd)	max. 10
Calcium (Ca)	max. 20
Cerium (Ce)	max. 10
Cesium (Cs)	max. 10
Chromium (Cr)	max. 20
Cobalt (Co)	max. 10
Copper (Cu)	max. 20
Dysprosium (Dy)	max. 1
Erbium (Er)	max. 1
Europium (Eu)	max. 1
Gadolinium (Gd)	max. 1
Gallium (Ga)	max. 10
Germanium (Ge)	max. 10
Gold (Au)	max. 20
Hafnium (Hf)	max. 10
Holmium (Ho)	max. 1
Indium (In)	max. 1
Iron (Fe)	max. 20
Lanthanum (La)	max. 1
Lead (Pb)	max. 10
Lithium (Li)	max. 10
Lutetium (Lu)	max. 1
Magnesium (Mg)	max. 10
Manganese (Mn)	max. 10
Mercury (Hg)	max. 100
Molybdenum (Mo)	max. 10
Neodymium (Nd)	max. 1
Nickel (Ni)	max. 50
Niobium (Nb)	max. 1
Palladium (Pd)	max. 20
Platinum (Pt)	max. 20

Certificate Provided Reporting Actual Lot Analysis

Potassium (K)	max. 10
Praseodymium (Pr)	max. 1
Rhenium (Re)	max. 10
Rhodium (Rh)	max. 10
Rubidium (Rb)	max. 10
Ruthenium (Ru)	max. 20
Samarium (Sm)	max. 1
Scandium (Sc)	max. 10
Selenium (Se)	act. value reported
Silver (Ag)	max. 10
Sodium (Na)	max. 10
Strontium (Sr)	max. 10
Tantalum (Ta)	act. value reported
Tellurium (Te)	max. 1
Terbium (Tb)	max. 1
Thallium (TI)	max. 10
Thorium (Th)	max. 1
Thulium (Tm)	max. 1
Tin (Sn)	max. 20
Titanium (Ti)	max. 10
Tungsten (W)	max. 10
Uranium (U)	max. 1
Vanadium (V)	max. 10
Ytterbium (Yb)	max. 1
Yttrium (Y)	max. 1
Zinc (Zn)	max. 20
Zirconium (Zr)	max. 10

PRODUCT	PACKING	CONT.
NO.		вох
6901.0500	500 ml Fluoropolymer, pre-leached	
6901.1000	1 l Fluoropolymer, pre-leached	
6901.2000	2 Fluoropolymer, pre-leached	



Nitric Acid 70% CMOS, Finyte Grade

See for detailed information section Microelectronic materials page 32 or the website www.jtbaker.com/micro

Mallinckrodt Baker's cGMP Manufactured Chemicals for the Biopharmaceutical industry are a necessity for uncomplicated scale-up.

See chapter 6 of this catalogue.

9598

Nitric Acid

69.0-70.0 $\%\,$ / 'BAKER INSTRA-ANALYZED' / for Trace Metal Analysis / ACS

M = 63.01 g/mol 1 I = 1.42 kg CAS NO. 7697-37-2 EINECS 231-714-2 NC CODE 2808 00 00 EC NO. 7 004 00 1 UN/ID NO. 2031 ADR/RID 8 CO1 IMDG 8/II

R: 35

S: 23-26-36-45 C corrosive

Meets ACS Specifications	
Assay	69.0-70.0%
Appearance	passes test
Color (APHA)	max. 10
Residue after Ignition	max. 2 ppm
Specific Gravity at 60°/60°F	1.416-1.420
Trace Impurities (in ppb):	
Aluminium (AI)	max. 30
Arsenic and Antimony (as As)	max. 5
Barium (Ba)	max. 1
Beryllium (Be)	max. 1
Bismuth (Bi)	max. 1
Boron (B)	max. 4
Cadmium (Cd)	max. 1
Calcium (Ca)	max. 50
Chromium (Cr)	max. 10
Cobalt (Co)	max. 1
Copper (Cu)	max. 1
Gallium (Ga)	max. 20
Germanium (Ge)	max. 4
Gold (Au)	max. 4
Heavy Metals (as Pb)	max. 100
Iron (Fe)	max. 10
Lead (Pb)	max. 0.5
Lithium (Li)	max. 1
Magnesium (Mg)	max. 7
Manganese (Mn)	max. 1
Mercury (Hg)	max. 0.5
Molybdenum (Mo)	max. 5
Nickel (Ni)	max. 1
Niobium (Nb)	max. 1
Potassium (K)	max. 5
Silicon (Si)	max. 20
Silver (Ag)	max. 1
Sodium (Na)	max. 200

max. 1

max. 2

max. 5

max. 5

max. 1 max. 5

Zircomum (2	<u>~1)</u>	IIIdX. I
Trace Impu	ırities (in ppm):	
Chloride (CI)		max. 0.04
Phosphate (PO ₄)	max. 0.1
Sulfate (SO ₄	1)	max. 0.4
PRODUCT	PACKING	CONT
NO.		ВОХ
9598.0500	500 ml	ϵ
9598.2500	2.5	4

Use J.T.Baker Ultrex II and BAKER INSTRA-ANALYZED acids for low level trace element analysis.

Strontium (Sr)

Tantalum (Ta)

Vanadium (V)

Thallium (TI)

Tin (Sn)

Zinc (Zn)

See chapter 3 of this catalogue for more details.



Nitri

Nitric Acid

6019 69.0-70.0 % / 'BAKER ANALYZED' / ACS

M = 63.01 g/mol
1 I = 1.42 kg
CAS NO. 7697-37-2
EINECS 231-714-2
NC CODE 2808 00 00
EC NO. 7 004 00 1
UN/ID NO. 2031
ADR/RID 8 CO1
IMDG 8/II
R: 35
S: 23-26-36-45

eagent		
Specifications for testing USP/NF monographs		
69.0-70.0%		
passes test		
max. 10		
max. 4 ppm		
1.416 - 1.420		
max. 100		
max. 4		
max. 50		
max. 200		
max. 100		
max. 50		
max. 200		
max. 100		
max. 5		
max. 50		
max. 300		
max. 200		
max. 200		
max. 200		
max. 0.1		
max. 0.2		

max. 0.5

PRODUCT NO.	DUCT PACKING	CONT
		B0
6019.2500	2.5	
6019.2500PE	2.5 I HDPE	



Nitric Acid, 69.0-70.0%

see for qualities meeting the requirements of cGMP, USP/NF, BP/Ph.Eur, JP or FCC, page 36

Sulfate (SO₄)

Calibrate and standardise your analytical methods and equipment with J.T.Baker Volumetric and Buffer solutions.

Refer to the Analytical applications section of this catalogue for more details.



Nitric Acid

65% (max. 5 ppb Hg) / 'BAKER ANALYZED'

Assay

6080

CONT.

вох

6

► HNO ₃	
M =	63.01 g/mol
1 i =	1.40 kg
CAS NO.	7697-37-2
EINECS	231-714-2
NC CODE	2808 00 00
EC NO.	7 004 00 1
UN/ID NO.	2031
ADR/RID	8 CO1
IMDG	8/II
R:	35
S:	23-26-36/37/39-45
	С

Chloride (CI)	max. 0.5 ppm
Dithizone determination	passes test
Mercury (Hg)	max. 5 ppb
Residue after Ignition	max. 5 ppm
Sulfate (SO ₄)	max. 1 ppm
Trace Impurities (in ppm):	
Aluminium (AI)	max. 0.05
Arsenic (As)	max. 0.05
Barium (Ba)	max. 0.02
Beryllium (Be)	max. 0.01
Cadmium (Cd)	max. 0.02
Calcium (Ca)	max. 0.5
Chromium (Cr)	max. 0.1
Cobalt (Co)	max. 0.01
Copper (Cu)	max. 0.01
Germanium (Ge)	max. 0.05
Iron (Fe)	max. 0.1
Lead (Pb)	max. 0.05
Lithium (Li)	max. 0.01
Magnesium (Mg)	max. 0.1
Manganese (Mn)	max. 0.01
Molybdenum (Mo)	max. 0.02
Nickel (Ni)	max. 0.05
Potassium (K)	max. 0.1
Silver (Ag)	max. 0.01
Sodium (Na)	max. 0.5
Strontium (Sr)	max. 0.01
Thallium (TI)	max. 0.05
Titanium (Ti)	max. 0.1
Vanadium (V)	max. 0.01
Zinc (Zn)	max. 0.05
· · · · · · · · · · · · · · · · · ·	

6080.9025

max. 0.1

PRODUCT

6080.1000

6080.2500

6080.2500PE

NO.

min. 65%

PACKING

11

2.5 l

2.5 | HDPE

Nitric Acid 65% CMOS, Finyte Grade

Nitric Acid

R: 34 **S:** 26-36/39-45 14

ctronic materials page 32 or the website www.jtbaker.com/micro

Zirconium (Zr)



2 mol/l / 'BAKER ANALYZED' 7119 1.995 - 2.005 PACKING CONT. ► HNO₃ Titer (mol/l) **PRODUCT M** = 63.01 g/mol BOX **CAS NO.** 7697-37-2 7119.1000 11 6 **EINECS** 231-714-2 7119.5000 5 I NC CODE 2808 00 00 20 I 7119.9020 **UN/ID NO.** 2031 ADR/RID 8 CO1 Volumetric Solution, ready for use.

Find more Chromatography information at www.jtbaker.com/chromatography



Nitri

Nitric Acid

7622 1 mol/l / 'BAKER ANALYZED'

► HNO ₃		Titer (mol/l)	0.995-1.005	PRODUCT	PACKING	CONT.
M =	63.01 g/mol			NO.		вох
CAS NO.	7697-37-2			7622.5000	5 HDPE	
EINECS	231-714-2			7622.9020	20 l Polycube	
NC CODE	2808 00 00					
UN/ID NO.	2031			Volumetric Sc	olution, ready for use.	
ADR/RID	8 CO1					
IMDG	8/11					
R:	34					
S:	26-36/39-45					
	C Corrosive					

Nitric Acid

4868 1 mol/l / DILUT

 PRODUCT NO.
 PACKING PACKING
 CONT. BOX

 4868
 1 amp.
 6

Volumetric Concentrate, for dilution to 1 l.

Nitric Acid

4712 0.1 mol/l / DILUT-IT

HNO ₃		PRODUCT	PACKING	CONT.
M =	63.01 g/mol	NO.		BOX
CAS NO.	7697-37-2	4712	1 amp.	
EINECS	231-714-2			
NC CODE	2808 00 00	Volumetric C	oncentrate, for dilution to 1 l.	
EC NO.	7 004 00 1			
UN/ID NO.	2031			
ADR/RID	8 CO1			
IMDG	8/II			
R:	34			
S:	23-26-36-45			
	C corrosive			

p-Nitrophenylphosphate, Disodium Salt

3351 'BAKER ANALYZED

► NO ₂ C ₆ H ₄ OPO ₃ Na ₂ .6H ₂ O	Assay	about 70%	PRODUCT	PACKING	CONT.
M = 371.15 g/mol	Free p-nitrophenol	max. 0.07%	NO.		вох
CAS NO. 4264-83-9			3351.0005	5 g	
EINECS 224-246-5			3351.0025	25 g Glass	
NC CODE 2919 00 90					

Find our up-to-date Product Literature at www.jtbaker.com/europe

вох

8159

2818

CONT.

n-Nonane

8280 'BAKER' **PRODUCT PACKING** CONT. Assay (by GC) min. 99%

150-151°C

0.717-0.718

min. 98%

125-126°C

max. 0.02

max. 0.02

PRODUCT

8159.0500

NO.

PACKING

500 ml

NO.

8280.2500

2.5 l

Boiling Point

Assay (by GC)

Boiling Point

Density (g/ml) at 20°C



M = 128.26 g/mol

1 l = 0.72 kg

FLASHPOINT 31 °C

CAS NO. 111-84-2

EINECS 203-913-4

NC CODE 2901 10 00

UN/ID NO. 1920

ADR/RID 3 F1

IMDG 3/III

R: 10-51/53

S: 16-61



environment

Octadecanoic Acid

See Stearic Acid

cis-9-Octadecenoic Acid

See Oleic Acid

iso-Octane

See 2,2,4-Trimethylpentane

n-Octane

'BAKER ANALYZED'

► CH₃(CH₂)₆CH₃ **M** = 114.23 g/mol

1 l = 0.70 kg

FLASHPOINT 12 °C

CAS NO. 111-65-9

EINECS 203-892-1

NC CODE 2901 10 00

EC NO. 601 009 00 8

UN/ID NO. 1262

ADR/RID 3F1

IMDG 3/II

R: 11-38-50/53-65-67

S: 16-29-33-60-61-62-9







highly flammable

1-Octanesulfonic Acid Sodium Salt

► CH₃(CH₂)₇SO₃Na

at 250 nm

at 260 nm

M = 216.27 g/mol

CAS NO. 5324-84-5 EINECS 226-195-4

NC CODE 2904 10 00

R: 36/38

S: 26



For Ion-Pair Chromatography	
Assay (acidimetric)	min. 98.0%
UV Absorbance of 0.25M Solution:	
at 200 nm	max. 0.10
at 210 nm	max. 0.05
at 220 nm	max. 0.04
at 230 nm	max. 0.03
at 240 nm	max. 0.02

PRODUCT	PACKING	CONT.
NO.		вох
2818.0025	25 g Glass	

Octan

7060

1-Octanol 'BAKER ANALYZED

► CH₃(CH₂)₆CH₂OH

M = 130.23 g/mol

1 l = 0.82 kg

FLASHPOINT 81 °C

CAS NO. 111-87-5

EINECS 203-917-6

NC CODE 2905 1680

R: 36/38

S: 23

Assay (by GC) min. 99% Boiling Range 193-195°C Density (g/ml) at 20°C 0.824-0.825 Water (H₂O) max. 0.1%

178-181°C

0.812-0.818

PRODUCT

7059.0100

7059.1000

PACKING

100 ml

11

CONT.

BOX

CONT. **PRODUCT PACKING** NO. вох 7060.1000 11 7060.9025 25 I For safe handling of 25 l tin cans, see Self-closing tap.

2-Octanol 7059

► CH₃(CH₂)₅CHOHCH₃

M = 130.23 g/mol

1 l = 0.82 kg

FLASHPOINT 74 °C

CAS NO. 123-96-6

EINECS 204-667-0 NC CODE 2905 16 20

R: 36/38 **S:** 23





Oil of Cedarwood

See Cedarwood Oil

6035

Oleic Acid

► CH₃(CH₂)₇CH:CH(CH₂)₇COOH

M = 282.47 g/mol

1 l = 0.89 kg


FLASHPOINT 189 °C

CAS NO. 112-80-1

EINECS 204-007-1

NC CODE 2916 15 00

Acid Value	196-204
Congealing Temperature	max. 10°C
lodine Value	85-95
Mineral Acids	passes test
Neutral Fat or Mineral Oil	passes test
Residue on Ignition	max. 0.01%
Specific Gravity	0.889-0.895





Oleic Acid

See for qualities meeting the requirements of cGMP, USP/NF, BP/Ph.Eur, JP or FCC, page 36

Boiling Point

Density (g/ml) at 25°C



Oleic Acid Sodium Salt

4425

Orange G

'BAKER ANALYZED' / Certified Stain

C₁₆H₁₀N₂Na₂O₇S₂

M = 452.38 g/mol **CAS NO.** 1936-15-8

EINECS 217-705-6

NC CODE 3204 12 00

Certified for Use in Histology (C.I. 16230) Absorbance at maximum (3.0 mg/200 ml in H₂O, 1-cm path) act. value reported Absorbance Maximum, nm act. value reported Biological Test passes test Total Dye Content act. value reported

PRODUCT	PACKING	CONT.
NO.		вох
4425.0025	25 g Glass	
Certified by	the Biological Stain Commission.	



Orange III

See Methyl Orange Sodium Salt



Orthoperiodic Acid

dic Acid

Orthophosphoric Acid

Oxalic Acid

0.05 mol/l 0.1N / 'BAKER ANALYZED'

7141

EINECS 205-634-3

Titer (mol/l) NC CODE 2917 11 00

PRODUCT PACKING NO.

CONT. BOX

7141.1000 11

0.0497-0.0503

7141.9020 20 l Polycube

Volumetric Solution, ready for use.

Each lot of this product is standardized potentiometrically against NIST tracable reference standard.

Oxalic Acid

 $0.05 \text{ mol/I} / \text{DILUT-IT} / {}^{1}/{}_{10} \text{ equiv.} = 4.502g; 0.1 N$

4665

▶ нососоон

M = 90.04 g/mol

EINECS 205-634-3 NC CODE 2917 11 00 **PRODUCT PACKING** CONT. 4665 1 amp.

Volumetric Concentrate, for dilution to 1 l.

Oxalic Acid Dihydrate 'BAKER ANALYZED' / ACS

0187

► HOCOCOOH.2H₂O

M = 126.07 g/mol **CAS NO.** 6153-56-6

EINECS 205-634-3

NC CODE 2917 11 00

EC NO. 607 006 00 8

R: 21/22 **S:** 24/25

Exceeds ACS Specifications	
Assay	99.5-102.5%
Insoluble Matter	max. 0.005%
Nitrogen Compounds (as N)	max. 0.001%
Residue after Ignition	max. 0.01%
Substances Darkened by Hot H ₂ SO ₄	passes test
Sulfate (SO ₄)	max. 0.002%
Trace Impurities (in ppm):	
Calcium (Ca)	max. 5

rrace impurities (in ppin).	
Calcium (Ca)	max. 5
Chloride (CI)	max. 5
Heavy Metals (as Pb)	max. 5
ron (Fe)	max. 2

PRODUCT	PACKING	CONT
NO.		B0
0187.1000	1 kg	
0187.9025	25 kg	

Oxidizing (ABI)

'BAKER ANALYZED' / for DNA/RNA synthesis

9488

FLASHPOINT -14 °C NC CODE 3822 00 00 **UN/ID NO.** 1993

ADR/RID 3 F1 IMDG 3/II

R: 11-19-20/21/22-36/37

S: 16-23-25-36/37-9





Sultable for Oligonacieolade Synthesis		
Assay Pyridine	18.5-21.0% (v/v)	
Molarity (lodine)	0.018-0.022	

PRODUCT	PACKING	CONT.
NO.		вох
9488.0200	200 ml Glass	
9488.0450	450 ml Glass	
9488.2000	2 I Glass	
9488.2500	2.5 Glass	

Oxidizing (ABI)

for DNA/RNA synth

9519

FLASHPOINT −21 °C NC CODE 3822 00 00 **UN/ID NO.** 1993

ADR/RID 3,3b

IMDG 3.1/II

R: 11-19-36/37

S: 16-26-33-9





Suitable for Oligonucleotide Synthesis	•	PRO
Assay lodine (g/l)	3.87 - 4.73	NO.
Assay Iodine (mol/I)	0.015 - 0.019	951
		951

PRODUCT	PACKING	CONT
NO.		B0
9519.0200	200 ml Glass	
9519.0450	450 ml Glass	
9519.0900	900 ml Glass	

Oxidi

Oxidizing reagents for use in DNA synthesis

ation section Reagents for DNA/RNA Synthe



Oxidizing reagents for use in DNA/RNA synthesis

d information section Reagents for DNA/RNA Synth

Boiling Point

2,2'-Oxydiethanol

► HOCH₂CH₂OCH₂CH₂OH

M = 106.12 g/mol

1 l = 1.11 kg FLASHPOINT 124 °C

CAS NO. 111-46-6

EINECS 203-872-2 NC CODE 2909 41 00

Packaging Options

d information section Packaging Options, page 12

Palladium 1000 μg/ml

5772 (Matrix: 2% nitric acid) / 'BAKER INSTRA-ANALYZED' / Plasma Standard

► Pd

M = 106.42 g/mol

NC CODE 3822 00 00

R: 36/38

S: 26



Certificate Provided Reporting Actual Lot Analysis Palladium (Pd)

998-1002 μg/ml

242-250°C

PRODUCT

8564.2500

PRODUCT PACKING NO. 5772.0100 100 ml

PACKING

2.5 [

CONT.

CONT.

CONT.

CONT.

BOX

BOX

вох

BOX

Prepared from the highest purity raw material available, generally greater than 99.999% spectral purity. The content of the solution is confirmed to be accurate to within \pm 0.2%. Typically 1000 μ g/ml. The certificate of analysis provided reports actual lot

analysis. The certificate also lists the trace impurities.

PACKING

Palladium 1000 μg/ml

(Matrix: 2% nitric acid) / 'BAKER INSTRA-ANALYZED' / Atomic Absorption Standard

► Pd

6953

M = 106.42 g/mol

NC CODE 3822 00 00

R: 36/38

S: 26



998-1002 μg/ml Palladium (Pd)

NO. 6953.0100 6953.0500 500 ml

PRODUCT

Prepared by dissolution of high purity raw materials (min. 99.99% spectral purity). Assays are verified by ICP against standards traceable to NIST. Standard Reference Material numbers (SRM) are printed on each label.

5739

Palladium 10000 μg/ml

(Matrix: 2% nitric acid) / 'BAKER INSTRA-ANALYZED' / Plasma Standard

► Pd

M = 106.42 g/mol

NC CODE 3822 00 00 R: 36/38

S: 26-37



Certificate Provided Reporting Actual Lot Analysis 9980-10020 $\mu g/ml$ Palladium (Pd)

NO.

PRODUCT PACKING 5739.0100 100 ml

Prepared from the highest purity raw material available, generally greater than 99.999% spectral purity. The content of the solution is confirmed to be accurate to within 0.2 %. Typically 10000 µg/ml.

The certificate of analysis provided reports actual lot analysis. The certificate also lists the trace impurities.

1156

der / 'BAKER ANALYZED'

► C₅H₄NN:NC₁₀H₆OH

M = 249.28 g/mol CAS NO. 85-85-8

EINECS 201-637-9 NC CODE 2933 39 95 Sensitivity as Indicator

passes test

PACKING	CON
	ВО
1 g	

286



Pancreatin

eeting the requirements of cGMP, USP/NF, BP/Ph.Eur, JP or FCC, page 36 ee for qualities m



3864

Papanicolaou solution 2a, Orange G solution OG6

FLASHPOINT 12 °C

NC CODE 2207 10 00

UN/ID NO. 1992 ADR/RID 3 FT1 IMDG 3/II

R: 11-20/21/22-68/20/21/22

S: 16-36/37-7/9





Staining solution used in the Papanicolaou procedure. For cytological cancer and cycle diagnosis

CONT. **PRODUCT PACKING** BOX 3864.1000 1 | Glass 3864.2500 2.5 | Glass

Papanicolaou 2a Orange G

ailed information www.jtbaker.com and select Clinical



Papanicolaou solution 2b, Orange II solution

NC CODE 2207 10 00

UN/ID NO. 1992

FLASHPOINT 12 °C

ADR/RID 3 FT1

IMDG 3/II

R: 11-20/21/22-68/20/21/22

S: 16-36/37-7/9





Staining solution used in the Papanicolaou procedure. For cytological cancer and cycle diagnosis

PRODUCT	PACKING	CONT		
NO.		ВО		
3865.1000	1 Glass			
3865.2500	2.5 Glass			
Excludina ex	cise.			

Papanicolaou 2b Orange II

See for detailed information www.jtbaker.com and select Clinical



3866

CONT.

BOX

Papanicolaou solution 3b, Polychrome solution EA 50

FLASHPOINT 12 °C NC CODE 2207 10 00 **UN/ID NO.** 1992

> ADR/RID 3 FT1 IMDG 3/II

> > R: 11-23/24/25-39/23/24/25

S: 16-36/37-45-7/9

Paraffin

FLASE





Staining solution used in the Papanicolaou procedure. For cytological cancer and cycle diagnosis

NO. 3866.1000 1 | Glass 3866.2500 2.5 | Glass

PACKING

PRODUCT



Papanicolaou 3b Polychrome Solution



	Liquid / BAKEN					7100
1 l =	0.88 kg	Color (APHA)	max. 20	PRODUCT	PACKING	CONT.
LASHPOINT	199 °C	Identification (by IR)	passes test	NO.		вох
CAS NO.	8012-95-1			7160.1000	11	6
EINECS	232-384-2			7160.5000	5 EcoTainer	
NC CODE	2712 20 10			7160.9200	200 l	

Paraf

Paraffin

9513

ellets / 'BAKER

CAS NO. 8002-74-2 52-54°C **PRODUCT PACKING** CONT. Melting Point **EINECS** 232-315-6 NO. вох NC CODE 2712 20 10 9513.1000 1 kg 9513.9025 25 kg

PRODUCT

3451.1500

NO.

min. 95%

34-37°C

PRODUCT

8239.2500

8239.9200

NO.

PACKING

2.5 I

200 l

PACKING

125 ml glass x 12

CONT.

CONT.

BOX

BOX

Paraffin Cleaner

3451 HISTO GRADE

FLASHPOINT 57 °C

NC CODE 3822 00 00 **UN/ID NO.** 1993 ADR/RID 3 F1 IMDG 3/III

> **R:** 41-65-67 **S:** 16-23-26-39



Paraffin for Histology

See UltraPar (54-56°C)



Paraffin pellets for Histology

n and select Clinical d information www.jtbaker.co

Assay

Boiling Range

1157

Paraformaldehyde

► (CH₂O)n **M** = 30.03 g/mol FLASHPOINT 70 °C **CAS NO.** 30525-89-4

EINECS 200-001-8 NC CODE 2912 60 00

UN/ID NO. 2213 **ADR/RID** 4.1 F1 IMDG 4.1/III

R: 20/22-36/37/38-40-43

S: 22-26-36/37



Reaction passes test Residue after Ignition max. 0.1% Solubility in Ammonium Hydroxide passes test

PRODUCT	PACKING	CONT.
NO.		вох
1157.1000	1 kg	

Pentane

/ 'BAKER'

► C₅H₁₂

8239

M = 72.15 g/mol

1 l = 0.62 kg

FLASHPOINT −40 °C **CAS NO.** 109-66-0

EINECS 203-692-4 NC CODE 2901 10 00

EC NO. 601 006 00 1 **UN/ID NO.** 1265

ADR/RID 3 F1 IMDG 3/I

R: 12-51/53-65-66-67

S: 16-29-33-61-62-9







dangerous for the extremely flammable environment



iso-Pentane

9333

n-Pentane

'BAKER ULTRA RESI-ANALYZED' / for Organic Residue Analysis

► CH₃(CH₂)₃CH₃

M = 72.15 g/mol

1 l = 0.62 kg

FLASHPOINT $-40\,^{\circ}\mathrm{C}$

CAS NO. 109-66-0

EINECS 203-692-4

NC CODE 2901 10 00 EC NO. 601 006 00 1

UN/ID NO. 1265

ADR/RID 3F1

IMDG 3/I


R: 12-51/53-65-66-67

S: 16-29-33-61-62-9









min. 99.0%

min. 98.0%

max. 2

min. 99%

max. 5 ppm

Assay (by GC) (corrected for water)

PRODUCT NO.	PACKING	CONT. BOX
9333.1000	11	6
9333.4000	4 l Glass	4

n-Pentane

'BAKER HPLC ANALYZED' / for use in High Performance Liquid Chromatography

Assay (by GC)

Assay (by GC)

Residue after Evaporation

Residue after Evaporation (in ppm)

9331

► CH₃(CH₂)₃CH₃

M = 72.15 g/mol

1 l = 0.62 kg

FLASHPOINT −40 °C

CAS NO. 109-66-0

EINECS 203-692-4

NC CODE 2901 10 00

EC NO. 601 006 00 1

UN/ID NO. 1265

ADR/RID 3 F1

IMDG 3/I

R: 12-51/53-65-66-67

S: 16-29-33-61-62-9







Substances Darkened by H ₂ SO ₄	passes test
Water (H ₂ O)	max. 0.01%
Fluorescence Trace Impurities (as qui	inine base),
opb:	
Measured at 450 nm	max. 0.5
Measured at Emission Maximum for	
Solvent Impurities	max. 1.5
Physical Data (not specifications):	
Density (g/ml) at 20°C	0.626

Ultraviolet Absorbance (1.00-cm path vs water):	
at 210 nm	max. 0.4
at 220 nm	max. 0.10
at 254 nm	max. 0.01
at 280 nm	max. 0.01
LIV Cut-off_nm	max 190

PRODUCT	PACKING	CONT. BOX
9331.1000	11	BUX
9331.2500	2.5	4

Filtered through a 0.2 micron filter. Packaged under Nitrogen.

n-Pentane

'BAKER ANALYZED' / Ultraviolet Spectrophotometry

8287

► CH₃(CH₂)₃CH₃

M = 72.15 g/mol

1 l = 0.62 kg

FLASHPOINT −40 °C

CAS NO. 109-66-0

EINECS 203-692-4

NC CODE 2901 10 00

EC NO. 601 006 00 1

UN/ID NO. 1265

ADR/RID 3 F1

IMDG 3/I

R: 12-51/53-65-66-67

S: 16-29-33-61-62-9







Water (H ₂ O)	max. 0.01%
Ultraviolet Absorbance (1.00-cm path vs water):	
at 205 nm	max. 1.00
at 210 nm	max. 0.30
at 220 nm	max. 0.10
at 245-400 nm	max. 0.01

PRODUCT NO.	PACKING	CONT. BOX



Penta

n-Pentane 8114 'BAKER ANALYZED

► CH₃(CH₂)₃CH₃

M = 72.15 g/mol

1 l = 0.62 kg

FLASHPOINT $-40\,^{\circ}\mathrm{C}$

CAS NO. 109-66-0 **EINECS** 203-692-4

NC CODE 2901 10 00

EC NO. 601 006 00 1

UN/ID NO. 1265

ADR/RID 3F1 IMDG 3/I

R: 12-51/53-65-66-67

S: 16-29-33-61-62-9







Assay (by GC)	min. 99%
Boiling Range	35-37°C
Density (g/ml) at 20°C	0.624-0.626
Residue after Evaporation	max. 0.001%
Sulfur Compounds (as S)	max. 0.005%
Water (H ₂ O)	max. 0.01%
Trace Impurities (in ppm):	
Aluminium (AI)	max. 0.5
Barium (Ba)	max. 0.1
Boron (B)	max. 0.02
Cadmium (Cd)	max. 0.05
Calcium (Ca)	max. 0.5
Chromium (Cr)	max. 0.02
Cobalt (Co)	max. 0.02
Copper (Cu)	max. 0.02
Iron (Fe)	max. 0.1
Lead (Pb)	max. 0.1
Magnesium (Mg)	max. 0.1
Manganese (Mn)	max. 0.02
Nickel (Ni)	max. 0.02
Tin (Sn)	max. 0.1
Zinc (Zn)	max. 0.1

NO.	вох
8114.1000	11 6
8114.2500	2.5 4
8114.5000	5 EcoTainer
8114.9025	25 l
8114.9200	200 l
EcoTainer, the r	metal solvent can for more safety in the lab.
For safe handlir	ng of 25 l tin cans, see Self-closing tap.

8685

► CH₃(CH₂)₃CH₃

n-Pentane

M = 72.15 g/mol

11= 0.62 kg

FLASHPOINT −40 °C

CAS NO. 109-66-0 **EINECS** 203-692-4

NC CODE 2901 10 00

EC NO. 601 006 00 1

UN/ID NO. 1265

ADR/RID 3 F1

IMDG 3/I

R: 12-51/53-65-66-67

S: 16-29-33-61-62-9









PRODUCT PACKING CONT. BOX 8685.2500 2.5 [8685.9200 200 l

1-Pentanesulfonic Acid Sodium Salt Monohydrate

at 230 nm

at 240 nm

at 250 nm

Assay (by GC)

Residue after Evaporation

'BAKER HPLC ANALYZED

CH₃(CH₂)₆SO₃Na.H₂O

2841

M = 220.26 g/mol

CAS NO. 22767-49-3

EINECS 245-208-4 NC CODE 2904 10 00

For Ion-Pair Chromatography of Basic Compounds min. 98.0% Assay (acidimetric) UV Absorbance of 0.25M Solution: at 200 nm at 210 nm at 220 nm

max. 0.2 max. 0.08 max. 0.06 max. 0.05 max. 0.05 max. 0.05

min. 95%

max. 0.001%

PRODUCT	PACKING	CONT.
NO.		вох
2841.0025	25 g Glass	
2841.0100	100 g	



1-Pentanol

CONT. BOX

Perchloric Acid

70% / ULTREX II Ultrapure Reagent

4806

► HClO₄

M = 100.46 g/mol **1 l** = 1.67 kg

CAS NO. 7601-90-3 **EINECS** 231-512-4

NC CODE 2811 1980 **EC NO.** 17 006 00 4

UN/ID NO. 1873

ADR/RID 5.1 CO1

IMDG 5.1/I

R: 35-5-8 **S:** 23-26-36-45





Assay	65-71%
Trace Impurities (in ppt) (p	g/g):
Aluminium (AI)	max. 100
Antimony (Sb)	max. 100
Arsenic (As)	max. 100
Barium (Ba)	max. 100
Beryllium (Be)	max. 100
Bismuth (Bi)	max. 10
Cadmium (Cd)	max. 10
Calcium (Ca)	max. 100
Cerium (Ce)	max. 10
Cesium (Cs)	max. 10
Cobalt (Co)	max. 100
Copper (Cu)	max. 100
Dysprosium (Dy)	max. 10
Erbium (Er)	max. 10
Europium (Eu)	max. 10
Gadolinium (Gd)	max. 10
Gallium (Ga)	max. 10
Hafnium (Hf)	act. value reported
Holmium (Ho)	max. 10
Indium (In)	max. 10
Iron (Fe)	max. 100
Lanthanum (La)	max. 10
Lead (Pb)	max. 10
Lithium (Li)	max. 100
Lutetium (Lu)	max. 10
Magnesium (Mg)	max. 100
Manganese (Mn)	max. 100
Molybdenum (Mo)	max. 100
Neodymium (Nd)	max. 10
Nickel (Ni)	max. 100
Niobium (Nb)	act. value reported
Palladium (Pd)	max. 10
Platinum (Pt)	max. 100
Potassium (K)	max. 100
Praseodymium (Pr)	max. 10
Rhodium (Rh)	max. 10
Rubidium (Rb)	max. 10
, ,	

Scandium (Sc)	max. 100
Silver (Ag)	max. 100
Sodium (Na)	max. 100
Strontium (Sr)	max. 100
Tantalum (Ta)	act. value reported
Tellurium (Te)	max. 10
Terbium (Tb)	max. 10
Thallium (TI)	max. 100
Thorium (Th)	max. 10
Thulium (Tm)	max. 10
Tin (Sn)	max. 100
Titanium (Ti)	max. 100
Tungsten (W)	act. value reported
Uranium (U)	max. 10
Vanadium (V)	max. 100
Ytterbium (Yb)	max. 10
Yttrium (Y)	max. 10
Zinc (Zn)	max. 100
Zirconium (Zr)	max. 100

PRODUCT

4806.0500

max. 10

PACKING

500 ml

The J.T.Baker CYCLE-TAINER High Purity Solvent Delivery System, preserves purity and protects people.

Samarium (Sm)

See chapter 3 of this catalogue for product details.



Perch

Perchloric Acid

6022 70-72% / 'BAKER ANALYZED' / ACS

M = 100.46 g/mol
11 = 1.67 kg
CAS NO. 7601-90-3
EINECS 231-512-4
NC CODE 2811 19 80
EC NO. 17 006 00 4
UN/ID NO. 1873
ADR/RID 5.1 O C1
IMDG 5.1/I
R: 35-5-8



S: 23-26-36-45



Exceeds ACS Specifications Assay 70.0-72.0% max. 0.0003% Chloride (CI) Color (APHA) max. 10 Free Chlorine max. 5 ppm Insoluble in Ethanol max. 0.001% Nitrogen Compounds (as N) max. 0.001% Residue after Ignition max 0.003% Silicate and Phosphate (as SiO₂) max. 5 ppm max. 0.001% Sulfate (SO₄) Trace Impurities (in ppm): Aluminium (AI) max. 0.05 Arsenic (As) max. 0.05 Barium (Ba) max. 0.02 Beryllium (Be max. 0.01 Cadmium (Cd) max. 0.02 Calcium (Ca) max. 0.5 Chromium (Cr) max. 0.3 max. 0.01 Cobalt (Co) Copper (Cu) max. 0.01 Germanium (Ge) max. 0.05 Heavy Metals (as Pb) max. 1 Iron (Fe) max. 1 Lead (Pb) max. 0.05 Lithium (Li) max. 0.01 Magnesium (Mg) max. 0.1 max. 0.01 Manganese (Mn) max. 0.02 Molybdenum (Mo) Nickel (Ni) max. 0.05 Potassium (K) max. 0.1 Silver (Ag) max. 0.10 Sodium (Na) max. 0.5 max. 0.01 Strontium (Sr) Thallium (TI) max. 0.05 Titanium (Ti) Vanadium (V) max. 0.01 max. 0.05 Zinc (Zn)

Zirconium (Zr)

PRODUCT	PACKING	CONT.
NO.		вох
6022.1000GL	1 Glass	6
6022.2500	2.5	4

Mallinckrodt Baker's cGMP Manufactured Chemicals for the Biopharmaceutical industry are a necessity for uncomplicated scale-up.

max. 0.1

See chapter 6 of this catalogue.



CONT. вох

Perchloric Acid

69-72 % / 'BAKER INSTRA-ANALYZED' / for Trace Metal Analysis / ACS

Meets ACS Specifications

9653

► HCIO₄	
M =	100.46 g/mol
1 =	1.67 kg
CAS NO.	7601-90-3
EINECS	231-512-4
NC CODE	2811 19 80
EC NO.	17 006 00 4
UN/ID NO.	1873
ADR/RID	5.1 CO1
IMDG	5.1/I
R:	35-5-8
S:	23-26-36-45
	C

ivieets ACS Specifications	
Assay	69.0-72.0%
Color (APHA)	max. 10
Manganese (Mn) (by AAS)	max. 0.005 ppm
Residue after Ignition	max. 10 ppm
Trace Impurities (in ppm):	
Aluminium (AI)	max. 0.05
Barium (Ba)	max. 0.1
Cadmium (Cd)	max. 0.005
Calcium (Ca)	max. 1
Chloride (CI)	max. 10
Chromium (Cr)	max. 0.05
Cobalt (Co)	max. 0.005
Copper (Cu)	max. 0.05
Heavy Metals (as Pb)	max. 0.1
Iron (Fe)	max. 0.05
Lead (Pb)	max. 0.01
Lithium (Li)	max. 0.1
Magnesium (Mg)	max. 0.05
Mercury (Hg)	max. 0.005
Nickel (Ni)	max. 0.005
Nitrogen Compounds (as N)	max. 10
Potassium (K)	max. 0.5
Silicate and Phosphate (as SiO ₂)	max. 5
Silicon (Si)	max. 0.1
Silver (Ag)	max. 0.005
Sodium (Na)	max. 1
Strontium (Sr)	max. 0.02
Sulfate (SO ₄)	max. 10

Manganese (Mn) (by AAS)	max. 0.005 ppm
Residue after Ignition	max. 10 ppm
Trace Impurities (in ppm):	
Aluminium (AI)	max. 0.05
Barium (Ba)	max. 0.1
Cadmium (Cd)	max. 0.005
Calcium (Ca)	max. 1
Chloride (CI)	max. 10
Chromium (Cr)	max. 0.05
Cobalt (Co)	max. 0.005
Copper (Cu)	max. 0.05
Heavy Metals (as Pb)	max. 0.1
Iron (Fe)	max. 0.05
Lead (Pb)	max. 0.01
Lithium (Li)	max. 0.1
Magnesium (Mg)	max. 0.05
Mercury (Hg)	max. 0.005
Nickel (Ni)	max. 0.005
Nitrogen Compounds (as N)	max. 10
Potassium (K)	max. 0.5
Silicate and Phosphate (as SiO ₂)	max. 5
Silicon (Si)	max. 0.1
Silver (Ag)	max. 0.005
Sodium (Na)	max. 1
Strontium (Sr)	max. 0.02
Sulfate (SO ₄)	max. 10
Tin (Sn)	max. 0.01
Zinc (Zn)	max. 0.05

PRODUCT

NO. 9653.0500 PACKING

500 ml

Use J.T.Baker Ultrex II and BAKER INSTRA-ANALYZED acids for low level trace element analysis.

See chapter 3 of this catalogue for more details.



Perch

Perchloric Acid

6023 60-62% / 'BAKER ANALYZED' / ACS

► HCIO₄

1 l = 1.54 kg **CAS NO.** 7601-90-3 **EINECS** 231-512-4

M = 100.46 g/mol

NC CODE 2811 19 80

EC NO. 17 006 00 4 **UN/ID NO.** 1873

ADR/RID 5.1 CO1

IMDG 5.1/I R: 35-5-8

S: 23-26-36-45





Exceeds ACS Specifications	
Assay	60.0-62.0%
Chloride (CI)	max. 0.001%
Color (APHA)	max. 10
Nitrogen Compounds (as N)	max. 0.001%
Residue after Ignition	max. 0.003%
Silicate and Phosphate (as SiO ₂)	max. 5 ppm
Sulfate (SO ₄)	max. 0.001%
Trace Impurities (in ppm):	
Aluminium (AI)	max. 0.05
Arsenic (As)	max. 0.05
Barium (Ba)	max. 0.02
Beryllium (Be)	max. 0.01
Cadmium (Cd)	max. 0.02
Calcium (Ca)	max. 0.5
Chromium (Cr)	max. 0.3
Cobalt (Co)	max. 0.01
Copper (Cu)	max. 0.01
Germanium (Ge)	max. 0.05
Heavy Metals (as Pb)	max. 1
Iron (Fe)	max. 1
Lead (Pb)	max. 0.05
Lithium (Li)	max. 0.01
Magnesium (Mg)	max. 0.1
Manganese (Mn)	max. 0.01
Molybdenum (Mo)	max. 0.02
Nickel (Ni)	max. 0.05
Potassium (K)	max. 0.1
Silver (Ag)	max. 0.10
Sodium (Na)	max. 0.5
Strontium (Sr)	max. 0.01
Thallium (TI)	max. 0.05
Titanium (Ti)	max. 0.1
Vanadium (V)	max. 0.01
Zinc (Zn)	max. 0.05
Zirconium (Zr)	max. 0.1

Perchloric acid in anhydrous acetic acid

221 0.5 mol/l / 'BAKER ANALYZED

NC CODE 3822 00 00 UN/ID NO. 2789 ADR/RID 8 CF1 IMDG 8/II

R: 10-35 **S:** 23-26

corrosive

Titer (mol/l) 0.49-0.51

 PRODUCT
 PACKING
 CONT.

 NO.
 BOX

 7221.1000
 1 |

Volumetric Solution, ready for use.

Perchloric acid in anhydrous acetic acid

Titer (mol/l)

0.1 mol/I / 'BAKER ANALYZED'

NC CODE 3822 00 00 UN/ID NO. 2789 8 CF1 IMDG 8/II R: 10-35

7142

S: 23-26



0.098-0.102 **PRODUCT NO.**7142.1000

PRODUCT

6023.1000

6023.2500

NO.

PACKING

11

2.5 l

CONT.

вох

6

 DUCT
 PACKING
 CONT.

 8.00
 BOX

 1.1000
 1 l
 6

Volumetric Solution, ready for use.



Perchloroethylene

See Tetrachloroethylene

1366

9272

9265

Periodic Acid

'BAKER ANALYZED

► H₅IO₆ **M** = 227.94 g/mol

CAS NO. 10450-60-9

EINECS 233-937-0

NC CODE 2811 1980

UN/ID NO. 3085

ADR/RID 5.1 OC2

IMDG 5.1/I

R: 34-8

S: 26-36/37/39-45





Assay	min. 99.5%
Insoluble Matter	max. 0.005%
lodide (I)	max. 0.001%
Residue after Ignition	max. 0.5%
Sulfate (SO ₄)	max. 0.01%

PRODUCT NO.	PACKING	CONT. BOX
1366.0025	25 g Glass	
1366.0100	100 g	

Petroleum Ether

20-40°C / 'BAKER ANALYZED'

1 l = 0.59 kg FLASHPOINT −50 °C

CAS NO. 8032-32-4

EINECS 232-453-7

NC CODE 2710 11 21

UN/ID NO. 1268

ADR/RID 3 F1

IMDG 3/II

R: 12-51/53-65 **S:** 16-29-33-9









Boiling Range 20-40°C Heavy Oils and Fats passes test Odor (faint) passes test Residue after Evaporation max. 0.002%

passes test

passes test

30-60°C

PRODUCT

9265.1000

9265,2500

251

NO.

PRODUCT PACKING CONT. вох NO. 9272.40005 4 | Safetainer

Petroleum Ether

30-60°C / 'BAKER ULTRA RESI-ANALYZED' / for Organic Residue Analysis

Acidity

Appearance and Color

PACKING CONT. BOX

FLASHPOINT −57 °C **CAS NO.** 8032-32-4 EINECS 232-453-7 NC CODE 2710 11 21

1 l = 0.64 kg

UN/ID NO. 1268

ADR/RID 3 F1

IMDG 3/II

R: 11-38-51/53-65-67

S: 16-23-24-29-33-61-62-9









Color (APHA) max. 10 Residue after Evaporation max. 5 ppm Water (by KF, coulometric) max. 0.05% ECD Sensitive Impurities (as Heptachlor Epoxide): Single Impurity Peak (pg/ml) FID-Sensitive Impurities (as 2-Octanol): Single Impurity Peak (ng/ml) max. 5 Neat solvent front characterization: ECD-Sensitive Impurities (as Ethylene Dibromide):

Boiling Range (initial to dry point)

Single Impurity Peak (ng/ml) max. 5

Petroleum Ether

30-60°C / 'BAKER ANALYZED' / For Hydrocarbon Oil Index determination

9266

CONT.

BOX

FLASHPOINT $-57~^{\circ}\text{C}$ CAS NO. 8032-32-4 **EINECS** 232-453-7 NC CODE 2710 1121 **UN/ID NO.** 1268

11 = 0.64 kg

ADR/RID 3 F1 IMDG 3/II

R: 11-38-51/53-65-67

S: 16-23-24-29-33-61-62-9









Color (APHA) max. 10 Residue after Evaporation max. 5 ppm Water (H₂O) max. 0.05% Hydrocarbon oil index concentration (as RIVM oil reference standard): Total peaks between n-decane $(C_{10}H_{22})$ and n-tetracontane (C₄₀H₈₂)

max. 0.5 mg/l

PRODUCT

9266,2500

NO.

Suitable for determination of Hydrocarbon Oil Index according to ISO 9377-2 and NEN 5733.

PACKING

2.5 |

Petro

Petroleum Ether

35-60°C / 'BAKER ANALYZED' / Ultraviolet Spectrophotometry

Acidity

Boiling Range

1 I = 0.64 kg
FLASHPOINT -57 °C
CAS NO. 8032-32-4
EINECS 232-453-7
NC CODE 2710 11 21
UN/ID NO. 1268
ADR/RID 3 F1
IMDG 3/II

R: 11-38-51/53-65-67 S: 16-23-24-29-33-61-62-9

N harmful

 Color (APHA)
 max. 10

 Residue after Evaporation
 max. 0.001%

 Ultraviolet Absorbance (1.00-cm path vs water):
 at 220 nm
 max. 1.0

 at 230 nm
 max. 0.20

 at 250 nm
 max. 0.05

 at 270 nm
 max. 0.01

 at 400 nm
 max. 0.01

passes test

35-60°C

PRODUCT	PACKING	CONT.
NO.		вох
9270.1000	11	
9270.4000	4 Glass	

Petroleum Ether

15 40-65°C / 'BAKER ANALYZED

1 I = 0.64 kg

FLASHPOINT −60 °C

CAS NO. 8032-32-4

EINECS 232-453-7

NC CODE 2710 11 21

UN/ID NO. 1268

ADR/RID 3 F1

IMMOG 3/II

R: 11-38-51/53-65-67

S: 16-23-24-29-33-61-62-9

dangerous for the

environment





highly

Boiling Range Residue after Evaporation max. 0.001% max. 0.02% Water (H₂O) Trace Impurities (in ppm): Aluminium (AI) max. 0.5 Barium (Ba) max. 0.1 Boron (B) max. 0.02 Cadmium (Cd) max. 0.05 Calcium (Ca) max. 0.5 Chromium (Cr) max. 0.02 Cobalt (Co) max. 0.02 max. 0.02 Copper (Cu) Iron (Fe) max. 0.1 Lead (Pb) max. 0.1 Magnesium (Mg) max. 0.1 Manganese (Mn) max. 0.02 Nickel (Ni) max. 0.02 Tin (Sn) max. 0.1 Zinc (Zn) max. 0.1

PRODUCT	PACKING CONT
NO.	B0
8115.1000	11
8115.5000	5 EcoTainer
8115.9025	25
8115.9030RC	30 Returnable Container
8115.9200	200
EcoTainer, the	metal solvent can for more safety in the lab
For safe hand	ling of 25 l tin cans, see Self-closing tap.
	, , , , , , , , , , , , , , , , , , , ,

140

Petroleum Ether

10-65°C / 'BAKER'

1 I = 0.64 kg FLASHPOINT −60 °C CAS NO. 8032-32-4 EINECS 232-453-7 NC CODE 2710 11 21 UN/ID NO. 1268

ADR/RID 3 F1
IMDG 3/II

R: 11-38-51/53-65-67

S: 16-23-24-29-33-61-62-9

dangerous





Appearance	passes test
Residue after Evaporation	max. 0.002%
Appearance at 0°C	passes test
Distillation Range	40-65°C

PRODUCT	PACKING C	ONT.
NO.		вох
8240.1000	11	
8240.5000	5 EcoTainer	
8240.9025	25 l	
8240.9200	200 l	
EcoTainer, th	e metal solvent can for more safety in the	lab.
For safe han	dling of 25 l tin cans, see Self-closing tap.	

www.jtbaker.com/europe

8116

8241

8242

Petroleum Ether

60-100°C / 'BAKER ANALYZED'

1 l =	0.67 kg
FLASHPOINT	−32 °C
CAS NO.	8032-32-4
EINECS	232-453-7
NC CODE	2710 11 21
UN/ID NO.	1268

UN/ID NO. 1268 ADR/RID 3 F1 IMDG 3/II

R: 11-38-51/53-65-67

S: 16-23-24-29-33-57-60-9







Boiling Range	60-100°C
Residue after Evaporation	max. 0.001%
Water (H ₂ O)	max. 0.02%
Trace Impurities (in ppm):	
Aluminium (AI)	max. 0.5
Barium (Ba)	max. 0.1
Boron (B)	max. 0.02
Cadmium (Cd)	max. 0.05
Calcium (Ca)	max. 0.5
Chromium (Cr)	max. 0.02
Cobalt (Co)	max. 0.02
Copper (Cu)	max. 0.02
Iron (Fe)	max. 0.1
Lead (Pb)	max. 0.1
Magnesium (Mg)	max. 0.1
Manganese (Mn)	max. 0.02
Nickel (Ni)	max. 0.02
Tin (Sn)	max. 0.1
Zinc (Zn)	max. 0.1

PRODUCT	PACKING CO	ONT.
NO.		BOX
8116.1000	11	
8116.5000	5 EcoTainer	
8116.9025	25 l	
EcoTainer, the	e metal solvent can for more safety in the	lab.

For safe handling of 25 l tin cans, see Self-closing tap.

Petroleum Ether

60-100°C / 'BAKER'

1 = 0.67 kgFLASHPOINT $-32 ^{\circ}\text{C}$ CAS NO. 8032-32-4

EINECS 232-453-7 **NC CODE** 2710 11 21

UN/ID NO. 1268
ADR/RID 3 F1
IMDG 3/II

R: 11-38-51/53-65-67

S: 16-23-24-29-33-57-60-9









Boiling Range

Residue after Evaporation

Boiling Range 60-100°C
Residue after Evaporation max. 0.002%

 PRODUCT NO.
 PACKING SOME SEASON
 CONT. BOX

 8241.9025
 25 l

 8241.9200
 200 l

For safe handling of 25 l tin cans, see Self-closing tap.

Petroleum Ether

80-110°C / 'BAKER'

1 I = 0.70 kgFLASHPOINT $-5 \,^{\circ}\text{C}$

CAS NO. 8032-32-4

EINECS 232-453-7 **NC CODE** 2710 11 21

UN/ID NO. 1268 ADR/RID 3 F1

IMDG 3/II

R: 11-38-51/53-65-67

S: 16-23-24-29-33-57-60-61-9



dangerous for the environment





80-110°C PRODUCT PACKING CONT.

max. 0.002% No. BOX
8242.1000 1 |
8242.9025 25 |

Innovation is principal to our business.



Petro

Petroleum Ether

8243 100-140°C / 'BAKER'

1 l = 0.73 kg FLASHPOINT 0°C **CAS NO.** 8032-32-4

EINECS 232-453-7 NC CODE 2710 1121 **UN/ID NO.** 1268

ADR/RID 3F1 IMDG 3/II

R: 11-38-51/53-65-67

S: 16-23-24-29-33-57-60-9



CONT. **PRODUCT PACKING** NO. вох 8243.1000 11 8243.5000 5 | EcoTainer 8243.9025 25 I 8243.9200 EcoTainer, the metal solvent can for more safety in the lab. For safe handling of 25 l tin cans, see Self-closing tap.

100-140°C

PRODUCT

1158.0005

1158.0025

passes test

PACKING

25 g Glass

CONT.

вох

max. 0.002%



Pharmaceutical Products

Boiling Range

Residue after Evaporation

Meets ACS Specifications

Suitability as Redox Indicator

Suitability for determining Iron

1158

1,10-Phenanthroline monohydrate

'BAKER ANALYZED' / ACS

C₁₂H₈N₂.H₂O

M = 198.23 g/mol

CAS NO. 66-71-7 **EINECS** 200-629-2 NC CODE 2933 99 90

EC NO. 613 092 00 8

UN/ID NO. 2811 **ADR/RID** 6.1 T2

IMDG 6.1/III

R: 25-50/53 **S:** 45-60-61

dangerous for the

environment









1,10-Phenanthroline-Iron(II) Salt Solution

1368

1,10-Phenanthrolinium Chloride 1aq

► C₁₂H₈N₂HCl.H₂O

M = 234.69 g/mol CAS NO. 3829-86-5 **EINECS** 223-325-1

NC CODE 2933 99 90

UN/ID NO. 2811 **ADR/RID** 6.1 T2

IMDG 6.1/III

R: 25 **S:** 45



Assay	min. 99%
Residue after Ignition	max. 0.05%
Suitability as Redox Indicator	passes test
Suitability for determining Iron	passes test
Water (H ₂ O)	about 9%

PRODUCT	PACKING	CONT.
NO.		вох
1368.0005	5 g	
1368.0025	25 g Glass	

www.jtbaker.com/europe

Phenol

white fused crystal / 'BAKER ULTRAPURE BIOREAGENT'

4056

► C ₆ H ₅ OH	
M =	94.11 g/mol
FLASHPOINT	79 °C
CAS NO.	108-95-2
EINECS	203-632-7
NC CODE	2907 11 00
EC NO.	604 001 00 2
UN/ID NO.	1671

ADR/RID 6.1 T2 **IMDG** 6.1/II **R:** 24/25-34 **S:** 28-45

Contains no preservative		
Suitable for Extraction of Nu	cleic Acids:	
Appearance	passes test	
Assay	min. 99%	
DNase Activity	none detected	
Heavy Metals (as Pb)	max. 0.001%	
Melting Point	39-41°C.	
pH (saturated)	4.5-6.0	
RNase Activity	none detected	
Water (H ₂ O)	max. 0.4%	

PRODUCT	PACKING	CONT.
NO.		вох
4056.0100	100 g	
4056.0500	500 g	

Phenol

crystal / 'BAKER ANALYZED' / ACS

1122

0505

EC NO. 604 001 00 2 UN/ID NO. 1671 ADR/RID 6.1 T2 IMDG 6.1/II

R: 23/24/25-34-48/20/21/22-68

S: 24/25-26-28C-36/37/39-45





Exceeds ACS Specifications	
Assay (by GC)	min. 99.0%
Clarity of Solution	passes test
Freezing Point	min. 40.5°C, dry basis
Residue after Evaporation	max. 0.03%
Water (H ₂ O)	max. 0.2%

PRODUCT NO.	PACKING	CONT. BOX
0188.0500	500 g	
0188.1000	1 kg	6

Phenol

crystal / 'BAKER'

► C_6H_5OH **M** = 94.11 g/mol **FLASHPOINT** 79 °C

> **CAS NO.** 108-95-2 **EINECS** 203-632-7 **NC CODE** 2907 11 00

EC NO. 604 001 00 2 **UN/ID NO.** 1671

ADR/RID 6.1 T2 **IMDG** 6.1/II

R: 23/24/25-34-48/20/21/22-68

S: 24/25-26-28C-36/37/39-45





Assay	99.0-100.5%
Acidity	passes test
Appearance of solution	passes test
Freezing Point	min. 39.5°C
dentification	passes test
Residue on Evaporation	max. 0.05%

PRODUCT	DUCT PACKING	CONT.
NO.		BOX
0505.1000	1 kg	6
0505.9025	25 kg	

Stored in an airtight container. Stored protected from light.

Mallinckrodt Baker's chemistry is Part of a pure process TM .

Pheno

Phenol

1480 quefied / 'BAKER ANALYZED'

► C₆H₅OH

M = 94.11 g/mol

1 l = 1.10 kg

FLASHPOINT $79\,^{\circ}\mathrm{C}$

CAS NO. 108-95-2

EINECS 203-632-7

NC CODE 2907 11 00

EC NO. 604 001 00 2

UN/ID NO. 2821

ADR/RID 6.1 T2

IMDG 6.1/II

R: 23/24/25-34-48/20/21/22-68

S: 24/25-26-28-36/37/39-45







Phenol

e for qualities meeting the requirements of cGMP, USP/NF, BP/Ph.Eur, JP or FCC, page 36

Assay (C₆H₅OH)

Residue after Evaporation

Suitability for Amino Acid Determination

Solubility and Reaction

Appearance

Water (H₂O)

1159

Phenolphthalein

'BAKER ANALYZED' / ACS

► Ç₆H₄COOÇ(C₆H₄-OH)₂

M = 318.33 g/mol

CAS NO. 77-09-8 **EINECS** 201-004-7

NC CODE 2932 29 10

Meets ACS Specifications Clarity of alcohol solution passes test Visual Transition Interval: pH 10.0 red pH 8.0 colorless **PRODUCT PACKING** CONT. вох 1159.0100 100 g

CONT.

CONT.

BOX

6

вох



7143

Phenolphthalein

See for qualities meeting the requirements of cGMP, USP/NF, BP/Ph.Eur, JP or FCC, page 36

Phenolphtalein

1% solution in ethanol / 'BAKER ANALYZED'

11 = 0.80 kg

FLASHPOINT 11 °C

NC CODE 2207 10 00

UN/ID NO. 1170

ADR/RID 3F1

IMDG 3/II

R: 11

S: 16-7



Visual Transition Interval: pH 10.0

pH 8.0

red colorless

min. 89.0%

passes test

max. 0.05%

passes test

passes test

act. value reported

PRODUCT

1480.0500

NO.

PACKING

500 ml

Contains about 20 ppm Citric Acid as a Preservative.

Preserve in tight, light-resistant glass containers.

NO. 7143.0250 250 ml Excluding excise.

PACKING

Solution, ready for use.

PRODUCT



L-Phenylalanine

See for qualities meeting the requirements of cGMP, USP/NF, BP/Ph.Eur, JP or FCC, page 36



N-Phenylaniline



(Phenylazo) formic Acid 2-Phenylhydrazide



Phenylmethanol

See Benzyl Alcohol



Phloxine B

'BAKER ANALYZED' / Certified Stain

1468

CAS NO. 18472-87-2 **EINECS** 242-355-6 NC CODE 3204 12 00

nl	
in 50% C_2H_5OH and 0.01% Na_2CO_3 , 1-cm	
act. value reported	
act. value reported	
passes test	
act. value reported	

PRODUCT NO.	PACKING	CONT.
		вох
1468.0025	25 g Glass	

Certified by the Biological Stain Commission.

Phosphate Reagent Solution

7491

NC CODE 3822 00 00

PRODUCT	PACKING	CONT.
NO.		вох
7491.9020	20 l Polycube	

Volumetric Solution, ready for use. Contains: di-Sodiumhydrogen Phosphate Dodecahydrate, tri-Sodiumphosphate Dodecahydrate, Water (max. $1\mu S$).

meta-Phosphoric Acid

0676

CAS NO.	37267-86-0
EINECS	253-433-4
NC CODE	2809 20 00
UN/ID NO.	3260
ADR/RID	8 C2
IMDG	8/111
R:	34
S:	26-36/37/39-45

Meets ACS Specifications	
Assay (HPO ₃)	33.5-36.5%
Chloride (CI)	max. 0.001%
Heavy Metals (as Pb)	max. 0.005%
Iron (Fe)	max. 0.005%
Nitrate (NO ₃)	max. 0.001%
Stabiliser (NaPO ₃)	57.0-63.0%
Substances Reducing KMnO ₄ (as H ₃ PO ₄)	max. 0.02%
Sulfate (SO ₄)	max. 0.005%
Trace Impurities (in ppm):	
Arsenic (As)	max. 1

PRODUCT	PACKING	CONT
NO.		BO)
0676.0500	500 g	
0676.2500	2.5 kg	

Vitreous Sodium Acid Metaphosphate. This product contains as a stabilizer a somewhat greater proportion of sodiummetaphosphate than that corresponding to the formula NaH(PO₃)₂

meta-Phosphoric Acid

CAS NO. 37267-86-0 **EINECS** 253-433-4

NC CODE 2809 20 00 **UN/ID NO.** 3260 ADR/RID 8 C2

> IMDG 8/III **R:** 34

> > **S:** 26-36/37/39-45



Assay of (HPO ₃)n	40-44%
Assay of (NaPO ₃)n	56-60%
Arsenic (As)	max. 5 ppm
Chloride (CI)	max. 0.002%
Heavy Metals (as Pb)	max. 0.001%
Iron (Fe)	max. 0.01%

вох
6

Phosphoric Acid

ualities meeting the requirements of cGMP, USP/NF, BP/Ph.Eur, JP or FCC, page 36



0190

Questions or suggestions, please contact us at jtbaker.nl@emea.tycohealthcare.com



Phosp

Phosphoric Acid

6908 85% / ULTREX II Ultrapure Reagent

Certificate Frovided Heporting A	Ctuar Lot Ariarysis
Actual Lot Analysis Lot No. Y095	84
Assay	86.5% (w/w)
Appearance	passes test
Density (g/ml) at 25°C	1.69
Trace Impurities (in ppb):	
Aluminium (AI)	5
Antimony (Sb)	< 0.1
Barium (Ba)	< 0.3
Beryllium (Be)	< 0.1
Bismuth (Bi)	< 0.1
Cadmium (Cd)	< 0.1
Calcium (Ca)	< 5
Cerium (Ce)	< 0.1
Cesium (Cs)	< 0.1
Chromium (Cr)	1
Cobalt (Co)	< 0.1
Copper (Cu)	< 5
Dysprosium (Dy)	< 0.1
Erbium (Er)	< 0.1
Europium (Eu)	< 0.1
Gadolinium (Gd)	< 0.1
Gallium (Ga)	< 0.2
Germanium (Ge)	< 0.1
Gold (Au)	< 0.1
Holmium (Ho)	< 0.1
Indium (In)	< 0.1
Iridium (Ir)	< 0.1
Iron (Fe)	< 0.1
Lanthanum (La)	< 0.1
Lead (Pb)	< 0.1
Lithium (Li)	< 0.1
Lutetium (Lu)	< 0.1
Magnesium (Mg)	3
Manganese (Mn)	< 0.1
Mercury (Hg)	0.04
Molybdenum (Mo)	< 0.1
Neodymium (Nd)	< 0.1
Nickel (Ni)	< 1
Niobium (Nb)	< 0.1
Palladium (Pd)	< 0.1

Certificate Provided Reporting Actual Lot Analysis

Platinum (Pt)		< 0.1
Potassium (K)	< 5
Praseodymiu	m (Pr)	< 0.1
Rhodium (Rh)	< 0.1
Rubidium (Rb	o)	< 0.1
Ruthenium (F	Ru)	< 0.1
Samarium (S	m)	< 0.1
Scandium (S	c)	< 0.1
Silver (Ag)		1.4
Sodium (Na)		< 5
Strontium (Sr	·)	< 0.1
Tantalum (Ta)	< 0.1
Terbium (Tb)		< 0.1
Thallium (TI)		< 0.1
Thorium (Th)		< 0.1
Thulium (Tm)		< 0.1
Tin (Sn)		< 0.1
Titanium (Ti)		< 1
Tungsten (W		< 0.1
Uranium (U)		< 0.1
Vanadium (V		< 0.1
Ytterbium (Yt	o)	< 0.1
Yttrium (Y)		< 0.1
Zinc (Zn)		< 5
Zirconium (Zi	-)	< 0.1
Trace Impui	ities (in ppm):	
Arsenic (As)		0.002
Boron (B)		< 0.0005
Nitrate (NO ₃)		2
Silicon (Si)		< 0.02
Sulfate (SO ₄)		12
PRODUCT	PACKING	CONT.
NO.		вох
6908.0050	50 q	

Calibrate and standardise your analytical methods and equipment with J.T.Baker Volumetric and Buffer solutions.

Refer to the Analytical applications section of this catalogue for more details.



Phosphoric Acid 85% / 'BAKER ANALYZED' / ACS

6024

	00 % / DAKE
► H ₃ PO ₄	
M =	98.00 g/mol
1 =	1.71 kg
CAS NO.	7664-38-2
EINECS	231-633-2
NC CODE	2809 20 00
EC NO.	15 011 00 6
UN/ID NO.	1805
ADR/RID	8 C2
IMDG	8/111
R:	34
S:	26-45
	C corrosive

Exceeds ACS Specifications	
Assay	min. 85.0%
Alkali phosphates	passes test
Chloride (CI)	max. 1 ppm
Color (APHA)	max. 10
Heavy Metals (as Pb)	max. 0.001%
Insoluble Matter	max. 0.001%
Nitrate (NO ₃)	max. 5 ppm
Phosphorous or Hypophosphorous Acid	passes test
Reducing Substances	passes test
Sulfate (SO ₄)	max. 0.003%
Volatile Acids (as CH ₃ COOH)	max. 0.001%
Trace Impurities (in ppm):	
Aluminium (AI)	max. 10
Antimony (Sb)	max. 4
Arsenic (As)	max. 0.5
Cadmium (Cd)	max. 5
Calcium (Ca)	max. 20
Copper (Cu)	max. 5
Iron (Fe)	max. 10
Lead (Pb)	max. 5
Magnesium (Mg)	max. 10
Manganese (Mn)	max. 0.5
Nickel (Ni)	max. 5
Potassium (K)	max. 5
Sodium (Na)	max. 10
Zinc (Zn)	max. 10

PRODUCT	PACKING	CONT.
NO.		вох
6024.1000	11	6
6024.2500	2.5	4
6024.9025	25 l	
6024.9200	200 l	

Phosphoric Acid 85% MOS, VLSI Grade

See for detailed information section Microelectronic materials page 32 or the website www.jtbaker.com/micro



Phosphoric Acid Solutions

See for qualities meeting the requirements of cGMP, USP/NF, BP/Ph.Eur, JP or FCC, page 36



Phosphorus, red

Assay



CONT.

вох

PRODUCT

9008.0250

min. 97%

PACKING

250 g

▶ P **M** = 30.97 g/mol **CAS NO.** 7723-14-0 **EINECS** 231-768-7 NC CODE 2804 70 00 **EC NO.** 15 002 00 7 **UN/ID NO.** 1338 **ADR/RID** 4.1 F3 IMDG 4.1/III **R:** 11-16-52/53 **S:** 43d-61-7



Phosphorus(V) Oxide See Phosphorus Pentoxide



Phosphotungstic Acid

See Tungstophosphoric Acid



Photoresist Strippers





Phosp

Phosphorus Pentoxide Powder / 'BAKER ANALYZED' / ACS

0192

▶ P₂O₅ **M** = 141.94 g/mol **CAS NO.** 1314-56-3 **EINECS** 215-236-1 NC CODE 2809 10 00 **EC NO.** 15 010 00 0 **UN/ID NO.** 1807

ADR/RID 8 C2 IMDG 8/II **R:** 35 **S:** 22-26-45

Exceeds ACS Specifications	
Assay	min. 98.0%
Ammonium (NH ₄)	max. 0.01%
Heavy Metals (as Pb)	max. 0.01%
Insoluble Matter	max. 0.01%
Phosphorus Trioxide (P ₂ O ₃)	passes test

00 g	BO)
00 g	

Phthalic Anhydride

0194 'BAKER ANALYZED' / ACS

► Ç₆H₄-1,2-COOÇO **M** = 148.12 g/mol **CAS NO.** 85-44-9 **EINECS** 201-607-5 NC CODE 2917 35 00 **EC NO.** 607 009 00 4 **UN/ID NO.** 2214 ADR/RID 8 C4 IMDG 8/III

> **R:** 22-37/38-41-42/43 **S:** 23-24/25-26-37/39-46

99.7-100.2% Assay white flaky crystals Appearance Chloride (CI) max. 0.001% Melting Range (includes 131°C) max. 3°C Residue after Ignition max. 0.005% Sulfate (SO₄) max. 0.003% Trace Impurities (in ppm): Heavy Metals (as Pb) max. 5 Iron (Fe) max. 5

Exceeds ACS Specifications

DOV		
ВОХ		NO.
	1 kg	0194.1000
	1 kg	194.1000

Phthalic Anhydride in Pyridine

160 g/l / 'BAKER ANALYZED

FLASHPOINT 17 °C NC CODE 3822 00 00 UN/ID NO. 2924 ADR/RID 3 FC IMDG 3/II


7526

R: 11-20/21/22-41-42/43

S: 23-26-36/37/39-45







Volumetric Solution, ready for use.



Pilocarpine Nitrate

See for qualities meeting the requirements of cGMP, USP/NF, BP/Ph.Eur, JP or FCC, page 36

Color (APHA)

Piperazine-N,N'-bis (2-ethanesulfonic acid)

4265 BAKER ULTRAPURE BIOREAGENT

C₈H₁₈N₂O₆S₂ M = 302.37 g/mol **CAS NO.** 5625-37-6 NC CODE 2933 59 95

Assay min. 99.0% Appearance passes test DNase Activity none detected Heavy Metals (as Pb) max. 5 ppm Insoluble Matter max. 0.005% Protease Activity none detected Residue after Ignition max. 0.3% RNase Activity none detected Product Information (not specifications): pKa at 20°C 6.80

NO.	
	BO
4265.0025 25 g	Glass
4265.0100 100	g
4265.1000 1 kg	




Piperazine-N,N'-bis(2-ethanesulfonic acid) Sodium Salt 'BAKER ULTRAPURE BIOREAGENT'

4266



M = 335.85 g/mol **CAS NO.** 100037-69-2 NC CODE 2933 59 95



PRODUCT	PACKING	CONT.
NO.		вох
4266.0025	25 g Glass	
4266.0100	100 g	
4266.1000	1 kg	

PIPES

See Piperazine-N,N'-bis (2-ethanesulfonic acid)



PIPES Sodium Salt



Platinum 1000 μg/ml

(Matrix: 5% hydrochloric acid) / 'BAKER INSTRA-ANALYZED' / Plasma Standard



► Pt

M = 195.08 g/mol NC CODE 3822 00 00

Certificate Provided Reporting Actual Lot Analysis 998-1002 μg/ml Platinum (Pt)

PRODUCT	PACKING	CONT.
NO.		вох
5773.0100	100 ml	

Prepared from the highest purity raw material available, generally greater than 99.999% spectral purity. The content of the solution is confirmed to be accurate to within \pm 0.2%. Typically 1000 μ g/ml. The certificate of analysis provided reports actual lot analysis. The certificate also lists the trace impurities.

Platinum 1000 μg/ml

BAKER INSTRA-ANALYZED

6951

M = 195.08 g/mol NC CODE 3822 00 00

Platinum (Pt)	998-1002 μg/ml

PRODUCT	PACKING	CONT.
NO.		вох
6951.0100	100 ml	
6951.0500	500 ml	

Prepared by dissolution of high purity raw materials (min. 99.99% spectral purity). Assays are verified by ICP against standards traceable to NIST. Standard Reference Material numbers (SRM) are printed on each label.

Platinum 10000 µg/ml

(Matrix: 10% hydrochloric acid) / 'BAKER INSTRA-ANALYZED' / Plasma Standard

5740

► Pt

M = 195.08 g/mol NC CODE 3822 00 00

R: 36/38

S: 26

Certificate Provided Re	eporting Actual Lot Analysis
Platinum (Pt)	9980-10020 ug/m

PRODUCT	PACKING	CONT.
NO.		вох
5740 0100	100 ml	

Prepared from the highest purity raw material available, generally greater than 99.999% spectral purity. The content of the solution is confirmed to be accurate to within 0.2 %. Typically 10000 μ g/ml.

The certificate of analysis provided reports actual lot analysis. The certificate also lists the trace impurities.

Platinum Chloride





Polya

5014

Polyamide 6
'BAKER-FLEX' / Flexible TLC Sheets, 20 X 20 cm

NC CODE 3908 10 00

PRODUCT **PACKING** CONT. NO. вох 5014 25 sheets

A flexible sheet coated with high purity polyamide 6

Polyamide 6 5046

'BAKER TLC' / for Thin Layer Chromatography

NC CODE 3908 10 00

Particle Size (<50 μ) min. 80% Suitability for TLC passes test PRODUCT PACKING CONT. вох NO. 5046.0100 100 g

Polyamide 6-F

5015 'BAKER-FLEX' / Flexible TLC Sheets, 20 X 20 cm

NC CODE 3908 10 00

PRODUCT **PACKING** CONT. вох 5015 25 sheets

A flexible sheet coated with high purity polyamide 6 powder containing a fluorscent indicator (activated at 2540 Å).

Polyethylene Glycol 7176

► H(OCH₂CH₂)_nOH

1 l = 1.12 kg CAS NO. 25322-68-3

EINECS 203-473-3

NC CODE 3907 20 11

Identity (by IR) passes test

PRODUCT NO.

PACKING 7176,1000 11

CONT. BOX

7178

Polyethylene Glycol

► H(OCH₂CH₂)_nOH **1 l** = 1.12 kg

CAS NO. 25322-68-3

EINECS 203-473-3

NC CODE 3907 20 11

4-8°C

PRODUCT

7178.1000

1613.1000

PACKING 11

CONT. вох

1613

Polyethylene Glycol

1000 / 'BAKER

► H(OCH₂CH₂)_nOH

EINECS 203-473-3

NC CODE 3907 20 11

Freezing Point

Freezing Point

37 - 40°C

PACKING PRODUCT NO.

1 kg

CONT. вох

1615

Polyethylene Glycol

► H(OCH₂CH₂)_nOH

EINECS 203-473-3 NC CODE 3907 20 11 Melting Point

55 - 60°C

PRODUCT PACKING NO. 1615.1000 1 kg

CONT. вох

Polyethylene Glycol

► H(OCH₂CH₂)_nOH

EINECS 203-473-3 NC CODE 3907 20 11 Melting Point

60 - 66°C

PRODUCT

1616.1000

PACKING CONT. вох 1 kg

Polypropylene Glycol

7179

► HOCH₂[CH(CH₃)OCH₂]_nCHOHCH₃

1 l = 1.01 kg **CAS NO.** 25322-69-4 **EINECS** 200-338-0 NC CODE 3907 20 29

PRODUCT PACKING CONT. NO. вох 7179.1000 11

Polysorbate 20

See for qualities meeting the requirements of cGMP, USP/NF, BP/Ph.Eur, JP or FCC, page 36



Polysorbate 80

ee for qualities meeting the requirements of cGMP, USP/NF, BP/Ph.Eur, JP or FCC, page 36



Potassium 1000 µg/ml (Matrix: 2% nitric acid) / 'BAKER INSTRA-ANALYZED' / Plasma Standard



▶ K

M = 39.10 g/mol

NC CODE 3822 00 00 R: 36/38

S: 26



Certificate Provided Reporting Actual Lot Analysis Potassium (K)

998-1002 μg/ml

PRODUCT	PACKING	CONT.
NO.		вох
5774.0100	100 ml	

Prepared from the highest purity raw material available, generally greater than 99.999% spectral purity. The content of the solution is confirmed to be accurate to within \pm 0.2%. Typically 1000 μ g/ml.

The certificate of analysis provided reports actual lot analysis. The certificate also lists the trace impurities.

Potassium 1000 μg/ml (Matrix: 2% nitric acid) / 'BAKER INSTRA-ANALYZED' / Atomic Absorption Standard

6950

CONT.

▶ K

M = 39.10 g/mol NC CODE 3822 00 00

R: 36/38

S: 26



Potassium (K) 998-1002 μg/ml

NO. 6950.0100 100 ml

PACKING

PRODUCT

PRODUCT

PRODUCT

6950.0500 500 ml

Prepared by dissolution of high purity raw materials (min. 99.99% spectral purity). Assays are verified by ICP against standards traceable to NIST. Standard Reference Material numbers (SRM) are printed on each label.

Potassium 1000 μ g/ml

'BAKER ANALYZED' / Atomic Absorption Standard

6820

CONT.

BOX

▶ K

M = 39.10 g/mol NC CODE 3822 00 00

Potassium (K)

998-1002 μg/ml

0.		
820.0100	100 ml	
820.0500	500 ml	

PACKING

Potassium nitrate in water.

Potassium 10000 μ g/ml

(Matrix: 2% nitric acid) / 'BAKER INSTRA-ANALYZED' / Plasma Standard

5741

K

M = 39.10 g/mol

NC CODE 3822 00 00

R: 36/38 **S**: 26



Certificate Provided Reporting Actual Lot Analysis 9980-10020 μg/ml Potassium (K)

NO.

PACKING 5741.0100 100 ml

CONT.

Prepared from the highest purity raw material available, generally greater than 99.999% spectral purity. The content of the solution is confirmed to be accurate to within 0.2 %. Typically 10000 $\mu g/ml$.

The certificate of analysis provided reports actual lot analysis. The certificate also lists the trace impurities.



Potassium Acetate 'BAKER ANALYZED' / ACS

0195

► CH₃COOK

M = 98.15 g/mol **CAS NO.** 127-08-2

EINECS 204-822-2 NC CODE 2915 29 00

Meets ACS Specifications	
Assay (by non-aqueous titrn.)	min. 99.0%
Calcium (Ca)	max. 0.005%
Chloride (CI)	max. 0.003%
Insoluble Matter	max. 0.005%
Magnesium (Mg)	max. 0.002%
pH of 5% Solution at 25°C	6.5-9.0
Phosphate (PO ₄)	max. 0.001%
Sodium (Na)	max. 0.03%
Sulfate (SO ₄)	max. 0.002%
Trace Impurities (in ppm):	
Heavy Metals (as Pb)	max. 5
Iron (Fe)	max. 5

PRODUCT CONT. **PACKING** вох 0195.1000

Potassium Acetate

See for qualities meeting the requirements of cGMP, USP/NF, BP/Ph.Eur, JP or FCC, page 36

Potassium Alum Powder

See for qualities meeting the requirements of cGMP, USP/NF, BP/Ph.Eur, JP or FCC, page 36

Potassium Bicarbonate

See Potassium Hydrogen Carbonate

Potassium Bicarbonate

ee for qualities meeting the requirements of cGMP, USP/NF, BP/Ph.Eur, JP or FCC, page 36

Potassium Bichromate

Potassium Biphosphate

Potassium Biphthalate

Potassium Bisulfate

See Potassium Hydrogen Sulfate

Potassium Bisulfate, fused

Potassium Bromate

1/60 mol/l / 0.1 N / 'BAKER ANALYZED'

M = 167.00 g/mol **EINECS** 231-829-8

NC CODE 2829 90 40 R: 45

S: 36-45-53

Titer (N)

0.0995- 0.1005

7133.1000

PRODUCT

7133.9020 20 l Polycube

CONT.

Volumetric Solution, ready for use.

PACKING

Innovation is principal to our business.

Potassium Bromate

4667 mol/L / 0.1 N / DILUT-IT

► KBrO₃

M = 167.00 g/mol **CAS NO.** 7758-01-2

EINECS 231-829-8 NC CODE 2829 90 40

UN/ID NO. 1484

ADR/RID 5.1 O2

IMDG 5.1/II **R:** 45

S: 36-45-53

PRODUCT PACKING CONT. NO. вох 4667 1 amp.

Volumetric Concentrate, for dilution to 1 l.

Potassium Bromide

'BAKER INSTRA-ANALYZED' / For IR Analysis

0373

0202

0218

► KBr

M = 119.01 g/mol **CAS NO.** 7758-02-3 **EINECS** 231-830-3 NC CODE 2827 51 00

Assay (argentometric titrn.) min. 99.0% Insoluble Matter max. 0.005% Loss on Drying at 105°C max. 0.05% pH of 5% Solution at 25°C 5.0-7.5 Suitability for Infrared Analysis passes test max. 0.002% Sulfate (SO₄) Trace Impurities (in ppm): Nitrogen Compounds (as N) max. 5 **PRODUCT PACKING** CONT. вох NO. 0373.0100 100 g

Potassium Bromide

'BAKER ANALYZED' / ACS

► KBr

► KBr

M = 119.01 g/mol **CAS NO.** 7758-02-3 **EINECS** 231-830-3 NC CODE 2827 51 00

Meets ACS Specifications. Meets Reagent Specifications for testing USP/NF monographs

Assay (argentometric titrn.) min. 99.0% Barium (Ba) max. 0.002% Bromate (BrO₃) max. 0.001% Calcium (Ca) max. 0.002% Chloride (CI) max. 0.2% Insoluble Matter max. 0.005% lodate (IO₃) max. 0.001% lodide (I) max. 0.001% Magnesium (Mg) max. 0.001% Nitrogen Compounds (as N) max. 0.005% pH of 5% Solution at 25°C 5.0-8.8 Sodium (Na) max. 0.02% Sulfate (SO₄) max. 0.005% Trace Impurities (in ppm): Heavy Metals (as Pb) max. 5 Iron (Fe) max. 5

PRODUCT	PACKING	CONT.
NO.		вох
0202.0250	250 g	6
0202.1000	1 kg	

Potassium Bromide

Assay (dried basis)

M = 119.01 g/mol Acidity or Alkalinity

CAS NO. 7758-02-3 **EINECS** 231-830-3 NC CODE 2827 51 00

98.0-100.5% passes test Appearance of solution passes test Bromates (as BrO₃) passes test Chlorides (as CI) max. 0.6% Heavy Metals (as Pb) max. 10 ppm Identification passes test lodides (as I) passes test Iron (Fe) max. 20 ppm Loss on Drying max. 1.0% Magnesium and Alkaline-earth Metals (as Ca) max. 200 ppm Sulfates (as SO₄) max. 100 ppm

PRODUCT	PACKING	CONT.
NO.		вох
0218.1000	1 kg	
0218.9025	25 kg	



Potassium Bromide - Borax Solution

pH 8.6 / 'BAKER ANALYZED' / For coulometric amm

NC CODE 3822 00 00

7486

PRODUCT **PACKING** CONT. NO. вох 7486.9020 20 l Polycube

CONT.

BOX

CONT.

BOX

6

Volumetric Solution, ready for use.

PACKING

1 kg

25 kg

PACKING

1 kg

50 kg

PRODUCT

0204.1000

0204.9025

PRODUCT

0205.1000

0205.9050

NO.

NO.

min. 99%

Potassium Carbonate Anhydrous

0204 'BAKER ANALYZED' / ACS

► K₂CO₃

M = 138.21 g/mol

CAS NO. 584-08-7

EINECS 209-529-3

NC CODE 2836 40 00

R: 22-36/37/38

S: 22-26



Exceeds ACS Specifications. Meets Reagents		
Specifications for testing USP/NF monographs		
Assay (by acid-base titration)	min. 99.0%	
Calcium (Ca)	max. 0.002%	
Chloride (CI)	max. 0.003%	
Copper (Cu)	max. 100 ppb	
Insoluble Matter	max. 0.01%	
Magnesium (Mg)	max. 0.002%	
Phosphate (PO ₄)	max. 0.001%	
Silica (SiO ₂)	max. 0.005%	
Sodium (Na)	max. 0.02%	
Sulfur Compounds (as SO ₄)	max. 0.004%	
Trace Impurities (in ppm):		
Heavy Metals (as Pb)	max. 5	
Iron (Fe)	max. 5	

Potassium Carbonate Anhydrous

► K₂CO₃

M = 138.21 g/mol

CAS NO. 584-08-7 EINECS 209-529-3

NC CODE 2836 40 00

R: 22-36/37/38

S: 22-26



Potassium Carbonate, Anhydrous

See for qualities meeting the requirements of cGMP, USP/NF, BP/Ph.Eur, JP or FCC, page 36

Assay

Potassium Chlorate

0207 'BAKER'

► KCIO₃

M = 122.55 g/mol

CAS NO. 3811-04-9

EINECS 223-289-7

NC CODE 2829 19 00

EC NO. 17 004 00 3

UN/ID NO. 1485

ADR/RID 5.1 O2

IMDG 5.1/II

R: 20/22-9 **S:** 13-16-27



Assay	min. 99.0%
Heavy Metals (as Pb)	max. 0.002%
Insoluble Matter	max. 0.01%
Iron (Fe)	max. 0.005%

PRODUCT	PACKING	CONT.
NO.		вох
0207.1000	1 kg HDPE	
0207.9050	50 kg	

Potassium Chloride

crystal / 'BAKER ULTRAPURE BIOREAGENT'

4001

KCI
NU.

M = 74.56 g/mol **CAS NO.** 7447-40-7 **EINECS** 231-211-8 NC CODE 3104 20 90

99.0-100.5% Assay (argentometric titrn.) max. 0.001% Barium (Ba) Bromide (Br) max. 0.1% Calcium, Magnesium and R₂O₃ Precipitate max. 0.005% DNase Activity none detected Insoluble Matter max. 0.005% lodide (I) max. 0.002% max. 1.0% Loss on Drying at 105°C pH of 5% Solution at 25°C 5.4-8.6 Protease Activity none detected RNase Activity none detected Sodium (Na) max. 0.005% Trace Impurities (in ppm): Arsenic (As) max. 1 Heavy Metals (as Pb) max. 5 Iron (Fe) max. 2

max. 1

PRODUCT	PACKING	CONT
NO.		ВОХ
4001.0500	500 g	
4001.2500	2.5 kg	

Potassium Chloride

'BAKER ANALYZED' / ACS

0208

► KCI

M = 74.56 g/mol **CAS NO.** 7447-40-7 **EINECS** 231-211-8 NC CODE 3104 20 90

Meets ACS Specifications. Meets Reagent Specifications for testing USP/NF monographs

Nickel (Ni)

Assay (argentometric titrn.)	99.0-100.5%
Barium (Ba)	passes test
Bromide (Br)	max. 0.01%
Calcium (Ca)	max. 0.002%
Chlorate and Nitrate (as NO ₃)	max. 0.003%
Insoluble Matter	max. 0.005%
lodide (I)	max. 0.002%
Loss on Drying at 105°C	max. 1.0%
Magnesium (Mg)	max. 0.001%
Nitrogen Compounds (as N)	max. 0.001%
pH of 5% Solution at 25°C	5.4-8.6
Sodium (Na)	max. 0.005%
Sulfate (SO ₄)	max. 0.001%
Trace Impurities (in ppm):	
Arsenic (As)	max. 1
Heavy Metals (as Pb)	max. 5
Iron (Fe)	max. 2
Nickel (Ni)	max. 1
Phosphate (PO ₄)	max. 5

PRODUCT	PACKING	CONT.
NO.		вох
0208.0250	250 g	
0208.1000	1 kg	6
0208.5000	5 kg	
0208.7100	100 lbs	

Potassium Chloride

0209

► KCI

M = 74.56 g/mol **CAS NO.** 7447-40-7 **EINECS** 231-211-8 NC CODE 3104 20 90

Assay	99.0-100.5%
Acid or alkaline reacting substances	passes test
Appearance of solution	passes test
Barium (Ba)	passes test
Bromide (Br)	max. 0.1%
Calcium and Magnesium	passes test
Heavy Metals (as Pb)	max. 10 ppm
Identification	passes test
lodide (I)	passes test
Loss on Drying at 105°C	max. 1.0%
Magnesium, earth-alkaline metals (as Ca)	max. 200 ppm
Organic Volatile Impurities	passes test
Sulfate (SO ₄)	max. 300 ppm
Trace Impurities (in ppm):	
Iron (Fe)	max. 20 ppm
Sodium (Na)	max. 0.1%

PRODUCT	PACKING	CONT.
NO.		вох
0209.1000	1 kg	6
0209.5000	5 kg	
0209.9050	50 kg	

D E F

Potas

Potassium Chloride

0509

► KCI

M = 74.56 g/mol **CAS NO.** 7447-40-7 **EINECS** 231-211-8

NC CODE 3104 20 90

min. 99% Assay Acidity and Alkalinity passes test Heavy Metals (as Pb) max. 0.001% Iron (Fe) max. 0.005% Sodium (Na) max. 0.1% Sulfate (SO₄) max. 0.01%

min. 99.0%

max. 0.005%

max. 0.005%

max. 0.005%

max. 0.02%

max. 0.03%

8.6-9.8

PRODUCT	PACKING	CONT.
NO.		вох
0509.0500	500 g	6
0509.1000	1 kg	6
0509.5000	5 kg	



Potassium Chloride

See for qualities meeting the requirements of cGMP, USP/NF, BP/Ph.Eur, JP or FCC, page 36

Assay

Calcium (Ca)

Chloride (CI)

Sulfate (SO₄)

Insoluble Matter

Insoluble Matter

pH of 5% Solution at 25°C

Sodium (Na) (by AAS)

0210

Potassium Chromate

'BAKER ANALYZED'

K₂CrO₄

M = 194.20 g/mol

CAS NO. 7789-00-6 **EINECS** 232-140-5

NC CODE 2841 50 00 EC NO. 24 006 00 8

R: 36/37/38-43-46-49-50/53

S: 45-53-60-61





dangerous for the environment

PRODUCT PACKING CONT. BOX NO. 0210.0500 500 g

Potassium Chromate

► K₂CrO₄

0211

M = 194.20 g/mol

CAS NO. 7789-00-6

EINECS 232-140-5

NC CODE 2841 50 00

EC NO. 24 006 00 8

R: 36/37/38-43-46-49-50/53

S: 45-53-60-61







CONT.

BOX

Potassium Citrate, Monohydrate

ee for qualities meeting the requirements of cGMP, USP/NF, BP/Ph.Eur, JP or FCC, page 36

Potassium Cyanide

0213 'BAKER ANALYZED' / ACS

► KCN

M = 65.12 g/mol

CAS NO. 151-50-8 **EINECS** 205-792-3

NC CODE 2837 19 00

EC NO. 6 007 00 5

UN/ID NO. 1680

ADR/RID 6.1 T5 IMDG 6.1/I

R: 26/27/28-32-50/53

S: 28-29-45-60-61-7



environment



Meets ACS Specifications	
Assay (argentometric titrn.)	min. 96.0%
Chloride (CI)	max. 0.5%
Iron, total (as Fe)	max. 0.03%
Phosphate (PO ₄)	max. 0.005%
Sodium (Na)	max. 0.5%
Sulfate (SO ₄)	max. 0.04%
Sulfide (S)	max. 0.003%
Thiocyanate (SCN)	max. 0.02%
Trace Impurities (in ppm):	
Lead (Pb)	max. 2

PRODUCT	PACKING	CONT.
NO.		вох
0213.0100	100 g	
0213.9050	50 kg	

Potassium Dichromate

Primary Standard / ACS

0214

	crystal / BAKER ANALYZED /
► K ₂ Cr ₂ O ₇	
M =	294.19 g/mol
CAS NO.	7778-50-9
EINECS	231-906-6
NC CODE	2841 50 00
EC NO.	24 002 00 6
UN/ID NO.	3288
ADR/RID	6.1 T5
IMDG	6.1/III
R:	21-25-26-34-42/43-45-46-48/23
	50/53-60-61-8
S:	45-53-60-61

Exceeds ACS Specifications	
Assay (dried basis)	99.95-100.05%
Calcium (Ca)	max. 0.001%
Chloride (CI)	max. 0.001%
Insoluble Matter	max. 0.005%
Iron (Fe)	max. 0.001%
Loss on Drying at 105°C	max. 0.02%
pH of 5% Solution at 25°C	3.5-4.0
Sodium (Na)	max. 0.01%
Sulfate (SO ₄)	max. 0.005%

PACKING	CONT
	BO)
500 g	







environment

Potassium Dichromate

'BAKER ANALYZED' / ACS

0215

0216

7250

CONT.

вох

	$K_2Cr_2O_7$
--	--------------

M = 294.19 g/mol **CAS NO.** 7778-50-9 **EINECS** 231-906-6 NC CODE 2841 50 00 **EC NO.** 24 002 00 6

UN/ID NO. 3288

ADR/RID 6.1 T5 IMDG 6.1/III

R: 21-25-26-34-42/43-45-46-48/23-50/53-60-61-8

S: 45-53-60-61



environment





Sodium (Na)

Sulfate (SO₄)

Assay

Exceeds ACS Specifications Assay min. 99.0% Calcium (Ca) max. 0.002% Chloride (CI) max. 0.001% Insoluble Matter max. 0.005% Iron (Fe) max. 0.001% Loss on Drying at 105°C max. 0.05%

max. 0.02%

max. 0.005%

min. 99%

PRODUCT

0216.9050

NO.

PACKING

50 kg

CONT. **PRODUCT PACKING** вох NO. 0215.0500 500 g 0215.1000 1 kg

Potassium Dichromate

M = 294.19 g/mol

CAS NO. 7778-50-9

EINECS 231-906-6

NC CODE 2841 50 00

EC NO. 24 002 00 6

UN/ID NO. 3288

K₂Cr₂O₇

ADR/RID 6.1 T5

IMDG 6.1/III

R: 21-25-26-34-42/43-45-46-48/23-50/53-60-61-8

S: 45-53-60-61







dangerous for the environment

Potassium Dichromate 1/24 mol/l / 'BAKER ANALYZED'

NC CODE 2841 50 00

R: 23-42/43-45-46-48/20-52/53-60-

S: 23-24-36/37-45-53-60-61



Titer (mol/l)	0.0415-0.0420

PRODUCT	PACKING	CONT.
NO.		вох
7250.1000	11	6

313



Potassium Dichromate

0.04 mol/I / 'BAKER ANALYZED' / For COD determination according AFNOR NFT90-101

Molarity (M)

UN/ID NO. 3264 ADR/RID 8 C1

7631

R: 23-35-42/43-45-46-48-52/53-60-61

S: 23-26-36/37/39-45-53-60-61



CONT. **PRODUCT PACKING** NO. вох 7631.1000

Potassium Dichromate / Sulfuric Acid Solution 0.04 mol/l.

CONT.

BOX

Potassium Dichromate

7134 0.02 mol/l / 'BAKER ANALYZED' / For COD determination according DIN 38409 vol.41

Titer (mol/l)

M = 294.19 g/mol **CAS NO.** 7778-50-9 **EINECS** 231-906-6

NC CODE 2841 50 00

R: 26/27/28-33-42/43-45-46-50/53-60-61

S: 29-36/37-45-53





dangerous for the environment

Volumetric Solution, ready for use

11

PACKING

PRODUCT

7134,1000

NO.

0.0396 - 0.0404

0.0198 - 0.0202

0.0048-0.0052

Potassium Dichromate

7249 0.005 mol/l / 'BAKER ANALYZED' / For COD determination according DIN 38409 vol.41

Titer (mol/l)

NC CODE 2841 50 00

R: 20-45-46 **S:** 36-45-53

PRODUCT **PACKING** CONT. BOX 7249.1000

Potassium Dichromate

4671 $^{1}/_{60}$ mol/l / $^{1}/_{10}$ equiv. = 4.903g; 0.1N / DILUT-IT

► K₂Cr₂O₇

M = 294.19 g/mol **CAS NO.** 7778-50-9

EINECS 231-906-6

NC CODE 2841 50 00 **EC NO.** 24 002 00 6

> R: 22-23-42/43-45-46-48/20-51/53-60-61

S: 29-36/37-45-53

environment



PRODUCT PACKING CONT. NO. BOX 1 amp.

Volumetric Concentrate, for dilution to 1 l.

J.T.Baker: over 100 years of experience.

See chapter 1 of this catalogue.



Potassium Dihydrogen Phosphate

4921

► KH₂PO₄

M = 136.09 g/mol **CAS NO.** 7778-77-0 **EINECS** 231-913-4 NC CODE 2835 24 00

Certificate Provided Reporting Actual Lot Analysis Actual Lot Analysis Lot No. A49470

Assay (dried basis)	99.7%	
1 5 1 15000 (6 0.11)	- 0.000/	
Loss on Drying at 150°C (for 2 Hours)	< 0.02%	
Particulate Matter	< 0.001%	
pH of 5% Solution at 25°C	4.3	
Metallic Impurities in parts per million (μg/g):		
Aluminium (AI)	0.4	

pH of 5% Solution at 25°C	4.3
Metallic Impurities in parts per n	nillion (μg/g):
Aluminium (AI)	0.4
Barium (Ba)	0.3
Bismuth (Bi)	0.2
Cadmium (Cd)	0.4
Calcium (Ca)	0.4
Chromium (Cr)	0.4
Cobalt (Co)	0.7
Copper (Cu)	0.4
Iron (Fe)	0.5
Lead (Pb)	0.4
Magnesium (Mg)	0.3
Manganese (Mn)	0.4
Mercury (Hg)	0.0002
Molybdenum (Mo)	0.2
Nickel (Ni)	0.3
Niobium (Nb)	0.4

Zirconium (Zr) Non-Metallic Impurities in parts per million ($\mu g/g$): Arsenic (As)

0.3

< 0.2 0.4

< 0.2

0.3

0.4

0.4

Silver (Ag)

Sodium (Na)

Strontium (Sr) Tin (Sn)

Titanium (Ti)

Vanadium (V)

Zinc (Zn)

Fluoride (F)	2
Halide (as CI)	< 10
Nitrogen Compounds (as N)	< 10
Silicon (Si)	0.2
Sulfate (SO ₄)	< 20

Ultraviolet Absorbance (1M aqueous Solution; 1.00-cm path vs water):

at 2 is nm	0.04
at 220 nm	0.04
at 254 nm	0.04
at 280 nm	0.03
at 400 nm	< 0.01

PRODUCT PACE	PACKING	CONT.
NO.		вох
4921.0100	100 g Glass	
4921.1000	1 kg Glass	

For Laboratory, Research or Manufacturing Use.

Potassium Dihydrogen Phosphate crystal / 'BAKER ULTRAPURE BIOREAGENT'

KΗ	2PO4	

M = 136.09 g/mol **CAS NO.** 7778-77-0 **EINECS** 231-913-4 NC CODE 2835 24 00

For Liquid Chromatography and Molecular Biology applications

Assay	min. 99.0%
DNAase Activity	none detected
Heavy Metals (as Pb)	max. 0.001%
Insoluble Matter, Calcium and NH ₄ OH	
Precipitate	max. 0.01%
Iron (Fe)	max. 0.002%
Loss on Drying over H ₂ SO ₄	max. 0.2%
pH of 5% Solution at 25°C	4.1-4.5
Protease Activity	none detected
RNAase Activity	none detected
Trace Impurities (in ppm):	

PRODUCT PA	PACKING	CONT
NO.		BOX
4008.0500	500 g	
4000 2500	2 E ka	

www.jtbaker.com/europe

Arsenic (As)

315



Potassium Dihydrogen Phosphate 'BAKER ANALYZED' / ACS

0240

► KH₂PO₄

M = 136.09 g/mol **CAS NO.** 7778-77-0 **EINECS** 231-913-4 NC CODE 2835 24 00

Exceeds ACS Specifications Assay min. 99.0% Chloride (CI) max. 0.001% Fluoride (F) max. 0.001% Heavy Metals (as Pb) max. 0.001% Insoluble Matter max. 0.01% Iron (Fe) max. 0.002% Lead (Pb) max. 0.001% Loss on Drying at 105°C max. 0.2% Nitrogen Compounds (as N) max. 0.001% pH of 5% Solution at 25°C 4.1-4.5 max. 0.005% Sodium (Na) max. 0.003% Sulfate (SO₄) Trace Impurities (in ppm): Arsenic (As) max. 3

PRODUCT NO.	PACKING	CONT. BOX
0240.1000	1 kg	6

Potassium Dihydrogen Phosphate

1737

► KH₂PO₄

M = 136.09 g/mol **CAS NO.** 7778-77-0

EINECS 231-913-4 NC CODE 2835 24 00 Assay min. 99% Loss on Drying at 105°C max. 0.1% **PRODUCT PACKING** BOX 1737,1000 1 kg 6



Potassium Dihydrogen Phosphate

See for qualities meeting the requirements of cGMP, USP/NF, BP/Ph.Eur, JP or FCC, page 36

0199

Potassium Disulfite

'BAKER ANALYZED

► K₂S₂O₅

M = 222.33 g/mol **CAS NO.** 16731-55-8

EINECS 240-795-3

NC CODE 2833 29 90

R: 31-37-41

S: 26-39



Assay	min. 95.0%
Chloride (CI)	max. 0.005%
Heavy Metals (as Pb)	max. 0.001%
Insoluble Matter	max. 0.005%
Iron (Fe)	max. 0.001%
Trace Impurities (in ppm):	
Arsenic (As)	max. 1

NO. вох 0199.1000 1 kg

CONT.

PACKING

PRODUCT



Potassium Ferricyanide

See Potassium Hexacyanoferrate(III)



Potassium Ferrocyanide

ee Potassium Hexacyanoferrate(II)

Potassium Fluoride Anhydrous

'BAKER ANALYZED'

► KF **M** = 58.10 g/mol **CAS NO.** 7789-23-3 **EINECS** 232-151-5 NC CODE 2826 19 00 **EC NO.** 9 005 00 2 **UN/ID NO.** 1812 **ADR/RID** 6.1 T5 IMDG 6.1/III

R: 23/24/25 **S:** 26-45

Assay	min. 99.0%
Chloride (CI)	max. 0.005%
Free Acid (as HF)	max. 0.05%
Heavy Metals (as Pb)	max. 0.002%
Insoluble Matter	max. 0.01%
Iron (Fe)	max. 0.001%
Potassium Fluorosilicate (K ₂ SiF ₆)	max. 0.05%
Sodium (Na)	max. 0.05%
Sulfate (SO ₄)	max. 0.005%
Sulfite (SO ₃)	max. 0.005%

PRODUCT PACKING CONT. вох 0225.0250 250 g 0225.1000 1 kg

1442

0220

0219

2940

Potassium Hexachloroplatinate(IV)

► K₂PtCl₆ PRODUCT PACKING CONT. Assay (Pt) min. 40% M = 486.01 g/mol Suitable for the preparation of APHA NO. вох **CAS NO.** 16921-30-5 1442.0001 passes test 1 g **EINECS** 240-979-3 NC CODE 2843 90 90 **UN/ID NO.** 1759

S: 26-28A-36/37/39

EINECS 237-722-2

NC CODE 2837 20 00

R: 36/37/38-42/43

ADR/RID 8 C10 IMDG 8/III

Potassium Hexacyanoferrate(II) crystal / 'BAKER ANALYZED' / ACS

max. 0.005%

max. 0.005%

max 0.01%

0220.1000

1 kg

Potassium Hexacyanoferrate(III)
'BAKER ANALYZED' / ACS

► K₃Fe(CN)₆ **Exceeds ACS Specifications PRODUCT PACKING** CONT. M = 329.26 g/mol Assay (by lodometry) min. 99.0% BOX CAS NO. 13746-66-2 Appearance passes test 0219.0250 **EINECS** 237-323-3 Chloride (CI) max. 0.01% 0219.1000 1 ka NC CODE 2837 20 00 Ferro Compounds as [Fe(CN)₆]⁴ max. 0.05% Insoluble Matter max 0.005%

Potassium Hydrogen Carbonate

'BAKER ANALYZED' / ACS

Sulfate (SO₄)

Insoluble Matter

Sulfate (SO₄)

► KHCO₃ Meets ACS Specifications. Meets Reagent **PRODUCT PACKING** CONT. **M** = 100.12 g/mol Specifications for testing USP/NF monographs BOX CAS NO. 298-14-6 Assay (dried basis) 99.7-100.5% 2940.0500 500 g **EINECS** 206-059-0 Calcium (Ca) max. 0.002% NC CODE 2836 40 00 Chloride (CI) max. 0.001% Insoluble Matter max. 0.01% Magnesium (Mg) max. 0.001% Sodium (Na) max. 0.03% Sulfur Compounds (as SO₄) max 0.003% Trace Impurities (in ppm): Ammonium (NH₄) max. 5 Heavy Metals (as Pb) max. 5 Iron (Fe) max. 5 Phosphate (PO₄) max. 5

Certificates of Analysis are available at www.jtbaker.com/europe

Potassium Hydrogen Phosphate Anhydrous

vder, 'BAKER ULTRAPURE BIOREAGE

▶ K₂HPO₄

M = 174.18 g/mol **CAS NO.** 7758-11-4 **EINECS** 231-834-5 NC CODE 2835 24 00

For Liquid Chromatography and Molecular Biology applications Assay min. 99.0% Chloride (CI) max. 0.001% DNAase Activity none detected Insoluble Matter max. 0.01% Iron (Fe) max. 0.001% Loss on Drying at 105°C max. 1.0% pH of 5% Solution at 25°C 8.5-9.6 Protease Activity none detected RNAase Activity none detected Trace Impurities (in ppm): Heavy Metals (as Pb) max. 5

CONT.
BOX

Potassium Hydrogen Phosphate Anhydrous

0241 'BAKER ANALYZED' / ACS

► K₂HPO₄

M = 174.18 g/mol **CAS NO.** 7758-11-4 **EINECS** 231-834-5 NC CODE 2835 24 00

Exceeds ACS Specifications. Meets Reagents Specifications for testing USP/NF monographs Assay Chloride (CI) max. 0.003% Fluoride (F) max. 0.001% Insoluble Matter max. 0.01% Iron (Fe) max. 0.001% Loss on Drying at 105°C max. 1.0% Nitrogen Compounds (as N) max. 0.001% pH of 5% Solution at 25°C Sodium (Na) max. 0.05% Sulfate (SO₄) max. 0.005% Trace Impurities (in ppm): Arsenic (As) max. 1 Heavy Metals (as Pb) max. 5 Lead (Pb) max. 5

PRODUCT NO.	PACKING	CONT.
		вох
0241.0500	500 g	6
0241.1000	1 kg	6

0242

Potassium Hydrogen Phosphate Anhydrous

► K₂HPO₄

M = 174.18 g/mol **CAS NO.** 7758-11-4 **EINECS** 231-834-5 NC CODE 2835 24 00

Assay 98.0-100.5% Arsenic (As) max. 3 ppm Carbonate (CO₃) passes test Chloride (CI) max. 0.03% Heavy Metals (as Pb) max. 0.001% Identification passes test Insoluble Substances max. 0.2% Iron (Fe) max. 0.003% Limit of fluoride max. 0.001% Limit of monobasic or tribasic salt passes test Loss on Drying at 105°C max. 1.0% pH of 5% Solution at 25°C 8.5 - 9.6 Sodium (Na) passes test Sulfate (SO₄) max. 0.1%

PRODUCI	PACKING	CON1.	
NO.		вох	
0242.1000	1 kg		



Potassium Hydrogen Phosphate Anhydrous

Innovation is principal to our business.

Potassium Hydrogen Phthalate

4889

► 1-KOCOC₆H₄-2-COOH

M = 204.23 g/mol **CAS NO.** 877-24-7 **EINECS** 212-889-4 NC CODE 2917 39 80

Certificate Provided Reporting Actual Lot Analysis Actual Lot Analysis Lot No. B02472

Assay (dried basis)	99.97%
Identification (by IR)	passes test
Particulate Matter	< 0.0001%
pH of 0.05 M Solution at 25°C	4.01

pH of 0.05 M Solution at 25°C	4.0
Metallic Impurities in parts per million $(\mu g/g)$:	
Aluminium (AI)	< 0.2
Barium (Ba)	< 0.2
Bismuth (Bi)	< 2.2
Cadmium (Cd)	< 0.
Calcium (Ca)	< 0.
Chromium (Cr)	< 0.2
Cobalt (Co)	< 0.2
Copper (Cu)	< 0.
Iron (Fe)	< 0.2
Lead (Pb)	< 0.2
Magnesium (Mg)	< 0.
Manganese (Mn)	< 0.
Mercury (Hg)	0.0002
Molybdenum (Mo)	< 0.
Nickel (Ni)	< 0.2
Niobium (Nb)	< 0.
Silver (Ag)	< 0.3
Sodium (Na)	< 0.
Strontium (Sr)	< 0.2
Tin (Sn)	< 0.2
Titanium (Ti)	< 0.2
Vanadium (V)	< 0.
Zinc (Zn)	< 0.2
Zirconium (Zr)	< 0.2
Non-Metallic Impurities in parts per million (µ	g/g):

PRODUCT	PACKING	CONT
NO.		B0
4889.0025	25 g Glass	

For Primary Standard and Research Applications. For Laboratory, Research or Manufacturing Use.

Potassium Hydrogen Phthalate 'BAKER ANALYZED' / Primary Standard / ACS

► 1-KOCOC₆H₄-2-COOH

M = 204.23 g/mol **CAS NO.** 877-24-7 **EINECS** 212-889-4 NC CODE 2917 39 80

Meets ACS Specifications	
Assay (dried basis)	99.95 - 100.05%
Chlorine Compounds (as CI)	max. 0.003%
Insoluble Matter	max. 0.005%
pH of 0.05 M Solution at 25°C	4.00 - 4.02
Sodium (Na)	max. 0.005%
Sulfur Compounds (as S)	max. 0.002%

< 30

< 0.1

< 5

Trace Impurities (in ppm):

Chloride Compounds (as CI)

Sulfur Compounds (as S)

Silicon (Si)

Heavy Metals (as Pb)	max. 5
Iron (Fe)	max. 5

		017

PRODUCI	PACKING	CONI
NO.		ВОХ
0197.0250	250 g	ϵ

This reagent is satisfactory for use as a pH standard. For use as an acidimetric standard this material should be lightly crushed and dried for 2 hours at 120°C to remove any absorbed moisture.

7525

BOX

Potassium Hydrogen Phthalate

BAKER ANALYZED

Assay (titrimetric) 1.04-1.08 g/l

PRODUCT	PACKING	
NO.		
7525.9025	25	

► 1-KOCOC₆H₄-2-COOH **M =** 204.23 g/mol **CAS NO.** 877-24-7 **EINECS** 212-889-4 NC CODE 2917 39 80



Potassium Hydrogen Sulfate

198 'BAKER ANALYZED

S: 26-36/37/39-45

M = 136.17 g/mol

CAS NO. 7646-93-7

EINECS 231-594-1

NC CODE 2833 29 90

EC NO. 16 056 00 4

UN/ID NO. 2509

ADR/RID 8 C2

IMDG 8/II

R: 34-37

35.0-37.0% Acidity (as H2SO4) Calcium and Magnesium Precipitate max. 0.005% max. 0.001% Chloride (CI) Heavy Metals (as Pb) max. 0.001% Insoluble Matter and NH₄OH Precipitate max. 0.01% Iron (Fe) max. 0.001% Phosphate (PO₄) max. 0.001% Trace Impurities (in ppm): Arsenic (As) max. 1

 PRODUCT
 PACKING
 CONT.

 NO.
 BOX

 0198.0500
 500 g

Potassium Hydroxide

0222 Pellets / 'BAKER ANALYZED' / ACS

M = 56.11 g/mol
CAS NO. 1310-58-3
EINECS 215-181-3
NC CODE 2815 20 10
EC NO. 19 002 00 8
UN/ID NO. 1813
ADR/RID 8 C6
IMDG 8/II
R: 22-35
S: 26-36/37/39-45

1.0

Exceeds ACS Specifications min. 86.0% Ammonium Hydroxide Precipitate max. 0.02% Chloride (CI) max. 0.001% max. 0.001% Heavy Metals (as Ag) Insoluble Matter max. 30 ppm Potassium Carbonate (K₂CO₃) max. 0.5% max. 0.05% Sodium (Na) Trace Impurities (in ppm): Iron (Fe) max. 3 Mercury (Hg) max. 0.1 Nickel (Ni) max 2 Nitrogen Compounds (as N) max. 3 Phosphate (PO₄) Sulfate (SO₄) max. 5

PRODUCT	PACKING	CONT
NO.		BOX
0222.1000	1 kg	ϵ
0222.5000	5 kg	

Potassium Hydroxide

0385 Pellets / 'BAKER ANALYZED'

M = 56.11 g/mol

CAS NO. 1310-58-3

EINECS 215-181-3

NC CODE 2815 20 10

EC NO. 19 002 00 8

UN/ID NO. 1813

ADR/RID 8 C6

IMDG 8/II

R: 22-35

5: 26-36/37/39-45

Assay min. 85% Aluminium (AI) max. 0.001% max. 0.001% Calcium (Ca) Potassium Carbonate (K₂CO₃) max. 1.0% Silica (SiO₂) max. 0.005% Sodium (Na) max. 0.5% Trace Impurities (in ppm): Chloride (CI) max. 5 Heavy Metals (as Pb) max. 5 Iron (Fe) max. 5 Nickel (Ni) max. 5 Nitrogen Compounds (as N) Phosphate (PO₄) max. 5 Sulfate (SO₄) max. 5

 PRODUCT NO.
 PACKING
 CONT.

 NO.
 BOX

 0385.0500
 500 g
 6

 0385.1000
 1 kg
 6

 0385.5000
 5 kg
 4

 0385.9050
 50 kg
 4

Find more Chromatography information at www.jtbaker.com/chromatography



Potassium Hydroxide Pellets / 'BAKER'

0224

► KOH

M = 56.11 g/mol **CAS NO.** 1310-58-3 **EINECS** 215-181-3

NC CODE 2815 20 10 **EC NO.** 19 002 00 8

UN/ID NO. 1813

ADR/RID 8 C6

IMDG 8/II

R: 22-35

S: 26-36/37/39-45



Assay (KOH)	min. 85.0%
Heavy Metals (as Pb)	max. 0.003%
Identification	passes test
Insoluble Substances	passes test
Potassium Carbonate (K ₂ CO ₃)	max. 3.5%

PRODUCT	PACKING	CONT.
NO.		вох
0224.1000	1 kg	
0224.5000	5 kg	
0224.9050	50 kg	

Potassium Hydroxide

0387

▶ КОН

M = 56.11 g/mol **CAS NO.** 1310-58-3 **EINECS** 215-181-3 NC CODE 2815 20 10 **EC NO.** 19 002 00 8

UN/ID NO. 1813 ADR/RID 8 C6

IMDG 8/II **R:** 22-35

S: 26-36/37/39-45



Assay	min. 84%
Aluminium (AI)	max. 0.002%
Carbonate (as K ₂ CO ₃)	max. 2%
Chloride (CI)	max. 0.01%
Heavy Metals (as Pb)	max. 0.002%
Iron (Fe)	max. 0.002%
Sulfate (SO ₄)	max. 0.01%

min. 45.0%

max. 0.5%

max. 0.02%

max. 5

max. 2

max. 5

max. 5

max. 2

max. 0.005% max. 0.005%

PRODUCT	PACKING	CONT.
NO.		вох
0387.9050	50 kg	

Potassium Hydroxide, Pellets
See for qualities meeting the requirements of cGMP, USP/NF, BP/Ph.Eur, JP or FCC, page 36

Chloride (CI)

Ammonium Hydroxide Precipitate

Potassium Carbonate (K₂CO₃)



Potassium Hydroxide

-/	ш	n	-5

► KOH	
M =	56.11 g/mol
1 i =	1.46 kg
CAS NO.	1310-58-3
EINECS	215-181-3
NC CODE	2815 20 90
EC NO.	19 002 00 8
UN/ID NO.	1814
ADR/RID	8 C5
IMDG	8/11

Sodium (Na) (by AAS) Sulfate (SO₄) max. 0.001% Trace Impurities (in ppm): Heavy Metals (as Ag) Iron (Fe) Nickel (Ni) **R:** 22-35 Nitrogen Compounds (as N) **S:** 26-36/37/39-45 Phosphate (PO₄)

FRODUCI	FACKING	COIVI.
NO.		вох
7063.1000	11	6
7063.9025	25	
For safe han	dling of 25 l tin cans, see Self-cl	osing tap.

www.jtbaker.com/europe



Potassium Hydroxide 1 mol/I / 'BAKER ANALYZED' 7129

M = 56.11 g/mol 0.995-1.005 PRODUCT PACKING Titer (mol/l) **1 l** = 1.05 kg NO. **CAS NO.** 1310-58-3 7129.1000 11 **EINECS** 215-181-3 7129.9020 20 l Polycube NC CODE 2815 20 90 **UN/ID NO.** 1814 Volumetric Solution, ready for use. ADR/RID 8 C5 Each lot of this product is standardized potentiometrically IMDG 8/II against NIST tracable reference standard. **R:** 35 **S:** 24-26-36/37/39-45

CONT.

вох

Potassium Hydroxide 0.5 mol/l / 'BAKER ANALYZED'

7203

▶ кон	Titer (mol/l)	0.4975-0.5025	PRODUCT	PACKING	CONT.
M = 56.11 g/mol			NO.		вох
1 l = 1.02 kg			7203.5000	51	
CAS NO. 1310-58-3					
EINECS 215-181-3			Volumetric .	Solution, ready for use.	
NC CODE 2815 20 90					
EC NO. 19 002 00 8					
UN/ID NO. 1814					
ADR/RID 8 C5					
IMDG 8/II					
R: 34					
S: 26-36/37/39-45					
C					

Potassium Hydroxide 0.23 mol/l / 'BAKER ANALYZED'

7639


► KOH	Molarity (M)	0.228-0.232	PRODUCT	PACKING	CONT.
M = 56.11 g/mol			NO.		вох
CAS NO. 1310-58-3			7639.9010	10 l	
EINECS 215-181-3					
NC CODE 2815 20 90			Volumetric S	olution, ready for use.	
EC NO. 19 002 00 8					
UN/ID NO. 1814					
ADR/RID 8 C5					
IMDG 8/II					
R: 36/38					
S: 26					
irritant Xi					

Potassium Hydroxide

7 150	O.THIOI/T / DAKEIT ANALTZED					
M =	56.11 g/mol	Titer (mol/l)	0.0995-0.1005	PRODUCT	PACKING	CONT.
1 =	1.01 kg			NO.		вох
CAS NO.	1310-58-3			7130.1000	11	6
EINECS	215-181-3			7130.9020	20 l Polycube	
NC CODE	2815 20 90					
UN/ID NO.	1814			Volumetric S	olution, ready for use.	
ADR/RID	8 C5			Each lot of th	nis product is standardized poten	tiometrically
IMDG	8/II			against NIST	tracable reference standard.	
R:	36/38					
S:	24-26					
	Xi					

Potassium Hydroxide 1 mol/I / DILUT-IT / 1 equiv. = 56.11g; 1N

4674



Potassium Hydroxide

0.5 mol/l / DILUT-IT

▶ КОН		PRODUCT	PACKING	CONT.
NC CODE	2815 20 90	NO.		вох
UN/ID NO.	1814	4866	1 amp.	
ADR/RID	8 C5			
IMDG	8/II	Volumetric Co	ncentrate, for dilution to 1 l.	
R:	35			
S:	26-36/37/39-45			
	C corrosive			

Potassium Hydroxide

0.1 mol/I / DILUT-IT / $^{1}/_{10}$ equiv. = 5.611g; 0.1N

4673

4866

► KOH	PRODUCT	PACKING	CONT.
M = 56.11 g/mol	NO.		вох
CAS NO. 1310-58-3	4673	1 amp.	
EINECS 215-181-3			
NC CODE 2815 20 90	Volumetric Co	oncentrate, for dilution to 1 l.	
EC NO. 19 002 00 8			
UN/ID NO. 1814			
ADR/RID 8 C5			
IMDG 8/II			
R: 35			
S: 26-27-37/39-45			
corrosive			

0.099 - 0.101

Potassium Hydroxide in tert-Butanol

Titer (mol/l)

0.1 mol/l / 'BAKER ANALYZED'

PRODUCT PACKING CONT. вох NO. 7113.1000 11 6

M =	56.11 g/mol
1 i =	0.77 kg
FLASHPOINT	11 °C
NC CODE	2905 14 10
UN/ID NO.	1120
ADR/RID	3 F1
IMDG	3/II
R:	11-20/21/22-

-36/38-68/20/21/22

S: 16-26-36/37



highly flammable

Volumetric Solution, ready for use. KEEP under NITROGEN.



Potassium Hydroxide in Ethanol 0.5 mol/I / 'BAKER ANALYZED'

Titer (mol/l)

Titer (mol/l)

Titer (mol/l)

7131

M = 56.11 g/mol **11=** 0.81 kg FLASHPOINT 11 °C NC CODE 2207 10 00 UN/ID NO. 2924 ADR/RID 3 FC IMDG 3/II

R: 11-34

S: 26-36/37/39-45





PRODUCT PACKING CONT. NO. вох 7131.1000 11 6 7131.2500 2.5 l 7131.9010 10 l 7131.9025 25 I For safe handling of 25 l tin cans, see Self-closing tap.

Volumetric Solution, ready for use.

A light red color can occur in time due to product characteristics.

Each lot of this product is standardized potentiometrically against NIST tracable reference standard.

Potassium Hydroxide in Ethanol

7132 0.1 mol/l / 'BAKER ANALYZED'

M = 56.11 g/mol**1 l** = 0.81 kg FLASHPOINT 11 °C NC CODE 2207 10 00 UN/ID NO. 2924

ADR/RID 3 FC IMDG 3/II

> **R:** 11-36/38 **S:** 16-26-7

> > highly flammable



0.0995-0.1005

0.4975-0.5025

PRODUCT PACKING CONT. вох NO. 7132.1000 11 6 7132.2500 2.5 | 4 7132.9025 25 I For safe handling of 25 l tin cans, see Self-closing tap.

Volumetric Solution, ready for use.

A light red color can occur in time due to product characteristics.

Each lot of this product is standardized potentiometrically against NIST tracable reference standard.

Potassium Hydroxide in IPA

7206 0.1 mol/l / 'BAKER ANALYZED'

FLASHPOINT 12 °C NC CODE 3822 00 00 UN/ID NO. 2924 ADR/RID 3 FC IMDG 3/II

> **R:** 11-36/38 **S:** 16-24/25-26-37-7/9

> > 0 highly flammable



Titer (mol/l) 0.0995-0.1005 **PRODUCT PACKING** CONT. NO. BOX 7206.1000 11

Volumetric Solution, ready for use.

Each lot of this product is standardized potentiometrically against Potassium Hydrogen Phtalate (NIST tracable reference standard).

Potassium Hydroxide in Methanol 7121

0.5 mol/l / 'BAKER ANALYZED'

M = 56.11 g/mol11 = 0.79 kg

► KOH

FLASHPOINT 11 °C

NC CODE 3822 00 00 UN/ID NO. 2924

ADR/RID 3 FC IMDG 3/II

R: 11-23/24/25-34-39/23/24/25

S: 26-36/37/39-45





0.4975-0.5025 PRODUCT

PACKING CONT. вох NO. 7121.1000 7121.9025 25 I For safe handling of 25 l tin cans, see Self-closing tap.

Volumetric Solution, ready for use.

Each lot of this product is standardized potentiometrically against NIST tracable reference standard.

Potassium Hydroxide in Methanol

0.1 mol/l / 'BAKER ANALYZED'



M = 56.11 g/mol

1 l = 0.80 kg

FLASHPOINT 11 $^{\circ}\mathrm{C}$

NC CODE 3822 00 00

UN/ID NO. 2924

ADR/RID 3 FC

IMDG 3/II

R: 11-23/24/25-34-39/23/24/25

S: 36/37/39-45





Titer (mol/l)	0.0997-0.1003
Chloride (CI)	max. 10 ppm
Heavy Metals (as Pb)	max. 1 ppm
Iron (Fe)	max. 0.5 ppm
Nitrogen Compounds (as N)	max. 1 ppm

PRODUCT	PACKING	CONT.
NO.		вох
7045.1000	11	6
7045.2500	2.5	4
7045.9025	25 l	

For safe handling of 25 l tin cans, see Self-closing tap.

Volumetric Solution, ready for use.

PACKING

PRODUCT

PRODUCT

4676

Each lot of this product is standardized potentiometrically against NIST tracable reference standard.

Potassium Iodate

CONT.

7045

► KIO₃ **M =** 214.00 g/mol CAS NO. 7758-05-6 **EINECS** 231-831-9 NC CODE 2829 90 80 **UN/ID NO.** 1479 **ADR/RID** 5.1 O2 IMDG 5.1/II **R:** 8 **S:** 17-35

Exceeds ACS Specifications 99.4-100.4% Chloride and Bromide (as CI) max. 0.01% Insoluble Matter max. 0.005% lodide (I) max. 0.001% Iron (Fe) max. 0.001% max. 0.002% Nitrogen Compounds (as N) pH of 5% Solution at 25°C 5.0-8.0 max. 0.005% Sodium (Na) Sulfate (SO₄) max. 0.005% Trace Impurities (in ppm): Heavy Metals (as Pb) max. 5

NO.		вох
0226.0250	250 g	

Potassium Iodate

 $^{1}/_{60}$ mol/l / $^{1}/_{10}$ equiv. = 3.567g; 0.1N / DILUT-IT

4676

CONT.

► KIO₃

► KI

M = 214.00 g/mol **CAS NO.** 7758-05-6 **EINECS** 231-831-9 NC CODE 2829 90 80

Volumetric Concentrate, for dilution to 1 l.

PACKING

1 amp.

Potassium Iodide

'BAKER ANALYZED' / ACS

0227

M = 166.01 g/mol **CAS NO.** 7681-11-0 **EINECS** 231-659-4

NC CODE 2827 60 00

Meets ACS Specifications min. 99.0% Assay Barium (Ba) max. 0.002% Calcium (Ca) max. 0.002% Chloride and Bromide (as CI) max. 0.01% Insoluble Matter max. 0.005% Loss on Drying at 150°C max. 0.2% max. 0.001% Magnesium (Mg) pH of 5% Solution at 25°C 6.0-9.2 Phosphate (PO₄) max. 0.001% Sodium (Na) max. 0.005% max. 0.005% Sulfate (SO₄) Trace Impurities (in ppm): Heavy Metals (as Pb) max. 5 lodate (IO₃) max. 3 Iron (Fe) max. 3

PRODUCT	PACKING	CONT.
NO.		вох
0227.0250	250 g	6
0227.0500	500 g	6
0227.1000	1 kg	6
0227.9050	50 kg	



Potas

Potassium Iodide

0228

M = 166.01 g/mol **CAS NO.** 7681-11-0

EINECS 231-659-4 NC CODE 2827 60 00

Assay (KI)(dried basis)	99.0 - 101.5%
Alkalinity	passes tes
Heavy Metals (as Pb)	max. 0.001%
Identification	passes tes
lodate (IO ₃)	max. 4 ppm
Loss on Drying at 105°C	max. 1.0%
Nitrate, nitrite and ammonia	passes tes
Organic Volatile Impurities	passes tes
Thiosulfate and Barium	passes tes

PRODUCT	PACKING	CONT.
NO.		вох
0228.1000	1 kg	6
0228.5000	5 kg	

Potassium Iodide

0230

M = 166.01 g/mol **CAS NO.** 7681-11-0

EINECS 231-659-4

NC CODE 2827 60 00

Assay (dried basis) 99.0-100.5% Alkalinity passes test Appearance of solution Heavy Metals (as Pb) Identification lodates (IO₃) passes test Iron (Fe) max. 20 ppm Loss on Drying max. 1.0% Sulfates (as SO₄) max. 150 ppm Thiosulfates (as S₂O₃) passes test

PRODUCT PACKING CONT. вох passes test 0230.0250 max. 10 ppm 0230.1000 passes test 0230.5000 5 kg

Stored protected from light.



Potassium Iodide

ee for qualities meeting the requirements of cGMP, USP/NF, BP/Ph.Eur, JP or FCC, page 36



Potassium Metabisulfite

See Potassium Disulfite



Potassium Metabisulfite

e for qualities meeting the requirements of cGMP, USP/NF, BP/Ph.Eur, JP or FCC, page 36



Potassium Metaperiodate

See Potassium Tetroxoiodate(VII)

Potassium Nitrate

0231 'BAKER ANALYZED' / ACS

KNO ₃	
M =	101.11 g/mol
CAS NO.	7757-79-1
EINECS	231-818-8
NC CODE	2834 21 00
UN/ID NO.	1486
ADR/RID	5.1 O2
IMDG	5.1/III
R:	8
S:	17-41
	Oxidizing

min. 99.0%
max. 0.005%
max. 0.002%
max. 0.005%
max. 0.002%
max. 0.001%
4.5-7.0
max. 0.005%
max. 0.003%
max. 5
max. 5
max. 2
max. 3

PRODUCT	PACKING	CONT.
NO.		вох
0231.0100	100 g	
0231.1000	1 kg	6

Innovation is principal to our business.

Potassium Nitrate

► KNO₃ min. 99% PRODUCT PACKING CONT. Assay **M** = 101.11 g/mol Insoluble Matter max. 0.05% NO. вох **CAS NO.** 7757-79-1 0232.1000 1 kg **EINECS** 231-818-8 NC CODE 2834 21 00 **UN/ID NO.** 1486

S: 17-41

ADR/RID 5.1 O2 **IMDG** 5.1/III **R:** 8

Potassium Nitrate

See for qualities meeting the requirements of cGMP, USP/NF, BP/Ph.Eur, JP or FCC, page 36




0233

0232

Potassium Nitrite 'BAKER ANALYZED' / ACS

- - •

M = 85.11 g/mol CAS NO. 7758-09-0 EINECS 231-832-4 NC CODE 2834 10 00 EC NO. 7 011 00 0 UN/ID NO. 1488 ADR/RID 5.1 02 IMDG 5.1/II R: 25-50-8 5: 45-61



PRODUCI	PACKING	CONI
NO.		BO
0233.0250	250 g	







Potassium Oxalate Monohydrate

'BAKER ANALYZED' / ACS

CONT.

► KOCOCOOK.	H ₂ O
M =	184.24 g/mol
CAS NO.	6487-48-5
EINECS	209-506-8
NC CODE	2917 11 00
EC NO.	607 007 00 3
UN/ID NO.	3282
ADR/RID	6.1 T3
IMDG	6.1/III
R:	21/22
S:	24/25
	Xn

Exceeds ACS Specifications	
Assay	98.5-101.0%
Ammonium (NH ₄)	max. 0.002%
Chloride (CI)	max. 0.002%
Heavy Metals (as Pb)	max. 0.002%
Insoluble Matter	max. 0.01%
Iron (Fe)	max. 0.001%
Neutrality	passes test
pH of 5% Solution at 25°C	7.0-8.5
Sodium (Na)	max. 0.02%
Substances Darkened by Hot H ₂ SO ₄	passes test
Sulfate (SO ₄)	max. 0.01%

NO.		B0)
0235.9050	50 kg	

PACKING

PRODUCT

Find our up-to-date Product Literature at www.jtbaker.com/europe



Potas

Potassium Oxalate Monohydrate

BAKER

► KOCOCOOK.H₂O

M = 184.24 g/mol

CAS NO. 583-52-8 **EINECS** 209-506-8

NC CODE 2917 11 00

EC NO. 607 007 00 3

UN/ID NO. 3282

ADR/RID 6.1 T3

IMDG 6.1/III

R: 21/22

S: 24/25

harmful

Assay	min. 99%
Chloride (CI)	max. 0.002%
Heavy Metals (as Pb)	max. 0.003%
Iron (Fe)	max. 0.001%
Sulfate (SO ₄)	max. 0.02%

PRODUCT	PACKING	CONT.
NO.		вох
1973.1000	1 kg	

Potassium Permanganate

3227 (max. 0.05 ppm Hg) / 'BAKER ANALYZED' / Suitable for Mercury Determination / ACS

Assay

Meets ACS Specifications

► KMnO₄

M = 158.04 g/mol

CAS NO. 7722-64-7

EINECS 231-760-3

NC CODE 2841 61 00

EC NO. 25 002 00 9

UN/ID NO. 1490 **ADR/RID** 5.1 O2

IMDG 5.1/II

R: 22-50/53-8

S: 60-61

N N





min. 99.0%

PRODUCT	PACKING	CONT.
NO.		вох
3227.0500	500 g	

Potassium Permanganate

0237 'BAKER ANALYZED' / ACS

► KIVInO₄

M = 158.04 g/mol

CAS NO. 7722-64-7

EINECS 231-760-3

NC CODE 28416100 EC NO. 25 002 00 9

EL NU. 25 002 (

UN/ID NO. 1490

ADR/RID 5.1 O2 IMDG 5.1/II

R: 22-50/53-8

S: 60-61









Meets ACS Specifications	
Assay	min. 99.0%
Chloride and Chlorate (as CI)	max. 0.005%
Insoluble Matter	max. 0.2%
Sulfate (SO ₄)	max. 0.02%

PRODUCT	PACKING	CONT.
NO.		вох
0237.0250	250 g	
0237.1000	1 kg	



Potassium Permanganate

See for qualities meeting the requirements of cGMP, USP/NF, BP/Ph.Eur, JP or FCC, page 36

www.jtbaker.com/europe

7247

CONT.

CONT.

BOX

6

вох

6

Potassium Permanganate

Titer (N)

0.2 mol/l (1.0N) / 'BAKER ANALYZED

M = 158.04 g/mol **CAS NO.** 7722-64-7

► KMnO₄

EINECS 231-760-3 NC CODE 2841 61 00

R: 51/53

S: 57

for the

Potassium Permanganate 7057

0.995-1.005

0.02 mol/l 0.1N / 'BAKER ANALYZED

► KMnO₄ **M** = 158.04 g/mol CAS NO. 7722-64-7 **EINECS** 231-760-3 NC CODE 2841 61 00

R: 52/53 **S:** 61

Titer (N) 0.0997-0.1003 Trace Impurities (in ppm): Chloride and Chlorate (as CI) max. 10 ppm Sulfate (SO₄) max. 20 ppm

PRODUCT NO. 7057.1000 7057.2500

PRODUCT

NO.

4677

PRODUCT

7247.1000

NO.

PACKING

11

Volumetric Solution, ready for use

PACKING

11

2.5 l

PACKING

1 amp.

Volumetric Concentrate, for dilution to 1 l.

Volumetric Solution, ready for use. Each lot of this product is standardized at 25°C against Sodium Oxalate (NIST tracable reference standard).

Potassium Permanganate

 $0.02 \text{ mol/l} / \frac{1}{10} \text{ equiv.} = 3.161g; 0.1N / DILUT-IT}$

4677 CONT.

BOX

► KMnO₄

M = 158.04 g/mol

CAS NO. 7722-64-7

EINECS 231-760-3

NC CODE 2841 61 00

R: 51/53

S: 57



environment

Potassium Permanganate

0.002 mol/I 0.01N / DILUT-IT

NC CODE 2841 61 00 R: 52/53

S: 61

PRODUCT CONT. **PACKING** BOX NO. 4873 1 amp.

Volumetric Concentrate, for dilution to 1 l.

Potassium Peroxodisulfate

For use in Kjeldahl Nitrogen, Phosphate and Mercury determination. / ACS

0377

4873

K₂S₂O₈

M = 270.33 g/mol

CAS NO. 7727-21-1

EINECS 231-781-8

NC CODE 2833 40 00

UN/ID NO. 1492 ADR/RID 5.1 O2

IMDG 5.1/III

R: 22-36/37/38-42/43-8

S: 22-24-26-37



Meets ACS Specifications. Meets Reagent Specifications for testing USP/NF monographs Assav min 99.0% max. 0.001% Chlorine Compounds (as CI) Heavy Metals (as Pb) max. 0.001% max. 0.005% Insoluble Matter Nitrogen Compounds (as N) max. 0.001% Trace Impurities (in ppm): Iron (Fe) max. 5

Manganese (Mn) max. 2 Mercury (Hg) act.value reported Phosphate (PO₄) max. 0.5

PRODUCT NO.	PACKING	CONT. BOX
0377.0500	500 g	6
0377.2500	2.5 kg	

Potas

Potassium Peroxodisulfate


BAKER ANALYZED

► K₂S₂O₈ **M** = 270.33 g/mol **CAS NO.** 7727-21-1 **EINECS** 231-781-8 NC CODE 2833 40 00 **UN/ID NO.** 1492

ADR/RID 5.1 O2 IMDG 5.1/III

> **R:** 22-36/37/38-42/43-8 **S:** 22-24-26-37





вох
6
4



Potassium Persulfate



Potassium Phosphate, Dibasic



Potassium Phosphate, Dibasic

See for qualities meeting the requirements of cGMP, USP/NF, BP/Ph.Eur, JP or FCC, page 36

Potassium Phosphate n-Hydrate

0243 'BAKER ANALYZED'

K₃PO₄.nH₂O **M** = 212.56 g/mol **CAS NO.** 7778-53-2 **EINECS** 231-907-1 NC CODE 2835 24 00 R: 36/38

Chloride (CI) max. 0.005% Dibasic Salt (as K₂HPO₄) max. 1.5% Excess Alkali (as KOH) max. 0.5% Heavy Metals (as Pb) max. 0.001% Insoluble Matter max. 0.01% max. 0.001% Iron (Fe) Nitrogen Compounds (as N) max. 0.002% Sulfate (SO₄) max. 0.005% Trace Impurities (in ppm): Arsenic (As) max. 1

PRODUCT **PACKING** CONT. NO. BOX 0243.0500 500 a

Potassium Phosphate, Monobasic



Potassium Phosphate, Monobasic

ee for qualities meeting the requirements of cGMP, USP/NF, BP/Ph.Eur, JP or FCC, page 36



Potassium Phosphate, Tribasic

See Potassium Phosphate n-Hydrat



Potassium Platinum(IV) Chloride

ee Potassium Hexachloroplatinate(IV)

Are you looking for a specialist in solvent applications?

See our solvent application program in chapter 3 of this catalogue.

0378

Potassium Pyrosulfate 'BAKER INSTRA-ANALYZED' / ACS

► K₂O₇S₂ **M** = 254.32 g/mol **CAS NO.** 7790-62-7 **EINECS** 232-216-8

NC CODE 2833 29 90

Meets ACS Specifications. Meets Reagent Specifications for testing USP/NF monographs Acidity (as H₂SO₄) 37.5-38.6% max. 0.002% Calcium (Ca) max. 0.001% Chloride (CI) Heavy Metals (as Pb) max. 0.001% Insoluble Matter max. 0.01% Iron (Fe) max. 0.002% Magnesium (Mg) max. 0.001% NH4OH ppt max. 0.005% max. 0.001% Phosphate (PO₄) max. 0.01% Sodium (Na) Suitability as Acid Flux passes test Water (H₂O) max. 2.5% Trace Impurities (in ppm):

PACKING	CONT.
	вох
500 g	

Flux Grade.

max. 1

Potassium Pyrosulfate

'BAKER ANALYZED' / ACS

Arsenic (As)

CAS NO. 7790-62-7 **EINECS** 232-216-8 NC CODE 2833 29 90

M = 254.32 g/mol

K₂S₂O₇

Exceeds ACS Specifications	
Acidity (as H ₂ SO ₄)	37.5-38.6%
Arsenic (As)	max. 0.0005%
Calcium (Ca)	max. 0.002%
Chloride (CI)	max. 0.001%
Heavy Metals (as Pb)	max. 0.001%
Insoluble Matter	max. 0.01%
Iron (Fe)	max. 0.002%
Magnesium (Mg)	max. 0.0005%
Phosphate (PO ₄)	max. 0.001%
Sodium (Na)	max. 0.01%
Water (H ₂ O)	max. 2.5%

PRODUCT	PACKING	CONT.
NO.		вох
0245 0100	100 a	

Potassium Pyrosulfite

See Potassium Disulfite



0246

1971

0245

Potassium Rhodanide

See Potassium Thiocyanate

Potassium Sodium Tartrate 4aq

'BAKER ANALYZED' / ACS

KNaC₄H₄O₆.4H₂O **M** = 282.23 g/mol **CAS NO.** 6381-59-5 **EINECS** 206-156-8 NC CODE 2918 13 00

Exceeds ACS Specifications	
Assay	99.0-102.0%
Ammonium (NH ₄)	max. 0.002%
Calcium (Ca)	max. 0.005%
Chloride (CI)	max. 0.001%
Insoluble Matter	max. 0.005%
pH of 5% Solution at 25°C	6.0-8.5
Phosphate (PO ₄)	max. 0.001%
Sulfate (SO ₄)	max. 0.005%
Trace Impurities:	
Heavy Metals (as Pb)	max. 5
Iron (Fo)	may 5

NO.	PACKING	BOX
0246.1000	1 kg	6
0246.9050	50 kg	
0246.9050	50 kg	

Potassium Sodium Tartrate 4aq

KNaC₄H₄O₆.4H₂O

M = 282.23 g/mol CAS NO. 6381-59-5 **EINECS** 206-156-8 NC CODE 2918 13 00

Assay min. 99.0% Appearance passes test Chloride (CI) max. 0.002% Sulfate (SO₄) max. 0.005% Trace Impurities: Heavy Metals (as Pb) max. 0.002% Iron (Fe) max. 0.002%

PRODUCT	PACKING	CONT.
NO.		вох
1971.1000	1 kg	

Potas



Potassium Sodium Tartrate, 4-Hydrate

e for qualities meeting the requirements of cGMP, USP/NF, BP/Ph.Eur, JP or FCC, page 36



Potassium Sorbate

ee for qualities meeting the requirements of cGMP, USP/NF, BP/Ph.Eur, JP or FCC, page 36

Potassium Sulfate

crystal / 'BAKER ANALYZED' / ACS

► K₂SO₄

M = 174.27 g/mol **CAS NO.** 7778-80-5

EINECS 231-915-5 NC CODE 3104 30 00

Meets ACS Specifications min. 99.0% max. 0.01% Calcium (Ca) Chloride (CI) max. 0.001% Insoluble Matter max. 0.01% Magnesium (Mg) max. 0.005% 5.5-8.5 pH of 5% Solution at 25°C Sodium (Na) max. 0.02%

Trace Impurities (in ppm): Heavy Metals (as Pb) max. 5 Iron (Fe) max. 5 Nitrogen Compounds (as N) max. 5

PRODUCT CONT. **PACKING** BOX 0247.1000 1 kg


Potassium Sulfate

0512 crystal / 'BAKER'

► K₂SO₄

M = 174.27 g/mol **CAS NO.** 7778-80-5

EINECS 231-915-5 NC CODE 3104 30 00





1169

Potassium Sulfate

See for qualities meeting the requirements of cGMP, USP/NF, BP/Ph.Eur, JP or FCC, page 36

Assav

Appearance

Potassium Tetroxoiodate(VII)

crystal / 'BAKER ANALYZED' / ACS

► KIO₄

M = 230.00 g/mol CAS NO. 7790-21-8

EINECS 232-196-0

NC CODE 2829 90 80

UN/ID NO. 1479 **ADR/RID** 5.1 O2

IMDG 5.1/I

R: 8

S: 24/25



Exceeds ACS Specifications Assay (dried basis) 99.8-100.3% max. 0.01% Other Halogens (as CI) Trace Impurities (in ppm): Manganese (Mn) max. 0.5

1169,9050 50 kg

PACKING

PRODUCT

PRODUCT

NO. BOX

CONT.

CONT.

Potassium Thiocyanate

'BAKER ANALYZED' / AC

► KSCN

M = 97.18 g/mol

CAS NO. 333-20-0 **EINECS** 206-370-1

NC CODE 2838 00 00 **EC NO.** 615 004 00 3

R: 20/21/22-32

S: 13



Exceeds ACS Specifications	
Assay	min. 98.5%
Ammonium (NH ₄)	max. 0.003%
Appearance	passes test
Chloride (CI)	max. 0.005%
Heavy Metals (as Pb)	max. 5 ppm
Insoluble in Water	max. 0.005%
lodine-consuming Substances (meq/g)	max. 0.004
Iron (Fe)	max. 2 ppm
pH of 5% Solution at 25°C	5.3-8.7
Sodium (Na)	max. 0.005%
Sulfate (SO ₄)	max. 0.005%

вох 0250.0100 100 g 0250.1000 1 kg

PACKING



Primary Drinking Water Standard I (Matrix: 2% nitric acid) / 'BAKER INSTRA-ANALYZED' / Plasma Standard

6021-01

NC CODE	
R:	36/38
S:	26-37
	X

Elemental Concent (µg/ml):	
Arsenic (As)	10
Barium (Ba)	100
Cadmium (Cd)	5
Chromium (Cr)	10
Lead (Pb)	10
Selenium (Se)	5
Silver (Ag)	10

PRODUCT	PACKING	CONT
NO.		B0)
6021-01	100 ml	

For use in EPA SW-846 Methods 6010 and 200.7. Traceable to NIST.

Primary Standards



L-Proline

See for qualities meeting the requirements of cGMP, USP/NF, BP/Ph.Eur, JP or FCC, page 36



Propanal



Propanedioic Acid



1,2-Propanediol

See Propylene Glycol



1,2,3-Propanetriol



Propanoic Acid

See Propionic Acid



1-Propanol

8066 BAKER ANALYZED



M = 60.10 g/mol

11= 0.80 kg

FLASHPOINT 15 °C

CAS NO. 71-23-8

EINECS 200-746-9

NC CODE 2905 12 00

EC NO. 603 003 00 0

UN/ID NO. 1274

ADR/RID 3 F1

IMDG 3/II

R: 11-41-67 **S:** 16-24-26-39-7





Assay (by GC)	min. 99%
2-Propanol (CH ₃ CHOHCH ₃) (by GC)	max. 0.05%
Boiling Range	max. 2.0°C
Color (APHA)	max. 10
Density (g/ml) at 25°C	0.800-0.804
Ethanol (C ₂ H ₅ OH) (by GC)	max. 0.05%
Free Acid (as C ₂ H ₅ COOH)	max. 0.001%
Methanol (by GC) (CH ₃ OH)	max. 0.05%
Recorded Boiling Point	97.2°C
Residue after Evaporation	max. 0.005%
Substances Darkened by H ₂ SO ₄	passes test
Substances Reducing KMnO ₄ (as O)	max. 0.0003%
Water (H ₂ O)	max. 0.1%
Trace Impurities (in ppm):	
Aluminium (AI)	max. 0.5
Barium (Ba)	max. 0.1
Boron (B)	may 0.02

Substances Darkened by H ₂ SO ₄	passes test
Substances Reducing KMnO ₄ (as O)	max. 0.0003%
Water (H ₂ O)	max. 0.1%
Trace Impurities (in ppm):	
Aluminium (AI)	max. 0.5
Barium (Ba)	max. 0.1
Boron (B)	max. 0.02
Cadmium (Cd)	max. 0.05
Calcium (Ca)	max. 0.5
Chromium (Cr)	max. 0.02
Cobalt (Co)	max. 0.02
Copper (Cu)	max. 0.02
Iron (Fe)	max. 0.1
Lead (Pb)	max. 0.1
Magnesium (Mg)	max. 0.1
Manganese (Mn)	max. 0.02
Nickel (Ni)	max. 0.02
Tin (Sn)	max. 0.1
Zinc (Zn)	max. 0.1

PRODUCT	PACKING	CONT.
NO.		вох
8066.1000	11	6
8066.2500	2.5	4



Propa

2-Propanol

9334 'BAKER ULTRA RESI-ANALYZED' / for Organic Residue Analysis

► CH₃CHOHCH₃

M = 60.10 g/mol

1 l = 0.78 kg

FLASHPOINT 12 $^{\circ}\mathrm{C}$

CAS NO. 67-63-0

EINECS 200-661-7

NC CODE 2905 12 00

EC NO. 603 003 00 4

UN/ID NO. 1219

ADR/RID 3F1

IMDG 3/II

R: 11-36-37

S: 16-24/25-26-7





Assay (by GC) (corrected for water)	min. 99.7%
Color (APHA)	max. 10
Residue after Evaporation	max. 1 ppm
Water (H ₂ O)	max. 0.1%

ECD Sensitive Impurities (as Heptachlor Epoxide): Single Impurity Peak (pg/ml) max. 1

FID-Sensitive Impurities (as 2-Octanol):

Single Impurity Peak (ng/ml) max. 5

PRODUCT	PACKING	CONT.
NO.		вох
9334.1000	11	6
9334.4000	4 Glass	

9830

2-Propanol

BAKER ANALYZED LC-MS Reagent

► CH₃CHOHCH₃

M = 60.10 g/mol

1 I = 0.78 kg

FLASHPOINT 12 $^{\circ}\mathrm{C}$

CAS NO. 67-63-0

EINECS 200-661-7 **NC CODE** 2905 12 00

EC NO. 603 117 00 0

UN/ID NO. 1219

ADR/RID 3F1

IMDG 3/II

R: 11-36-67 **S:** 16-24/25-26-7







Certificate Provided Reporting Ac	tual Lot Analysis
Assay (by GC)	min. 99.5%
Residue after Evaporation	max. 1 ppm
Water (H ₂ O)	max. 0.02%
LC-Gradient-Diode Array Detectio	n (a.u.), test
solution is modified with 0.1% (v/s	v) formic acid:
at 220 nm	max. 0.005
at 254 nm	max. 0.001
LC-MS Gradient Suitability Test (1	TIC, 100 to 2000
m/z), test solution is modified with	h 0.1% (v/v)
formic acid:	
Positive ESI-MS Sensitive Impurities (a	as
Reserpine)	max. 50 ng/m
Product Information (not specifica	ntions):
Density (g/ml) at 20°C	0.78
Trace Impurities (in ppb):	
Aluminium (AI)	max. 50

 Aluminium (AI)
 max. 50

 Calcium (Ca)
 max. 50

 Iron (Fe)
 max. 50

 Magnesium (Mg)
 max. 50

 Potassium (K)
 max. 50

 Sodium (Na)
 max. 50

 Ultraviolet Absorbance (1.00-cm path vs water):
 at 220 nm

 at 254 nm
 max. 0.30

 at 254 nm
 max. 0.02

max. 0.01

 PRODUCT
 PACKING
 CONT.

 NO.
 BOX

 9830.1000GL
 1 | Glass
 6

Element concentrations are at time of lot release.

2-Propanol

8175 'BAKER HPLC ANALYZED' / for use in High Performance Liquid Chromatography

at 280 nm

► CH₃CHOHCH₃

M = 60.10 g/mol

1 l = 0.78 kg

FLASHPOINT 12 °C **CAS NO.** 67-63-0

EINECS 200-661-7

NC CODE 2905 12 00

EC NO. 603 117 00 0

UN/ID NO. 1219

ADR/RID 3 F1

IMDG 3/II

R: 11-36-67

S: 16-24/25-26-7





Assay (by GC) (corrected for water)	mın. 99./%
Residue after Evaporation (in ppm)	max. 2
Water (H ₂ O)	max. 0.03%
Product Information (not specification	ıs):
Density (g/ml) at 20°C	0.785
Ultraviolet Absorbance (1.00-cm path	vs water):
at 225 nm	max. 0.16
at 254 nm	max. 0.02
at 280 nm	max. 0.01
at 350 nm	max. 0.01
UV Cut-off, nm	max. 205

PRODUCT	PACKING	CONT.
NO.		вох
8175.1000	11	6
8175.2500	2.5	4
8175.5000	5 EcoTainer	
EcoTainer, th	e metal solvent can for more	safety in the lab.

Filtered through a 0.2 micron filter. Packaged under Nitrogen.



The reagents program for histopathology contains reagents for histology and cytology.

Major fields of application for histology and cytology:

Hospital

We offer a complete product line:

Histology

- Fixatives: formaldehyde buffered: pH 7.0 and 40mM phosphate
- Waste disposal: formalin neutraliser to neutralise formaldehyde waste
- Xylene/toluene replacement: UltraClear, an iso-paraffin clearing agent
- Paraffin: UltraPar 54-56°C, our top quality paraffin
- Mounting medium: UltraKitt, our xylene/toluene free cover slipping reagent. Patent for this product is granted and filed
- H&E stains: hematoxyline, is available according to Mayer or modified according to Gill II / Harris. In addition we also offer a Scotch solution, which enables a stable and optimised bluing of the cell nucleus. Eosine Y stains are offered as an alcoholic

- or as an aqueous solution. Our Eosine Y counter stains are optimised in combination with our hematoxyline stains
- Giemsa: we offer a high quality Giemsa for staining of tissue sections, such as gastric sections
- Solvents: wide range of alcohols, but also Histo grade xylene and toluene

Cytology

- Fixatives: Cervix Spray fixative
- Papanicolaou staining: we offer Papanicolaou 1
 (hematoxyline according to Gill II / Harris),
 Papanicolaou 2A (Orange G6), Papanicolaou 2B
 (Orange II) and Papanicolaou 3B (EA 50)
- May Grünwald and Giemsa, Wright and Leishman:
 we offer a broad range of eosine-methylene blue
 stains for cytology smears and bone marrow
- Other products: see also UltraClear, UltraKitt and our Histo grade Solvent range

For more details and recommended staining procedures, product information bulletins are available at our website www.jtbaker.com



Propa

2-Propanol

8067

'BAKER ANALYZED' / ACS

► CH₃CHOHCH₃

M = 60.10 g/mol **1 l** = 0.78 kg

FLASHPOINT 12 $^{\circ}\mathrm{C}$

CAS NO. 67-63-0

EINECS 200-661-7

NC CODE 2905 12 00 EC NO. 603 117 00 0

UN/ID NO. 1219

ADR/RID 3F1

IMDG 3/II

R: 11-36-67

S: 16-24/25-26-7





Assay	min. 99.5%
Color (APHA)	max. 10
Residue after Evaporation	max. 0.001%
Solubility in Water	passes tes
Titrable Acid or Base (meq/g)	max. 0.000
Water (H ₂ O)	max. 0.1%
Trace Impurities (in ppm):	
Aluminium (AI)	max. 0.5
Barium (Ba)	max. 0.
Boron (B)	max. 0.02
Cadmium (Cd)	max. 0.09
Calcium (Ca)	max. 0.
Chromium (Cr)	max. 0.02
Cobalt (Co)	max. 0.02
Copper (Cu)	max. 0.02
Iron (Fe)	max. 0.
Lead (Pb)	max. 0.
Magnesium (Mg)	max. 0.
Manganese (Mn)	max. 0.02
Nickel (Ni)	max. 0.02
Tin (Sn)	max. 0.
Zinc (Zn)	max. 0.

Exceeds ACS Specifications

PRODUCT	PACKING	CONT.
NO.		вох
8067.1000	11	6
8067.2500	2.5	4
8067.5000	5 EcoTainer	4
8067.9020	20	
8067.9025	25	
8067.9200	200	
EcoTainer, th	e metal solvent can for more saf	ety in the lab.

For safe handling of 25 l tin cans, see Self-closing tap.

2-Propanol

'BAKER ANALYZED' / Ultraviolet Spectrophotometry / ACS

at 300 nm

at 330-400 nm

► CH₃CHOHCH₃

8068

M = 60.10 g/mol

1 l = 0.78 kg

FLASHPOINT 12 $^{\circ}\mathrm{C}$

CAS NO. 67-63-0 **EINECS** 200-661-7

NC CODE 2905 12 00

EC NO. 603 117 00 0

UN/ID NO. 1219

ADR/RID 3 F1

IMDG 3/II

R: 11-36-67 **S:** 16-24/25-26-7





Exceeds ACS Specifications Assay min. 99.5% Color (APHA) max. 10 Residue after Evaporation max. 5 ppm Solubility in Water passes test Titrable Acid or Base (meq/g) max. 0.0001 Water (H₂O) max. 0.05% Ultraviolet Absorbance (1.00-cm path vs water): max. 1.00 at 210 nm at 220 nm max. 0.40 at 230 nm max. 0.20 at 245 nm max. 0.08 at 260 nm max. 0.04 at 275 nm max. 0.03

max. 0.02

max. 0.01

PRODUCT	PACKING	CONT.
NO.		вох
8068.1000	11	6
8068.2500	2.5	4

2-Propanol 'BAKER

8119

► CH₃CHOHCH₃

M = 60.10 g/mol

1 l = 0.78 kg **FLASHPOINT** 12 °C

CAS NO. 67-63-0

EINECS 200-661-7 NC CODE 2905 12 00

EC NO. 603 117 00 0

UN/ID NO. 1219

ADR/RID 3 F1

IMDG 3/II

R: 11-36-67

S: 16-24/25-26-7





Assay (by GC)	min. 99%
Acidity (as CH ₃ COOH)	max. 0.002%
Alkalinity (as NH ₃)	max. 0.002%
Boiling Range	max. 1.0°C
Density (g/ml) at 25°C	0.781-0.783
Recorded Boiling Point	82°C
Residue after Evaporation	max. 0.005%
Water (H ₂ O)	max. 0.2%

PRODUCI	PACKING	CONT.
NO.		вох
8119.1000	11	6
8119.2500	2.5	4
8119.5000	5 EcoTainer	4
8119.9025	25 l	
8119.9200	200 l	
EcoTainer, th	e metal solvent can for more	e safety in the lab.
- / /	00 710	





The four cornerstones of doing business with us: Innovation, Collaboration, Support and Quality.



Propi

Propionic Acid 6034 'BAKER ANAI YZED

Assay (acidimetric)

Density (g/ml) at 25°C

Identification (by IR)

Residue after Evaporation

Boiling Range

Water (H₂O)

► CH₃CH₂COOH

M = 74.08 g/mol

1 l = 0.99 kg

 $\textbf{FLASHPOINT} \ \ 54\ ^{\circ}\text{C}$

CAS NO. 79-09-4

EINECS 201-176-3

NC CODE 2915 50 00 EC NO. 607 089 00 0

UN/ID NO. 1848

ADR/RID 8 C3

IMDG 8/III

R: 34

S: 23-36-45





iso-Propyl Alcohol



n-Propyl Alcohol



iso-Propyl Ether

Propylene Carbonate

7439 ▶ QCH(CH₃)CH₂OÇO

M = 102.09 g/mol

11 = 1.20 kg

FLASHPOINT 135 °C

CAS NO. 108-32-7

EINECS 203-572-1

NC CODE 2920 90 10

EC NO. 607 194 00 1

R: 36

S: 26



Propylene Glycol

7064

► CH₃CHOHCH₂OH

M = 76.10 g/mol **1 l** = 1.03 kg

FLASHPOINT 99 °C

CAS NO. 57-55-6

NC CODE 2905 32 00

EINECS 200-338-0

max. 0.005% Acidity (as CH₃COOH) Boiling Range 185-189°C Chloride (CI) max. 0.001% Density (g/ml) at 25°C 1 032-1 034 Residue after Ignition max. 0.005% Water (H₂O) max. 0.2% Trace Impurities (in ppm): Arsenic (As) max. 1

7064.0100 100 ml 7064 1000 11

PACKING

PRODUCT

max. 5

PRODUCT

6034.0500

PACKING

PACKING

11

7439.1000

min. 99.5%

139-142°C

0.984-0.989

max. 0.1%

passes test

max. 0.005%

CONT.

вох

BOX

CONT.

Innovation is principal to our business.

Heavy Metals (as Pb)

Propylene Glycol

7169

min. 99%

33-35°C

max. 10

max. 0.01%

PRODUCT

8822.0100

8822.0500

8822.1000

8822.9200

PACKING

100 ml

500 ml

11

200 I



M = 76.10 g/mol

1 l = 1.03 kg

 $\textbf{FLASHPOINT} \ 99\ ^{\circ}\text{C}$

CAS NO. 57-55-6 **EINECS** 200-338-0

NC CODE 2905 32 00

Assay	min. 99.5%
Acidity	passes test
Appearance of solution	passes test
Chloride (CI)	max. 0.007%
Density (g/ml) at 20°C	1.035-1.040
Identification	passes test
Organic Volatile Impurities	passes test
Oxidizing substances	passes test
Reducing Substances	passes test
Refractive Index at 25°C, n ²⁵	1.431-1.433
Residue on Ignition (as SO ₄)	max. 0.007%
Specific Gravity 25°C/25°C	1.035-1.037
Sulfate (SO ₄)	max. 0.006%
Water (H ₂ O)	max. 0.2%
Trace Impurities (in ppm):	
Heavy Metals (as Pb)	max. 5 ppm

PRODUCT NO.	PACKING	CONT. BOX
7169.2500	2.5	
7169.9025	25 l	
7169.9200	200 l	

Propylene Glycol

qualities meeting the requirements of cGMP, USP/NF, BP/Ph.Eur, JP or FCC, page 36

Assay (by GC)

Boiling Point

Color (APHA)

Water (H₂O)



8822

CONT.

BOX

Propylene Oxide

M = 58.08 g/mol 11= 0.83 kg

FLASHPOINT −37 °C

CAS NO. 75-56-9

EINECS 200-879-2

NC CODE 2910 20 00

EC NO. 603 055 00 4

UN/ID NO. 1280

► CH₃ÇHCH₂Q

ADR/RID 3 F1 **IMDG** 3/I

R: 12-20/21/22-36/37/38-45-46

S: 45-53





Propylparaben

alities meeting the requirements of cGMP, USP/NF, BP/Ph.Eur, JP or FCC, page 36



PRS-1000 Positive Resist Stripper



PRS-2000 Positive Resist Stripper

nic materials page 32 or the website www.jtbaker.com/micro



PRS-3000 Positive Resist Stripper

 $See for \ detailed \ information \ section \ Microelectronic \ materials \ page \ 32 \ or \ the \ website \ www.jtbaker.com/microelectronic \ materials \ page \ 32 \ or \ the \ website \ www.jtbaker.com/microelectronic \ materials \ page \ 32 \ or \ the \ website \ www.jtbaker.com/microelectronic \ materials \ page \ 32 \ or \ the \ website \ www.jtbaker.com/microelectronic \ materials \ page \ 32 \ or \ the \ website \ www.jtbaker.com/microelectronic \ materials \ page \ 32 \ or \ the \ website \ www.jtbaker.com/microelectronic \ materials \ page \ 32 \ or \ the \ website \ www.jtbaker.com/microelectronic \ materials \ page \ 32 \ or \ the \ website \ www.jtbaker.com/microelectronic \ materials \ page \ 32 \ or \ the \ website \ www.jtbaker.com/microelectronic \ page \ 32 \ or \ the \ website \ www.jtbaker.com/microelectronic \ page \ 32 \ or \ the \ website \ page \ 32 \ or \ the \ website \ page \ 32 \ or \ the \ page \ 32$



Pumice Stone

Water Soluble Substances



granular 1.5 mm / 'BAKER'

NC CODE 2513 19 00 Acid Soluble Substances max. 6.0% **PRODUCT PACKING** CONT. Iron (Fe) BOX passes test NO.

max. 0.20%

1926.1000

1 kg



1927 granular 3 mm / 'BAKER' NC CODE 2513 19 00 Acid Soluble Substances max. 6.0% PACKING CONT. PRODUCT Iron (Fe) вох passes test NO. Water Soluble Substances max. 0.20% 1927.1000 1 kg

Purge

9393

Purge and Trap Solvents

See for detailed information section Solvent applications, page 2

Pyridine

Low Water / 'BAKER HPLC ANALYZED' / for Use in Liquid Chromatography and Spectrophotometry

at 350 nm at 375 nm

at 400 nm

N:CHCH:CHCH:CH

M = 79.10 g/mol

1 I = 0.98 kg

FLASHPOINT 17 °C

CAS NO. 110-86-1

EINECS 203-809-9
NC CODE 2933 31 00
EC NO. 613 002 00 7

UN/ID NO. 1282 ADR/RID 3 F1 IMDG 3/II

R: 11-20/21/22

S: 26-28

Xn armful



Assay (by GC) (corrected for water)	min. 99.0%
Ammonia (NH ₃)	max. 0.002%
Color (APHA)	max. 10
Reducing Substances	passes test
Residue after Evaporation	max. 5 ppm
Solubility in Water	passes test
Sulfate (SO ₄)	max. 0.001%
Water (by KF, volumetric)	max. 0.01%
Trace Impurities (in ppm):	
Chloride (CI)	max. 5
Copper (Cu)	max. 5
Ultraviolet Absorbance (1.00-cm pat	h vs water):
at 330 nm	max. 1.00
at 340 nm	max 0.10

max. 0.01

max. 0.005

max. 0.1

max. 0.02

max. 0.02

max. 0.1

max. 1.0

)	PRODUCT	PACKING	CONT. BOX
	9393.1000	11	
t			
1			
t			
)			
)			
5			
))			
)			
)			

8073

Pyridine

'BAKER ANALYZED' / ACS

► C₅H₅N

M = 79.10 g/mol **1** l = 0.98 kg

FLASHPOINT 17 °C CAS NO. 110-86-1 EINECS 203-809-9

NC CODE 2933 31 00 EC NO. 613 002 00 7

UN/ID NO. 1282 **ADR/RID** 3 F1 **IMDG** 3/II

R: 11-20/21/22

S: 26-28





Exceeds ACS Specifications	
Assay (by GC)	min. 99.0%
Ammonia (NH ₃)	max. 0.002%
Chloride (CI)	max. 0.001%
Reducing Substances	passes test
Residue after Evaporation	max. 0.002%
Solubility in Water	passes test
Suitability for OH groups determ. by	
phthalisation	passes test
Sulfate (SO ₄)	max. 0.001%
Water (H ₂ O)	max. 0.1%
Trace Impurities (in ppm):	
Aluminium (AI)	max. 0.5
Barium (Ba)	max. 0.1
Boron (B)	max. 0.1
Cadmium (Cd)	max. 0.05
Calcium (Ca)	max. 0.5
Chromium (Cr)	max. 0.02
Cobalt (Co)	max. 0.02
Copper (Cu)	max. 0.02
Iron (Fe)	max. 0.1
Lead (Pb)	max. 0.1

PRODUCT	PACKING	CONT
NO.		B0)
8073.0500	500 ml	(
8073.1000	11	(
8073.2500	2.5	4

EcoTainer, the metal solvent can for more safety in the lab.

Mallinckrodt Baker's chemistry is Part of a pure process $^{\text{TM}}$.

Magnesium (Mg)

Manganese (Mn)

Nickel (Ni)

Tin (Sn)

Zinc (Zn)

8074

CONT.

вох

Pyridine

► C₅H₅N

M = 79.10 g/mol **1 l** = 0.98 kg

FLASHPOINT 17 °C

CAS NO. 110-86-1 **EINECS** 203-809-9

NC CODE 2933 31 00 EC NO. 613 002 00 7

UN/ID NO. 1282 ADR/RID 3F1

IMDG 3/II

R: 11-20/21/22

S: 26-28





Pyridoxine Hydrochloride

e for qualities meeting the requirements of cGMP, USP/NF, BP/Ph.Eur, JP or FCC, page 36

Chloride (CI)

Water (H₂O)

Residue after Evaporation



1-(2-Pyridylazo)-2-naphthol



Pyrogallol

'BAKER ANALYZED' / ACS

1171

► 1,2,3-(HO)₃C₆H₃

M = 126.11 g/mol

CAS NO. 87-66-1

EINECS 201-762-9

NC CODE 2907 29 00

EC NO. 604 009 00 6

UN/ID NO. 2811

ADR/RID 6.1 T2

IMDG 6.1/III

R: 20/21/22-40-52/53

S: 36/37-61



Meets ACS Specifications	
Chloride (CI)	max. 0.001%
Iron (Fe)	max. 0.001%
Melting Point	131.0-135.0°C
Residue after Ignition	max. 0.005%
Sulfate (SO ₄)	max. 0.005%
Trace Impurities (in ppm):	

Heavy Metals (as Pb) max. 5 **PRODUCT PACKING** CONT. NO. вох 1171.0100 100 g 1171.9050 50 kg

Pyruvic Acid, Sodium Salt

e Sodium Pyruvate



1983

1984

BOX

Quartzsand

NC CODE 2506 10 00	Chloride (CI)	max. 0.01%
	Iron (Fe)	max. 0.02%
	Loss after Ignition	max. 0.2%

Soluble in HCI

ax. 0.01%	PR
ax. 0.02%	NO
nax. 0.2%	19
nax 0.5%	

max. 0.005%

max. 0.005%

max. 0.25%

PRODUCT

8074.2500

8074.9025

NO.

PACKING

2.5 l

25 I

RODUCT **PACKING** 983.1000 1 kg

Quartzsand

Granular 1-2 mm / 'BAKER'

Chloride (CI)	max. 0.01%
Iron (Fe)	max. 0.01%
Loss after Ignition	max. 0.5%
Soluble in HCl	max. 0.05%

CONT.	
BOX	

Quinol

NC CODE 2506 10 00

See Hydroquinone



RBS35

Reagent for silicic acid determination

7488

CAS NO. 7664-93-9 **EINECS** 231-163-9 NC CODE 3822 00 00 EC NO. 16 020 00 8

UN/ID NO. 2796 ADR/RID 8 C1 IMDG 8/II

R: 36/38 **S:** 26

PRODUCT PACKING CONT. BOX 7488.9020 20 l Polycube

Volumetric Solution, ready for use. Contains: Silicate-free water, Sulfuric acid, According to Bran and Lübbe.

Reagent for silicic acid determination

7534 'BAKER ANALYZED' / Reduction solution (Reagent 294)

NC CODE 3822 00 00 R: 43-52/53

S: 24-26-37-61

PRODUCT PACKING CONT. BOX NO. 7534.9020 20 I

Contains: Silicate-free water, 4-(Methylamino)phenolsulfate, Potassium Disulfite.

Reagent for silicic acid determination

7535 'BAKER ANALYZED' / Molybdate solution (Reagent 293)

NC CODE 3822 00 00

PRODUCT PACKING CONT. вох NO. 7535.9020 20 l Polycube

Contains: Silicate-free water, Ammonium heptamolybdate tetrahydrate, Silicate-free Ammonium Hydroxide.

Reinecke Salt 'BAKER ANALYZED' / ACS

1388

► NH₄[Cr(SCN)₄(NH₃)₂].H₂O **M** = 354.44 g/mol **CAS NO.** 13573-16-5 **EINECS** 237-003-3 NC CODE 2842 90 90

Meets ACS Specifications Assay min. 93.0% Insoluble in Dilute HCI max. 0.05% Sensitivity passes test

PRODUCT PACKING CONT. вох 1388.0025 25 g Glass

Residue Analysis Solvents

nation section Solvent applications, page 22



Resorcinol

e for qualities meeting the requirements of cGMP, USP/NF, BP/Ph.Eur, JP or FCC, page 36

REZI-38 Positive Resist Stripper

See for detailed information section Microelectronic materials page 32 or the website www.jtbaker.com/micro

0252

5776

6

4

Rhodamine B

1389 C₂₈H₃₁CIN₂O₃ PRODUCT PACKING CONT. Sensitivity as Indicator passes test **M** = 479.02 g/mol Solubility NO. вох passes test **CAS NO.** 81-88-9 1389.0100 100 g **EINECS** 201-383-9 NC CODE 3204 19 00 **R:** 41-52/53 **S:** 22-26-39-61 Saccharin Sodium See for qualities meeting the requirements of cGMP, USP/NF, BP/Ph.Eur, JP or FCC, page 36 **Safety Products** See for detailed information section Safety, page 400 Salicylaldehyde 7066 ► HOC₆H₄CHO min. 98.0% PRODUCT PACKING CONT. Assay (by GC) **M** = 122.12 g/mol Density (g/ml) at 20/4°C 1.164-1.166 вох **1 l** = 1.16 kg Freezing Point 1-2°C 7066.0100 100 ml FLASHPOINT 78 °C 7066 1000 11 **CAS NO.** 90-02-8 **EINECS** 201-961-0 NC CODE 2912 49 00 **UN/ID NO.** 2810 ADR/RID 6.1 T1 IMDG 6.1/III **R:** 21/22 **S:** 36/37 × **Salicylic Acid** See for qualities meeting the requirements of cGMP, USP/NF, BP/Ph.Eur, JP or FCC, page 36 Sand See Quartzsand

Sand

washed and ignited / 'BAKER'

CONT. NC CODE 2505 90 00 **PACKING** Appearance passes test **PRODUCT** вох NO. 0252.1000 1 kg 0252.5000 5 kg 0252.9050 50 kg

Scandium 1000 µg/ml

(Matrix: 2% nitric acid) / 'BAKER INSTRA-ANALYZED' / Plasma Standard

► Sc Certificate Provided Reporting Actual Lot Analysis **PRODUCT PACKING** CONT. **M** = 44.96 g/mol 998-1002 μg/ml Scandium (Sc) NO. BOX NC CODE 3822 00 00 5776.0100 100 ml **R:** 36/38 **S:** 26-37

Prepared from the highest purity raw material available, generally greater than 99.999% spectral purity. The content of the solution is confirmed to be accurate to within \pm 0.2%. Typically 1000 $\mu\text{g/ml.}$ The certificate of analysis provided reports actual lot analysis. The certificate also lists the trace impurities.

Scand

Scandium 10000 $\mu g/ml$ (Matrix: 7% nitric acid) / 'BAKER INSTRA-ANALYZED' / Plasma Standard 5742

► Sc

M = 44.96 g/mol

NC CODE 3822 00 00

UN/ID NO. 2031

ADR/RID 8 CO1

IMDG 8/II

R: 34

S: 26-36/39-45

Certificate Provided Reporting Actual Lot Analysis Scandium (Sc)

9980-10020 μg/ml

min. 97%

PRODUCT

6064.0100

PRODUCT CONT. **PACKING** NO. вох 5742.0100 100 ml

Prepared from the highest purity raw material available, generally greater than 99.999% spectral purity. The content of the solution is confirmed to be accurate to within 0.2 %. Typically 10000 $\mu\text{g/ml.}$

The certificate of analysis provided reports actual lot analysis. The certificate also lists the trace impurities.

Scotch Buffer

3872 HISTO GRADE / Bluing solution

NC CODE 382 20 00

Ready to use Bluing Solution for use in HE Staining

Assay (as Chlorine)

PRODUCT CONT. **PACKING** BOX

CONT.

BOX

PACKING

100 ml

Sebacoyl Chloride

► CICO(CH₂)₈COCI

6064

M = 239.15 g/mol

11 = 1.13 kg

FLASHPOINT 66 °C

CAS NO. 111-19-3

EINECS 203-843-4 NC CODE 2917 19 90

UN/ID NO. 3265

ADR/RID 8 C3

IMDG 8/II

R: 34

S: 26-36/37/39-45



Selenium

2036 Granular / 'BAKER ANALYZED'

M = 78.96 g/mol

CAS NO. 7782-49-2

EINECS 231-957-4

NC CODE 2804 90 00 EC NO. 34 001 00 2

UN/ID NO. 3283

ADR/RID 6.1 T5

IMDG 6.1/III

R: 23/25-33-53

S: 20/21-28-45-61

Heavy Metals (as Pb) max. 0.01% Iron (Fe) max. 0.05% Nitrogen (N) max. 0.1% Residue on Ignition max. 1.0%

PRODUCT PACKING CONT. BOX 2036.0100 100 g

www.jtbaker.com/europe



Selfc

Selenium 1000 μ g/ml

(Matrix: 2% nitric acid) / 'BAKER INSTRA-ANALYZED' / Plasma Standard

5777

•	Se	

M = 78.96 g/mol

NC CODE 3822 00 00

EC NO. 34 002 00 8

R: 36/37/38 **S:** 26-37

Certificate Provided Reporting Actua	al Lot Analysis
Colonium (Co)	000 1000 110/2

Selenium (Se)

PRODUCT PACKING CONT. NO. вох 5777.0100 100 ml

Prepared from the highest purity raw material available, generally greater than 99.999% spectral purity. The content of the solution is confirmed to be accurate to within \pm 0.2%. Typically 1000 $\mu\text{g/ml.}$

The certificate of analysis provided reports actual lot analysis. The certificate also lists the trace impurities.

PACKING

100 ml

numbers (SRM) are printed on each label.

PRODUCT

6938.0100


NO.

Selenium 1000 µg/ml (Matrix: 2% nitric acid) / 'BAKER INSTRA-ANALYZED' / Atomic Absorption Standard

6938

CONT.

BOX



Selenium (Se) 998-1002 μg/ml

6938.0500 500 ml Prepared by dissolution of high purity raw materials (min. 99.99% spectral purity). Assays are verified by ICP against standards traceable to NIST. Standard Reference Material

Selenium 10000 μg/ml (Matrix: 2% nitric acid) / 'BAKER INSTRA-ANALYZED' / Plasma Standard

5743



Certificate Provided Reporting Actual Lot Analysis

Selenium (Se) 9980-10020 μg/ml

CONT. PRODUCT PACKING NO. вох 5743.0100 100 ml

M = 78.96 g/molNC CODE 3822 00 00 **EC NO.** 34 002 00 8 R: 33-36/38 **S:** 26-36

Prepared from the highest purity raw material available, generally greater than 99.999% spectral purity. The content of the solution is confirmed to be accurate to within 0.2 %. Typically 10000 $\mu g/ml$.

The certificate of analysis provided reports actual lot analysis. The certificate also lists the trace impurities.

Self-closing tap adaptor

4532

PRODUCT	PACKING	CONT.
NO.		вох
4532	1 unit	

Self-closing tap for flammable solvents in 10 l drums

4535

PRODUCT	PACKING	CONT.
NO.		вох
4535	1 unit	

Self-closing tap for flammable solvents in 25 l tin-cans

4531

PRODUCT	PACKING	CONT.
NO.		вох
4531	1 unit	

Self-closing tap for flammable solvents in 200 l drums

4529

PRODUCT	PACKING	CONT.
NO.		вох
4529	1 unit	

Serin

L-Serine

See for qualities meeting the requirements of cGMP, USP/NF, BP/Ph.Eur, JP or FCC, page 36

Shorr's stain for hormonal cytodiagnosis

Shorr Staining Solution

3876 Cytolog

1 I = 0.90 kgFLASHPOINT 13 °C

NC CODE 3822 00 00

UN/ID NO. 1992

ADR/RID 3 FT1 IMDG 3/II

R: 11-20/21/22-68/20/21/22

S: 16-33-36/37-7/9





Silica Gel

0253 0.063 - 0.200 mm / 'BAKER ANALYZED' / for Chromatography

EINECS 215-683-2	Alcohol-Ether Soluble Substances	max. 0.05%	PRODUCT	PACKING	CONT.
NC CODE 2811 22 00	Loss on Drying	max. 6.0%	NO.		BOX
			0253.0500	500 g	
			0253.1000	1 kg	
			0253.5000	5 kg	
			0253.9025	25 kg	

Silica Gel

with Indicator: Cobalt(II)chloride / 3-6 mm / 'BAKER'

EINECS	215-683-2	Appearance	passes test	PRODUCT	PACKING	CONT.
NC CODE	2811 22 00			NO.		BOX
R:	49			1827.1000	1 kg	6
S:	45-53					
	₩					

Silica Gel 7

5056 Powder / 'BAKER TLC' / for Thin Layer Chromatography

CAS NO. 1343-98-2	Recommended Slurry proporti	ons: 25g/in 60 ml H₂O	PRODUCT	PACKING	CONT.
EINECS 215-683-2	Chloride (CI)	max. 0.02%	NO.		ВОХ
NC CODE 2811 22 00	Heavy Metals (as Pb)	max. 0.004%	5056.0500	500 g	
	Iron (Fe)	max. 0.01%			
	Loss on Ignition	max. 20.0%			
	Nonvolatile with HF	max. 0.3%			
	Particle Size (< 40 µm)	min. 95%			
	pH (10% slurry)	7.0-8.0			
	Suitability for TLC	passes test			

Silica Gel IB-2

5003

'BAKER-FLEX' / Flexible TLC Sheets, 20 X 20 cm

	BOX
25 sheets	
	25 sheets

containing an inert binder.

The binder does not char on sulfuric acid-heat treatment.



Silica Gel IB2-F

'BAKER-FLEX' / Flexible TLC Sheets, 20 X 20 cm

5004

PRODUCT	PACKING	CONT.
NO.		вох
5004	25 sheets	

A flexible sheet coated with high purity silica gel containing an inert binder and a fluorescent indicator (activated at 2540 Å). Neithter the binder nor the indicator chars on sulfuric acid-heat treatment.

Silicic Acid

0324-01

$\sim SiO_2.nH_2O$	
CAS NO.	1343-98-2
EINECS	215-683-2
NC CODE	2811 22 00

Assay (as SiO ₂)	min. 84%
Chloride (CI)	max. 0.01%
Heavy Metals (as Pb)	max. 0.002%
Iron (Fe)	max. 0.002%
Loss on Ignition (as H ₂ O)	max. 16%
Nonvolatile with HF	max. 0.2%
Sulfate (SO ₄)	max 0.02%

RODUCT	PACKING	CONT.
10.		вох
324-01	500 g HDPE	

Silicon 1000 μ g/ml

s a trace of hydrofluoric acid) / 'BAKER INSTRA-ANALYZED' / Plasma Standard (Matrix: 2% nitric acid p

5778

▶ Si

M = 28.09 g/mol

1 l = 1.03 kg

NC CODE 3822 00 00

EC NO. 9 012 00 0

R: 36/38

S: 26-37



Certificate Provided Reporting Actual Lot Analysis Silicon (Si) 998-1002 μg/ml

PRODUCT PACKING CONT. вох NO. 5778.0100 100 ml

Prepared from the highest purity raw material available, generally greater than 99.999% spectral purity. The content of the solution is confirmed to be accurate to within \pm 0.2%. Typically 1000 $\mu\text{g/ml.}$

The certificate of analysis provided reports actual lot analysis. The certificate also lists the trace impurities

Silicon 1000 μg/ml

us a trace of hydrofluoric acid) / 'BAKER INSTRA-ANALYZED' / Atomic Absorption Standard (Matrix: 2% nitric acid p

6939

CONT.

► Si

M = 28.09 g/mol

11 = 1.00 kg

NC CODE 3822 00 00

EC NO. 9 012 00 0

R: 36/38

S: 26-37



Silicon (Si) $998\text{-}1002~\mu\text{g/ml}$

PRODUCT NO.

N 0

6939.0100 100 ml 6939.0500 500 ml

PACKING

Prepared by dissolution of high purity raw materials (min. 99.99% spectral purity). Assays are verified by ICP against standards traceable to NIST. Standard Reference Material numbers (SRM) are printed on each label.

Silicon 10000 µg/ml

a trace of hydrofluoric acid) / 'BAKER INSTRA-ANALYZED' / Plasma Standard

5744

► Si

M = 28.09 g/mol

NC CODE 3822 00 00

EC NO. 9 012 00 0

R: 36/38 **S:** 26-37



Certificate Provided Reporting Actual Lot Analysis 9980-10020 μg/ml Silicon (Si)

CONT. **PRODUCT** PACKING вох NO. 5744.0100 100 ml

Prepared from the highest purity raw material available, generally greater than 99.999% spectral purity. The content of the solution is confirmed to be accurate to within 0.2 %. Typically 10000 μ g/ml.

The certificate of analysis provided reports actual lot analysis. The certificate also lists the trace impurities.

Silicone Emulsion



Silve

Silver 1000 μ g/ml

(Matrix: 2% nitric acid) / 'BAKER INSTRA-ANALYZED' / Plasma Standard 5779

► Ag **M** = 107.87 g/mol

NC CODE 3822 00 00

R: 36/38

S: 26



Certificate Provided Reporting Actual Lot Analysis

Silver (Ag) 998-1002 μg/ml

998-1002 μg/ml

 $998\text{-}1002~\mu\text{g/ml}$

PRODUCT CONT. **PACKING** NO. вох 5779.0100 100 ml

Prepared from the highest purity raw material available, generally greater than 99.999% spectral purity. The content of the solution is confirmed to be accurate to within \pm 0.2%. Typically 1000 μ g/ml.

The certificate of analysis provided reports actual lot analysis. The certificate also lists the trace impurities.

Silver 1000 μ g/ml (Matrix: 2% nitric acid) / 'BAKER INSTRA-ANALYZED' / Atomic Absorption Standard 6940

Silver (Ag)

Silver (Ag)

► Ag **M** = 107.87 g/mol NC CODE 3822 00 00

R: 36/38

S: 26



PRODUCT PACKING CONT. NO. BOX 6940.0100 100 ml 6940.0500 500 ml

Prepared by dissolution of high purity raw materials (min. 99.99% spectral purity). Assays are verified by ICP against standards traceable to NIST. Standard Reference Material numbers (SRM) are printed on each label.

Silver 1000 μg/ml
'BAKER ANALYZED' / Atomic Absorption Standard 6821

► Ag

M = 107.87 g/mol NC CODE 3822 00 00

R: 36/38

S: 26-37



PRODUCT PACKING CONT. NO. вох 6821.0100 100 ml 6821.0500 500 ml

Silver nitrate in nitric acid 0.5 mol/l.

Silver 10000 μg/ml

(Matrix: 2% nitric acid) / 'BAKER INSTRA-ANALYZED' / Plasma Standard

M = 107.87 g/mol

► Ag

NC CODE 3822 00 00

R: 36/38

S: 26



Certificate Provided Reporting Actual Lot Analysis 9980-10020 μg/ml Silver (Ag)

PRODUCT PACKING CONT. вох NO. 5745.0100 100 ml

Prepared from the highest purity raw material available, generally greater than 99.999% spectral purity. The content of the solution is confirmed to be accurate to within 0.2 %. Typically 10000 µg/ml.

The certificate of analysis provided reports actual lot analysis. The certificate also lists the trace impurities.

Silver Chloride

BAKER ANALYZED

► AgCl **M** = 143.32 g/mol CAS NO. 7783-90-6 **EINECS** 232-033-3

NC CODE 2843 29 00


1180

Assay min. 99.0% Heavy Metals (as Pb) max. 0.005% Iron (Fe) max. 0.003% Nitrate (NO₃) max. 0.001% Trace Impurities (in ppm): Copper (Cu) max. 5

PACKING CONT. **PRODUCT** BOX 1180.0025 25 g Glass 1180.0100 100 q 1180.1000 1 kg

Silver Nitrate

'BAKER ANALYZED' / ACS



EC NO. 47 001 00 2 **UN/ID NO.** 1493 **ADR/RID** 5.1 O2 IMDG 5.1/II

> **R:** 34-50/53 **S:** 26-45-60-61

corrosive

Exceeds ACS Specifications	
Assay	min. 99.9%
Clarity of Solution	passes test
Free Acid	passes test
Lead (Pb)	max. 0.001%
Substances not Precipitated by HCI	max. 0.01%
Sulfate (SO ₄)	max. 0.002%
Trace Impurities (in ppm):	
Chloride (CI)	max. 5
Copper (Cu)	max. 2
Iron (Fe)	max. 2

0.995-1.005

PRODUCT NO.	PACKING	CONT. BOX
1182.0025	25 g Glass	6
1182.0100	100 g	6
1182.0250	250 g	6
1182.1000	1 kg	

Silver Nitrate

See for qualities meeting the requirements of cGMP, USP/NF, BP/Ph.Eur, JP or FCC, page 36

Titer (mol/l)



1182

Silver Nitrate

1 ml = 1 mg Chloride / 'BAKER ANALYZED' / Standard Solution

7095

► AgNO₃ **M** = 169.87 g/mol **CAS NO.** 7761-88-8

EINECS 231-853-9 NC CODE 2843 21 00 **R:** 52/53

S: 61

Titer 0.995-1.005mg Cl/ml Trace Impurities (in ppm): Chloride (CI) max. 1 Sulfate (SO₄) max. 1 **PRODUCT PACKING** BOX 11 7095.1000 7095.9020 20 l Polycube

Volumetric Solution, ready for use. Each lot of this product is standardized potentiometrically against NIST tracable reference standard.

Silver Nitrate

1 mol/l / 'BAKER ANALYZED'

7246

► AgNO₃ **M** = 169.87 g/mol

1 l = 1.14 kg

CAS NO. 7761-88-8

EINECS 231-853-9

NC CODE 2843 21 00

UN/ID NO. 1760

ADR/RID 8 C9

IMDG 8/II

R: 34-51/53

S: 26-36/37/39-45-57





PRODUCT	PACKING	CONT.
NO.		вох
7246.1000	11	6

Volumetric Solution, ready for use.

Silver Nitrate

INTERNATION

7096	-	~~	
	- /	W4	ъ
		~~	

0.111101/1 / DAKEITAIVAE	TELD				1070
► AgNO ₃	Titer (mol/l)	0.0997-0.1003	PRODUCT	PACKING	CONT.
M = 169.87 g/mol	Trace Impurities (in ppm):		NO.		вох
1 = 1.01 kg	Chloride (CI)	max. 1	7096.1000	11	6
CAS NO. 7761-88-8	Sulfate (SO ₄)	max. 1	7096.5000	5 l Polycube	
EINECS 231-853-9			7096.9010	10 l Polycube	
NC CODE 2843 21 00			7096.9020	20 l Polycube	
R: 52/53					
S: 61			Volumetric S	Solution, ready for use.	

Silve

Silver Nitrate

0.05 mol/l / 'BAKER ANALYZED' 7245

Titer (mol/l)

► AgNO₃ **M** = 169.87 g/mol **1 l** = 1.01 kg

> **CAS NO.** 7761-88-8 **EINECS** 231-853-9 NC CODE 2843 21 00

> > **R:** 52/53 **S:** 61

PRODUCT PACKING CONT. NO. вох 7245.1000 11

Volumetric Solution, ready for use.

0.0495-0.0505

Silver Nitrate

4682

► AgNO₃

D F G

M = 169.87 g/mol **CAS NO.** 7761-88-8 **EINECS** 231-853-9

NC CODE 2843 21 00

UN/ID NO. 1493 **ADR/RID** 5.1 O2 IMDG 5.1/II

> **R:** 34-50/53 **S:** 26-37/39-45-57



environment

PRODUCT PACKING CONT. NO. BOX 4682 1 amp.

Volumetric Concentrate, for dilution to 1 l.

Silver Nitrate 4681

0.1 mol/I / DILUT-IT

► AgNO₃

M = 169.87 g/mol **CAS NO.** 7761-88-8

EINECS 231-853-9 NC CODE 2843 21 00

UN/ID NO. 1493 **ADR/RID** 5.1 O2 IMDG 5.1/II

R: 34-50/53

S: 26-37/39-45-57





PRODUCT PACKING CONT. вох NO. 1 amp. 6

Volumetric Concentrate, for dilution to 1 l.

Silver Sulfate

Low in Nitrate / 'BAKER ANALYZED' / ACS

► Ag₂SO₄

M = 311.80 g/mol **CAS NO.** 10294-26-5

EINECS 233-653-7

NC CODE 2843 29 00 **R:** 41

S: 22-26-39



Exceeds ACS Specifications	
Assay	min. 99.5%
Insoluble Matter and Silver Chloride	max. 0.02%
Iron (Fe)	max. 0.001%
Nitrate (NO ₃)	max. 0.001%
Substances not Precipitated by HCI	may 0.03%

PRODUCT	PACKING	CONT.
NO.		вох
1183.0025	25 g Glass	6
1183.0100	100 g	

Silver Sulfate

'BAKER ANALYZED'

M = 311.80 g/mol **M =** 311.80 g/mol **CAS NO.** 10294-26-5

EINECS 233-653-7 **NC CODE** 2843 29 00

R: 41

S: 22-26-39



Assay	min. 99.5%
Insoluble Matter and Chloride	max. 0.02%
Iron (Fe)	max. 0.001%
Nitrate (NO ₃)	max. 0.01%
Substances not Precipitated by HCI	max. 0.03%

PACKING	CONT.
	вох
25 g Glass	6
100 g	
250 g	
	25 g Glass 100 g

Silver Sulfate-Sulfuric acid

 $6.6\ g/I\ /\ 'BAKER\ ANALYZED'\ /\ for\ the\ COD\ determiniation\ according\ AFNOR\ NF\ T\ 90-101$

7641

1536

1 I = 1.84 kg NC CODE 3822 00 00 UN/ID NO. 1830 ADR/RID 8 C1 IMDG 8/II R: 35

S: 26-36/37/39-45



PRODUCT	PACKING	CONT.
NO.		вох
7641.2500	2.5	

Volumetric Solution, ready for use.

Silver Sulfate-Sulfuric acid

10g Ag₂SO₄ + 35ml H₂O + 965ml H₂SO₄ 96% / 'BAKER ANALYZED' / For COD determination according DIN 38409 vol. 41

7147

1 I = 1.84 kg NC CODE 3822 00 00 UN/ID NO. 1830 ADR/RID 8 C1 IMDG 8/II

R: 35 **S:** 26-30-36/37/39-45

C C

PRODUCT	PACKING	CONT.
NO.		вох
7147.1000	11	6
7147.2500	2.5	4

Volumetric Solution, ready for use.

Silver Sulfate-Sulfuric acid

60g Ag₂SO₄ in Sulfuric Acid 96% / 'BAKER ANALYZED' / For COD determination according DIN 38409 vol.41

7148

1 I = 1.84 kg
NC CODE 3822 00 00
ADR/RID 8 C1
IMDG 8/II
R: 35

S: 26-30-36/37/39-45



PRODUCT	PACKING	CONT.
NO.		вох
7148.1000	11	
7148.2500	2.5	

Volumetric Solution, ready for use.

Silver Sulfate-Sulfuric acid

80g Ag₂SO₄ in Sulfuric Acid 96% / 'BAKER ANALYZED' / For COD determination according DIN 38409-H41.

7149

1 l = 1.84 kg NC CODE 3822 00 00 UN/ID NO. 1830 ADR/RID 8 C1 IMDG 8/II R: 35

S: 26-30-36/37/39-45

C corrosive

PRODUCT	PACKING	CONT
NO.		ВОХ
7149.1000	11	
7149.2500	2.5	

Volumetric Solution, ready for use.

SLSI chemicals

See for detailed information section Microelectronic materials page 32 or the website www.jtbaker.com/microelectronic materials page 32 or the website was 32 or the webs





Soda

Soda Lime

1783

11 = 0.94 kg CAS NO. 8006-28-8 NC CODE 3822 00 00

UN/ID NO. 1907 ADR/RID 8 C6

> IMDG 8/III **R:** 35

> > **S:** 26-36/37/39-45



Identification passes test Loss on Drying at 105°C 12-19%

CONT. **PRODUCT PACKING** NO. вох 1783.1000 1 kg 1783.9050 50 kg

Sodium 1000 μ g/ml

5780 (Matrix: 2% nitric acid) / 'BAKER INSTRA-ANALYZED' / Plasma Standard

► Na

M = 22.99 g/mol

NC CODE 3822 00 00

R: 36/38 **S:** 26-37

Certificate Provided Reporting Actual Lot Analysis Sodium (Na) 998-1002 μg/ml **PRODUCT PACKING** CONT. NO. вох 5780.0100 100 ml

Prepared from the highest purity raw material available, generally greater than 99.999% spectral purity. The content of the solution is confirmed to be accurate to within \pm 0.2%. Typically 1000 $\mu\text{g/ml.}$ The certificate of analysis provided reports actual lot analysis. The certificate also lists the trace impurities.

Sodium 1000 µg/ml

6941 (Matrix: 2% nitric acid) / 'BAKER INSTRA-ANALYZED' / Atomic Absorption Standard

Sodium (Na)

► Na

M = 22.99 g/mol

NC CODE 3822 00 00

R: 36/38

S: 26-37



998-1002 µg/ml **PRODUCT**

6941.0100 100 ml 6941 0500 500 ml

CONT.

вох

CONT.

BOX

PACKING

Prepared by dissolution of high purity raw materials (min. 99.99% spectral purity). Assays are verified by ICP against standards traceable to NIST. Standard Reference Material numbers (SRM) are printed on each label.

Sodium 1000 μg/ml 6822

'BAKER ANALYZED' / Atomic Absorption Standard

► Na

M = 22.99 g/mol NC CODE 3822 00 00

Sodium (Na) 998-1002 μg/ml

PRODUCT NO.

6822.0100 100 ml

6822 0500 500 ml

PACKING

Sodium nitrate in water

5746

 $\begin{array}{lll} \textbf{Sodium 10000} & \mu \underline{g} / \underline{m} \\ \textbf{(Matrix: } 2\% & \text{nitric acid)} & \text{'BAKER INSTRA-ANALYZED'} & \text{Plasma Standard} \\ \end{array}$

► Na

M = 22.99 g/mol NC CODE 3822 00 00

R: 36/38

S: 26



Certificate Provided Reporting Actual Lot Analysis Sodium (Na) 9980-10020 µg/ml PRODUCT PACKING CONT. NO. BOX 5746.0100 100 ml

Prepared from the highest purity raw material available, generally greater than 99.999% spectral purity. The content of the solution is confirmed to be accurate to within 0.2 %. Typically 10000 μ g/ml.

The certificate of analysis provided reports actual lot analysis. The certificate also lists the trace impurities.



Sodium Acetate Anhydrous

'BAKER ANALYZED' / ACS

0258

CH ₃ COONa

M = 82.03 g/molCAS NO. 127-09-3EINECS 204-823-8NC CODE 2915 22 00

Meets ACS Specifications	
Assay	min. 99.0%
Calcium (Ca)	max. 0.005%
Chloride (CI)	max. 0.002%
Heavy Metals (as Pb)	max. 0.001%
Insoluble Matter	max. 0.01%
Iron (Fe)	max. 0.001%
Loss on Drying at 120°C	max. 1.0%
Magnesium (Mg)	max. 0.002%
pH of 5% Solution at 25°C	7.0-9.2
Phosphate (PO ₄)	max. 0.001%
Sulfate (SO ₄)	max. 0.003%

PRODUCT NO.	PACKING	CONT. BOX
0258.1000	1 kg	6

Sodium Acetate Anhydrous

'BAKER'

► CH₃COONa

M = 82.03 g/molCAS NO. 127-09-3EINECS 204-823-8NC CODE 2915 22 00

Assay (by Perchloric Acid titrn.)	min. 99%
Chloride (CI)	max. 0.2%
Heavy Metals (as Pb)	max. 0.001%
Insoluble Matter	max. 0.01%
Iron (Fe)	max. 0.002%
Loss on Drying at 120°C	max. 1%
nH of 5% Solution at 25°C	7.0-9.3

PRODUCT NO.	PACKING	CONT. BOX
0259.1000	1 kg	
0259.9025	25 kg	

Sodium Acetate, Anhydrous

See for qualities meeting the requirements of cGMP, USP/NF, BP/Ph.Eur, JP or FCC, page 36

at 254 nm



0259

Sodium Acetate Trihydrate

crystal / 'BAKER ULTRAPURE BIOREAGENT'

CH₃COONa.3H₂	0
--	---

M = 136.08 g/molCAS NO. 6131-90-4EINECS 204-823-8NC CODE 2915 22 00

Chloride (CI) max. 0.001% DNAase Activity none detected Insoluble Matter max. 0.005% pH of 5% Solution at 25°C 7.5-9.2 Potassium (K) max. 0.005% Protease Activity none detected RNAase Activity none detected Substances Reducing KMnO₄ passes test max. 0.002% Sulfate (SO₄) Absorbance of a 1 M Solution (1-cm path vs water):

at 280 nm at 350 nm Trace Impurities (in ppm):

Heavy Metals (as Pb) max. 5
Iron (Fe) max. 5
Phosphate (PO₄) max. 5

max. 0.01

max. 0.01

4	nn	١٥	
-4	υv	"	

 PRODUCT
 PACKING
 CONT.

 NO.
 BOX

 4009.1000
 1 kg

 4009.5000
 5 kg

Questions or suggestions, please contact us at jtbaker.nl@emea.tycohealthcare.com



Sodiu

Sodium Acetate Trihydrate

BAKER HPLC ANALYZED' / for use in High Performance Liquid Chromatography

0393 ► CH₃COONa.3H₂O

M = 136.08 g/mol **CAS NO.** 6131-90-4 **EINECS** 204-823-8 NC CODE 2915 22 00

Assay (by non-aqueous titrn.) 99.0-101.0% Calcium (Ca) max. 0.005% Chloride (CI) max. 0.001% Insoluble Matter max. 0.005% Magnesium (Mg) max. 0.002%

pH of 5% Solution at 25°C 7.5-9.2

Potassium (K) (by FES) max. 0.005% Substances Reducing KMnO₄ passes test Sulfate (SO₄) max. 0.002% Trace Impurities (in ppm): Heavy Metals (as Pb) max. 5 Iron (Fe) max. 5 Phosphate (PO₄) max. 5 Ultraviolet Absorbance (1.00-cm path vs water): at 254 nm max. 0.02 at 280 nm max. 0.01 at 350 nm max. 0.01

PRODUCT

0393.1000

PRODUCT

PACKING

1 kg

Sodium Acetate Trihydrate

0256 'BAKER ANALYZED' / ACS

► CH₃COONa.3H₂O

M = 136.08 g/mol **CAS NO.** 6131-90-4 **EINECS** 204-823-8 NC CODE 2915 22 00

Meets ACS Specifications 99.0-101.0% Assay Calcium (Ca) max. 0.005% Chloride (CI) max. 0.001% Insoluble Matter max. 0.005% Magnesium (Mg) max. 0.002% pH of 5% Solution at 25°C 7.5-9.2 Potassium (K) max. 0.005% Substances Reducing KMnO₄ passes test Sulfate (SO₄) max. 0.002% Trace Impurities (in ppm): Heavy Metals (as Pb) max. 5 Iron (Fe) max. 5

NO. BOX 0256.0500 500 g 6 0256.1000 1 kg 6 0256.5000 5 kg 4

PACKING

CONT.

вох

CONT.

0517

Sodium Acetate Trihydrate

► CH₃COONa.3H₂O **M** = 136.08 g/mol

CAS NO. 6131-90-4 **EINECS** 204-823-8 NC CODE 2915 22 00 Assay (by Perchloric Acid titrn.)

min. 98%

0.0495-0.0505

max. 5

PRODUCT PACKING CONT. вох 0517.1000 1 kg 0517.9025 25 kg 0517.9050 50 kg



Sodium Acetate, TrihydrateSee for qualities meeting the requirements of cGMP, USP/NF, BP/Ph.Eur, JP or FCC, page 36

Titer (mol/l)

Phosphate (PO₄)

Sodium Arsenite

7253 0.05 mol/l / 'BAKER ANALYZED'

► NaAsO₂ **M** = 129.91 g/mol

> 11 = 1.01 kg NC CODE 2842 90 90 UN/ID NO. 1686

ADR/RID 6.1 T4 IMDG 6.1/II

> **R:** 45-52/53 **S:** 20-36/39-53

PRODUCT PACKING CONT. вох 7253.1000 11

Sodium Azide

9099 Assay min. 99%

max. 0.1%

max. 0.001%

max. 0.05% max. 0.1%

•	NaN ₃	
	M =	65.01 g/mol
	CAS NO.	26628-22-8
	EINECS	247-852-1
	NC CODE	2850 00 50

EC NO. 11 004 00 7 **UN/ID NO.** 1687 **ADR/RID** 6.1 T5

> IMDG 6.1/II **R:** 28-32-50/53 **S:** 28-45-60-61



environment



Insoluble Matter	
Loss on Drying	

Alkalinity

Heavy Metals (as Pb)

PRODUCT	PACKING	CONT.
NO.		вох
9099.0100	100 g	6
9099.1000	1 kg	

Sodium Benzoate

2021 CONT. ► C₆H₅COONa Assay 99.0-100.5% **PRODUCT PACKING M** = 144.11 g/mol Alkalinity вох passes test **CAS NO.** 532-32-1 Heavy Metals (as Pb) max. 0.001% 2021.1000 1 kg **EINECS** 208-534-8 Identification 2021.9025 25 kg passes test NC CODE 2916 31 00 Organic Volatile Impurities passes test 2021.9050 50 kg Water (H₂O) max. 1.5%

Sodium Benzoate

See for qualities meeting the requirements of cGMP, USP/NF, BP/Ph.Eur, JP or FCC, page 36

Sodium Bicarbonate ee Sodium Hydrogen Carbonate



Sodium Bicarbonate

See for qualities meeting the requirements of cGMP, USP/NF, BP/Ph.Eur, JP or FCC, page 36

See for qualities meeting the requirements of cGMP, USP/NF, BP/Ph.Eur, JP or FCC, page 36



Sodium Bisulfite See Sodium Hydrogen Sulfite

Sodium Bisulfite



Sodium Borate



See Disodium Tetraborate Decahydrate

Sodium Borate, 10-Hydrate See for qualities meeting the requirements of cGMP, USP/NF, BP/Ph.Eur, JP or FCC, page 36

Mercury (Hg)

Selenium (Se)



Sodium Borohydride

'BAKER ANALYZED' / for cold vapour analysis



CONT.

BOX

► NaBH₄	
M =	37.83 g/mol
CAS NO.	16940-66-2
EINECS	241-004-4
NC CODE	2850 00 20
UN/ID NO.	1426
ADR/RID	4.3 W2
IMDG	4.3/I
R:	15-25-34

S:	26-36-43a-7/8				
	B				
	highly flammable	toxic			

min. 97%
max. 0.5%
max. 0.001%
max. 0.001%
max. 0.005%
max. 0.2
max. 0.2
max. 0.2

9161.0100	100 g	6

PACKING

PRODUCT

NO.

max 0.001 max. 0.2



Sodiu

Sodium Borohydride

Assay

► NaBH₄

9098

M = 37.83 g/mol **CAS NO.** 16940-66-2 **EINECS** 241-004-4

NC CODE 2850 00 20 **UN/ID NO.** 1426

ADR/RID 4.3 W2

IMDG 4.3/I

R: 15-25-34 **S:** 26-36-43a-7/8



PRODUCT **PACKING** min. 97% NO. 9098.0025 25 g Glass

100 g

9098.0100

CONT.

вох

Sodium Bromide

0271

M = 102.89 g/mol CAS NO. 7647-15-6 **EINECS** 231-599-9

NC CODE 2827 51 00

min. 99.0% Assay Alkalinity (as Na₂CO₃) max. 0.015% Barium (Ba) max. 0.002% Bromate (BrO₃) max. 0.001% Calcium and Magnesium (as Ca) max. 0.005% Chloride (CI) max. 0.2% Insoluble Matter max. 0.005% pH of 5% Solution at 25°C 5.5-7.5 max. 0.002% Sulfate (SO₄)

Trace Impurities (in ppm): Heavy Metals (as Pb) max. 5 Iron (Fe) max. 5 Nitrogen Compounds (as N) max. 5

PRODUCT PACKING CONT. BOX 0271.0250 250 g 0271.1000 1 kg 0271.9050 50 kg

2022

Sodium Bromide

► NaBr

M = 102.89 g/mol **CAS NO.** 7647-15-6 EINECS 231-599-9

NC CODE 2827 51 00

Assay (dried basis) 98.0-100.5% Acidity or Alkalinity passes test Appearance of solution passes test Bromates (as BrO₃) passes test Chlorides (as CI) max. 0.6% Identification passes test lodides (as I) passes test Loss on Drying at 105°C max. 3.0% Magnesium and Alkaline-earth Metals (as Ca) max. 200 ppm Sulfates (as SO₄) max. 100 ppm Trace Impurities (in ppm): Heavy Metals (as Pb) max. 10 ppm

PRODUCT PACKING CONT. NO. BOX 1 kg 2022.1000

Stored in an airtight container

Sodium Carbonate

7138 0.5 mol/l / 'BAKER ANALYZED'

EINECS 207-838-8 NC CODE 2836 20 00

0.4975-0.5025 Titer (mol/l)

max. 20 ppm

Iron (Fe)

PRODUCT PACKING CONT. вох NO. 7138.1000 11 7138.9010 10 l Polycube

Volumetric Solution, ready for use. Each lot of this product is standardized potentiometrically against NIST tracable reference standard.



CONT.

вох

Sodium Carbonate

0.05 mol/l / 'BAKER ANALYZED'

Titer (mol/l)

7139

EINECS 207-838-8 NC CODE 2836 20 00
> 7139.1000 1 l 7139.9010 10 l Polycube

Volumetric Solution, ready for use.

Each lot of this product is standardized potentiometrically against NIST tracable reference standard.

Sodium Carbonate Anhydrous

ULTREX Ultrapure Reagen

4923

Na ₂ CO ₃	
M =	105.99 g/mol
CAS NO.	497-19-8
EINECS	207-838-8
NC CODE	2836 20 00
EC NO.	11 005 00 2
R:	36
S:	22-26
	Xi

Actual Analysis Lot. No. V04622	
Assay (dried basis)	99.6%
Ammonium Hydroxide Precipitate	0.001%
Arsenic (As)	< 0.1
Iron (Fe)	< 1
Loss on Drying at 285°C	0.005%
Particulate Matter	0.0002%
Metallic Impurities (in ppm):	
Aluminium (AI)	< 1
Barium (Ba)	< 10
Bismuth (Bi)	< 10
Cadmium (Cd)	< 5
Calcium (Ca)	26
Chromium (Cr)	< 1
Cobalt (Co)	< 1
Copper (Cu)	< 1
Lead (Pb)	< 10
Lithium (Li)	< 100
Magnesium (Mg)	< '
Manganese (Mn)	< 1
Mercury (Hg)	< 0.001
Molybdenum (Mo)	1
Nickel (Ni)	< 5
Niobium (Nb)	< 5
Potassium (K)	< 100
Silver (Ag)	< 1
Strontium (Sr)	< 1
Tin (Sn)	< 10
Titanium (Ti)	< 1
Vanadium (V)	< '
Zinc (Zn)	< 10
Zirconium (Zr)	< 1
Non Metallic Impurities (in ppm):	
Halide (as CI)	1
Nitrogen Compounds (as N)	1
Phosphate (PO ₄)	0.5
Silicon (Si)	< 1

PRODUCT	PACKING	CONT.
NO.		вох
4923.0025	25 g Glass	
4923.0100	100 g Glass	

J.T.Baker: over 100 years of experience.

Sulfur Compounds (as SO₄)

See chapter 1 of this catalogue.



Sodiu

Sodium Carbonate Anhydrous

'BAKER ANALYZED' / Primary Stand 1400

► Na₂CO₃ **M** = 105.99 g/mol **CAS NO.** 497-19-8 **EINECS** 207-838-8 NC CODE 2836 20 00 **EC NO.** 11 005 00 2 **R:** 36 **S:** 22-26



Meets ACS Specifications Assay (after 2 hr at 285°C) 99.95-100.05 Calcium (Ca) max. 0.02 Chloride (CI) max. 0.001 Heavy Metals (as Pb) max. 0.0005 Insoluble Matter max. 0.01 Iron (Fe) max. 0.0005 Loss on Heating at 285°C max. 1.0⁰ Magnesium (Mg) max. 0.004 Phosphate (as PO₄) max. 0.001 max. 0.005 Potassium (K) Silica (SiO₂) max. 0.005 Solubility (8 in 50) passes te Sulfur Compounds (as SO₄) max. 0.003 Product Information (not specifications): Appearance white granules

PRODUCT	PACKING	CONT
NO.		ВОХ
1400.0250	250 g	

Sodium Carbonate Anhydrous 0274

'BAKER ANALYZED' / ACS

► Na₂CO₃ **M** = 105.99 g/mol **CAS NO.** 497-19-8 **EINECS** 207-838-8 NC CODE 2836 20 00 EC NO. 11 005 00 2 **R:** 36

22-26		
Xi		
irritant		

Meets ACS Specifications. Meets Reagent		
Specifications for testing USP/NF monographs		
Assay (dried basis) (by acidimetry)	min. 99.59	
Calcium (Ca)	max. 0.039	
Chloride (CI)	max. 0.0019	
nsoluble Matter	max. 0.019	
oss on Heating at 285°C	max. 1.09	
Magnesium (Mg)	max. 0.0059	
Phosphate (PO ₄)	max. 0.0019	
Potassium (K)	max. 0.0059	
Silica (SiO ₂)	max. 0.0059	
2 1/ 00)	0.000	

Magnesium (Mg)	max. 0.005%
Phosphate (PO ₄)	max. 0.001%
Potassium (K)	max. 0.005%
Silica (SiO ₂)	max. 0.005%
Sulfur Compounds (as SO ₄)	max. 0.003%
Trace Impurities (in ppm):	
Heavy Metals (as Pb)	max. 5
Iron (Fe)	max. 5

PRODUCT	PACKING	CONT.
NO.		вох
0274.0500	500 g	6
0274.1000	1 kg	6
0274.5000	5 kg	4

Sodium Carbonate Anhydrous

► Na₂CO₃

2024

M = 105.99 g/mol **CAS NO.** 497-19-8 **EINECS** 207-838-8 NC CODE 2836 20 00 **EC NO.** 11 005 00 2 R: S:

36	
22-26	
▲ Xi	
X	
irritant	

Assay	99.5-100.5%
Alkali Hydroxides and Bicarbonates	passes test
Appearance of solution	passes test
Arsenic (As)	max. 5 ppm
Chloride (CI)	max.125 ppm
Heavy Metals (as Pb)	max. 50 ppm
Identification	passes test
Iron (Fe)	max. 50 ppm
Loss on Drying	max. 1.0%
Sulfate (SO ₄)	max. 250 ppm

PRODUCT	PACKING	CONT.
NO.		вох
2024.1000	1 kg	6
2024.5000	5 kg	4

Stored in an airtight container.

Sodium Carbonate Anhydrous

1994

► Na₂CO₃ **M** = 105.99 g/mol **CAS NO.** 497-19-8 **EINECS** 207-838-8 NC CODE 2836 20 00

> **EC NO.** 11 005 00 2 **R:** 36 **S:** 22-26



min. 99.5% Assay

PRODUCT PACKING CONT. NO. BOX 1994.9050 50 kg



Sodium Carbonate, Anhydrous

See for qualities meeting the requirements of cGMP, USP/NF, BP/Ph.Eur, JP or FCC, page 36



Sodium Carbonate, Monohydrate

See for qualities meeting the requirements of cGMP, USP/NF, BP/Ph.Eur, JP or FCC, page 36



Sodium Chloride

ULTREX Ultrapure Reagent

4924	

► NaCl	
M =	58.44 g/mol
CAS NO.	7647-14-5
EINECS	231-598-3
NC CODE	2501 00 99

Certificate Provided Reporting Actual Lot Analysis Actual Lot Analysis Lot No. B21582		
Loss on drying at 650°C	0.003%	
Particulate Matter	0.003%	
Metallic Impurities in parts per m	illion (μg/g):	
Aluminium (AI)	0.7	
Barium (Ba)	< 0.2	
Bismuth (Bi)	< 2.2	
Cadmium (Cd)	< 0.1	
Calcium (Ca)	1.2	
Chromium (Cr)	< 0.2	
Cobalt (Co)	< 0.2	
Copper (Cu)	0.5	
Iron (Fe)	0.2	
Lead (Pb)	< 0.2	
Lithium (Li)	< 0.2	
Magnesium (Mg)	0.8	
Manganese (Mn)	< 0.1	
Mercury (Hg)	0.4	
Molybdenum (Mo)	< 0.1	
Nickel (Ni)	< 0.2	
Potassium (K)	2.2	
Silver (Ag)	< 0.3	
Strontium (Sr)	0.2	
Tin (Sn)	< 0.2	
Titanium (Ti)	0.2	
Vanadium (V)	< 0.1	
Zinc (Zn)	< 0.2	
Zirconium (Zr)	< 0.2	
Non-Metallic Impurities in parts p	per million (μg/g):	
Arsenic (As)	< 0.8	
Bromide (Br)	0.0008	
lodide (I)	< 0.005	

Nitrogen Compounds (as N) Phosphate (PO₄)

Sulfur Compounds (as SO₄)

PRODUCT	PACKING	CONT.
NO.		вох
4924.0100	100 g Glass	

For Laboratory, Research or Manufacturing Use.

Mallinckrodt Baker's cGMP Manufactured Chemicals for the Biopharmaceutical industry are a necessity for uncomplicated scale-up.

0.8

< 0.2

See chapter 6 of this catalogue.



Sodiu

Sodium Chloride

4058 stal / 'BAKER ULTRAPURE BIOREAGENT'

► NaCl

M = 58.44 g/mol **CAS NO.** 7647-14-5 **EINECS** 231-598-3 NC CODE 2501 00 99

For Liquid Chromatography and Molecular Biology applications

min. 99.0% Assay **DNAase Activity** none detected pH of 5% Solution at 25°C 5.5-7.5 Protease Activity none detected RNAase Activity none detected

Maximum Limits of Impurities:

Iron (Fe)

Phosphate (PO₄)

Nitrogen Compounds (as N)

Barium (Ba) 0.001% Bromide (Br) 0.01% Calcium, Magnesium and R₂O₃ Precipitate 0.005% Chlorate and Nitrate (as NO₃) 0.003% Insoluble Matter 0.005% lodide (I) 0.002% Potassium (K) 0.005% Sulfate (SO₄) 0.003% Trace Impurities (in ppm): Heavy Metals (as Pb) max. 2

max. 1

max. 5

max. 5

max. 5

PRODUCT **PACKING** CONT. NO. вох 4058.0500 500 g 4058.2500 2.5 kg 4058.9012 12 kg

0277

Sodium Chloride

crystal / 'BAKER ANALYZED' / ACS

► NaCl

M = 58.44 g/mol **CAS NO.** 7647-14-5 **EINECS** 231-598-3 NC CODE 2501 00 99

Meets ACS Specifications Assay min. 99.0% Barium (Ba) max. 0.002% Bromide (Br) max. 0.01% Calcium (Ca) max. 0.002% Chlorate and Nitrate (as NO₃) max. 0.003% max. 0.005% Insoluble Matter lodide (I) max. 0.002% Magnesium (Mg) max. 0.001% pH of 5% Solution at 25°C 5.0-9.0

max 0.005% Potassium (K) Sulfate (SO₄) max. 0.004% Trace Impurities (in ppm): max. 5

Heavy Metals (as Pb)

Iron (Fe) max. 2 Phosphate (PO₄) max. 5

PRODUCT NO.	PACKING	CONT. BOX
0277.1000	1 kg	6
0277.5000	5 kg	4
0277.9010	10 kg	
0277.9050	50 kg	

0278

Sodium Chloride

crystal / 'BAKER ANALYZED'

► NaCl

M = 58.44 g/mol **CAS NO.** 7647-14-5 FINECS 231-598-3

NC CODE 2501 00 99

Assay min. 99.5% Barium (Ba) max. 0.002% Bromide (Br) max. 0.01% Calcium (Ca) max. 0.002% Insoluble Matter max. 0.005% lodide (I) max. 0.005% Magnesium (Mg) max. 0.001% 5.0-8.0 pH of 5% Solution at 25°C Potassium (K) max. 0.01% Sulfate (SO₄) max. 0.01% Trace Impurities (in ppm):

Arsenic (As) max. 3 Heavy Metals (as Pb) max. 5 Iron (Fe) max. 3

PRODUCT NO.	PACKING	CONT. BOX
0278.1000	1 kg	6
0278.5000	5 kg	4
0278.9050	50 kg	



Sodium Chloride

e for qualities meeting the requirements of cGMP, USP/NF, BP/Ph.Eur, JP or FCC, page 36

Phosphate (PO₄)

Sodium Chloride

0.1 mol/L / 'BAKER ANALYZED'

EINECS 231-598-3 0.099-0.101 PRODUCT **PACKING** CONT. Titer (mol/l) NC CODE 2501 00 91 NO. вох 7140.1000 11 7140.9010 10 l Polycube Volumetric Solution, ready for use.

Each lot of this product is standardized potentiometrically against NIST tracable reference standard.

Sodium Chloride

4684

7140

N	aCl
---	-----

M = 58.44 g/mol **CAS NO.** 7647-14-5 **EINECS** 231-598-3 NC CODE 2501 00 91

PRODUCT	PACKING	CONT.
NO.		вох
4684	1 amp.	6

Volumetric Concentrate, for dilution to 1 l.

Sodium Citrate Dihydrate crystal / 'BAKER ANALYZED'

0280

M = 294.10 g/mol **CAS NO.** 6132-04-3 **EINECS** 200-675-3 NC CODE 2918 15 00

Trace Impurities (in ppm):	
Sulfate (SO ₄)	max. 0.002%
pH of 5% Solution at 25°C	7.0-9.0
Insoluble Matter	max. 0.005%
Chloride (CI)	max. 0.001%
Calcium (Ca)	max. 0.005%
Assay	min. 99.0%

max. 5

max. 5

PRODUCT	PACKING	CONT.
NO.		вох
0280.0500	500 g	6
0280.1000	1 kg	6
0280.5000	5 kg	4

Sodium Citrate Dihydrate crystal / 'BAKER'

1758

► HOC(COONa)(CH₂COONa)₂.2H₂O **M** = 294.10 g/mol

CAS NO. 6132-04-3 **EINECS** 200-675-3 NC CODE 2918 15 00

Assay	99.0-100.5%
Acid or alkaline reacting substances	passes test
Appearance of solution	passes test
Chloride (CI)	max. 0.003%
Heavy Metals (as Pb)	max. 10 ppm
Identification	passes test
Oxalate (C ₂ O ₄)	max. 300 ppm
Substances Darkened by H ₂ SO ₄	passes test
Sulfate (SO ₄)	max. 0.005%
Tartrate	passes test
Water (H ₂ O)	11.0-13.0%

PRODUCI	PACKING	CONI.
NO.		вох
1758.1000	1 kg	6
1758.9050	50 kg	

Sodium Citrate Dihydrate

See for qualities meeting the requirements of cGMP, USP/NF, BP/Ph.Eur, JP or FCC, page 36

Heavy Metals (as Pb)

Iron (Fe)



Sodium Cyanide 'BAKER ANALYZED' / ACS

0281

► NaCN	
--------	--

M = 49.01 g/mol

CAS NO. 143-33-9 EINECS 205-599-4 NC CODE 2837 11 00

EC NO. 6 007 00 5 **UN/ID NO.** 1689

ADR/RID 6.1 T5

IMDG 6.1/I **R:** 26/27/28-32-50/53

S: 28-29-45-60-61-7



environment



Meets ACS Specifications	
Assay (argentometric titrn.)	min. 95.0%
Chloride (CI)	max. 0.15%
Iron, total (as Fe)	max. 0.005%
Phosphate (PO ₄)	max. 0.02%
Sulfate (SO ₄)	max. 0.05%
Sulfide (S)	max. 0.005%
Thiocyanate (SCN)	max. 0.02%
Trace Impurities (in ppm):	
Lead (Pb)	max. 5

PRODUCT	PACKING	CONT.
NO.		вох
0281.0100	100 g	
0281.9050	50 kg	



Sodium Dichromate Dihydrate

0282 'BAKER ANALYZED' / ACS

► Na₂Cr₂O₇.2H₂O

M = 298.00 g/mol

CAS NO. 7789-12-0

EINECS 234-190-3

NC CODE 2841 30 00

EC NO. 24 004 00 7

UN/ID NO. 3288

ADR/RID 6.1 T5

IMDG 6.1/III

R: 21-25-26-34-42/43-45-46-48/23-

Specifications for testing USP/NF monographs 99.5-100.5% Assay Aluminium (AI) max. 0.002% Calcium (Ca) max. 0.002% Chloride (CI) max. 0.005% Insoluble Matter max. 0.005% Magnesium (Mg) max. 0.005% pH of 5% Solution at 25°C 3.4-4.8 Potassium (K) max. 0.01% max. 0.01% Sulfate (SO₄)

Meets ACS Specifications. Meets Reagent

PRODUCT	PACKING	CONT.
NO.		вох
0282.0500	500 g	
0282.9050	50 kg	

S: 45-53-60-61

50/53-60-61-8





environment

Sodium Diphosphate Decahydrate

0309 'BAKER ANALYZED' / ACS

Na₄P₂O₇.10H₂O **M** = 446.06 g/mol **CAS NO.** 13472-36-1 **EINECS** 231-767-1 **NC CODE** 2835 39 00

Meets ACS Specifications	
Assay	99.0-103.0%
Chloride (CI)	max. 0.002%
Heavy Metals (as Pb)	max. 0.001%
Insoluble Matter	max. 0.01%
Iron (Fe)	max. 0.001%
Nitrogen Compounds (as N)	max. 0.001%
pH of 5% Solution at 25°C	9.5-10.5
Sulfate (SO ₄)	max. 0.005%

PRODUCT	PACKING	CONT.
NO.		вох
0309.1000	1 kg	6
0309.9050	50 kg	

1766

Sodium Diphosphate Decahydrate

► Na₄P₂O₇.10H₂O **M** = 446.06 g/mol **CAS NO.** 13472-36-1 **EINECS** 231-767-1

NC CODE 2835 39 00

Assay min. 99% Heavy Metals (as Pb) max. 0.001%
 PRODUCT NO.
 PACKING CONT.

 1766.1000
 1 kg

1768

Sodium Dihydrogen Phosphate Dihydrate

1768 'BAK

► NaH₂PO₄.2H₂O **M** = 156.01 g/mol **CAS NO.** 13472-35-0 **EINECS** 231-449-2 **NC CODE** 2835 22 00 Assay min. 98%

Trace Impurities (in ppm):

Arsenic (As) max. 2

 PRODUCT NO.
 PACKING CONT. BOX

 1768.1000
 1 kg

 1768.9050
 50 kg

4011

Sodium Dihydrogen Phosphate Monohydrate

ystal / 'BAKER ULTRAPURE BIOREAGENT'

► NaH₂PO₄.H₂O **M** = 137.99 g/mol **CAS NO.** 10049-21-5 **EINECS** 231-449-2 **NC CODE** 2835 22 00

For Liquid Chromatography and Molecular Biology 99.0-102.0% Assay DNase Activity none detected Heavy Metals (as Pb) max. 0.001% Insoluble Matter, Calcium and NH₄OH max. 0.01% Iron (Fe) max. 0.001% pH of 5% Solution at 25°C 4.1-4.5 Protease Activity none detected RNase Activity none detected

B0)

J.T.Baker

Sodium Dihydrogen Phosphate Monohydrate 'BAKER ANALYZED' / ACS

► NaH₂PO₄.H₂O

M = 137.99 g/mol **CAS NO.** 10049-21-5 **EINECS** 231-449-2 NC CODE 2835 22 00

Meets ACS Specifications	
Assay	98.0-102.0%
Calcium (Ca)	max. 0.005%
Heavy Metals (as Pb)	max. 0.001%
Insoluble Matter	max. 0.01%
Iron (Fe)	max. 0.001%
pH of 5% Solution at 25°C	4.1-4.5
Potassium (K)	max. 0.01%
Sulfate (SO ₄)	max. 0.003%
Trace Impurities (in ppm):	
Chloride (CI)	max. 5

PRODUCT	PACKING	CONT.
NO.		вох
0303.0500	500 g	6
0303.1000	1 kg	6
0303.9025	25 kg	

Sodium Dihydrogen Phosphate, Monohydrate

See for qualities meeting the requirements of cGMP, USP/NF, BP/Ph.Eur, JP or FCC, page 36



0267

1767

0303

Sodium Disulfite

'BAKER ANALYZED' / ACS

► Na₂S₂O₅ **M** = 190.10 g/mol

> **CAS NO.** 7681-57-4 **EINECS** 231-673-0 NC CODE 2832 10 00 **EC NO.** 16 063 00 2

> > **R:** 22-31-41 **S:** 26-39-46



Exceeds ACS Specifications	
Assay	min. 97.0%
Assay (SO ₂)	min. 58.5%
Chloride (CI)	max. 0.005%
Heavy Metals (as Pb)	max. 0.001%
Insoluble Matter	max. 0.005%
Thiosulfate (S ₂ O ₃)	max. 0.05%
Trace Impurities (in ppm):	
Iron (Fe)	max. 5

Iron (Fe)

PRODUCT	PACKING	CONT.
NO.		вох
0267.0500	500 g	6
0267.1000	1 kg	6
0267.9050	50 kg	

Contains Sodium Disulfite and Sodium meta-Bisulfite.

Sodium Disulfite

M = 190.10 g/mol **CAS NO.** 7681-57-4

EINECS 231-673-0 NC CODE 2832 10 00

EC NO. 16 063 00 2 **R:** 22-31-41

S: 26-39-46



Assay min. 97%

1767.1000 1 kg

PRODUCT PACKING CONT. BOX

Sodium Dithionite

2071

► Na₂S₂O₄

► Na₂S₂O₅

M = 174.14 g/mol

CAS NO. 7775-14-6

EINECS 231-890-0

NC CODE 2831 10 00 **EC NO.** 16 028 00 1

UN/ID NO. 1384

ADR/RID 4.2 S4

IMDG 4.2/II **R:** 22-31-7

S: 26-28-43a-7/8



Assay	min. 85%
Chloride (CI)	max. 0.8%
Iron (Fe)	max. 0.005%

PRODUCT	PACKING	CONT.
NO.		вох
2071.1000	1 kg	



Sodium Dodecyl Sulfate

BAKER ULTRAPURE BIOREAGENT / For Electrophoresis and Protein Solubilization

► CH₃(CH₂)₁₁OSO₃Na **M** = 288.38 g/mol **CAS NO.** 151-21-3 **EINECS** 205-788-1 NC CODE 3402 19 00 **UN/ID NO.** 1325

ADR/RID 4.1 F1 IMDG 4.1/III

R: 11-20/21-36/37/38

S: 26-36/37





Assay (as C₁₂H₂₅OSO₃Na) min. 99% Appearance passes test Chloride (CI) max. 0.03% DNase Activity none detected pH of 1% Solution at 25°C 5.0-8.0 Protease Activity none detected RNase Activity Absorbance of a 3% (w/v) Aqueous Solution (1-cm

path vs water): at 280 nm Trace Impurities (in ppm): Heavy Metals (as Pb) max. 2 Iron (Fe) max. 1 Nitrogen (N) Phosphate (PO₄) max. 1

NO.	0.01
	B0)
1095.0025 2	5 g Glass
1095.0100 1	00 g
1095.1000 1	kg
1095.2000 2	kg
1095.9012 1	2 kg

Sodium Dodecyl Sulfate

2811

► CH₃(CH₂)₁₁OSO₃Na M = 288.38 g/mol

CAS NO. 151-21-3 **EINECS** 205-788-1


NC CODE 3402 19 00 UN/ID NO. 1325 ADR/RID 4.1 F1 IMDG 4.1/III

R: 11-20/21-36/37/38

S: 26-36/37







Sodium Ethylenediaminetetraacetate

Assay

Sodium Fluoride

0285 'BAKER ANALYZED' / ACS

► NaF M = 41.99 g/mol**CAS NO.** 7681-49-4 **EINECS** 231-667-8 NC CODE 2826 11 00 **EC NO.** 9 004 00 7 **UN/ID NO.** 1690 **ADR/RID** 6.1 T5 IMDG 6 1/III **R:** 25-32-36/38 **S:** 22-36-45

Meets ACS Specifications. Meets Reagent Specifications for testing USP/NF monographs min. 99.0% max. 0.005% Heavy Metals (as Pb) max 0.002% Insoluble Matter max. 0.02% Iron (Fe) max. 0.003% Loss on Drying at 150°C max. 0.2% max. 0.02% Potassium (K) Sodium Fluosilicate (Na₂SiF₆) max. 0.1% Sulfate (SO₄) max. 0.02% Sulfite (SO₂) max. 0.005% Titrable Acid (meq/g) max. 0.03 Titrable Base (meq/g) max. 0.01

BOX 0285.0500 500 q

CONT.

PACKING

PRODUCT

Innovation is principal to our business.

CONT.

вох

Sodium Fluoride

2016

98.5-100.5%

passes test

passes test

passes test

passes test

max. 0.5%

max. 200 ppm

max. 200 ppm

PRODUCT

2016.1000

PRODUCT

0287.0500

0287.9025

0287.9050

PRODUCT

0394.2500

NO.

NO.

PACKING

1 kg

► NaF

M = 41.99 g/mol **CAS NO.** 7681-49-4

EINECS 231-667-8

NC CODE 2826 11 00

EC NO. 9 004 00 7

UN/ID NO. 1690

ADR/RID 6.1 T5

IMDG 6.1/III

R: 25-32-36/38

S: 22-36-45



Sodium Fluoride See for qualities meeting the requirements of cGMP, USP/NF, BP/Ph.Eur, JP or FCC, page 36

Sulfate (SO₄)

Assay

Acidity or Alkalinity

Chlorides (as CI)

Fluosilicate

Identification

Loss on Drying

Sulfates (as SO₄)

Appearance of solution



Sodium Formate

'BAKER ANALYZED' / ACS

0287

► HCOONa

CAS NO. 141-53-7

EINECS 205-488-0

M = 68.01 g/mol

NC CODE 2915 12 00

R: 36 **S**: 26

Exceeds ACS Specifications Assay

Calcium (Ca) max. 0.001% Chloride (CI) max. 0.001% Insoluble Matter max. 0.005% pH of 5% Solution at 25°C 6.5-8.5

Trace Impurities (in ppm):

Heavy Metals (as Pb) max. 5 Iron (Fe) max. 5



вох

CONT.

Sodium Hexametaphosphate

► (NaPO₃)₆

M = 611.77 g/mol **CAS NO.** 10124-56-8 **EINECS** 233-343-1 NC CODE 2835 39 00

P₂O₅ Content

68.0-71.0%

99.7-100.3%

min. 99.0%

max. 0.001%

PRODUCT PACKING NO. 0908.3000 3 kg

PACKING

2.5 kg

PACKING

500 g

25 kg

50 kg

0908 CONT.

BOX

BOX

Sodium Hydrogen Carbonate

der / 'BAKER HPLC ANALYZED' / for use in High Performance Liquid Chromatography

0394

► NaHCO₃

M = 84.01 g/mol

CAS NO. 144-55-8 **EINECS** 205-633-8

NC CODE 2836 30 00

Calcium, Magnesium and $\ensuremath{R_2O_3}$ Precipitate max. 0.02% max. 0.003% Chloride (CI) Insoluble Matter max. 0.015% Iron (Fe) max. 0.001% Phosphate (PO₄) max 0.001% Potassium (K) max. 0.005% Sulfur Compounds (as SO₄) max. 0.003%

Trace Impurities (in ppm):

Ammonium (NH₄) max. 5 Heavy Metals (as Pb) max. 5

Ultraviolet Absorbance (1.00-cm path vs water):

at 254 nm max. 0.05 at 280 nm max. 0.02 at 350 nm max. 0.01 CONT.

365



Sodium Hydrogen Carbonate

0263 der / 'BAKER ANALYZED' / AC

▶ NaHCO₃ **M** = 84.01 g/mol **CAS NO.** 144-55-8 **EINECS** 205-633-8 NC CODE 2836 30 00

Meets ACS Specifications Assay (dried basis) 99.7-100.3% Calcium (Ca) max. 0.02% Chloride (CI) max. 0.003% Insoluble Matter max. 0.015% Iron (Fe) max. 0.001% Magnesium (Mg) max. 0.005% Phosphate (PO₄) max. 0.001% Potassium (K) max. 0.005% Sulfur Compounds (as SO₄) max. 0.003% Trace Impurities (in ppm): Ammonium (NH₄) max. 5 Heavy Metals (as Pb) max. 5

PRODUCT	PACKING	CONT.
NO.		вох
0263.0500	500 g	6
0263.1000	1 kg	6
0263.5000	5 kg	4

Sodium Hydrogen Carbonate

1780

► NaHCO₃

M = 84.01 g/mol**CAS NO.** 144-55-8 EINECS 205-633-8 NC CODE 2836 30 00

99.0-101.0% Assay Ammonium (NH₄) max. 20 ppm Appearance of solution passes test Arsenic (As) max. 2 ppm Calcium (Ca) max. 100 ppm Carbonate (CO₃) passes test Chloride (CI) max. 150 ppm Heavy Metals (as Pb) max. 10 ppm Identification passes test Iron (Fe) max. 20 ppm Sulfate (SO₄) max. 150 ppm

PRODUCT	PACKING	CONT.
NO.		вох
1780.1000	1 kg	6
1780.5000	5 kg	
1780.9050	50 kg	

4062

Sodium Hydrogen Phosphate Anhydrous

► Na₂HPO₄

M = 141.96 g/mol **CAS NO.** 7558-79-4 **EINECS** 231-448-7 NC CODE 2835 22 00

For Liquid Chromatography and Molecular Biology applications Assay min. 99.0% **DNAase Activity** none detected Heavy Metals (as Pb) max. 0.001% Insoluble Matter max.0.01% Iron (Fe) max. 0.001% Loss on Drying at 105°C max. 0.2% pH of 5% Solution at 25°C 8.7-9.3 Protease Activity none detected RNAase Activity none detected

NO.		
		ВОХ
4062.0500	500 g	
4062.2500	2.5 kg	

0306

Sodium Hydrogen Phosphate Anhydrous

'BAKER ANALYZED' / ACS

► Na₂HPO₄

M = 141.96 g/mol **CAS NO.** 7558-79-4 **EINECS** 231-448-7 NC CODE 2835 22 00

Meets ACS Specifications Assav min. 99.0% Chloride (CI) max. 0.002% Heavy Metals (as Pb) max. 0.001% Insoluble Matter max. 0.01% Iron (Fe) max. 0.002% Loss on Drying at 105°C max. 0.2% pH of 5% Solution at 25°C 8.7-9.3 max. 0.005%

PRODUCI	PACKING	CONI.
NO.		вох
0306.0500	500 g	6
0306.1000	1 kg	6
0306.5000	5 kg	



Sodium Hydrogen Phosphate, Anhydrous

for qualities meeting the requirements of cGMP, USP/NF, BP/Ph.Eur, JP or FCC, page 36

Sulfate (SO₄)



Sodium Hydrogen Phosphate Dihydrate 'BAKER ANALYZED'

0326

► Na₂HPO₄.2H₂O **M** = 177.99 g/mol **CAS NO.** 10028-24-7 **EINECS** 231-448-7

NC CODE 2835 22 00

Assay	min. 99%
Chloride (CI)	max. 0.001%
Heavy Metals (as Pb)	max. 0.001%
Insoluble Matter	max. 0.005%
Iron (Fe)	max. 0.001%
pH of 0.1 M Solution at 25°C	9.0-9.5
Sulfate (SQ.)	max 0.005%

PRODUCT	PACKING	CONT.
NO.		вох
0326.1000	1 kg	6
0326.9050	50 kg	

Sodium Hydrogen Phosphate Dodecahydrate

0304

Na₂HPO₄.12H₂O **M** = 358.14 g/mol **CAS NO.** 10039-32-4 **EINECS** 231-448-7

NC CODE 2835 22 00

Chloride (CI)	max. 0.001%
Heavy Metals (as Pb)	max. 0.001%
Insoluble Matter	max. 0.005%
Iron (Fe)	max. 0.001%
pH of 5% Solution at 25°C	8.7-9.3
Sulfate (SO ₄)	max. 0.02%
Trace Impurities (in ppm):	
Nitrogen Compounds (as N)	max. 5

PRODUCT	PACKING	CONT.
NO.		вох
0304.1000	1 kg	6
0304.9050	50 kg	

Sodium Hydrogen Phosphate Heptahydrate

crystal / 'BAKER ANALYZED' / ACS

0305

► Na₂HPO₄.7H₂O

M = 268.07 g/mol **CAS NO.** 7782-85-6 **EINECS** 231-448-7 NC CODE 2835 22 00

Meets ACS Specifications. Meets Reagent		
Specifications for testing USP/NF monographs		
Assay	98.0-102.0%	
Chloride (CI)	max. 0.001%	
Heavy Metals (as Pb)	max. 0.001%	
Insoluble Matter	max. 0.005%	
Iron (Fe)	max. 0.001%	
Loss on Drying at 105°C	43-50%	
pH of 5% Solution at 25°C	8.7-9.3	
Sulfate (SO ₄)	max. 0.002%	

PRODUCT	PACKING	CONT
NO.		ВОХ
0305.0500	500 g	
0305.7100	100 lbs	

Sodium Hydrogen Phosphate, 7-HydrateSee for qualities meeting the requirements of cGMP, USP/NF, BP/Ph.Eur, JP or FCC, page 36



0264

Sodium Hydrogen Sulfate Monohydrate

'BAKER ANALYZED'

► NaHSO₄.H₂O

M = 138.07 g/mol **CAS NO.** 10034-88-5 **EINECS** 240-778-0

NC CODE 2833 19 00 **EC NO.** 16 046 00 0

UN/ID NO. 3260 ADR/RID 8 C2 IMDG 8/III

> **R:** 41 **S:** 24-26



Assay (as H ₂ SO ₄)	35.0-36.5%	
Calcium and Magnesium Precipitate	max. 0.005%	
Chloride (CI)	max. 0.001%	
Insoluble Matter and NH ₄ OH Precipitate	max. 0.005%	
Iron (Fe)	max. 0.001%	
Phosphate (PO ₄)	max. 0.001%	
Trace Impurities (in ppm):		
A:- (A -)	1	

Trace Impurities (in ppm):	
Arsenic (As)	max. 1
Heavy Metals (as Pb)	max. 5

PRODUCT	PACKING	CONT.
NO.		вох
0264.1000	1 kg	
02/4 0050	EO ka	

Sodium Hydrogen Sulfite

'BAKER ANALYZED' / ACS

0266

CAS NO.	7631-90-5
EINECS	231-548-0
NC CODE	2832 10 00
EC NO.	16 064 00 8
R:	22-31
S:	25-46
	Xn

Meets ACS Specifications	
Assay (as SO ₂) (by iodometry)	min. 58.5%
Chloride (CI)	max. 0.02 %
Heavy Metals (as Pb)	max. 0.001%
Insoluble Matter	max. 0.005%
Iron (Fe)	max. 0.002%
Thiosulfate (S ₂ O ₃)	act. value reported

PRODUCT	PACKING	CONT.
NO.		вох
0266.0500	500 g	6
0266.1000	1 kg	
0266.9050	50 kg	

Sodium Hydrosulfite

See Sodium Dithionit

0403

Sodium Hydroxide

Pellets / 'BAKER ANALYZED

M = 40.00 g/mol
CAS NO. 1310-73-2
EINECS 215-185-5
NC CODE 2815 11 00
EC NO. 11 002 00 6

NC CODE 2815 11 00 EC NO. 11 002 00 UN/ID NO. 1823 ADR/RID 8 C6 IMDG 8/II R: 35

RED 8 C6

DG 8/II

R: 35

S: 26-37/39-45

Corrosive

min. 99.0% Ammonium Hydroxide Precipitate max. 0.02% Chloride (CI) max. 5 ppm Heavy Metals (as Ag) max. 5 ppm Heavy Metals (as Pb) max. 5 ppm Na₂CO₃ max. 1.0% Potassium (K) max. 0.01% Trace Impurities (in ppm): Iron (Fe) max. 3 Mercury (Hg) (by AAS) max. 0.1 Nickel (Ni) max. 5 Nitrogen Compounds (as N) max. 3

max. 10

max. 5

Phosphate (as PO₄)

Sulfate (SO₄)

 PRODUCT NO.
 PACKING
 CONT.

 0403.5000
 5 kg
 4

 0403.9050
 50 kg
 4

0402

Sodium Hydroxide

Pellets / 'BAKER ANALYZED

► NaOH

M = 40.00 g/mol CAS NO. 1310-73-2 EINECS 215-185-5 NC CODE 2815 11 00 EC NO. 11 002 00 6 UN/ID NO. 1823

ADR/RID 8 C6 **IMDG** 8/II **R:** 35

S: 26-37/39-45 C corrosive

Assay min. 98% Carbonate (as Na₂CO₃) max. 1% Potassium (K) max. 0.05% Silica (SiO₂) max. 0.001% Trace Impurities (in ppm): Aluminium (AI) max. 5 Calcium (Ca) Chloride (CI) max. 5 max. 5 Heavy Metals (as Pb) Iron (Fe) max. 5 Nickel (Ni) max. 5 Phosphate (PO₄) max. 5 Sulfate (SO₄) max. 5 Total Nitrogen (N) max. 3

 PRODUCT NO.
 PACKING
 CONT.

 0402.0500
 500 g
 6

 0402.1000
 1 kg
 6

 0402.5000
 5 kg
 4

 0402.9050
 50 kg

Sodium Hydroxide

Pellets / 'BAKER ANALYZED' / ACS

► NaOH

M = 40.00 g/molCAS NO. 1310-73-2EINECS 215-185-5NC CODE 2815 11 00

EC NO. 11 002 00 6 UN/ID NO. 1823

ADR/RID 8 C6 IMDG 8/II R: 35

S: 26-37/39-45



Specifications for testing USP/NF monographs		
Assay	min. 97.0%	
Calcium (Ca)	max. 0.005%	
Chloride (CI)	max. 0.005%	
Heavy Metals (as Ag)	max. 0.001%	
Magnesium (Mg)	max. 0.002%	
Potassium (K)	max. 0.01%	
Sodium Carbonate (Na ₂ CO ₃)	max. 0.5%	
Trace Impurities (in ppm):		

Exceeds ACS Specifications. Meets Reagents

Trace Impurities (in ppm):		
Iron (Fe)	max. 3	
Mercury (Hg)	max. 0.1	
Nickel (Ni)	max. 5	
Nitrogen Compounds (as N)	max. 3	
Phosphate (PO ₄)	max. 2	
Sulfate (SO ₄)	max. 5	

PRODUCT	PACKING	CONT.
NO.		вох
0003.0500	E00 -	

Low in carbonate.

www.jtbaker.com/europe



Sodium Hydroxide

Pellets / 'BAKER ANALYZED' / ACS

0288

	reliets / DAIL
► NaOH	
M =	40.00 g/mol
CAS NO.	1310-73-2
EINECS	215-185-5
NC CODE	2815 11 00
EC NO.	11 002 00 6
UN/ID NO.	1823
ADR/RID	8 C6
IMDG	8/11
R:	35
S:	26-37/39-45
	C C

Meets ACS Specifications. Meets Re	eagent
Specifications for testing USP/NF me	onographs
Assay	min. 97.0%
Calcium (Ca)	max. 0.005%
Chloride (CI)	max. 0.005%
Heavy Metals (as Ag)	max. 0.001%
Magnesium (Mg)	max. 0.002%
Potassium (K)	max. 0.01%
Sodium Carbonate (Na ₂ CO ₃)	max. 1.0%
Trace Impurities (in ppm):	
Iron (Fe)	max. 3
Mercury (Hg)	max. 0.1
Nickel (Ni)	max. 5
Nitrogen Compounds (as N)	max. 3
Phosphate (PO ₄)	max. 2

max. 5

PRODUCT	PACKING	CONT.
NO.		вох
0288.1000	1 kg	6
0288.5000	5 kg	4
0288.9025	25 kg	
0288.9050	50 kg	

Sodium Hydroxide

Pellets / 'BAKFR'

corrosive

0404

► NaOH	
M =	40.00 g/mol
CAS NO.	1310-73-2
EINECS	215-185-5
NC CODE	2815 11 00
EC NO.	11 002 00 6
UN/ID NO.	1823
ADR/RID	8 C6
IMDG	8/11
R:	35

S: 26-37/39-45

Assay	min. 97%
Aluminium (AI)	max. 0.002%
Carbonate (as Na ₂ CO ₃)	max. 1%
Chloride (CI)	max. 0.01%
Heavy Metals (as Pb)	max. 0.002%
Iron (Fe)	max. 0.002%
Sulfate (SO ₄)	max. 0.01%

Sulfate (SO₄)

PRODUCT	PACKING	CONT
NO.		ВОХ
0404.1000	1 kg	ϵ

Sodium Hydroxide



Appearance of solution	passes test
Assay	97.0-100.5%
Carbonate (as Na ₂ CO ₃)	max. 2.0%
Chlorides (as CI)	max. 50 ppm
Iron (Fe)	max. 10 ppm
Heavy Metals (as Pb)	max. 20 ppm
Identification	passes test
Potassium (K)	passes test
Sulfates (as SO ₄)	max. 50 ppm

PRODUCT NO.	PACKING	CONT. BOX
0405.5000	5 kg	4
0405.9025	25 kg	

Stored in an airtight non-metallic container.

Sodium Hydroxide

See for qualities meeting the requirements of cGMP, USP/NF, BP/Ph.Eur, JP or FCC, page 36



0405

Innovation is principal to our business.



Sodium Hydroxide

/ 'BAKER ANALYZED

► NaOH **M** = 40.00 g/mol **1 l** = 1.54 kg **CAS NO.** 1310-73-2 **EINECS** 215-185-5 NC CODE 2815 12 00 **EC NO.** 11 002 00 6 UN/ID NO. 1824 ADR/RID 8 C5 IMDG 8/II R: 35

50-52% Ammonium Hydroxide Precipitate max. 0.01% max. 0.002% Chloride (CI) Heavy Metals (as Ag) max. 0.001% Potassium (K) max. 0.01% Sodium Carbonate (Na₂CO₃) max. 0.1% Sulfate (SO₄) max. 0.001% Trace Impurities (in ppm): Iron (Fe) max. 5 Nickel (Ni) max. 5 Nitrogen Compounds (as N) max. 5 Phosphate (PO₄) max. 5

PRODUCT	PACKING	CONT.
NO.		вох
7067.1000	11	6
7067.2500	2.5	4
7067.9025	25 l	

Sodium Hydroxide

6165

S: 26-37/39-45

► NaOH **M** = 40.00 g/mol

1 I = 1.38 kg **CAS NO.** 1310-73-2 **EINECS** 215-185-5 NC CODE 2815 12 00 EC NO. 11 002 00 6 **UN/ID NO.** 1824 ADR/RID 8 C5

IMDG 8/II R: 35 **S:** 26-37/39-45

min. 40.0% Density (g/ml) at 20°C 1.43 Trace Impurities (in ppm): Iron (Fe) max. 5 Nitrogen Compounds (as N) max. 5

PRODUCT PACKING CONT. 6165.9010 10 l Polycube

Sodium Hydroxide

7109

► NaOH **M** = 40.00 g/mol 11 = 1.35 kg **CAS NO.** 1310-73-2 **EINECS** 215-185-5 NC CODE 2815 12 00 EC NO. 11 002 00 6 **UN/ID NO.** 1824 ADR/RID 8 C5

IMDG 8/II **R:** 35

S: 26-36/37/39-45

Assay (acidimetric) min. 32% Chloride (CI) max. 0.002% Silica (SiO₂) max. 0.003% Sodium Carbonate (acidimetric) max. 1% Sulfate (SO₄) max. 0.001% Trace Impurities (in ppm): Aluminium (AI) max. 5 Heavy Metals (as Pb) max. 5 Iron (Fe) max. 5 Phosphate (PO₄) max 5 Total Nitrogen (N) max. 5

PACKING PRODUCT BOX 7109.2500 2.5 [7109.5000 5 I HDPE 10 l Polycube 7109.9025 25 I

For Nitrogen determination.

Certificates of Analysis are available at www.jtbaker.com/europe

7620

Volumetric Solution, ready for use.

Sodium Hydroxide 20% / 'BAKER ANALYZED'

6164 ► NaOH min. 20% PRODUCT PACKING CONT. Assay **CAS NO.** 1310-73-2 NO. вох **EINECS** 215-185-5 6164.5000 NC CODE 2815 12 00 6164.9010 10 l Polycube **EC NO.** 11 002 00 6 **UN/ID NO.** 1824 ADR/RID 8 C5 IMDG 8/II **R:** 35

Sodium Hydroxide

S: 26-36/37/39-45

5.95-6.05 ► NaOH Molarity (M) **PRODUCT PACKING** CONT. **CAS NO.** 1310-73-2 вох **EINECS** 215-185-5 7620.5000 5 I NC CODE 2815 12 00

UN/ID NO. 1824 ADR/RID 8 C5 IMDG 8/II

R: 35

S: 26-36/37/39-45



EC NO. 11 002 00 6

Sodium Hydroxide

7598

M = 40.00 g/mol Titer (mol/l) 4.90-5.10 **PRODUCT PACKING CAS NO.** 1310-73-2 **EINECS** 215-185-5 7598.5000 5 I NC CODE 2815 12 00 **UN/ID NO.** 1824 Volumetric Solution, ready for use. ADR/RID 8 C5 IMDG 8/II **R:** 35 **S:** 26-34-37/39-45

Sodium Hydroxide

7202 4 mol/l / 'BAKER ANALYZED'

► NaOH 3.98-4.02 **PACKING** CONT. Titer (mol/l) **PRODUCT M** = 40.00 g/mol вох NO. **1 l** = 1.18 kg 7202.5000 5 I **CAS NO.** 1310-73-2 **EINECS** 215-185-5 Volumetric Solution, ready for use. NC CODE 2815 12 00

UN/ID NO. 1824 ADR/RID 8 C5 IMDG 8/II

R: 35

S: 26-34-37/39-45



Sodium Hydroxide 3 mol/I / 'BAKER ANALYZED

7621

► NaOH **CAS NO.** 1310-73-2 **EINECS** 215-185-5

> NC CODE 2815 12 00 **EC NO.** 11 002 00 6

UN/ID NO. 1824 ADR/RID 8 C5

> IMDG 8/II **R:** 35

> > **S:** 26-36/37/39-45

2.98-3.02 Molarity (M)

PRODUCT PACKING CONT. NO. вох 7621.5000 5 I

Volumetric Solution, ready for use.

Sodium Hydroxide

7036

► NaOH **M** = 40.00 g/mol

1 l = 1.09 kg **CAS NO.** 1310-73-2 **EINECS** 215-185-5

NC CODE 2815 12 00 **UN/ID NO.** 1824

ADR/RID 8 C5 IMDG 8/II

S: 26-37/39-45

R: 35

Titer (mol/l) 1.995-2.005 Trace Impurities (in ppm): Chloride (CI) max. 5 Heavy Metals (as Pb) max. 1 Iron (Fe) max. 0.5

PRODUCT	PACKING	CONT.
NO.		вох
7036.1000	11	6
7036.5000	5 l	4

Volumetric Solution, ready for use.

Each lot of this product is standardized potentiometrically against a NIST tracable reference standard.

Sodium Hydroxide

7557

► NaOH **CAS NO.** 1310-73-2

EINECS 215-185-5 NC CODE 2815 12 00

UN/ID NO. 1824 ADR/RID 8 C5 IMDG 8/II

R: 35

S: 26-37/39-45



Titer (mol/l) 1.495 - 1.505

7557.9010

PRODUCT

BOX 10 I

CONT.

Volumetric Solution, ready for use.

PACKING

Sodium Hydroxide

7097 1 mol/l / 'BAKER ANALYZED

► NaOH **M** = 40.00 g/mol **1 l** = 1.05 kg

CAS NO. 1310-73-2 **EINECS** 215-185-5 NC CODE 2815 12 00 **EC NO.** 11 002 00 6 UN/ID NO. 1824

ADR/RID 8 C5 IMDG 8/II

> **R:** 34 **S:** 26-36/37/39-45

Titer (mol/l)	0.997-1.003
Trace Impurities (in ppm):	
Chloride (CI)	max. 5
Heavy Metals (as Pb)	max. 1
Iron (Fe)	max. 0.5

PRODUCT	PACKING	CONT.
NO.		вох
7097.1000	11	6
7097.9010	10 l Polycube	4
7097 9020	20	

Volumetric Solution, ready for use.

Each lot of this product is standardized potentiometrically against NIST tracable reference standard.

J.T.Baker

Sodium Hydroxide 0.5 mol/l / 'BAKER ANALYZED'

Titer (mol/l)

► NaOH **M** = 40.00 g/mol

> **1 l** = 1.02 kg **CAS NO.** 1310-73-2 **EINECS** 215-185-5

NC CODE 2815 12 00

UN/ID NO. 1824

ADR/RID 8 C5

IMDG 8/II

R: 36/38 **S:** 26

0.4975-0.5025 Titer (mol/l) Trace Impurities (in ppm): Chloride (CI) max. 5 Heavy Metals (as Pb) max. 1 Iron (Fe) max. 0.5

PRODUCT	PACKING	CONT
NO.		ВОХ
7105.1000	11	
7105.5000	5 l Polycube	4
7105.9010	10 l Polycube	4
7105.9020	20 l Polycube	

Volumetric Solution, ready for use.

Each lot of this product is standardized potentiometrically against NIST tracable reference standard.

Sodium Hydroxide

0.357 mol/l 1/2.8 mol/l / 'BAKER ANALYZED'

7214

7105

► NaOH **M** = 40.00 g/mol

1 l = 1.02 kg

CAS NO. 1310-73-2

EINECS 215-185-5

NC CODE 2815 12 00

UN/ID NO. 1824 ADR/RID 8 C5

IMDG 8/II

R: 36/38

S: 26

PACKING CONT. PRODUCT вох NO. 7214.1000 7214.9010 10 l Polycube

Volumetric Solution, ready for use.

Sodium Hydroxide

0.33 mol/l / 'BAKER ANALYZED'

► NaOH

M = 40 **1 | =** 1.0

CAS NO. 1310-73-2

EINECS 215-185-5

NC CODE 2815 12 00

UN/ID NO. 1824

ADR/RID 8 C5

IMDG 8/II

R: 36/38 **S:** 26

	Titer (mol/l)	0.325-0.335	PRODUCT	PACKING	CONT.
0.00 g/mol			NO.		вох
.02 kg			7201.1000	11	
310-73-2			7201 9010	10	

0.355 - 0.359

Volumetric Solution, ready for use.

Sodium Hydroxide

0.25 mol/l / 'BAKER ANALYZED'

7112

7201

► NaOH	
M =	40.00 g/mol
1 l =	1.01 kg
CAS NO.	1310-73-2
EINECS	215-185-5
NC CODE	2815 12 00
UN/ID NO.	1824

ADR/RID 8 C5

IMDG 8/II

0.2485-0.2515 Titer (mol/l) Trace Impurities (in ppm): Chloride (CI) max. 5 Heavy Metals (as Pb) max. 1 max. 0.5 Iron (Fe)

PRODUCT	PACKING	CONT.
NO.		вох
7112.1000	11	
7112.9010	10 l Polycube	

Volumetric Solution, ready for use. Each lot of this product is standardized potentiometrically against NIST tracable reference standard.



Sodium Hydroxide 7216

0.2 mol/l / 'BAKER ANALYZED

1 l = 1.01 kg **CAS NO.** 1310-73-2 **EINECS** 215-185-5 NC CODE 2815 12 00 **UN/ID NO.** 1824 ADR/RID 8 C5 IMDG 8/II R: 36/38 **S**: 26

0.199-0.201 Titer (mol/l)

PRODUCT CONT. **PACKING** NO. вох 7216.1000 11 7216.9020 20 l Polycube

Volumetric Solution, ready for use.

Each lot of this product is standardized potentiometrically against NIST tracable reference standard.

Sodium Hydroxide

7098

► NaOH **M** = 40.00 g/mol **1 l** = 1.00 kg **CAS NO.** 1310-73-2 **EINECS** 215-185-5 NC CODE 2815 12 00 UN/ID NO. 1824 ADR/RID 8 C5 IMDG 8/II

Titer (mol/l) 0.0997-0.1003 Trace Impurities (in ppm): Chloride (CI) max. 5 Heavy Metals (as Pb) max. 1 Iron (Fe) max. 0.5

PRODUCT PACKING CONT. вох 11 7098.1000 7098.5000 5 | Jerrycan 7098.9010 10 l Polycube 7098.9020 20 l Polycube

Volumetric Solution, ready for use. Each lot of this product is standardized potentiometrically against NIST tracable reference standard.

Sodium Hydroxide

7480

► NaOH M = 40.00 g/mol **CAS NO.** 1310-73-2 **EINECS** 215-185-5 NC CODE 2815 12 00 **UN/ID NO.** 1824 ADR/RID 8 C5

IMDG 8/III

Titer (mol/l) 0.0495-0.0505 **PRODUCT PACKING** CONT. BOX 7480 1000 11

Volumetric Solution, ready for use. Each lot of this product is standardized potentiometrically against NIST tracable reference standard.

Sodium Hydroxide 7475

0.0357 mol/l / 'BAKER ANALYZED'

M = 40.00 g/mol **CAS NO.** 1310-73-2 **EINECS** 215-185-5 NC CODE 2815 12 00 **UN/ID NO.** 1824 ADR/RID 8 C5

IMDG 8/III

Titer (mol/l) 0.035 - 0.037 **PRODUCT PACKING** CONT. BOX 7475.9010 10 l Polycube

Volumetric Solution, ready for use.

Sodium Hydroxide 7099

0.02 mol/l / 'BAKER ANALYZED'

► NaOH **M** = 40.00 g/mol **CAS NO.** 1310-73-2 **EINECS** 215-185-5 NC CODE 2815 12 00 **UN/ID NO.** 1824 ADR/RID 8 C5 IMDG 8/III

0.0195-0.0205 Titer (mol/l)

PRODUCT PACKING CONT. NO. BOX 7099.0500 500 ml 20 l Polycube 7099.9020

Volumetric Solution, ready for use.

Each lot of this product is standardized potentiometrically against NIST tracable reference standard.

Sodium Hydroxide 5 mol/I / DILUT-IT / 5 equiv. = 200.0g; 5N

4690

► NaOH PRODUCT PACKING CONT. **M** = 40.00 g/mol вох NO. **CAS NO.** 1310-73-2 4690 500 ml Sealed Bottle **EINECS** 215-185-5 NC CODE 2815 12 00 Volumetric Concentrate, for dilution to 1 l. **EC NO.** 11 002 00 6 **UN/ID NO.** 1824 ADR/RID 8 C5 IMDG 8/II **R:** 35 **S:** 26-36/37/39-45

Sodium Hydroxide

= 40.00g, 1N / DILUT-IT

	2.7.7						
NaOH			PR	RODUCT	PACKING		CON
M =	40.00 g/mol		NO) .			ВО
CAS NO.	1310-73-2		46	89	1 amp.		
EINECS	215-185-5						
NC CODE	2815 12 00		Vo	olumetric (Concentrate, for dilu	ition to 1 l.	
UN/ID NO.	1824						
ADR/RID	8 C5						
IMDG	8/II						
R:	35						
S:	26-36/37/39-45						
	C C C C Corrosive						

Sodium Hydroxide 0.5 mol/I / DILUT-IT

4691

► NaOH **PRODUCT PACKING** CONT. **M** = 40.00 g/mol вох **CAS NO.** 1310-73-2 4691 1 amp. **EINECS** 215-185-5 NC CODE 2815 12 00 Volumetric Concentrate, for dilution to 1 l. **UN/ID NO.** 1824 ADR/RID 8 C5 IMDG 8/II **R:** 35 **S:** 26-36/37/39-45

Sodium Hydroxide

4865

► NaOH **M =** 40.00 g/mol NC CODE 2815 12 00 **UN/ID NO.** 1824 ADR/RID 8 C5 IMDG 8/II **R:** 35 **S:** 26-36/37/39-45

PRODUCT	PACKING	CONT
NO.		ВОХ
4865	1 amp.	ϵ

Volumetric Concentrate, for dilution to 1 l.



Sodium Hydroxide

4687 0.1 mol/l, $\frac{1}{10}$ equiv. = 4.000g, 0.1 N / DILUT-IT

M = 40.00 g/mol

CAS NO. 1310-73-2

EINECS 215-185-5

NC CODE 2815 12 00

UN/ID NO. 1824

ADR/RID 8 C5

IMDG 8/II R: 35

S: 26-34-37/39-45



Sodium Hydroxide

4715 0.02 mol/l / DILUT-IT / 1/50 equiv. = 0.800g; 0.02N

► NaOH **M** = 40.00 g/mol **CAS NO.** 1310-73-2

EINECS 215-185-5 NC CODE 2815 12 00 UN/ID NO. 1824

ADR/RID 8 C5 IMDG 8/II R: 36/38

S: 26

irritant)

 PRODUCT
 PACKING
 CONT.

 NO.
 BOX

 4715
 1 amp.
 6

Volumetric Concentrate, for dilution to 1 l.

PRODUCT

NO.

4687

PACKING

1 amp.

Volumetric Concentrate, for dilution to 1 l.

CONT.

вох

Sodium Hydroxide

4864 0.01 mol/l / DILUT-IT

■ NaOH

M = 40.00 g/mol

NC CODE 2815 12 00

UN/ID NO. 1824

ADR/RID 8 C5

IMDG 8/II

R: 36/38

S: 26

40.00 g/mol 2815 12 00 1824 3 C5 3/II
 PRODUCT NO.
 PACKING BOX

 4864
 1 amp.
 6

Volumetric Concentrate, for dilution to 1 l.

Sodium Hydroxide Solutions

See for qualities meeting the requirements of cGMP, USP/NF, BP/Ph.Eur, JP or FCC, page 36

Sodium Hydroxide on support

Granular 1.6-3 mm / 'BAKER ANALYZED

EINECS 215-185-5 NC CODE 3822 00 00 EC NO. 11 002 00 6

EC NO. 11 002 00 6 UN/ID NO. 1823

ADR/RID 8 C6 **IMDG** 8/II **R:** 35

S: 26-37/39-45



 ${\rm CO_2}$ absorption capacity min. 50% Loss on Ignition 18-25%

 PRODUCT
 PACKING

 NO.
 0458.0500

 500 g
 500 g

CONT.

вох

Sodium Hypochlorite 5% Solution, Sterile

See for qualities meeting the requirements of cGMP, USP/NF, BP/Ph.Eur, JP or FCC, page 36

0293

0294



XC

Sodium Iodide 'RAKER ANAI YZED

► Nal Assay (dried basis) min. 99.0% **M** = 149.89 g/mol Barium (Ba) max. 0.002% **CAS NO.** 7681-82-5 Chloride and Bromide (as CI) max. 0.01% **EINECS** 231-679-3 Heavy Metals (as Pb) max. 5 ppm NC CODE 2827 60 00 Insoluble Matter max. 0.005%

lodate (IO₃) max. 3 ppm Iron (Fe) max. 5 ppm Nitrogen Compounds (as N) max. 0.001% pH of 5% Solution at 25°C 7.5-9.5 Phosphate (PO₄) max. 0.002% Sulfate (SO₄) max. 0.005%

ВО

Sodium Iodide

► Nal

99.0-100.5% Assay **M** = 149.89 g/mol Alkalinity passes test CAS NO. 7681-82-5 Appearance of solution **EINECS** 231-679-3 Heavy Metals (as Pb) NC CODE 2827 60 00 Identification

passes test max. 10 ppm passes test lodate (IO₃) passes test Iron (Fe) max. 20 ppm Limit of nitrate, nitrite and ammonia passes test Loss on Drying max. 3.0% Potassium (K) passes test Sulfate (SO₄) max. 150 ppm Thiosulfates (as S₂O₃) passes test Water (H₂O) max. 2.0%

PACKING	CONT
	ВОХ
250 g	
500 g	ϵ
1 kg	ϵ
50 kg	
	250 g 500 g 1 kg

Preserve in tight containers Stored protected from light.

Sodium Iodide

e requirements of cGMP, USP/NF, BP/Ph.Eur, JP or FCC, page 36



Sodium Lactate

	60 %	/	syrup	/	'BAKER
► CH₃CHOHCO	ONa				

M = 112.06 g/m**1 l** = 1.28 kg **CAS NO.** 72-17-3 **EINECS** 200-772-0 NC CODE 2918 11 00

	Assay	59-61%	PRODUCT	PACKING	CONT.
nol			NO.		BOX
			7391.0500	500 ml	
			7391.2500	2.5	
			7391.9025	25 l	

Sodium meta-Bisulfite



Sodium meta-Bisulfite

See for qualities meeting the requirements of cGMP, USP/NF, BP/Ph.Eur, JP or FCC, page 36



Sodium Metaperiodate



Find more Chromatography information at www.jtbaker.com/chromatography



Sodium Metasilicate Pentahydrate

0302 'BAKER ANALYZED'

Na₂SiO₃·5H₂O

M = 212.20 g/mol

CAS NO. 6834-92-0

EINECS 229-912-9

NC CODE 2839 11 00

EC NO. 14 010 00 8

UN/ID NO. 3253 ADR/RID 8 C6

> IMDG 8/III R: 34-37

> > **S:** 13-24/25-36/37/39-45

corrosive

Assay	min. 93%
Appearance of 3% Solution	clear
Chloride (CI)	max. 0.05%
Heavy Metals (as Pb)	max. 0.001%
Iron (Fe)	max. 0.02%
pH of 3% Solution at 25°C	min. 12.3
Sulfate (SO ₄)	max. 0.01%

0302.0250 250 g 0302.5000 5 kg 4	PRODUCT	PACKING	CONT. BOX
0302.5000 5 kg	0302.0250	250 g	
	0302.5000	5 kg	4

Sodium Molybdate Dihydrate 'BAKER ANALYZED' / ACS

► Na₂MoO₄.2H₂O

M = 241.95 g/molCAS NO. 10102-40-6EINECS 231-551-7NC CODE 2841 70 00

Meets ACS Specifications Assay 99.5-103.0% Ammonium (NH₄) max. 0.001% Chloride (CI) max. 0.005% Insoluble Matter max. 0.005% max. 0.001% Iron (Fe) pH of 5% Solution at 25°C 7.0-10.5 Sulfate (SO₄) max. 0.015% Trace Impurities (in ppm): Heavy Metals (as Pb) max. 5 Phosphate (PO₄) max. 5

	PACKING	CONT.
١0.		вох
)295.0250	250 g	
295.1000	1 kg	
0295.9050	50 kg	

Sodium Nitrate

0296 'BAKER ANALYZED' / ACS

► NaNO₃

M = 84.99 g/mol **CAS NO.** 7631-99-4

EINECS 231-554-3

NC CODE 3102 50 90

UN/ID NO. 1498

ADR/RID 5.1 O2

IMDG 5.1/III

R: 8 **S:** 17-41

oxidizing

Calcium (Ca)	max. 0.005%
Chloride (CI)	max. 0.001%
Insoluble Matter	max. 0.005%
Magnesium (Mg)	max. 0.002%
Nitrite (NO ₂)	max. 0.001%
pH of 5% Solution at 25°C	5.5-8.3
Sulfate (SO ₄)	max. 0.002%
Trace Impurities (in ppm):	
Heavy Metals (as Pb)	max. 5
lodate (IO ₃)	max. 5
Iron (Fe)	max. 2

min. 99.0%

max. 5

max. 1.0%

Exceeds ACS Specifications

	ВОХ
00 g	6
kg	
0 kg	
	00 g kg

2017

Sodium Nitrate

BAKER'

► NaNO₃

M = 84.99 g/mol

CAS NO. 7631-99-4

EINECS 231-554-3 **NC CODE** 3102 50 90

UN/ID NO. 1498

ADR/RID 5.1 O2

IMDG 5.1/III

R: 8

S: 17-41

oxidizing

PRODUCT	PACKING	CONT.
NO.		вох
2017.1000	1 kg	
2017.9050	50 kg	



Sodium Nitrate

See for qualities meeting the requirements of cGMP, USP/NF, BP/Ph.Eur, JP or FCC, page 36

Phosphate (PO₄)

Loss on Drying

0297

0298

J.T.Baker

Sodium Nitrite

'BAKER ANALYZED' / ACS

► NaNO₂

M = 69.00 g/mol

CAS NO. 7632-00-0

EINECS 231-555-9

NC CODE 2834 10 00

EC NO. 7 010 00 4

UN/ID NO. 1500

ADR/RID 5.1 OT2

IMDG 5.1/III

R: 25-50-8

S: 45-61







Specifications for testing USP/NF monographs	
Assay	min. 97.0%
Calcium (Ca)	max. 0.01%
Chloride (CI)	max. 0.005%
Heavy Metals (as Pb)	max. 0.001%
Insoluble Matter	max. 0.01%
Potassium (K)	max. 0.005%
Sulfate (SO ₄)	max. 0.01%
Trace Impurities (in ppm):	
Arsenic (As)	max. 3
Iron (Fe)	max. 5

Exceeds ACS Specifications. Meets Reagents

PRODUCT	PACKING	CONT.
NO.		вох
0297.0100	100 g	
0297.1000	1 kg	
0297.9025	25 kg	
0297.9050	50 kg	

Sodium Nitrite

► NaNO₂

M = 69.00 g/mol

CAS NO. 7632-00-0

EINECS 231-555-9

NC CODE 2834 10 00

EC NO. 7 010 00 4

UN/ID NO. 1500 ADR/RID 5.1 OT2

IMDG 5.1/III

R: 25-50-8

S: 45-61







PRODUCT PACKING CONT. вох NO. 0298.1000 1 kg 0298.9050 50 kg

Stored in an airtight container.

Sodium Nitrite

ting the requirements of cGMP, USP/NF, BP/Ph.Eur, JP or FCC, page 36

Assay (dried basis)

Acidity or Alkalinity

Chloride (CI)

Identification

Loss on Drying

Sulfates (as SO₄)

Appearance of solution

Heavy Metals (as Pb)



Sodium Nitroferrricyanide

See Sodium Pentacyanonitrosylferrate(III) Dihydrate



1916 'BAKER'

► CH₃(CH₂)₇CH:CH(CH₂)₇COONa

M = 304.45 g/mol **CAS NO.** 143-19-1

EINECS 205-591-0 NC CODE 2916 15 00 lodine number of free acid 85.0-95.0

98.5-100.5%

passes test

passes test

max. 50 ppm

max. 20 ppm

passes test

max. 0.25%

max. 200 ppm

NO. 1916.0250 250 g

PACKING

PRODUCT

CONT. BOX

0301

Sodium Oxalate

'BAKER ANALYZED' / ACS

► NaOCOCOONa

M = 134.00 g/mol **CAS NO.** 62-76-0

EINECS 200-550-3

NC CODE 2917 11 00

EC NO. 607 007 00 3

UN/ID NO. 3282

ADR/RID 6.1 T3

IMDG 6.1/III

R: 21/22 **S:** 24/25

Xn

Exceeds ACS Specifications Assay (by KMnO ₄ titrn.)	min. 99.8%
Ammonium (NH ₄)	max. 0.002%
Chloride (CI)	max. 0.002%
Heavy Metals (as Pb)	max. 0.002%
Insoluble Matter	max. 0.005%
Iron (Fe)	max. 0.001%
Loss on Drying at 105°C	max. 0.01%
Neutrality	passes test
Potassium (K)	max. 0.005%
Substances Darkened by Hot H ₂ SO ₄	passes test
Sulfate (SO ₄)	max. 0.002%

PRODUCT	PACKING	CONT
NO.		BOX
0301.0250	250 g	
0301.9050	50 kg	



Sodium Oxalate

'BAKER'

2018 'B
➤ NaOCOCOONa

M = 134.00 g/mol **CAS NO.** 62-76-0

EINECS 200-550-3

NC CODE 2917 11 00

EC NO. 607 007 00 3

UN/ID NO. 3282

ADR/RID 6.1 T3

IMDG 6.1/III

R: 21/22 **S:** 24/25

X

Assay	min. 99.0%
Chloride (CI)	max. 0.005%
Heavy Metals (as Pb)	max. 0.01%
Iron (Fe)	max. 0.005%
Sulfate (SO.)	max 0.05%

PRODUCT	PACKING	CONT.
NO.		вох
2018.1000	1 kg	

Sodium Pentacyanonitrosylferrate(III) Dihydrate

1190 'BAKER ANALYZED' / ACS

Na₂Fe(CN)₅NO.2H₂O **M** = 297.95 g/mol **CAS NO.** 13755-38-9

EINECS 238-373-9 **NC CODE** 2837 20 00

UN/ID NO. 1588 **ADR/RID** 6.1 T5

IMDG 6.1/II R: 25

S: 22-37-45

toxic T

Exceeds ACS Specifications	
Assay	99.0-102.0%
Chloride (CI)	max. 0.02%
Insoluble Matter	max. 0.01%
pH of 5% Solution at 25°C	4.0-7.0
Sulfate (SO ₄)	max. 0.01%

PRODUCT	PACKING	CONT.
NO.		вох
1190.0100	100 g	

Sodium Perborate, 4-Hydrate

See for qualities meeting the requirements of cGMP, USP/NF, BP/Ph.Eur, JP or FCC, page 36

Sodium Perchlorate Monohydrate

2815 'BAKER ANALYZED' / ACS

► NaClO₄·H₂O

M = 140.46 g/mol

CAS NO. 7791-07-3

EINECS 231-511-9

NC CODE 2829 90 10


EC NO. 17 010 00 6

UN/ID NO. 1502

ADR/RID 5.1 O2 **IMDG** 5.1/II **R:** 22-9

S: 13-22-27

harmful Xn



PRODUCI	PACKING	CONI.
NO.		вох
2815.0100	100 g	6
2815.1000	1 kg	6



Sodiummeta-Periodate

See Sodium Tetroxoiodate(VII)

www.jtbaker.com/europe

Sodium Peroxide

Granular 0.3-2.0 mm / 'BAKER ANALYZED' / ACS

9015

► Na ₂ O ₂	
M =	77.98 g/mol
CAS NO.	1313-60-6
EINECS	215-209-4
NC CODE	2815 30 00
EC NO.	11 003 00 1
UN/ID NO.	1504
ADR/RID	5.1 02
IMDG	5.1/I
R:	35-8
S:	27-39-45-8
	c O

Meets ACS Specifications	
Assay	min. 93.0%
Chloride (CI)	max. 0.002%
Heavy Metals (as Pb)	max. 0.002%
Iron (Fe)	max. 0.005%
Sulfate (SO ₄)	max. 0.001%
Trace Impurities (in ppm):	
Phosphate (PO ₄)	max. 5

IG CONT.
вох
6

Sodium Phosphate Dodecahydrate

'BAKER ANALYZED' / ACS

0307

la ₃ PO ₄ . 12H ₂ 0)
M =	380.12 g/mol
CAS NO.	10101-89-0
EINECS	231-509-8
NC CODE	2835 23 00
R:	36/38
S:	26
	Xi

► N

Meets ACS Specifications	
Assay	98.0-102.0%
Chloride (CI)	max. 0.001%
Excess Alkali (as NaOH)	max. 2.5%
Heavy Metals (as Pb)	max. 0.001%
Insoluble Matter	max. 0.01%
Iron (Fe)	max. 0.001%
Sulfate (SO ₄)	max. 0.01%

min. 98%

PRODUCT

BO 2
000 1 kg
000 I Kg
050 50 kg
050 50 kg

Sodium Phosphate Dodecahydrate

Na₃PO₄.12H₂O

M = 380.12 g/mol **CAS NO.** 10101-89-0

EINECS 231-509-8

NC CODE 2835 23 00

R: 36/38

S: 26



NO.		вох
0308.1000	1 kg	6
0308.5000	5 kg	
0308.9050	50 kg	

PACKING

Sodium Phosphate, Dibasic, Anhydrous

See for qualities meeting the requirements of cGMP, USP/NF, BP/Ph.Eur, JP or FCC, page 36

Assay



0308

CONT.

Sodium Phosphate, Dibasic, 7-Hydrate

See for qualities meeting the requirements of cGMP, USP/NF, BP/Ph.Eur, JP or FCC, page 36



Sodium Phosphate, Monobasic, Monohydrate

See for qualities meeting the requirements of cGMP, USP/NF, BP/Ph.Eur, JP or FCC, page 36



Sodium Pyrophosphate See Sodium Diphosphate Decahydrate



Sodium Pyrosulfite

See Sodium Disulfite

NC CODE 2918 30 00



2456

Sodium Pyruvate 'BAKER'

 ► CH₃COCOONa
 Assay (by Perchloric Acid titrn.)
 min. 98%
 PRODUCT
 PACKING
 CONT.

 M = 110.05 g/mol
 Water (H₂O)
 max. 2%
 NO.
 BOX

 CAS NO. 113-24-6
 2456.0100
 100 g

 EINECS 204-024-4
 2456.0100
 100 g

Sodium Rhodanide

ee Sodium Thiocvanate

2028

Sodium Salicylate

'BAKER ANALYZED

► 2-HOC₆H₄COONa

M = 160.11 g/mol

CAS NO. 54-21-7 **EINECS** 200-198-0

NC CODE 2918 21 00

R: 22

S: 24/25



 Assay (after drying)
 min. 99.5%

 Chloride (CI)
 max. 0.002%

 Heavy Metals (as Pb)
 max. 0.001%

 Iron (Fe)
 max. 0.001%

 Sulfate (SO₄)
 max. 0.02%

 Water (H₂O)
 max. 0.2%

 PRODUCT NO.
 PACKING PACKING
 CONT.

 2028.0250
 250 g
 6



Sodium Salicylate

See for qualities meeting the requirements of cGMP, USP/NF, BP/Ph.Eur, JP or FCC, page 36

0620 0620

Sodium Succinate Hexahydrate

620 'BAKE

► NaOCOCH₂CH₂COONa.6H₂O

M = 270.15 g/mol **CAS NO.** 150-90-3

FINECS 205-778-7

NC CODE 2917 19 90

Water (H₂O)

38-42%

max. 0.5%

max. 5

max. 5

max. 0.005%

PRODUCT PACKING NO.

0620.0500 500 g 0620.5000 5 kg CONT. BOX

S 12

Sodium Sulfate Anhydrous

12-60 mesh / 'BAKER ULTRA RESI-ANALYZED' / for Organic Residue Analysis / ACS

Loss on Ignition

Magnesium (Mg)

► Na₂SO₄

M = 142.04 g/mol

CAS NO. 7757-82-6

EINECS 231-820-9

NC CODE 2833 11 00

Meets ACS Specifications	
Assay	min. 99.0%
Calcium (Ca)	max. 0.01%
Chloride (CI)	max. 0.001%
Extraction-Concentration Suitability	passes test
Insoluble Matter	max. 0.01%
Iron (Fe)	max. 0.001%

pH of 5% Solution at 25°C 5.2-9.2 Phosphate (as PO $_4$) max. 0.001% Potassium (K) max. 0.002%

Trace Impurities (in ppm): Heavy Metals (as Pb)

Nitrogen Compounds (as N)

PRODUCT	PACKING	CONT.
NO.		вох
3375.0500	500 g	
3375.1000	1 kg	

Sodium Sulfate Anhydrous

3377-00 12-60 mesh / 'BAKER ANALYZED' / For Hydrocarbon Oil Index determination

► Na₂SO₄

M = 142.04 g/mol

CAS NO. 7757-82-6 EINECS 231-820-9 NC CODE 2833 11 00 Assay min. 99.0% Insoluble Matter max. 0.01% Loss on Ignition max. 0.5%

 PRODUCT
 PACKING
 CONT.

 NO.
 BOX

 3377-00
 100 g

Suitable for determination of Hydrocarbon Oil Index according to ISO 9377-2 and NEN 5733.

After opening: Store in desiccator.

Innovation is principal to our business.

Sodium Sulfate Anhydrous 'BAKER ANALYZED'

0313

► Na₂SO₄ **M** = 142.04 g/mol

CAS NO. 7757-82-6 **EINECS** 231-820-9 NC CODE 2833 11 00

Assay	min. 99.0%
Calcium, Magnesium and R ₂ O ₃ Precipitate	max. 0.01%
Chloride (CI)	max. 0.001%
Insoluble Matter	max. 0.01%
Loss on Ignition	max. 0.5%
pH of 5% Solution at 25°C	5.2-9.2
Potassium (K) (by AAS)	max. 0.01%
Trace Impurities (in ppm):	
Arsenic (As)	max. 0.5
Heavy Metals (as Pb)	max. 5
Iron (Fe)	max. 5
Nitrogen Compounds (as N)	max. 5

PRODUCT NO.	PACKING	CONT. BOX
0313.0500	500 g	
0313.1000	1 kg	6
0313.5000	5 kg	4
0313.9025	25 kg	
0313.9050	50 kg	

Sodium Sulfate Anhydrous

'BAKER ANALYZED' / ACS

► Na₂SO₄ **M** = 142.04 g/mol CAS NO. 7757-82-6 **EINECS** 231-820-9 NC CODE 2833 11 00

Meets ACS Specifications	
Assay	min. 99.0%
Calcium (Ca)	max. 0.01%
Chloride (CI)	max. 0.001%
Insoluble Matter	max. 0.01%
Iron (Fe)	max. 0.001%
Loss on Ignition	max. 0.5%
Magnesium (Mg)	max. 0.005%
pH of 5% Solution at 25°C	5.2-9.2
Phosphate (as PO ₄)	max. 0.001%
Potassium (K)	max. 0.002%
Trace Impurities (in ppm):	
Heavy Metals (as Pb)	max. 5
Nitrogen Compounds (as N)	max. 5

PRODUCT	PACKING	CONT. BOX
0312.0100	100 g	
0312.1000	1 kg	6
0312.9050	50 kg	

Sodium Sulfate, Anhydrous

See for qualities meeting the requirements of cGMP, USP/NF, BP/Ph.Eur, JP or FCC, page 36

Iron (Fe)

Nitrogen Compounds (as N)



0311

0312

Sodium Sulfate Decahydrate

'BAKER ANALYZED'

Na₂SO₄. 10H₂O **M =** 322.19 g/mol **CAS NO.** 7727-73-3

EINECS 231-820-9 NC CODE 2833 11 00

Assay	min. 98.0%
Calcium and Magnesium (as Ca)	max. 0.005%
Chloride (CI)	max. 0.001%
nsoluble Matter	max. 0.005%
oH of 5% Solution at 25°C	4.5-6.5
Trace Impurities (in ppm):	
Arsenic (As)	max. 0.5
Heavy Metals (as Pb)	max. 5

111ax. 0.005%
4.5-6.5
max. 0.5
max. 5
max. 5
max. 5

PRODUCT	PACKING	CONT.
NO.		вох
0311.1000	1 ka	

Sodium Sulfate Decahydrate 'BAKER'

Na₂SO₄.10H₂O

M = 322.19 g/mol **CAS NO.** 7727-73-3 **EINECS** 231-820-9 NC CODE 2833 11 00

Assay (dried basis) 98.5-101.0% Acidity or Alkalinity passes test Appearance of solution passes test Chloride (CI) max. 200 ppm Heavy Metals (as Pb) max. 20 ppm Identification passe test Loss on Drying 52.0-57.0%

PRODUCT	PACKING	CONT.
NO.		BOX
1774.5000	5 kg	
1774.9050	50 kg	

Sodium Sulfate, 10-Hydrate

See for qualities meeting the requirements of cGMP, USP/NF, BP/Ph.Eur, JP or FCC, page 36



1774



Sodium Sulfite Anhydrous 'BAKER ANALYZED' / ACS

0316

Na₂SO₃ **M** = 126.04 g/mol **CAS NO.** 7757-83-7 **EINECS** 231-821-4

NC CODE 2832 10 00

Meets ACS Specifications. Meets Reagent Specifications for testing USP/NF monographs min. 98.0% Assay (by iodometry) Appearance passes test Chloride (CI) max. 0.02% Free Acid passes test Heavy Metals (as Pb) max. 0.001% Insoluble Matter max. 0.005% Titrable Free Base (meq/g) max. 0.03 Trace Impurities (in ppm): Iron (Fe) max. 5

PRODUCT	PACKING	CONT.
NO.		вох
0316.0250	250 g	6
0316.1000	1 kg	6
0316.9012PE	12 kg PE Pail	

Sodium Sulfite Anhydrous

1777

► Na₂SO₃

M = 126.04 g/mol **CAS NO.** 7757-83-7 **EINECS** 231-821-4 NC CODE 2832 10 00

Assay min. 95% Chloride (CI) max. 0.02% max. 0.005% Heavy Metals (as Pb) max. 0.002% Iron (Fe) Trace Impurities (in ppm): Arsenic (As) max. 2

PRODUCT	PACKING	CONT.
NO.		вох
1777.1000	1 kg	6

0317

Sodium Tartrate Dihydrate

► NaOCO(CHOH)₂COONa.2H₂O M = 230.08 g/mol **CAS NO.** 6106-24-7 **EINECS** 212-773-3 NC CODE 2918 13 00

Exceeds ACS Specifications Assay 99.0-101.0% Ammonium (NH₄) max. 0.002% Calcium (Ca) max. 0.01% Insoluble Matter max. 0.005% 15.61-15.71% Loss on Drying at 150°C pH of 5% Solution at 25°C 7.0-9.0 Sulfate (SO₄) max. 0.005% Trace Impurities (in ppm): Chloride (CI) max. 5 Heavy Metals (as Pb) max. 5 Iron (Fe) max. 5 Phosphate (PO₄) max. 5

PACKING	CONT.
	вох
250 g	
1 kg	
50 kg	
	250 g 1 kg

This reagent is suitable for standardization of Karl Fischer Reagent as used for the determination of water.



Sodium Tetraborate



Sodium Tetrahydroborate

See Sodium Borohydride

Sodium Tetraphenylborate

1194 'BAKER ANALYZED' / ACS

NaB(C₆H₅)₄ **M** = 342.23 g/mol **CAS NO.** 143-66-8 **EINECS** 205-605-5 NC CODE 2931 00 95

R: 22

ivieets ACS Specifications	
Assay	min. 99.5%
Clarity of Solution	passes test
Loss on Drying at 105°C	max. 0.5%
,	

PRODUCT	PACKING	CONT.
NO.		вох
1194.0025	25 g Glass	
1194.0100	100 g	



Sodium Tetraphenylboron

e Sodium Tetraphenylbora

Sodium Tetroxoiodate(VII)

BAKER ANALYZED 1192

► NalO₄

M = 213.89 g/mol

CAS NO. 7790-28-5 **EINECS** 232-197-6

NC CODE 2829 90 80

UN/ID NO. 1479

ADR/RID 5.1 02

IMDG 5.1/I

R: 8

S: 17



Assay	min. 99.8%
Other Halogens (as CI)	max. 0.01%
Sulfate (SO ₄)	max. 0.005%
Trace Impurities (in ppm):	
Manganese (Mn)	max. 1

PRODUCT	PACKING	CONT.
NO.		вох
1192.0100	100 g	

Sodium Thiocyanate

'BAKER ANALYZED'

Iron (Fe)

Sulfide (S)

► NaSCN

M = 81.07 g/mol

CAS NO. 540-72-7

EINECS 208-754-4 NC CODE 2838 00 00

EC NO. 615 004 00 3

R: 20/21/22-32

S: 13



Exceeds ACS Specifications	
Assay	min. 98.0%
Ammonium (NH ₄)	max. 0.002%
Carbonate (as Na ₂ CO ₃)	max. 0.2%
Chloride (CI)	max. 0.005%
Insoluble Matter	max. 0.005%
Sulfate (SO ₄)	max. 0.01%
Sulfide (S)	max. 0.001%
Trace Impurities (in ppm):	
Heavy Metals (as Pb)	max. 5

max. 2

max. 1

passes test

32.0-37.0%

PRODUCT	PACKING	CONT.
NO.		вох
0318.0500	500 g	
0318.7110	110 lbs	

0318

0320

0319

1779

CONT.

BOX

Sodium Thiosulfate Anhydrous

BAKER ANALYZED

► Na₂S₂O₃

M = 158.11 g/mol

CAS NO. 7772-98-7

EINECS 231-867-5

NC CODE 2832 30 00

Assay (by lodometry) min 98.0% max. 0.01% Calcium (Ca) Heavy Metals (as Pb) max. 0.001% max. 0.01% Insoluble Matter Trace Impurities (in ppm):

PACKING CONT. PRODUCT вох NO. 0320.0500 500 g 50 kg 0320.9050

Sodium Thiosulfate Pentahydrate

'BAKER ANALYZED' / ACS

Na₂S₂O₃.5H₂O

M = 248.18 g/mol

NC CODE 2832 30 00

CAS NO. 10102-17-7 FINECS 231-867-5

Meets ACS Specifications		
Assay	99.5-101.0%	
Insoluble Matter	max. 0.005%	
Nitrogen Compounds (as N)	max. 0.002%	
pH of 5% Solution at 25°C	6.0-8.4	
Sulfate and Sulfite (as SO ₄)	max. 0.1%	
Sulfide (S)	may 1 nnm	

PRODUCT	PACKING	CONT.
NO.		вох
0319.0500	500 g	6
0319.1000	1 kg	6
0319.9050	50 kg	

Sodium Thiosulfate Pentahydrate

Na₂S₂O₃.5H₂O

M = 248.18 g/mol

CAS NO. 10102-17-7

EINECS 231-867-5

NC CODE 2832 30 00

99.0-101.0% Appearance of solution passes test Calcium (Ca) passes test Heavy Metals (as Pb) max. 10 ppm Identification passes test 6.0-8.4 Sulfates and Sulfites max. 0.2%

PRODUCT PACKING 1 kg 1779.1000 1779.5000 5 kg 50 kg 1779.9050 Preserve in tight containers.

Stored in an airtight container.

Sodium Thiosulfate, 5-Hydrate

ee for qualities meeting the requirements of cGMP, USP/NF, BP/Ph.Eur, JP or FCC, page 36

Water (H₂O)

Sulfides





Sodium Thiosulfate

7100 1 mol/L / 'BAKER ANALYZED

Na₂S₂O₃ **M** = 158.11 g/mol **CAS NO.** 10102-17-7 **EINECS** 231-867-5 NC CODE 2832 30 00

0.995-1.005 Titer (mol/l) Trace Impurities (in ppm): Heavy Metals (as Pb) max. 1 Sulfate and Sulfite (as SO₄) max. 50 PRODUCT CONT. **PACKING** NO. вох 7100.1000 11 7100.9020 20 l Polycube

Volumetric Solution, ready for use. Each lot of this product is standardized potentiometrically against NIST tracable reference standard.

Sodium Thiosulfate 7101

Na₂S₂O₃ **M** = 158.11 g/mol **CAS NO.** 1010-21-7 **EINECS** 231-867-5 NC CODE 2832 30 00

0.0997-0.1003 Titer (mol/l) Trace Impurities (in ppm): Heavy Metals (as Pb) max. 1 Sulfate and Sulfite (as SO₄) max. 50 **PRODUCT PACKING** CONT. BOX 7101.1000 11 5 l Polycube 7101.5000 7101.9010 10 l Polycube 20 l Polycube 7101.9020

Volumetric Solution, ready for use. Each lot of this product is standardized potentiometrically against NIST tracable reference standard.

Sodium Thiosulfate 4695 0.1 mol/l / DILUT-IT

► Na₂S₂O₃

M = 158.11 g/mol **CAS NO.** 10102-17-7 **EINECS** 231-867-5 NC CODE 2832 30 00

PRODUCT PACKING CONT. NO. BOX 4695 1 amp. 6

Volumetric Concentrate, for dilution to 1 l.

Volumetric Concentrate, for dilution to 1 l.

Sodium Thiosulfate 4693

Na₂S₂O₃

M = 158.11 g/mol CAS NO. 10102-17-7 **EINECS** 231-867-5 NC CODE 2832 30 00

PRODUCT PACKING CONT. NO. BOX 4693 1 amp. 6

Sodium Tungstate Dihydrate 1195

'BAKER ANALYZED' / ACS

► Na₂WO₄.2H₂O

M = 329.86 g/mol **CAS NO.** 10213-10-2 FINECS 236-743-4 NC CODE 2841 80 00

99.0-101.0%
max. 0.005%
max. 0.001%
max. 0.01%
max. 0.001%
max. 0.01%
max. 0.02 meq/g

PRODUCT PACKING CONT. вох 1195.0100 100 g 1195,9050 50 kg

Sodium Tungstate Dihydrate

► Na₂WO₄.2H₂O

M = 329.86 g/mol **CAS NO.** 10213-10-2 **EINECS** 236-743-4 NC CODE 2841 80 00 R: 22



Assay min. 99.0% Alkalinity (as Na₂CO₃) max. 0.2% Arsenic (As) max. 0.002% Chloride (CI) max. 0.01% Heavy Metals and Iron (as Pb) max. 0.001% max. 0.02% Insoluble Matter max. 0.05% Molybdenum (Mo) pH of 5% Solution at 25°C 8.0-10.0 Sulfate (SO₄) max. 0.02%

PRODUCT PACKING CONT. BOX 0321.0500 500 g 0321.9050 50 kg



Sodium Wolframate

See Sodium Tungstate Dihydrate



Solid Phase Extraction

See for detailed information section Chromatography products, page 448



2039

D(-)-Sorbitol

► HOCH₂(CHOH)₄CH₂OH.H₂O **M** = 200.19 a/mol

M = 200.19 g/molCAS NO. 50-70-4EINECS 200-061-5NC CODE 2905 44 91


98.0-101.0% Assay Appearance of solution passes test Chloride (CI) max. 50 ppm Identification passes test Lead (Pb) max. 0.5 ppm Nickel (Ni) max. 1 ppm Reducing Sugars (as Glucose) passes test Residue after Ignition (as SO₄) max. 0.1% Specific rotation +4.0 - +7.0Sulfate (SO₄) max. 100 ppm Water (H₂O) max. 1.5%

PRODUCT	PACKING	CONT
NO.		BOX
2039.1000	1 kg	

Sørenson buffer HEMATOLOGY/HISTOLOGY

Phosphate Buffer Concentrate (20x). For use in

combination with staining solutions



Speedisk

See for detailed information section Chromatography products, page 448



3716

Spill Cleanup Products

See for detailed information section Safety, page 400



4450

Spill Kit Accessories: Cleanup Centre

•

NC CODE 3822 00 00 PRODUCT PACKING CONT.

NO. BOX
4450 1 unit

Contains: Acid, Caustic and Solvent Spill Cleanup Kit plus a storage cabinet.

Spill Kit Accessories: Safety equipment kit

4509

 NC CODE
 3824 90 95
 PRODUCT
 PACKING
 CONT.

 NO.
 BOX

 4509
 1 unit

Spill Kit Accessories: Storage Cabinet for 3 kits

4435

NC CODE 7326 90 97

PRODUCT PACKING CONT.

NO.

4435 1 unit

Spill Kit Acids

4442

NC CODE 3822 00 00 PRODUCT PACKING CONT.

NO. BOX
4442 1 unit

Contains: 3.2 kg Neutrasorb, goggles, gloves, scoop, brush, waste bags, sponge.



Spill

	5 P				
4510	Spill Kit Acids Neutraso	rb			
NC CODE	3822 00 00		PRODUCT NO.	PACKING	CONT. BOX
			4510	3.2 kg	
4456	Spill Kit Acids Neutraso	rb			
NC CODE	3822 00 00		PRODUCT	PACKING	CONT. BOX
			4456.9045	45 kg	
1555-02	Spill Kit Acids TEAM Lo	w Na⁺			
NC CODE	3822 00 00	Especially recommended for use in clean rooms	PRODUCT NO.	PACKING	CONT. BOX
			4555-02	1 unit	ВОХ
			Contains: 6	x 946 ml; 1x handspraye	r.
4441	Spill Kit Caustics				
NC CODE	3822 00 00		PRODUCT	PACKING	CONT. BOX
			4441	1 unit	
				2 kg Neutrasorb, goggle e bags, sponge.	s, gloves, scoop,
4460	Spill Kit Caustics Neutra	cit			
NC CODE	3822 00 00		PRODUCT NO.	PACKING	CONT. BOX
			4460	23 kg	
	Spill Kit Caustics Neutra	ocit			

PRODUCT	PACKING	CONT.
NO.		вох
4512	1.2 kg	

Spill Kit Mercury

	Spill Kit Mercury Resisorb
4455	For Mercury vapor

CONT.	PACKING	PRODUCT
BOX		NO.
	900 g	4455.0900E
	11.3 kg	4455.7025E



Spill Kit Mercury: Mercury Sponge

J.T.Baker

						4444	
NC CODE	7323 10 00			PRODUCT NO.	PACKING	CON1 BO	
				4444	1 unit		
				Contains: ac	ctivator for the absorptic cury.	on of up to 10 g of	
	Spill Kit Solvents					4437	
NC CODE	2022 00 00			PROPUET	DACKING		
NC CODE	3822 00 00			PRODUCT NO.	PACKING	CONT BO	
				4437	1 unit	ВО	^
					1 kg Solusorb; goggles e bags; sponge.	; gloves; scoop;	
	Spill Kit Solvents Sol	usorb					
						4458	
NC CODE	3822 00 00			PRODUCT	PACKING	CONT	
				NO.	40 kg	ВО	X
				4458.7040	18 kg		
	Spill Kit Solvents Sol	usorb					
	•					4511	
NC CODE	3822 00 00			PRODUCT	PACKING	CONT	T.
				NO.		ВО	X
				4511	1.1 kg		
	a. 1 1 6 7 .						
	Standards for Trace A See for detailed information Stan						
	Stannous Chloride Di	hydrate					
	See Tin(II) Chloride Dihydrate						
	•						
	Starch 'BAKER ANALYZED' / Soluble p	otato for indometry / ACS				0327	
CAC NO	9005-25-8			PRODUCT	DACKING		_
	232-679-6	Meets ACS Specifications pH of 2% Solution at 25°C	5.0-7.0	NO.	PACKING	CONT BO	
	3505 10 90	Residue after Ignition	max. 0.4%	0327.0250	250 g		6
		Sensitivity	passes test	0327.2500	2.5 kg		Ē
		Solubility	passes test				
	Starch (Lintner) soluble / 'BAKER ANALYZED'					1130	
CAS NO.	9005-25-8	Apparent Viscosity, 1 hour	max. 2.0	PRODUCT	PACKING	CONT	г.
	232-686-4	Apparent Viscosity, 24 hours	max. 4.0	NO.		ВО	
	3505 10 90	Dextrin	passes test	1130.0500	500 g		
		Loss on Drying	max. 12%				
		pH of 2% Solution at 25°C	4.5-7.0				
		Dad	O 7E0/				

max. 0.75%

max. 0.5%

passes test

Reducing Substances (as maltose)

Residue after Ignition

Solubility



Stear

Stearic Acid

0329

► CH₃(CH₂)₁₆COOH M = 284.49 g/mol

> **CAS NO.** 57-11-4 **EINECS** 200-313-4 NC CODE 2915 70 25

Congealing Temperature	min. 54°C
Heavy Metals (as Pb)	max. 0.001%
lodine Value	max. 4
Mineral Acids	passes test
Neutral fat or Paraffin	passes test
Residue on Ignition	max. 0.1%

		CONT
NO.		BO)
0329.1000	1 kg	



Stearic Acid

See for qualities meeting the requirements of cGMP, USP/NF, BP/Ph.Eur, JP or FCC, page 36



Stearic Acid, 92%

See for qualities meeting the requirements of cGMP, USP/NF, BP/Ph.Eur, JP or FCC, page 36

Strontium 1000 μ g/ml

(Matrix: 2% nitric acid) / 'BAKER INSTRA-ANALYZED' / Plasma Standard 5781

► Sr

M = 87.62 g/mol NC CODE 3822 00 00

R: 36/38

S: 26



Certificate Provided Reporting Actual Lot Analysis Strontium (Sr) 998-1002 μg/ml

PRODUCT PACKING CONT. NO. вох 5781.0100 100 ml

Prepared from the highest purity raw material available, generally greater than 99.999% spectral purity. The content of the solution is confirmed to be accurate to within \pm 0.2%. Typically 1000 μ g/ml. The certificate of analysis provided reports actual lot analysis. The certificate also lists the trace impurities.

 $\begin{array}{lll} \textbf{Strontium 1000 } \ \mu \textbf{g/ml} \\ \textbf{(Matrix: } 2\% \ nitric \ acid) \ / \ \ \textbf{'BAKER INSTRA-ANALYZED'} \ / \ \ \textbf{Atomic Absorption Standard} \\ \end{array}$ 6942

► Sr

M = 87.62 g/mol

NC CODE 3822 00 00

R: 36/38 **S**: 26

Strontium (Sr)	998-1002 μ g/m

6942.0100 100 ml

PACKING

PRODUCT

PRODUCT

NO. BOX 6942.0500 500 ml

CONT.

CONT.

CONT.

BOX

BOX

Prepared by dissolution of high purity raw materials (min. 99.99% spectral purity). Assays are verified by ICP against standards traceable to NIST. Standard Reference Material numbers (SRM) are printed on each label.

Strontium 1000 µg/ml

6823 'BAKER ANALYZED' / Atomic Absorption Standard

► Sr

M = 87.62 g/mol

NC CODE 3822 00 00

R: 36/38 **S:** 26-37

Strontium (Sr)

998-1002 µg/ml

NO. 6823.0100 100 ml 6823.0500 500 ml

Strontium chloride in nitric acid 0.5 mol/l.

PACKING

Strontium 10000 µg/ml

(Matrix: 2% nitric acid) / 'BAKER INSTRA-ANALYZED' / Plasma Standard

► Sr

5747

M = 87.62 g/mol NC CODE 3822 00 00

R: 36/38

S: 26-37



Certificate Provided Reporting Actual Lot Analysis Strontium (Sr) 9980-10020 $\mu g/ml$

PRODUCT PACKING NO. 5747.0100 100 ml

Prepared from the highest purity raw material available, generally greater than 99.999% spectral purity. The content of the solution is confirmed to be accurate to within 0.2 %. Typically 10000 μ g/ml.

The certificate of analysis provided reports actual lot analysis. The certificate also lists the trace impurities.



Strontium Chloride Hexahydrate

'BAKER ANALYZED' / ACS 0331

➤ SrCl₂.6H₂O **M** = 266.62 g/mol **CAS NO.** 10025-70-4 **EINECS** 233-971-6 NC CODE 2827 39 90

R: 20/21/22

S: 36-39

Meets ACS Specifications Assay 99.0-103.0% Barium (Ba) max. 0.05% Calcium (Ca) max. 0.05% Insoluble Matter max. 0.005% pH of 5% Solution at 25°C 5.0-7.0 Sulfate (SO₄) max. 0.001% Trace Impurities (in ppm): Heavy Metals (as Pb) max. 5 Iron (Fe) max. 5 Magnesium (Mg) max. 2

PRODUCT	PACKING	CONT.
NO.		вох
0331.0500	500 g	
0331.9050	50 kg	

Styrene

8133

max. 10

max. 5

► C₆H₅CH:CH₂

M = 104.15 g/mol **1 l** = 0.91 kg

FLASHPOINT 32 °C

CAS NO. 100-42-5 EINECS 202-851-5

NC CODE 2902 50 00

EC NO. 601 026 00 0 **UN/ID NO.** 2055

ADR/RID 3 F1 IMDG 3/III

R: 10-20-36/38

S: 23

Density (g/ml) at 20°C	0.905-0.907

Color (APHA)

PRODUCT PACKING CONT. NO. BOX 8133.0100 100 ml 8133.1000 11

Stabilized with tert-Butylpyrocatechol.

Succinic Acid

'BAKER ANALYZED' / ACS

M = 118.09 g/mol **CAS NO.** 110-15-6

► (CH₂COOH)₂

EINECS 203-740-4 NC CODE 2917 19 90

> **R:** 36 **S**: 26



Meets ACS Specifications	
Assay	min. 99.0%
Chloride (CI)	max. 0.001%
Insoluble Matter	max. 0.01%
Melting Point	185.0-191.0°C
Nitrogen Compounds (as N)	max. 0.001%
Phosphate (PO ₄)	max. 0.001%
Residue after Ignition	max. 0.02%
Sulfate (SO ₄)	max. 0.003%
Trace Impurities (in ppm):	
Heavy Metals (as Pb)	max. 5

PRODUCT NO.	PACKING	CONT BO
0333.0250	250 g	

Succinic Acid Disodium Salt



0333

MSDS (Material Safety Data Sheets) are available in 16 languages at www.jtbaker.com/europe

Iron (Fe)

Sucro

Sucrose

ULTREX Ultrapure Reagent

► C₁₂H₂₂O₁₁

M = 342.30 g/mol **CAS NO.** 57-50-1 **EINECS** 200-334-9 NC CODE 1701 99 10

Certificate Provided Reporting Actual Lot Analysis Actual Analysis Lot.No. L37479

Ash (sulfated) Chloride (CI) 0.003 Homogeneity (by GLC) passes test Homogenity (by TLC) no extraneous spots Invert Sugar < 0.05% Iron (Fe) < 1 Loss on Drying at 105°C 0.02% Particulate Matter 0.003% Specific Rotation $[\alpha]^{25}_{\overline{D}}$ 66.4° Sulfate and Sulfite (as SO₄) 0.005 Titrable Acid (µeq/g)

Metallic Impurities in parts per million ($\mu g/g$):

Aluminium (AI) Bismuth (Bi) < 1 Cadmium (Cd) Calcium (Ca) < 1 Chromium (Cr) < 1 Cobalt (Co) < 10 Copper (Cu) Lead (Pb) < 10 Magnesium (Mg) < 0.5 < 10 Manganese (Mn) Mercury (Hg) < 0.0001 Molybdenum (Mo) < 1 Nickel (Ni) Silver (Ag) Sodium (Na) Tin (Sn) < 10 Titanium (Ti) Vanadium (V) < 1 Zinc (Zn) Non-Metallic Impurities in parts per million ($\mu g/g$):

PRODUCT NO.	PACKING	CONT. BOX
4890.0100	100 g Glass	
4890.0500	500 g Glass	

Sucrose

'BAKER ULTRAPURE BIOREAGENT'

► C₁₂H₂₂O₁₁

M = 342.30 g/mol **CAS NO.** 57-50-1 **EINECS** 200-334-9 NC CODE 1701 99 10

For Density Gradient Centrifugation Assay (by HPLC) min. 99.9% Clarity and Color of Solution passes test Clarity of Solution passes test DNase Activity none detected max. 0.1% Protease Activity none detected RNase Activity none detected

< 0.5

< 2

Trace Impurities (in ppm):

Arsenic (As) Boron (B)

Silicon (Si)

Nitrogen Compounds (as N)

Total of Copper (Cu), Iron (Fe) and Lead

PACKING	CONT.	
	ВОХ	
1 kg		
5 kg		
	1 kg	

Find our up-to-date Product Literature at www.jtbaker.com/europe



Sucrose

0334 'BAKER ANALYZED' / ACS

► C₁₂H₂₂O₁₁

M = 342.30 g/mol **CAS NO.** 57-50-1 **EINECS** 200-334-9 NC CODE 1701 99 10

Meets ACS Specifications Chloride (CI) max. 0.005% Insoluble Matter max. 0.005% Invert Sugar max. 0.05% Iron (Fe) max. 5 ppm Loss on Drying at 105°C max. 0.03% Residue after Ignition max. 0.01% Specific Rotation [α] $\frac{26}{5}$ +66.3°-+66.8° Sulfate and Sulfite (as SO₄) max. 0.005% Titrable Acid (meq/g) max. 0.0008 Trace Impurities (in ppm): Heavy Metals (as Pb) max. 5

PRODUCT	PACKING	CONT.
NO.		вох
0334.1000	1 kg	6
0334.5000	5 kg	
0334.9050	50 kg	

Sucrose

See for qualities meeting the requirements of cGMP, USP/NF, BP/Ph.Eur, JP or FCC, page 36

Iron (Fe)



0707

1828

Sudan III

► C₆H₂N:NC₆H₄N:NC₁₀H₆OH **M =** 352.40 g/mol

> **CAS NO.** 85-86-9 **EINECS** 201-638-4 NC CODE 2927 00 00

PRODUCT PACKING CONT. NO. BOX 0707.0025 25 g Glass

C.I. 26100.

Sulfamic Acid

'BAKER ANALYZED

M = 97.09 g/mol **CAS NO.** 5329-14-6 **EINECS** 226-218-8

► NH₂SO₃H

NC CODE 2811 1980 **EC NO.** 16 026 00 0 **UN/ID NO.** 2967

ADR/RID 8 C2 IMDG 8/III

R: 36/38-52/53 **S:** 26-28-61



min. 99% Assav Chloride (CI) max. 0.001% Heavy Metals (as Pb) max. 0.001% Insoluble Matter max. 0.005% Residue after Ignition max. 0.01% Sulfate (SO₄) max. 0.02%

Trace Impurities (in ppm): max. 5

PRODUCT PACKING CONT. 1828.0250 250 g 1828,1000 1828.5000 5 kg

Sulfamic Acid, Ammonium Salt



0709

Sulfanilamide

4-NH₂C₆H₄SO₂NH₂ **M** = 172.21 g/mol

CAS NO. 63-74-1 **EINECS** 200-563-4 NC CODE 2935 00 90 Melting Point

164-166°C

PRODUCT PACKING NO. 0709.1000 1 kg

CONT. BOX

www.jtbaker.com/europe



Sulfa

Sulfanilic Acid Anhydrous

1197 BAKER ANALYZED' / ACS

► 4-NH₂C₆H₄SO₃H **M** = 173.19 g/mol **CAS NO.** 121-57-3 **EINECS** 204-482-5 NC CODE 2921 42 10 **EC NO.** 612 014 00 0 **R:** 36/38-43 **S:** 24-37

Exceeds ACS Specifications 98.0-102.09 Assay max. 0.0029 Chloride (CI) Insoluble in Sodium Carbonate Solution max. 0.019 Residue after Ignition max. 0.019 Suitability for Nitrogen Dioxide Determination(ASTM D1607) passes tes Sulfate (SO₄) max. 0.005% Trace Impurities (in ppm): max. 0.5 Nitrite (NO₂)

PRODUCT	PACKING	CONT
NO.		B0)
1197.0100	100 g	

Sulfosalicylic Acid Dihydrate

► 2-HOC₆H₃-1-COOH-5-SO₃H.2H₂O **M** = 254.22 g/mol **CAS NO.** 5965-83-3 **EINECS** 202-555-6 NC CODE 2918 29 10 **UN/ID NO.** 3261 ADR/RID 8 C4 IMDG 8/III R: 36/38 **S:** 26

Meets ACS Specifications Assay (acidimetric) 99.0-101.0% Chloride (CI) max. 0.001% Heavy Metals (as Pb) max. 0.002% Insoluble Matter max. 0.02% max. 0.001% Iron (Fe) Residue after Ignition max. 0.1% Salicylic Acid (HOC₆H₄COOH) max. 0.04% Sulfate (SO₄) max. 0.02%

PRODUCT NO.	PACKING	CONT.	
		вох	
1198.0100	100 g		
1198.1000	1 kg		

Sulfur

0335 Precipitated / 'BAKER'

S M = 32.06 g/mol **CAS NO.** 7704-34-9 **EINECS** 231-722-6 NC CODE 2802 00 00 UN/ID NO. 1350 **ADR/RID** 4.1 F3

IMDG 4.1/III

passes test Appearance Assay (anhydrous basis) 99.5 - 100.5% Identification passes test Other forms of sulfur passes test Reaction passes test Residue on Ignition max. 0.3% Water (H₂O) max. 0.5%

PRODUCT NO.	PACKING	CONT.	
		вох	
0335.1000	1 kg	6	



See for qualities meeting the requirements of cGMP, USP/NF, BP/Ph.Eur, JP or FCC, page 36

Use J.T.Baker Ultrex II and BAKER INSTRA-ANALYZED acids for low level trace element analysis.

See chapter 3 of this catalogue for more details.



Sulfuric Acid

95% / ULTREX II Ultrapure Reagent

6902

► H ₂ SO ₄	
M =	98.08 g/mol
1 =	1.84 kg
CAS NO.	7664-93-9
EINECS	231-639-5
NC CODE	2807 00 10
EC NO.	16 020 00 8
UN/ID NO.	1830
ADR/RID	8 C1
IMDG	8/II
R:	35
S:	26-30-45
	C S

Certificate Provided Reporting Ac	tual Lot Analysis
Assay	93 -98% (w/w)
Trace Impurities (in ppt) (pg/g):	
Aluminium (AI)	max. 50
Antimony (Sb)	max. 50
Arsenic (As)	max. 500
Barium (Ba)	max. 10
Beryllium (Be)	max. 10
Bismuth (Bi)	max. 10
Cadmium (Cd)	max. 10
Calcium (Ca)	max. 100
Cerium (Ce)	max. 10
Cesium (Cs)	max. 10
Chromium (Cr)	max. 10
Cobalt (Co)	max. 10
Copper (Cu)	max. 10
Dysprosium (Dy)	max. 10
Erbium (Er)	max. 10
Europium (Eu)	max. 10
Gadolinium (Gd)	max. 10
Gallium (Ga)	max. 10
Germanium (Ge)	max. 100
Hafnium (Hf)	max. 10
Holmium (Ho)	max. 10
Indium (In)	max. 10
Iron (Fe)	max. 50
Lanthanum (La)	max. 10
Lead (Pb)	max. 10
Lithium (Li)	max. 10
Lutetium (Lu)	max. 10
Magnesium (Mg)	max. 50
Manganese (Mn)	max. 10
Mercury (Hg)	max. 100
Molybdenum (Mo)	max. 10
Neodymium (Nd)	max. 10
Nickel (Ni)	max. 50
Niobium (Nb)	max. 10
Palladium (Pd)	act. value reported
Platinum (Pt)	act. value reported
Potassium (K)	max. 50
Praseodymium (Pr)	max. 10

Rhodium (Rh)	max. 50
Rubidium (Rb)	max. 10
Samarium (Sm)	max. 10
Scandium (Sc)	max. 10
Selenium (Se)	max. 500
Silver (Ag)	max. 50
Sodium (Na)	max. 50
Strontium (Sr)	max. 10
Tantalum (Ta)	act. value reported
Tellurium (Te)	max. 100
Terbium (Tb)	max. 10
Thallium (TI)	max. 10
Thorium (Th)	max. 10
Thulium (Tm)	max. 10
Tin (Sn)	max. 50
Titanium (Ti)	max. 50
Tungsten (W)	max. 10
Uranium (U)	max. 10
Vanadium (V)	max. 10
Ytterbium (Yb)	max. 10
Yttrium (Y)	max. 10
Zinc (Zn)	max. 50
Zirconium (Zr)	max. 10
PRODUCT PACKING	CONT.

PRODUCI	PACKING	CONT.
NO.		вох
6902.0500	500 ml Fluoropolymer, pre-leached	
6902.1000	1 Fluoropolymer, pre-leached	

Calibrate and standardise your analytical methods and equipment with J.T.Baker Volumetric and Buffer solutions.

Refer to the Analytical applications section of this catalogue for more details.



Sulfu

Sulfuric Acid

95-98% / 'BAKER INSTRA-ANALYZED' / for Trace Metal Analysis

M = 98.08 g/mol

1 I = 1.84 kg

CAS NO. 7664-93-9

EINECS 231-639-5

NC CODE 2807 00 10

EC NO. 16 020 00 8

UN/ID NO. 1830

ADR/RID 8 C1

IMDG 8/II

R: 35

S: 26-30-45

20-30-	40
I.e.	С
corrosiv	е

Certificate Provided Reports Actual Lot Analysis		
(About 36.0 N)		
Assay (by acid-base titration)	95.0-98.0%	
Color (APHA)	max. 10	
Trace Impurities (in ppm):		
Aluminium (AI)	max. 0.03	
Ammonium (NH ₄)	max. 1	
Arsenic and Antimony (as As)	max. 0.004	
Barium (Ba)	max. 0.02	
Beryllium (Be)	max. 0.01	
Bismuth (Bi)	max. 0.05	
Boron (B)	max. 0.01	
Cadmium (Cd)	max. 0.002	
Calcium (Ca)	max. 0.05	
Chloride (CI)	max. 0.1	
Chromium (Cr)	max. 0.006	
Cobalt (Co)	max. 0.0005	
Copper (Cu)	max. 0.001	
Gallium (Ga)	max. 0.02	
Germanium (Ge)	max. 0.10	
Gold (Au)	max. 0.05	
Heavy Metals (as Pb)	max. 0.5	
Iron (Fe)	max. 0.05	
Lead (Pb)	max. 0.0005	
Lithium (Li)	max. 0.02	
Magnesium (Mg)	max. 0.007	
Manganese (Mn)	max. 0.001	
Mercury (Hg)	max. 0.0005	
Molybdenum (Mo)	max. 0.05	
Nickel (Ni)	max. 0.002	
Niobium (Nb)	max. 0.05	
Nitrate (NO ₃)	max. 0.2	
Phosphate (PO ₄)	max. 0.5	
Potassium (K)	max. 0.5	
Residue after Ignition	max. 3	
Selenium (Se)	value on label	
Silicon (Si)	max. 0.1	
Silver (Ag)	max. 0.001	
Sodium (Na)	max. 0.5	
Strontium (Sr)	0.005	
	0.000	

DDODUCT DACKING	CONT
Zirconium (Zr)	max. 0.01
Zinc (Zn)	max. 0.005
Vanadium (V)	max. 0.01
Titanium (Ti)	max. 0.01
Tin (Sn)	0.005
Thallium (TI)	max. 0.08

PRODUCT NO.	PACKING	CONT. BOX
6163.2500	2.5	4

The J.T.Baker CYCLE-TAINER
High Purity Solvent Delivery System,
preserves purity and protects people.

Tantalum (Ta)

See chapter 3 of this catalogue for product details.



Sulfuric Acid

95-98% (max. 5 ppb Hg) / 'BAKER ANALYZED' / ACS

6027

► H ₂ SO ₄	
M =	98.08 g/mol
1 =	1.84 kg
CAS NO.	7664-93-9
EINECS	231-639-5
NC CODE	2807 00 10
EC NO.	16 020 00 8
UN/ID NO.	1830
ADR/RID	8 C1
IMDG	8/11
R:	35
S:	26-30-45
	C corrosive

Assay	95.0-98.0%
Ammonium (NH ₄)	max. 1 ppm
Appearance	passes test
Color (APHA)	max. 10
Mercury (Hg)	max. 5 ppb
Residue after Ignition	max. 4 ppm
Substances Reducing KMnO ₄ (as SO ₂)	max. 2 ppm
Trace Impurities (in ppm):	
Aluminium (AI)	max. 0.05
Arsenic (As)	max. 0.01
Barium (Ba)	max. 0.02
Beryllium (Be)	max. 0.01
Cadmium (Cd)	max. 0.02
Calcium (Ca)	max. 0.5
Chloride (CI)	max. 0.1
Chromium (Cr)	max. 0.1
Cobalt (Co)	max. 0.01
Copper (Cu)	max. 0.01
Germanium (Ge)	max. 0.05
Heavy Metals (as Pb)	max. 1
Iron (Fe)	max. 0.1
Lead (Pb)	max. 0.05
Lithium (Li)	max. 0.01
Magnesium (Mg)	max. 0.1
Manganese (Mn)	max. 0.01
Molybdenum (Mo)	max. 0.02
Nickel (Ni)	max. 0.05
Nitrate (NO ₃)	max. 0.5
Potassium (K)	max. 0.1
Silver (Ag)	max. 0.01
Sodium (Na)	max. 0.5
Strontium (Sr)	max. 0.01
Thallium (TI)	max. 0.05
Titanium (Ti)	max. 0.1
Vanadium (V)	max. 0.01
Zinc (Zn)	max. 0.05

Zirconium (Zr)

PRODUCT	PACKING	CONT.
NO.		вох
6027.1000	11	
6027.2500	2.5	4

Mallinckrodt Baker's cGMP Manufactured Chemicals for the Biopharmaceutical industry are a necessity for uncomplicated scale-up.

max. 0.1

See chapter 6 of this catalogue.

R S



Sulfu

Sulfuric Acid

95-97% / 'BAKER ANALYZED' 6057

► H₂SO₄ **M** = 98.08 g/mol **1 l** = 1.84 kg **CAS NO.** 7664-93-9 **EINECS** 231-639-5 NC CODE 2807 00 10 **EC NO.** 16 020 00 8 **UN/ID NO.** 1830 ADR/RID 8 C1 IMDG 8/II **R:** 35 **S:** 26-30-45



95-97% Assay Ammonium (NH₄) max. 1 ppm Color (APHA)

Color (APHA)	max. 10
Heavy Metals (as Pb)	max. 5 ppm
Nitrate (NO ₃)	max. 1 ppm
Residue after Ignition	max. 5 ppm
Substances Reducing KMnO ₄ (as SO ₂)	max. 5 ppm
Trace Impurities (in ppm):	
Aluminium (AI)	max. 0.05
Arsenic (As)	max. 0.05
Barium (Ba)	max. 0.02
Beryllium (Be)	max. 0.01
Cadmium (Cd)	max. 0.02
Calcium (Ca)	max. 0.5
Chloride (CI)	max. 0.5
Chromium (Cr)	max. 0.1
Cobalt (Co)	max. 0.01
Copper (Cu)	max. 0.01
Germanium (Ge)	max. 0.05
ron (Fe)	max. 0.1
Lead (Pb)	max. 0.05
Lithium (Li)	max. 0.01
Magnesium (Mg)	max. 0.1
Manganese (Mn)	max. 0.01
Molybdenum (Mo)	max. 0.02
Nickel (Ni)	max. 0.05
Potassium (K)	max. 0.1
Silver (Ag)	max. 0.01
Sodium (Na)	max. 0.5
Strontium (Sr)	max. 0.01
Thallium (TI)	max. 0.05
Titanium (Ti)	max. 0.1

PRODUCT	PACKING	CONT.
NO.		вох
6057.1000	11	6
6057.2500	2.5	4
6057.2500PE	2.5 HDPE	4
6057.9025	25 l	
6057.9200	200	

CONT.

вох



Sulfuric Acid 96% CMOS, Finyte, Finyte-1, Ultryte Grade

Vanadium (V) Zinc (Zn)

Zirconium (Zr)

See for detailed information section Microelectronic materials page 32 or the website www.jtbaker.com/micro

Sulfuric Acid

6147 90-91% $\,/\,$ 'BAKER ANALYZED' $\,/\,$ For fat determination according to Gerber

► H₂SO₄ **M** = 98.08 g/mol **1 l** = 1.82 kg **CAS NO.** 7664-93-9 **EINECS** 231-639-5 NC CODE 2807 00 10 **EC NO.** 16 020 00 8

UN/ID NO. 1830 ADR/RID 8 C1 IMDG 8/II

R: 35

S: 26-30-45



Assay (acidimetric)	90-91%	PRODUCT	PACKING
Density (g/ml)	1.815-1.821	NO.	
Nitrate (NO ₃)	max. 0.2 ppm	6147.2500	2.5
Suitability for Gerber test	passes test	6147.9025	25 l


max. 0.01

max. 0.05

max. 0.1

7242

Sulfuric Acid
5 mol/l / 'BAKER ANALYZED'



Sulfuric Acid

7610 4 mol/l / 'BAKER ANALYZED'

► H ₂ SO ₄		Titer (mol/l)	3.95-4.05	PRODUCT	PACKING	CONT.
M =	98.08 g/mol			NO.		BOX
CAS NO.	7664-93-9			7610.1000	11	6
EINECS	231-639-5					
NC CODE	2807 00 10			Volumetric Sc	olution, ready for use.	
EC NO.	16 020 00 8					
UN/ID NO.	2796					
ADR/RID	8 C1					
IMDG	8/II					
R:	35					
S:	26-36/37/39-45					
	C corrosive					

Sulfuric Acid 2.6 mol/l / 'BAKER ANALYZED'

► H ₂ SO ₄	Titer (mol/l)	2.55-2.65	PRODUCT	PACKING	CONT.
M = 98.08 g/mol			NO.		вох
EINECS 231-639-5			7471.9010	10 l Polycube	
NC CODE 2807 00 10			7471.9025	25 l	
UN/ID NO. 2796			7471.9200	200 l	
ADR/RID 8 C1					
IMDG 8/II			Volumetric S	Solution, ready for use.	
R: 35			Each lot of the	his product is standardized	d potentiometrically
S: 26-30-36/37/39-45			against NIST	tracable reference standa	ard.
Corrosive					

Sulfuric Acid

7624 2 mol/l / 'BAKER ANALYZED'

► H ₂ SO ₄		Titer (mol/l)	1.99 - 2.01	PRODUCT	PACKING	CONT.
CAS NO.	7664-93-9			NO.		вох
EINECS	231-639-5			7624.5000	51	
NC CODE	2807 00 10			7624.9010	10 l Polycube	4
UN/ID NO.	2796					
ADR/RID	8 C1			Volumetric Sc	olution, ready for use.	
IMDG	8/II					
R:	35					
S:	26-30-36/37/39-45					
	C corrosive					



J.T.Baker Spill Clean-up Products

Spills and hazardous waste, such as acids and flammable solvents, are becoming a bigger concern to all users of chemicals. Those are also the most hazardous to handle and spills of these products can lead to fire, explosion, burns and exposure to toxic substances. It is self-evident to all users that it is necessary to be prepared for accidental chemical spills.

The J.T.Baker SAF-T-Spill Kits do not just absorb; they also neutralise the hazard or reduce the flammable or toxic vapours. We have added in other features for your safety and convenience.

Besides the SAF-T-Spill Kits for laboratory use, products for larger spills are also available.

Built-in excellence	Making life simpler
HDPE bottle	Easy-grip and easy-pour
Colour coded	Maximum safety
Tamper-evident carton	Guarantees unused packaging
Step-by-step instructions	Trouble-free working

Safety products we offer

Spill clean-up products for laboratory use

Code	Description
4450	Spill Cleanup Centre (contains an Acid, Caustic and Solvent Spill Cleanup Kit plus a storage cabinet)
4442	Acid Spill Cleanup Kit (contains 3.2 kg Neutrasorb)
4441	Caustic Spill Cleanup Kit (contains 1.2 kg Neutracit)
4437	Solvent Spill Cleanup Kit (contains 1.1 kg Solusorb)
4439	Mercury Spill Cleanup Kit (contains 0.9 kg Resisorb, 250 g Cinnasorb Base, 20 g Cinnasorb Activator)
4435	Storage Cabinet for 3 kits

Speciality safety products

Code	Description
4444	Mercury Sponge with activator for the absorption of up to 10 g of spilled mercury

Spill clean-up products for larger spills

Code	Description
4510	1 box Neutrasorb of 3.2 kg for Acids
4512	1 box Neutracit of 1.2 kg for Caustics
4511	1 box Solusorb of 1.1 kg for Solvents
4455.0900	1 box Resisorb of 900 g for Mercury Vapour
4455.7025	1 box Resisorb of 11.3 kg for Mercury Vapour
4509	Safety equipment kit

Spill clean-up products in bulk containers for pilot plants or warehouse

Code	Description
4456.9045	45 kg Neutrasorb for Acids, in fibre drum with scoop
4460	23 kg Neutracit for Caustics, in fibre drum with scoop
4458.7040	18 kg Solusorb for Solvents, in fibre drum with scoop

Cleaning capacity Spill clean-up products

Veutrasorb	for	Acids	32	ka
VCULIUSCIS	101	, was,	0.2	119

Neutrasorb for Acids, 5.2 kg	
Acetic Acid (98%)	0.47 L
Hydrochloric Acid (38%)	1.94 L
Nitric Acid (70%)	1.51 L
Phosphoric Acid (87%)	0.47 L
Sulphuric Acid (98%)	0.66 L
Neutracit for Caustics, 1.2 kg	
Ammonium Hydroxide (56%)	0.56 L
Potassium Hydroxide (45%)	0.75 L
Sodium Hydroxide (50%)	0.47 L
Solusorb for Solvents, 1.1 kg	
Flammable solvents (Toluene, Ether etc.)	0.38 L
Toxic and Noxious liquids (Chloroform, Pyridine, THF, etc)	0.38 L
Halogens-Corrosive hazards (Bromine, Iodine, etc.)	0.45 kg



Sulfu

Sulfuric Acid

7102 0.5 mol/l 1N / 'BAKER ANALYZED'

M = 98.08 g/mol
CAS NO. 7664-93-9
EINECS 231-639-5
NC CODE 2807 00 10

Titer (N)	0.997-1.003
Trace Impurities (in ppm):	
Chloride (CI)	max. 1
Heavy Metals (as Pb)	max. 0.3
Iron (Fe)	max. 0.5
Phosphate (PO ₄)	max. 1

PRODUCT	PACKING	CONT.
NO.		вох
7102.1000	11	6
7102.9010	10 l Polycube	
7102.9020	20 l Polycube	

Volumetric Solution, ready for use.

Each lot of this product is standardized potentiometrically against Sodium Carbonate (NIST tracable reference standard).

Sulfuric Acid

7110 0.25 mol/I 0.5N / 'BAKER ANALYZED'

M = 98.08 g/mol
CAS NO. 7664-93-9
EINECS 231-639-5
NC CODE 2807 00 10

Titer (N)	0.4975-0.5025
Trace Impurities (in ppm):	
Chloride (CI)	max. 1
Heavy Metals (as Pb)	max. 0.3
Iron (Fe)	max. 0.5
Phosphate (PO ₄)	max. 1

PRODUCT NO.	PACKING	CONT. BOX
7110.1000	11	6
7110.9010	10 l Polycube	

Volumetric Solution, ready for use.

Each lot of this product is standardized potentiometrically against Sodium Carbonate (NIST tracable reference standard).

Sulfuric Acid

7638 0.13 mol/l / 'BAKER ANALYZED'

► H₂SO₄ 1 I = 1.00 kg CAS NO. 7664-93-9 EINECS 231-639-5 NC CODE 2807 00 10
 0.129 - 0.131
 PRODUCT
 PACKING
 CONT.

 NO.
 BOX

 7638.9010
 10 l Polycube

Volumetric Solution, ready for use.

7215

Sulfuric Acid 0.1 mol/10.2 N / 'BAKER ANALYZED'

M = 98.08 g/mol
CAS NO. 7664-93-9
EINECS 231-639-5
NC CODE 2807 00 10

Titer (N) 0.199-0.201

Titer (mol/l)

 PRODUCT
 PACKING
 CONT.

 NO.
 BOX

 7215.9010
 10 | Polycube
 4

Volumetric Solution, ready for use.

Each lot of this product is standardized potentiometrically against Sodium Carbonate (NIST tracable reference standard).

7102

Sulfuric Acid

0.05 mol/l 0.1N / 'BAKER ANALYZED'

M = 98.08 g/mol

CAS NO. 7664-93-9

EINECS 231-639-5

NC CODE 2807 00 10

Titer (N)	0.0997-0.1003
Trace Impurities (in ppm):	
Chloride (CI)	max. 1
Heavy Metals (as Pb)	max. 0.3
Iron (Fe)	max. 0.5
Phosphate (PO ₄)	max. 1

PRODUCT	PACKING	CONT.
NO.		вох
7103.1000	11	6
7103.9010	10 l Polycube	
7103.9020	20 l Polycube	

Volumetric Solution, ready for use.

Each lot of this product is standardized potentiometrically against Sodium Carbonate (NIST tracable reference standard).

PRODUCT

PRODUCT

PRODUCT

NO.

4699

NO.

4700

NO.

4701

PACKING

Volumetric Concentrate, for dilution to 1 l.

PACKING

1 amp.

PACKING

1 amp.

Volumetric Concentrate, for dilution to 1 l.

500 ml Sealed Bottle

4701

CONT.

вох

Sulfuric Acid

2.5 mol/l / DILUT-IT

 H_2SO_4 **M** = 98.08 g/mol

CAS NO. 7664-93-9 **EINECS** 231-639-5

NC CODE 2807 00 10

EC NO. 16 020 00 8

UN/ID NO. 1830

ADR/RID 8 C1

IMDG 8/I

R: 35

S: 26-30-36/37/39-45



Sulfuric Acid

0.5 mol/l / DILUT-IT

► H₂SO₄

M = 98.08 g/mol

1 l = 1.29 kg

CAS NO. 7664-93-9

EINECS 231-639-5

NC CODE 2807 00 10

EC NO. 16 020 00 8

UN/ID NO. 2796 **ADR/RID** 8 C1

IMDG 8/II

R: 35

S: 26-30-36/37/39-45



Sulfuric Acid

0.05 mol/l / DILUT-IT

► H₂SO₄

M = 98.08 g/mol

 $\textbf{CAS NO.} \ \ 7664\text{-}93\text{-}9$

EINECS 231-639-5

NC CODE 2807 00 10

EC NO. 16 020 01 5

UN/ID NO. 2796

ADR/RID 8 C1

IMDG 8/II

R: 36/38

S: 26

irritant

Sulfuric Acid

0.01 mol/l / DILUT-IT

► H₂SO₄

M = 98.08 g/mol

CAS NO. 7664-93-9

EINECS 231-639-5 **NC CODE** 2807 00 10

EC NO. 16 020 01 5

PRODUCT PACKING

NO. 4704 1 amp.

Volumetric Concentrate, for dilution to 1 l.

4699

CONT.

вох

вох

CONT.

4700

Volumetric Concentrate, for dilution to 1 l.

CONT.

BOX

Sulfu

Sulfuric Acid

4869 0.005 mol/I 0.01N / DILUT-IT

NC CODE 2807 00 10

 PRODUCT
 PACKING
 CONT.

 NO.
 BOX

 4869
 1 amp.

Volumetric Concentrate, for dilution to 1 l.



Sulfuric Acid

See for qualities meeting the requirements of cGMP, USP/NF, BP/Ph.Eur, JP or FCC, page 36



Sulfuric Acid Solutions

See for qualities meeting the requirements of cGMP, USP/NF, BP/Ph.Eur, JP or FCC, page 36

Sulfurous Acid

6050 (solution of SO₂ in water) / 'BAKER ANALYZED' / ACS

► H₂SO₃

M = 82.08 g/mol

11= 1.03 kg

CAS NO. 7782-99-2 **NC CODE** 2811 19 80

UN/ID NO. 1833

ADR/RID 8 C1

IMDG 8/11

R: 20-34

S: 26-36/37/39-45



Meets ACS Specifications. Meets Reagent Specifications for testing USP/NF monographs

 $\begin{array}{ll} \mbox{Assay (as SO_2) (by iodometry)} & 6.0 - 9.0\% \\ \mbox{Residue after Ignition} & \mbox{max. } 0.005\% \\ \end{array}$

Trace Impurities (in ppm):

 Chloride (CI)
 max. 5

 Heavy Metals (as Pb)
 max. 2

 Iron (Fe)
 max. 5

PRODUCT	PACKING	CONT.
NO.		вох
	man Lol	

0339

Talc

Powder / 'BAKER'

 CAS NO.
 14807-96-6
 Loss on drying (at 180°C)
 max. 1.0%
 PRODUCT PACKING
 PACKING
 CONT.

 EINECS
 238-877-9
 Organic Substances
 passes test
 NO.
 BOX

 NC CODE
 2526 20 00
 0339.2500
 2.5 kg



Talc

See for qualities meeting the requirements of cGMP, USP/NF, BP/Ph.Eur, JP or FCC, page 36



Tannic Acid

See Tannin



Tannic Acid

See for qualities meeting the requirements of cGMP, USP/NF, BP/Ph.Eur, JP or FCC, page 36

Tannin

99 'BAKER ANALYZED'

► C₇₆H₅₂O₄₆

M = 1701.23 g/mol **CAS NO.** 1401-55-4

EINECS 215-753-2 **NC CODE** 3201 90 90

Loss on Drying at 100°C max. 12%
Residue after Ignition max. 0.1%

 PRODUCT
 PACKING

 NO.
 250 g

CONT.

вох

Innovation is principal to our business.

5782

► Ta

M = 180.95 g/mol

NC CODE 3822 00 00 **UN/ID NO.** 1790

ADR/RID 8 CT1 IMDG 8/II

R: 20/21/22-36

S: 26-36/37



Certificate Provided Reporting Actual Lot Analysis

Tantalum (Ta) 998-1002 μg/ml

PRODUCT	PACKING	CONT
NO.		ВОХ
5782.0100	100 ml	

Prepared from the highest purity raw material available, generally greater than 99.999% spectral purity. The content of the solution is confirmed to be accurate to within \pm 0.2%. Typically 1000 $\mu\text{g/ml.}$

The certificate of analysis provided reports actual lot analysis. The certificate also lists the trace impurities.

Tantalum 10000 μg/ml

'BAKER INSTRA-ANALYZED' / Plasma Standard (Matrix: 1% hydrofluoric acid) /

5748

0340

► Ta

M = 180.95 g/mol

NC CODE 3822 00 00 **UN/ID NO.** 1790 ADR/RID 8 CT1

IMDG 8/II

R: 23/24/25-34

S: 20-26-28-36/37/39-4-45-9





Certificate Provided Reporting Actual Lot Analysis Tantalum (Ta) 9980-10020 μg/ml

CONT. **PRODUCT PACKING** NO. 5748.0100 100 ml

Prepared from the highest purity raw material available, generally greater than 99.999% spectral purity. The content of the solution is confirmed to be accurate to within 0.2 %. Typically 10000 $\mu g/ml$.

The certificate of analysis provided reports actual lot analysis. The certificate also lists the trace impurities.

L(+)-Tartaric Acid

'BAKER ANALYZED' / ACS

► HOOC(CHOH),COOH M = 150.09 g/mol

CAS NO. 87-69-4 **EINECS** 201-766-0

NC CODE 2918 12 00

R: 36 **S:** 24/25



Meets ACS Specifications	
Assay	min. 99.0%
Chloride (CI)	max. 0.001%
Heavy Metals (as Pb)	max. 5 ppm
Insoluble Matter	max. 0.005%
Oxalate (C ₂ O ₄)	passes test
Phosphate (PO ₄)	max. 0.001%
Residue after Ignition	max. 0.02%
Sulfur Compounds (as S)	max. 0.002%
Trace Impurities (in ppm):	
Iron (Fe)	max. 5

PRODUCI	PACKING	CONT.
NO.		вох
0340.1000	1 kg	6

L(+)-Tartaric Acid

2054

► HOOC(CHOH)₂COOH

M = 150.09 g/mol

CAS NO. 87-69-4 **EINECS** 201-766-0

NC CODE 2918 12 00

R: 36

S: 24/25



Assay	99.5 - 101.0%
Appearance of solution	passes test
Calcium (Ca)	max. 200 ppm
Chlorides (as CI)	max. 100 ppm
Heavy Metals (as Pb)	max. 10 ppm
Identification	passes test
Loss on Drying	max. 0.2%
Oxalic acid	max. 350 ppm
Specific Optical Rotation	+ 12.0 - + 12.8
Sulfated Ash	max. 0.1%
Sulfates (as SO ₄)	max. 150 ppm

PRODUCT	PACKING	CONT.
NO.		вох
2054.1000	1 kg	
2054.9025	25 kg	
2054.9050	50 kg	

Tartaric Acid

See for qualities meeting the requirements of cGMP, USP/NF, BP/Ph.Eur, JP or FCC, page 36



Tellu

Tellurium 1000 μ g/ml

5783 (Matrix: 10% hydrochloric aci I) / 'BAKER INSTRA-ANALYZED' / Plasma Standard

► Te

M = 127.60 g/mol

NC CODE 3822 00 00

R: 36/37/38 **S:** 26-45

Certificate Provided Reporting Actual Lot Analysis

Tellurium (Te) 998-1002 μg/ml

PRODUCT CONT. **PACKING** NO. вох 5783.0100 100 ml

Prepared from the highest purity raw material available, generally greater than 99.999% spectral purity. The content of the solution is confirmed to be accurate to within \pm 0.2%. Typically 1000 $\mu\text{g/ml.}$

The certificate of analysis provided reports actual lot analysis. The certificate also lists the trace impurities.

Tellurium 10000 μg/ml /Μοστίν: 40% hydrochloric acid) / 'BAKER INSTRA-ANALYZED' / Plasma Standard 5749

► Te

M = 127.60 g/mol NC CODE 3822 00 00

UN/ID NO. 1789

ADR/RID 8 C1

IMDG 8/II **R:** 34-37

S: 20-26-36/37/39-45

Certificate Provided Reporting Actual Lot Analysis Tellurium (Te) 9980-10020 µg/ml

NO. 5749.0100 100 ml

PACKING

PRODUCT

CONT. BOX

Prepared from the highest purity raw material available, generally greater than 99.999% spectral purity. The content of the solution is confirmed to be accurate to within 0.2 %. Typically 10000 μ g/ml.

The certificate of analysis provided reports actual lot analysis. The certificate also lists the trace impurities.



TEMED

See N,N,N',N'-Tetramethylethylene diamine

7069

1,1,2,2-Tetrabromoethane

'BAKER ANALYZED

► Br₂CHCHBr₂

M = 345.65 g/mol

11 = 2.96 kg

CAS NO. 79-27-6

FINECS 201-191-5

NC CODE 2903 30 36

EC NO. 602 016 00 9

UN/ID NO. 2504

ADR/RID 6.1 T1

IMDG 6.1/III

R: 26-36-52/53

S: 24-27-45-61



min. 98% Assay Color (APHA) max. 75 Density (g/ml) at 25°C 2.946-2.956 Free Halogens passes test Inorganic halides (as CI) max. 0.0015%

CONT. PRODUCT PACKING вох NO. 7069.0500

2',4',5',7'-Tetrabromofluorescein

3',3,5',5-Tetrabromophenolsulfonphthalein



Tetrabutylammonium Dihydrogen Phosphate

Tetrabutylammonium Hydrogen Sulfate BAKER HPLC ANALYZED

► (CH₃CH₂CH₂CH₂)₄NHSO₄ **M** = 339.54 g/mol **CAS NO.** 32503-27-8 EINECS 251-068-5

NC CODE 2923 90 00

For Ion-Pair Chromatography Assay min. 98% Melting Range 167 - 173°C

PRODUCT PACKING CONT. BOX 2846.0500 500 g

Tetrabutylammonium Hydroxide 0.4M in water / Titrant / 'BAKER HPLC ANALYZED'

2843

► (CH ₃ CH ₂ CH ₂ CH ₂) ₄ NOH			
M = 259.48 g/mol			
CAS NO.	5574-97-0		
EINECS 218-147-6			
NC CODE 2923 90 00			
UN/ID NO.	3267		
ADR/RID	8 C7		

IMDG 8/II **R:** 34

Absorbance (neat) at 420 nm	max. 0.15
Appearance	passes test
Bromide (Br)	max. 0.05%
Molarity	0.37-0.43

PRODUCT	PACKING	CONT. BOX
2843.0500	500 g	ВОХ

S: 26-36/37/39-45

Tetrabutylammonium Phosphate

'BAKER HPLC ANALYZED

	CONT

► (CH ₃ CH ₂ CH ₂ CH ₂) ₄ NH ₂ PO ₄		
M =	339.46 g/mol	
CAS NO.	5574-97-0	
EINECS	226-947-1	
NC CODE	2923 90 00	

For Ion-Pair Chromatography		
Assay	min. 97.0%	
Appearance (white to off-white powder)	passes test	
Ultraviolet Absorbance (0.1% solution in water):		
at 196 nm	max. 1.00	
at 200 nm	max. 0.21	
at 210 nm	max. 0.15	
at 215-400 nm	max. 0.01	

BOX			NO.
		25 g Glass	2842.0025

PACKING

PACKING

4 | Glass

PRODUCT

PRODUCT

9360,4000

Tetrachloroethylene

'BAKER ULTRA RESI-ANALYZED' / For Trace Hydrocarbon Analysis by IR Spectrophotometry

9360

CONT.

BOX

► Cl₂C:CCl₂

M = 165.85 g/mol

1 l = 1.62 kg **CAS NO.** 127-18-4

EINECS 204-825-9

NC CODE 2903 23 00

EC NO. 602 028 00 4

UN/ID NO. 1897

ADR/RID 6.1 T1

IMDG 6.1/III

R: 40-51/53

S: 23-36/37-61





700 cm⁻¹,
max. 0.005%
max. 0.0005
max. 10 ppm
passes test
max. 1 ppm
min. 99.8%

(Hexadecane:isooctane:chlorobenzene standard

Hydrocarbons (excl. of stabilizer) max. 5 ppm

Tetrachloroethylene



M = 165.85 g/mol

1 l = 1.62 kg

CAS NO. 127-18-4

EINECS 204-825-9

NC CODE 2903 23 00 EC NO. 602 028 00 4

UN/ID NO. 1897

ADR/RID 6.1 T1

IMDG 6.1/III **R:** 40-51/53

S: 23-36/37-61





Assay (by GC)	min. 99%
Chloride (CI)	max. 0.001%
Residue after Evaporation	max. 0.005%

PRODUCT	PACKING	CONT.
NO.		вох
7165.1000	11	
7165.2500	2.5	4
7165.9025	25	
For safe han	dling of 25 l tin cans, see Se	lf-closing tap.

Stabilized with 0.015% 4-Methylmorpholine.



Tetra

Tetrachloromethane

7106 'BAKER HPLC ANALYZED'

M = 153.82 g/mol 1 I = 1.59 kg CAS NO. 56-23-5 EINECS 200-262-8

> NC CODE 2903 14 00 EC NO. 602 008 00 5

UN/ID NO. 1846 ADR/RID 6.1 T1 IMDG 6.1/II

R: 23/24/25-40-48/23-52/53-59

S: 23-36/37-45-59-61

N dangerous

dangerous toxic for the environment

 Assay (by GC)
 min. 99.7%

 Chloride (CI) (in ppm)
 max. 10

 Residue after Evaporation (in ppm)
 max. 3

 Titrable Acid (meq/g)
 max. 0.0005

 Water (H₂O)
 max. 0.01%

Fluorescence Trace Impurities (as quinine base), ppb:

Measured at 450 nm act. value reported
Measured at Emission Maximum for
Solvent Impurities act. value reported

Product Information (not specifications):

Density (g/ml) at 20°C

Ultraviolet Absorbance (1.00-cm path vs water):
at 280 nm

max. 0.06
at 330 nm

max. 0.01

max. 265

PRODUCT

 PRODUCT
 PACKING
 CONT.

 NO.
 BOX

 7106.2500
 2.5 l
 4

Filtered through a 0.2 micron filter. Packaged under Nitrogen.

PACKING

Tetrachloromethane

7306 'BAKER INSTRA-ANALYZED' / GC-Spectrophotometric quality

UV Cut-off, nm

M = 153.82 g/mol 1 I = 1.59 kg CAS NO. 56-23-5 EINECS 200-262-8 NC CODE 2903 14 00 EC NO. 602 008 00 5

EC NO. 602 008 00 UN/ID NO. 1846 ADR/RID 6.1 T1

> IMDG 6.1/II R: 23/24/25-40-48/23-52/53-59

S: 23-36/37-45-59-61



environment

gerous toxic

Assay (by GC) (corrected for water) min. 99.9% Color (APHA) max. 10 Free Chlorine passes test Infrared Transmittance (10 mm path vs air) 2800 - 3200 cm⁻¹ min. 80% lodine-consuming Substances passes test Residue after Evaporation max. 0.001% Substances Darkened by H₂SO₄ passes test Suitability for Use in Dithizone Test passes test Sulfur Compounds (as S) passes test Water Soluble Titrable Acid (meq/g) max. 0.0005 Ultraviolet Absorbance (1.00-cm path vs water): at 265 nm

 Ultraviolet Absorbance (1.00-cm path vs water):

 at 265 nm
 max. 1.00

 at 270 nm
 max. 0.35

 at 280 nm
 max. 0.10

 at 290 nm
 max. 0.05

 at 300 nm
 max. 0.02

 at 330 nm
 max. 0.01

 at 400 nm
 max. 0.01

NO. BOX 7306.1000 1 I 6

CONT.

Tetrachloromethane

7016 'BAKER ANALYZED' / Ultraviolet Spectrophotometry

► CCI₄

M = 153.82 g/mol

1 I = 1.59 kg CAS NO. 56-23-5 EINECS 200-262-8

NC CODE 2903 14 00

EC NO. 602 008 00 5

UN/ID NO. 1846 **ADR/RID** 6.1 T1

IMDG 6.1/II

R: 23/24/25-40-48/23-52/53-59

S: 23-36/37-45-59-61



ngerous toxic the

min. 99.9% Assay (by GC) Color (APHA) max. 10 Free Chlorine passes test Iodine-consuming Substances passes test Residue after Evaporation max. 0.001% Substances Darkened by H2SO4 passes test Suitability for Use in Dithizone Test passes test Sulfur Compounds (as S) passes test Water Soluble Titrable Acid (meq/g) max. 0.0005 Ultraviolet Absorbance (1.00-cm path vs water):

at 265 nm max. 1.00
at 270 nm max. 0.35
at 280 nm max. 0.10
at 290 nm max. 0.05
at 300 nm max. 0.02
at 330 nm max. 0.01
at 400 nm max. 0.01

PRODUCT	PACKING	CONT.
NO.		вох
7016.1000	11	6
7016.2500	2.5	4

Suitable for dithizone.



Tetrachloromethane

7420 'BAKER ANALYZED' / For Dithize

► CCI₄

M = 153.82 g/mol

1 l = 1.59 kg

CAS NO. 56-23-5

EINECS 200-262-8

NC CODE 2903 14 00

EC NO. 602 008 00 5

UN/ID NO. 1846

ADR/RID 6.1 T1

IMDG 6.1/II

R: 23/24/25-40-48/23-52/53-59

S: 23-36/37-45-59-61





environment

Assay (by GC)	min. 99.9%
Color (APHA)	max. 10
Free Chlorine	passes test
lodine-consuming Substances	passes test
Residue after Evaporation	max. 0.001%
Substances Darkened by H ₂ SO ₄	passes test
Suitability for Use in Dithizone Test	passes test
Sulfur Compounds (as S)	max. 0.005%
Water Soluble Titrable Acid (meq/g)	max. 0.0005
Trace Impurities (in ppm):	
Aluminium (AI)	max. 0.5
Barium (Ba)	max. 0.1
Boron (B)	max. 0.02
Cadmium (Cd)	max. 0.05
Calcium (Ca)	max. 0.5
Chromium (Cr)	max. 0.02
Cobalt (Co)	max. 0.02
Copper (Cu)	max. 0.02
Iron (Fe)	max. 0.1
Lead (Pb)	max. 0.1
Magnesium (Mg)	max. 0.1
Manganese (Mn)	max. 0.02
Nickel (Ni)	max. 0.02
Tin (Sn)	max. 0.1
Zinc (Zn)	max. 0.1

PRODUCT	PACKING	CONT.
NO.		вох
7420.1000	11	6
7420.2500	2.5	4

n-Tetradecane

8192

CONT.

► CH₃(CH₂)₁2CH₃

M = 198.40 g/mol

1 l = 0.76 kg

FLASHPOINT 100 °C

CAS NO. 629-59-4

EINECS 211-096-0

NC CODE 2901 10 00

ADR/RID 3,32c

Tetraethylrhodamine



Tetrahydrofuran

' 'BAKER HPLC ANALYZED' / for use in Liquid Chromatography, Organic and Biosynthesis

Residue after Evaporation

Color (APHA)

Peroxide (as H₂O₂)

Assay (by GC) (corrected for water)

Assay (by GC)

Freezing Point

9439

▶ QCH₂CH₂CH₂CH₂

M = 72.11 g/mol

1 | = 0.88 kg

FLASHPOINT −20 °C **CAS NO.** 109-99-9

EINECS 203-726-8

NC CODE 2932 11 00

EC NO. 603 025 00 0

UN/ID NO. 2056

ADR/RID 3T1

IMDG 3/II

R: 11-19-36-37

S: 16-29-33





Water (by KF, coulometric) max. 50 ppm Ultraviolet Absorbance (1.00-cm path vs water): at 230 nm max. 0.40 at 254 nm max. 0.15 at 280 nm max. 0.030 at 320-400 nm max. 0.005 UV Cut-off, nm max. 212

PRODUCT NO.	PACKING	CONT. BOX

Contains no preservative.

min. 98%

min. 99.8%

max. 10 ppm

max. 2 ppm

max. 10

5-7 °C

PRODUCT

8192.1000

NO.

PACKING

11



Tetra

Tetrahydrofuran

BAKER HPLC ANALYZED' / for Use in Liquid Chromatography and Spectrophotometry

1) At time of lot release

9440 ▶ QCH₂CH₂CH₂CH₂CH₂

M = 72.11 g/mol

1 l = 0.88 kg

FLASHPOINT $-20 \, ^{\circ}\text{C}$

CAS NO. 109-99-9

EINECS 203-726-8

NC CODE 2932 11 00 EC NO. 603 025 00 0

UN/ID NO. 2056

ADR/RID 3T1

IMDG 3/II

R: 11-19-36-37

S: 16-29-33





Assay (by GC) (corrected for water)	min. 99.8%
Peroxide (as H ₂ O ₂) ¹)	max. 0.005%
Residue after Evaporation	max. 2 ppm
Water (by KF, coulometric)	max. 0.02%
Ultraviolet Absorbance (1.00-cm path vs water):	
at 220 nm	may 0.45

at 230 nm max 0.45 at 254 nm max. 0.15 at 280 nm max. 0.030 at 320-400 nm max. 0.01 UV Cut-off, nm max. 212

PRODUCT PACKING CONT. NO. вох 9440.1000 11

Contains no preservative.

Tetrahydrofuran

BAKER HPLC ANALYZED' / for use in High Performance Liquid Chromatography

Assay (by GC)

1) At time of lot release

9441 ▶ QCH₂CH₂CH₂CH₂

M = 72.11 g/mol

1 l = 0.88 kg

FLASHPOINT −20 °C

CAS NO. 109-99-9

EINECS 203-726-8

NC CODE 2932 11 00

EC NO. 603 025 00 0

UN/ID NO. 2056

ADR/RID 3T1 IMDG 3/II

R: 11-19-36-37

S: 16-29-33

0



Peroxide (as H₂O₂)¹) max. 0.010% Residue after Evaporation max. 3 ppm max. 0.02% Water (H₂O) Product Information (not specifications): Density (g/ml) at 20°C 0.889 Ultraviolet Absorbance (1.00-cm path vs water): at 230 nm max. 1.0 at 260 nm max. 0.15 at 270 nm max. 0.08 at 310 nm max. 0.01

min. 99.8%

PRODUCT PACKING CONT. BOX 9441.1000 11 9441.2500 2.5 [

Filtered through a 0.2 micron filter. Packaged under Nitrogen.

Tetrahydrofuran

9446 BakerDRY / Low Water Solvent / ACS

► QCH₂CH₂CH₂CH₂

M = 72.11 g/mol

1 l = 0.88 kg

FLASHPOINT $-20 \, ^{\circ}\text{C}$ **CAS NO.** 109-99-9

EINECS 203-726-8

NC CODE 2932 11 00

EC NO. 603 025 00 0

UN/ID NO. 2056

ADR/RID 3T1

IMDG 3/II

R: 11-19-36-37

S: 16-29-33





Meets ACS Specifications Assay (by GC) (corrected for water) min. 99.8% Color (APHA) max. 10 Peroxide (as H₂O₂) max. 10 ppm Residue after Evaporation max. 1 ppm Water (by KF, coulometric) max. 10 ppm

PRODUCT PACKING CONT. BOX NO. 9446.1000 11

Packaged under Argon.

8117

Tetrahydrofuran

'BAKER ANALYZED'



M = 72.11 g/mol

1 l = 0.88 kg

FLASHPOINT $-20 \,^{\circ}\text{C}$ CAS NO. 109-99-9

EINECS 203-726-8

NC CODE 2932 11 00

EC NO. 603 025 00 0

UN/ID NO. 2056

ADR/RID 3 F1 IMDG 3/II

R: 11-19-36/37

S: 16-29-33





Assay (by GC)	min. 99%
Boiling Range	65-67°C
Color (APHA)	passes test
Density (g/ml) at 25°C	0.880-0.883
Peroxide (as H ₂ O ₂)	max. 0.02%
Residue after Evaporation	max. 0.002%
Water (H ₂ O)	max. 0.02%
Trace Impurities (in ppm):	
Aluminium (AI)	max. 0.5
Barium (Ba)	max. 0.1
Boron (B)	max. 0.02
Cadmium (Cd)	max. 0.05
Calcium (Ca)	max. 0.5
Chromium (Cr)	max. 0.02
Cobalt (Co)	max. 0.02
Copper (Cu)	max. 0.02
Iron (Fe)	max. 0.1
Lead (Pb)	max. 0.1
Magnesium (Mg)	max. 0.1
Manganese (Mn)	max. 0.02
Nickel (Ni)	max. 0.02
Tin (Sn)	max. 0.1
Zinc (Zn)	max. 0.1

PRODUCT	PACKING	CONT.
NO.		вох
8075.1000	11	6
8075.1000S	1 EcoTainer	
8075.2500	2.5	4
8075.5000	5 EcoTainer	
8075.9025	25 l	
market and the		and the state of the state of

EcoTainer, the metal solvent can for more safety in the lab. For safe handling of 25 l tin cans, see Self-closing tap.

Stabilized with 0.025% BHT.

Tetrahydrofuran

'BAKER'

QCH₂CH₂CH₂CH₂

M = 72.11 g/mol

1 I = 0.88 kg

 $\textbf{FLASHPOINT} \quad -20~^{\circ}\text{C}$

CAS NO. 109-99-9 **EINECS** 203-726-8

NC CODE 2932 11 00

EC NO. 603 025 00 0

UN/ID NO. 2056

ADR/RID 3 F1

IMDG 3/II

R: 11-19-36/37

S: 16-29-33





 Assay
 min. 99%

 Boiling Range
 65-68°C

 Water (H₂O)
 max. 0.2%

 PRODUCT NO.
 PACKING
 CONT. C

Stabilized with 0.025% BHT.

6131.0025

25 g Glass

Tetramethylammonium Hydroxide Pentahydrate

► (CH₃)₄NOH.5H₂O Assay min. 98% PRODUCT PACKING CONT.

M = 181.23 g/mol NO. BOX

M = 181.23 g/mol CAS NO. 75-59-2 EINECS 200-882-9

NC CODE 2923 90 00 UN/ID NO. 1835

ADR/RID 8 C7

IMDG 8/II

R: 25-34-52/53

S: 26-36/37/39-45-61



Tetramethylammonium Hydroxide (TMAH), Finyte Grade

See for detailed information section Microelectronic materials page 32 or the website www.jtbaker.com/microelectronic materials page 32 or the website was 32 or the website



6131

Tetramethylene Oxide

See Tetrahydrofura

Tetra

N,N,N',N'-Tetramethylethylene diamine

BAKER ULTRAPURE BIOREAGENT' / Catalyst for polyacrylamide gels

Assay

DNase Activity

RNase Activity

Protease Activity

Refractive Index at 25

► (CH₃)₂NCH₂CH₂N(CH₃)₂

4098

M = 116.21 g/mol

1 l = 0.78 kg

 $\textbf{FLASHPOINT} \quad 17~^{\circ}\text{C}$

CAS NO. 110-18-9

EINECS 203-744-6

NC CODE 2921 29 00

EC NO. 612 103 00 3

UN/ID NO. 2372

ADR/RID 3 F1

IMDG 3/II

R: 11-20/22-34

S: 16-26-36/37/39-45





	act. value reported
	none detected
	none detected
°C, n ²⁶	1.414-1.420
	none detected

PRODUCT NO.	PACKING	CONT. BOX
4098.0005	5 ml	
4098.0025	25 ml	
4098.0050	50 ml	

5784

Thallium 1000 μg/ml

(Matrix: 2% nitric acid) / Pl

M = 204.38 g/mol

NC CODE 3822 00 00

R: 20/22-36/38 **S:** 26-36/37



Certificate Provided Reporting Actual Lot Analysis 998-1002 $\mu g/ml$

PRODUCT PACKING CONT. NO. BOX 5784.0100 100 ml

Prepared from the highest purity raw material available, generally greater than 99.999% spectral purity. The content of the solution is confirmed to be accurate to within \pm 0.2%. Typically 1000 μ g/ml. The certificate of analysis provided reports actual lot analysis. The certificate also lists the trace impurities.

Thallium 10000 μg/ml

5761 (Matrix: 2% nitric acid) / 'BAKER ANALYZED' / Plasma Standard

► TI

M = 204.38 g/mol

NC CODE 3822 00 00

R: 23/25-33-36/38

S: 26-28-36/37/39-45



Certificate Provided Reporting Actual Lot Analysis Thallium (TI)

9980-10020 μg/ml

PRODUCT PACKING CONT. NO. BOX 5761.0100 100 ml

Prepared from the highest purity raw material available, generally greater than 99.999% spectral purity. The content of the solution is confirmed to be accurate to within 0.2 %. Typically 10000 μ g/ml.

The certificate of analysis provided reports actual lot analysis. The certificate also lists the trace impurities.



Thiamine Hydrochloride

ee for qualities meeting the requirements of cGMP, USP/NF, BP/Ph.Eur, JP or FCC, page 36

1200

Thioacetamide

crystal / 'BAKER ANALYZED'

► CH₃CSNH₂

M = 75.13 g/mol

CAS NO. 62-55-5

EINECS 200-541-4

NC CODE 2930 90 70 **EC NO.** 616 026 00 6

R: 22-36/38-45-52/53

S: 45-53-61



Assay	min. 99.0%
Melting Point	111.0-114.0°C.
Residue after Ignition	max. 0.05%

PRODUCT	PACKING	CONT.
NO.		вох
1200.0125	125 g	



Thiocarbamide

Thioethylene Glycol

Thiourea BAKER ANALYZED

0342

► (NH ₂) ₂ CS

M = 76.12 g/mol

CAS NO. 62-56-6

EINECS 200-543-5

NC CODE 2930 90 70

EC NO. 612 082 00 0

UN/ID NO. 2811

ADR/RID 6.1 T2

IMDG 6.1/III

R: 22-40-51/53-63

S: 22-36/37-61







Iron (Fe) max. 0.001% Loss on Drying max. 0.5% 174-177°C Melting Range Residue after Ignition max. 0.05% Sensitivity for Bismuth passes test Solubility in Water passes test

min. 99.0%

PRODUCT PACKING CONT. BOX NO. 0342.5000 5 kg 0342.9050 50 kg

Thorium 1000 μg/ml

(Matrix: 2% nitric acid) / 'BAKER INSTRA-ANALYZED' / Plasma Standard

Assay

5785

► Th

M = 232.04 g/mol

NC CODE 3822 00 00

R: 36/38

S: 26



Certificate Provided Reporting Actual Lot Analysis

Thorium (Th) 998-1002 µg/ml

PRODUCT	PACKING	CONT
NO.		вох
5785.0100	100 ml	

Prepared from the highest purity raw material available, generally greater than 99.999% spectral purity. The content of the solution is confirmed to be accurate to within \pm 0.2%. Typically 1000 μ g/ml.

The certificate of analysis provided reports actual lot analysis. The certificate also lists the trace impurities.

Thorium 1000 µg/ml

(Matrix: 2% nitric acid) / 'BAKER INSTRA-ANALYZED' / Atomic Absorption Standard

6954

BOX

► Th

M = 232.04 g/mol

NC CODE 3822 00 00

R: 36/38

S: 26-37

Thorium (Th)	998-1002 μg/ml

NO.	
6954.0100	100 ml
(054.0500	E00 ml

PACKING

PRODUCT

Prepared by dissolution of high purity raw materials (min. 99.99% spectral purity). Assays are verified by ICP against standards traceable to NIST. Standard Reference Material numbers (SRM) are printed on each label.

Thorium 10000 µg/ml

(Matrix: 5% nitric acid) / 'BAKER INSTRA-ANALYZED' / Plasma Standard

5750

► Th

M = 232.04 g/mol

NC CODE 3822 00 00

R: 34



S: 23-26-36-45

Certificate Provided Reporting Actual Lot Analysis 9980-10020 $\mu g/ml$

PRODUCT PACKING BOX 5750 0100 100 ml

Prepared from the highest purity raw material available, generally greater than 99.999% spectral purity. The content of the solution is confirmed to be accurate to within 0.2 %. Typically 10000 μ g/ml.

The certificate of analysis provided reports actual lot analysis. The certificate also lists the trace impurities.

L-Threonine

ting the requirements of cGMP, USP/NF, BP/Ph.Eur, JP or FCC, page 36



Thymol

See for qualities meeting the requirements of cGMP, USP/NF, BP/Ph.Eur, JP or FCC, page 36





Thymo

Thymolphthalein 1206 'BAKER ANALYZED' / ACS

> **CAS NO.** 125-20-2 **EINECS** 204-729-7 **NC CODE** 3204 19 00

Meets ACS Specifications	
Clarity of Solution	passes test
Visual Transition Interval:	
from pH 8.8(colourless) to pH 10.5(blue)	passes test

PRODUCT	PACKING	CONT.
NO.		вох
1206.0005	5 g	
1206.1000	1 kg	

Tin

9509 Shot / 'BAKER ANALYZED' / ACS

M = 118.69 g/mol
CAS NO. 7440-31-5
EINECS 231-141-8
NC CODE 8005 00 00

Exceeds ACS Specifications	
Assay (by EDTA titrn.)	min. 99.9%
Antimony (Sb)	max. 0.02%
Arsenic (As)	max. 1 ppm
Copper (Cu)	max. 0.002%
Iron (Fe)	max. 0.01%
Lead (Pb)	max. 0.005%
Zinc (Zn)	max. 0.005%

NO.		B0)
9509.0500 5	00 g	

Tin 1000 μg/ml

5786 (Matrix: 2% nitric acid plus a trace of hydrofluoric acid) / 'BAKER INSTRA-ANALYZED' / Plasma Standard

► Sn

M = 118.71 g/mol

NC CODE 3822 00 00

R: 36/38

6943



 $\begin{array}{c} \textit{Certificate Provided Reporting Actual Lot Analysis} \\ \text{Tin (Sn)} & 998\text{-}1002 \ \mu\text{g/ml} \end{array}$

 PRODUCT
 PACKING
 CONT.

 NO.
 BOX

 5786.0100
 100 ml

Prepared from the highest purity raw material available, generally greater than 99.999% spectral purity. The content of the solution is confirmed to be accurate to within \pm 0.2%. Typically 1000 $\mu g/ml$. The certificate of analysis provided reports actual lot analysis. The certificate also lists the trace impurities.

Tin 1000 μg/ml

0.10% (w/v) / (Matrix: 2% nitric acid plus a trace of hydrofluoric acid) / 'BAKER INSTRA-ANALYZED' / Atomic Absorption Standard

998-1002 µg/ml

Sn Tin (Sn)

M = 118.71 g/mol

NC CODE 3822 00 00

R: 20/21/22-36/38

 PRODUCT NO.
 PACKING BOX

 6943.0100
 100 ml

 6943.0500
 500 ml

S: 26-36/37



Prepared by dissolution of high purity raw materials (min. 99.99% spectral purity). Assays are verified by ICP against standards traceable to NIST. Standard Reference Material numbers (SRM) are printed on each label.

Tin 1000 μg/ml

'BAKER ANALYZED' / Atomic Absorption Standard

M = 118.71 g/mol NC CODE 3822 00 00 R: 36/37/38

► Sn

S: 26-37-45



Tin (Sn)	995-1005 μg/ml

PRODUCT	PACKING	CONT.
NO.		вох
6824.0100	100 ml	
6824.0500	500 ml	

Tin(II)choride in hydrochloric acid 5 mol/l.

www.jtbaker.com/europe

Tin 10000 μ g/ml

(Matrix: 2% nitric acid plus a trace of hydrofluoric acid) / 'BAKER ANALYZED' / Plasma Standard

5751

► Sn

M = 118.71 g/mol

NC CODE 3822 00 00

R: 36/38

S: 23-26-37



Certificate Provided Reporting Actual Lot Analysis

Tin (Sn) 9980-10020 μg/ml

PRODUCT	PACKING	CONT.
NO.		вох
5751 0100	100 ml	

Prepared from the highest purity raw material available, generally greater than 99.999% spectral purity. The content of the solution is confirmed to be accurate to within 0.2 %. Typically 10000 $\mu\text{g/ml.}$

The certificate of analysis provided reports actual lot analysis. The certificate also lists the trace impurities.

Tin(II) Chloride Dihydrate

0323

► SnCl₂.2H₂O

M = 225.63 g/mol

CAS NO. 10025-69-1

EINECS 231-868-0 NC CODE 2827 39 10

R: 22-36/37/38

S: 24/25-26



min. 98% Assav Ammonium (NH₄) max. 0.002% Arsenic (As) max. 0.0001% Copper (Cu) max. 0.001% max. 0.002% Iron (Fe) Lead (Pb) max. 0.005% Mercury (Hg) max. 0.01 ppm Substances not Precipitated by H₂S (as max. 0.05% SO₄)

Sulfate (SO₄) max. 0.002% **PRODUCT** PACKING CONT.

Tin(II) Chloride Dihydrate

'BAKER ANALYZED' / ACS

0325

0520

➤ SnCl₂.2H₂O

M = 225.63 g/mol

CAS NO. 10025-69-1

EINECS 231-868-0 NC CODE 2827 39 10

R: 22-36/37/38

S: 24/25-26



Exceeds ACS Specifications	
Assay	98.0-103.0%
Calcium (Ca)	max. 0.005%
Iron (Fe)	max. 0.003%
Lead (Pb)	max. 0.01%
Potassium (K)	max. 0.005%
Sodium (Na)	max. 0.01%
Solubility in HCI	passes test
Sulfate (SO ₄)	max. 0.002%

PRODUCI	PACKING	CONT.
NO.		вох
0325.0100	100 g	6
0325.0250	250 g	6
0325.1000	1 kg	6
0325.9050	50 kg	

Tin(II) Chloride Dihydrate

► SnCl₂.2H₂O

M = 225.63 g/mol

CAS NO. 10025-69-1

FINECS 231-868-0

NC CODE 2827 39 10 **R:** 22-36/37/38

S: 24/25-26



Assay	min. 95%
Iron (Fe)	max. 0.01%
Sulfate (SO ₄)	max. 0.02%
Trace Impurities (in ppm):	
Arsenic (As)	max 2

PRODUCT	PACKING	CONT.
NO.		вох
0520.1000	1 kg	
0520.9050	50 kg	

Mallinckrodt Baker's chemistry is Part of a pure process™.



TinIV

Tin(IV) Chloride Pentahydrate Lump / 'BAKER ANALYZED'

0322

➤ SnCl₄.5H₂O M = 350.58 g/mol **CAS NO.** 10026-06-9

EINECS 231-588-9 NC CODE 2827 39 10

EC NO. 50 001 00 5

UN/ID NO. 2440

ADR/RID 8 C2 IMDG 8/III

R: 34-52/53

S: 26-45-61-7/8

corrosive

Assay	99.0-102.0%
Arsenic (As)	max. 0.002%
Insoluble in HCl	max. 0.01%
Iron (Fe)	max. 0.003%
Sulfate (SO ₄)	max. 0.002%

PRODUCT NO.	PACKING	CONT. BOX
0322.0500	500 g	
0322.9050	50 kg	

Titanium 1000 μg/ml

5787 e of hydrofluoric acid) / 'BAKER INSTRA-ANALYZED'

▶ Ti **M** = 47.88 g/mol NC CODE 3822 00 00

> **R:** 36/38 S: 26

Certificate Provided Reporting Actual Lot Analysis Plasma Standard

998-1002 μg/ml Titanium (Ti)

PACKING **PRODUCT** CONT. вох NO. 5787.0100 100 ml

Prepared from the highest purity raw material available, generally greater than 99.999% spectral purity. The content of the solution is confirmed to be accurate to within \pm 0.2%. Typically 1000 $\mu g/ml$ The certificate of analysis provided reports actual lot

analysis. The certificate also lists the trace impurities.

Titanium 1000 μg/ml

(Matrix: 2% nitric acid plus a trace of hydrofluoric acid) / 'BAKER INSTRA-ANALYZED' / Atomic Absorption Standard

► Ti **M** = 47.88 g/mol NC CODE 3822 00 00

6964

R: 36/38 **S:** 26-36/37

998-1002 μg/ml Titanium (Ti)

PRODUCT PACKING CONT. NO. вох 6964.0100 100 ml 6964.0500 500 ml

Prepared by dissolution of high purity raw materials (min. 99.99% spectral purity). Assays are verified by ICP against standards traceable to NIST. Standard Reference Material numbers (SRM) are printed on each label.

CONT.

BOX

Titanium 1000 μg/ml

6825 'BAKER ANALYZED' / Ator rotion Standard

M = 47.88 g/molNC CODE 3822 00 00

► Ti

R: 36/37/38

S: 26-37-45

 $995\text{-}1005~\mu\text{g/ml}$ Titanium (Ti)

PRODUCT **PACKING** NO. 6825.0100 100 ml 6825.0500 500 ml

Titanium(III)chloride in hydrochloric acid 5 mol/l.

Titanium 10000 μg/ml

5752 (Matrix: 2% nitric acid plus a trace of hydrofluoric acid) / 'BAKER INSTRA-ANALYZED' / Plasma Standard

► Ti **M** = 47.88 g/mol NC CODE 3822 00 00 R: 36/38

S: 26-36/37



Certificate Provided Reporting Actual Lot Analysis Titanium (Ti) 9980-10020 μg/ml

PRODUCT PACKING CONT. NO. вох 5752.0100 100 ml

Prepared from the highest purity raw material available, generally greater than 99.999% spectral purity. The content of the solution is confirmed to be accurate to within 0.2 %. Typically 10000 µg/ml.

The certificate of analysis provided reports actual lot analysis. The certificate also lists the trace impurities.

TMAH



J.T.Baker

Titanium(IV) Chloride

► TiCl₄ min. 99% PRODUCT PACKING CONT. Assay (TiCl₄) max. 10 ppm **M** = 189.71 g/mol Iron (Fe) NO. вох **1 l** = 1.73 kg 9019.0500 500 ml **CAS NO.** 7550-45-0

EINECS 231-441-9 NC CODE 2827 39 80 **EC NO.** 22 001 00 5 UN/ID NO. 1838

R: 14-34

ADR/RID 8 C1 IMDG 8/II

S: 26-36/37/39-45-7/8



Titanium Dioxide

See Titanium(IV) Oxide



9019

Titanium(IV) Oxide

4962

$11O_2$		
	M =	79.90 g/mol
	CAS NO.	13463-67-7
	EINECS	236-675-5
	NC CODE	2823 00 00

Certificate Provided Reporting Actual Lot Analysis	
Actual Lot Analysis Lot No. B21421	
Assay (based on loss on ignition)	99.9%
Identification	passes test
Loss on Ignition at 1150°C	0.08%
Metallic Impurities in parts per million	on (μg/g):
Aluminium (AI)	< 0.2
Calcium (Ca)	2
Chromium (Cr)	9
Copper (Cu)	8
Iron (Fe)	1
Lead (Pb)	< 0.1
Magnesium (Mg)	< 0.1
Manganese (Mn)	< 0.1
Molybdenum (Mo)	< 0.1
Nickel (Ni)	< 0.2
Niobium (Nb)	< 0.2

Loss on ignition at 1100 C	0.0070
Metallic Impurities in parts per millio	n (μg/g):
Aluminium (AI)	< 0.2
Calcium (Ca)	2
Chromium (Cr)	9
Copper (Cu)	8
Iron (Fe)	1
Lead (Pb)	< 0.1
Magnesium (Mg)	< 0.1
Manganese (Mn)	< 0.1
Molybdenum (Mo)	< 0.1
Nickel (Ni)	< 0.2
Niobium (Nb)	< 0.2
Silver (Ag)	13
Strontium (Sr)	< 0.2
Vanadium (V)	4
Zirconium (Zr)	< 0.2
Non-Metallic Impurities in parts per r	million (μg/g):
Chloride (CI)	< 10

PRODUCT	PACKING	CONT.
NO.		ВОХ
4962.0100	100 g Glass	

For Laboratory, Research or Manufacturing Use.

Titanium Tetrachloride

See Titanium(IV) Chloride



TMAH

Silicon (Si)



Questions or suggestions, please contact us at jtbaker.nl@emea.tycohealthcare.com



Tolue

Toluene

9336 'BAKER ULTRA RESI-ANALYZED' / for Organic Residue Analysis

► C₆H₅CH₃

M = 92.14 g/mol

1 l = 0.86 kg

FLASHPOINT $4\,^\circ\mathrm{C}$

CAS NO. 108-88-3

EINECS 203-625-9 NC CODE 2902 30 00

EC NO. 601 021 00 3

UN/ID NO. 1294

ADR/RID 3F1

IMDG 3/II R: 11-20

S: 16-25-29-33





FCD Sensitive Impurities (as Hentac	hlor Fnoxide)
Water (by KF, coulometric)	max. 0.03%
Substances Darkened by H ₂ SO ₄	passes test
Residue after Evaporation	max. 1 ppm
Color (APHA)	max. 10
Assay (by GC) (corrected for water)	min. 99.7%

Single Impurity Peak (pg/ml)

FID-Sensitive Impurities (as 2-Octanol):

Single Impurity Peak (ng/ml) max. 10

PRODUCT	PACKING	CONT.
NO.		вох
9336.1000	11	
9336.2500	2.5	4

Toluene

'BAKER HPLC ANALYZED' / for use in High Performance Liquid Chromatography

Assay (by GC)

9351 ► C₆H₅CH₃

M = 92.14 g/mol

1 l = 0.86 kg

FLASHPOINT 4 °C

CAS NO. 108-88-3

EINECS 203-625-9

NC CODE 2902 30 00

EC NO. 601 021 00 3

UN/ID NO. 1294

ADR/RID 3 F1

IMDG 3/II

R: 11-38-48/20-63-65-67

S: 36/37-46-62





Substances Darkened by H₂SO₄ passes test Maximum Limits of Impurities: Residue after Evaporation (in ppm) Water (H₂O) 0.02% Physical Data (not specifications): Density (g/ml) at 20°C 0.867 Ultraviolet Absorbance (1.00-cm path vs water): at 288 nm max. 0.40 at 300 nm max. 0.10 at 350 nm max. 0.01 UV Cut-off, nm max. 285

min. 99.7%

passes test

3

PRODUCT PACKING BOX 9351.1000 11 6 9351.5000 5 | EcoTainer EcoTainer, the metal solvent can for more safety in the lab.

Filtered through a 0.2 micron filter. Packaged under Nitrogen.

Toluene

PHOTREX Reagent / For Spectrophotmetry and Liquid Scintillation Counting / ACS

► C₆H₅CH₃

8714

M = 92.14 g/mol

1 l = 0.86 kg

FLASHPOINT $4\,^\circ\mathrm{C}$

CAS NO. 108-88-3

EINECS 203-625-9

NC CODE 2902 30 00

EC NO. 601 021 00 3

UN/ID NO. 1294

ADR/RID 3 F1

IMDG 3/II **R:** 11-20

S: 16-25-29-33





Assay (by GC)	min. 99.5%
Color (APHA)	max. 5
Counting Efficiency for ³ H in Prepared	
'Cocktail'	act. value reported
Residue after Evaporation	max. 0.0005%
Substances Darkened by H ₂ SO ₄	passes test
Sulfur Compounds (as S)	max. 0.003%
Water (by KF, coulometric)	max. 0.02%
Optical Absorbance (1-cm path vs	water):
at 286 nm	max. 1.00
at 288 nm	max. 0.50
at 293 nm	max. 0.20
at 300 nm	max. 0.10
at 310 nm	max. 0.05
at 335 nm	max. 0.02
at 350 nm	max. 0.01
at 380-400 nm	max. 0.005
Product Information (not specifica	tions):
Boiling Point (typical)	110.6°C

Windows of Infrared Transmittance (0.1-mm path, 50-100%T), μ m:

PRODUCT	PACKING	CONT
NO.		BO
8714.0500	500 ml	

9364

Toluene

BakerDRY / Low Water Solvent / ACS

► C₆H₅CH₃

M = 92.14 g/mol

1 l = 0.86 kg

FLASHPOINT 4°C

CAS NO. 108-88-3

EINECS 203-625-9 NC CODE 2902 30 00

EC NO. 601 021 00 3

UN/ID NO. 1294

ADR/RID 3F1

IMDG 3/II

R: 11-38-48/20-63-65-67

S: 36/37-46-62





Meets ACS Specifications	
Assay (C ₆ H ₅ CH ₃)(by GC,corrected for	
water)	min. 99.5%
Boiling Point (typical)	110.6°C
Color (APHA)	max. 10
Residue after Evaporation	max. 1 ppm
Substances Darkened by H ₂ SO ₄	passes test
Sulfur Compounds (as S)	max. 0.003%
Water (by KF, coulometric)	max. 10 ppm

PRODUCT	PACKING	CONT.
NO.		вох
9364.1000	11	

Toluene

'BAKER ANALYZED' / ACS

8077

8078

► C₆H₅CH₃ **M =** 92.14 g/mol

1 l = 0.86 kg

FLASHPOINT 4 °C

CAS NO. 108-88-3 EINECS 203-625-9

NC CODE 2902 30 00

EC NO. 601 021 00 3

UN/ID NO. 1294

ADR/RID 3 F1

IMDG 3/II

R: 11-38-48/20-63-65-67

S: 36/37-46-62





Assay	min. 99.5%
Benzene	max. 0.1%
Color (APHA)	max. 10
Residue after Evaporation	max. 0.001%
Substances Darkened by H ₂ SO ₄	passes test
Sulfur Compounds (as S)	max. 0.003%
Water (H ₂ O)	max. 0.03%
Trace Impurities (in ppm):	
Aluminium (AI)	max. 0.5
Barium (Ba)	max. 0.1
Boron (B)	max. 0.02
Cadmium (Cd)	max. 0.05
Calcium (Ca)	max. 0.5
Chromium (Cr)	max. 0.02
Cobalt (Co)	max. 0.02
Copper (Cu)	max. 0.02
Iron (Fe)	max. 0.1
Lead (Pb)	max. 0.1
Magnesium (Mg)	max. 0.1
Manganese (Mn)	max. 0.02
Nickel (Ni)	max. 0.02
Tin (Sn)	max. 0.1

Exceeds ACS Specifications

Zinc (Zn)

PRODUCT	PACKING	CONT
NO.		BO)
8077.1000	11	(
8077.2500	2.5	4
8077.5000	5 EcoTainer	4
8077.9025	25 l	
8077.9200	200 l	
e e : 0	and the second	Carried Life

EcoTainer, the metal solvent can for more safety in the lab. For safe handling of 25 l tin cans, see Self-closing tap.

Toluene

M = 92.14 g/mol **1 l** = 0.86 kg

FLASHPOINT 4 °C

► C₆H₅CH₃

CAS NO. 108-88-3

EINECS 203-625-9

NC CODE 2902 30 00

EC NO. 601 021 00 3

UN/ID NO. 1294

ADR/RID 3 F1

IMDG 3/II

R: 11-38-48/20-63-65-67

S: 36/37-46-62





Assay (by GC)	min. 98%
Boiling Range	109-112°C
Residue after Evaporation	max. 0.002%
Water (H ₂ O)	max. 0.05%

max. 0.1

PRODUCT	PACKING	CONT
NO.		B0)
8078.2500	2.5	4
8078.5000	5 EcoTainer	4
8078.9025	25 l	
8078.9200	200 l	

EcoTainer, the metal solvent can for more safety in the lab. For safe handling of 25 l tin cans, see Self-closing tap.

Tolue

Toluene

3411 99% / HISTO GRADE

► C₆H₅CH₃

M = 92.14 g/mol

1 l = 0.86 kg

 $\textbf{FLASHPOINT} \hspace{0.2cm} 4 \hspace{0.1cm} ^{\circ} \! \mathbb{C}$

CAS NO. 108-88-3 **EINECS** 203-625-9

NC CODE 2902 30 00

EC NO. 601 021 00 3

UN/ID NO. 1294

ADR/RID 3 F1

IMDG 3/II

R: 11-38-48/20-63-65-67

S: 36/37-46-62







Histo-Grade implicates that this reagent is specially tested and therefore solely intended for use in histo-pathology applications. This reagent is of an analytical quality.



Toluene MOS, VLSI Grade

See for detailed information section Microelectronic materials page 32 or the website www.jtbaker.com/micro

Benzene

Toluidine Blue 0

0716 'BAKER

► (CH₃)₂NC₆H₃N:C₆H₂(CH₃)(NH₂):ŞCI

M = 305.83 g/mol

CAS NO. 92-31-9 **EINECS** 202-146-2

NC CODE 2934 30 90

Identification (by IR)

passes test

max. 0.02%

PRODUCT PACKING NO.0716.0025 25 g Glas

PACKING CONT.
BOX
25 g Glass

CONT.

вох

C.I. 52040

PRODUCT

NO.

Trace Metal Standard I

(Matrix: 5% nitric acid) / 'BAKER INSTRA-ANALYZED' / Plasma Standard

NC CODE 3822 00 00 UN/ID NO. 3264

ADR/RID 8,17c

6031-01

IMDG 8/III

R: 34 **S:** 23-26-36/37/39-45



Elemental Concent (µg/ml):		
Aluminium (Al)	500	
Arsenic (As)	100	
Beryllium (Be)	100	
Cadmium (Cd)	25	
Chromium (Cr)	100	
Cobalt (Co)	100	
Copper (Cu)	100	
Iron (Fe)	100	
Lead (Pb)	100	
Manganese (Mn)	100	
Mercury (Hg)	5	
Nickel (Ni)	100	
Selenium (Se)	25	
Vanadium (V)	250	
Zinc (Zn)	100	

6031-01	100 ml	
For use in I	PA SW-846 Methods 6010	and 200.7

For use in EPA SW-846 Methods 6010 and 200.2 Traceable to NIST.

PACKING

Trace Metal Standard III

(Matrix: 2% nitric acid) / 'BAKER INSTRA-ANALYZED' / Plasma Standard

NC CODE 3822 00 00

R: 36/38

S: 26-37



Elemental Concent (µg/ml):	
500	
500	
100	
500	
100	
500	

PRODUCT	PACKING	CONT.
NO.		вох
6033-01	100 ml	

For use in EPA SW-846 Methods 6010 and 200.7. Traceable to NIST.



6033-01

Tribromomethane

See Bromoform



Tributyl Phosphate

► (CH₃CH₂CH₂CH₂O)₃PO **PRODUCT** PACKING CONT. Appearance passes test **M** = 266.32 g/mol вох Identification (by IR) NO. passes test **1 l** = 0.98 kg 7074.1000 **FLASHPOINT** 146 °C

CAS NO. 126-73-8 **EINECS** 204-800-2 NC CODE 2919 00 10 **EC NO.** 15 014 00 2 **R:** 22 **S**: 25



Tricalcium Phosphate

See for qualities meeting the requirements of cGMP, USP/NF, BP/Ph.Eur, JP or FCC, page 36



7074

Trichloroacetaldehyde Hydrate



Trichloroacetic Acid



'BAKE	R ANALYZED' / For Tui	bidimetry / ACS				0344
► Cl ₃ CCOOH		Exceeds ACS Specifications		PRODUCT	PACKING	CONT.
M = 163.3	g/mol	Assay	min. 99.0%	NO.		вох
CAS NO. 76-03	.9	Chloride (CI)	max. 0.001%	0344.0250	250 g	6
EINECS 200-9	27-2	Heavy Metals (as Pb)	max. 0.002%	0344.1000	1 kg	6
NC CODE 2915	00 04	Insoluble Matter	max. 0.005%	0344.9050	50 kg	
EC NO. 607 0	04 00 7	Iron (Fe)	max. 0.001%			
UN/ID NO. 1839		Nitrate (NO ₃)	max. 0.002%			
ADR/RID 8 C4		Residue after Ignition	max. 0.02%			
IMDG 8/II		Substances Darkened by H ₂ SO ₄	passes test			
R: 35-50	/53	Sulfate (SO ₄)	max. 0.02%			
S: 26-36	/37/39-45-60-61	Trace Impurities (in ppm):				
	C N	Phosphate (PO ₄)	max. 5			
1	1/2					
corrosiv						
	for the environment					

1,2,4-Trichlorobenzene

_	-	-		
ч	л	л	7	

	BAKER ANALYZED / FOR U	ise in Gel Permeation Chromatography				9444
► C ₆ H ₃ Cl ₃		Assay (by GC)	min. 99%	PRODUCT	PACKING	CONT.
M =	181.45 g/mol	Maximum Limits of Impurities:		NO.		вох
CAS NO.	120-82-1	Residue after Evaporation (in ppm)	5	9444.2500	2.5	4
EINECS	204-428-0	Water (H ₂ O)	0.01%			
NC CODE	2903 69 90	Physical Data (not specifications):				
UN/ID NO.	2321	Boiling Point	213.5°C			
ADR/RID	6.1 T1	Density (g/ml) at 20°C	1.454			
IMDG	6.1/III	Refractive Index n ²⁰	1.572			
R:	22-36/37/38-51/53	Vapor Pressure at 100°C	20 mm Hg			
S:	26-61	Ultraviolet Absorbance (1.00-cm pati	vs water):			
	₩ N W Xn	at 350 nm	max. 0.15			
		at 375 nm	max. 0.05			
	dangerous harmful	at 400 nm	max. 0.01			
	for the environment	UV Cut-off, nm	max. 310			

Innovation is principal to our business.



Trich

1,2,4-Trichlorobenzene

Assay (by GC)

8430

► C₆H₃Cl₃ **M** = 181.45 g/mol

> **CAS NO.** 120-82-1 **EINECS** 204-428-0

NC CODE 2903 69 90 **EC NO.** 602 087 00 6

UN/ID NO. 2321

ADR/RID 6.1 T1

IMDG 6.1/III

R: 22-50/53

S: 23-37/39-60-61





dangerous for the environment **PRODUCT PACKING** CONT. NO. вох 8430.0100 100 ml 8430.1000 11 8430.2500 2.5 |

8430.9025

min. 99%

For safe handling of 25 l tin cans, see Self-closing tap.

1,1,1-Trichloroethane

7309

► CH₃CCI₃

M = 133.41 g/mol

1 l = 1.34 kg

CAS NO. 71-55-6

EINECS 200-756-3

NC CODE 2903 1980

EC NO. 602 013 00 2

UN/ID NO. 2831 **ADR/RID** 6.1 T1

IMDG 6.1/III

R: 20-59

S: 24/25-59-61



dangerous for the environment

Boiling Range (1 ml-95 ml) 70-88 °C Color (APHA) max. 15

PRODUCT PACKING CONT. NO. BOX 7309.1000 11 6 7309.2500 2.5 l

Certified by the Biological Stain Commission. Contains 4% stabilizer

PACKING

CONT.

вох

PRODUCT

NO.

max. 0.1 max. 0.02

max. 0.02

max. 0.1

max. 0.1

Trichloroethylene

► CICH:CCI₂

7075

M = 131.39 g/mol

11= 1.46 kg

FLASHPOINT 32 °C **CAS NO.** 79-01-6

EINECS 201-167-4

NC CODE 2903 22 00

EC NO. 602 027 00 9

UN/ID NO. 1710 ADR/RID 6.1 T1

IMDG 6.1/III

R: 36/38-45-52/53-67

S: 45-53-61



Assay (by GC)	min. 99%
Boiling Range (initial to dry point)	max. 1.0°C
Color (APHA)	max. 10
Density (g/ml) at 25°C	1.455-1.460
Free Halogens	passes test
Recorded Boiling Point	87.1°C
Residue after Evaporation	max. 0.002%
Titrable Acid (meq/g)	max. 0.0001
Titrable Base (meq/g)	max. 0.0003
Water (H ₂ O)	max. 0.01%
Trace Impurities (in ppm):	
Aluminium (AI)	max. 0.5
Barium (Ba)	max. 0.1
Boron (B)	max. 0.02
Cadmium (Cd)	max. 0.05
Calcium (Ca)	max. 0.5
Chromium (Cr)	max. 0.02
Cobalt (Co)	max. 0.02
Copper (Cu)	max. 0.02
Iron (Fe)	max. 0.1
Lead (Pb)	max. 0.1

Magnesium (Mg)

Manganese (Mn) Nickel (Ni)

Tin (Sn)

Zinc (Zn)

NU.		BU
7075.2500	2.5	
7075.9200	200 l	



Trichloroethylene MOS, VLSI Grade

See for detailed information section Microelectronic materials page 32 or the website www.jtbaker.com/micro



Trichloromethane



1,1,2-Trichloro-1,2,2-trifluorethane

9132

8691

CONT.

► Cl₂FCCF₂CI

M = 187.40 g/mol **1 l** = 1.57 kg

CAS NO. 76-13-1 **EINECS** 200-936-1

NC CODE 2903 43 00

R: 59

S: 61



dangerous for the environment

Assay (by GC)	min. 99.5%
Residue after Evaporation	max. 5 ppm
Water (H ₂ O)	max. 0.05%

FT-Infrared Absorbance (4-cm path vs air)

measured from 3125 to 2777 cm ⁻¹ T	77 cm ⁻¹ Transmission (T):	
at 2924 cm ⁻¹ (-CH ₂ -)	min. 45%	
at 2958 cm ⁻¹ (-CH ₃)	min. 40%	
at 3030 cm ⁻¹ (-C ₆ H ₆)	min. 30%	

PRODUCT	PACKING	CONT.
NO.		вох
9132.2500	2.5	4

The E.U. Comittee has forbidden the use of FREON 113 for oil in water analysis.

Tricresyl Phosphate

1 l = 1.16 kg **CAS NO.** 78-30-8

EINECS 201-103-5 NC CODE 2919 00 10 EC NO. 15 015 00 8

UN/ID NO. 2574 **ADR/RID** 6.1 T1


IMDG 6.1/II

R: 39/23/24/25-51/53

S: 20-21-28-45-61



dangerous for the



PRODUCT

8691.1000 11

From mixed Cresols.

PACKING

Triethanolamine

See for qualities meeting the requirements of cGMP, USP/NF, BP/Ph.Eur, JP or FCC, page 36



9111

CONT.

BOX

Triethylamine

BAKER BIO-ANALYZED

► (C₂H₅)₃N

M = 101.19 g/mol **11** = 0.73 kg

FLASHPOINT − 17 °C

CAS NO. 121-44-8

EINECS 204-469-4

NC CODE 2921 19 10 **EC NO.** 612 004 00 5

UN/ID NO. 1296

ADR/RID 3 FC

IMDG 3/II

R: 11-20/21/22-35

S: 16-26-29-3-36/37/39-45





Assay (by GC) (corrected for water)	min. 99.0%
Color (APHA)	max. 10
Water (H ₂ O)	max. 0.20%

i koboci	1 ACKING
NO.	
9111.0500GL	500 ml Glass

Triet

Triethylamine

8079 'BAKER

► (C₂H₅)₃N

M = 101.19 g/mol

1 l = 0.73 kg

 $\textbf{FLASHPOINT} \ -17~^{\circ}\text{C}$

CAS NO. 121-44-8

EINECS 204-469-4

NC CODE 2921 19 10

EC NO. 612 004 00 5

UN/ID NO. 1296

ADR/RID 3 FC

IMDG 3/II

R: 11-20/21/22-35

S: 16-26-29-3-36/37/39-45





Assay	min. 98.0%
Boiling Point	88-90°C
Color (APHA)	max. 15
Density (g/ml) at 20°C	0.725-0.734
Water (H _o O)	max. 0.2%

PRODUCT	PACKING	CONT.
NO.		вох
8079.1000	11	6
8079.2500	2.5	4

8689

Triethylene Glycol

'BAKE

► HO(CH₂CH₂O)₂CH₂CH₂OH

M = 150.18 g/mol

1 l = 1.12 kg

FLASHPOINT $166~^{\circ}\text{C}$

CAS NO. 112-27-6

EINECS 203-953-2

NC CODE 2909 49 19

S: 24/25

Color (APHA) max. 25 Identification (by IR) passes test

PRODUCT	PACKING	CONT.
NO.		вох
8689.1000	11	

9470 T

Trifluoroacetic Acid

'BAKER HPLC ANALYZED'

► F₃CCOOH

M = 114.03 g/mol

1 l = 1.49 kg

CAS NO. 76-05-1 **EINECS** 200-929-3

EINECS 200-929-3

NC CODE 2915 90 80

EC NO. 607 091 00 1

UN/ID NO. 2699

ADR/RID 8 C3

IMDG 8/I

R: 20-35-52/53

S: 26-27-28-45-61-9



Assay min. 99.5%
Heavy Metals (as Pb) max. 10 ppm
Residue after Evaporation max. 10 ppm

Gradient Elution Test (a.u.):
at 215 nm max. 0.05

Trace Impurities (in ppm):

Iron (Fe) max. 5

Ultraviolet Absorbance (1.00-cm path vs water):
at 230 nm max. 0.15
at 254 nm max. 0.01

max. 210

UV Cut-off, nm

 PRODUCT NO.
 PACKING PACKING
 CONT. BOX

 9470.0010
 1 ml x 10
 10 ml

 9470.0070
 70 ml
 10 ml

 9470.1000
 1 l
 10 ml

Gradient tested for HPLC Peptide and Protein Analysis.

Trifluoroacetic Acid

'BAKER ANALYZED'

► F₃CCOOH

7079

M = 114.03 g/mol

1 l = 1.49 kg

CAS NO. 76-05-1

EINECS 200-929-3

NC CODE 2915 90 80

EC NO. 607 091 00 1

UN/ID NO. 2699

UN/IU NU. 2098

ADR/RID 8 C3

IMDG 8/1

R: 20-35-52/53

S: 26-27-28-45-61-9



Assay	min. 99.0%
Heavy Metals (as Pb)	max. 0.001%
Residue after Evaporation	max. 0.01%
Trace Impurities (in ppm):	
Iron (Fe)	max. 1

PRODUCT	PACKING	CONT.
NO.		вох
7079.0500	500 ml	



Triglycol

1,2,3-Trihydroxybenzene



Trimethylcarbinol



2,2,4-Trimethylpentane

'BAKER RESI-ANALYZED / for Organic Residue Analysis

Assay (by GC)

► (CH ₃) ₃ CCH ₂ CH(CH ₃) ₂			
M = 114.23 g/mol			
1 l = 0.70 kg			
FLASHPOINT − 12 °C			

CAS NO. 540-84-1 EINECS 208-759-1 NC CODE 2901 10 00 EC NO. 601 009 00 8 **UN/ID NO.** 1262

ADR/RID 3 F1 IMDG 3/II R: 11-38-50/53-65-67

S: 16-29-33-60-61-62-9

environment





Color (APHA) max. 10 Residue after Evaporation max. 1 ppm Substances Darkened by H₂SO₄ passes test Water (H₂O) (by Coulometry) max. 0.03% ECD Sensitive Impurities (as Heptachlor Epoxide): Single Impurity Peak (pg/ml) FID-Sensitive Impurities (as 2-Octanol): Single Impurity Peak (ng/ml) max. 5 Neat solvent front characterization: ECD-Sensitive Impurities (as Ethylene Dibromide): Single Impurity Peak (ng/ml) max. 5

min. 99.8%

min. 99.8%

max. 0.01%

max. 2

PRODUCT

PRODUCT NO.	PACKING	CONT. BOX
9335.4000	4 Glass	4

2,2,4-Trimethylpentane

'BAKER HPLC ANALYZED' / for use in High Performance Liquid Chromatography

Assay (by GC)

Water (H₂O)

Residue after Evaporation (in ppm)

9480

M =	114.23 g/mol
1 =	0.69 kg
FLASHPOINT	− 12 °C
CAS NO.	540-84-1
EINECS	208-759-1

(CH₃)₃CCH₂CH(CH₃)₂

NC CODE 2901 10 00 EC NO. 601 009 00 8

UN/ID NO. 1262 ADR/RID 3 F1 IMDG 3/II

R: 11-38-50/53-65-67

S: 16-29-33-60-61-62-9







Water Soluble Titrable Acid (meq/g)	max. 0.0003
Physical Data (not specifications):	
Density (g/ml) at 20°C	0.690
Ultraviolet Absorbance (1.00-cm path vs water):	
at 225 nm	max. 0.10
at 254 nm	max. 0.015
at 280 nm	max. 0.01
at 350 nm	max. 0.01
UV Cut-off, nm	max. 205

BOX 9480.1000 9480.2500 2.5 [

Filtered through a 0.2 micron filter Packaged under Nitrogen.

PACKING

Analytical applications are available in our Technical Library at www.jtbaker.com/europe

Trime

2,2,4-Trimethylpentane 'BAKER ANALYZED' / Ultraviolet Spe

Exceeds ACS Specifications

Water Soluble Titrable Acid (meq/g)

Exceeds ACS Specifications

Ultraviolet Absorbance (1.00-cm path vs water):

Residue after Evaporation

Sulfur Compounds (as S)

min. 99.0%

max. 5 ppm

max. 0.005%

max. 0.0003

max. 1.0

max. 0.20

max. 0.10

max. 0.04

max. 0.01

max. 10

Assay (by GC)

Color (APHA)

at 210 nm

at 220 nm

at 230 nm

at 240 nm

at 250 - 400nm

CH ₃) ₃ CCH ₂ CH(CH ₃) ₂			
M = 114.23 g/mol			
1 l =	0.69 kg		
FLASHPOINT − 12 °C			
CAS NO.	540-84-1		
1 I = 0.69 kg FLASHPOINT - 12 °C CAS NO. 540-84-1			

EINECS 208-759-1 NC CODE 2901 10 00 EC NO. 601 009 00 8 **UN/ID NO.** 1262

ADR/RID 3F1 IMDG 3/II

R: 11-38-50/53-65-67

S: 16-29-33-60-61-62-9





PRODUCT	PACKING	CONT.
NO.		вох
8715.1000	11	6
8715.2500	2.5	4

2,2,4-Trimethylpentane

'BAKER ANALYZED' / ACS

► (CH₃)₃CCH₂CH(CH₃)₂

M = 114.23 g/mol **1 l** = 0.69 kg

FLASHPOINT −12 °C

CAS NO. 540-84-1 **EINECS** 208-759-1

NC CODE 2901 10 00 EC NO. 601 009 00 8

UN/ID NO. 1262 ADR/RID 3 F1

IMDG 3/II

R: 11-38-50/53-65-67

S: 16-29-33-60-61-62-9







Assay (by GC)	mın. 99.0%
Color (APHA)	max. 10
Density (g/ml) at 25°C	max. 0.690
Residue after Evaporation	max. 0.001%
Sulfur Compounds (as S)	max. 0.005%
Water Soluble Titrable Acid (meq/g)	max. 0.0003
Trace Impurities (in ppm):	
Aluminium (AI)	max. 0.5
Barium (Ba)	max. 0.1
Boron (B)	max. 0.02
Cadmium (Cd)	max. 0.05
Calcium (Ca)	max. 0.5
Chromium (Cr)	max. 0.02
Cobalt (Co)	max. 0.02
Copper (Cu)	max. 0.02
Iron (Fe)	max. 0.1
Lead (Pb)	max. 0.1
Magnesium (Mg)	max. 0.1
Manganese (Mn)	max. 0.02
Nickel (Ni)	max. 0.02
Tin (Sn)	max. 0.1
Zinc (Zn)	max. 0.1

PRODUCT	PACKING	CONT.
NO.		вох
8217.1000	11	6
8217.2500	2.5	4

2,4,6-Trimethylpyridine

Boiling Point

7081

► N:C(CH₃)CH:C(CH₃)CH:CCH₃

M = 121.18 g/mol

1 l = 0.91 kg FLASHPOINT 57 °C

CAS NO. 108-75-8

EINECS 203-613-3

NC CODE 2933 39 99

UN/ID NO. 1992

ADR/RID 3 FT1

IMDG 3/III

R: 10-22-36/38

S: 23



BOX 7081.0100 100 ml

PRODUCT

169-171°C

PACKING

CONT.



1,3,7-Trimethylxanthine



CONT.

Trioctylphosphine Oxide

1462

► [CH₃(CH₂)₇]₃PO

M = 386.65 g/mol

CAS NO. 78-50-2 **EINECS** 201-121-3

NC CODE 2931 00 95

R: 34-50/53 **S:** 26-29-36/39-45





50-55°C.

PRODUCT **PACKING** NO.

вох 1462.0100 100 g

Tripotassium Phosphate



TRIS Hydrochloride

See for qualities meeting the requirements of cGMP, USP/NF, BP/Ph.Eur, JP or FCC, page 36

Melting Point



Tris(hydroxymethyl)aminomethane 'BAKER ULTRAPURE BIOREAGENT'

4109

NH₂C(CH₂OH)₃

M = 121.14 g/mol

CAS NO. 77-86-1

EINECS 201-064-4

NC CODE 2922 1980

R: 36/38 **S:** 24/25-26



ror iviolecular Biology Applications		
Assay	99.0-101.0%	
Assay (dried basis)	min. 99.9%	
Appearance (white crystals)	passes test	
DNase Activity	none detected	
Identification	passes test	
Insoluble Matter	max. 0.005%	
pH of 0.1 M Solution at 25°C	10.0-11.0	
Protease Activity	none detected	
RNase Activity	none detected	
Water (H ₂ O)	max. 0.3%	

Absorbance of a 1 M Solution, Maximum (1-cm

•	•	
at 260 nm		0.06
at 280 nm		0.06
at 400 nm		0.01
Trace Imp	purities (in ppm):	
Arsenic (A	s)	max. 1
Calcium (C	Ca)	max. 1
Copper (C	u)	max. 1

Manganese (Mn)	max. 1
Description	passes test
Heavy Metals (as Pb)	max. 0.001%
Loss on Drying	max. 1.0%
Melting Range	168-172°C
Organic Volatile Impurities (1)	conforms to USP/NF

max. 1

10.0-11.5

max. 0.1%

PRODUCT	PACKING	CONT.
NO.		вох
4109.0500	500 g	
4109.1000	1 kg	
4109.5000	5 kg	
4109.9012	12 kg	

Material also meets the chemical specifications of Tromethamine, U.S.P.

www.jtbaker.com/europe

Iron (Fe) Lead (Pb) Magnesium (Mg)

pH (1 in 20)

Residue after Ignition

Trish

Tris(hydroxymethyl)aminomethane 'BAKER ANALYZED' / ACS

1414 NH₂C(CH₂OH)₃

M = 121.14 g/mol

CAS NO. 77-86-1 **EINECS** 201-064-4

NC CODE 2922 1980

R: 36/38

S: 24/25-26

Meets ACS Specifications	
Assay (dry basis)	99.8-100.1%
Absorbance	passes test
Heavy Metals (as Pb)	max. 5 ppm
Insoluble Matter	max. 0.005%
Iron (Fe)	max. 5 ppm
Water (H ₂ O)	max. 2%

min. 99.5%

PRODUCT

1533.1000

1533.9025

PRODUCT

4103.0500

4103.1000

4103.5000

NO.

PACKING

PACKING

500 g

1 kg

5 kg

1 kg

PRODUCT PACKING NO.	CONT.	
		вох
1414.0500	500 g	6
1414.1000	1 kg	6
1414.9050	50 kg	

CONT.

вох

CONT.

BOX

Tris(hydroxymethyl)aminomethane

NH₂C(CH₂OH)₃

M = 121.14 g/mol **CAS NO.** 77-86-1

EINECS 201-064-4 NC CODE 2922 1980

R: 36/38

S: 24/25-26



Tris (hydroxymethyl)aminomethane

ee for qualities meeting the requirements of cGMP, USP/NF, BP/Ph.Eur, JP or FCC, page 36

Assay

Tris(hydroxymethyl)aminomethane hydrochloride

4103 'BAKER ULTRAPURE BIOREAGENT'

► NH₂C(CH₂OH)₃.HCI

M = 157.60 g/mol **CAS NO.** 1185-53-1

EINECS 201-064-4

NC CODE 2922 19 90

R: 36/38

S: 24/25-26



For Liquid Chromatography	and Molecular Biology
applications	
Assay (dried basis)	min. 99.0%
Appearance	passes test
DNase Activity	none detected
Protease Activity	none detected
RNase Activity	none detected
Water (H ₂ O)	max. 0.5%
Absorbance of a 1 M Solution	on, Maximum (1-cm
path vs water):	
at 260 nm	0.06
000	

at 260 nm	0.06
at 280 nm	0.06
at 400 nm	0.01
Frace Impurities (in ppm):	
Calcium (Ca)	max. 1
Copper (Cu)	max. 1
ron (Fe)	max. 1
Lead (Pb)	max. 1
Magnesium (Mg)	max. 1

Triton X100

2838

M = 646.87 g/mol

1 l = 1.00 kg

CAS NO. 9002-93-1

NC CODE 3402 13 00

R: 22-41

S: 24-26-39



(Alkylaryl Polyether Alcohol)	
Color (APHA)	max. 60
Identification (by IR)	passes test

PRODUCT	PACKING	CONT.
NO.		вох
2838.0500	500 ml	
2838.5000	51	

Trolamine

See for qualities meeting the requirements of cGMP, USP/NF, BP/Ph.Eur, JP or FCC, page 36



Tromethamine, U.S.P.

See for qualities meeting the requirements of cGMP, USP/NF, BP/Ph.Eur, JP or FCC, page 36



L-Tryptophan

See for qualities meeting the requirements of cGMP, USP/NF, BP/Ph.Eur, JP or FCC, page 36



Tungstophosphoric Acid n-Hydrate

crystal / 'BAKER ANALYZED'



P ₂ O ₅ .24WO ₃ .nH ₂ O		
CAS NO.	12501-23-4	
EINECS	215-682-7	
NC CODE	2811 19 80	
UN/ID NO.	3260	
ADR/RID	8 C2	
IMDG	8/111	
R:	34	

Meets Reagent Specifications for testing USP/NF	
monographs	
Ammonium (NH ₄)	max. 0.005%
Chloride (CI)	max. 0.03%
Heavy Metals (as Pb)	max. 0.005%
Insoluble Matter	max. 0.01%
Iron (Fe)	max. 0.003%
Nitrate (NO ₃)	max. 0.005%
Sulfate (SO ₄)	max. 0.02%

PRODUCT	PACKING	CONT.
NO.		вох
1165.0125	125 g	

Tween 20 'BAKER'

S: 26-36/37/39-45

corrosive

1 l = 1.11 kg CAS NO. 9005-64-5 NC CODE 3402 13 00

PRODUCT	PACKING	CONT.
NO.		вох
7374.1000	11	

Tween 80 'BAKER'

1 l = 1.08 kg CAS NO. 9005-65-6 NC CODE 3402 13 00

PRODUCT	PACKING	CONT.
NO.		вох
7394.0500	500 ml	

L-Tyrosine

See for qualities meeting the requirements of cGMP, USP/NF, BP/Ph.Eur, JP or FCC, page 36



7394

ULSI chemicals

See for detailed information section Microelectronic materials page 32 or the website www.jtbaker.com/micro



3905

CONT.

UltraClear **HISTO GRADE**

FLASHPOINT > 62 °C

CAS NO. 90622-57-4 **EINECS** 292-459-0 NC CODE 2710 11 25

R: 65

Clearing agent for use in Histology and Cytology Aromates (as Xylene) max. 60 ppm **PRODUCT PACKING** NO. 3905.2500

BOX 2.5 | Glass 3905.5000 5 I 3905.9010 10 l 3905.9200 200 l



A biodegradable paraffin clearing agent. Contains Hydrocarbons.

Mounting medium for coverslipping of slides.

UltraKitt

HISTO GRADE / Mounting Medium

3921

NC CODE 3822 00 00	For use in Histology and Cytology	PRODUCT	PACKING	CONT.
		NO.		вох
		3921.0500	500 ml Glass	
		3921.0600	100 ml Glass x 6	

Ultra

UltraPar (54-56°C) 3925

HISTO GRADE

CAS NO. 8002-74-2 NC CODE 271 22 01

PRODUCT **PACKING** CONT. вох 3925.5000 5 kg

BOX

CONT.

BOX



Ultra pure acids

e for detailed information section Reagents for trace element analysis, page 21



Ultryte chemicals

detailed information section Microelectronic materials page 32 or the website www.jtbaker.com/micro

Assay (by GC)

Boiling Point

8271

n-Undecane

► CH₃(CH₂)₉CH₃

M = 156.31 g/mol

11 = 0.74 kg

FLASHPOINT 66 °C

CAS NO. 1120-21-4 EINECS 214-300-6

NC CODE 2901 10 00 **UN/ID NO.** 2330

ADR/RID 3 F1 IMDG 3/III

Uranin

See Fluorescein Disodium Salt

Uranium 1000 μg/ml

(Matrix: 1% nitric acid) / 'BAKER INSTRA-ANALYZED' / Plasma Standard

M = 238.03 g/mol

NC CODE 3822 00 00

R: 20/22-36/38



Certificate Provided Reporting Actual Lot Analysis

998-1002 μg/ml

min. 99%

195-196°C

PRODUCT

8271.0500

PACKING

500 ml

PRODUCT PACKING 5788.0100 100 ml

Prepared from the highest purity raw material available, generally greater than 99.999% spectral purity. The content of the solution is confirmed to be accurate to within \pm 0.2%. Typically 1000 μ g/ml. The certificate of analysis provided reports actual lot

analysis. The certificate also lists the trace impurities.

5753

Uranium 10000 μg/ml

(Matrix: 1% nitric acid) / 'BAKER INSTRA-ANALYZED' / Plasma Standard

► U

M = 238.03 g/mol

NC CODE 3822 00 00

R: 23/25-33-36/38 **S:** 26-36/39-45

Certificate Provided Reporting Actual Lot Analysis Uranium (U) 9980-10020 $\mu g/ml$

BOX 5753.0100 100 ml

Prepared from the highest purity raw material available, generally greater than 99.999% spectral purity. The content of the solution is confirmed to be accurate to within 0.2 %. Typically 10000 µg/ml.

The certificate of analysis provided reports actual lot analysis. The certificate also lists the trace impurities.

Innovation is principal to our business.

4111

Urea 'BAKER ULTRAPURE BIOREAGENT' / For Protein Solubilization and Denaturation

► NH₂CONH₂ **M** = 60.06 g/mol

CAS NO. 57-13-6 **EINECS** 200-315-5 NC CODE 3102 10 10

Assay (NH2CONH2) min. 99.5% Appearance passes test max. 0.01% Conductivity of 8.5M Solution, μ mho/cm max. 30 Cyanate (CNO) none detected **DNase Activity** none detected Insoluble Matter max. 0.005% Protease Activity none detected RNase Activity none detected Water (H₂O) max. 0.5% Absorbance of a 5M Solution (1-cm path vs water): at 260 nm max. 0.05 at 280 nm max. 0.03 Trace Impurities (in ppm): Chloride (CI) max. 1 Copper (Cu) max. 1 Iron (Fe) max. 1 Lead (Pb) max. 1

PRODUCT NO.	PACKING	CONT. BOX
4111.0500	500 g	
4111.2500	2.5 kg	
4111.9012	12 kg	

Urea

BAKER ANALYZED' / ACS

0345

► NH₂CONH₂

M = 60.06 g/mol **CAS NO.** 57-13-6 **EINECS** 200-315-5 NC CODE 3102 10 10

Meets ACS Specifications Heavy Metals (as Pb) max. 0.001% Insoluble Matter max. 0.01% Iron (Fe) max. 0.001% Melting Point 132-135°C Residue after Ignition max. 0.01% Sulfate (SO₄) max. 0.001% Trace Impurities (in ppm): Chloride (CI)

max. 5

PRODUCT	PACKING	CONT.
NO.		вох
0345.1000	1 kg	6
0345.9025	25 kg	

Urea 'BAKER

NH₂CONH₂

M = 60.06 g/mol **CAS NO.** 57-13-6 **EINECS** 200-315-5 NC CODE 3102 10 10

max. 0.007% Chloride (CI) Heavy Metals (as Pb) max. 0.002% Insoluble in Alcohol max. 0.04% Melting Range 132-135°C max. 0.1% Residue on Ignition Sulfate (SO₄) max. 0.010%

CONT. **PRODUCT PACKING** NO. вох 0346.1000 1 kg 6 0346.9025 25 ka

See for qualities meeting the requirements of cGMP, USP/NF, BP/Ph.Eur, JP or FCC, page 36



0346

L-Valine

See for qualities meeting the requirements of cGMP, USP/NF, BP/Ph.Eur, JP or FCC, page 36



Vanadium 1000 μg/ml

(Matrix: 2% nitric acid) / 'BAKER INSTRA-ANALYZED' / Plasma Standard

5789

V

M = 50.94 g/mol

1 l = 1.03 kg

NC CODE 3822 00 00 **EC NO.** 7 004 00 1

R: 36/38

S: 26



Certificate Provided Reporting Actual Lot Analysis

998-1002 µg/ml Vanadium (V)

PRODUCT PACKING CONT. BOX NO. 5789.0100 100 ml

Prepared from the highest purity raw material available, generally greater than 99.999% spectral purity. The content of the solution is confirmed to be accurate to within \pm 0.2%. Typically 1000 $\mu g/\text{ml}.$

The certificate of analysis provided reports actual lot analysis. The certificate also lists the trace impurities.

Vanad

Vanadium 1000 μg/ml

6945 (Matrix: 0.5 M sulfuric acid) / 'BAKER INSTRA-ANALYZED' / Atomic Absorption Standard

Vanadium (V)

M = 50.94 g/mol NC CODE 3822 00 00

> **R:** 36/38 **S**: 26

CONT. **PRODUCT PACKING** вох 6945.0100 100 ml 6945.0500 500 ml

Prepared by dissolution of high purity raw materials (min. 99.99% spectral purity). Assays are verified by ICP against standards traceable to NIST. Standard Reference Material numbers (SRM) are printed on each label.

Vanadium 1000 μg/ml

6826 'BAKER ANALYZED' / Atom

M = 50.94 g/molNC CODE 3822 00 00

998-1002 μg/ml Vanadium (V)

998-1002 μg/ml

min. 99.5%

max. 0.03%

max. 0.05%

max. 0.02%

CONT **PRODUCT PACKING** BOX 6826.0100 100 ml 6826.0500 500 ml

Vanadium(V)oxide in sulfuric acid 0.5 mol/l.

Vanadium 10000 μg/ml

5754 (Matrix: 5% nitric acid) / 'BAKER INSTRA-ANALYZED' / Plasma Standard

Assay

Iron (Fe)

Potassium (K)

Sodium (Na)

▶ V **M** = 50.94 g/mol

1 l = 1.03 kg NC CODE 3822 00 00 **EC NO.** 7 004 00 1

R: 34

S: 23-26-36-45



Certificate Provided Reporting Actual Lot Analysis 9980-10020 μg/ml Vanadium (V)

PRODUCT PACKING CONT. вох NO. 100 ml

Prepared from the highest purity raw material available, generally greater than 99.999% spectral purity. The content of the solution is confirmed to be accurate to within 0.2 %. Typically 10000 $\mu g/ml$.

The certificate of analysis provided reports actual lot analysis. The certificate also lists the trace impurities.

Vanadium(V) Oxide

1217 'BAKER ANALYZED

V₂O₅

M = 181.88 g/mol **CAS NO.** 1314-62-1 **EINECS** 215-239-8

NC CODE 2825 30 00

EC NO. 23 001 00 8

UN/ID NO. 2862

ADR/RID 6.1 T5 IMDG 6.1/II

R: 20/22-37-48/23-51/53-63-68

S: 36/37-38-45-61





PRODUCT PACKING CONT. вох 1217.0100

J.T.Baker: over 100 years of experience.

See chapter 1 of this catalogue.

Vanadium Pentoxide

See Vanadium(V) Oxide



See Malachite Green Oxalat

Vinylbenzene

See Styrene

Vitamin C

See L(+)-Ascorbic Acid



See for detailed information section Microelectronic materials page 32 or the website www.jtbaker.com/micro

Water

'BAKER ULTRA RESI-ANALYZED'

4417	
	-

H_2O	

M = 18.01 g/molCAS NO. 7732-18-5EINECS 231-791-2NC CODE 2851 00 10

Suitable for environmental inorganic and organic		
trace analysis		
Filtration Test	passes test	
Residue after Evaporation	max. 2 ppm	
For Use in EPA Methods:		
Drinking Water	500 series	
Inorganic Non-Metals	300 series	
Metals	200 series	
Solid Waste	846 series	
Waste Water	600 series	
EPA Contract Lab Program:		
Inorganic Target Analytes (below CQRL)	passes test	
Volatile Organic Trace Analysis: Gas		
Chromatography with Purge and Trap		

Chromatography with Purge and Trap		
concentration (EPA Contract Required Quantitation		
Limit-CRQL):		
Electroconductivity Detection (ELCD),		
Below CRQL	passes test	
Photoionization Detection (PID)	passes test	
Total Organic Carbon (100 ppb max.)	passes test	

PRODUCT	PACKING	CONT.
NO.		вох
4219.4000	4 Glass	4

Filtered through a 0.2 micron filter.

Resistivity during production is 10 megohm-cm minimum

The J.T.Baker CYCLE-TAINER High Purity Solvent Delivery System, preserves purity and protects people.

See chapter 3 of this catalogue for product details.



Water

Water

6906 **ULTREX II Ultrapure Reagent**

► H₂O

M = 18.01 g/mol **CAS NO.** 7732-18-5 **EINECS** 231-791-2 NC CODE 2851 00 10

Certificate Provided Reporting Actual Lot Analysis Trace Impurities (in ppt) (pg/g):

rrace impurities (in ppr) (pg/g).	
Aluminium (AI)	max. 20
Antimony (Sb)	max. 10
Arsenic (As)	max. 10
Barium (Ba)	max. 10
Beryllium (Be)	max. 10
Bismuth (Bi)	max. 10
Boron (B)	max. 50
Cadmium (Cd)	max. 10
Calcium (Ca)	max. 20
Cerium (Ce)	max. 10
Cesium (Cs)	max. 10
Chromium (Cr)	max. 10
Cobalt (Co)	max. 10
Copper (Cu)	max. 10
Dysprosium (Dy)	max. 1
Erbium (Er)	max. 1
Europium (Eu)	max. 1
Gadolinium (Gd)	max. 1
Gallium (Ga)	max. 10
Germanium (Ge)	max. 10
Gold (Au)	max. 10
Hafnium (Hf)	max. 1
Holmium (Ho)	max. 1
Indium (In)	max. 1
Iron (Fe)	max. 20
Lanthanum (La)	max. 1
Lead (Pb)	max. 10
Lithium (Li)	max. 10
Lutetium (Lu)	max. 1
Magnesium (Mg)	max. 10
Manganese (Mn)	max. 10
Mercury (Hg)	max. 20
Molybdenum (Mo)	max. 10
Neodymium (Nd)	max. 1
Nickel (Ni)	max. 10
Niobium (Nb)	max. 10
Palladium (Pd)	max. 10
Platinum (Pt)	max. 10
Potassium (K)	max. 10

NO.		вох
PRODUCT	PACKING	CONT.
Sulfate (SO ₄)		max. 200
Phosphate (P	O ₄)	max. 100
Chloride (CI)		max. 100
Trace Impui	ities in ppb (ng/g):	
Zirconium (Zi	-)	max. 10
Zinc (Zn)		max. 10
Yttrium (Y)		max. 1
Ytterbium (Yt	o)	max. 10
Vanadium (V		max. 10
Tungsten (W		max. 10
Titanium (Ti)		max. 10
Tin (Sn)		max. 50
Thulium (Tm)		max. 10
Thorium (Th)		max. 10
Thallium (TI)		max. 10
Terbium (Tb)	·	max. 10
Tellurium (Te	<u> </u>	max. 1
Tantalum (Ta	<i>'</i>	max. 10
Strontium (Sr	-)	max. 10
Sodium (Na)		max. 10
Silver (Ag)	1	max. 10
Selenium (Se	,	max. 50
Scandium (S	,	max. 10
Samarium (S	,	max. 10
Ruthenium (F	/	max. 10
Rhodium (Rh Rubidium (Rb	,	max. 10
Rhenium (Re	·	max. 10
DI . /D	\	10

Water

9823

BAKER ANALYZED LC-MS Reagent

► H₂O

M = 18.01 g/mol **CAS NO.** 7732-18-5 **EINECS** 231-791-2 NC CODE 2851 00 10

Certificate Provided Reporting Actual Lot Analysis Residue after Evaporation

max. 1 ppm LC-Gradient-Diode Array Detection (a.u.), test

solution is modified with 0.1% (v/v) formic acid: max. 0.001 at 254 nm

LC-MS Gradient Suitability Test (TIC, 100 to 2000

m/z), test solution is modified with 0.1% (v/v) formic acid:

Positive ESI-MS Sensitive Impurities (as

Reserpine) max. 50 ng/ml Product Information (not specifications):

1.00

Density (g/ml) at 20°C

Trace Impurities (in ppb):

Aluminium (AI) max. 50 Calcium (Ca) max. 50 Iron (Fe) max. 50 Magnesium (Mg) max. 50 Potassium (K) max. 50 Sodium (Na) max. 50

вох
6

Element concentrations are at time of lot release.

6906.1000

1 | PE

Water

'BAKER HPLC ANALYZED' / HPLC Gradient Grade

4218

4217

► H₂O

M = 18.01 g/mol **1** | 1 kg **CAS NO.** 7732-18-5 **EINECS** 231-791-2

NC CODE 2851 00 10

max. 0.0001% Particulate Matter Polarity Index 9.0 Residue after Evaporation max. 2 ppm Solvent Group 9

Ultraviolet Absorbance of largest eluted peak: max. 0.005 au

max. 0.001 au

PACKING	CONT
	BOX
11	(
2.5	4
	11

Water

BAKER ANALYZED



M = 18.01 g/mol **CAS NO.** 7732-18-5 **EINECS** 231-791-2 NC CODE 2851 00 10

Ammonium (NH₄) max. 0.3 ppm Conductivity at 20°C max. 5 μS/cm Heavy Metals (as Pb) passes test pH at 25°C 5.0-7.0 Residue after Evaporation max. 10 ppm Substances Reducing KMnO₄ passes test Trace Impurities in ppb (ng/g): Chloride (CI) max. 1 ppm Sulfate (SO₄) passes test

PRODUCT	PACKING	CONT.
NO.		вох
4217.9020	20 l Polycube	

Water

See for qualities meeting the requirements of cGMP, USP/NF, BP/Ph.Eur, JP or FCC, page 36

at 220 nm

at 254 nm



Water Determination

See for detailed information section Karl Fischer reagents, page 195




Water Standard I & II

 $\textbf{Matrix I: 5}\% \ nitric \ acid \ - \ \textbf{Matrix II: 20}\% \ \ \textbf{HCI} \ \ / \ \ \textbf{'BAKER INSTRA-ANALYZED'} \ \ / \ \ \textbf{Plasma Standard}$

6122-01

NC CODE 3822 00 00 **UN/ID NO.** 2031 ADR/RID 8,2b IMDG 8/II **R:** 34

S: 23-26-36/37/39-45 corrosive



PRODUCT NO.	PACKING	CONT. BOX
6122-01	100 ml x 2	

For EPA Contract Laboratory Program (CLP). Traceable to NIST.

Wijs' Solution

See Iodine Solution



Wood's Alloy

sticks / 'BAKER ANALYZED'

1219

CAS NO. 8049-22-7 NC CODE 8106 00 90

Bismuth (Bi)	49.0-51.0%
Cadmium (Cd)	12.0-13.0%
Lead (Pb)	24.5-25.5%
Melting Point	70.0-74.0°C
Tin (Sn)	12.0-13.0%

PRODUCT	PACKING	CONT.
NO.		ВОХ
1219.0500	500 g	



Wrigh

Wright

3816 r / Hematology / Oxidized Methylene Blue/Eosine Y Stain, suitable for staining blood smears

CAS NO. 68988-92-1

R: 22-36

S: 22-26-36



PRODUCT CONT. **PACKING** NO. вох 3816.0025 25 g Glass

Wright

3878

1 l = 0.79 kg FLASHPOINT 11C °C **UN/ID NO.** 1230 IMDG 3.2/II

R: 11-23/24/25-36-39/23/24/25

S: 16-20-4-7/9





	Т
toxic	

PRODUCT	PACKING	CONT.
NO.		вох
3878.1000	11	
3878.2500	2.5	

Xylene

'BAKER ANALYZED' / Ultraviolet Spectrophotometry

C₆H₄(CH₃)₂

9516

M = 106.17 g/mol

1 l = 0.86 kg FLASHPOINT 24 °C

CAS NO. 1330-20-7

EINECS 215-535-7

NC CODE 2902 44 00

EC NO. 601 022 00 9

UN/ID NO. 1307 ADR/RID 3 F1

IMDG 3/III

R: 10-20/21-38

S: 25



137-144°C Boiling Range (initial to dry point) Color (APHA) max. 10 Free Acid, μ eq/g max. 0.12 max. 0.001% Residue after Evaporation Substances Darkened by H₂SO₄ passes test Sulfur Compounds (as S) max. 0.003% Water (H₂O) max. 0.02% Ultraviolet Absorbance (1.00-cm path vs water):

Wrights stain for differetial blood picture staining

Wrights stain for differetial blood picture staining

at 295 nm max. 1.00 at 300 nm max. 0.30 max. 0.05 at 330 nm

max. 0.01

max. 0.02

max. 0.1

max. 0.1

at 380-400 nm

Nickel (Ni)

Tin (Sn)

Zinc (Zn)

PRODUCT PACKING CONT. NO. BOX 9516.1000 11

Contains Ethylbenzene.

Xylene

'BAKER ANALYZED' / ACS

C₆H₄(CH₃)₂

8080

M = 106.17 g/mol

1 l = 0.86 kg

 $\textbf{FLASHPOINT} \ \ 24\ ^{\circ}\text{C}$

CAS NO. 1330-20-7

EINECS 215-535-7 NC CODE 2902 44 00

EC NO. 601 022 00 9

UN/ID NO. 1307

ADR/RID 3 F1

IMDG 3/III

R: 10-20/21-38

S: 25



Assay	min. 98.5%
Color (APHA)	max. 10
Ethylbenzene	max. 4%
Residue after Evaporation	max. 0.002%
Substances Darkened by H ₂ SO ₄	passes test
Sulfur Compounds (as S)	max. 0.003%
Water (H ₂ O)	max. 0.03%
Trace Impurities (in ppm):	
Aluminium (AI)	max. 0.5
Barium (Ba)	max. 0.1
Boron (B)	max. 0.02
Cadmium (Cd)	max. 0.05
Calcium (Ca)	max. 0.5
Chromium (Cr)	max. 0.02
Cobalt (Co)	max. 0.02
Copper (Cu)	max. 0.02
Iron (Fe)	max. 0.1
Lead (Pb)	max. 0.1
Magnesium (Mg)	max. 0.1
Manganese (Mn)	max. 0.02

PRODUCT	PACKING	CONT
NO.		ВОХ
8080.1000	11	ϵ
8080.2500	2.5	4
8080.5000	5 EcoTainer	4
8080.9025	25 l	4
8080.9200	200 l	
EcoTainer, th	e metal solvent can for mo	re safety in the lab.
For safe hand	lling of 25 l tin cans, see Se	elf-closing tap.
Contains Eth		en-closing tap.
CONTAINS EU	ylborizorio.	
Contains Eth	ylbenzene.	

8118

3410

8146

Xylene

M = 106.17 g/mol

1 l = 0.86 kg

FLASHPOINT 24 °C

► C₆H₄(CH₃)₂

CAS NO. 1330-20-7

EINECS 215-535-7

NC CODE 2902 44 00

EC NO. 601 022 00 9

UN/ID NO. 1307

ADR/RID 3 F1

IMDG 3/III

R: 10-20/21-38

S: 25

Assay min. 98% Appearance passes test Color

max. 10

PRODUCT PACKING CONT. NO. вох 8118.1000 11 6 8118.2500 2.5 l 8118.5000 5 | EcoTainer 8118.9025 8118.9200

EcoTainer, the metal solvent can for more safety in the lab. For safe handling of 25 l tin cans, see Self-closing tap.

Xylene

M = 106.17 g/mol

11= 0.86 kg

FLASHPOINT 24 °C

C₆H₄(CH₃)₂

CAS NO. 1330-20-7

FINECS 215-535-7 NC CODE 2902 44 00

EC NO. 601 022 00 9

UN/ID NO. 1307

ADR/RID 3 F1

IMDG 3/III

R: 10-20/21-38

S: 25



Benzene

max. 0.02%

min 99%

138-140°C

PRODUCT PACKING CONT. NO. BOX 2.5 | HDPE 3410.2500PE 3410.5000PE 5 I HDPE 3410.9010 10 l 3410.9025 25 I

Histo-Grade implicates that this reagent is specially tested and therefore solely intended for use in histo-pathology applications. This reagent is of an analytical quality.

m-Xylene

C₆H₄(CH₃)₂ **M** = 106.17 g/mol

1 l = 0.86 kg

FLASHPOINT 25 °C

CAS NO. 108-38-3

EINECS 203-576-3

NC CODE 2902 42 00

FC NO. 601 022 00 9

UN/ID NO. 1307

ADR/RID 3 F1 IMDG 3/III

R: 10-20/21-38

×

Assay (by GC)

Boiling Point

PACKING CONT. **PRODUCT**

www.jtbaker.com/europe

Xylen

o-Xylene

8147

 $C_6H_4(CH_3)_2$ **M** = 106.17 g/mol

1 l = 0.88 kg

FLASHPOINT 30 $^{\circ}\mathrm{C}$

CAS NO. 95-47-6

EINECS 202-422-2

NC CODE 2902 41 00

EC NO. 601 038 00 6

UN/ID NO. 1307

ADR/RID 3 F1

IMDG 3/II

R: 10-20/21-38

S: 25





📘 Xylene MOS, VLSI Grade

See for detailed information section Microelectronic materials page 32 or the website www.jtbaker.com/micro

Assay (by GC)

Boiling Point



Xylene replacement for Histology

ee UltraClear



Xylene replacement for Histology

See for detailed information www.jtbaker.com and select Clinic

1421

Xylenol Orange Tetrasodium Salt

'BAKER ANALYZED' / ACS

C₃₁H₃₂N₂O₁₃S (asAcid)

M = 672.67 g/mol

CAS NO. 1611-35-4

EINECS 216-553-8

NC CODE 2934 30 90

Meets ACS Specifications Clarity of Solution Suitability for Zinc titration

 PRODUCT
 PACKING

 passes test
 NO.

 passes test
 1421.0001
 1 g

PRODUCT

8147.1000

NO.

min. 98%

144-145°C

PACKING

11

CONT. BOX

CONT.

вох

1565

D(+)-Xylose

'BAKER ANALYZED' Biochemical

► QCH₂(CHOH)₃CHOH

M = 150.13 g/mol

CAS NO. 58-86-6

EINECS 200-400-7

NC CODE 2940 00 00

Meets NAS/NRC Specifications and Criteria for Biochemical Compounds

Heavy Metals (as Cu)max. 0.001%Loss on Drying at 105°Cmax. 0.15%Residue after Ignitionmax. 0.05%

in water) $+ 18.2^{\circ}$ to $+ 19.4^{\circ}$

Trace Impurities (in ppm):

Arsenic (As) max. 0.5

PRODUCI	PACKING	CONI.
NO.		вох
1565.0100	100 g	

Yttrium 1000 µg/ml

(Matrix: 1% nitric acid) / 'BAKER INSTRA-ANALYZED' / Plasma Standard

Y

M = 88.91 g/mol

NC CODE 3822 00 00

R: 36/38



Certificate Provided Reporting Actual Lot AnalysisYttrium (Υ) 998-1002 μg/ml

 PRODUCT
 PACKING
 CONT.

 NO.
 BOX

 5790.0100
 100 ml

Prepared from the highest purity raw material available, generally greater than 99.999% spectral purity. The content of the solution is confirmed to be accurate to within \pm 0.2%. Typically 1000 $\mu g/ml$.

The certificate of analysis provided reports actual lot analysis. The certificate also lists the trace impurities.

Yttrium 10000 µg/ml (Matrix: 1% nitric acid) / 'BAKER INSTRA-ANALYZED' / Plasma Standard

5755

M = 88.91 g/mol

1 l = 1.03 kg

NC CODE 3822 00 00

R: 36/38 **S:** 26-37

Certificate Provided Reporting Actual Lot Analysis

Yttrium (Y) 9980-10020 μg/ml PRODUCT PACKING CONT. NO. вох 5755.0100 100 ml

Prepared from the highest purity raw material available, generally greater than 99.999% spectral purity. The content of the solution is confirmed to be accurate to within \pm 0.2%. Typically 1000 μ g/ml.

The certificate of analysis provided reports actual lot analysis. The certificate also lists the trace impurities.

Zinc

Granular 2 mm / 'BAKER ANALYZED' / ACS

0351

▶ 7n

M = 65.37 g/mol **CAS NO.** 7440-66-6 **EINECS** 231-175-3 NC CODE 7901 12 30 EC NO. 30 001 00 1

Exceeds ACS Specifications min. 99.8% Assay Iron (Fe) (by AAS) max. 0.01% Lead (Pb) (by AAS) max. 0.01% Suitability for determination of arsenic passes test Trace Impurities (in ppm): Arsenic (As) max. 0.1

PRODUCT **PACKING** CONT. BOX 0351.0500 500 g

Zinc

ılar 0.84 mm / 'BAKER ANALYZED' / ACS

0352

► Zn

M = 65.37 g/mol CAS NO. 7440-66-6 **EINECS** 231-175-3 NC CODE 7901 12 30 EC NO. 30 001 00 1

Exceeds ACS Specifications Assay min. 99.8% Iron (Fe) (by AAS) max. 0.01% Lead (Pb) (by AAS) max. 0.01% Suitability for determination of arsenic passes test Trace Impurities (in ppm): Arsenic (As) max. 0.1

PACKING CONT. BOX 500 g 0352.0500 0352.9025 25 ka

Zinc

Granular 0.42 mm (40 Mesh) / 'BAKER ANALYZED' / ACS

0354

Zn

M = 65.37 g/mol **CAS NO.** 7440-66-6 **EINECS** 231-175-3 NC CODE 7901 12 30 EC NO. 30 001 00 1

Meets ACS Specifications. Meets Reagent Specifications for testing USP/NF monographs Assay (by EDTA titrn.) min. 99.8% Iron (Fe) max. 0.01% max. 0.01% Lead (Pb)

Suitability for Arsenic Determination passes test Mesh: On U.S. No. 40 Sieve max. 5%

max. 30%

max. 0.1

Trace Impurities (in ppm):

Thru U.S. Sieve No. 70

Arsenic (As)

PRODUCT PACKING CONT. NO. BOX 0354.0500 500 g

Zinc

Powder / 'BAKER'

0356

▶ 7n

M = 65.37 g/mol

CAS NO. 7440-66-6 **EINECS** 231-175-3

NC CODE 7901 12 30

EC NO. 30 001 00 1

UN/ID NO. 1436 ADR/RID 4.3 WS

IMDG 4.3/II

R: 15-17-50/53

S: 43A-46-60-61





Assay (by EDTA titrn.)	min. 98%
Nitrogen Compounds (as N)	max. 0.002%
Trace Impurities (in ppm):	
Arsenic (As)	max. 0.5

PRODUCT	PACKING	CONT.
NO.		вох
0356.1000	1 kg	

Zinc

Zinc 1000 μ g/ml

5791 (Matrix: 2% nitric acid) / 'BAKER INSTRA-ANALYZED' / Plasma Standard

➤ Zn

M = 65.39 g/mol

NC CODE 3822 00 00

R: 36/38

S: 26-37-60



Certificate Provided Reporting Actual Lot Analysis

Zinc (Zn) 998-1002 μg/ml

998-1002 µg/ml

 $998\text{-}1002~\mu\text{g/ml}$

PRODUCT CONT. **PACKING** NO. вох 5791.0100 100 ml

Prepared from the highest purity raw material available, generally greater than 99.999% spectral purity. The content of the solution is confirmed to be accurate to within \pm 0.2%. Typically 1000 $\mu g/ml$

The certificate of analysis provided reports actual lot analysis. The certificate also lists the trace impurities.

Zinc 1000 μ g/ml (Matrix: 2% nitric acid) / 'BAKER INSTRA-ANALYZED' / Atomic Absorption Standard 6946

Zinc (Zn)

➤ Zn **M** = 65.39 g/mol NC CODE 3822 00 00

R: 36/38

S: 26



PRODUCT PACKING CONT. NO. BOX 6946.0100 100 ml 6946.0500 500 ml

Prepared by dissolution of high purity raw materials (min. 99.99% spectral purity). Assays are verified by ICP against standards traceable to NIST. Standard Reference Material numbers (SRM) are printed on each label.

Zinc 1000 μg/ml
'BAKER ANALYZED' / Atomic Absorption Standard 6827

➤ Zn Zinc (Zn)

M = 65.39 g/molNC CODE 3822 00 00

R: 36/38

S: 26-37



PRODUCT PACKING CONT. NO. вох 6827.0100 100 ml 6827.0500 500 ml

Zinc nitrate in nitric acid 0.5 mol/l.

Zinc 10000 μ g/ml

(Matrix: 2% nitric acid) / 'BAKER INSTRA-ANALYZED' / Plasma Standard

► Zn **M** = 65.39 g/mol

NC CODE 3822 00 00 R: 36/38

S: 26



Certificate Provided Reporting Actual Lot Analysis 9980-10020 μg/ml Zinc (Zn)

PRODUCT PACKING CONT. вох 5756.0100 100 ml

Prepared from the highest purity raw material available, generally greater than 99.999% spectral purity. The content of the solution is confirmed to be accurate to within 0.2 %. Typically 10000 μ g/ml.

The certificate of analysis provided reports actual lot analysis. The certificate also lists the trace impurities.

Zinc Acetate Dihydrate

'BAKER ANALYZED' / AC

► Zn(CH₃COO)₂.2H₂O

0357

M = 219.49 g/mol

CAS NO. 5970-45-6 EINECS 209-170-2

NC CODE 2915 29 00

R: 22

ĸ.	22
S:	25
	Xn harmful

Exceeds ACS Specifications	
Assay	99.0-101.0%
Calcium (Ca)	max. 0.005%
Chloride (CI)	max. 5 ppm
Insoluble Matter	max. 0.005%
Iron (Fe)	max. 5 ppm
Lead (Pb)	max. 0.002%
Magnesium (Mg)	max. 0.005%
pH of 5% Solution at 25°C	6.0-7.0
Potassium (K)	max. 0.01%
Sodium (Na)	max. 0.05%
Sulfate (SO ₄)	max. 0.002%

NO.		CONT BO)
0357.0500	500 g	Ć
0357.9050	50 kg	

Z

Zinc Bromide

granular / 'BAKER ANALYZED' 6123

► ZnBr₂

M = 225.18 g/molCAS NO. 7699-45-8EINECS 231-718-4

NC CODE 2827 59 00

UN/ID NO. 3077 **ADR/RID** 9 M7

, IMDG 9/III

R: 34

S: 26-36/37/39-45-7/8



69.5-71.0%
max. 0.005%
max. 0.001%
max. 0.005%
max. 0.01%

PRODUCT	PACKING	CONT
NO.		BOX
6123.0500	500 g	

Zinc Carbonate

powder / 'BAKER ANALYZED'

4312

➤ Zn(CO₃)₂

M = 185.39 g/molCAS NO. 3486-35-9NC CODE 2836 99 18

min. 70.0% Average Particle Diameter, μm (APD) (by Sedigraph)(typical) act value reported Bulk Density (g/cc)(typical) act value reported Calcium (Ca) act value reported Chloride (CI) max. 0.002% Insoluble in H₂SO₄ max. 0.02% Iron (Fe) max. 0.002% max. 0.005% Lead (Pb) Nitrate (NO₃) max. 0.005% Silicon (Si) act value reported Sodium (Na) act value reported Specific Surface Area, m^2/g (typical) act value reported Substances not Precipitated by (NH₄)₂S (as max. 0.40% Sulfate (SO₄) max. 0.01% Mesh (Wet Screen Analysis): On U.S. No. 325 Sieve act value reported

PRODUCT	PACKING	CONT.
NO.		вох
4312.0500	500 g	

Zinc Chloride

Lump / 'BAKER ANALYZED'

0359

➤ ZnCl₂

M = 136.28 g/mol CAS NO. 7646-85-7 EINECS 231-592-0

NC CODE 2827 36 00 EC NO. 30 003 00 2

UN/ID NO. 2331 **ADR/RID** 8 C2 **IMDG** 8/III

R: 22-34-50/53

S: 26-36/37/39-45-60-61





97.0-100.5%
max. 1.0%
passes test
passes test
max. 0.005%
passes test
passes test
passes test
max. 0.03%

 PRODUCT NO.
 PACKING SONT.

 0359.0500
 500 g

Innovation is principal to our business.

D

Zinc

Zinc Chloride

0360

► ZnCl₂ **M** = 136.28 g/mol

> **CAS NO.** 7646-85-7 **EINECS** 231-592-0

NC CODE 2827 36 00

EC NO. 30 003 00 2

UN/ID NO. 2331

ADR/RID 8 C2

IMDG 8/III

R: 34 **S:** 28-45-7/8

corrosive

Assay	min. 97.0%
Insoluble Matter	max. 0.1%
Lead (Pb)	max. 0.01%
pH of 5% Solution at 25°C	5.5-7.0

PRODUCT	PACKING	CONT.
NO.		вох
0360.1000	1 kg	
0360.9050	50 kg	



Zinc Chloride

See for qualities meeting the requirements of cGMP, USP/NF, BP/Ph.Eur, JP or FCC, page 36

Iron (Fe)

Assay

0362

Zinc Nitrate Hexahydrate

BAKER ANALYZED

► Zn(NO₃)₂.6H₂O

M = 297.47 g/mol

CAS NO. 10196-18-6

EINECS 231-943-8

NC CODE 2834 29 80

UN/ID NO. 1514

ADR/RID 5.1 O2 IMDG 5.1/II

R: 22-36/37/38-8

S: 24/25-26



Assay (by EDTA titrn.) 99.0-101.0% Chloride (CI) max. 0.002% Insoluble Matter max. 0.005% Lead (Pb) ma pH of 5% Solution at 25°C Sulfate (SO₄) ma Trace Impurities (in ppm): Copper (Cu)

ıx. 0.005%	03
3.5-5.5	03
x. 0.005%	
max. 5	
max. 5	

99-101%

PRODUCT PACKING CONT. BOX 0362.1000 1 kg 362.7100 100 lbs 62.9025 25 kg

1795

► Zn(NO₃)₂.6H₂O

Zinc Nitrate Hexahydrate

M = 297.47 g/mol

CAS NO. 10196-18-6

EINECS 231-943-8

NC CODE 2834 29 80

UN/ID NO. 1514

ADR/RID 5.1 O2 IMDG 5.1/II

R: 22-36/37/38-8

S: 24/25-26





PACKING CONT. **PRODUCT** вох 1795.1000 1 kg

Zinc Nitrate solution

10% in H₂O / 'BAKER ANALYZED'

➤ Zn(NO₃)₂.6H₂O **CAS NO.** 7779-88-6

7476

EINECS 231-943-8 NC CODE 2834 29 90 Assay

9.7 - 10.7%

PRODUCT **PACKING** NO.

CONT. вох 7476.9020 20 l Polycube

Volumetric Solution, ready for use.

Zinc Oxide

0363 'BAKER ANALYZED' / ACS

➤ ZnO

M = 81.37 g/mol

CAS NO. 1314-13-2 **EINECS** 215-222-5

NC CODE 2817 00 00

R: 50/53

S: 60-61

dangerous for the environment

Exceeds ACS Specifications	
Assay	min. 99.0%
Alkalinity	passes test
Arsenic (As)	max. 2 ppm
Calcium (Ca)	max. 0.005%
Chloride (CI)	max. 0.001%
Insoluble in Dilute H ₂ SO ₄	max. 0.01%
Iron (Fe)	max. 0.001%
Lead (Pb)	max. 0.005%
Magnesium (Mg)	max. 0.005%
Nitrate (NO ₃)	max. 0.003%
Potassium (K)	max. 0.01%
Sodium (Na)	max. 0.05%
Sulfur Compounds (as SO ₄)	max. 0.01%
Trace Impurities (in ppm):	
Manganese (Mn)	max. 5

PRODUCT NO.	PACKING	CONT. BOX
0363.1000	1 kg	
0363.9050	50 kg	

Zinc Oxide

'BAKER'

► ZnO

M = 81.37 g/mol

CAS NO. 1314-13-2

EINECS 215-222-5 NC CODE 2817 00 00

R: 50/53

S: 60-61



99.0-100.5% Assay Alkalinity passes test Arsenic (As) max. 5 ppm Cadmium (Cd) max. 10 ppm Carbonates and Substances Insoluble in Acids Identification

passes test passes test Iron (Fe) max. 200 ppm Lead (Pb) max. 50 ppm max. 1.0% Loss on Ignition

CONT. PRODUCT **PACKING** вох NO. 0364.1000 1 kg

Zinc Oxide

See for qualities meeting the requirements of cGMP, USP/NF, BP/Ph.Eur, JP or FCC, page 36



0364

Zinc Stearate

See for qualities meeting the requirements of cGMP, USP/NF, BP/Ph.Eur, JP or FCC, page 36

Iron (Fe)

Manganese (Mn)



0365

Zinc Sulfate Heptahydrate

'BAKER ANALYZED' / ACS

M = 287.54 g/mol **CAS NO.** 7446-20-0 **EINECS** 231-793-3 NC CODE 2833 26 00

➤ ZnSO₄.7H₂O

EC NO. 30 006 00 9

R: 36/38-50/53 **S:** 22-25-60-61





Exceeds ACS Specifications	
Assay	99.0-103.0%
Ammonium (NH ₄)	max. 0.001%
Calcium (Ca)	max. 0.005%
Insoluble Matter	max. 0.01%
Lead (Pb)	max. 0.003%
Magnesium (Mg)	max. 0.005%
Nitrate (NO ₃)	max. 0.002%
pH of 5% Solution at 25°C	4.4-6.0
Potassium (K)	max. 0.01%
Sodium (Na)	max. 0.05%
Trace Impurities (in ppm):	
Chloride (CI)	max. 5

max. 5

max. 3

PRODUCI	PACKING	CONI.
NO.		вох
0365.0100	100 g	
0365.1000	1 kg	
0365.9050	50 kg	

Zinc

Zinc Sulfate Heptahydrate

1796 'BAKER

➤ ZnSO₄.7H₂O

M = 287.54 g/mol **CAS NO.** 7446-20-0

EINECS 231-793-3

NC CODE 2833 26 00

EC NO. 30 006 00 9

R: 36/38-50/53

S: 22-25-60-61



environment



Assay	99.0-104.0%
Acidity	passes test
Alkalies and Alkaline Earths	max. 0.9%
Appearance of solution	passes test
Arsenic (As)	passes test
Chloride (CI)	max. 300 ppm
Iron (Fe)	max. 100 ppm
Identification	passes test
Lead (Pb)	max. 0.002%
рН	4.4-5.6

PRODUCT	PACKING	CONT.
NO.		вох
1796.1000	1 kg	

Preserve in tight containers.



Zinc Sulfate, 7-Hydrate

See for qualities meeting the requirements of cGMP, USP/NF, BP/Ph.Eur, JP or FCC, page 36

7254

Zinc Sulfate

0.1 mol/l / 'BAKER ANALYZED'

➤ ZnSO₄.7H₂O

M = 287.54 g/mol

1 l = 1.02 kg

CAS NO. 7446-20-0

EINECS 231-793-3

NC CODE 2833 26 00

R: 51/53

S: 57



dangerous for the environment Titer (mol/l)

0.0995-0.1005

NO.	
7254.1000	11
7254.9010	10 l

PACKING

PRODUCT

Zinc Sulfate

0.05 mol/l / 'BAKER ANALYZED'

CAS NO. 7733-02-0 **EINECS** 231-793-3 **NC CODE** 2833 26 00

Titer (mol/l) 0.0495 - 0.0505

 PRODUCT
 PACKING

 NO.
 20 l Polycube

Volumetric Solution, ready for use.

4874

7218

Zinc Sulfate

0.1 mol/l / DILUT-IT

1 l = 1.16 kg **CAS NO.** 7446-20-0 **NC CODE** 2833 26 00
 PRODUCT
 PACKING

 NO.
 4874

 1 amp.

Volumetric Concentrate, for dilution to 1 l.

5792 (

Zirconium 1000 μg/ml

(Matrix: 1% nitric acid) / 'BAKER INSTRA-ANALYZED' / Plasma Standard

➤ Zr

M = 91.22 g/mol **NC CODE** 3822 00 00

 Certificate Provided Reporting Actual Lot Analysis

 Zirconium (Zr)
 998-1002 μg/ml

μg/ml NO.

 PRODUCT
 PACKING

 NO.
 5792.0100

 100 ml

CONT. BOX

CONT.

BOX

CONT.

CONT.

BOX

BOX

6

Prepared from the highest purity raw material available, generally greater than 99.999% spectral purity. The content of the solution is confirmed to be accurate to within \pm 0.2%. Typically 1000 $\mu g/ml$. The certificate of analysis provided reports actual lot

The certificate of analysis provided reports actual lot analysis. The certificate also lists the trace impurities.



Zirconium 10000 µg/ml (Matrix: 1% nitric acid) / 'BAKER INSTRA-ANALYZED' / Plasma Standard

5757

▶ 7r

M = 91.22 g/mol

NC CODE 3822 00 00

R: 36/38 **S**: 26

Certificate Provided Reporting Actual Lot Analysis

Zirconium (Zr) 9980-10020 μg/ml

PRODUCT	PACKING	CONT
NO.		B0
5757.0100	100 ml	

Prepared from the highest purity raw material available, generally greater than 99.999% spectral purity. The content of the solution is confirmed to be accurate to within 0.2 %. Typically 10000 μ g/ml.

The certificate of analysis provided reports actual lot analysis. The certificate also lists the trace impurities.

Mallinckrodt Baker's cGMP Manufactured Chemicals for the Biopharmaceutical industry are a necessity for uncomplicated scale-up.

See chapter 6 of this catalogue.

Z



Liquid Chromatography (HPLC)

Fields of application for I.T.Baker HPLC solvents

- Pharmaceutical / vitamins
- Food / drinks / beverages / flavours / wine
- Water / environment
- Petrochemistry
- Biology / biochemistry
- Cosmetics
- Paints / lacquers / solvents
- Other chemical industries

We offer a complete product line:

- LC/MS solvents and reagents
- HPLC solvents
- HPLC acids
- HPLC buffer salts
- HPLC ion pair reagents

Quality definition of HPLC Acetonitrile, Methanol and Water

In analytical and preparative HPLC, acetonitrile, methanol and water are the most common and most critical mobile phases. Mallinckrodt Baker has an extensive HPLC solvent and reagent program. As an example, acetonitrile and water qualities versus the most common applications are represented on the next page.

Acetonitrile HPLC qualities versus applications	Preparative HPLC	Isocratic Grade HPLC	Far UV HPLC	Gradient elution 254 nm	Gradient elution 210 nm	Protein profiling	Pesticide Analysis	PAH analysis	Proteomics and LC/MS application
9821 LC/MS Grade									
9017 Ultra Gradient Grade									
9012 Far UV/Gradient Grade									
8275 Isocratic Grade									

Methanol HPLC qualities versus applications	Preparative HPLC	Isocratic Grade HPLC	Gradient elution 254 nm	Proteomics and LC/MS application
9822 LC/MS Grade				
8402 Gradient Grade				
8404 Isocratic Grade				

Water HPLC qualities versus applications	Preparative HPLC	Isocratic Grade HPLC	Gradient elution 254 nm	Gradient elution 210 nm	Proteomics and LC/MS application
9823 LC/MS Grade					
4218 Gradient Grade					

Chromatography products:

Solid-Phase Extraction Products

BAKERBOND Columns and Media
Biochromatography Products and Media



Solid-Phase Extraction Products



As a key to successful sample preparation, we offer a wide range of silica and polymer based Solid-Phase Extraction (SPE) products – traditional BAKERBOND spe columns and innovative BAKERBOND Speedisk columns and disks.

Following the demands in research and development for purity, as well as for improved detection and quantification limits in analytical techniques, we provide the market with a variety of J.T.Baker SPE products with guaranteed uniformity and reproducibility.

Silica and Polymer based SPE products

Silica based products

Mallinckrodt Baker carefully defines and controls critical surface chemistry parameters to ensure performance consistency. Our knowledge and experience have led to the development of a wide range of silica products - endcapped offering high hydrolytic stability and non-endcaped, used in extraction of more polar analytes. SPE silica-based sorbents provide predictable and consistent extractions for discrete subsets of a broad range of sample types. Mixed-mode sorbents are recommended for the extraction of compounds from more complex matrices. Our product range includes weak and strong cation and anion exchange sorbents.

Polymer based products

Among the wide range of SPE sorbents, Mallinckrodt Baker offers you polymeric sorbents that will improve the recovery of your sample preparation:

 BAKERBOND spe SDB phases that have large surface area, are highly rigid and stable over the whole pH range and involve surface-modified polymeric sorbents with different functionalities (ion-exchange moiety). BAKERBOND Speedisk polymer columns that feature BAKER polymer resins which are the product of J.T.Baker ultra-clean polymer micro particle technology. These patented polymeric sorbents are available in hydrophilic and hydrophobic forms as well as mixed mode ion exchangers. They are stable over pH 1-14 and highly recommended when advanced detection methods have to be used, offering the unique properties of a polymer.

Silica and polymer based SPE products are applicable in different chromatographic modes:

- reverse phase
- normal phase
- ion exchange
- mixed-mode sorbents
- adsorption sorbents
- · SPE columns for proteins and polynucleotides
- application related SPE columns

The family of SPE products includes:

- BAKERBOND spe Columns
- BAKERBOND Speedisk Columns
- BAKERBOND Speedisk 96 Columns
- BAKERBOND Speedisk 96-Well Plate
- BAKERBOND Speedisk Extraction Disks

Application related BAKERBOND spe and BAKERBOND Speedisk columns:

Extraction columns for drugs of abuse testing

Narc-1 columns are specifically formulated for rapid, reproducible extraction of Δ^9 THC-carboxylic acid from urine using a unique, patented carboxy ester bonded phase.

Narc-1 has a high selectivity for Δ^9 THC-carboxylic acid and provides highly consistent recoveries without co-extracting many other commonly used drugs. Narc-2 columns are mixed mode SPE phases (C_8 /Strong Cation Exchanger) for the extraction of basic compounds, for instance in important drug screening extractions of opiates, LSD, phencyclidine, amine-based drugs, cocaine, et cetera. Narc-2 can be used for basic drug screening, as well as acidic/neutral drugs.



For the extraction and clean-up of Policyclic Aromatic Hydrocarbons (PAH's), including the 16-priority pollutant EPA PAHs the J.T.Baker brand offers several applications using different SPE columns with combination (double) phases:

- for example BAKERBOND spe PAH SOIL (500 mg Cyano/1000mg SiOH) is designed for the clean-up of PAHs in soil extracts
- for example BAKERBOND spe PAH AQUA (1000 mg C₁₈ / 500 mg NH₂) is designed for the extraction of PAHs from water (DIN 38407)

Extraction columns for PCB applications

For the extraction and clean up of PCBs Mallinckrodt Baker designed two different columns:

- for example BAKERBOND spe PCB-N, (500mg Ar-SO₃ / 500 mg SiOH) for the extraction of PCBs from oil (DIN 51527, part 1)
- for example BAKERBOND spe PCB-A (500 mg SiOH/H₂SO₄ / 500 mg Ar-SO₃) in combination with BAKERBOND spe Silica Gel, for the extraction of PCBs in oil (dirty samples)


Extraction column for Acrylamide applications

For the extraction and purification of organic compounds such as acrylamide and other polar compounds we developed a spherical activated carbon:

 for example BAKERBOND spe Carbon, (6 ml/1000 mg)

Extraction columns for Pesticides applications

For the extraction of pesticides from water, we designed the following double phase column:





- for example BAKERBOND spe (500 mg C₁₈ / 200 mg SDB-1)
- for example BAKERBOND spe (250 mg C₁₈ Polar Plus / 100 mg SDB-1)

Extraction columns for Mineral Oil Index application

For the clean up of extracts of mineral oil Mallinckrodt Baker has exclusively developed several new products which are tested according to ISO 9377-2 including ready to use SPE glass columns (packed with Florisil at bottom and Sodium Sulfate on top)

 for example BAKERBOND spe Clean-up Column (2000 mg Na₂SO₄ / 2000 mg Florisil)

More comprehensive application related information can be found at: http://www.jtbaker.com/europe/ techlib/default.asp

BAKERBOND spe and BAKERBOND Speedisk

BAKERBOND spe

With BAKERBOND spe columns, silica and polymer based, you can choose the solid phase extraction column that best fits your sample size, performance requirements and equipment.


Use standard BAKERBOND spe 1, 3, 6 and 8 ml columns round rimmed and ear shaped, in different types of column housing (ultraclean polypropylene and glass). Selected standard columns and wide-mouth columns are available in a foil-lined fibre drum called the Jumbo Pack. The columns are loose packed in the foil liner to eliminate the inconvenience of opening individual pouches. In addition, Jumbo Packs provide cost savings compared to the standard package.

BAKERBOND Speedisk - High Performance SPE

The new *Speedisk* column for solid phase extraction features BAKERBOND HPLC particle chemistry that permits rapid, efficient separation without sacrificing capacity.

Speedisk columns:

- run 9 times faster than traditional SPE columns
- eliminate or shorten pre-filtration and evaporation steps.



Featuring a unique laminar configuration, Speedisk columns:

- operate with smaller solvent volumes and
- have higher capacity per milligram sorbent than conventional SPE columns.



The Speedisk design:

- shortens analysis times
- increases capacity and
- may eliminate pre-filtration and evaporation steps.

Columns come in 1 ml, 3 ml or 6 ml sizes, and are filled to levels from 10 mg to 200 mg with a choice of over 20 BAKERBOND sorbents.

Speedisk 96 Columns



Speedisk 96 Columns raise high throughput SPE to one level higher performance. 96 "rimless" Speedisk Columns are inserted in a Speedisk 96 Column Holder and the disposable assembly is ready out of the box for placement on a microplate processor such as the Speedisk 96 Processor.

Speedisk 96-Well Plate



Once method development is complete, use the *Speedisk* 96-well plate, a one-piece, moulded plate that is pre-assembled with the sorbent of your choice. Our new *Speedisk* 96-well plate was designed with standard geometries to *adapt to most popular automated liquid handling systems*. Plates are pre-assembled with the silica or polymer sorbent of your choice.

Table 7a Comparison between BAKERBOND spe columns and BAKERBOND Speedisk columns

Sample preparation step	BAKERBOND spe columns	BAKERBOND Speedisk columns
Column Size / Sorbent	1 ml / 100 mg	1 ml /20 mg
Particle Size	40 μm	25 μm
Sample Volume	2 ml	1 ml
Column conditioning	2 ml (20-40 sec)	0.5 ml (5-10 sec)
Sample addition	2 ml (100 sec)	50µl -0,5 ml (50 sec)
Washing	1.5 ml (15-20 sec)	0.4 ml (2-5 sec)
Elution	1-2 ml (15-20 sec)	0.3-0.6 ml (2-5 sec)
Sample concentration/ evaporation	3-10 minutes	reduced or eliminated

BAKERBOND Speedisk Extraction Disks

The correct choice for samples from 200 ml to 2 L is our patented BAKERBOND *Speedisk* Extraction *Disks*, pre-assembled for use in preparing aqueous samples for analysis.

BAKERBOND

Speedisk

Products are

Protected by U.S.

Patent No.

5,595,653



BAKERBOND *Speedisk* Extraction Disks shorten extraction time to less than one hour. Versatile, silica and polymer based, BAKERBOND *Speedisk* Extraction Disks can be used for clear and particle-laden samples

BAKERBOND Speedisk Extraction Disks:

- Ensure rapid run completion even with dirty samples
- Reduce solvent consumption and hazardous waste
- Improve precision with its optimised flow path design
- Provide additional technology options to meet EPA requirements

The patented BAKERBOND *Speedisk* extraction disk is neither cartridge nor membrane.

A thin bed of micro-particles of BAKERBOND sorbent is supported in a laminar structure to maintain speed and capacity and enhance reproducibility of adsorption.

The laminar configuration provides filtration capacity and inlet characteristics that maximise access of analyte molecules to the microparticulate sorbent. BAKERBOND Speedisk Extraction Disk design resists clogging and

- d
 e
 f
 BAKERBOND Speedisk
 C₁₈ 8055-06
- Retaining ring Screen
- c Glass fiber filter
- d Sorbent
- e Glass fiber filter
- f Screen

ensures high throughput rates even when samples contain solids.

Capacity, recovery, and precision are high due to the unique disk configuration and performance of BAKERBOND sorbent.

With BAKERBOND *Speedisk* Extraction Disks sample contamination is virtually eliminated:

- Your hands never touch the wetted parts of the pre-assembled disk
- The sorbent and disk housing are pre-cleaned
- Polyester packaging provides a barrier that repels moisture and eliminates the risk of contamination by plastic additives (e.g. phthalates)

BAKERBOND *Speedisk* Extraction Disks are compatible with J.T.Baker standard *vacuum processors and Speedisk Extraction Stations*.


Selection Guide for spe sorbents and solvents

Selection Guide for SPE Sorbents and Solvents

In order to help you in choosing the right sorbent and solvent for your solid phase extraction method, Mallinckrodt Baker provides you with the general guidance - the *Sorbent Selection Guide*.

This Sorbent Selection Guide is a systematic guide classifying samples according to polarity, ionisability and solubility in water or organic solvents. This information is useful in selecting the necessary components of a preliminary extraction method.

Table 7b Organic Samples MW< 2000 (in solution)



SDB: styrene divinyl benzene

DVB: divinyl benzene

 $\rm H_2O\text{-}Phobic\,WA\,DVB:\,Weak}$ anion exchanger $\rm H_2O\text{-}Phobic\,SC\,DVB:\,Strong\,cation\,exchanger$ $\rm H_2O\text{-}Philic\,SA\,DVB\,Strong\,anion\,exchanger$ $\rm H_2O\text{-}Philic\,SC\,DVB:\,Strong\,cation\,exchanger$

3, 4 Solvents:

9262 Hexane, ULTRA RESI-ANALYZED

9257 Chloroform, ULTRA RESI-ANALYZED

9264 Dichloromethane, ULTRA RESI-ANALYZED

9260 Ethyl acetate, ULTRA RESI-ANALYZED

9254 Acetone, ULTRA RESI-ANALYZED

9255 Acetonitrile, ULTRA RESI-ANALYZED

9077 Methanol , ULTRA RESI-ANALYZED

4219 Water, ULTRA RESI-ANALYZED



More Information

This Sorbent Selection Guide is also available on:


http://www.jtbaker.com/chromatography/SolidPhaseExtraction.asp

Processor and Extraction stations

Standard Vacuum Processors for Extraction Columns

J.T.Baker standard vacuum processors offer the flexibility of processing SPE devices of different heights, diameters, or formats during the same experiment. The vacuum box design is familiar

throughout the industry, and it supports all devices and accessories with luer-type fittings such as BAKERBOND spe, *Speedisk* columns as well as *Speedisk* extraction disks.




BAKER spe-12G Column Processor (PTFE Design) (7018-94)

Parts and Accessories

The BAKER spe-12G Column Processor (PTFE Design) complete includes: 1 borosilicate glass vacuum chamber, 1 white-coloured polyamide lid including 12 luer PTFE connectors, plugs (12x) for lid, 1 polyethylene gasket, 12 Luer PP Stopcocks, 1 PTFE sample collection rack set including height adjustable shelves, 1 vacuum gauge/PTFE controller assembly, 1 instruction sheet, 1 BAKERBOND spe Application Notes Manual, 1 Solid Phase Extraction for Sample Preparation Manual.

Description	Quantity per box	Product Numbe
Extraction Column Processor replacement parts for BAKER spe-12G Co	lumn Processor (PTFE desi	gn)
Borosilicate Glass Vacuum Chamber	1	7512-00
Polyamide Lid (with 12 PTFE luer connectors)	1	7513-00
PTFE Sample Collection Rack Set	1	7516-00
Polyethylene Gasket Seals	2	7430-00
Neoprene Gasket Seals	2	7433-00
Plugs for lid that fits into luer PTFE connectors (4586)	30	7517-00
Luer PTFE Stopcocks	12	7514-00
Luer PTFE Connectors (one-piece)	12	4586
Polypropylene Stopcocks	10	7241-00
Vacuum Gauge/PTFE Controller Assembly	1 Assembly	7515-00
Accessories for BAKER spe-12g Column Processor (PTFE design)		
Autosampler Plate	1	7516-01
Drying Top, polyamide for drying (columns) or evaporation (eluates) purpose	1	4581
PTFE Lined Luer locks (taps), Stainless Steel	12	4514
Stop Cocks, polypropylene	10	7241-00



BAKER spe-12G Column Processor (7018-00)

Parts and Accessories

BAKER spe-12G Column Processor complete includes:


1 glass vacuum chamber, 1 blue-coloured nylon cover
with 12 luer connectors and 1 polyethylene gasket,

12 individual flow control polypropylene stopcocks,

12 plugs for lid, 12 stainless steel needles, 1 sample
collection rack with 3 support posts, 3 height
adjustable shelves, 9 shelf support clips, 1 vacuum
gauge/polypropylene controller assembly,

1 instruction sheet, 1 BAKERBOND spe Application
Notes Manual, 1 Solid Phase Extraction for Sample
Preparation Manual.

Description	Quantity per box	Product Number
Extraction Column Processor replacement parts for BAKER s	spe-12G Column Processor	
Glass Vacuum Chamber	1	7421-00
Nylon (blue-coloured) Lid (with 12 luer connectors)	1	7424-00
Polypropylene Rack Set	1	7427-00
Polyethylene Gasket Seals	2	7430-00
Neoprene Gasket Seals	2	7433-00
Polypropylene Replacement Needles	12	7436-00
Polypropylene Rack Shelf Clips	12	7438-00
Vacuum Gauge / Controller Assembly	1 assembly	7437-00
Plugs for Lid	30	7327-00
Polypropylene Stopcocks	10	7241-00
Polypropylene Luer Connector Female	12	2120-02
Polypropylene Luer Connector Male	12	2121-20
Stainless steel needles	12	7323-00
Extraction Column Processor replacement parts for BAKER s	pe-12G Column Processor	
Inert Flow Control Valves	12	7425-00
Inert Flow Control Valves	150	7425-01
Luer PTFE Connectors (one-piece)	12	4586
Luer Locks (taps), Stainless Steel	12	4505
PTFE Lined Luer Locks (taps), Stainless Steel	12	4514



BAKER spe-24G Column Processor (7208-00)

Parts and Accessories

The BAKER spe-24G Column Processor complete includes: 1 glass vacuum chamber, 1 blue-coloured nylon cover with luer fitting connectors and 1 polyethylene gasket, 24 individual flow control polypropylene stopcocks, 24 stainless steel needles, 1 sample collection rack with 3 support posts, 3 height adjustable shelves, 9 shelf support clips, 1 vacuum gauge/polypropylene controller assembly, 1 instruction sheet, 1 BAKERBOND spe Application Notes Manual, 1 Solid Phase Extraction for Sample Preparation Manual.

Description	Quantity per box	Product Number
Extraction Column Processor replacement parts for BAKER	spe-24G Column Processor	
Glass Vacuum Chamber	1	7423-00
Nylon (blue-coloured) Lid (with luer connectors)	1	7426-00
Polypropylene Rack Set	1	7429-00
Polyethylene Gasket Seals	2	7432-00
Neoprene Gasket Seals	2	7435-00
Polypropylene Replacement Needles	12	7436-00
Polypropylene Rack Shelf Clips	12	7438-00
Vacuum Gauge / Controller Assembly	1 assembly	7437-00
Plugs for Lid	30	7327-00
Polypropylene Stopcocks	10	7241-00
Polypropylene Luer Connector Female	12	2120-02
Polypropylene Luer Connector Male	12	2121-20
Stainless steel needles	12	7323-00
Accessories for BAKER spe-24G Column Processor		
Inert Flow Control Valves	12	7425-00
Inert Flow Control Valves	150	7425-01
Luer PTFE Connectors (one-piece)	12	4586
Luer Locks (taps), Stainless Steel	12	4505
PTFE Lined Luer Locks (taps), Stainless Steel	12	4514

BAKERBOND Solid Phase Extraction Automation

High productivity and reproducibility in solid phase extraction as a sample preparation technique can be achieved using automation. As a manufacturer of standard SPE columns, extraction disks as well as 96-well plates, J.T.Baker offers you a possibility to use our products with automated processors.

With the sorbents we offer, you will have lot-to-lot flow reproducibility and low backpressures that will enable you to use automation to achieve higher productivity and reliable results.



- Twice the capacity of a typical vacuum processor
- Exceptional precision, control, and reliability
- Perfect for processing 1, 3 and 6 ml columns in batches of 1–48 samples
- Distributes N₂ or compressed air to displace liquid from SPE columns
- Easily enables a high level of control and precision



Speedisk 96 Positive Pressure Processor (8129-00)

- The Speedisk 96 Pressure Processor provides the means to batch process up to 96 solid phase extraction columns assembled to conform to the 96-well micro-plate footprint
- Higher pressure differentials are achievable to overcome resistance to flow caused by high viscosity
- The processor produces the same flow rate in each of 96 columns arrayed in the 8 x 12 micro-plate pattern

Description	Quantity per box	Product Number
Accessories for <i>Speedisk</i> 48 Positive Pressure Processor ¹		
Rack for 1 ml SPE Columns	1	8122-01
Rack for 3 ml SPE Columns	1	8123-01
Rack for 6 ml SPE Columns	1	8124-01
Collection Tube Rack, 12 x 75 mm Tubes	1	8119-01
Collection Tube Rack, 13 x 100 mm Tubes	1	8120-01
Collection Tube Rack, 16 x 100 mm Tubes	1	8121-01
Collection Vial Rack, 12 x 32 mm Auto-Sampler	1	8125-01
Waste Bin	1	8126-01
48-Column Sealing Gasket	1	8127-01
Gas Supply Adapter ²	1	8128-01

¹ Complete—Includes 3 ml SPE column rack, collection tube rack for 12 x 75 mm tubes, waste bin, gas supply adapter and tubing, 48-sealing gasket

 $^{^2}$ 6 ft. length of 1 D4" tubing and 1 D8"-1 D4" adapter fittings.

Description	Quantity per box	Product Number
Accessories for <i>Speedisk</i> 96 ¹ Positive Pressure Prod	cessor	
10 ml x 24 Well Collection Tray	1	8197-24
1 ml x 96 Well Collection Tray	1	8188-96
2 ml x 96 Well Collection Tray	1	8131-96
Speedisk 96 Column Holder	1	8150-00
96-Column Sealing Gasket	1	8130-01
Gas Supply Adapter ²	1	8128-01

¹ Complete—includes 1 ml x 96 collection tray, 2 ml x 96 collection tray, 96-column sealing gasket, Speedisk 96 column holder, gas supply adapter and tubing.

 $^{^2}$ 6 ft. length of 1 D4" tubing and 1 D8"-1 D4" adapter fittings

Speedisk Extraction Stations

Whatever your space and sample loading requirements, we have a vacuum extraction disk processor to meet your needs. Choose from models such as our:

- Speedisk Single Extraction Station
- Expanded Extraction Station used in reservoir, inverted, or remote sample feed modes





Speedisk Single Extraction Station (8093-01)

The Speedisk Single Extraction Station includes a single-port vacuum stand and accessories to support the extraction of analytes using BAKERBOND Speedisk laminar extraction disks. Supports any sample loading technique.

Speedisk Expanded Extraction Station (8095-06)

The Speedisk Expanded Extraction Station includes a six-port vacuum manifold and the accessories needed to support the extraction of analyte by BAKERBOND Speedisk laminar extraction disks.

The manifold has a rectangular footprint and inter-port spacing to accommodate six, side-by-side, 1 liter sample reservoirs. Each vacuum port has an individual open/close valve.

Description	Quantity per box	Product Numbe
Speedisk Extraction Stations		
Expanded extraction station		
Six-port processing system (six side-by-side 1L sample reservoirs)	1	8095-06
Single extraction station single-port processing system	1	8093-01
Speedisk Adapters		
Remote sample adapter		
For transfer of sample from remote container to <i>Speedisk</i> Disk	6	8099-06
Flask Adapter		
Single port with #8 stopper.		
Connection of disks to vacuum flask and many other vacuum manifolds		
(such as Nalge and Contes). Accepts disks or collection chamber	1	8070-01
Adapter ring		
For 40-35 tapered outer joint. Accepts disks or collection chamber	6	8100-06
Speedisk Reservoirs		
185 ml Reservoir		
Holds inverted 1L reservoir or 185 ml sample	6	8097-06
1L Glass reservoir		
1L sample reservoir, fits into a <i>Speedisk</i> extraction disk	1	8098-01
Collection chamber (includes sample vial)	2	8096-02
Collection vials (Sample vials)	100	8102-01
Sample tray		
Holds up to four 1L bottle at a tilt to ensure complete sample uptake		
by remote sample adapter suction tube	1	8101-01
70mm/ Mason Jar Adapter		
Enables inverted feed directly to extraction disk from sample jar	4	8102-04

BAKERBOND Columns and Media

BAKERBOND Columns and Media are designed to deliver optimum performance,


reproducibility, and easy scale-up.

- BAKERBOND bonded phases are synthesised using trifunctional silane synthesis chemistry
- The resulting bonded phases provide increased resistance to hydrolysis, greater stability at extremes of pH, and reduced silanol interactions for higher recoveries and less "tailing"
- Trifunctional silane chemistry provides a more consistent ligand density for increased resolution
- These high-quality spherical and irregular "polymeric" particles increase column lifetime versus monofunctional-bonded silicas
- BAKERBOND bonded phases with a given functionality are synthesised using the same chemistry regardless of particle size. Purification scale-up from analytical to process scale is straightforward
- BAKERBOND bonded phases are available as reverse, normal and ion exchange phases

BAKERBOND PAH-16 PLUS HPLC Column

The separation of PAHs is greatly influenced by the type of synthesis that is used for the preparation of C_{18} . In reverse phase chromatography two types of silane modification can be distinguished that resulted in either monomeric or polymeric bond linkages. The polymeric phases provided significant enhanced

selectivity characteristics, particularly for isomeric PAHs. BAKERBOND PAH-16 Plus column contains a wide pore C_{18} bonded phase synthesised under polymer-forming conditions. Specially designed for the separation of the sixteen PAH's as defined by the U.S. Environmental Protection Agency (EPA).



SRM	M 1647d, diluted in Acetonitrile/Water (40/60v/v)	μg/ml
1	Naphthalene	4.03
2	Acenaphthylene	3.10
3	Acenaphthene	4.15
4	Fluorene	0.95
5	Phenanthene	0.68
6	Anthracene	0.16
7	Fluoranthene	1.53
8	Pyrene	1.69
9	Benz[a]anthracene	0.82
10	Chrysene	0.73
11	Benzo[b]fluoranthene	0.83
12	Benzo[k]fluoranthene	0.94
13	Benzo[a]pyrene	0.98
14	Dibenz[a,h]anthracene	0.71
15	Benzo[qhi]perylene	0.74
16	Indeno[1,2,3-cd]pyrene	0.86

BAKERBOND PAH 16 Plus, $5\mu m$, $250 \times 3.0 \ mm$ + guard

(7504-00)


Mobile Phase: A=Water
B=Acetonitrile
Gradient: 50% B for 5 min
50 - 100% B in 15 min
100% B for 25 min

Flow Rate: 0.5 mL/min
Temperature: 25 °C
Detection: 254 nm

BAKERBOND BDC C₁₈ HPLC Column

BAKERBOND BDC C_{18} receives a special base deactivation treatment during silica material manufacture to reduce silanophilic activity. BAKERBOND BDC C_{18} is endcapped and suitable for chromatography of basic compounds. BAKERBOND BDC C_{18} phases are ideal for the determination of pharmaceutical active components.

- Base Deactivated Column
- Sharp peaks and improved symmetry for basic compounds
- Value and performance from a globally referenced column
- Long column lifetime with excellent columnto-column and batch-to-batch reproducibility
- Excellent peak symmetry



BAKERBOND BDC C_{18} , $5\mu m$, $150 \times 4.6 mm$ (7441-01)

Gradient: 60% Methanol 40% 0.05 M Phosphate buffer


Flow Rate: 1.0 m L/min Detection: 254 nm

Specifications BAKERBOND BDC C	18:
Available Particle Sizes (µm):	3,5
Pore Diameter (Å):	130
Specific Surface (m2/g):	170
Carbon (%):	11.0
Endcapped:	Yes
USP Listing:	L1

BAKERBOND ENV HPLC Column

BAKERBOND ENV is an alkyl-bonded phase, tailored for the analysis of polar pollutants commonly found in industrial and agricultural residues. The bonded phase has a uniform alkyl coverage providing a high shielding of the silica surface and a high degree of reproducibility and repeatable analysis, as well as low batch-to-batch variation.

Triazine Herbicides



BAKERBOND ENV, $5\mu m$, $150 \times 4.6 mm$ (7446-01) Gradient: A: 65% 0.1M Sodium acetate, pH 6;

B: 35% Acetonitrile

Flow Rate: 1.5 m L/min Detection: 220 nm

Specifications BAKERBOND ENV	
Available Particle Sizes (µm):	5
Pore Diameter (Å):	120
Specific Surface (m2/g):	170
Carbon (%):	13.0
Endcapped:	Yes
Pore Volume (cc/g):	0.65

BAKERBOND Chiral Columns

BAKERBOND Chiral Phase HPLC columns offer economic alternatives for the separation of racemic mixtures, including isomers of pharmaceuticals and other biologically active compounds.



Biochromatography Products and Media

Biomolecule Purification

Biochromatography products and media are specifically designed for purification and analysis of proteins, peptides and polynucleotides.

Combining chromatographic properties of wide-pore silica with a proprietary surface chemistry has created a versatile family of products for a wide range of large molecule applications. Scale-up from analytical to process scale is facilitated by providing the same surface chemistry on all three-particle sizes.

The BAKEBROND Biochromatography product line.

The BAKERBOND Biochromatography product line consists of the following types of media:

- Reverse Phase
- Ion Exchange and Antibody Purification
- Hydrophobic Interaction
- Affinity Chromatography

BAKERBOND ABx (antibody exchanger) is a unique mixed-mode ion exchanger designed specifically for the purification and analysis of antibodies. The "affinity-like" properties of ABx allow it to preferentially bind antibodies yielding purities of 80–99% in a single step, binding all antibody classes and subclasses. Since ABx does not bind most albumins, transferrins, proteases or phenol red, most of the column capacity is used for binding antibody. This provides capacities of up to 50 mg/ml.

BAKERBOND CARBOXY-SULFON directly binds/ captures, concentrates and purifies neutral and basic target proteins from neat cell culture fluids without a pre-processing step in the purification of proteins from cell culture fluids and fermentation broth (buffer exchange, dialysis, dilution). BAKERBOND CARBOXY-SULFON is a weak-strong cation exchanger providing unique binding selectivity as well as extremely high ligand densities which results in the tight binding of proteins.





	Extraction Column Size (ml)	Quantity per Box	Sorbent Weight (mg)	Product Number	
ADSORPTION					
Activated Spherical Carbon					
	6	30	500	7575-06	
	6	30	1000	7575-07	
Alumina, Neutral (Al ₂ O ₃)					
	3	50	500	7214-03	
	3	50	1000	7214-04	
	6	30	1000	7214-07	
Jumbo Pack	6	250	1000	7214-27	
Florisil (Mg ₂ SiO ₃)					
V 92 37	3	50	500	7213-03	
	6	30	1000	7213-07	
	6	30	2000	7213-08	
Jumbo Pack	3	250	500	7213-23	
Jumbo Pack	6	250	1000	7213-27	
Glass	8	30	500	7420-06	
Glass	8	30	1000	7420-07	
Silica Gel (SiOH)					
Sinta dei (Sioti)	1	100	50	7086-00	
	1	100	100	7086-01	
	3	50	200	7086-02	
	3	50	500	7086-03	
	6	30	500	7086-06	
	6	30	1000	7086-07	
	6	30	2000	7086-08	
Jumbo Pack	3	400	200	7086-22	
Jumbo Pack	3	400	500	7086-23	
Jumbo Pack	6	250	500	7086-26	
Jumbo Pack	6	250	1000	7086-28	
	15	20	2000	7086-40	
	45	20	5000	7086-41	
	70	10	10000	7086-42	
Glass	1	100	100	7337-01	
Glass	8	30	1000	7337-07	
Glass	8	30	2000	7337-08	
ANION EXCHANGE					
Diamino (NH/NH ₂)	3	50	500	7089-03	
	J	30	230	, 555 66	
Quaternary Amine (N°)		4	-	700	
	1	100	50	7091-00	
	1	100	100	7091-01	
	3	50	500	7091-03	



	Extraction Column Size (ml)	Quantity per Box	Sorbent Weight (mg)	Product Number	
BIOCHROMATOGRAPHY					
Sephadex G-25	6	30	500	7219-07	
WP Butyl (C ₄)	6	30	500	7216-06	
WP CBX (COOH)	6	30	500	7217-06	
WP HI-Propyl (C ₃)	6	30	500	7238-06	
WP PEI (NH)	1	100	150	7218-02	
CATION EXCHANGE					
Aromatic Sulfonic Acid (C _c H ₅ -SO ₃ H)	1 1 3	100 100 50	50 100 500	7090-00 7090-01 7090-03	
Jumbo Pack	6	30 400	1000 500	7090-07 7090-29	
Carboxylic Acid (COOH) Propyl Sulfonic Acid (CH3CH2CH2SO3H)	1 1 3	100 100 50	50 100 500	7211-00 7211-01 7211-03	
METHOD DEVELOPMENT	3	50	500	7155-03	
Method development kit, 5 each:(octadecyl, octyl, phenyl, si gel, cyano, carboxylic acid, amino, diol, diamino, quaternary aromatic sulfonic acid)		55	500	7096-00	
protein sorbent selection kit, 3 each (CBX, PEI, HI-propyl, but sephadex G-25)		15	500	7239-09	
NORMAL PHASE					
Amino (NH ₂)	1 1 3 3 6	100 100 50 50 30	50 100 200 500 2000	7088-00 7088-01 7088-02 7088-03 7088-09	
wide mouth	19	50	500	7088-13	



	Extraction Column Size (ml)	Quantity per Box	Sorbent Weight (mg)	Product Number	
NORMAL PHASE					
Cyano (CN)	1 1 3 3 6	100 100 50 50 30	50 100 100 500 1000	7021-00 7021-01 7021-02 7021-03 7021-07	
Diol (COHCOH)	1 1 3 6	100 100 50 30	50 100 500 500	7094-00 7094-01 7094-03 7094-06	
REVERSED PHASE					
Cyclohexyl(C _c H ₁₁) Glass Glass	1 1 3 8 8	100 100 50 30 30	50 100 500 500 1000	7212-00 7212-01 7212-03 7419-06 7419-07	
Ethyl (C ₂)	1 1 3 3	100 100 50 50	50 100 200 500	7273-00 7273-01 7273-02 7273-03	
Wide-Mouth Wide-Mouth Jumbo Pack Jumbo Pack, also available as ear shaped: PN 7020-22T Jumbo Pack Jimbo Pack Wide-Mouth	1 1 3 3 6 6 6 19 19 1 1 3 6 6 19 15 45 70 1 3	100 100 50 50 30 30 30 50 50 400 1000 400 250 250 250 20 10 100 50	50 100 200 500 500 1000 2000 100 500 100 500 500 2000 5000 10000 1000 1	7020-00 7020-01 7020-02 7020-03 7020-06 7020-07 7020-08 7020-11 7020-13 7020-21 7020-22 7020-23 7020-26 7020-27 7020-33 7020-40 7020-41 7020-42 7334-01 7334-03	
Glass Glass Glass Glass	3 8 8 8	50 30 30 30	1000 500 1000 2000	7334-04 7334-06 7334-07 7334-08	



	Extraction Column Size (ml)	Quantity per Box	Sorbent Weight (mg)	Product Number	
REVERSED PHASE					
Octadecyl (C ₁₈), Light Load					
Also available as ear shaped: PN 7189-01T	1	100	100	7189-01	
	3	50	200	7189-02	
	3	50	500	7189-03	
	3	50	1000	7189-04	
	6	30	500	7189-06	
Octadecyl (C ₁₈), Polar Plus					
	1	100	50	7466-00	
	1	100	100	7466-01	
Also available as ear shaped: PN 7466-03T	3	50	500	7466-03	
Also available as ear shaped: PN 7466-04T	3	50	1000	7466-04	
	6	30	500	7466-06	
	6	30	1000	7466-07	
	6	30	2000	7466-08	
Also available as ear shaped: PN 7466-10T	6	30	3000	7466-10	
	1	100	200	7466-12	
Jumbo Pack, also available as ear shaped: PN 7466-22T	3	400	200	7466-22	
Jumbo Pack	6	250	1000	7466-27	
Octyl (C ₈)					
	1	100	50	7087-00	
	1	100	100	7087-01	
	3	50	200	7087-02	
	3	50	500	7087-03	
	6	30	500	7087-06	
Jumbo Pack	6	250	500	7087-26	
Glass	8	30	500	7418-06	
Glass	8	30	1000	7418-07	
Phenyl(C ₆ H ₅)					
	1	100	50	7095-00	
	1	100	100	7095-01	
	3	50	500	7095-03	
	6	30	500	7095-06	
	6	30	2500	7095-09	
SDB-1 Polymer Phase					
	1	100	50	7519-00	
	3	50	100	7519-01	
Also available as ear shaped: PN 7519-02T	3	50	200	7519-02	
	6	30	200	7519-05	
Jumbo Pack, also available as ear shaped: PN 7519-22T	3	400	200	7519-22	
Jumbo Pack	6	250	200	7519-25	
Glass	3	50	100	7609-01	
Glass	3	50	200	7609-02	
SDB-2 Polymer Phase					
	3	50	50	7523-00	
Also available as ear shaped: PN 7523-02T	3	50	200	7523-02	
Also available as ear shaped: PN 7523-05T	6	30	200	7523-05	
Jumbo Pack, also available as ear shaped: PN 7523-22T	3	400	200	7523-22	



	Extraction Column Size (ml)	Quantity per Box	Sorbent Weight (mg)	Product Number	
SPECIAL APPLICATIONS					
Clean-up column of extracts of Hydrocarbons from Water. Suitable for ISO-9377-2					
Contains 0.2 g Anhydrous Sodium Sulfate (top) and 0.2 g Florisil, separated by a PTFE frit	3	50	400	7495-04	
Contains 2 g Anhydrous Sodium Sulfate and 2 g Florisil, separated by a PTFE frit	8	15	4000	7495-18	
Extraction of Pesticides from water Contains 500 mg $\rm C_{18}$ (top) and 200 mg SDB-1, separated by a PE frit. Also available as ear shaped: PN 7650-07T	6	30	700	7650-07	
Extraction of Polar Pesticides from water Contains 250 mg $\rm C_{18}$ Polar Plus (top) and 100 mg SDB-1, separated by a PE frit	6	30	350	7704-06	
Narc-1 (Δ°-carboxy THC)	3	50	500	7221-03	
Narc-2 (cocaine, benzoylecgonine)	3 6 6	50 30 30	125 250 500	7225-04 7225-05 7225-06	
Jumbo Pack	3	400	125	7225-00	
PAH Aqua Contains 200 mg Amino (top) and 500 mg C_{18} separated by a PE frit Contains 500 mg Amino (top) and 1000 mg C_{18} , separated by a PE frit	3 6	50 30	700 1500	7490-07 7490-08	
PAH Soil Contains 500 mg Cyano (top) and 1000 mg Silica gel, separated by a PE frit	6	30	1500	7518-08	
PCB-A Contains 500 mg Sulfuric Acid treated Silica gel (top) and 500 mg Aromatic Sulfonic Acid, separated by a PE frit Contains 1000 mg Sulfuric Acid treated Silica gel (top) and 700 mg	3	50	1000	7511-04	
Aromatic Sulfonic Acid, separated by a PE frit	3	50	1700	7511-06	
PCB-N Contains 500 mg Aromatic Sulfonic Acid (top) and 500 mg Silica gel, separated by a PE frit	3	50	1000	7524-04	
spe-500 (EPA approved for organochloropesticides)	6	30	500	7222-06	



BAKERBOND spe Glass Columns with PTFE Frits

	Extraction Column Size (ml)	Quantity per Box	Sorbent Weight (mg)	Product Number	
EMPTY spe GLASS COLUMNS WITH PTFE FRITS					
spe GLASS COLUMNS WITH PTFE Frits	1 3 8	100 50 30		7328-01 7328-03 7328-06	
PTFE FRITS FOR spe GLASS COLUMNS					
PTFE Frits for 1 ml spe Glass Columns PTFE Frits for 3 ml spe Glass Columns PTFE Frits for 8 ml spe Glass Columns		250 250 250		7329-01 7329-03 7329-06	



BAKERBOND speColumn Processors Parts and Accessories

Product Number

ACCESSORIES FOR BAKER spe-12G COLUMN PROCESSOR (PTFE design)

PTFE Lined Luer Locks (taps), Stainless Steel (pkg. of 12)	4514	
Drying Top, polyamide for drying purpose (columns) or evaporation purpose (eluates)	4581	
Luer PTFE Stopcocks (pkg of 12)	7514-00	
Vacuum Gauge/PTFE Controller Assembly	7515-00	
Autosampler Plate	7516-01	
Plugs for lid that fits into Luer PTFE Connectors (PN 4586) (pkg. of 30)	7517-00	

EXTRACTION COLUMN PROCESSOR REPLACEMENT PARTS FOR BAKER spe-12G COLUMN PROCESSOR

for both designs of BAKER spe-12G Column Processor		
Glass vacuum chamber, suitable for PN 7018-00 BAKER spe-12G Column Processor	7421-00	
Nylon (blue-colored) Lid (with luer connectors), suitable for PN 7018-00 BAKER spe-12G Column Processor	7424-00	
Polypropylene Rack Set, suitable for PN 7018-00 BAKER spe-12G Column Processor	7427-00	
Polyethylene Gasket Seals (pkg of 2), suitable for both designs of BAKER spe-12G Column Processors PN 7018-00 and PN		
7018-94)	7430-00	
Neoprene Gasket Seals (pkg of 2), suitable for both designs of BAKER spe-12G Column Processor PN 7018-00 and PN 7018-94	7433-00	
Glass Vacuum Chamber, suitable for PN 7018-94 BAKER spe-12G Column Processor (PTFE Design)	7512-00	
Polyamide Lid (with 12 PTFE luer connections), suitable for PN 7018-94 BAKER spe-12G Column Processor (PTFE Design)	7513-00	
PTFE Rack Set, suitable for PN 7018-94 BAKER spe-12G Column Processor (PTFE Design)	7516-00	

EXTRACTION COLUMN PROCESSOR REPLACEMENTS PARTS FOR BAKER spe-24G PROCESSOR

7423-00
7426-00
7429-00
7432-00
7435-00

EXTRACTION COLUMN PROCESSORS

BAKER spe-12G Column Processor	7018-00	
BAKER spe-12G Column Processor (PTFE Design))	7018-94	
BAKER spe-24G Column Processor	7208-00	

GENERAL ACCESSORIES FOR BAKER spe-12G and BAKER-24G

General SPE Accessories		
Polypropylene Luer Connector Female (pkg. of 12)	2120-02	
Polypropylene Luer Connector Male (pkg. of 12)	2121-20	
Luer Locks (taps), Stainless Steel (pkg. of 12)	4505	
10 ml Volumetric Flasks (pkg. of 10)	4515	
BAKERBOND spe Application Notes Manual	4520	
Luer PTFE Connectors (one-piece)	4586	
Polypropylene Stopcocks (pkg. of 10)	7241-00	
Stainless Steel Replacement Needles (pkg. of 12)	7323-00	
Plugs for Female Luer (pkg. of 30)	7327-00	
Inert Flow Control valves (pkg. of 12)	7425-00	
Inert Flow Control Valves (pkg. of 150)	7425-01	
Polypropylene Replacement Needles (pkg. of 12)	7436-00	
Vacuum Gauge / controller assembly	7437-00	



BAKERBOND speColumn Processors Parts and Accessories

Product Number

GENERAL SPE ACCESSORIES	
Reservoirs For Use With 1, 3 and 6 ml spe Columns, 15 ml (pkg. of 10)	7119-01
Reservoirs For Use With 3 and 6 ml Columns, 75 ml (pkg. of 10)	7120-03
Filtration Columns, Packed With Dual 20 μm Frits, 1 ml(pkg.of 100)	7121-01
Filtration spe Columns, Packed With Dual 20 μm Frits, 3 ml (pkg. of 50), also available as ear shaped: PN 7121-0	7121-03
Filtration spe Columns, Packed With Dual 20 μm Frits, 3 ml (pkg. of 100)	7121-04
Frits for 3ml Polypropylene spe Columns (pkg. of 50)	7121-05
Filtration spe Columns, Packed With Dual 20 μm Frits, 6 ml (pkg. of 30)	7121-06
Filtration spe Columns, Packed With 20 μm Frits, 6 ml (pkg.of 100)	7121-08
Adaptors	
Adaptor (white-colored) PTFE for Glass spe Columns (pkg. of 12)	4528
Adaptor (blue-colored) for Attaching Reservoir or Luer Tip (pkg. of 12)	7122-00
3 ml spe Column adaptor for use with Gilson Aspec (pkg. of 250)	8137-01
3 ml spe Column adaptor for use with Gilson Aspec (pkg. of 1000)	8137-10
1 ml spe Column adaptor for use with Gilson Aspec (pkg. of 250)	8138-01
1 ml spe Column adaptor for use with Gilson Aspec (pkg. of 1000)	8138-10



	Extraction Column Size (ml)	Quantity per Box	Sorbent Weight (mg)	Product Number	
ADSORPTION					
Silica	1 1 1 3 3 3 3 6 6	100 100 100 50 50 50 50 30	10 20 35 35 50 100 50	8163-00 8163-01 8163-02 8163-03 8163-04 8163-06 8163-07 8163-08	
	6 1 (rimless)	30 96	200 20	8163-09 8163-11	
ANION EXCHANGE					
Amino (NH ₂) Quarternary Amine (N*)	1 1 3 3 6 6 6 1 (rimless)	100 100 100 50 50 30 30 30 96 100 100 50 50 50 30 30 30 30	10 20 35 35 35 100 50 100 200 20 10 20 35 35 50 100 50 100 200 20	8165-00 8165-01 8165-02 8165-03 8165-06 8165-07 8165-08 8165-09 8165-11 8168-00 8168-01 8168-02 8168-03 8168-04 8168-06 8168-07 8168-09 8168-09 8168-11	
CATION EXCHANGE					
Aromatic Sulfonic Acid (C ₆ H ₅ -SO ₃ H)	1 1 1 3 3 3 6 6 6 6 6 1 (rimless)	100 100 100 50 50 50 30 30 30 96	10 20 35 35 50 100 50 100 200 20	8170-00 8170-01 8170-02 8170-03 8170-04 8170-06 8170-07 8170-08 8170-09 8170-11	



	Extraction Column Size (ml)	Quantity per Box	Sorbent Weight (mg)	Product Number	
CATION EXCHANGE					
Carboxylic Acid (COOH)					
	1	100	10	8172-00	
	1	100	20	8172-01	
	1	100	35	8172-02	
	3	50	35	8172-03	
	3	50	100	8172-06	
	6	30	50	8172-07	
	6	30	100	8172-08	
	6 1 (rimless)	30 96	200 20	8172-09 8172-11	
	i (iiiiiicss)	00	20	0172-11	
NORMAL PHASE					
Cyano(CN)					
	1	100	10	8159-00	
	3	50	35	8159-03	
	6	30	50	8159-07	
	6	30	200	8159-09	
Diol (COHCOH)					
	1	100	10	8167-00	
	1	100	35	8167-02	
	3	50	35	8167-03	
	6	30	50	8167-07	
	1 (rimless)	96	20	8167-11	
POLYMER COLUMN-HYDROPHOBIC					
H ₂ O-Phobic DVB					
	1	100	10	8109-00	
	1	100	20	8109-01	
	1	100	35	8109-02	
	3	50	35	8109-03	
	3	50	50	8109-04	
	3	50	100	8109-06	
	6	30	50	8109-07	
	6	30	100	8109-08	
	6	30	200	8109-09	
	1 (rimless)	96	20	8109-11	
H ₂ O-Phobic SC-DVB (SO ₃)	4	100	10	0100.00	
	1	100	10 20	8196-00 9196-01	
	1	100	20 25	8196-01	
	1 3	100 50	35 35	8196-02 8196-03	
	3	50 50	35 50	8196-03 8196-04	
	3	50 50	100	8196-04 8196-06	
	6	30	50	8196-06	
	6	30	100	8196-07	
	6	30	200	8196-08	
			200	8196-09	
	1 (rimless)	96	-7(1	8146-11	



	Extraction Column Size (ml)	Quantity per Box	Sorbent Weight (mg)	Product Number	
POLYMER COLUMN-HYDROPHOBIC					
POLYMER COLUMNS-HYDROPHILIC	1 1 1 3 3 3 6 6 6 6 1 (rimless)	100 100 100 50 50 50 30 30 30 96	10 20 35 35 50 100 50 100 200 20	8115-00 8115-01 8115-02 8115-03 8115-04 8115-06 8115-07 8115-08 8115-09 8115-11	
H ₂ o-Philic DVB					
	1 1 3 3 3 6 6 6 6 1 (rimless)	100 100 100 50 50 50 30 30 30 30	10 20 35 35 50 100 50 100 200 20	8108-00 8108-01 8108-02 8108-03 8108-04 8108-06 8108-07 8108-08 8108-09 8108-11	
H ₂ O-Philic SA-DVB	1 1 1 3 3 3 6 6 6 6 1 (rimless)	100 100 100 50 50 50 30 30 30 96	10 20 35 35 50 100 50 100 200 20	8113-00 8113-01 8113-02 8113-03 8113-04 8113-06 8113-07 8113-08 8113-09 8113-11	
H ₂ O-Philic SC-DVB (SO ₃)	1 1 1 3 3 3 6 6 6 6 1 (rimless)	100 100 100 50 50 50 30 30 30 30	10 20 35 35 50 100 50 100 200 20	8111-00 8111-01 8111-02 8111-03 8111-04 8111-06 8111-07 8111-08 8111-09	



	Extraction Column Size (ml)	Quantity per Box	Sorbent Weight (mg)	Product Number	
REVERSED PHASE					
Butyl (C ₄)					
	1	100	10	8184-00	
	1	100	20	8184-01	
	1	100	35	8184-02	
	3	50	35	8184-03	
	3	50	50	8184-04	
	3	50	100	8184-06	
	6	30	50	8184-07	
	6	30	100	8184-08	
	6	30	200	8184-09	
	1 (rimless)	96	20	8184-11	
Ethyl (C ₂)		4	-		
	1	100	35	8157-02	
	3	50	35	8157-03	
	6	30	50	8157-07	
	6	30	100	8157-08	
	6	30	200	8157-09	
Octadecyl (C ₁₈)					
	1	100	20	7606-01	
	1	100	35	7606-02	
	3	50	35	7606-03	
	3	50	50	7606-04	
	3	50	100	7606-06	
	6	30	50	7606-07	
	6	30	100	7606-08	
	6	30	200	7606-09	
	1 (rimless)	96	20	7606-11	
Octadecyl Light Load (C ₁₈)					
	1	100	10	8151-00	
	1	100	20	8151-01	
	1	100	35	8151-02	
	3	50	35 	8151-03	
	3	50	50	8151-04	
	3	50	100	8151-06	
	6	30	50	8151-07	
	6	30	100	8151-08	
	6	30	200	8151-09	
	1 (rimless)	96	20	8151-11	
Octyl (C _s)					
	1	100	10	8154-00	
	1	100	20	8154-01	
	1	100	35	8154-02	
	3	50	35	8154-03	
	3	50	50	8154-04	
	3	50	100	8154-06	
	6	30	50	8154-07	
	6	30	100	8154-08	
	6	30	200	8154-09	
	1 (rimless)	96	20	8154-11	



	Extraction Column Size (ml)	Quantity per Box	Sorbent Weight (mg)	Product Number	
REVERSED PHASE					
Phenyl (C ₆ H ₅)	1 1	100 100	10 35	8160-00 8160-02	
	3 3 6	50 50 30	50 100 50	8160-04 8160-06 8160-07	
PolarPlus Octadecyl (C ₁₈)	6	30 100	200	8160-09 8153-00	
	1 1 3 3	100 100 50	20 35 35 50	8153-01 8153-02 8153-03 8153-04	
	3 6 6 6 1 (rimless)	50 50 30 30 30 96	100 50 100 200	8153-04 8153-06 8153-07 8153-08 8153-09 8153-11	
PolarPlus Octyl (C ₈)	1 1 1 3	100 100 100 50	10 20 35 35	8156-00 8156-01 8156-02 8156-03	
	3 6 1 (rimless)	50 30 96	50 50 20	8156-04 8156-07 8156-11	
SPECIAL APPLICATIONS					
Narc-1 (Δ° - Carboxy THC)	1 3 3 6 6	100 50 50 30 30	10 50 100 50 100	8174-00 8174-04 8174-06 8174-07 8174-08	
Narc-2 (Basic Drugs, Cocaine, Benzoylecgonine)	1 1 1 3 3 3 6 6 6 6 1 (rimless)	100 100 100 50 50 50 30 30 30 30 96	10 20 35 35 50 100 50 100 200 20	8175-00 8175-01 8175-02 8175-03 8175-04 8175-06 8175-07 8175-08 8175-11	



SPEEDISK 96 POLYMER COLUMNS (incl. holder)			
H ₂ O-Philic SA-DVB 1 (rimless) 96	6 20	8113-31	
H ₂ O-Phobic DVB 1 (rimless) 96	6 20	8109-31	
H ₂ O-Phobic SC-DVB 1 (rimless) 96	6 20	8196-31	
H ₂ O-Phobic WA-DVB 1 (rimless) 96	6 20	8115-31	
H ₂ O-Philic DVB 1 (rimless) 96	6 20	8108-31	
H ₂ O-Philic SC-DVB 1 (rimless) 96	6 20	8111-31	
SPEEDISK 96 SILICA COLUMNS (incl. holder)			
Amino (NH ₂) 1 (rimless) 96	6 20	8165-31	
Aromatic Sulfonic Acid (C ₆ H ₅ -SO ₃ H) 1 (rimless) 96	6 20	8170-31	
Butyl (C ₄) 1 (rimless) 96	6 20	8184-31	
Butyl WP 1 (rimless) 96	6 20	8180-31	
Carboxylic Acid 1 (rimless) 96	6 20	8172-31	
CBx WP	6 20	8176-31	
Ethyl (C ₂) 1 (rimless) 96	6 20	8157-31	
HI Propyl (C ₃) 1 (rimless) 96		8182-31	
Narc-2 1 (rimless) 96		8175-31	
Octadecyl (C ₁₈) 1 (rimless) 96		7606-31	
Octadecyl LightLoad (C ₁₈) 1 (rimless) 96		8151-31	
Octyl (C _s) 1 (rimless) 96 1 (rimless)		8154-31	



	Extraction Column Size (ml)	Quantity per Box	Sorbent Weight (mg)	Product Number	
SPEEDISK 96 SILICA COLUMNS (incl. holder)					
PEI WP	1 (rimless)	96	20	8178-31	
Phenyl (C₀H₅)	1 (rimless)	96	20	8160-31	
PolarPlus Octadecyl (C ₁₈)	1 (rimless)	96	20	8153-31	
PolarPlus Octyl (C _s)	1 (rimless)	96	20	8156-31	
Quarternary Amine (N ⁻)	1 (rimless)	96	20	8168-31	
Silica	1 (rimless)	96	20	8163-31	



BAKERBOND Speedisk extraction products

	Quantity per box	Product Number	
BAKERBOND SPEEDISKS EXTRACTION DISKS			
C ₁₈ (octadecyl) 50 mm disks for normal water-samples 50 mm disks for water samples, High capacity	20 20	8055-06 8055-07	
C ₁₈ Polar Plus 50 mm disks for water-samples containing slightly polar to non-polar analytes	20	8061-06	
C ₁₈ XF (Extra filter) 50 mm disks for crude and dirty samples.	20	8056-06	
C ₈ (octyl) 50 mm disks for diquat/paraquat.	20	8057-06	
H ₂ O Philic DVB (DiVinylBenzene) 50 mm disks for non-polar to polar analytes 50 mm disks, High capacity	20 20	8072-06 8072-07	
H ₂ O-Phobic DVB (DiVinylBenzene) 50 mm disks for non-polar to moderately polar analytes 50 mm disks, High capacity	20 20	8068-06 8068-07	
Oil & Grease 50 mm disks for hydrocarbons/Oil & Grease.	20	8060-06	
SAX (Strong Anion Exchanger) 50 mm disks for haloacetic acids/Dalapon.	20	8058-06	
SPEEDISK ADAPTERS	'		
Adapters Connection of disks to vacuum flask and many other vacuum manifolds (such as Nalge and Contes). Adapter ring for 40-35 tapered outer joint	1 6	8070-01 8100-06	
Remote sample adapter For transfer of sample from remote container.	6	8099-06	
SPEEDISK EXTRACTION STATIONS			
Compact extraction station Six-port processing system. (three side-by-side 1L sample reservoirs).	1	8094-06	
Expanded extraction station Six-port processing system. (six side-by-side 1L sample reservoirs).	1	8095-06	
Single extraction station single-port processing system.	1	8093-01	



BAKERBOND Speedisk extraction products

	Quantity per box	Product Number	
SPEEDISK RESERVOIRS			
185 ml Reservoir Holds inverted 1L reservoir or 185 ml sample.	6	8097-06	
1L Glass reservoir 1L sample reservoir, fits into a Speedisk extraction disk.	1	8098-01	
70 mm/Mason Jar Adaptor Enables inverted geed directly to extraction disk from sample jar	4	8102-04	
Collection chamber (includes sample vial)	2	8096-02	
Collection Vials (Sample vials)	100	8102-01	
Sample tray Holds up to four 1L bottles at a tilt to ensure complete sample uptake by remote sample adapter suction tube.	1	8101-01	



BAKERBOND Speedisk Processing Systems System Accessories

	Product Number	
BAKERBOND Speedisk PROCESSING SYSTEMS		
Speedisk 48 Positive Pressure Processor	8118-00	
Speedisk 96 Positive Pressure Processor	8129-00	
GENERAL ACCESSORIES for Speedisk 48 and Speedisk 96		
Gas Supply Adapter, 6 ft. length of 1/4 tubing and 1/8 -1/4 adapter fittings	8128-01	
Speedisk 48 ACCESSORIES		
Collection Tube Rack, 12 x 75 mm Tubes Collection Tube Rack, 13 x 100 mm Tubes Collection Tube Rack, 16 x 100 mm Tubes Rack for 1 ml SPE Columns Rack for 3 ml SPE Columns Rack for 6 ml SPE Columns Collection Vial Rack, 12 x 32 mm Auto-Sampler Waste Bin 48-Column Sealing Gasket	8119-01 8120-01 8121-01 8122-01 8123-01 8124-01 8125-01 8126-01 8127-01	
Speedisk 96 ACCESSORIES		
96-Column Sealing Gasket 2 ml x 96 Well Collection Tray Speedisk 96 Column Holder 1 ml x 96 Well Collection Tray 10 ml x 24 Well Collection Tray	8130-01 8131-96 8150-00 8188-96 8197-24	



HPLC Columns

	Å	column size (mm)	micron	Product Number	
ADSORPTION HPLC COLUMNS					
Silica Gel (SiOH)	120 120	4.6 x 150 4.6 x 50	5 3	7097-01 7161-05	
NORMAL PHASE HPLC COLUMNS					
Cyano (CN)	120	4.6 x 250	5	7111-00	
REVERSED PHASE HPLC COLUMNS					
octadecyl (C ₁₈ BDC) with guard with guard with guard Replacement guard cartridges	120 120 120 120 120 120 120	4.6 x 250 4.6 x 150 4.6 x 250 4.6 x 150 4.6 x 150	5 5 5 5 3 5	7441-00 7441-01 7441-02 7441-03 7441-08 7442-05	
Octadecyl (C ₁₈)	120 120 120	4.6 x 250 4.6 x 150 4.6 x 50	5 5 3	7098-00 7098-01 7160-05	
Octyl (C ₈)	120	4.6 x 250	5	7109-00	
SPECIAL APPLICATION HPLC COLUMNS					
Bakerbond ENV Replacement column for Bakerbond ENV Bakerbond ENV with guard (4.6x10 mm) Replacement guard cartridges (2) for Bakerbond ENV	120 120 120	4.6 x 150 4.6 x 150 4.6 x 10	5 5 5	7445-01 7446-01 7448-01	
Bakerbond PAH 16 plus Bakerbond PAH 16 plus including guard column, manufit connection system and certificate of analysis Replacement guard cartridges for Bakerbond PAH 16 plus (3) Replacement column for Bakerbond PAH 16 plus including certificate of analysis	120 120 120	3 x 250 3 x 20 3 x 250	5 5 5	7504-00 7505-00 7506-00	



	Å	micron	Pkg Size	Product Number	
ADSORPTION					
Silica Gel (SiOH)	60 60 60 60 60 150	40 40 40 40 40 10	5 kg 1 kg 25 kg 500 g 125 g 10 g	7024-00 7024-01 7024-02 7024-05 7024-09 7048-01	
ION EXCHANGE					
Aromatic Sulfonic Acid (C _o H ₂ SO ₃ H)	60	40	100 g	7046-00	
Carboxylic Acid (COOH)	60	40	100 g	7044-00	
Diamino (NH/NH ₂)	60	40	100 g	7042-00	
Quaternary Amine (N°)	60	40	100 g	7043-00	
Sulfonic Acid (SO ₃ H)	60	40	100 g	7045-00	
NORMAL PHASE					
Amino (NH ₂)	60 60 150 120	40 40 10 5	100 g 1 kg 10 g 10 g	7028-00 7028-01 7052-01 7070-01	
Diol (COHCOH)	60 60	40 40	100 g 1 kg	7047-00 7047-01	
REVERSED PHASE					
Butyl (C ₄)	60 150	40 10	100 g 10 g	7037-00 7055-01	
Cyano (CN)	60 150 120	40 10 5	100 g 10 g 10 g	7027-00 7051-01 7069-01	
Cyclohexyl (C _c H ₁₁)	150	10	10 g	7057-01	
Ethyl (C ₂)	60 60	40 40	100 g 1 kg	7199-00 7199-01	



	Å	micron	Pkg Size	Product Number	
REVERSED PHASE					
Octadecyl (C ₁₈)					
	60	40	100 g	7025-00	
	60	40	1 kg	7025-01	
	150	10	10 g	7049-01	
Octadecyl (C ₁₈), High Hydrophobic					
	60	40	100 g	7243-00	
Octadecyl (C ₁₈), Polar Plus					
7 (10)	60	40	100 g	7465-00	
	60	40	1 kg	7465-01	
	60	40	20 g	7465-02	
Octyl (C ₈)					
	60	40	100 g	7026-00	
	150	10	10 g	7050-01	
	120	5	10 g	7068-01	
Phenyl(C ₆ H ₅)					
	60	40	100 g	7040-00	
	150	10	10 g	7058-01	
	120	5	10 g	7076-01	
Styrene divinyl benzene (SDB 1)					
			100 g	7530-00	
SPECIAL APPLICATIONS					
Floricii CDT andrast for along up anlungo Guitable for					
Florisil. SPE sorbent for clean-up columns. Suitable for determination of Hydrocarbon Oil Index according ISO-9377-2 and					
NEN 5733					
			100 g	7061-00	
Magnesium Sulfate Heptahydrate. Suitable for the determination of					
Hydrocarbon Oil Index according ISO-9377-2 and NEN 5377. Avoids					
formation of emulsions					
			100g-500g-1kg	0168	
Sodium Sulfate Anhydrous. SPE sorbent for clean-up columns.					
Suitable for determination of Hydrocarbon Oil Index according					
ISO-9377-2 and NEN 5733					
the state of the s			100 g	3377-00	



HPLC Columns HPLC Columns Accessories

	Å	micron	Size (mm)	Product Number	
GUARD COLUMNS (empty)					
Guard Columns Recommended for use with columns ID 10.0 mm Recommended for use with columns ID 21.2 mm Recommended for use with columns ID 50.8 mm			4.6 x 50 10.0 x 50 21.2 x 100	7128-02 7128-04 7128-08	
HYDROPHOBIC INTERACTION					
WP HI-Propyl (C ₃) STANDARD VERSA-TEN SEMI-PREP	300 300 300	5 5 15	4.6 x 250 7.75 x 100 10 x 250	7276-00 7276-06 7277-00	
ION EXCHANGE HPLC COLUMNS					
ABx (antibody exchanger) SEMI-PREP PREP	300 300	15 15	10 x 250 21.2 x 250	7272-00 7272-43	
WP CARBOXY-SULFON (weak-strong cation exchanger) STANDARD SCOUT SEMI-PREP	300 300 300	5 5 15	4.6 x 250 4.6 x 50 10 x 250	7159-00 7159-05 7274-00	
WP CBX (weak cation exchanger) STANDARD VERSA-TEN SEMI-PREP	300 300 300	5 5 15	4.6 × 250 7.75 × 100 10 × 250	7114-00 7114-06 7279-00	
WP DEAM (weak anion exchanger) VERSA-TEN SEMI-PREP	300 300	5 15	7.75 x 100 10.0 x 250	7474-06 7475-39	
WP PEI (weak anion exchanger) SEMI-PREP	300	15	10 x 250	7278-00	
WP QUAT (strong anion exchanger) SCOUT VERSA-TEN SEMI-PREP	300 300 300	5 5 15	4.6 x 50 7.75 x 100 10.0 x 250	7158-05 7158-06 7275-00	
REVERSED PHASE HPLC COLUMNS					
WP Butyl (C4) STANDARD VERSA-TEN SEMI-PREP SEMI-PREP PREP	300 300 300 300 300	5 5 5 15 15	4.6 x 250 7.75 x 100 10 x 250 10 x 250 21.2 x 250	7116-00 7116-06 7116-39 7280-00 7280-43	
WP Diphenyl ((C ₆ H ₅)) ₂ SEMI-PREP	300	15	10 x 250	7298-39	



HPLC Columns HPLC Columns Accessories

	Å	micron	Size (mm)	Product Number	
REVERSED PHASE HPLC COLUMNS					
WP Octadecyl (C ₁₈)					
STANDARD	300	5	4.6 x 250	7104-00	
STANDARD-TEN	300	5	4.6 x 100	7104-01	
SEMI PREP	300	5	10.0 x 250	7104-39	
SEMI-PREP	300	15	10.0 x 250	7297-00	
PREP	300	15	21.2 x 250	7297-43	
PREP	300	15	50.8 x 250	7297-47	
WP Octyl (C ₈)					
STANDARD-TEN	300	5	4.6 x 100	7105-01	
SEMI-PREP	300	15	10 x 250	7296-39	



	Å	micron	Pkg Size	Product Number	
AFFINITY					
Glutaraldehyde-P					
didialdenyde r	300	40	100 g	7567-00	
	300	40	10 g	7567-02	
	300	40	500 g	7567-05	
Chitaraldahuda D					
Glutaraldehyde-P	500	40	1 kg	7568-01	
	500	40	10 g	7568-02	
ION EXCHANGE					
ABx (antibody exchanger)					
	300	5	1 g	7154-01	
	300	15	100 g	7157-00	
	300	15	10 g	7157-02	
	300	15	500 g	7157-05	
	275	40	100 g	7269-00	
	275	40	10 g	7269-02	
	275	40	500 g	7269-05	
ABx (antibody exchanger) 500 Å					
	300	15	100 g	7182-00	
	300	15	10 g	7182-02	
	300	15	500 g	7182-05	
	275	40	100 g	7285-00	
	275	40	10 g	7285-02	
	275	40	500 g	7285-05	
	300	5	1 g	7291-01	
ABx Plus (antibody exchanger)					
	275	40	100 g	7254-00	
	275	40	10 g	7254-02	
	275	40	500 g	7254-05	
WP CARBOXY-SULFON (weak-strong cation exchanger)					
	300	15	10 g	7184-02	
	275	40	100 g	7252-00	
	275	40	10 g	7252-02	
	275	40	500 g	7252-05	
WP CBX (weak cation exchanger)					
	300	5	1 g	7136-01	
	300	15	100 g	7181-00	
	300	15	10 g	7181-02	
	300	15	500 g	7181-05	
	275	40	100 g	7263-00	
	275 275	40 40	10 g	7263-02 7263-05	
	2/0	40	500 g	7203-05	
WP DEAM (weak anion exchanger)	000	4-	40	7.470.00	
	300	15	10 g	7472-02	
	275	40	100 g	7473-00	
		4(1)	10 g	7473-02	
	275 275	40 40	500 g	7473-05	



	Å	mirron	Pkg Size	Product Number	
	А	micron	PKG SIZE	Nullibei	
ION EXCHANGE					
WP PEI (weak anion exchanger)					
	300	5	1 g	7134-01	
	300 300	15 15	100 g	7180-00 7180-02	
	300	15	10 g 500 g	7180-02	
	275	40	100 g	7264-00	
	275	40	1 kg	7264-01	
	275	40	10 g	7264-02	
	275	40	500 g	7264-05	
WP QUAT (strong anion exchanger)					
	300	15	10 g	7183-02	
	275	40	100 g	7251-00	
	275 275	40 40	10 g 500 g	7251-02 7251-05	
	2/5	40	500 g	7251-05	
REVERSED PHASE CHROMATOGRAPHY MEDIA					
WP Butyl (C ₄)					
111 boty. (c ₄)	300	5	1 g	7137-01	
	300	15	100 g	7179-00	
	300	15	10 g	7179-02	
	300	15	500 g	7179-05	
	275	40	100 g	7283-00	
	275 275	40 40	10 g 500 g	7283-02 7283-05	
	2/5	40	500 g	7203-03	
WP Cyano (CN)	300	5	1 ~	7129-01	
	300	5	1 g	7129-01	
WP Diphenyl ((C ₆ H ₅) ₂)			400	7400.00	
	300 300	15 15	100 g 10 g	7192-00 7192-02	
	300	15	500 g	7192-02	
		_			
WP Octadecyl (C ₁₈)	300	5	1 ~	7133-01	
	300	5	1 g 50 g	7133-06	
	300	15	100 g	7191-00	
	300	15	10 g	7191-02	
	300	15	500 g	7191-05	
	300	15 15	100 g	7207-00	
	300	15 15	10 g	7207-02	
	300 275	15 40	500 g 100 g	7207-05 7248-00	
	275	40	100 g	7248-02	
	275	40	500 g	7248-05	
WP Octyl (C ₈)					
σειμ (εξ)	300	5	1 g	7132-01	
			Ü		



	Å	micron	Pkg Size	Product Number	
WIDE-PORE SILICA					
WP Silica	1000 1000	40 40	1 kg 10 g	7315-01 7315-02	
WP Silica	500 500	40 40	1 kg 10 g	7314-01 7314-02	



Flash Chromatography Columns Column Accessories

J.T. BAKER FLASH CHROMATOGRAPHY COLUMN A	Column Size	Joint Size	Single Pass Sample Size on Silica Gel	Product Number	
Flash Chromatography Accessories Silanized Glass Wool Reservoir for use with 100 ml column Reservoir for use with 250 ml and 450 ml column Reservoir for use with 650 ml column (includes an extra clamp) Flow Controller Replacement for use with 100 ml column Flow Controller Replacement for use with 250 ml and 450 ml column Flow Controller Replacement for use with 650 ml column	50 g 500 ml 500 ml 1 l			7084-05 7092-00 7092-04 7092-06 7118-01 7118-04 7118-06	
J.T.BAKER FLASH CHROMATOGRAPHY COLUMNS					
Flash Chromatography Columns 19 mm I.D. x 35.6 cm length (includes 20 g C ₁₈ bonded phase) 30 mm I.D. x 36.8 cm length (includes 45 g silica gel) 33 mm I.D. x 52.7 cm length (includes 80 g silica gel) 37 mm I.D. x 61.0 cm length (includes 125 g silica gel)	100 ml 250 ml 450 ml 650 ml	35/20 50/30 50/30 65/40	0-350 mg 0-600 mg 0-750 mg 0-1000 mg	7022-01 7022-02 7022-04 7022-06	



Chiral Chromatography Summary

HPLC Columns

	Column Size (mm)	micron	Product Number	
PIRKLE-TYPE HPLC COLUMNS AND PACKINGS				
α-Naphtyl Urea	4.6 x 250	5	7309-00	
(R) DNBPG (Covalent)	4.6 x 250	5	7113-00	
(R) DNBPG(ionic)	4.6 x 250	5	7103-00	
PROTEIN-TYPE HPLC COLUMN				
Guard Column Guard Column Semi-Prep Semi-Prep Semi-Prep Semi-Prep Semi-Prep CHIRAL-CBH Guard Column Guard Column	4.0 x 100 3.0 x 100 3.0 x 150 2.0 x 100 2.0 x 150 4.0 x 50 4.0 x 150 3 x 10 2 x 10 10.0 x 100 10.0 x 150 3 x 10 2 x 10 4.0 x 150 3 x 10 2 x 10 10.0 x 150 3 x 10 2 x 10 4.0 x 150 3 x 10 2 x 10 4.0 x 50 3.0 x 100	5 5 5 5 5 5 5 5 5 5 5 5 5 5 5 5 5 5 5	7165-00 7165-01 7165-02 7165-03 7165-04 7165-26 7165-30 7400-00 7400-01 7404-07 7404-08 7404-37 7404-38	
CHIRAL-HSA Guard Column	3 × 10 3.0 × 100 4.0 × 150	5 5 5 5 5	7496-29 7496-30 7401-00 7401-01 7401-29 7401-30	
General Accessories Guard Column Holder Guard Column Coupler			7402-00 CON-2	



Thin Layer Chromatography Summary

BAKER Pre-Coated Hard Layer TLC Plates

	Plate Size (cm)	Quantity per Box	Product Number			
SILICA GEL PLATES 250 μm Analytical Layer						
TLC Plate With Fluorescent Indicator BAKER Si250F BAKER Si250F BAKER Si250F	5 x 20 10 x 20 20 x 20	100 50 25	7001-00 7001-03 7001-04			
TLC Plate With Fluorescent Indicator BAKER Si250F-PA BAKER Si250F (19C) BAKER Si250F-PA	20 x 20 20 x 20 20 x 20	25 25 25	7004-04 7007-04 7010-04			
TLC Plate Without Fluorescent Indicator BAKER 5i250 BAKER 5i250	10 x 20 20 x 20	50 25	7000-03 7000-04			
TLC Plate Without Fluorescent Indicator BAKER Si250-PA BAKER Si250 (19C) BAKER Si250-PA (19)	20 x 20 20 x 20 20 x 20	25 25 25	7003-04 7005-04 7009-04			
SILICA GEL PLATES 500 μm Preparative Layer						
TLC Plate With Fluorescent Indicator BAKER Si500F	20 x 20	20	7002-05			
SILICA GEL PLATES High Performance TLC, 200 μm Analytic	SILICA GEL PLATES High Performance TLC, 200 μm Analytical Layer					
TLC Plate With Fluorescent Indicator BAKER SiHPF	10 x 20	50	7011-03			
SILICA GEL PLATES - Reversed Phase, Octadecyl (C18), 200	SILICA GEL PLATES - Reversed Phase, Octadecyl (C18), 200 μm Analytical Layer					
TLC Plate With Fluorescent Indicator BAKER Si-C ₁₈ F	20 x 20	25	7013-04			
TLC Plate Without Fluorescent Indicator BAKER Si-C ₁₈	20 x 20	25	7012-04			



Trademarks

The following list of trademarks is accurate to the best of our knowledge at the time of printing. Any errors brought to our attention will be corrected in the next issue of this catalogue. Please contact the trademark owners or public resources for specific information about these trademarks.

Trademarks of Mallinckrodt Baker, Inc.*

BAKER REZI

AGILE GRANUSIC
ANHYDRONE HYDRA-POINT

BAKER narc

BAKER ALEG

BAKER ANALYZED

narc

NEUTRACIT

NEUTRASORB

BAKER BIO-ANALYZED pCMP
BAKERBOND PolarPlus
BAKERBOND ABx QUAT
BAKERBOND spe RESISORB

BAKERBOND Speedisk Right From the Start
BAKERBOND Wide-Pore SAFETAINER
BakerDRY SOLUSORB
BakerClean Speedisk

Baker-flex ULTRA RESI-ANALYZED

BAKER INSTRA-ANALYZED ULTREX
BAKER PRS-1000 VLSI

BAKER PRS-3000 Part of a pure process

CARBOXY-SULFON

Trademarks of Mallinckrodt Baker B.V. Deventer
CYCLE-TAINER

CyMet

DESICCHLORA DetectoTerge
DILUT-IT Diluid
DISKMATE Lyzerglobin
DUAL-TINT ProClean
Flowmor Rapid Stat
GelTwin UltraClear

Global Pharma UltraKitt UltraPar

The following Registered Trademarks of the companies indicated appear in this catalogue

Brij ICI Americas, Inc.
Celite Celite Corporation

CELLOSOLVE Union Carbide Corporation
COOMASSIE BASF Aktiengesellschaft Inc.
Coulter Coulter Electronics Inc.
Dowex Dow Chemical Company

Florisil Floridin Company

Freon E.I.duPont deNemours Company

Hyamine Lonza, Inc.

Responsible Care Chemical Manufacturers Association

Sephadex Pharmacia Biotech

Sysmex Toa Medical Electronics Co., Ltd.
Teflon E.I.duPont deNemours Company
Tergitol Union Carbide Corporation
Triton Union Carbide Corporation

Tween ICI Americas, Inc.

The following Trademarks of the companies indicated appear in this catalogue

CHIRAL-AGP Chrom Tech AB
CHIRAL-HAS Chrom Tech AB

Service marks of Mallinckrodt Baker, Inc.

Solv-IT

Solv-IT Center

Design marks of Mallinckrodt Baker, Inc.







Risk and Safety phrases

List of Risk phrases

D.1	5 1 1 1	D40/40	M
R1	Explosive when dry.	R42/43	May cause sensitization by inhalation and skin contact.
R2 R3	Risk of explosion by shock, friction, fire or other sources of ignition.	R43	May cause sensitization by skin contact.
R4	Extreme risk of explosion by shock, friction, fire or other sources of ignition.	R44 R45	Risk of explosion if heated under confinement.
R5	Forms very sensitive explosive metallic compounds. Heating may cause an explosion.	R46	May cause cancer.
R6	Explosive with or without contact with air.	R48	May cause heritable genetic damage. Danger of serious damage to health by prolonged exposure.
R7	May cause fire.		
R8	Contact with combustible material may cause fire.	R48/20	Harmful: danger of serious damage to health by prolonged exposure through inhalation.
R9	Explosive when mixed with combustible material.	R48/20/21	Harmful: danger of serious damage to health by prolonged exposure through
R10	Flammable.	1140/20/21	inhalation and in contact with skin.
R11	Highly flammable.	R48/20/21/22	Harmful: danger of serious damage to health by prolonged exposure through
R12	Extremely flammable.	1140/20/21/22	inhalation, in contact with skin and if swallowed.
R14	Reacts violently with water.	R48/20/22	Harmful: danger of serious damage to health by prolonged exposure through
R14/15	Reacts violently with water, liberating extremely flammable gases.	1140/20/22	inhalation and if swallowed.
R15	Contact with water liberates extremely flammable gases.	R48/21	Harmful: danger of serious damage to health by prolonged exposure in
R15/29	Contact with water liberates extremely flammable gases.	1140/21	contact with skin.
R16	Explosive when mixed with oxidizing substances.	R48/21/22	Harmful: danger of serious damage to health by prolonged exposure in
R17	Spontaneously flammable in air.	1140/21/22	contact with skin and if swallowed.
R18	In use, may form flammable/explosive vapour-air mixture.	R48/22	Harmful: danger of serious damage to health by prolonged exposure if
R19	May form explosive peroxides.	1140/22	swallowed.
R20	Harmful by inhalation.	R48/23	Toxic: danger of serious damage to health by prolonged exposure through
R20/21	Harmful by inhalation and in contact with skin.	1140/20	inhalation.
R20/21/22	Harmful by inhalation, in contact with skin and if swallowed.	R48/23/24	Toxic: danger of serious damage to health by prolonged exposure through
R20/22	Harmful by inhalation and if swallowed.	1110/20/21	inhalation and in contact with skin.
R21	Harmful in contact with skin.	R48/23/24/25	Toxic: danger of serious damage to health by prolonged exposure through
R21/22	Harmful in contact with skin and if swallowed.	1110/20/21/20	inhalation, in contact with skin and if swallowed.
R22	Harmful if swallowed.	R48/23/25	Toxic: danger of serious damage to health by prolonged exposure through
R23	Toxic by inhalation.	1110/20/20	inhalation and if swallowed.
R23/24	Toxic by inhalation and in contact with skin.	R48/24	Toxic: danger of serious damage to health by prolonged exposure in contact
R23/24/25	Toxic by inhalation, in contact with skin and if swallowed.	,	with skin.
R23/24/25E	Also toxic by inhalation, in contact with skin and if swallowed.	R48/24/25	Toxic: danger of serious damage to health by prolonged exposure in contact
R23/25	Toxic by inhalation and if swallowed.		with skin and if swallowed.
R23/25E	Also toxic by inhalation and if swallowed.	R48/25	Toxic: danger of serious damage to health by prolonged exposure if
R24	Toxic in contact with skin.	, ,	swallowed.
R24/25	Toxic in contact with skin and if swallowed.	R49	May cause cancer by inhalation.
R25	Toxic if swallowed.	R50	Very toxic to aquatic organisms.
R26	Very toxic by inhalation.	R50/53	Very toxic to aquatic organisms, may cause long-term adverse effects in the
R26/27	Very toxic by inhalation and in contact with skin.	· ·	aquatic environment.
R26/27/28	Very toxic by inhalation, in contact with skin and if swallowed.	R51	Toxic to aquatic organisms.
R26/28	Very toxic by inhalation and if swallowed.	R51/53	Toxic to aquatic organisms, may cause long-term adverse effects in the
R27	Very toxic in contact with skin.		aquatic environment.
R27/28	Very toxic in contact with skin and if swallowed.	R52	Harmful to aquatic organisms.
R28	Very toxic if swallowed.	R52/53	Harmful to aquatic organisms, may cause long-term adverse effects in the
R29	Contact with water liberates toxic gas.		aquatic environment.
R30	Can become highly flammable in use.	R53	May cause long-term adverse effects in the aquatic environment.
R31	Contact with acids liberates toxic gas.	R54	Toxic to flora.
R32	Contact with acids liberates very toxic gas.	R55	Toxic to fauna.
R33	Danger of cumulative effects.	R56	Toxic to soil organisms.
R34	Causes burns.	R57	Toxic to bees.
R35	Causes severe burns.	R58	May cause long-term adverse effects in the environment.
R36	Irritating to eyes.	R59	Dangerous for the ozone layer.
R36/37	Irritating to eyes and respiratory system.	R60	May impair fertility.
R36/37/38	Irritating to eyes, respiratory system and skin.	R61	May cause harm to the unborn child.
R36/38	Irritating to eyes and skin.	R62	Possible risk of impaired fertility.
R37	Irritating to respiratory system.	R63	Possible risk of harm to the unborn child.
R37/38	Irritating to respiratory system and skin.	R64	May cause harm to breastfed babies
R38	Irritating to skin.	R65	Harmful: may cause lung damage if swallowed.
R39	Danger of very serious irreversible effects.	R66	Repeated exposure may cause skin dryness or cracking.
R39/23	Toxic: danger of very serious irreversible effects through inhalation.	R67	Vapours may cause drowsiness and dizziness.
R39/23/24	Toxic: danger of very serious irreversible effects through inhalation and in	R68	Possible risks of irreversible effects.
			Harmful: passible risk of irroversible offeets through inhalation
B30/33/34/3E	contact with skin. Toxic: danger of year serious irreversible effects through inhelation, in contact	R68/20 R68/20/21	Harmful: possible risk of irreversible effects through inhalation.
R39/23/24/25	Toxic: danger of very serious irreversible effects through inhalation, in contact	R68/20/21	Harmful: possible risk of irreversible effects through inhalation and in contact
	Toxic: danger of very serious irreversible effects through inhalation, in contact with skin and if swallowed.	R68/20/21	Harmful: possible risk of irreversible effects through inhalation and in contact with skin.
R39/23/24/25 R39/23/25	Toxic: danger of very serious irreversible effects through inhalation, in contact with skin and if swallowed. Toxic: danger of very serious irreversible effects through inhalation and if		Harmful: possible risk of irreversible effects through inhalation and in contact with skin. Harmful: possible risk of irreversible effects through inhalation, in contact with
R39/23/25	Toxic: danger of very serious irreversible effects through inhalation, in contact with skin and if swallowed. Toxic: danger of very serious irreversible effects through inhalation and if swallowed.	R68/20/21 R68/20/21/22	Harmful: possible risk of irreversible effects through inhalation and in contact with skin. Harmful: possible risk of irreversible effects through inhalation, in contact with skin and if swallowed.
R39/23/25 R39/24	Toxic: danger of very serious irreversible effects through inhalation, in contact with skin and if swallowed. Toxic: danger of very serious irreversible effects through inhalation and if swallowed. Toxic: danger of very serious irreversible effects in contact with skin.	R68/20/21	Harmful: possible risk of irreversible effects through inhalation and in contact with skin. Harmful: possible risk of irreversible effects through inhalation, in contact with skin and if swallowed. Harmful: possible risk of irreversible effects through inhalation and if
R39/23/25	Toxic: danger of very serious irreversible effects through inhalation, in contact with skin and if swallowed. Toxic: danger of very serious irreversible effects through inhalation and if swallowed. Toxic: danger of very serious irreversible effects in contact with skin. Toxic: danger of very serious irreversible effects in contact with skin and if	R68/20/21 R68/20/21/22 R68/20/22	Harmful: possible risk of irreversible effects through inhalation and in contact with skin. Harmful: possible risk of irreversible effects through inhalation, in contact with skin and if swallowed. Harmful: possible risk of irreversible effects through inhalation and if swallowed.
R39/23/25 R39/24 R39/24/25	Toxic: danger of very serious irreversible effects through inhalation, in contact with skin and if swallowed. Toxic: danger of very serious irreversible effects through inhalation and if swallowed. Toxic: danger of very serious irreversible effects in contact with skin. Toxic: danger of very serious irreversible effects in contact with skin and if swallowed.	R68/20/21 R68/20/21/22 R68/20/22 R68/21	Harmful: possible risk of irreversible effects through inhalation and in contact with skin. Harmful: possible risk of irreversible effects through inhalation, in contact with skin and if swallowed. Harmful: possible risk of irreversible effects through inhalation and if swallowed. Harmful: possible risk of irreversible effects in contact with skin.
R39/23/25 R39/24 R39/24/25 R39/25	Toxic: danger of very serious irreversible effects through inhalation, in contact with skin and if swallowed. Toxic: danger of very serious irreversible effects through inhalation and if swallowed. Toxic: danger of very serious irreversible effects in contact with skin. Toxic: danger of very serious irreversible effects in contact with skin and if swallowed. Toxic: danger of very serious irreversible effects if swallowed.	R68/20/21 R68/20/21/22 R68/20/22	Harmful: possible risk of irreversible effects through inhalation and in contact with skin. Harmful: possible risk of irreversible effects through inhalation, in contact with skin and if swallowed. Harmful: possible risk of irreversible effects through inhalation and if swallowed. Harmful: possible risk of irreversible effects in contact with skin. Harmful: possible risk of irreversible effects in contact with skin and if
R39/23/25 R39/24 R39/24/25 R39/25 R39/26	Toxic: danger of very serious irreversible effects through inhalation, in contact with skin and if swallowed. Toxic: danger of very serious irreversible effects through inhalation and if swallowed. Toxic: danger of very serious irreversible effects in contact with skin. Toxic: danger of very serious irreversible effects in contact with skin and if swallowed. Toxic: danger of very serious irreversible effects if swallowed. Very toxic: danger of very serious irreversible effects through inhalation.	R68/20/21 R68/20/21/22 R68/20/22 R68/21 R68/21/22	Harmful: possible risk of irreversible effects through inhalation and in contact with skin. Harmful: possible risk of irreversible effects through inhalation, in contact with skin and if swallowed. Harmful: possible risk of irreversible effects through inhalation and if swallowed. Harmful: possible risk of irreversible effects in contact with skin. Harmful: possible risk of irreversible effects in contact with skin and if swallowed.
R39/23/25 R39/24 R39/24/25 R39/25	Toxic: danger of very serious irreversible effects through inhalation, in contact with skin and if swallowed. Toxic: danger of very serious irreversible effects through inhalation and if swallowed. Toxic: danger of very serious irreversible effects in contact with skin. Toxic: danger of very serious irreversible effects in contact with skin and if swallowed. Toxic: danger of very serious irreversible effects if swallowed. Very toxic: danger of very serious irreversible effects through inhalation. Very toxic: danger of very serious irreversible effects through inhalation and in	R68/20/21 R68/20/21/22 R68/20/22 R68/21	Harmful: possible risk of irreversible effects through inhalation and in contact with skin. Harmful: possible risk of irreversible effects through inhalation, in contact with skin and if swallowed. Harmful: possible risk of irreversible effects through inhalation and if swallowed. Harmful: possible risk of irreversible effects in contact with skin. Harmful: possible risk of irreversible effects in contact with skin and if
R39/23/25 R39/24 R39/24/25 R39/25 R39/26 R39/26/27	Toxic: danger of very serious irreversible effects through inhalation, in contact with skin and if swallowed. Toxic: danger of very serious irreversible effects through inhalation and if swallowed. Toxic: danger of very serious irreversible effects in contact with skin. Toxic: danger of very serious irreversible effects in contact with skin and if swallowed. Toxic: danger of very serious irreversible effects if swallowed. Very toxic: danger of very serious irreversible effects through inhalation. Very toxic: danger of very serious irreversible effects through inhalation and in contact with skin.	R68/20/21 R68/20/21/22 R68/20/22 R68/21 R68/21/22	Harmful: possible risk of irreversible effects through inhalation and in contact with skin. Harmful: possible risk of irreversible effects through inhalation, in contact with skin and if swallowed. Harmful: possible risk of irreversible effects through inhalation and if swallowed. Harmful: possible risk of irreversible effects in contact with skin. Harmful: possible risk of irreversible effects in contact with skin and if swallowed.
R39/23/25 R39/24 R39/24/25 R39/25 R39/26	Toxic: danger of very serious irreversible effects through inhalation, in contact with skin and if swallowed. Toxic: danger of very serious irreversible effects through inhalation and if swallowed. Toxic: danger of very serious irreversible effects in contact with skin. Toxic: danger of very serious irreversible effects in contact with skin and if swallowed. Toxic: danger of very serious irreversible effects if swallowed. Very toxic: danger of very serious irreversible effects through inhalation. Very toxic: danger of very serious irreversible effects through inhalation and in contact with skin. Very toxic: danger of very serious irreversible effects through inhalation, in	R68/20/21 R68/20/21/22 R68/20/22 R68/21 R68/21/22	Harmful: possible risk of irreversible effects through inhalation and in contact with skin. Harmful: possible risk of irreversible effects through inhalation, in contact with skin and if swallowed. Harmful: possible risk of irreversible effects through inhalation and if swallowed. Harmful: possible risk of irreversible effects in contact with skin. Harmful: possible risk of irreversible effects in contact with skin and if swallowed.
R39/23/25 R39/24 R39/24/25 R39/25 R39/26 R39/26/27	Toxic: danger of very serious irreversible effects through inhalation, in contact with skin and if swallowed. Toxic: danger of very serious irreversible effects through inhalation and if swallowed. Toxic: danger of very serious irreversible effects in contact with skin. Toxic: danger of very serious irreversible effects in contact with skin and if swallowed. Toxic: danger of very serious irreversible effects if swallowed. Very toxic: danger of very serious irreversible effects through inhalation. Very toxic: danger of very serious irreversible effects through inhalation and in contact with skin. Very toxic: danger of very serious irreversible effects through inhalation, in contact with skin and if swallowed.	R68/20/21 R68/20/21/22 R68/20/22 R68/21 R68/21/22	Harmful: possible risk of irreversible effects through inhalation and in contact with skin. Harmful: possible risk of irreversible effects through inhalation, in contact with skin and if swallowed. Harmful: possible risk of irreversible effects through inhalation and if swallowed. Harmful: possible risk of irreversible effects in contact with skin. Harmful: possible risk of irreversible effects in contact with skin and if swallowed.
R39/23/25 R39/24 R39/24/25 R39/25 R39/26 R39/26/27	Toxic: danger of very serious irreversible effects through inhalation, in contact with skin and if swallowed. Toxic: danger of very serious irreversible effects through inhalation and if swallowed. Toxic: danger of very serious irreversible effects in contact with skin. Toxic: danger of very serious irreversible effects in contact with skin and if swallowed. Toxic: danger of very serious irreversible effects if swallowed. Very toxic: danger of very serious irreversible effects through inhalation. Very toxic: danger of very serious irreversible effects through inhalation and in contact with skin. Very toxic: danger of very serious irreversible effects through inhalation, in contact with skin and if swallowed. Very toxic: danger of very serious irreversible effects through inhalation, in contact with skin and if swallowed.	R68/20/21 R68/20/21/22 R68/20/22 R68/21 R68/21/22	Harmful: possible risk of irreversible effects through inhalation and in contact with skin. Harmful: possible risk of irreversible effects through inhalation, in contact with skin and if swallowed. Harmful: possible risk of irreversible effects through inhalation and if swallowed. Harmful: possible risk of irreversible effects in contact with skin. Harmful: possible risk of irreversible effects in contact with skin and if swallowed.
R39/23/25 R39/24 R39/24/25 R39/25 R39/26 R39/26/27 R39/26/27/28 R39/26/28	Toxic: danger of very serious irreversible effects through inhalation, in contact with skin and if swallowed. Toxic: danger of very serious irreversible effects through inhalation and if swallowed. Toxic: danger of very serious irreversible effects in contact with skin. Toxic: danger of very serious irreversible effects in contact with skin and if swallowed. Toxic: danger of very serious irreversible effects if swallowed. Very toxic: danger of very serious irreversible effects through inhalation. Very toxic: danger of very serious irreversible effects through inhalation and in contact with skin. Very toxic: danger of very serious irreversible effects through inhalation, in contact with skin and if swallowed. Very toxic: danger of very serious irreversible effects through inhalation and if swallowed.	R68/20/21 R68/20/21/22 R68/20/22 R68/21 R68/21/22	Harmful: possible risk of irreversible effects through inhalation and in contact with skin. Harmful: possible risk of irreversible effects through inhalation, in contact with skin and if swallowed. Harmful: possible risk of irreversible effects through inhalation and if swallowed. Harmful: possible risk of irreversible effects in contact with skin. Harmful: possible risk of irreversible effects in contact with skin and if swallowed.
R39/23/25 R39/24 R39/24/25 R39/25 R39/26 R39/26/27 R39/26/27/28 R39/26/28	Toxic: danger of very serious irreversible effects through inhalation, in contact with skin and if swallowed. Toxic: danger of very serious irreversible effects through inhalation and if swallowed. Toxic: danger of very serious irreversible effects in contact with skin. Toxic: danger of very serious irreversible effects in contact with skin and if swallowed. Toxic: danger of very serious irreversible effects if swallowed. Very toxic: danger of very serious irreversible effects through inhalation. Very toxic: danger of very serious irreversible effects through inhalation and in contact with skin. Very toxic: danger of very serious irreversible effects through inhalation, in contact with skin and if swallowed. Very toxic: danger of very serious irreversible effects through inhalation and if swallowed. Very toxic: danger of very serious irreversible effects through inhalation and if swallowed.	R68/20/21 R68/20/21/22 R68/20/22 R68/21 R68/21/22	Harmful: possible risk of irreversible effects through inhalation and in contact with skin. Harmful: possible risk of irreversible effects through inhalation, in contact with skin and if swallowed. Harmful: possible risk of irreversible effects through inhalation and if swallowed. Harmful: possible risk of irreversible effects in contact with skin. Harmful: possible risk of irreversible effects in contact with skin and if swallowed.
R39/23/25 R39/24 R39/24/25 R39/25 R39/26 R39/26/27 R39/26/27/28 R39/26/28	Toxic: danger of very serious irreversible effects through inhalation, in contact with skin and if swallowed. Toxic: danger of very serious irreversible effects through inhalation and if swallowed. Toxic: danger of very serious irreversible effects in contact with skin. Toxic: danger of very serious irreversible effects in contact with skin and if swallowed. Toxic: danger of very serious irreversible effects if swallowed. Very toxic: danger of very serious irreversible effects through inhalation. Very toxic: danger of very serious irreversible effects through inhalation and in contact with skin. Very toxic: danger of very serious irreversible effects through inhalation, in contact with skin and if swallowed. Very toxic: danger of very serious irreversible effects through inhalation and if swallowed. Very toxic: danger of very serious irreversible effects in contact with skin. Very toxic: danger of very serious irreversible effects in contact with skin.	R68/20/21 R68/20/21/22 R68/20/22 R68/21 R68/21/22	Harmful: possible risk of irreversible effects through inhalation and in contact with skin. Harmful: possible risk of irreversible effects through inhalation, in contact with skin and if swallowed. Harmful: possible risk of irreversible effects through inhalation and if swallowed. Harmful: possible risk of irreversible effects in contact with skin. Harmful: possible risk of irreversible effects in contact with skin and if swallowed.
R39/23/25 R39/24 R39/24/25 R39/25 R39/26 R39/26/27 R39/26/27/28 R39/26/28	Toxic: danger of very serious irreversible effects through inhalation, in contact with skin and if swallowed. Toxic: danger of very serious irreversible effects through inhalation and if swallowed. Toxic: danger of very serious irreversible effects in contact with skin. Toxic: danger of very serious irreversible effects in contact with skin and if swallowed. Toxic: danger of very serious irreversible effects if swallowed. Very toxic: danger of very serious irreversible effects through inhalation. Very toxic: danger of very serious irreversible effects through inhalation and in contact with skin. Very toxic: danger of very serious irreversible effects through inhalation, in contact with skin and if swallowed. Very toxic: danger of very serious irreversible effects through inhalation and if swallowed. Very toxic: danger of very serious irreversible effects in contact with skin. Very toxic: danger of very serious irreversible effects in contact with skin and if swallowed.	R68/20/21 R68/20/21/22 R68/20/22 R68/21 R68/21/22	Harmful: possible risk of irreversible effects through inhalation and in contact with skin. Harmful: possible risk of irreversible effects through inhalation, in contact with skin and if swallowed. Harmful: possible risk of irreversible effects through inhalation and if swallowed. Harmful: possible risk of irreversible effects in contact with skin. Harmful: possible risk of irreversible effects in contact with skin and if swallowed.
R39/23/25 R39/24 R39/24/25 R39/25 R39/26 R39/26/27 R39/26/27/28 R39/26/28 R39/27 R39/27/28	Toxic: danger of very serious irreversible effects through inhalation, in contact with skin and if swallowed. Toxic: danger of very serious irreversible effects through inhalation and if swallowed. Toxic: danger of very serious irreversible effects in contact with skin. Toxic: danger of very serious irreversible effects in contact with skin and if swallowed. Toxic: danger of very serious irreversible effects if swallowed. Very toxic: danger of very serious irreversible effects through inhalation. Very toxic: danger of very serious irreversible effects through inhalation and in contact with skin. Very toxic: danger of very serious irreversible effects through inhalation, in contact with skin and if swallowed. Very toxic: danger of very serious irreversible effects through inhalation and if swallowed. Very toxic: danger of very serious irreversible effects in contact with skin. Very toxic: danger of very serious irreversible effects in contact with skin.	R68/20/21 R68/20/21/22 R68/20/22 R68/21 R68/21/22	Harmful: possible risk of irreversible effects through inhalation and in contact with skin. Harmful: possible risk of irreversible effects through inhalation, in contact with skin and if swallowed. Harmful: possible risk of irreversible effects through inhalation and if swallowed. Harmful: possible risk of irreversible effects in contact with skin. Harmful: possible risk of irreversible effects in contact with skin and if swallowed.
R39/23/25 R39/24 R39/24/25 R39/25 R39/26 R39/26/27 R39/26/27/28 R39/26/28 R39/27/28 R39/27/28	Toxic: danger of very serious irreversible effects through inhalation, in contact with skin and if swallowed. Toxic: danger of very serious irreversible effects through inhalation and if swallowed. Toxic: danger of very serious irreversible effects in contact with skin. Toxic: danger of very serious irreversible effects in contact with skin and if swallowed. Toxic: danger of very serious irreversible effects if swallowed. Very toxic: danger of very serious irreversible effects through inhalation. Very toxic: danger of very serious irreversible effects through inhalation and in contact with skin. Very toxic: danger of very serious irreversible effects through inhalation, in contact with skin and if swallowed. Very toxic: danger of very serious irreversible effects through inhalation and if swallowed. Very toxic: danger of very serious irreversible effects in contact with skin. Very toxic: danger of very serious irreversible effects in contact with skin and if swallowed. Very toxic: danger of very serious irreversible effects in contact with skin and if swallowed.	R68/20/21 R68/20/21/22 R68/20/22 R68/21 R68/21/22	Harmful: possible risk of irreversible effects through inhalation and in contact with skin. Harmful: possible risk of irreversible effects through inhalation, in contact with skin and if swallowed. Harmful: possible risk of irreversible effects through inhalation and if swallowed. Harmful: possible risk of irreversible effects in contact with skin. Harmful: possible risk of irreversible effects in contact with skin and if swallowed.
R39/23/25 R39/24 R39/24/25 R39/25 R39/26 R39/26/27 R39/26/27/28 R39/26/28 R39/27 R39/27/28 R39/27	Toxic: danger of very serious irreversible effects through inhalation, in contact with skin and if swallowed. Toxic: danger of very serious irreversible effects through inhalation and if swallowed. Toxic: danger of very serious irreversible effects in contact with skin. Toxic: danger of very serious irreversible effects in contact with skin and if swallowed. Toxic: danger of very serious irreversible effects if swallowed. Very toxic: danger of very serious irreversible effects through inhalation. Very toxic: danger of very serious irreversible effects through inhalation and in contact with skin. Very toxic: danger of very serious irreversible effects through inhalation, in contact with skin and if swallowed. Very toxic: danger of very serious irreversible effects through inhalation and if swallowed. Very toxic: danger of very serious irreversible effects in contact with skin. Very toxic: danger of very serious irreversible effects in contact with skin and if swallowed. Very toxic: danger of very serious irreversible effects in contact with skin and if swallowed. Very toxic: danger of very serious irreversible effects in contact with skin and if swallowed. Very toxic: danger of very serious irreversible effects in contact with skin and if swallowed.	R68/20/21 R68/20/21/22 R68/20/22 R68/21 R68/21/22	Harmful: possible risk of irreversible effects through inhalation and in contact with skin. Harmful: possible risk of irreversible effects through inhalation, in contact with skin and if swallowed. Harmful: possible risk of irreversible effects through inhalation and if swallowed. Harmful: possible risk of irreversible effects in contact with skin. Harmful: possible risk of irreversible effects in contact with skin and if swallowed.



List of Safety phrases

S1	Koon looked up
	Keep locked up.
S1/2 S2	Keep locked up and out of reach of children.
	Keep out of reach of children.
S3	Keep in a cool place.
S3/7	Keep container tightly closed in a cool place.
S3/9/14	Keep in a cool, well ventilated place away from(incompatible materials to be
00/0/44/40	indicated by the manufacturer)
S3/9/14/49	Keep only in the original container in a cool well ventilated place away from
	(incompatible materials to be indicated by the manufacturer)
S3/9/49	Keep only in the original container in a cool well ventilated place.
S3/14	Keep in a cool place away from (incompatible materials to be indicated by
	the manufacturer).
S4	Keep away from living quarters.
S5	Keep contents under water.
S6	Keep under (inert gas to be specified by the manufacturer)
S7	Keep container tightly closed.
S7/8	Keep container tightly closed and dry.
S7/9	Keep container tightly closed and in a well-ventilated place.
S7/47	Keep container tightly closed and at a temperature not exceeding C.
S8	Keep container dry.
S9	Keep container in a well-ventilated place.
S12	Do not keep the container sealed.
S13	Keep away from food, drink and animal feedingstuffs.
S14	Keep away from highly flammable substances.
S15	Keep away from heat.
S16	Keep away from combustible material
S17	Keep away from combustible material.
S18	Handle and open container with care.
S20	When using do not eat or drink.
S20/21	When using do not eat, drink or smoke.
S21	When using do not smoke.
S22	Do not breathe dust.
S23	Do not breathe gas/fumes/vapour/spray
S23A	Do not breathe vapour.
S24	Avoid contact with skin.
S24/25	Avoid contact with skin and eyes.
S25	Avoid contact with eyes.
S26	In case of contact with eyes, rinse immediately with plenty of water and seek
	medical advice.
S27	Take off immediately all contaminated clothing.
S27/28	After contact with skin, take off immediately all contaminated clothing and
	wash immediately with plenty of(to be specified by the manufacturer)
S28	After contact with skin wash immediately with plenty of water.
S28A	After contact with skin wash skin immediately with water and soap.
S28B	After contact with skin, wash immediately with plenty of copper sulfate
	solution 2%.
S28C	After contact with skin, wash immediately with plenty of polyethylene glycol
	400, then rinse with plenty of water.
S29	Do not empty into drains.
S29/35	Do not empty into drains; dispose of this material and its container in a safe
	way.
S29/56	Do not empty into drains, dispose of this material and its container at
	hazardous or special waste collection point.
S30	Never add water to this product.
S33	Take precautionary measures against static discharges.
S35	This material and its container must be disposed of in a safe way.
S36	Wear suitable protective clothing.
S36/37	Wear suitable protective clothing and gloves.
S36/37/39	Wear suitable protective clothing, gloves and eye/face protection.
S36/39	Wear suitable protective clothing and eye/face protection.
S37	Wear suitable gloves.
S37/39	Wear suitable gloves and eye/face protection.
S38	In case of insufficient ventilation, wear suitable respiratory equipment.
S39	Wear eye/face protection.
S40	To clean the floor and all objects contaminated by this material, use (to be
0.10	specified by the manufacturer).
S41	In case of fire and/or explosion do not breathe fumes.
S41A	In case of fire do not breathe fumes.
S42	During fumigation/spraying wear suitable respiratory equipment.
S43A	In case of fire, use special extinguish powder or dry sand NEVER use water.
S43B	In case of fire use large amounts of water, dry sand, special powder.
S43C	In case of fire use powder, waterspray, foam, halons, carbondioxide.
S43C S43D	
	In case of fire use large amounts of water.
S43h	In case of fire, use sand, carbondioxide or powdered extinguishing agent. Never use water.
S45	
340	In case of accident or if you feel unwell, seek medical advice immediately
246	(show the label where possible).
S46	If swallowed, seek medical advice immediately and show this container or
C47	label.
S47	Keep at temperature not exceeding C (to be specified by the manufaturer).
S47/49	Keep only in the original container at temperature not exceeding (to be
	specified by the manufacturer).

S48 Keep wetted with \ldots (appropriate material to be specified by the hanufatuer).

Keep only in the original container.

Do not mix with . . . (to be specified by the manufacturer).

Use only in well ventilated areas.

Not recommended for interior use on large surface areas.

Avoid exposure - obtain special instructions before use. S49 S50 S51 S52 S53 S56 Dispose of this material and its container at hazardous or special waste collection point. S57 S59 Use appropriate container to avoid environmental contamination.

Refer to manufacturer/supplier for information on recovery/recycling. S5A S5B Keep contents under paraffin oil. Keep contents under petroleum. S60 S61 This material and/or its container must be disposed of as hazardous waste Avoid release to the environment. Refer to special instructions/Safety data sheets.
If swallowed, do not induce vomiting: seek medical advice immediately and show this container or label.

In case of accident by inhalation: remove casualty to fresh air and keep at rest If swallowed, rinse mouth with water (only if the person is conscious) S63 S64

General conditions Mallinckrodt Baker B.V.1)

Deventer, The Netherlands (EU), June 2005

Article 1 General

In these general conditions, the following terms are understood to have the following meanings:

Agreement:

every agreement entered into between Mallinckrodt Baker and Customer, all amendments and/or supplements thereto, and all (legal) acts to prepare to and to carry out an Agreement;

Customer:

a party who places Orders with Mallinckrodt Baker, buys the Products of Mallinckrodt Baker and/or concludes an Agreement with Mallinckrodt Baker; *Mallinckrodt Baker*:

Mallinckrodt Baker B.V., a private company with limited liability incorporated under the laws of The Netherlands, having its registered office and principal place of business in Deventer;

Offer

an offer or tender of Mallinckrodt Baker, including but not limited to delivery-, quality-, and payment details and/or specifications inserted therein, whether oral or in writing, being an invitation to Customer to place an Order;

Order.

every order of Customer, whether oral or in writing;

Order Acknowledgment:

acceptance of an Order by Mallinckrodt Baker, including but not limited to delivery-, quality-, and payment details and/or specifications inserted therein, including in default of a specific Order Acknowledgement any invoice following the (actual) execution of an Order and so replacing such Order Acknowledgement;

Products

all products manufactured, marketed and/or distributed by Mallinckrodt Baker that are subject to an Agreement.

Article 2 Applicability

- 2.1 These general conditions are applicable to all Agreements and are also applicable to all Offers, Orders and Order Acknowledgements, however named. References made in an Agreement, Offer and/or Order Acknowledgement which are an apparent mistake, do not bind Mallinckrodt Baker and do not give Customer any action against Mallinckrodt Baker.
- 2.2 The applicability of any other general terms and conditions used by Customer or referred to by Customer in any way is hereby explicitly excluded and declared not to be applicable.
- 2.3 Apart from provisions of mandatory (rules of) law, the contents of an Agreement shall exclusively be formed by the contents of the Order Acknowledgement and/or invoice following an Order and these general conditions. Stipulations varying from these general conditions must be expressly agreed upon in writing in the Agreement.
- 2.4 If any provision of these general conditions is held to be illegal, void or unenforceable by a court of competent jurisdiction, the other provisions of these general conditions shall remain in full force and effect and the provision(s) that is/are held to be illegal, void or unenforceable shall be deemed to have been replaced by a provision which as closely as possible meets the intention of Mallinckrodt Baker when inserting the original provision.

Article 3 Offers, Orders, Order Acknowledgements, Agreements

3.1 All Offers shall be valid during thirty (30) days, unless a longer or shorter term is stipulated in the Offer and are always made without any

- obligation and are subject to the approval of Mallinckrodt Baker, unless explicitly stated otherwise. Offers are therefore not more than an invitation to Customer to place an Order.
- 3.2 In case Customer places an Order, regardless whether or not Mallinckrodt Baker made a prior Offer, an Agreement will only be realized after Mallinckrodt Baker has confirmed this Order, by means of an Order Acknowledgement, or as soon as Mallinckrodt Baker will actually carry out the Order. As from that moment, Customer and Mallinckrodt Baker shall be deemed to have entered into an Agreement, unless Customer immediately declares in writing that he does not accept the contents of the Order Acknowledgement.
- 3.3 Modifications and/or supplements of an Order Acknowledged and/or an Agreement can only be made by an instrument in writing duly signed by Mallinckrodt Baker and Customer.

Article 4 Establishment Weight and Quantity

- 4.1 Measurements and weighting as mentioned by Mallinckrodt Baker on the submitted measurement and weight certificate shall be decisive for the contents of the Agreement and therefore subsequently decisive for invoicing and payment.
- 4.2 With respect to the Products delivered by Mallinckrodt Baker, in terms of measurement and/or weight a five (5) percent tolerance (upwards or downwards) shall be applied. This means that Customer shall not have any right of reclamation with respect to measurement and/or weight of the Products or other rights in case of deviations within the aforesaid limit of five (5) percent.

Article 5 Delivery and (transfer of) risk

- .1 Without prejudice to Article 7 hereunder, the delivery ('aflevering') of the Products shall be 'ex works', as defined in the Inco Terms as drawn up by the International Chamber of Commerce as in force at the time of the conclusion of the Agreement, unless another delivery clause is included in the Agreement. In case another delivery clause is agreed upon in the Agreement, the agreed delivery clause shall, in the absence of rules included in the INCO TERMS, with regard to its interpretation, be interpreted in accordance with the general accepted terms in trade in The Netherlands.
- 5.2 The delivery period mentioned in the Order Acknowledgement or otherwise indicated by Mallinckrodt Baker, shall start from the date of the Order Acknowledgement. This delivery period is only indicative and shall not be considered as a deadline, unless expressly agreed otherwise in writing. In case Mallinckrodt Baker does not deliver at the deadline of the indicative delivery period, Customer has the possibility to give notice of default ('ingebrekestelling') to Mallinckrodt Baker. Customer has neither the right to dissolve or to terminate the Agreement, unless, without prejudice to article 11 of these general conditions. Mallinckrodt Baker has exceeded the period of delivery after a notice of default of Customer to such extent, that it can not reasonably be expected that Customer shall maintain the Agreement. without Mallinckrodt Baker being obliged to pay damages. Mallinckrodt Baker shall inform Customer as soon as possible, in case it becomes apparent that Mallinckrodt Baker shall exceed the indicative delivery
- 5.3 In case Customer does not take possession of the Products by the time the Products are ready for delivery including but not limited to forwarding and/or shipment to Customer and Mallinckrodt Baker has confirmed this to Customer, the Products shall be deemed to be delivered and as from that moment the risk with regard to the

- Products passes to Customer. Mallinckrodt Baker shall continue to hold the uncollected Products available for Customer and Mallinckrodt Baker shall store these Products at the expense and risk of Customer. Mallinckrodt Baker may always exercise its authority to sell the Products, as referred to in Article 90 of Book 6 of the Dutch Civil Code ('Burgerlijk Wetboek').
- 5.4 As soon as the Products are ready for delivery including but not limited to forwarding and/or shipment to Customer and Customer is informed thereof, or on the date on which the Agreement is terminated or dissolved pursuant to these conditions, the risk for all damage to be sustained by or on account of the Products shall be transferred to Customer, unless expressly agreed otherwise.
- 5.5 In case Mallinckrodt Baker organizes or arranges shipment of the Products to the Customer, Mallinckrodt Baker shall do so on behalf of the Customer and the Customer shall solely be responsible and liable for such shipment. Customer indemnifies and holds Mallinckrodt Baker harmless for all consequences of such shipment. Regardless of the foregoing, the previous paragraphs of this article shall at all times apply as regards the delivery of the Products.

Article 6 Transport and Packaging

- 6.1 In case means of transport and/or packaging is given on loan by Mallinckrodt Baker to Customer, the packaging conditions of the Association of Dutch Chemical Industry and the Federation of Traders in Chemical Products shall apply.
- 6.2 Loading or filling of means of transport and/or packaging made available by Customer shall be effected at Customer's risk, regardless whether such loading or filling and/or (other) means of transport or packaging were subject of advice by or on behalf of Mallinckrodt Baker, unless gross negligence or wilful intention of Mallinckrodt Baker.
- 6.3 If in the opinion of Mallinckrodt Baker the means of transport and/or the packaging do not meet the requirements imposed thereon for security and/or safety reasons, Mallinckrodt Baker shall be authorized to refuse to load or fill such means of transport and/or packaging. In that case Mallinckrodt Baker shall not be liable for whatever consequences, especially not for damage ensuing from delay. Article 5 Paragraph 3 shall be applicable mutatis mutandis.

Article 7 Retention of Title

- 7.1 Mallinckrodt Baker shall remain the owner of all Products sold to Customer, and the legal title of these Products shall as a result be retained, as long as Customer has not duly and adequately fulfilled all obligations resulting or arising from the Agreement (including but not limited to full payment). Delivery as referred to in article 5 above shall therefore only imply that Customer shall become the detentor ('houder') of the Products. Legal transfer of title ('juridische levering') shall be effected by and after full payment is made and all other obligations resulting or arising from the Agreement have been adequately and duly fulfilled. The foregoing also applies in case of partial deliveries.
- 7.2 Without explicit and written permission of Mallinckrodt Baker, Customer shall not be authorised to give third parties the use of the Products, to rent or pledge these, to transfer the title to the Products or to alienate or encumber the Products in any other way, until the ownership shall be transferred to Customer, even if that is part of the usual performance of his business or even if that is the usual destination of the Products. In case Customer contravenes these provisions, Customer shall be obliged to assign his claims against third parties relating to this alienation and/or encumbering of the Products to Mallinckrodt Baker up to the outstanding claims of Mallinckrodt Baker.
- 7.3 If Mallinckrodt Baker is unable to exercise its retention of title to the Products, as a result of the fact that Customer processes, forms or mixes delivered Products or has the same products processed, formed

- or mixed into (new) products, Customer shall have to act in accordance with the instructions of Mallinckrodt Baker as long as Customer has not fulfilled all obligations arising and/or relating to the Agreement and shall act as an agent of Mallinckrodt Baker under the obligation, insofar as necessary, to transfer all its rights to Mallinckrodt Baker, with respect to the (new formed) products, without prejudice to all (other) rights of Mallinckrodt Baker.
- 7.4 As long as Customer has not fulfilled all obligations arising and/or relating to the Agreement and until the ownership of the Products shall be transferred to Customer and/ or Mallinckrodt Baker has reasonable grounds to fear that the Customer will come in default. Mallinckrodt Baker has the right to immediately claim, respectively reclaim the (title to) the Products delivered at any moment. Customer hereby grants Mallinckrodt Baker its irrevocable and unconditional consent to enter its premises for this purpose. On first demand, Customer shall be obliged to return the Products immediately to Mallinckrodt Baker, carriage paid, without any court intervention, summons or notice of default being required. In that event the Agreement shall not be automatically dissolved as a consequence thereof, such without prejudice to the other rights of Mallinckrodt Baker with this respect. If Mallinckrodt Baker dissolves the Agreement after having taken back the Products, the amounts paid by Customer for the execution of his obligations shall be paid back to Customer, taking into account all claims of Mallinckrodt Baker against Customer.
- 7.5 As long as Customer is detentor and not the owner of the Products and has not the title to the Products, Customer is obliged to insure the Products properly, in any case against risks of theft, damage and loss.

Article 8 Prices

- 8.1 All quotations and prices, whether stated in an Offer, Order Acknowledgement, invoice or otherwise, shall be 'ex works' as defined in the IncoTerms referred to in Article 5 of these general conditions, unless agreed upon otherwise between Mallinckrodt Baker and Customer.
- 8.2 All prices are inclusive of interior packaging and cartons, but exclusive of special exterior packaging and exclusive of VAT, alcohol excise, unless agreed otherwise in writing. Special exterior packaging shall be charged at cost price and shall not be taken back, unless agreed otherwise in writing.
- 8.3 Mallinckrodt Baker shall be authorized to increase the agreed prices, after an Order Acknowledgement and/or the conclusion of an Agreement, as a result of cost increases, irrespective of their nature and reason (for example as a result of changes of the prices of auxiliary materials, raw materials, freight, import duties, levies, taxes and other charges such as wages, salaries and social premiums).
- 8.4 Excise duty. In case of delivery of products containing alcohol, the dedicated administrative documents have to be returned on due time to Mallinckrodt Baker, otherwise excise duty will be charged to the Customer

Article 9 Payment

- 9.1 Payment has to be made effectively in the currency provided in the Agreement, respectively in the currency mentioned on the invoices.
- 9.2 In case of partial deliveries, Mallinckrodt Baker is entitled to invoice each part separately, unless otherwise agreed in writing. Mallinckrodt Baker is also entitled to send advance invoices.
- 9.3 All payments have to be made within thirty (30) days after the date of the invoice, unless another term of payment has been agreed upon in writing. Deduction, setoff, discounts or suspension of the payment obligation for whatever reason is not allowed, unless agreed otherwise.
- 9.4 Payments shall be made by Customer by means of money transfer to a bank account designated by Mallinckrodt Baker.

- 9.5 Claims about invoices shall have to be submitted to Mallinckrodt Baker in writing within fourteen (14) days after the date of the respective invoice. If no claim has been submitted within the aforementioned term, the invoice shall be deemed as accepted and approved by
- 9.6 If Customer fails to pay within thirty (30) days after the date of invoice or another term of payment agreed in writing, the consequence thereof will be that:
 - (i) Customer is in default ipso jure, without any notice of default being required and Mallinckrodt Baker will be entitled to dissolve the Agreement wholly or partially, without court intervention being required or to suspend the execution of the Agreement wholly or partially, without any obligation to compensate damages;
 - (ii) All other claims of Mallinckrodt Baker against Customer shall become immediately due and payable;
 - (iii) Customer is due an interest rate of 2.5% points on top of the Dutch statutory interest ('wettelijke rente') on the outstanding amount, calculated as from the due date of the invoice until the day of full payment, without any notice of default being required and without prejudice to any other rights of Mallinckrodt Baker;
 - (iv) Customer is bound to compensate for all damages sustained and to be sustained by Mallinckrodt Baker as a consequence of his overdue payment, including expenses and extra judicial collection costs, inter alia consisting of costs of third parties contracted by Mallinckrodt Baker in order to determine and collect the debt. In case of overdue payment, the amount of the extra judicial costs shall be determined ('begroot') in accordance with the collection fee of the Dutch Bar ('incassotarief van de Nederlandse Orde van Advocaten').
- 9.7 Payments made by Customer always serve to settle all costs due primarily, then (accrued) interest and subsequently collectible invoices outstanding for the longest period, irrespective of whether Customer mentions that the settlement concerns a later invoice.
- 9.8 Mallinckrodt Baker reserves the right to require an adequate security or guarantee (such as letters of credit, bank guarantees etc.) for the fulfilment of Customer's obligations arising from the Agreement, especially in case Customer fails to pay an invoice timely. If Customer refuses to issue the requested security or guarantee, Mallinckrodt Baker shall be entitled to suspend further deliveries and/or to dissolve the Agreement wholly or partially. In that case Customer shall be obliged to compensate for the damage sustained by Mallinckrodt Baker including loss of profits.

Article 10 Failure to Fulfil Obligations

In case of non fulfilment or untimely fulfilment of any obligation arising from and/or relating to the Agreement by Customer, as well as in case Customer applies for moratorium or voluntary winding-up, obtains moratorium or a winding-up order, liquidation of Customer or similar circumstances, including attachment at his expense, all claims of Mallinckrodt Baker against Customer arising from the Agreement shall immediately be entirely due and payable and Customer shall be deemed to be in default by law with immediate effect. Moreover, Mallinckrodt Baker shall be entitled to dissolve the Agreement without court intervention and without any obligation to compensate damages and to take back the delivered unpaid Products in conformity with the provisions set out in Article 7 Paragraph 4, without prejudice to other rights of Mallinckrodt Baker (for example to claim damages, including loss of profits).

Article 11 Force maieure

- 11.1 If as a consequence of force majeure Mallinckrodt Baker is not able to deliver or to deliver timely, the obligation of delivery shall be suspended for a period equal to the duration of force majeure.
- 11.2 With respect to price increases during the period of force majeure, Mallinckrodt Baker shall be authorized to adapt the agreed price - as referred to in article 8 - and to invoice the adapted price to Customer.
- 11.3 If, as a consequence of a force majeure, the proper performance of the Agreement is permanently impossible or if the proper performance of the Agreement is impossible for a period totalling more than one month, Mallinckrodt Baker and Customer shall be entitled to dissolve the Agreement wholly or partially without court intervention. In the event of a force majeure Customer is not entitled to any compensation.
- Force majeure shall include any measure of any authority of whatever nature as a consequence of which delivery is (temporarily) limited or prevented, for example in the event importation and exportation are prohibited or prevented, physical or economic stagnations in transport, war, state of war or emergency, mobilization, danger of war, civil war, regardless of the fact whether The Netherlands are involved therein, riots, strike or lock-out, work-to-rule, picketing or entrance blockades and/or industrial interruptions, fire, flooding, frost, explosion, illness of personnel. both in the enterprise of Mallinckrodt Baker and in enterprises supplying raw materials or auxiliary materials, and in general all events and circumstances beyond the control of Mallinckrodt Baker and in general circumstances which interrupt the regular course of business of Mallinckrodt Baker, which are of such nature that these circumstances cannot reasonably be attributed to Mallinckrodt Baker and/or which were not reasonably foreseeable for Mallinckrodt Baker at the time the Agreement was concluded, and/or circumstances that are as such that Mallinckrodt Baker cannot reasonably be required to fulfill the obligations arising from the Agreement.

Article 12 Liability and Indemnifications

- 12.1 Customer is acquainted with the production process of the Products and knows that the Products are produced and intended for in vitro laboratory use only, unless the label explicitly indicates otherwise. In the event the label does not explicitly indicate that the Products could be used for other purposes than for in vitro laboratory use, Customer declares and undertakes vis à vis Mallinckrodt Baker that, in case Customer does not use the Products for in vitro laboratory purposes, such use is at his own risk and in full compliance and in accordance with all applicable public and/or private regulations. Customer assumes full liability for such use of the Products and Customer indemnifies and holds Mallinckrodt Baker harmless against any and all liability resulting from such use of the Products.
- 12.2 Notices given by Mallinckrodt Baker or in the name of Mallinckrodt Baker relating to inter alia the quality, composition, treatment in the broadest sense, possibilities of application and properties of the Products do not bind Mallinckrodt Baker, unless expressly made in writing in the form of a warranty ('garantie').
- 12.3 Customer shall check accurately the quantity and the quality of the delivered Products immediately upon delivery as referred to in Article 5 of these general conditions. Claims with respect to the quantity and the quality of the Products delivered have to be filed within fourteen (14) days after delivery, in writing, provided that defects with respect to not observable defects in the quality of the Products that could not be observed reasonably by Customer within the aforementioned term, have to be filed in writing immediately after Customer observed or could have observed this defect. If no claim is filed in writing within the aforementioned term, Mallinckrodt Baker shall be deemed to have duly and correctly fulfilled its obligations resulting or arising from the Agreement with regard to inter alia the quantity and quality.

- 12.4 Claims do not give Customer the right to suspend his payments or to appeal to setoff or discount.
- 12.5 Any partial delivery shall be considered to be a separate delivery.
- 12.6 Unless expressly provided otherwise Products shall exclusively be delivered for in vitro laboratory purposes only. Mallinckrodt Baker shall aim at taking care of appropriate labelling, but Mallinckrodt Baker does not accept any responsibility and/or liability therefore. Customer is liable for taking all measures to avoid damage to persons and/or goods. Customer shall be obliged to store the Products and to use them in their original packaging and Customer is not allowed to change the labelling or to remove instructions for a correct use of the Products from the Products and/or the packaging.
- 12.7 Customer shall not resell the Products. Customer declares and undertakes vis à vis Mallinckrodt Baker that, in case Customer resells the Products, this resale is at his own risk. Customer assumes full liability for any resale of the Products and Customer indemnifies and holds Mallinckrodt Baker harmless against any and all liability resulting from such resale of the Products.
- 12.8 Mallinckrodt Baker shall not be liable for consequential damage or loss of profits sustained by Customer, his subordinate or persons employed by him or the Customer's clients of whatever name or for damage to third parties, caused by whole or partial delivery of the Products, delayed or unsound delivery or failure of delivery of the Products. The aforementioned provision does not affect the possible liability of Mallinckrodt Baker pursuant to the mandatory legal provisions with regard to product liability.
- 12.9 The total liability of Mallinckrodt Baker shall never exceed the nett sales price or the nett invoice amount of the Products in question which the Customer is due in consideration of the Agreement, with deduction of possible amounts credited by Mallinckrodt Baker. In the event that Mallinckrodt Baker's liability should exceed the limitation set forth in the previous sentence, Mallinckrodt Baker shall in any case never be liable for more than the amount of the liability insurance payment under the insurance coverage of Mallinckrodt Baker.
- 12.10 Except for cases of gross negligence or intention of Mallinckrodt Baker, Mallinckrodt Baker shall not be liable for material and/or emotional damage including economic damage caused directly or indirectly by the release of any liquid or solid substance or any gas and/or noise pollution and/or vibrations, in general any nuisance to the whole of living and non living elements of the environment, including all moveable and immoveable goods, separately and in conjunction with each other, being: water, air, soil, persons, animals, plants and goods and the relations between them: ecosystem, nature and landscape and/or any other damage to the environment, by or on account of the delivered Products. Except for cases of gross negligence or intention of Mallinckrodt Baker, Customer will have to hold Mallinckrodt Baker harmless against claims of third parties with respect to the abovementioned environmental damage regardless of the form of those claims, and will have to be duly insured against such claims.
- 12.11 Customer shall be obliged to hold Mallinckrodt Baker harmless against any and all claims of third parties, regardless of their form and for whatever reason they are filed, which claims are related to the delivered Products which are marketed by Mallinckrodt Baker or on behalf of Mallinckrodt Baker. Customer is bound to be duly insured against such claims.

Article 13 Intellectual Property

- 13.1 Goods marketed by Mallinckrodt Baker cannot be offered or traded without Mallinckrodt Baker's prior written permission and on conditions to be imposed by Mallinckrodt Baker under the trademarks to which are entitled or mentioning the trade name of Mallinckrodt Baker.
- 13.2 The Products marketed by Mallinckrodt Baker in retail packaging may only be sold in that original retail packaging without any change or damage of the packaging or its contents. The Products may never be

- offered or delivered as bonus or jointly with other products not marketed by Mallinckrodt Baker at one price without the prior written authorization of Mallinckrodt Baker.
- 13.3 All intellectual property rights with respect to the Products, such as patents, royalties, trademarks, models and furthermore all industrial rights and licenses related to the preceding intellectual property rights shall remain property of Mallinckrodt Baker and/or the licensor. In case of a possible infringement of any of the intellectual property rights of Mallinckrodt Baker by third parties, Customer is obliged to inform Mallinckrodt Baker thereof immediately.
- 13.4 Customer shall be bound to impose the preceding obligations set out in this article to his own clients and Customer shall indemnify Purchaser against possible infringements by its clients of the intellectual property rights of Mallinckrodt Baker with respect to the Products.
- 13.5 For each violation or inaccurate fulfilment of any of the provisions of this article, Customer shall forfeit to Mallinckrodt Baker, without any demand, prior notice or judicial intervention being required, an immediately payable amount of EUR 15,000.- (in words: fifteen thousand euros) for such breach to be increased by a penalty of EUR 1,000 (in words: thousand euros) for each day or part of a day that any such breach shall continue, without Mallinckrodt Baker being required to prove any damage or loss and without prejudice to Mallinckrodt Baker's right to recover any damages or losses in excess of the penalties.

Article 14

Transfer of Rights and Obligations and Transition of Control

- 14.1 Without the prior written permission of Mallinckrodt Baker, Customer shall not be entitled to transfer all or part of his rights and obligations under the Agreement to a third party or to have an Agreement executed wholly or partially by a third party.
- 14.2 If after conclusion of an Agreement the control of all or a substantial part of Customer's enterprise is directly or indirectly transferred to others, Mallinckrodt Baker shall be entitled, within thirty (30) days after being informed by Customer of such transfer, to terminate the Agreement unilaterally by registered letter, subject to a notice period of fourteen (14) days, without court intervention, without Mallinckrodt Baker being obliged to pay compensation for any damages and/or losses sustained by Customer.

Article 15 Applicable Law and Jurisdiction

- 15.1 These general conditions and all Offers, Orders, Order Acknowledgements and Agreements to which these general conditions are applicable or obligations ensuing there from, shall be governed by and construed in accordance with Dutch law.
- 15.2 Parties explicitly agree that the United Nations Convention on Contracts for the International Sale of Goods ('Verdrag der Verenigde Naties inzake internationale koopovereenkomsten betreffende roerende zaken') is not applicable to these general conditions.
- 15.3 All disputes relating to or ensuing from Offers, Orders, Order Acknowledgements and/or Agreements to which these general conditions are applicable, shall exclusively be settled by the competent court in The Netherlands.



Catalogue Index

CAS	NO.	PROD. NO.	EINECS	PAGE	CAS NO.	PROD. NO.	EINECS	PAGE	CAS NO.	PROD. NO.	EINECS	PAGE
50.4		2052	000 004 0	470	00.74.7	4450	202 202 2	200	70.05.4	7070	202 202 2	404
	00-0	3858	200-001-8	178	66-71-7	1158	200-629-2	298	76-05-1	7079	200-929-3	424
	00-0	3859	200-001-8	177	67-52-7	1547	200-658-0	85	76-05-1	9470	200-929-3	424
	00-0	3933	200-001-8	179	67-56-1	2685	200-659-6	259	76-13-1	9132	200-936-1	423
	-00-0	7040	200-001-8	178	67-56-1	3400	200-659-6	260	76-59-5	1513	200-971-2	96
	00-0	7041	200-001-8	178	67-56-1	8045	200-659-6	259	76-60-8	1511	200-972-8	95
	00-0	7390	200-001-8	179	67-56-1	8046	200-659-6	257	77-09-8	1159	201-004-7	300
50-0	-00-0	7424	200-001-8	179	67-56-1	8047	200-659-6	260	77-73-6	8599	201-052-9	151
50-0	01-1	4045	200-002-3	189	67-56-1	8402	200-659-6	256	77-86-1	1414	201-064-4	428
50-0	01-1	4075	200-002-3	189	67-56-1	8404	200-659-6	257	77-86-1	1533	201-064-4	428
50-0	01-1	4168	200-002-3	189	67-56-1	8898	200-659-6	203	77-86-1	4109	201-064-4	427
	70.4	0000	000 004 5	007	07.50.4	0040	000 050 0	050	77.00.0	0000		400
	70-4 81-7	2039 1018	200-061-5 200-066-2	387 76	67-56-1 67-56-1	9049 9077	200-659-6 200-659-6	258 255	77-92-9 77-92-9	0090 1995	201-069-1 201-069-1	129 129
	81-7	1914	200-066-2	76	67-56-1	9091	200-659-6	258	78-30-8	8691	201-103-1	423
											201-103-3	
	.99-7 .99-7	0113	200-075-1 200-075-1	185 185	67-56-1 67-56-1	9097 9263	200-659-6 200-659-6	257 255	78-50-2 78-78-4	1462 8113	201-121-3	427 263
		0115										
	99-7	4893	200-075-1	145	67-56-1	9822	200-659-6	256	78-83-1	8056	201-148-0	232
	89-1	1576	200-157-7	142	67-63-0	3401	200-661-7	337	78-83-1	8210	201-148-0	232
	21-7	2028	200-198-0	382	67-63-0	8067	200-661-7	336	78-83-1	9048-03	201-148-0	231
	55-0	2794	200-237-1	262	67-63-0	8068	200-661-7	336	78-92-2	8103	201-158-5	107
56-2	23-5	7016	200-262-8	408	67-63-0	8119	200-661-7	336	78-93-3	8052	201-159-0	265
56-3	23-5	7106	200-262-8	408	67-63-0	8175	200-661-7	334	78-93-3	8053	201-159-0	265
	23-5	7306	200-262-8	408	67-63-0	9334	200-661-7	334	78-93-3	9214	201-159-0	265
	23-5	7420	200-262-8	409	67-63-0	9830	200-661-7	334	79-01-6	7075	201-167-4	422
	40-6	1504	200-202-8	187	67-64-1	3403	200-662-2	46	79-06-1	4081	201-107-4	51
	40-6	1933	200-272-2	187	67-64-1	8001	200-662-2	45	79-08-1	6034	201-173-7	338
	40-6	4059	200-272-2	186	67-64-1	8002	200-662-2	46	79-11-8	0177	201-178-4	123
	81-5	4043	200-289-5	186	67-64-1	8142	200-662-2	45	79-14-1	1930	201-180-5	187
	81-5	7044	200-289-5	186	67-64-1	9254	200-662-2	45	79-27-6	7069	201-191-5	406
	81-5	7158	200-289-5	186	67-66-3	7018	200-663-8	125	79-33-4	6017	201-196-2	234
57-	11-4	0329	200-313-4	390	67-66-3	7019	200-663-8	126	79-33-4	6060	201-196-2	235
57-	13-6	0345	200-315-5	431	67-66-3	7152	200-663-8	126	81-88-9	1389	201-383-9	343
	13-6	0346	200-315-5	431	67-66-3	7386	200-663-8	126	84-74-2	7024	201-557-4	147
	13-6	4111	200-315-5	431	67-66-3	9174	200-663-8	124	85-44-9	0194	201-607-5	304
	48-7	1563	200-333-3	183	67-66-3	9175	200-663-8	125	85-85-8	1156	201-637-9	286
	50-1	0334	200-333-3	393	67-66-3	9257	200-663-8	124	85-86-9	0707	201-638-4	393
	50-1	4097	200-334-9	392	67-68-5	7033	200-664-3	158	87-66-1	1171	201-762-9	341
		4890	200-334-9	392	67-68-5	7093	200-664-3	158	87-69-4	0340	201-766-0	405
	50-1											
	55-6	7064	200-338-0	338	67-68-5	9234	200-664-3	158	87-69-4	2054	201-766-0	405
	55-6 08-2	7169 1514	200-338-0 200-362-1	339 110	68-12-2 68-12-2	7032 7400	200-679-5 200-679-5	157 157	90-02-8 90-15-3	7066 1150	201-961-0 201-969-4	343 273
	00 2	1011	200 002 1		00 12 2		200 070 0	107	00 10 0	1100	201 000 1	270
	86-6	1565	200-400-7	438	68-12-2	9213	200-679-5	156	90-44-8	1231	201-994-0	74
	00-4	1109	200-448-9	162	68-12-2	9344	200-679-5	157	91-17-8	7080	202-046-9	144
60-2	24-2	4049	200-464-6	250	69-65-8	1125	200-711-8	249	92-31-9	0716	202-146-2	420
60-2	29-7	8033	200-467-2	152	69-65-8	1998	200-711-8	249	93-58-3	8123	202-259-7	262
60-2	29-7	8254	200-467-2	154	71-00-1	1603	200-745-3	201	95-14-7	2363	202-394-1	88
60-2	29-7	8401	200-467-2	154	71-23-8	8066	200-746-9	333	95-47-6	8147	202-422-2	438
60-2	29-7	9237	200-467-2	152	71-36-3	8017	200-751-6	106	95-48-7	7085	202-423-8	139
	29-7	9250	200-467-2	153	71-41-0	8009	200-752-1	72	95-50-1	7025	202-425-9	147
	29-7	9259	200-467-2	153	71-43-2	8014	200-753-7	88	95-50-1	9233	202-425-9	147
	73-4	1329	200-515-2	264	71-55-6	7309	200-756-3	422	100-10-7	1279	202-819-0	156
00.1	FO 0	7004	000 500 0	70	70 17 0	7004	000 770 0	077	100 10 5	0400	000 051 5	001
	53-3	7004	200-539-3	73	72-17-3	7391	200-772-0	377	100-42-5	8133	202-851-5	391
	53-3	7005	200-539-3	73	75-05-8	2680	200-835-2	50	100-51-6	3402	202-859-9	89
	54-4	0059	200-540-9	111	75-05-8	8004	200-835-2	50	100-51-6	7010	202-859-9	89
	55-5	1200	200-541-4	412	75-05-8	8134	200-835-2	49	100-66-3	8013	202-876-1	73
	56-6	0342	200-543-5	413	75-05-8	8144	200-835-2	50	100-97-0	1113	202-905-8	196
	76-0	0301	200-550-3	379	75-05-8	8257	200-835-2	49	100-97-0	3371	202-905-8	194
	76-0	2018	200-550-3	380	75-05-8	9012	200-835-2	48	105-53-3	7028	203-305-9	155
	74-1	0709	200-563-4	393	75-05-8	9017	200-835-2	48	106-93-4	7368	203-444-5	146
	17-5	3405	200-578-6	170	75-05-8	9035	200-835-2	49	107-06-2	8042	203-458-1	148
04-	17-5	3406	200-578-6	168	75-05-8	9255	200-835-2	47	107-06-2	8707	203-458-1	148
64-	17-5	3407	200-578-6	170	75-05-8	9478		52	107-15-3	8409	203-468-6	174
	17-5	3408	200-578-6	170	75-05-8	9821	200-835-2	47	107-21-1	7037	203-473-3	175
64-	17-5	8006	200-578-6	167	75-07-0	9020	200-836-8	42	107-22-2	M834-09	203-474-9	187
	17-5	8007	200-578-6	169	75-09-2	7053	200-838-9	150	108-10-1	4855	203-550-1	266
	17-5	8098	200-578-6	167	75-09-2	7153	200-838-9	151	108-10-1	8049	203-550-1	267
	17-5	8228	200-578-6	168	75-09-2	7305	200-838-9	150	108-10-1	8224	203-550-1	267
	17-5	8229	200-578-6	169	75-09-2	9264	200-838-9	149	108-20-3	8072	203-560-6	155
	17-5	8462	200-578-6	166	75-09-2	9295	200-838-9	149	108-24-7	6004	203-564-8	44
	18-6	6008	200-579-1	182	75-09-2	9315	200-838-9	149	108-32-7	7439	203-572-1	338
	18-6	6037	200-579-1	181	75-09-2	9316	200-838-9	150	108-38-3	8146	203-576-3	437
		2422			75 40 5	****	000 010 0	400	400 75 5	7004		*05
	18-6 18-6	6166 8699	200-579-1 200-579-1	181 181	75-12-7 75-12-7	4028 7042	200-842-0 200-842-0	180 180	108-75-8 108-87-2	7081 8050	203-613-3 203-624-3	426 264
	18-6	9820	200-579-1			7354			108-87-2	3411	203-625-9	420
			200-079-1	180	75-25-2		200-854-6	95				
	18-6	9832	200 500 7	182	75-31-0	8821	200-860-9	233	108-88-3	8077	203-625-9	419
	19-7	4706	200-580-7	44	75-56-9	8822	200-879-2	339	108-88-3	8078	203-625-9	419
	19-7	6052	200-580-7	44	75-59-2	6131	200-882-9	411	108-88-3	8714	203-625-9	418
64	19-7	6152	200-580-7	43	75-65-0	8019	200-889-7	107	108-88-3	9336	203-625-9	418
	10.7	6903	200-580-7	42	75-65-0	8212	200-889-7	107	108-88-3	9351	202 625 0	418
64-											203-625-9	
64-	19-7	9524	200-580-7	43	75-85-4	8012	200-889-7	72	108-88-3	9364	203-625-9	419



	CAS NO.	PROD. NO.	EINECS	PAGE	CAS NO.	PROD. NO.	EINECS	PAGE	CAS NO.	PROD. NO.	EINECS	PAGE
_	108-94-1	8028	203-806-2	141	127-18-4	9360	204-825-9	407	877-24-7	4889	212-889-4	319
	108-94-1	0188	203-632-7	299	127-10-4	7030	204-826-4	156	877-24-7	7525	212-889-4	319
	108-95-2	0505	203-632-7	299	130-22-3	1471	204-981-8	53	920-66-1	7108	213-059-4	194
	108-95-2	1480	203-632-7	300	139-33-3	4652	205-358-3	165	1010-21-7	7101	231-867-5	386
	108-95-2	4056	203-632-7	299	139-33-3	4653	205-358-3	164	1066-33-7	0013	213-911-5	62
	109-66-0	8114	203-692-4	290	139-33-3	7125	205-358-3	164	1120-21-4	8271	214-300-6	430
	109-66-0 109-66-0	8239 8287	203-692-4 203-692-4	288 289	139-33-3 140-22-7	7613 1295	205-358-3 205-403-7	164 160	1132-61-2 1135-40-6	4004 4118	214-478-5 214-492-1	272 118
	109-66-0	8685	203-692-4	290	141-43-5	7055	205-483-3	171	1185-53-1	4103	201-064-4	428
	109-66-0	9331	203-692-4	289	141-53-7	0287	205-488-0	365	1185-57-5	1110	214-686-6	67
	400.00.0		000 000 4	000	444.70.0	0007	005 500 4	470	4000 45 0	4007		474
	109-66-0 109-86-4	9333 8084	203-692-4 203-713-7	289 262	141-78-6 141-78-6	8037 8038	205-500-4 205-500-4	173 174	1239-45-8 1303-86-2	4007 1176	214-984-6 215-125-8	171 93
	109-89-7	8030	203-716-3	152	141-78-6	8204	205-500-4	173	1303-96-4	0268	215-540-4	161
	109-99-9	8075	203-726-8	411	141-78-6	9260	205-500-4	172	1303-96-4	0270	215-540-4	161
	109-99-9	8117	203-726-8	411	141-78-6	9276	205-500-4	173	1304-76-3	1032	215-134-7	91
	109-99-9	9439	203-726-8	409	141-78-6	9282	205-500-4	172	1305-62-0	0073	215-137-3	115
	109-99-9 109-99-9	9440 9441	203-726-8 203-726-8	410 410	141-78-6 141-82-2	9831 1124	205-500-4 205-503-0	172 247	1305-62-0 1305-78-8	0074 0078	215-137-3 215-138-9	115 116
	109-99-9	9446	203-726-8	410	142-82-5	8110	205-563-8	192	1305-78-8	0154	215-138-9	116
	110-15-6	0333	203-740-4	391	142-82-5	8111	205-563-8	193	1308-38-9	0083	215-160-9	128
	110-16-7	2790	203-742-5	247	142-82-5	8662	205-563-8	193	1309-37-1	0122	215-168-2	230
	110-18-9 110-26-9	4098 4031	203-744-6 203-750-9	412 264	142-82-5 142-82-5	9177 9185	205-563-8 205-563-8	191 192	1309-48-4 1310-58-3	1122 0222	215-171-9 215-181-3	245 320
	110-54-3	8044	203-750-9	199	142-82-5	9338	205-563-8	191	1310-58-3	0224	215-181-3	321
	110-54-3	8205	203-777-6	197	142-82-5	9365	205-563-8	192	1310-58-3	0385	215-181-3	320
	110-54-3	8668	203-777-6	196	143-19-1	1916	205-591-0	379	1310-58-3	0387	215-181-3	321
	110-54-3	8669	203-777-6	199	143-33-9	0281	205-599-4	361	1310-58-3	4673	215-181-3	323
	110-54-3 110-54-3	8710 9262	205-563-8 203-777-6	198 197	143-66-8 144-55-8	1194 0263	205-605-5 205-633-8	384 366	1310-58-3 1310-58-3	4674 7063	215-181-3 215-181-3	323 321
	110-54-3	9277	203-777-6	197	144-55-8	0394	205-633-8	365	1310-58-3	7129	215-181-3	322
		<u></u>										
	110-54-3	9304	203-777-6	198	144-55-8	1780	205-633-8	366	1310-58-3	7130	215-181-3	322
	110-80-5	8083	203-804-1	171	150-90-3	0620	205-778-7	382	1310-58-3	7203	215-181-3	322
	110-82-7 110-82-7	8026 8104	203-806-2 203-806-2	140 141	151-21-3 151-21-3	2811 4095	205-788-1 205-788-1	364 364	1310-58-3 1310-73-2	7639 0288	215-181-3 215-185-5	322 369
	110-82-7	8706	203-806-2	141	151-50-8	0213	205-792-3	312	1310-73-2	0402	215-185-5	368
	110-82-7	9258	203-806-2	140	288-32-4	1747	206-019-2	224	1310-73-2	0403	215-185-5	368
	110-82-7	9292	203-806-2	140	298-14-6	2940	206-059-0	317	1310-73-2	0404	215-185-5	369
	110-83-8	G045-01	203-807-8	142	302-17-0	1052	206-117-5	123	1310-73-2	0405	215-185-5	369
	110-86-1 110-86-1	8073 8074	203-809-9 203-809-9	340 341	333-20-0 357-57-3	0250 2459	206-370-1 206-614-7	332 96	1310-73-2 1310-73-2	0903 4687	215-185-5 215-185-5	368 376
	110 00 1	5074	200 000 0	011	007 07 0	2-100	200 014 7	- 00	1010 70 2	4007	210 100 0	070
	110-86-1	9393	203-809-9	340	373-02-4	0179	206-761-7	275	1310-73-2	4689	215-185-5	375
	111-19-3	6064	203-843-4	344	471-34-1	0061	207-439-9	113	1310-73-2	4690	215-185-5	375
	111-42-2	7026	203-868-0	151 286	471-34-1 471-34-1	0503 1464	207-439-9	113	1310-73-2	4691	215-185-5	375 376
	111-46-6 111-65-9	8564 8159	203-872-2 203-892-1	283	471-34-1	4922	207-439-9 207-439-9	112 112	1310-73-2 1310-73-2	4715 6164	215-185-5 215-185-5	371
	111-76-2	8414	203-905-0	105	493-52-7	1146	207-776-1	269	1310-73-2	6165	215-185-5	370
	111-84-2	8280	203-913-4	283	497-19-8	0274	207-838-8	358	1310-73-2	7036	215-185-5	372
	111-87-5	7060	203-917-6	284	497-19-8	1400	207-838-8	358	1310-73-2	7067	215-185-5	370
	111-92-2 112-27-6	8578 8689	203-921-8 203-953-2	146 424	497-19-8 497-19-8	1994 2024	207-838-8 207-838-8	358 358	1310-73-2 1310-73-2	7097 7098	215-185-5 215-185-5	372 374
	112-27-0	0000	200-000-2	727	407-10-0	2024	207-030-0	330	1310-73-2	7000	210-100-0	3/4
	112-34-5	8417	203-961-6	105	497-19-8	4923	207-838-8	357	1310-73-2	7099	215-185-5	374
	112-40-3	8172	203-967-9	162	518-47-8	1314	208-253-0	177	1310-73-2	7105	215-185-5	373
	112-53-8	L034-07	203-982-0	162	532-32-1	2021	208-534-8	355	1310-73-2	7109	215-185-5	370
	112-80-1 113-24-6	6035 2456	204-007-1 204-024-4	284 381	538-62-5 540-72-7	1296 0318	208-698-0 208-754-4	160 385	1310-73-2 1310-73-2	7112 7201	215-185-5 215-185-5	373 373
	115-39-9	1369	204-086-2	95	540-84-1	3904	200-734-4	272	1310-73-2	7202	215-185-5	371
	115-40-2	1330	204-087-8	95	540-84-1	8217	208-759-1	426	1310-73-2	7214	215-185-5	373
	119-26-6	1293	204-309-3	159	540-84-1	8715	208-759-1	426	1310-73-2	7216	215-185-5	374
	120-82-1 120-82-1	8430 9444	204-428-0 204-428-0	422 421	540-84-1 540-84-1	9335 9480	208-759-1 208-759-1	425 425	1310-73-2 1310-73-2	7475 7480	215-185-5 215-185-5	374 374
	120-02-1	9444	204-420-0	421	040-04-1	9400	200-759-1	420	1310-73-2	7400	2 10- 100-0	3/4
	121-44-8	8079	204-469-4	424	544-76-3	7194	208-878-9	194	1310-73-2	7557	215-185-5	372
	121-44-8	9111	204-469-4	423	547-58-0	1145	208-925-3	268	1310-73-2	7598	215-185-5	371
	121-54-0	7255	204-479-9	201	548-62-9	1065	208-953-6	139	1310-73-2	7620	215-185-5	371
	121-57-3	1197	204-482-5	394 160	554-13-2	0155	209-062-5	241	1310-73-2	7621	215-185-5	372 236
	122-39-4 123-31-9	1071 1115	204-539-4 204-617-8	221	573-58-0 583-52-8	1262 1973	209-358-4 209-506-8	134 328	1312-81-8 1313-13-9	4991 0170	215-200-5 215-202-6	236
	123-38-6	8819	204-623-0	337	584-08-7	0204	209-529-3	310	1313-60-6	9015	215-209-4	381
	123-51-3	8010	204-633-5	231	584-08-7	0205	209-529-3	310	1313-99-1	1153	215-215-7	276
	123-86-4	8018	204-658-1	106	598-62-9	0172	209-942-9	248	1314-13-2	0363	215-222-5	443
	123-91-1	8031	204-661-8	160	598-82-3	6069	201-196-2	234	1314-13-2	0364	215-222-5	443
	123-91-1	8434	204-661-8	159	617-45-8	1546	210-513-3	77	1314-56-3	0192	215-236-1	304
	123-96-6	7059	204-667-0	284	629-59-4	8192	211-096-0	409	1314-56-3	0506	215-236-1	188
	124-18-5	8168	204-686-4	144	631-61-8	0011	211-162-9	57	1314-62-1	1217	215-239-8	432
	124-68-5	3321	204-709-8	56	631-61-8	0390	211-162-9	57	1317-36-8	2338	215-267-0	239
	125-20-2 126-73-8	1206 7074	204-729-7 204-800-2	414 421	632-99-5 766-09-6	0619 8782	211-189-6 212-161-6	146 175	1317-38-0 1317-38-0	0101 0102	215-269-1 215-269-1	138 137
	120-73-8	0195	204-800-2	308	865-49-6	7413	212-101-0	127	1317-38-0	1878-01	215-269-1	137
	127-09-3	0258	204-823-8	353	872-50-4	9261	212-828-1	268	1330-20-7	3410	215-535-7	437
	127-09-3	0259	204-823-8	353	872-50-4	9345	212-828-1	268	1330-20-7	8080	215-535-7	436
	127-18-4	7165	204-825-9	407	877-24-7	0197	212-889-4	319	1330-20-7	8118	215-535-7	437



CAS NO.	PROD. NO.	EINECS	PAGE	CAS NO.	PROD. NO.	EINECS	PAGE	CAS NO.	PROD. NO.	EINECS	PAGE
1220 20 7	0F10	245 525 7	420	C407.40.F	0225	200 500 0	207	7047 14 5	4004	221 500 2	201
1330-20-7	9516	215-535-7 215-540-4	436 161	6487-48-5	0235 0302	209-506-8	327 378	7647-14-5	4684 4924	231-598-3 231-598-3	361 359
1330-43-4 1333-82-0	0367 0086	215-607-8	128	6834-92-0 7080-50-4	1255	229-912-9 204-854-7	123	7647-14-5 7647-15-6	0271	231-599-9	356
1335-32-6	0142	215-630-3	237	7119-22-7	4163	204-004-7	272	7647-15-6	2022	231-599-9	356
1336-21-6	4807	215-647-6	67	7296-20-1	2690	230-732-8	158	7647-17-8	1955	231-600-2	122
1336-21-6	6005	215-647-6	66	7365-45-9	4018	230-732-8	190	7647-17-8	2012	231-600-2	122
1336-21-6	6051	215-647-6	66	7429-90-5	1003	231-072-3	53	7647-17-8	4042	231-600-2	121
1336-21-6	6125	215-647-6	63	7429-90-5	9502	231-072-3	53	7664-38-2	6024	231-633-2	303
1336-21-6	6162	215-647-6	63	7439-89-6	0134	231-072-3	228	7664-38-2	6908	231-633-2	302
1341-49-7	0014	215-676-4	62	7439-89-6	0135	231-096-4	228	7664-39-3	6013	231-634-8	217
1343-88-0	0567	215-681-1	177	7439-95-4	1120	231-104-6	242	7664-39-3	6079	231-634-8	218
1343-88-0	7061-00	215-681-1	177	7439-97-6	9006	231-106-7	250	7664-39-3	6904	231-634-8	215
1343-98-2	5056	215-683-2	346	7440-31-5	9509	231-141-8	414	7664-39-3	9563	231-634-8	216
1343-98-2	0324-01	215-683-2	347	7440-44-0	1991	231-153-3	51	7664-93-9	4699	231-639-5	403
1344-28-1	0007	215-691-6	55	7440-44-0	E343	231-153-3	51	7664-93-9	4700	231-639-5	403
1344-28-1	8000	215-691-6	55	7440-50-8	1728	231-159-6	134	7664-93-9	4701	231-639-5	403
1344-28-1	0537	215-691-6	55	7440-50-8	9503	231-159-6	134	7664-93-9	4704	231-639-5	403
1344-28-1	1848	215-691-6	55	7440-50-8	9504	231-159-6	134	7664-93-9	6027	231-639-5	397
1344-28-1	1850	215-691-6	55	7440-66-6	0351	231-175-3	439	7664-93-9	6057	231-639-5	398
1401-55-4	1199	215-753-2	404	7440-66-6	0352	231-175-3	439	7664-93-9	6147	231-639-5	398
1465-25-4	1341	215-981-2	274	7440-66-6	0354	231-175-3	439	7664-93-9	6163	231-639-5	396
1600-27-7	1128	216-491-1	252	7440-66-6	0356	231-175-3	439	7664-93-9	6902	231-639-5	395
1611-35-4	1421	216-553-8	438	7446-20-0	0365	231-793-3	443	7664-93-9	7102	231-639-5	402
1634-04-4	9042	216-653-1	263	7446-20-0	1796	231-793-3	444	7664-93-9	7103	231-639-5	402
1634-04-4	9043	216-653-1	263	7446-20-0	4874		444	7664-93-9	7110	231-639-5	402
1662-01-7	1691	216-767-1	161	7446-20-0	7254	231-793-3	444	7664-93-9	7215	231-639-5	402
1762-95-4	0034	217-175-6	71	7447-40-7	0208	231-211-8	311	7664-93-9	7242	231-639-5	399
1762-95-4	2662	217-175-6	71	7447-40-7	0209	231-211-8	311	7664-93-9	7488	231-163-9	342
1762-95-4	4650	217-175-6	71	7447-40-7	0509	231-211-8	312	7664-93-9	7610	231-639-5	399
1762-95-4	7122	217-175-6	71	7447-40-7	4001	231-211-8	311	7664-93-9	7624	231-639-5	399
1787-61-7	1107	217-250-3	166	7447-41-8	0156	231-212-3	241	7664-93-9	7638	231-639-5	402
1936-15-8	4425	217-705-6	284	7447-41-8	0157	231-212-3	241	7681-11-0	0227	231-659-4	325
2353-45-9	1476	217 700 0	175	7447-41-8	4002	231-212-3	241	7681-11-0	0228	231-659-4	326
2437-29-8	1325	219-441-7	246	7487-88-9	0366	231-298-2	246	7681-11-0	0230	231-659-4	326
2832-45-3	2175	220-601-3	200	7487-94-7	1131	231-299-8	252	7681-49-4	0285	231-667-8	364
3012-65-5	0020	221-146-3	62	7550-45-0	9019	231-441-9	417	7681-49-4	2016	231-667-8	365
3051-09-0	1149	221-266-6	273	7553-56-2	0132	231-442-4	225	7681-57-4	0267	231-673-0	363
3076-05-9	4147	221-352-3	122	7553-56-2	1118	231-442-4	225	7681-57-4	1767	231-673-0	363
3244-88-0	1420	221-816-5	51	7553-56-2	4660	231-441-9	227	7681-82-5	0293	231-679-3	377
3251-23-8	0099	221-838-5	137	7553-56-2	4662	231-441-9	227	7681-82-5	0294	231-679-3	377
3333-67-3	0181	222-068-2	275	7553-56-2	7577	231-441-9	226	7697-37-2	4712	231-714-2	282
3486-35-9	4312	222 000 2	441	7558-79-4	0306	231-448-7	366	7697-37-2	6019	231-714-2	280
3811-04-9	0207	223-289-7	310	7558-79-4	4062	231-448-7	366	7697-37-2	6080	231-714-2	281
3829-86-5	1368	223-325-1	298	7601-90-3	4806	231-512-4	291	7697-37-2	6901	231-714-2	278
4264-83-9	3351	224-246-5	282	7601-90-3	6022	231-512-4	292	7697-37-2	7119	231-714-2	281
4432-31-9	4014	224-632-3	271	7601-90-3	6023	231-512-4	294	7697-37-2	7622	231-714-2	282
5324-84-5	2818	226-195-4	283	7601-90-3	9653	231-512-4	293	7697-37-2	9598	231-714-2	279
5329-14-6	1828	226-218-8	393	7631-90-5	0266	231-548-0	367	7699-45-8	6123	231-718-4	441
5470-11-1	1116	226-798-2	221	7631-99-4	0296	231-554-3	378	7704-34-9	0335	231-722-6	394
5470-11-1	1521	226-798-2	222	7631-99-4	2017	231-554-3	378	7705-08-0	9051	231-729-4	229
5470-11-1	2195	226-798-2	221	7632-00-0	0297	231-555-9	379	7722-64-7	0237	231-760-3	328
5574-97-0	2842	226-947-1	407	7632-00-0	0298	231-555-9	379	7722-64-7	3227	231-760-3	328
5574-97-0	2843	218-147-6	407	7646-85-7	0359	231-592-0	441	7722-64-7	4677	231-760-3	329
5625-37-6	4265		304	7646-85-7	0360	231-592-0	442	7722-64-7	7057	231-760-3	329
5808-22-0	1261	204-972-9	129	7646-93-7	0198	231-594-1	320	7722-64-7	7247	231-760-3	329
5949-29-1	0088	201-069-1	129	7647-01-0	4654	231-595-7	214	7722-76-1	0029	231-764-5	61
5949-29-1	2064	201-069-1	130	7647-01-0	4655	231-595-7	214	7722-76-1	0392	231-764-5	60
5965-83-3	1198	202-555-6	394	7647-01-0	4657	231-595-7	213	7722-76-1	4931	231-764-5	60
5970-45-6	0357	209-170-2	440	7647-01-0	6011	231-595-7	206	7722-84-1	5155	231-765-0	219
6046-93-1	0095	205-553-3	135	7647-01-0	6012	231-595-7	207	7722-84-1	7047	231-765-0	220
6080-56-4	0141	206-104-4	238	7647-01-0	6070	231-595-7	209	7722-84-1	7048	231-765-0	218
6104-58-1	1473	228-058-4	134	7647-01-0	6081	231-595-7	209	7723-14-0	9008	231-768-7	303
6104-59-2	1474	228-060-5	134	7647-01-0	6167	231-595-7	211	7726-95-6	9003	231-778-1	94
6106-24-7	0317	212-773-3	384	7647-01-0	6900	231-595-7	208	7726-95-6	9094	231-778-1	94
6131-90-4	0256	204-823-8	354	7647-01-0	7038	231-595-7	213	7727-21-1	0239	231-7781-8	330
6131-90-4	0393	204-823-8	354	7647-01-0	7088	231-595-7	213	7727-21-1	0239	231-781-8	329
6131-90-4	0517	204-823-8	354	7647-01-0	7111	231-595-7	212	7727-43-7	0053	231-784-4	87
6131-90-4	4009	204-823-8	353	7647-01-0	7211	231-595-7	213	7727-43-7	0028	231-786-5	69
6132-04-3	0280	200-675-3	361	7647-01-0	7212	231-595-7	212	7727-54-0	1899	231-786-5	69
6132-04-3	1758	200-675-3	361	7647-01-0	7219	231-595-7	212	7727-73-3	0311	231-820-9	383
6147 50 4	1057	200 755 0	100	7647.01.0	7000	221 EOF 7	212	7707 70 0	1774	221 020 0	202
6147-53-1 6153-56-6	1057 0187	200-755-8 205-634-3	132 285	7647-01-0 7647-01-0	7230 7477	231-595-7 231-595-7	212 213	7727-73-3 7732-18-5	1774 4217	231-820-9 231-791-2	383 435
6303-21-5	6016	228-601-5	222	7647-01-0	7477	231-595-7	213	7732-18-5	4218	231-791-2	435
6381-59-5	0246	206-156-8	331	7647-01-0	7651	231-595-7	213	7732-18-5	4219	231-791-2	433
6381-59-5	1971	206-156-8	331	7647-01-0	7655	231-595-7	211	7732-18-5	6906	231-791-2	434
6381-59-5	7479	205-698-2	176	7647-01-0	9530	231-595-7	207	7732-18-5	9823	231-791-2	434
6381-92-6	1073	205-358-3	163	7647-01-0	9600	231-595-7	210	7733-02-0	7218	231-793-3	444
6381-92-6	2060	205-358-3	163	7647-14-5	0277	231-598-3	360	7757-79-1	0231	231-818-8	326
6381-92-6	4040	205-358-3	163	7647-14-5	0278	231-598-3	360	7757-79-1	0232	231-818-8	327
6484-52-2	0025	229-347-8	69	7647-14-5	4058	231-598-3	360	7757-82-6	0312	231-820-9	383
1 3.0.022	0020	0 0			.000		000		00.2	5 0	000



CAS NO.	PROD. NO.	EINECS	PAGE	CAS NO.	PROD. NO.	EINECS	PAGE	CAS NO.	PROD. NO.	EINECS	PAGE
7757 00 0	0010	221 020 0	202	7700 04 2	1040	222 221 6	100	10125 12.0	0007	221 210 2	100
7757-82-6 7757-82-6	0313 3375	231-820-9 231-820-9	383 382	7790-84-3 7791-07-3	1043 2815	233-331-6 231-511-9	109 380	10125-13-0 10125-13-0	0097 0515	231-210-2 231-210-2	136 137
7757-82-6	3377-00	231-820-9	382	7791-13-1	1059	231-589-4	133	10125-13-0	0362	231-943-8	442
7757-83-7	0316	231-821-4	384	7791-18-6	0162	232-094-6	244	10196-18-6	1795	231-943-8	442
7757-83-7	1777	231-821-4	384	7791-18-6	0163	232-094-6	244	10213-10-2	0321	236-743-4	386
7758-01-2	4667	231-829-8	309	7791-18-6	4003	232-094-6	244	10213-10-2	1195	236-743-4	386
7758-02-3	0202	231-830-3	309	7791-20-0	0182	231-743-0	275	10277-43-7	0597	233-238-0	236
7758-02-3	0218	231-830-3	309	7791-20-0	0183	231-743-0	276	10294-26-5	1183	233-653-7	350
7758-02-3	0373	231-830-3	309	7803-55-6	1010	232-261-3	68	10294-26-5	1536	233-653-7	351
7758-05-6	0226	231-831-9	325	8002-74-2	3925		430	10294-42-5	2008	237-029-5	120
7758-05-6 7758-09-0	4676 0233	231-831-9 231-832-4	325 327	8002-74-2 8006-28-8	9513 1783	232-315-6	288 352	10294-42-5 10326-27-9	7114 0045	237-029-5 233-788-1	121 87
7758-11-4	0241	231-834-5	318	8012-95-1	7160	232-384-2	287	10326-27-9	7248	255-766-1	86
7758-11-4	0242	231-834-5	318	8032-32-4	8115	232-453-7	296	10326-27-9	7645		86
7758-11-4	4012	231-834-5	318	8032-32-4	8116	232-453-7	297	10361-29-2	0015	233-786-0	58
7758-89-6	0108	231-842-9	136	8032-32-4	8240	232-453-7	296	10361-29-2	0017	233-786-0	58
7758-95-4	0145	231-845-5	238	8032-32-4	8241	232-453-7	297	10361-29-2	0391	233-786-0	58
7758-98-7	0107	231-847-6	138	8032-32-4	8242	232-453-7	297	10450-60-9	1366	233-937-0	295
7758-98-7	7478	231-847-6	176	8032-32-4	8243	232-453-7	298	12007-60-2	1465	234-514-3	242
7758-99-8	0104	231-847-6	138	8032-32-4	9265	232-453-7	295	12054-85-2	0024	234-722-4	61
7758-99-8	0105	231-847-6	139	8032-32-4	9266	232-453-7	295	12054-85-2	1890	234-722-4	61
7761-88-8	1182	231-853-9	349	8032-32-4	9270	232-453-7	296	12125-01-8	0023	235-185-9	61
7761-88-8	4681	231-853-9	350	8032-32-4	9272	232-453-7	295	12125-02-9	0018	235-186-4	59
7761-88-8	4682	231-853-9	350	8049-22-7	1219	200 250 4	435	12125-02-9	0019	235-186-4	59
7761-88-8	7095	231-853-9	349	9002-18-0	1892	232-658-1	52	12125-28-9	0161	235-192-7	244
7761-88-8 7761-88-8	7096 7245	231-853-9 231-853-9	349 350	9002-92-0 9002-93-1	1240 2838		94 428	12135-76-1 12501-23-4	7003 1165	235-223-4 215-682-7	71 429
7761-88-8	7246	231-853-9	349	9002-93-1	0327	232-679-6	389	12627-53-1	3813	235-732-1	239
7772-98-7	0320	231-867-5	385	9005-25-8	1130	232-686-4	389	13291-61-7	1269	236-308-9	145
7773-06-0	1011	231-871-7	70	9005-64-5	7374	202 000 4	429	13446-18-9	0164	233-826-7	245
7774-29-0	1133	231-873-8	252	9005-65-6	7394		429	13446-34-9	0173	231-869-6	248
7775-14-6	2071	231-890-0	363	9012-36-6	A426	232-731-8	52	13453-69-5	0382	236-631-5	240
7778-18-9	0536	231-900-3	117	9048-46-8	0604	232-936-2	53	13463-67-7	4962	236-675-5	417
7778-18-9	1701	231-900-3	116	10022-31-8	0052	233-020-5	87	13472-35-0	1768	231-449-2	362
7778-18-9	1702	231-900-3	117	10025-69-1	0323	231-868-0	415	13472-36-1	0309	231-767-1	362
7778-18-9	1703	231-900-3	117	10025-69-1	0325	231-868-0	415	13472-36-1	1766	231-767-1	362
7778-50-9	0214	231-906-6	313	10025-69-1	0520	231-868-0	415	13477-34-4	0076	233-332-1	116
7778-50-9 7778-50-9	0215 0216	231-906-6 231-906-6	313 313	10025-70-4 10025-77-1	0331 0119	233-971-6 231-729-4	391 229	13478-00-7 13478-10-9	0184 0125	236-068-5 213-843-4	276 229
7778-50-9	4671	231-906-6	314	10025-91-9	0037	233-047-2	75	13573-16-5	1388	237-003-3	342
7778-50-9	7134	231-906-6	314	10026-06-9	0322	231-588-9	416	13746-66-2	0219	237-323-3	317
7778-53-2	0243	231-907-1	330	10026-22-9	0091	233-402-1	133	13755-38-9	1190	238-373-9	380
7778-77-0	0240	231-913-4	316	10026-24-1	0093	233-334-2	133	14459-95-1	0220	237-722-2	317
7778-77-0	1737	231-913-4	316	10028-24-7	0326	231-448-7	367	14634-91-4	7333	238-676-6	176
7778-77-0	4008	231-913-4	315	10034-81-8	0035	233-108-3	72	14807-96-6	0339	238-877-9	404
7778-77-0	4921	231-913-4	315	10034-88-5	0264	240-778-0	367	15086-94-9	1303	239-138-3	165
7778-80-5	0247	231-915-5	332	10034-96-5	0174	232-089-9	249	15244-10-7	0123	233-072-9	230
7778-80-5	0512	231-915-5	332	10034-99-8	0168	231-298-2	246	15708-41-5	1883	239-802-2	165
7779-88-6 7782-49-2	7476 2036	231-943-8 231-957-4	442 344	10034-99-8 10035-04-8	0169 0064	231-298-2 233-140-8	246 114	16674-78-5 16731-55-8	0160 0199	205-554-9 240-795-3	243 316
7782-61-8	0121	233-899-5	230	10035-04-8	0504	233-140-8	114	16774-21-3	1050	240-827-6	58
7782-63-0	0121	231-753-5	230	10035-04-8	1029	233-791-8	91	16921-30-5	1442	240-979-3	317
7782-63-0	0127	231-753-5	231	10035-10-6	6010	233-113-0	204	16940-66-2	9098	241-004-4	356
7782-85-6	0305	231-448-7	367	10039-26-6	0138	200-559-2	235	16940-66-2	9161	241-004-4	355
7782-91-4	0176	231-970-5	271	10039-32-4	0304	231-448-7	367	17927-65-0	1889	233-135-0	56
7782-99-2	6050		404	10043-35-3	0055	233-139-2	92	18472-87-2	1468	242-355-6	301
7783-20-2	0032	231-984-1	70	10043-35-3	0501	233-139-2	93	18497-13-7	2890-03	241-010-7	127
7783-20-2	4027	231-984-1	70	10043-35-3	4035	233-139-2	92	21645-51-2	0005	244-492-7	54
7783-28-0 7783-35-9	0031 1140	231-987-8 231-992-5	62 254	10043-35-3 10043-52-4	5168 0070	233-139-2 233-140-8	92 113	21908-53-2 22767-49-3	1136 2841	244-654-7 245-208-4	253 290
7783-35-9	7135	231-992-5	254	10043-52-4	7116	233-140-8	114	22767-50-6	2173	245-210-5	193
7783-35-9	7642	000 000 4	254	10043-52-4	7232	233-140-8	114	25322-68-3	7176	203-473-3	306
7783-83-7	0118 0124	233-382-4	68 68	10045-94-0	1134	233-152-3	253 363	25322-68-3	7178	203-473-3	306 307
7783-85-9 7783-90-6	1180	233-151-8 232-033-3	348	10049-21-5 10049-21-5	0303 4011	231-449-2 231-449-2	362	25322-69-4 26628-22-8	7179 9099	200-338-0 247-852-1	355
7784-27-2	0006	236-751-8	54	10043-21-3	1055	233-038-3	128	27565-41-9	F780	247-002-1	130
7784-31-8	0010	233-135-0	56	10099-58-8	0383	233-237-5	236	28300-74-5	6968		74
7789-00-6	0210	232-140-5	312	10099-74-8	0148	233-245-9	238	30525-89-4	1157	200-001-8	288
7789-00-6	0211	232-140-5	312	10101-41-4	0081	231-900-3	117	32503-27-8	2846	251-068-5	406
7789-09-5	0021	232-143-1	59	10101-89-0	0307	231-509-8	381	34725-61-6	1472	252-170-2	96
7789-12-0	0282	234-190-3	362	10101-89-0	0308	231-509-8	381	37267-86-0	0190	253-433-4	301
7789-23-3	0225	232-151-5	316	10101-97-0	0185	232-104-9	277	37267-86-0	0676	253-433-4	301
7789-45-9	1067	232-167-2	136	10102-17-7	0319	231-867-5	385	51811-82-6	3815	257-438-2	185
7789-75-5	1046	232-188-7	115	10102-17-7	1779	231-867-5	385	64742-49-0	8043		200
7790-21-8	1169	232-196-0	332	10102-17-7	4693	231-867-5	386	64742-49-0	9267		200
7790-28-5 7790-62-7	1192 0245	232-197-6 232-216-8	385 331	10102-17-7 10102-17-7	4695 7100	231-867-5 231-867-5	386 386	64742-49-0 64742-96-7	9305 8387	265-200-4	199 234
7790-62-7	0378	232-216-8	331	10102-17-7	0295	231-551-7	378	68855-54-9	3376	272-489-0	120
7790-69-4	0158	232-218-9	242	10102-40-0	1039	233-296-7	109	68988-92-1	3816	2,2 100 0	436
7790-78-5	1208-01	233-296-7	109	10124-56-8	0908	233-343-1	365	75277-39-3	4153		191



CAS NO.	PROD. NO.	EINECS	PAGE	CAS NO.	PROD. NO.	EINECS	PAGE	CAS NO	. PROD. NO.	EINECS	PAGE
75621-03-3 90622-57-4	4145 3905	292-459-0	122 429								
100037-69-2	4266		305								



Catalogue Index

EINECS	PROD. NO.	CAS NO.	PAGE	EINECS	PROD. NO.	CAS NO.	PAGE	EINECS	PROD. NO.	CAS NO.	PAGE
000 004 0	4457	00505 00 4	000	000 500 7	0450	04.40.7	40	000 070 0	0000	75 50 0	000
200-001-8	1157 3858	30525-89-4 50-00-0	288 178	200-580-7	6152	64-19-7 64-19-7	43 42	200-879-2	8822	75-56-9 75-59-2	339
200-001-8 200-001-8	3859	50-00-0	176	200-580-7 200-580-7	6903 9524	64-19-7	43	200-882-9 200-889-7	6131 8019	75-65-0	411 107
200-001-8	3933	50-00-0	177	200-618-2	1270	65-85-0	88	200-889-7	8212	75-65-0	107
200-001-8	7040	50-00-0	178	200-629-2	1158	66-71-7	298	200-803-7	8012	75-85-4	72
200-001-8	7041	50-00-0	178	200-658-0	1547	67-52-7	85	200-908-9	0344	76-03-9	421
200-001-8	7385	30-00-0	179	200-659-6	2685	67-56-1	259	200-929-3	7079	76-05-1	424
200-001-8	7390	50-00-0	179	200-659-6	3400	67-56-1	260	200-929-3	9470	76-05-1	424
200-001-8	7424	50-00-0	179	200-659-6	8045	67-56-1	259	200-929-3	9132	76-13-1	423
200-001-3	4045	50-01-1	189	200-659-6	8046	67-56-1	257	200-971-2	1513	76-59-5	96
200-002-3	4075	50-01-1	189	200-659-6	8047	67-56-1	260	200-972-8	1511	76-60-8	95
200-002-3	4168	50-01-1	189	200-659-6	8402	67-56-1	256	201-004-7	1159	77-09-8	300
200-061-5	2039	50-70-4	387	200-659-6	8404	67-56-1	257	201-052-9	8599	77-73-6	151
200-066-2	1018	50-81-7	76	200-659-6	8898	67-56-1	203	201-064-4	1414	77-86-1	428
200-066-2	1914	50-81-7	76	200-659-6	9049	67-56-1	258	201-064-4	1533	77-86-1	428
200-075-1	0113	50-99-7 50-99-7	185 185	200-659-6	9077	67-56-1	255	201-064-4	4103 4109	1185-53-1	428 427
200-075-1	0115			200-659-6	9091	67-56-1	258	201-064-4		77-86-1	
200-075-1	4893	50-99-7	145	200-659-6	9097	67-56-1	257	201-069-1	0088	5949-29-1	129
200-157-7 200-198-0	1576 2028	52-89-1 54-21-7	142 382	200-659-6 200-659-6	9263 9822	67-56-1 67-56-1	255 256	201-069-1 201-069-1	0090 1995	77-92-9 77-92-9	129 129
200-237-1	2794	55-55-0	262	200-661-7	3401	67-63-0	337	201-069-1	2064	5949-29-1	130
200-262-8	7016	56-23-5	408	200-661-7	8067	67-63-0	336	201-103-5	8691	78-30-8	423
200-262-8	7106	56-23-5	408	200-661-7	8068	67-63-0	336	201-121-3	1462	78-50-2	427
200-262-8	7306	56-23-5	408	200-661-7	8119	67-63-0	336	201-142-8	8113	78-78-4	263
200-262-8	7420	56-23-5	409	200-661-7	8175	67-63-0	334	201-148-0	8056	78-83-1	232
200-272-2	1504	56-40-6	187	200-661-7	9334	67-63-0	334	201-148-0	8210	78-83-1	232
200-272-2	1933	56-40-6	187	200-661-7	9830	67-63-0	334	201-148-0	9048-03	78-83-1	231
200-272-2	4059	56-40-6	186	200-662-2	3403	67-64-1	46	201-158-5	8103	78-92-2	107
200-289-5 200-289-5	4043 7044	56-81-5 56-81-5	186 186	200-662-2 200-662-2	8001 8002	67-64-1 67-64-1	45 46	201-159-0 201-159-0	8052 8053	78-93-3 78-93-3	265 265
200-289-5	7158	56-81-5	186	200-662-2	8142	67-64-1	45	201-159-0	9214	78-93-3	265
200-313-4	0329	57-11-4	390	200-662-2	9254	67-64-1	45	201-167-4	7075	79-01-6	422
200-315-5	0345	57-13-6	431	200-663-8	7018	67-66-3	125	201-173-7	4081	79-06-1	51
200-315-5	0346	57-13-6	431	200-663-8	7019	67-66-3	126	201-176-3	6034	79-09-4	338
200-315-5	4111	57-13-6	431	200-663-8	7152	67-66-3	126	201-178-4	0177	79-11-8	123
200-333-3	1563	57-48-7	183	200-663-8	7386	67-66-3	126	201-180-5	1930	79-14-1	187
200-334-9	0334	57-50-1	393	200-663-8	9174	67-66-3	124	201-191-5	7069	79-27-6	406
200-334-9	4097	57-50-1	392	200-663-8	9175	67-66-3	125	201-196-2	6017	79-33-4	234
200-334-9 200-338-0	4890 7064	57-50-1 57-55-6	392 338	200-663-8 200-664-3	9257 7033	67-66-3 67-68-5	124 158	201-196-2 201-196-2	6060 6069	79-33-4 598-82-3	235 234
200-338-0	7169	57-55-6	339	200-664-3	7093	67-68-5	158	201-383-9	1389	81-88-9	343
200-338-0	7179	25322-69-4	307	200-664-3	9234	67-68-5	158	201-557-4	7024	84-74-2	147
200-362-1	1514	58-08-2	110	200-675-3	0280	6132-04-3	361	201-607-5	0194	85-44-9	304
200-400-7	1565	58-86-6	438	200-675-3	1758	6132-04-3	361	201-637-9	1156	85-85-8	286
200-448-9	1109	60-00-4	162	200-679-5	7032	68-12-2	157	201-638-4	0707	85-86-9	393
200-464-6	4049	60-24-2	250	200-679-5	7400	68-12-2	157	201-762-9	1171	87-66-1	341
200-467-2	8033	60-29-7	152	200-679-5	9213	68-12-2	156	201-766-0	0340	87-69-4	405
200-467-2	8254	60-29-7	154	200-679-5	9344	68-12-2	157	201-766-0	2054	87-69-4	405
200-467-2	8401	60-29-7	154	200-711-8	1125	69-65-8	249	201-961-0	7066	90-02-8	343
200-467-2	9237	60-29-7	152	200-711-8	1998	69-65-8	249	201-969-4	1150	90-15-3	273
200-467-2	9250	60-29-7	153	200-745-3	1603	71-00-1	201	201-994-0	1231	90-44-8	74
200-467-2	9259	60-29-7	153	200-746-9	8066	71-23-8	333	202-046-9	7080	91-17-8	144
200-515-2	1329	61-73-4	264	200-751-6	8017	71-36-3	106	202-146-2	0716	92-31-9	420
200-539-3	7004	62-53-3	73	200-752-1	8009	71-41-0	72	202-259-7	8123	93-58-3	262
200-539-3	7005	62-53-3	73	200-753-7	8014	71-43-2	88	202-394-1	2363	95-14-7	88
200-540-9	0059	62-54-4	111	200-755-8	1057	6147-53-1	132	202-422-2	8147	95-47-6	438
200-540-9	7117	00.55.5	111	200-756-3	7309	71-55-6	422	202-423-8	7085	95-48-7	139
200-541-4	1200	62-55-5	412	200-772-0	7391	72-17-3	377	202-425-9	7025	95-50-1	147
200-543-5	0342	62-56-6	413	200-835-2	2680	75-05-8	50	202-425-9	9233	95-50-1	147
200-550-3	0301	62-76-0	379	200-835-2	8004	75-05-8	50	202-555-6	1198	5965-83-3	394
200-550-3	2018	62-76-0	380	200-835-2	8134	75-05-8	49	202-819-0	1279	100-10-7	156
200-559-2	0138	10039-26-6	235	200-835-2	8144	75-05-8	50	202-851-5	8133	100-42-5	391
200-563-4	0709	63-74-1	393	200-835-2	8257	75-05-8	49	202-859-9	3402	100-51-6	89
200-578-6	3405	64-17-5	170	200-835-2	9012	75-05-8	48	202-859-9	7010	100-51-6	89
200-578-6	3406	64-17-5	168	200-835-2	9017	75-05-8	48	202-876-1	8013	100-66-3	73
200-578-6	3407	64-17-5	170	200-835-2	9035	75-05-8	49	202-905-8	1113	100-97-0	196
200-578-6	3408	64-17-5	170	200-835-2	9255	75-05-8	47	202-905-8	3371	100-97-0	194
200-578-6	8006	64-17-5	167	200-835-2	9821	75-05-8	47	203-305-9	7028	105-53-3	155
200-578-6	8007	64-17-5	169	200-836-8	9020	75-07-0	42	203-444-5	7368	106-93-4	146
200-578-6	8098	64-17-5	167	200-838-9	7053	75-09-2	150	203-458-1	8042	107-06-2	148
200-578-6	8228	64-17-5	168	200-838-9	7153	75-09-2	151	203-458-1	8707	107-06-2	148
200-578-6	8229	64-17-5	169	200-838-9	7305	75-09-2	150	203-468-6	8409	107-15-3	174
200-578-6	8462	64-17-5	166	200-838-9	9264	75-09-2	149	203-473-3	1613		306
200-579-1	6008	64-18-6	182	200-838-9	9295	75-09-2	149	203-473-3	1615		306
200-579-1	6037	64-18-6	181	200-838-9	9315	75-09-2	149	203-473-3	1616		306
200-579-1	6166	64-18-6	181	200-838-9	9316	75-09-2	150	203-473-3	7037	107-21-1	175
200-579-1	8699	64-18-6	181	200-842-0	4028	75-12-7	180	203-473-3	7176	25322-68-3	306
200-579-1	9820	64-18-6	180	200-842-0	7042	75-12-7	180	203-473-3	7178	25322-68-3	306
200-579-1	4706	64-19-7	44	200-854-6	7354	75-25-2	95	203-474-9	M834-09	107-22-2	187
200-580-7	6052	64-19-7	44	200-860-9	8821	75-31-0	233	203-550-1	4855	108-10-1	266
250 000 7	0002	31 10 7		230 000 0	0021	, , , , , , ,	200	200 000 1	-1000	.00 10 1	200



EINECS	PROD. NO.	CAS NO.	PAGE	EINECS	PROD. NO.	CAS NO.	PAGE	EINECS	PROD. NO.	CAS NO.	PAGE
203-550-1	8049	108-10-1	267	204-617-8	1115	123-31-9	221	207-439-9	4922	471-34-1	112
203-550-1	8224	108-10-1	267	204-623-0	8819	123-38-6	337	207-776-1	1146	493-52-7	269
203-560-6	8072	108-20-3	155	204-633-5	8010	123-51-3	231	207-838-8	0274	497-19-8	358
203-564-8	6004	108-24-7	44	204-658-1	8018	123-86-4	106	207-838-8	1400	497-19-8	358
203-572-1	7439	108-32-7	338	204-661-8	8031	123-91-1	160	207-838-8	1994	497-19-8	358
203-576-3	8146	108-38-3	437	204-661-8	8434	123-91-1	159	207-838-8	2024	497-19-8	358
203-613-3	7081	108-75-8	426	204-667-0	7059	123-96-6	284	207-838-8	4923	497-19-8	357
203-624-3	8050	108-87-2	264	204-686-4	8168	124-18-5	144	207-838-8	7138		356
203-625-9	3411	108-88-3	420	204-709-8	3321	124-68-5	56	207-838-8	7139		357
203-625-9	8077	108-88-3	419	204-729-7	1206	125-20-2	414	208-253-0	1314	518-47-8	177
203-625-9	8078	108-88-3	419	204-800-2	7074	126-73-8	421	208-534-8	2021	532-32-1	355
203-625-9	8714	108-88-3	418	204-822-2	0195	127-08-2	308	208-698-0	1296	538-62-5	160
203-625-9	9336	108-88-3	418	204-823-8	0256	6131-90-4	354	208-754-4	0318	540-72-7	385
203-625-9	9351	108-88-3	418	204-823-8	0258	127-09-3	353	208-759-1	8217	540-84-1	426
203-625-9	9364	108-88-3	419	204-823-8	0259	127-09-3	353	208-759-1	8715	540-84-1	426
203-628-5	8023	108-90-7	124	204-823-8	0393	6131-90-4	354	208-759-1	9335	540-84-1	425
203-632-7	0188	108-95-2	299	204-823-8	0517	6131-90-4	354	208-759-1	9480	540-84-1	425
203-632-7	0505	108-95-2	299	204-823-8	4009	6131-90-4	353	208-878-9	7194	544-76-3	194
203-632-7	1480	108-95-2	300	204-825-9	7165	127-18-4	407	208-925-3	1145	547-58-0	268
203-632-7	4056	108-95-2	299	204-825-9	9360	127-18-4	407	208-953-6	1065	548-62-9	139
203-032-7	4000	100-90-2	299	204-020-9	9300	127-10-4	407	200-903-0	1000	040-02-9	139
203-692-4	8114	109-66-0	290	204-826-4	7030	127-19-5	156	209-062-5	0155	554-13-2	241
203-692-4		109-66-0	288	204-854-7	1255	7080-50-4	123		0357	5970-45-6	
	8239							209-170-2			440
203-692-4	8287	109-66-0	289	204-972-9	1261	5808-22-0	129	209-358-4	1262	573-58-0	134
203-692-4	8685	109-66-0	290	204-981-8	1471	130-22-3	53	209-506-8	0235	6487-48-5	327
203-692-4	9331	109-66-0	289	205-358-3	1073	6381-92-6	163	209-506-8	1973	583-52-8	328
203-692-4	9333	109-66-0	289	205-358-3	2060	6381-92-6	163	209-529-3	0204	584-08-7	310
203-713-7	8084	109-86-4	262	205-358-3	4040	6381-92-6	163	209-529-3	0205	584-08-7	310
203-716-3	8030	109-89-7	152	205-358-3	4652	139-33-3	165	209-942-9	0172	598-62-9	248
203-726-8	8075	109-99-9	411	205-358-3	4653	139-33-3	164	210-513-3	1546	617-45-8	77
203-726-8	8117	109-99-9	411	205-358-3	4870		165	211-096-0	8192	629-59-4	409
203-726-8	9439	109-99-9	409	205-358-3	7124		163	211-162-9	0011	631-61-8	57
203-726-8	9440	109-99-9	410	205-358-3	7125	139-33-3	164	211-162-9	0390	631-61-8	57
203-726-8	9441	109-99-9	410	205-358-3	7613	139-33-3	164	211-189-6	0619	632-99-5	146
203-726-8	9446	109-99-9	410	205-403-7	1295	140-22-7	160	212-161-6	8782	766-09-6	175
203-740-4	0333	110-15-6	391	205-483-3	7055	141-43-5	171	212-742-4	7413	865-49-6	127
203-742-5	2790	110-16-7	247	205-488-0	0287	141-53-7	365	212-773-3	0317	6106-24-7	384
203-744-6	4098	110-18-9	412	205-500-4	8037	141-78-6	173	212-828-1	9261	872-50-4	268
203-750-9	4031	110-26-9	264	205-500-4	8038	141-78-6	174	212-828-1	9345	872-50-4	268
203-777-6	8044	110-54-3	199	205-500-4	8204	141-78-6	173	212-889-4	0197	877-24-7	319
203-777-6	8205	110-54-3	197	205-500-4	9260	141-78-6	172	212-889-4	4889	877-24-7	319
203-777-6	8668	110-54-3	196	205-500-4	9276	141-78-6	173	212-889-4	7525	877-24-7	319
203-777-6	8669	110-54-3	199	205-500-4	9282	141-78-6	172	213-059-4	7108	920-66-1	194
203-777-6	9262	110-54-3	197	205-500-4	9831	141-78-6	172	213-843-4	0125	13478-10-9	229
203-777-6	9277	110-54-3	197	205-503-0	1124	141-82-2	247	213-911-5	0013	1066-33-7	62
203-777-6	9304	110-54-3	198	205-553-3	0095	6046-93-1	135	214-300-6	8271	1120-21-4	430
203-804-1	8083	110-80-5	171	205-554-9	0160	16674-78-5	243	214-478-5	4004	1132-61-2	272
203-806-2	8026	110-82-7	140	205-563-8	8110	142-82-5	192	214-492-1	4118	1135-40-6	118
203-806-2	8028	108-94-1	141	205-563-8	8111	142-82-5	193	214-686-6	1110	1185-57-5	67
203-806-2	8104	110-82-7	141	205-563-8	8662	142-82-5	193	214-984-6	4007	1239-45-8	171
203-806-2	8706	110-82-7	141	205-563-8	8710	110-54-3	198	215-125-8	1176	1303-86-2	93
203-806-2	9258	110-82-7	140	205-563-8	9177	142-82-5	191	215-134-7	1032	1304-76-3	91
203-806-2	9292	110-82-7	140	205-563-8	9185	142-82-5	192	215-137-3	0073	1305-62-0	115
203-807-8	G045-01	110-83-8	142	205-563-8	9338	142-82-5	191	215-137-3	0074	1305-62-0	115
203-809-9	8073	110-86-1	340	205-563-8	9365	142-82-5	192	215-138-9	0078	1305-78-8	116
203-809-9	8074	110-86-1	341	205-591-0	1916	143-19-1	379	215-138-9	0154	1305-78-8	116
203-809-9	9393	110-86-1	340	205-599-4	0281	143-33-9	361	215-160-9	0083	1308-38-9	128
203-843-4	6064	111-19-3	344	205-605-5	1194	143-66-8	384	215-168-2	0122	1309-37-1	230
203-868-0	7026	111-42-2	151	205-633-8	0263	144-55-8	366	215-171-9	1122	1309-48-4	245
203-872-2	8564	111-46-6	286	205-633-8	0394	144-55-8	365	215-181-3	0222	1310-58-3	320
203-892-1	8159	111-65-9	283	205-633-8	1780	144-55-8	366	215-181-3	0224	1310-58-3	321
203-905-0	8414	111-76-2	105	205-634-3	0187	6153-56-6	285	215-181-3	0385	1310-58-3	320
203-913-4	8280	111-84-2	283	205-634-3	4665		285	215-181-3	0387	1310-58-3	321
203-917-6	7060	111-87-5	284	205-634-3	7141		285	215-181-3	4673	1310-58-3	323
203-921-8	8578	111-92-2	146	205-698-2	7479	6381-59-5	176	215-181-3	4674	1310-58-3	323
203-953-2	8689	112-27-6	424	205-778-7	0620	150-90-3	382	215-181-3	7063	1310-58-3	321
203-961-6	8417	112-34-5	105	205-788-1	2811	151-21-3	364	215-181-3	7129	1310-58-3	322
203-967-9	8172	112-40-3	162	205-788-1	4095	151-21-3	364	215-181-3	7130	1310-58-3	322
203-982-0	L034-07	112-53-8	162	205-792-3	0213	151-50-8	312	215-181-3	7203	1310-58-3	322
204-007-1	6035	112-80-1	284	206-019-2	1747	288-32-4	224	215-181-3	7639	1310-58-3	322
204-024-4	2456	113-24-6	381	206-059-0	2940	298-14-6	317	215-185-5	0288	1310-73-2	369
204-086-2	1369	115-39-9	95	206-104-4	0141	6080-56-4	238	215-185-5	0402	1310-73-2	368
204-087-8	1330	115-40-2	95	206-117-5	1052	302-17-0	123	215-185-5	0403	1310-73-2	368
204-309-3	1293	119-26-6	159	206-156-8	0246	6381-59-5	331	215-185-5	0404	1310-73-2	369
204-428-0	8430	120-82-1	422	206-156-8	1971	6381-59-5	331	215-185-5	0405	1310-73-2	369
204-428-0	9444	120-82-1	421	206-370-1	0250	333-20-0	332	215-185-5	0458		376
204-469-4	8079	121-44-8	424	206-614-7	2459	357-57-3	96	215-185-5	0903	1310-73-2	368
204-469-4	9111	121-44-8	423	206-761-7	0179	373-02-4	275	215-185-5	4687	1310-73-2	376
204-479-9	7255	121-54-0	201	207-439-9	0061	471-34-1	113	215-185-5	4689	1310-73-2	375
204-473-5	1197	121-57-3	394	207-439-9	0503	471-34-1	113	215-185-5	4690	1310-73-2	375
204-539-4	1071	122-39-4	160	207-439-9	1464	471-34-1	112	215-185-5	4691	1310-73-2	375
20.000 7	137 1	50 -	.00	20, 400 0	. 101	., . 5 1	. 12	2.0 100 0	.001		570



EINECS	PROD. NO.	CAS NO.	PAGE	EINECS	PROD. NO.	CAS NO.	PAGE	EINECS	PROD. NO.	CAS NO.	PAGE
215-185-5	4715	1310-73-2	376	221-816-5	1420	3244-88-0	51	231-589-4	1059	7791-13-1	133
215-185-5	6164	1310-73-2	371	221-838-5	0099	3251-23-8	137	231-592-0	0359	7646-85-7	441
215-185-5	6165	1310-73-2	370	222-068-2	0181	3333-67-3	275	231-592-0	0360	7646-85-7	442
215-185-5	7036	1310-73-2	372	223-289-7	0207	3811-04-9	310	231-594-1	0198	7646-93-7	320
215-185-5	7067	1310-73-2	370	223-325-1	1368	3829-86-5	298	231-595-7	4654	7647-01-0	214
215-185-5	7097	1310-73-2	372	224-246-5	3351	4264-83-9	282	231-595-7	4655	7647-01-0 7647-01-0	214
215-185-5 215-185-5	7098 7099	1310-73-2 1310-73-2	374 374	224-632-3 226-195-4	4014 2818	4432-31-9 5324-84-5	271 283	231-595-7 231-595-7	4657 4867	7647-01-0	213 214
215-185-5	7105	1310-73-2	373	226-218-8	1828	5329-14-6	393	231-595-7	6011	7647-01-0	206
215-185-5	7109	1310-73-2	370	226-798-2	1116	5470-11-1	221	231-595-7	6012	7647-01-0	207
215-185-5	7112	1310-73-2	373	226-798-2	1521	5470-11-1	222	231-595-7	6070	7647-01-0	209
215-185-5	7201	1310-73-2	373	226-798-2	2195	5470-11-1	221	231-595-7	6081	7647-01-0	205
215-185-5	7202	1310-73-2	371	226-947-1	2842	5574-97-0	407	231-595-7	6167	7647-01-0	211
215-185-5 215-185-5	7214 7216	1310-73-2 1310-73-2	373 374	228-058-4 228-060-5	1473 1474	6104-58-1 6104-59-2	134 134	231-595-7 231-595-7	6900 7038	7647-01-0 7647-01-0	208 213
215-185-5	7475	1310-73-2	374	228-601-5	6016	6303-21-5	222	231-595-7	7088	7647-01-0	212
215-185-5	7480	1310-73-2	374	229-347-8	0025	6484-52-2	69	231-595-7	7111	7647-01-0	212
215-185-5	7557	1310-73-2	372	229-912-9	0302	6834-92-0	378	231-595-7	7211	7647-01-0	213
215-185-5	7598	1310-73-2	371	230-732-8	2690	7296-20-1	158	231-595-7	7212	7647-01-0	212
215-185-5	7620	1310-73-2	371	230-907-9	4018	7365-45-9	190	231-595-7	7219	7647-01-0	212
215-185-5	7621	1310-73-2	372	231-072-3	1003	7429-90-5	53	231-595-7	7230	7647-01-0	212
215-200-5	4991	1312-81-8 1313-13-9	236	231-072-3 231-096-4	9502	7429-90-5 7439-89-6	53	231-595-7	7477	7647-01-0	213
215-202-6 215-209-4	0170 9015	1313-13-9	248 381	231-096-4	0134 0135	7439-89-6	228 228	231-595-7 231-595-7	7493 7651	7647-01-0 7647-01-0	213 213
215-215-7	1153	1313-99-1	276	231-104-6	1120	7439-95-4	242	231-595-7	7655	7647-01-0	211
215-222-5	0363	1314-13-2	443	231-106-7	9006	7439-97-6	250	231-595-7	9530	7647-01-0	207
215-222-5	0364	1314-13-2	443	231-141-8	9509	7440-31-5	414	231-595-7	9600	7647-01-0	210
215-236-1	0192	1314-56-3	304	231-153-3	1991	7440-44-0	51	231-598-3	0277	7647-14-5	360
215-236-1 215-239-8	0506 1217	1314-56-3 1314-62-1	188 432	231-153-3 231-159-6	E343 1728	7440-44-0 7440-50-8	51 134	231-598-3 231-598-3	0278 4058	7647-14-5 7647-14-5	360 360
215 207 0	2220	1017.00.0	220	221 150 0	0500	7440 50 0	104	221 500 2	4004	7047 14 5	201
215-267-0 215-269-1	2338 0101	1317-36-8 1317-38-0	239 138	231-159-6 231-159-6	9503 9504	7440-50-8 7440-50-8	134 134	231-598-3 231-598-3	4684 4924	7647-14-5 7647-14-5	361 359
215-269-1	0102	1317-38-0	137	231-163-9	7488	7664-93-9	342	231-598-3	7140	7047-14-3	361
215-270-7	1878-01	1317-39-1	137	231-175-3	0351	7440-66-6	439	231-599-9	0271	7647-15-6	356
215-535-7	3410	1330-20-7	437	231-175-3	0352	7440-66-6	439	231-599-9	2022	7647-15-6	356
215-535-7 215-535-7	8080 8118	1330-20-7 1330-20-7	436 437	231-175-3 231-175-3	0354 0356	7440-66-6 7440-66-6	439 439	231-600-2 231-600-2	1955 2012	7647-17-8 7647-17-8	122 122
215-535-7	9516	1330-20-7	436	231-210-2	0097	10125-13-0	136	231-600-2	4042	7647-17-8	121
215-540-4	0268	1303-96-4	161	231-210-2	0515	10125-13-0	137	231-633-2	6024	7664-38-2	303
215-540-4	0270	1303-96-4	161	231-211-8	0208	7447-40-7	311	231-633-2	6908	7664-38-2	302
215-540-4	0367	1330-43-4	161	231-211-8	0209	7447-40-7	311	231-634-8	6013	7664-39-3	217
215-607-8	0086	1333-82-0	128	231-211-8	0509	7447-40-7	312	231-634-8	6079	7664-39-3	218
215-630-3	0142	1335-32-6	237	231-211-8	4001	7447-40-7	311	231-634-8	6904	7664-39-3	215
215-647-6 215-647-6	4807 6005	1336-21-6 1336-21-6	67 66	231-212-3 231-212-3	0156 0157	7447-41-8 7447-41-8	241 241	231-634-8 231-639-5	9563 4699	7664-39-3 7664-93-9	216 403
215-647-6	6051	1336-21-6	66	231-212-3	4002	7447-41-8	241	231-639-5	4700	7664-93-9	403
215-647-6	6125	1336-21-6	63	231-298-2	0168	10034-99-8	246	231-639-5	4701	7664-93-9	403
215-647-6	6162	1336-21-6	63	231-298-2	0169	10034-99-8	246	231-639-5	4704	7664-93-9	403
215-676-4	0014	1341-49-7	62	231-298-2	0366	7487-88-9	246	231-639-5	6027	7664-93-9	397
215-681-1	0567	1343-88-0	177	231-299-8	1131	7487-94-7	252	231-639-5	6057	7664-93-9	398
215-681-1	7061-00	1343-88-0	177	231-441-9	4660	7553-56-2	227	231-639-5	6147	7664-93-9	398
215-682-7 215-683-2	1165 0253	12501-23-4	429 346	231-441-9 231-441-9	4662 7128	7553-56-2	227 226	231-639-5 231-639-5	6163 6902	7664-93-9 7664-93-9	396 395
215-683-2	1827		346	231-441-9	7210		226	231-639-5	7102	7664-93-9	402
215-683-2	5056	1343-98-2	346	231-441-9	7577	7553-56-2	226	231-639-5	7103	7664-93-9	402
215-683-2	0324-01	1343-98-2	347	231-441-9	9019	7550-45-0	417	231-639-5	7110	7664-93-9	402
215-691-6	0007	1344-28-1	55	231-442-4	0132	7553-56-2	225	231-639-5	7215	7664-93-9	402
215-691-6 215-691-6	0008 0537	1344-28-1 1344-28-1	55 55	231-442-4 231-448-7	1118 0304	7553-56-2 10039-32-4	225 367	231-639-5 231-639-5	7242 7471	7664-93-9	399 399
215-691-6	1848	1344-28-1	55	231-448-7	0305	7782-85-6	367	231-639-5	7610	7664-93-9	399
215-691-6	1850	1344-28-1	55	231-448-7	0306	7558-79-4	366	231-639-5	7624	7664-93-9	399
215-753-2	1199	1401-55-4	404	231-448-7	0326	10028-24-7	367	231-639-5	7638	7664-93-9	402
215-981-2	1341	1465-25-4	274	231-448-7	4062	7558-79-4	366	231-659-4	0227	7681-11-0	325
216-491-1	1128	1600-27-7	252	231-449-2	0303	10049-21-5	363	231-659-4	0228	7681-11-0	326
216-553-8	1421	1611-35-4	438	231-449-2	1768	13472-35-0	362	231-659-4	0230	7681-11-0	326
216-653-1 216-653-1	9042 9043	1634-04-4 1634-04-4	263 263	231-449-2 231-509-8	4011 0307	10049-21-5 10101-89-0	362 381	231-659-4 231-667-8	4876 0285	7681-49-4	225 364
216-767-1	1691	1662-01-7	161	231-509-8	0308	10101-89-0	381	231-667-8	2016	7681-49-4	365
217-175-6	0034	1762-95-4	71	231-511-9	2815	7791-07-3	380	231-673-0	0267	7681-57-4	363
217-175-6	2662	1762-95-4	71	231-512-4	4806	7601-90-3	291	231-673-0	1767	7681-57-4	363
217-175-6	4650	1762-95-4	71	231-512-4	6022	7601-90-3	292	231-673-0	3934		178
217-175-6	7122	1762-95-4	71	231-512-4	6023	7601-90-3	294	231-679-3	0293	7681-82-5	377
217-250-3	1107	1787-61-7	166	231-512-4	9653	7601-90-3	293	231-679-3	0294	7681-82-5	377
217-705-6 218-147-6	4425 2843	1936-15-8 5574-97-0	284 407	231-548-0 231-551-7	0266 0295	7631-90-5 10102-40-6	367 378	231-714-2 231-714-2	4712 6019	7697-37-2 7697-37-2	282 280
219-441-7	1325	2437-29-8	246	231-554-3	0296	7631-99-4	378	231-714-2	6080	7697-37-2	281
220-601-3	2175	2832-45-3	200	231-554-3	2017	7631-99-4	378	231-714-2	6901	7697-37-2	278
221-146-3	0020	3012-65-5	62	231-555-9	0297	7632-00-0	379	231-714-2	7119	7697-37-2	281
221-266-6 221-352-3	1149 4147	3051-09-0 3076-05-9	273 122	231-555-9 231-588-9	0298 0322	7632-00-0 10026-06-9	379 416	231-714-2 231-714-2	7622 9598	7697-37-2 7697-37-2	282 279
ZZ 1-00Z-0	717/	5070-00-0	122	201-300-3	0022	10020-00-0	-110	201-714-2	3330	7007-07-2	210



EINECS	PROD. NO.	CAS NO.	PAGE	EINECS	PROD. NO.	CAS NO.	PAGE	EINECS	PROD. NO.	CAS NO.	PAGE
004.740.4	0400	7000 45 0		004.007.5	7404	1010 01 7	200	000 400 0	0504	10010 05 0	
231-718-4 231-722-6	6123 0335	7699-45-8 7704-34-9	441 394	231-867-5 231-868-0	7101 0323	1010-21-7 10025-69-1	386 415	233-139-2	0501 4035	10043-35-3 10043-35-3	93 92
231-729-4	0119	10025-77-1	229	231-868-0	0325	10025-69-1	415	233-139-2 233-139-2	5168	10043-35-3	92
231-729-4	9051	7705-08-0	229	231-868-0	0520	10025-69-1	415	233-140-8	0064	10045-35-3	114
231-743-0	0182	7791-20-0	275	231-869-6	0173	13446-34-9	248	233-140-8	0070	10033-04-8	113
231-743-0	0183	7791-20-0	276	231-871-7	1011	7773-06-0	70	233-140-8	0504	10043-52-4	114
231-743-0	0126	7782-63-0	230	231-873-8	1133	7774-29-0	252	233-140-8	7116	10033-04-8	114
231-753-5	0127	7782-63-0	230	231-890-0	2071	7775-14-6	363	233-140-8	7232	10043-52-4	114
231-760-3	0237	7722-64-7	328	231-900-3	0081	10101-41-4	117	233-151-8	0124	7783-85-9	68
231-760-3	3227	7722-64-7	328	231-900-3	0536	7778-18-9	117	233-151-8	1134	10045-94-0	253
231-760-3 231-760-3	4677	7722-64-7	329	231-900-3	1701	7778-18-9	116 117	233-237-5 233-238-0	0383	10099-58-8	236
231-760-3	7057 7247	7722-64-7 7722-64-7	329 329	231-900-3 231-900-3	1702 1703	7778-18-9 7778-18-9	117	233-245-9	0597 0148	10277-43-7 10099-74-8	236 238
231-764-5	0029	7722-76-1	61	231-906-6	0214	7778-50-9	313	233-296-7	1039	10108-64-2	109
231-764-5	0392	7722-76-1	60	231-906-6	0215	7778-50-9	313	233-296-7	1208-01	7790-78-5	109
231-764-5	4931	7722-76-1	60	231-906-6	0216	7778-50-9	313	233-331-6	1043	7790-84-3	109
231-765-0	5155	7722-84-1	219	231-906-6	4671	7778-50-9	314	233-332-1	0076	13477-34-4	116
231-765-0	7047	7722-84-1	220	231-906-6	7134	7778-50-9	314	233-334-2	0093	10026-24-1	133
231-765-0	7048	7722-84-1	218	231-907-1	0243	7778-53-2	330	233-343-1	0908	10124-56-8	365
231-767-1	0309	13472-36-1	362	231-913-4	0240	7778-77-0	316	233-382-4	0118	7783-83-7	68
231-767-1	1766	13472-36-1	362	231-913-4	1737	7778-77-0	316	233-402-1	0091	10026-22-9	133
231-768-7	9008	7723-14-0	303	231-913-4	4008	7778-77-0	315	233-653-7	1183	10294-26-5	350
231-778-1	9003	7726-95-6	94	231-913-4	4921	7778-77-0	315	233-653-7	1536	10294-26-5	351
231-778-1	9094	7726-95-6	94	231-915-5	0247	7778-80-5	332	233-786-0	0015	10361-29-2	58
231-781-8	0239	7727-21-1	330	231-915-5	0512	7778-80-5	332	233-786-0	0017	10361-29-2	58
231-781-8	0377	7727-21-1	329	231-943-8	0362	10196-18-6	442	233-786-0	0391	10361-29-2	58
231-784-4	0053	7727-43-7	87	231-943-8	1795	10196-18-6	442	233-788-1	0045	10326-27-9	87
231-786-5	0028	7727-54-0	69	231-943-8	7476	7779-88-6	442	233-791-8	1029	10035-06-0	91
231-786-5	1899	7727-54-0	69	231-957-4	2036	7782-49-2	344	233-826-7	0164	13446-18-9	245
231-791-2	4217	7732-18-5	435	231-970-5	0176	7782-91-4	271	233-899-5	0121	7782-61-8	230
231-791-2	4218	7732-18-5	435	231-984-1	0032	7783-20-2	70	233-937-0	1366	10450-60-9	295
231-791-2	4219	7732-18-5	433	231-984-1	4027	7783-20-2	70	233-971-6	0331	10025-70-4	391
231-791-2	6906	7732-18-5	434	231-987-8	0031	7783-28-0	62	234-190-3	0282	7789-12-0	362
231-791-2	7092		163	231-992-5	1140	7783-35-9	254	234-514-3	1465	12007-60-2	242
231-791-2	7481		164	231-992-5	7135	7783-35-9	254	234-722-4	0024	12054-85-2	61
231-791-2	9823	7732-18-5	434	231-992-5	7464		254	234-722-4	1890	12054-85-2	61
231-793-3	0365	7446-20-0	443	232-033-3	1180	7783-90-6	348	235-185-9	0023	12125-01-8	61
231-793-3	1796	7446-20-0	444	232-089-9	0174	10034-96-5	249	235-186-4	0018	12125-02-9	59
231-793-3	7218	7733-02-0	444	232-094-6	0162	7791-18-6	244	235-186-4	0019	12125-02-9	59
231-793-3	7254	7446-20-0	444	232-094-6	0163	7791-18-6	244	235-192-7	0161	12125-28-9	244
231-818-8	0231	7757-79-1	326	232-094-6	4003	7791-18-6	244	235-223-4	7003	12135-76-1	71
231-818-8	0232	7757-79-1	327	232-104-9	0185	10101-97-0	277	235-732-1	3813	12627-53-1	239
231-820-9	0311	7727-73-3	383	232-140-5	0210	7789-00-6	312	236-068-5	0184	13478-00-7	276
231-820-9	0312	7757-82-6	383	232-140-5	0211	7789-00-6	312	236-308-9	1269	13291-61-7	145
231-820-9	0313	7757-82-6	383	232-143-1	0021	7789-09-5	59	236-631-5	0382	13453-69-5	240
231-820-9	1774	7727-73-3	383	232-151-5	0225	7789-23-3	316	236-675-5	4962	13463-67-7	417
231-820-9	3375	7757-82-6	382	232-167-2	1067	7789-45-9	136	236-743-4	0321	10213-10-2	386
231-820-9	3377-00	7757-82-6	382	232-188-7	1046	7789-75-5	115	236-743-4	1195	10213-10-2	386
231-821-4	0316	7757-83-7	384	232-196-0	1169	7790-21-8	332	236-751-8	0006	7784-27-2	54
231-821-4	1777	7757-83-7	384	232-197-6	1192	7790-28-5	385	237-003-3	1388	13573-16-5	342
231-829-8	4667	7758-01-2	309	232-216-8	0245	7790-62-7	331	237-029-5	2008	10294-42-5	120
231-829-8	7133		308	232-216-8	0378	7790-62-7	331	237-029-5	7114	10294-42-5	121
231-830-3	0202	7758-02-3	309	232-218-9	0158	7790-69-4	242	237-323-3	0219	13746-66-2	317
231-830-3	0218	7758-02-3	309	232-261-3	1010	7803-55-6	68	237-722-2	0220	14459-95-1	317
231-830-3	0373	7758-02-3	309	232-315-6	9513	8002-74-2	288	238-373-9	1190	13755-38-9	380
231-831-9	0226	7758-05-6	325	232-384-2	7160	8012-95-1	287	238-676-6	7333	14634-91-4	176
231-831-9	4676	7758-05-6	325	232-453-7	8115	8032-32-4	296	238-877-9	0339	14807-96-6	404
231-832-4	0233	7758-09-0	327	232-453-7	8116	8032-32-4	297	239-138-3	1303	15086-94-9	165
231-834-5	0241	7758-11-4	318	232-453-7	8240	8032-32-4	296	239-802-2	1883	15708-41-5	165
231-834-5	0242	7758-11-4	318	232-453-7	8241	8032-32-4	297	240-778-0	0264	10034-88-5	367
231-834-5	4012	7758-11-4	318	232-453-7	8242	8032-32-4	297	240-795-3	0199	16731-55-8	316
231-842-9	0108	7758-89-6	136	232-453-7	8243	8032-32-4	298	240-827-6	1050	16774-21-3	58
231-845-5	0145	7758-95-4	238	232-453-7	9265	8032-32-4	295	240-979-3	1442	16921-30-5	317
231-847-6	0104	7758-99-8	138	232-453-7	9266	8032-32-4	295	241-004-4	9098	16940-66-2	356
231-847-6	0105	7758-99-8	139	232-453-7	9270	8032-32-4	296	241-004-4	9161	16940-66-2	355
231-847-6	0107	7758-98-7	138	232-453-7	9272	8032-32-4	295	241-010-7	2890-03	18497-13-7	127
231-847-6	7478	7758-98-7	176	232-658-1	1892	9002-18-0	52	242-355-6	1468	18472-87-2	301
231-853-9	1182	7761-88-8	349	232-679-6	0327	9005-25-8	389	244-492-7	0005	21645-51-2	54
231-853-9 231-853-9	4681 4682	7761-88-8 7761-88-8	350 350	232-686-4 232-731-8	1130 A426	9005-25-8 9012-36-6	389 52	244-654-7 245-208-4	1136 2841	21908-53-2 22767-49-3	253 290
231-000-9			300				ÜZ			22101-43-3	230
231-853-9	7095	7761-88-8	349	232-936-2	0604	9048-46-8	53	245-210-5	2173	22767-50-6	193
231-853-9	7096	7761-88-8	349	233-020-5	0052	10022-31-8	87	247-852-1	9099	26628-22-8	355
231-853-9	7245	7761-88-8	350	233-038-3	1055	10060-12-5	128	251-068-5	2846	32503-27-8	406
231-853-9	7246	7761-88-8	349	233-047-2	0037	10025-91-9	75	252-170-2	1472	34725-61-6	96
231-867-5	0319	10102-17-7	385	233-072-9	0123	15244-10-7	230	253-433-4	0190	37267-86-0	301
231-867-5	0320	7772-98-7	385	233-108-3	0035	10034-81-8	72	253-433-4	0676	37267-86-0	301
231-867-5	1779	10102-17-7	385	233-113-0	6010	10035-10-6	204	257-438-2	3815	51811-82-6	185
231-867-5	4693	10102-17-7	386	233-135-0	0010	7784-31-8	56 56	265-200-4	8387	64742-96-7	234
231-867-5	4695	10102-17-7	386	233-135-0	1889	17927-65-0	56	272-489-0	0902	COOEL L4 O	119
231-867-5	7100	10102-17-7	386	233-139-2	0055	10043-35-3	92	272-489-0	3376	68855-54-9	120



EINECS	PROD. NO.	CAS NO.	PAGE	EINECS	PROD. NO.	CAS NO.	PAGE	EINECS	PROD. NO.	CAS NO.	PAGE
292-459-0	3905	90622-57-4	429								



Catalogue Index

PROD. NO.	CAS NO.	EINECS	PAGE	PROD. NO.	CAS NO.	EINECS	PAGE	PROD. NO.	CAS NO.	EINECS	PAGE
2005	04045 54 0	044 400 7		0.400	7704 40 0	000 004 0	044	2000	7700 40 0	004 400 0	000
0005	21645-51-2	244-492-7	54	0162	7791-18-6	232-094-6	244	0282	7789-12-0	234-190-3	362
0006	7784-27-2	236-751-8	54	0163	7791-18-6	232-094-6	244	0285	7681-49-4	231-667-8	364
0007	1344-28-1	215-691-6	55	0164	13446-18-9	233-826-7	245	0287	141-53-7	205-488-0	365
0008	1344-28-1	215-691-6	55	0168	10001 00 0	204 200 2	482	0288	1310-73-2	215-185-5	369
0010	7784-31-8	233-135-0	56	0168	10034-99-8	231-298-2	246	0293	7681-82-5	231-679-3	377
0011	631-61-8	211-162-9	57	0169	10034-99-8	231-298-2	246	0294	7681-82-5	231-679-3	377
0013	1066-33-7	213-911-5	62	0170	1313-13-9	215-202-6	248	0295	10102-40-6	231-551-7	378
0014	1341-49-7	215-676-4	62	0172	598-62-9	209-942-9	248	0296	7631-99-4	231-554-3	378
0015	10361-29-2	233-786-0	58	0173	13446-34-9	231-869-6	248	0297	7632-00-0	231-555-9	379
0017	10361-29-2	233-786-0	58	0174	10034-96-5	232-089-9	249	0298	7632-00-0	231-555-9	379
0018	12125-02-9	235-186-4	59	0176	7782-91-4	231-970-5	271	0301	62-76-0	200-550-3	379
0019	12125-02-9	235-186-4	59	0177	79-11-8	201-178-4	123	0302	6834-92-0	229-912-9	378
0020	3012-65-5	221-146-3	62	0179	373-02-4	206-761-7	275	0303	10049-21-5	231-449-2	363
0021	7789-09-5	232-143-1	59	0181	3333-67-3	222-068-2	275	0304	10039-32-4	231-448-7	367
0023	12125-01-8	235-185-9	61	0182	7791-20-0	231-743-0	275	0305	7782-85-6	231-448-7	367
0024	12054-85-2	234-722-4	61	0183	7791-20-0	231-743-0	276	0306	7558-79-4	231-448-7	366
0025	6484-52-2	229-347-8	69	0184	13478-00-7	236-068-5	276	0307	10101-89-0	231-509-8	381
0028	7727-54-0	231-786-5	69	0185	10101-97-0	232-104-9	277	0308	10101-89-0	231-509-8	381
0029	7722-76-1	231-764-5	61	0187	6153-56-6	205-634-3	285	0309	13472-36-1	231-767-1	362
0031	7783-28-0	231-987-8	62	0188	108-95-2	203-632-7	299	0311	7727-73-3	231-820-9	383
0032	7783-20-2	231-984-1	70	0190	37267-86-0	253-433-4	301	0312	7757-82-6	231-820-9	383
0034	1762-95-4	217-175-6	71	0192	1314-56-3	215-236-1	304	0313	7757-82-6	231-820-9	383
0035	10034-81-8	233-108-3	72	0194	85-44-9	201-607-5	304	0316	7757-83-7	231-821-4	384
0037	10025-91-9	233-047-2	75	0195	127-08-2	204-822-2	308	0317	6106-24-7	212-773-3	384
0045	10326-27-9	233-788-1	87	0197	877-24-7	212-889-4	319	0318	540-72-7	208-754-4	385
0052	10022-31-8	233-020-5	87	0198	7646-93-7	231-594-1	320	0319	10102-17-7	231-867-5	385
0053	7727-43-7	231-784-4	87	0199	16731-55-8	240-795-3	316	0320	7772-98-7	231-867-5	385
0055	10043-35-3	233-139-2	92	0202	7758-02-3	231-830-3	309	0321	10213-10-2	236-743-4	386
0059	62-54-4	200-540-9	111	0204	584-08-7	209-529-3	310	0322	10026-06-9	231-588-9	416
0061	471-34-1	207-439-9	113	0205	584-08-7	209-529-3	310	0323	10025-69-1	231-868-0	415
0064	10035-04-8	233-140-8	114	0207	3811-04-9	223-289-7	310	0324-01	1343-98-2	215-683-2	347
0070	10033-04-8	233-140-8	113	0207	7447-40-7	231-211-8	311	0325	10025-69-1	231-868-0	415
0073	1305-62-0	215-137-3	115	0209	7447-40-7	231-211-8	311	0326	10025-05-1	231-448-7	367
0073		215-137-3	115	0209	7789-00-6	232-140-5	312	0327	9005-25-8	232-679-6	389
	1305-62-0					232-140-5					
0076	13477-34-4	233-332-1	116	0211	7789-00-6		312	0329	57-11-4	200-313-4	390
0078	1305-78-8	215-138-9	116	0213	151-50-8	205-792-3	312	0331	10025-70-4	233-971-6	391
0081	10101-41-4	231-900-3	117	0214	7778-50-9	231-906-6	313	0333	110-15-6	203-740-4	391
0083	1308-38-9	215-160-9	128	0215	7778-50-9	231-906-6	313	0334	57-50-1	200-334-9	393
0086	1333-82-0 5949-29-1	215-607-8 201-069-1	128 129	0216 0218	7778-50-9 7758-02-3	231-906-6 231-830-3	313 309	0335 0339	7704-34-9 14807-96-6	231-722-6 238-877-9	394 404
0090	77-92-9	201-069-1	129	0219	13746-66-2	237-323-3	317	0340	87-69-4	201-766-0	405
0091	10026-22-9	233-402-1	133	0220	14459-95-1	237-722-2	317	0342	62-56-6	200-543-5	413
0093	10026-24-1	233-334-2	133	0222	1310-58-3	215-181-3	320	0344	76-03-9	200-927-2	421
0095	6046-93-1	205-553-3	135	0224	1310-58-3	215-181-3	321	0345	57-13-6	200-315-5	431
0097	10125-13-0	231-210-2	136	0225	7789-23-3	232-151-5	316	0346	57-13-6	200-315-5	431
0099	3251-23-8	221-838-5	137	0226	7758-05-6	231-831-9	325	0351	7440-66-6	231-175-3	439
0101	1317-38-0	215-269-1	138	0227	7681-11-0	231-659-4	325	0352	7440-66-6	231-175-3	439
0102	1317-38-0	215-269-1	137	0228	7681-11-0	231-659-4	326	0354	7440-66-6	231-175-3	439
0104	7758-99-8	231-847-6	138	0230	7681-11-0	231-659-4	326	0356	7440-66-6	231-175-3	439
0105	7758-99-8	231-847-6	139	0231	7757-79-1	231-818-8	326	0357	5970-45-6	209-170-2	440
0107	7758-98-7	231-847-6	138	0232	7757-79-1	231-818-8	327	0359	7646-85-7	231-592-0	441
0108	7758-89-6	231-842-9	136	0233	7758-09-0	231-832-4	327	0360	7646-85-7	231-592-0	442
0112			144	0235	6487-48-5	209-506-8	327	0362	10196-18-6	231-943-8	442
0113	50-99-7	200-075-1	185	0237	7722-64-7	231-760-3	328	0363	1314-13-2	215-222-5	443
0115	50-99-7	200-075-1	185	0239	7727-21-1	231-781-8	330	0364	1314-13-2	215-222-5	443
0116			234	0240	7778-77-0	231-913-4	316	0365	7446-20-0	231-793-3	443
0118	7783-83-7	233-382-4	68	0241	7758-11-4	231-834-5	318	0366	7487-88-9	231-298-2	246
0119	10025-77-1	231-729-4	229	0242	7758-11-4	231-834-5	318	0367	1330-43-4	215-540-4	161
0121	7782-61-8	233-899-5	230	0243	7778-53-2	231-907-1	330	0373	7758-02-3	231-830-3	309
0122	1309-37-1	215-168-2	230	0245	7790-62-7	232-216-8	331	0377	7727-21-1	231-781-8	329
0123	15244-10-7	233-072-9	230	0246	6381-59-5	206-156-8	331	0378	7790-62-7	232-216-8	331
0124	7783-85-9	233-151-8	68	0247	7778-80-5	231-915-5	332	0382	13453-69-5	236-631-5	240
0125	13478-10-9	213-843-4	229	0250	333-20-0	206-370-1	332	0383	10099-58-8	233-237-5	236
0126	7782-63-0	231-753-5	230	0252			343	0385	1310-58-3	215-181-3	320
0127	7782-63-0	231-753-5	231	0253		215-683-2	346	0387	1310-58-3	215-181-3	321
0132	7553-56-2	231-442-4	225	0256	6131-90-4	204-823-8	354	0390	631-61-8	211-162-9	57
0134	7439-89-6	231-096-4	228	0258	127-09-3	204-823-8	353	0391	10361-29-2	233-786-0	58
0135	7439-89-6	231-096-4	228	0259	127-09-3	204-823-8	353	0392	7722-76-1	231-764-5	60
0138	10039-26-6	200-559-2	235	0263	144-55-8	205-633-8	366	0393	6131-90-4	204-823-8	354
0141	6080-56-4	206-104-4	238	0264	10034-88-5	240-778-0	367	0394	144-55-8	205-633-8	365
0142	1335-32-6	215-630-3	237	0266	7631-90-5	231-548-0	367	0402	1310-73-2	215-185-5	368
0145	7758-95-4	231-845-5	238	0267	7681-57-4	231-673-0	363	0403	1310-73-2	215-185-5	368
0148	10099-74-8	233-245-9	238	0268	1303-96-4	215-540-4	161	0404	1310-73-2	215-185-5	369
0154	1305-78-8	215-138-9	116	0270	1303-96-4	215-540-4	161	0405	1310-73-2	215-185-5	369
0155	554-13-2	209-062-5	241	0271	7647-15-6	231-599-9	356	0458		215-185-5	376
0156	7447-41-8	231-212-3	241	0274	497-19-8	207-838-8	358	0501	10043-35-3	233-139-2	93
0157	7447-41-8	231-212-3	241	0277	7647-14-5	231-598-3	360	0503	471-34-1	207-439-9	113
0158	7790-69-4	232-218-9	242	0278	7647-14-5	231-598-3	360	0504	10035-04-8	233-140-8	114
0160	16674-78-5	205-554-9	243	0280	6132-04-3	200-675-3	361	0505	108-95-2	203-632-7	299
0161	12125-28-9	235-192-7	244	0281	143-33-9	205-599-4	361	0506	1314-56-3	215-236-1	188
				0201	5 55 5						.00



PROD. NO.	CAS NO.	EINECS	PAGE	PROD. NO.	CAS NO.	EINECS	PAGE	PROD. NO.	CAS NO.	EINECS	PAGE
0509	7447-40-7	231-211-8	312	1198	5965-83-3	202-555-6	394	1828	5329-14-6	226-218-8	393
0512	7778-80-5	231-915-5	332	1199	1401-55-4	215-753-2	404	1848	1344-28-1	215-691-6	55
0515	10125-13-0	231-210-2	137	1200	62-55-5	200-541-4	412	1850	1344-28-1	215-691-6	55
0517	6131-90-4	204-823-8	354	1206	125-20-2	204-729-7	414	1878-01	1317-39-1	215-270-7	137
0520	10025-69-1	231-868-0	415	1208-01	7790-78-5	233-296-7	109	1883	15708-41-5	239-802-2	165
0536	7778-18-9	231-900-3	117	1217	1314-62-1	215-239-8	432	1889	17927-65-0	233-135-0	56
0537	1344-28-1	215-691-6	55	1219	8049-22-7		435	1890	12054-85-2	234-722-4	61
0567	1343-88-0	215-681-1	177	1231	90-44-8	201-994-0	74	1892	9002-18-0	232-658-1	52
0597	10277-43-7	233-238-0	236	1240	9002-92-0	204.054.7	94	1899	7727-54-0	231-786-5	69
0604	9048-46-8	232-936-2	53	1255	7080-50-4	204-854-7	123	1914	50-81-7	200-066-2	76
0619	632-99-5	211-189-6	146	1261	5808-22-0	204-972-9	129	1916	143-19-1	205-591-0	379
0620	150-90-3	205-778-7	382	1262	573-58-0	209-358-4	134	1926			339
0676	37267-86-0	253-433-4	301	1269	13291-61-7	236-308-9	145	1927			339
0707 0709	85-86-9	201-638-4	393 393	1270 1279	65-85-0	200-618-2 202-819-0	88	1930	79-14-1	201-180-5	187 187
0716	63-74-1 92-31-9	200-563-4 202-146-2	420	1293	100-10-7 119-26-6	204-309-3	156 159	1933 1955	56-40-6 7647-17-8	200-272-2 231-600-2	122
0810	32-31-3	202-140-2	270	1295	140-22-7	205-403-7	160	1971	6381-59-5	206-156-8	331
0902		272-489-0	119	1296	538-62-5	208-698-0	160	1973	583-52-8	209-506-8	328
0903	1310-73-2	215-185-5	368	1303	15086-94-9	239-138-3	165	1983			341
0908	10124-56-8	233-343-1	365	1314	518-47-8	208-253-0	177	1984			341
1002			52	1325	2437-29-8	219-441-7	246	1991	7440-44-0	231-153-3	51
1003	7429-90-5	231-072-3	53	1329	61-73-4	200-515-2	264	1994	497-19-8	207-838-8	358
1010	7803-55-6	232-261-3	68	1330	115-40-2	204-087-8	95	1995	77-92-9	201-069-1	129
1011	7773-06-0	231-871-7	70	1341	1465-25-4	215-981-2	274	1998	69-65-8	200-711-8	249
1018	50-81-7	200-066-2	76	1366	10450-60-9	233-937-0	295	2008	10294-42-5	237-029-5	120
1029	10035-06-0	233-791-8	91	1368	3829-86-5	223-325-1	298	2012	7647-17-8	231-600-2	122
1032	1304-76-3	215-134-7	91	1369	115-39-9	204-086-2	95	2016	7681-49-4	231-667-8	365
1039	10108-64-2	233-296-7	109	1388	13573-16-5	237-003-3	342	2017	7631-99-4	231-554-3	378
1043 1046	7790-84-3 7789-75-5	233-331-6 232-188-7	109 115	1389 1400	81-88-9 497-19-8	201-383-9 207-838-8	343 358	2018 2021	62-76-0 532-32-1	200-550-3 208-534-8	380 355
1050	16774-21-3	240-827-6	58	1414	77-86-1	201-064-4	428	2022	7647-15-6	231-599-9	356
1052	302-17-0	206-117-5	123	1420	3244-88-0	221-816-5	51	2024	497-19-8	207-838-8	358
1055 1057	10060-12-5 6147-53-1	233-038-3 200-755-8	128 132	1421 1442	1611-35-4 16921-30-5	216-553-8 240-979-3	438 317	2028 2036	54-21-7 7782-49-2	200-198-0 231-957-4	382 344
1059	7791-13-1	231-589-4	133	1462	78-50-2	201-121-3	427	2039	50-70-4	200-061-5	387
1065	548-62-9	208-953-6	139	1464	471-34-1	207-439-9	112	2054	87-69-4	201-766-0	405
1067	7789-45-9	232-167-2	136	1465	12007-60-2	234-514-3	242	2060	6381-92-6	205-358-3	163
1071	122-39-4	204-539-4	160	1468	18472-87-2	242-355-6	301	2064	5949-29-1	201-069-1	130
1073	6381-92-6	205-358-3	163	1471	130-22-3	204-981-8	53	2071	7775-14-6	231-890-0	363
1099			227	1472	34725-61-6	252-170-2	96	2120-02			468
1100			227	1473	6104-58-1	228-058-4	134	2121-20			468
1107	1787-61-7	217-250-3	166	1474	6104-59-2	228-060-5	134	2173	22767-50-6	245-210-5	193
1109	60-00-4	200-448-9	162	1476	2353-45-9		175	2175	2832-45-3	220-601-3	200
1110	1185-57-5	214-686-6	67	1480	108-95-2	203-632-7	300	2195	5470-11-1	226-798-2	221
1113	100-97-0	202-905-8	196	1504	56-40-6	200-272-2	187	2338	1317-36-8	215-267-0	239
1115 1116	123-31-9 5470-11-1	204-617-8 226-798-2	221 221	1511 1513	76-60-8 76-59-5	200-972-8 200-971-2	95 96	2363 2456	95-14-7 113-24-6	202-394-1 204-024-4	88 381
1118	7553-56-2	231-442-4	225	1514	58-08-2	200-371-2	110	2459	357-57-3	206-614-7	96
1120	7439-95-4	231-104-6	242	1521	5470-11-1	226-798-2	222	2662	1762-95-4	217-175-6	71
1122	1309-48-4	215-171-9	245	1533	77-86-1	201-064-4	428	2680	75-05-8	200-835-2	50
1124	141-82-2	205-503-0	247	1536	10294-26-5	233-653-7	351	2685	67-56-1	200-659-6	259
1125	69-65-8	200-711-8	249	1546	617-45-8	210-513-3	77	2690	7296-20-1	230-732-8	158
1128	1600-27-7	216-491-1	252	1547	67-52-7	200-658-0	85	2790	110-16-7	203-742-5	247
1130	9005-25-8	232-686-4	389	1563	57-48-7	200-333-3	183	2794	55-55-0	200-237-1	262
1131	7487-94-7	231-299-8	252	1565	58-86-6	200-400-7	438	2811	151-21-3	205-788-1	364
1133	7774-29-0	231-873-8	252	1569			270	2815	7791-07-3	231-511-9	380
1134	10045-94-0	233-152-3	253	1576	52-89-1	200-157-7	142	2818	5324-84-5	226-195-4	283
1136	21908-53-2	244-654-7	253	1603	71-00-1	200-745-3	201	2838	9002-93-1		428
1140	7783-35-9	231-992-5	254	1613		203-473-3 203-473-3	306	2841	22767-49-3	245-208-4	290
1145	547-58-0	208-925-3	268	1615		203-473-3	306	2842	5574-97-0	226-947-1	407
1146	493-52-7	207-776-1	269	1616		203-473-3	306	2843	5574-97-0	218-147-6	407
1149	3051-09-0	221-266-6	273	1691	1662-01-7	216-767-1	161	2846	32503-27-8	251-068-5	406
1150	90-15-3	201-969-4	273	1701	7778-18-9 7778-18-9	231-900-3	116	2867-01			77
1153 1156	1313-99-1 85-85-8	215-215-7 201-637-9	276 286	1702 1703	7778-18-9	231-900-3 231-900-3	117 117	2868-01			78 78
1157	30525-89-4	200-001-8	288	1703	7440-50-8	231-159-6	134	2869-01 2871-01			77
1158	66-71-7	200-629-2	298	1737	7778-77-0	231-913-4	316	2873-01			77
1159	77-09-8	201-004-7	300	1747	288-32-4	206-019-2	224	2874-01			78
1165	12501-23-4	215-682-7	429	1758	6132-04-3	200-675-3	361	2875-01			78
1169	7790-21-8	232-196-0	332	1766	13472-36-1	231-767-1	362	2876-01			77
1171	87-66-1	201-762-9	341	1767	7681-57-4	231-673-0	363	2878-01			77
1176	1303-86-2	215-125-8	93	1768	13472-35-0	231-449-2	362	2879-01			77
1180	7783-90-6	232-033-3	348	1774	7727-73-3	231-820-9	383	2880-01			78
1182	7761-88-8	231-853-9	349	1777	7757-83-7	231-821-4	384	2881-01			78
1183	10294-26-5	233-653-7	350	1779	10102-17-7	231-867-5	385	2890-03	18497-13-7	241-010-7	127
1190	13755-38-9	238-373-9	380	1780	144-55-8	205-633-8	366	2940	298-14-6	206-059-0	317
1192 1194	7790-28-5 143-66-8	232-197-6 205-605-5	385 384	1783 1795	8006-28-8 10196-18-6	231-943-8	352 442	3058 3227	7722-64-7	231-760-3	224 328
1194	10213-10-2	236-743-4	386	1796	7446-20-0	231-793-3	442	3321	124-68-5	204-709-8	328 56
1197	121-57-3	204-482-5	394	1827	7-10-20-0	215-683-2	346	3351	4264-83-9	224-246-5	282



PROD. NO.	CAS NO.	EINECS	PAGE	PROD. NO.	CAS NO.	EINECS	PAGE	PROD. NO.	CAS NO.	EINECS	PAGE
3369			177	4145	75621-03-3		122	4657	7647-01-0	231-595-7	213
3371	100-97-0	202-905-8	194	4147	3076-05-9	221-352-3	122	4660	7553-56-2	231-441-9	227
3375	7757-82-6	231-820-9	382	4149	0070 00 0	LL: OOL O	75	4662	7553-56-2	231-441-9	227
3376	68855-54-9	272-489-0	120	4153	75277-39-3		191	4665		205-634-3	285
3377-00			482	4163	7119-22-7		272	4667	7758-01-2	231-829-8	309
3377-00	7757-82-6	231-820-9	382	4168	50-01-1	200-002-3	189	4671	7778-50-9	231-906-6	314
3400	67-56-1	200-659-6	260	4217	7732-18-5	231-791-2	435	4673	1310-58-3	215-181-3	323
3401	67-63-0	200-661-7	337	4218	7732-18-5	231-791-2	435	4674	1310-58-3	215-181-3	323
3402	100-51-6	202-859-9	89	4219	7732-18-5	231-791-2	433	4676	7758-05-6	231-831-9	325
3403	67-64-1	200-662-2	46	4265	5625-37-6		304	4677	7722-64-7	231-760-3	329
3405	64-17-5	200-578-6	170	4266	100037-69-2		305	4681	7761-88-8	231-853-9	350
3406	64-17-5	200-578-6	168	4312	3486-35-9		441	4682	7761-88-8	231-853-9	350
3407	64-17-5	200-578-6	170	4390	3400-33-3		83	4684	7647-14-5	231-598-3	361
3408	64-17-5	200-578-6	170	4391			84	4687	1310-73-2	215-185-5	376
3410	1330-20-7	215-535-7	437	4392			84	4689	1310-73-2	215-185-5	375
3411	108-88-3	203-625-9	420	4393			84	4690	1310-73-2	215-185-5	375
3451			288	4394			84	4691	1310-73-2	215-185-5	375
3716			387	4395			84	4693	10102-17-7	231-867-5	386
3719			190	4396			85	4695	10102-17-7	231-867-5	386
3813	12627-53-1	235-732-1	239	4397			85	4699	7664-93-9	231-639-5	403
3814			249	4398			85	4700	7664-93-9	231-639-5	403
3815	51811-82-6	257-438-2	185	4399			82	4701	7664-93-9	231-639-5	403
3816	68988-92-1		436	4400			82	4704	7664-93-9	231-639-5	403
3855			250	4401			82	4706	64-19-7	200-580-7	44
3856			184	4402			82	4708			95
3858	50-00-0	200-001-8	178	4403			82	4712	7697-37-2	231-714-2	282
3859	50-00-0	200-001-8	177	4404			83	47 15	1310-73-2	215-185-5	376
3864			287	4405			83	4795			97
3865			287	4406			83	4796			98
3866			287	4407			83	4797			99
3869			121	4408			78	4806	7601-90-3	231-512-4	291
3870			190	4409			79	4807	1336-21-6	215-647-6	67
3871			165	4410			79	4850			96
3872			344	4411			79	4851			96
3873			190	4412			79	4852			97
3874 3876			165 346	4413 4414			79 80	4853 4854			97 97
3878			436	44 15			80	4855	108-10-1	203-550-1	266
3879			239	4416			80	4856	100-10-1	203-330-1	98
3880			94	4417			80	4857			99
3904	540-84-1	000 450 0	272	4418			81	4858			99
3905	90622-57-4	292-459-0	429	4419			81	4860			97
3921 3925	8002-74-2		429 430	4420 4421			81 80	4861 4862			98 99
3930	0002-74-2		144	4425	1936-15-8	217-705-6	284	4863			98
3933	50-00-0	200-001-8	179	4435	1330-13-0	217-700-0	387	4864			376
3934	30-00-0	231-673-0	178	4437			389	4865			375
4001	7447-40-7	231-211-8	311	4439			388	4866			323
4002	7447-41-8	231-212-3	241	4441			388	4867		231-595-7	214
4003	7791-18-6	232-094-6	244	4442			387	4868			282
4004	1132-61-2	214-478-5	272	4444			389	4869			404
4007	1239-45-8	214-984-6	171	4450			387	4870		205-358-3	165
4008	7778-77-0	231-913-4	315	4455			388	4871		200 000 0	164
4009	6131-90-4	204-823-8	353	4456			388	4873			329
4011	10049-21-5	231-449-2	362	4458			389	4874	7446-20-0		444
4012	7758-11-4	231-834-5	318	4460			388	4876		231-659-4	225
4014	4432-31-9	224-632-3	271	4505			468	4878			98
4018	7365-45-9	230-907-9	190	4509			387	4889	877-24-7	212-889-4	319
4027	7783-20-2	231-984-1	70	4510			388	4890	57-50-1	200-334-9	392
4028	75-12-7	200-842-0	180	4511			389	4893	50-99-7	200-075-1	145
4031	110-26-9	203-750-9	264	4512			388	4921	7778-77-0	231-913-4	315
4035	10043-35-3	233-139-2	92	4514			468	4922	471-34-1	207-439-9	112
4040	6381-92-6	205-358-3	163	4515			468	4923	497-19-8	207-838-8	357
4042	7647-17-8	231-600-2	121	4520			468	4924	7647-14-5	231-598-3	359
4043	56-81-5	200-289-5	186	4528			469	4931	7722-76-1	231-764-5	60
4045	50-01-1	200-002-3	189	4529			345	4962	13463-67-7	236-675-5	417
4049	60-24-2	200-464-6	250	4531			345	4991	1312-81-8	215-200-5	236
4056	108-95-2	203-632-7	299	4532			345	5003			346
4058	7647-14-5	231-598-3	360	4533			162	5004			347
4059	56-40-6	200-272-2	186	4535			345	5006			55
4062	7558-79-4	231-448-7	366	4555-02			388	5007			120
4075	50-01-1	200-002-3	189	4581			468	5012			120
4081	79-06-1	201-173-7	51	4586			468	5014			306
4095	151-21-3	205-788-1	364	4601			227	5015			306
4097	57-50-1	200-334-9	392	4631	1700 05 4	047 475 0	227	5016			120
4098	110-18-9 1185-53-1	203-744-6 201-064-4	412 428	4650 4652	1762-95-4	217-175-6 205-358-3	71 165	5036 5046			120 306
4103	77-86-1	201-064-4	428	4653	139-33-3 139-33-3	205-358-3	164	5056	1343-98-2	215-683-2	346
4111	57-13-6	200-315-5	431	4654	7647-01-0	231-595-7	214	5155	7722-84-1	231-765-0	219
4118	1135-40-6	214-492-1	118	4655	7647-01-0	231-595-7	214	5168	10043-35-3	233-139-2	92
5					, J. J			0.00			



PROD. NO.	CAS NO.	EINECS PAGE	PROD. NO.	CAS NO.	EINECS	PAGE	PROD. NO.	CAS NO.	EINECS	PAGE
5651		104	5779			348	6813			237
5652		104	5780			352	6814			240
5653		102	5781			390	6815			243
5654		101	5782			405	6816			247
5655		103	5783			406	6817			251
5656		102	5784			412	6818			270
5657		100	5785			413	6819			274
5701		53	5786			414	6820			307
5703		74	5787			416	6821			348
5704		75	5788			430	6822			352
0,01		, 0	0,00			.00	0022			002
5705		85	5789			431	6823			390
5706		89	5790			438	6824			414
5707		90	5791			440	6825			416
5708		93	5792			444	6826			432
5709		108	5793			247	6827			440
5710		110	6003-01			269	6900	7647-01-0	231-595-7	208
5711		127	6004	108-24-7	203-564-8	44	6901	7697-37-2	231-714-2	278
5712		132	6004-01			269	6902	7664-93-9	231-639-5	395
5713		135	6005	1336-21-6	215-647-6	66	6903	64-19-7	200-580-7	42
5714		183	6005-01			270	6904	7664-39-3	231-634-8	215
5716		54	6008	64-18-6	200-579-1	182	6906	7732-18-5	231-791-2	434
5717		75	6010	10035-10-6	233-113-0	204	6908	7664-38-2	231-633-2	302
5718		76	6011	7647-01-0	231-595-7	206	6917			54
5719		86	6011-01	7647.04.0	221 FOF 7	224	6919			76
5720		90	6012	7647-01-0	231-595-7	207	6920			86
5721		91 93	6012-01	7664-39-3	221 624 0	224	6921 6922			89 90
5722 5723		108	6013 6014-01	7004-39-3	231-634-8	217 225	6923			93
5724		111	6016	6303-21-5	228-601-5	222	6924			108
5727		128	6017	79-33-4	201-196-2	234	6925			110
3727		120	0017	70-00-4	201-100-2	254	0020			110
5728		132	6019	7697-37-2	231-714-2	280	6926			127
5729		135	6021-01	7007 07 2	2017142	333	6927			132
5730		188	6022	7601-90-3	231-512-4	292	6928			135
5731		229	6023	7601-90-3	231-512-4	294	6929			228
5732		237	6024	7664-38-2	231-633-2	303	6930			237
5733		240	6027	7664-93-9	231-639-5	397	6931			240
5734		243	6031-01			420	6932			243
5735		248	6033-01			420	6933			247
5736		251	6034	79-09-4	201-176-3	338	6934			251
5737		271	6035	112-80-1	204-007-1	284	6935			270
5738		275	6037	64-18-6	200-579-1	181	6936			274
5739		286	6050	7782-99-2	045 047 0	404	6938			345
5740		305	6051	1336-21-6	215-647-6	66	6939			347
5741		307	6052	64-19-7	200-580-7	200	6940			348
5742 5743		344 345	6057 6060	7664-93-9 79-33-4	231-639-5 201-196-2	398 235	6941 6942			352 390
5744		347	6064	111-19-3	203-843-4	344	6943			414
5745		348	6069	598-82-3	201-196-2	234	6945			432
5746		352	6070	7647-01-0	231-595-7	209	6946			440
5747		390	6079	7664-39-3	231-634-8	218	6947			235
07 17			0070	7001000	20.00.0	2.0	0017			200
5748		405	6080	7697-37-2	231-714-2	281	6948			235
5749		406	6081	7647-01-0	231-595-7	205	6949			236
5750		413	6103-01			223	6950			307
5751		415	6106-01			119	6951			305
5752		416	6119-01			223	6952			188
5753		430	6122-01			435	6953			286
5754		432	6123	7699-45-8	231-718-4	441	6954			413
5755		439	6124			226	6964			416
5756		440	6125	1336-21-6	215-647-6	63	6968	28300-74-5		74
5757		445	6131	75-59-2	200-882-9	411	7000-03			490
5750		404	0447	7004.00.0	004 000 5	000	7000.04			100
5758		184	6147	7664-93-9	231-639-5	398	7000-04			490
5759		184	6152	64-19-7	200-580-7	43	7001-00			490
5760		277	6162	1336-21-6	215-647-6	63	7001-03			490
5761 5762		412 184	6163	7664-93-9 1310-73-2	231-639-5	396	7001-04			490 490
5762 5763		188	6164 6165	1310-73-2	215-185-5 215-185-5	371 370	7002-05 7003	12135-76-1	235-223-4	71
5764		228	6166	64-18-6	200-579-1	181	7003-04	12 130-70-1	230-223-4	490
5765		237	6167	7647-01-0	231-595-7	211	7003-04	62-53-3	200-539-3	73
5766		237	6801	, 047-01-0	201-000-7	54	7004-04	02-33-3	200-000-0	490
5767		243	6802			74	7004-04	62-53-3	200-539-3	73
3707		210	0002				7,000	32 30 0		,,,
5768		251	6803			76	7005-04			490
5769		270	6804			86	7007-04			490
5770		274	6805			90	7009-04			490
5771		277	6806			91	7010	100-51-6	202-859-9	89
5772		286	6807			108	7010-04			490
5773		305	6808			110	7011-03			490
5774		307	6809			128	7012-04			490
5776		343	6810			132	7013-04			490
5777		345	6811			135	7016	56-23-5	200-262-8	408
5778		347	6812			228	7018	67-66-3	200-663-8	125



PAGE	EINECS	CAS NO.	PROD. NO.	PAGE	EINECS	CAS NO.	PROD. NO.	PAGE	EINECS	CAS NO.	PROD. NO.
402	231-639-5	7664-93-9	7103	482			7061-00	468			7018-00
489	201 000 0	700+000	7103-00	177	215-681-1	1343-88-0	7061-00	468			7018-94
484			7104-00	321	215-181-3	1310-58-3	7063	126	200-663-8	67-66-3	7019
484			7104-01	338	200-338-0	57-55-6	7064	464			7020-00
484			7104-39	343	201-961-0	90-02-8	7066	464			7020-01
373	215-185-5	1310-73-2	7105	370	215-185-5	1310-73-2	7067	464			7020-02
484			7105-01	482			7068-01	464			7020-03
408	200-262-8	56-23-5	7106	406	201-191-5	79-27-6	7069	464			7020-06
194	213-059-4	920-66-1	7108	481			7069-01	464			7020-07
370	215-185-5	1310-73-2	7109	481			7070-01	464			7020-08
480 402	231-639-5	7664-93-9	7109-00 7110	421 422	204-800-2 201-167-4	126-73-8 79-01-6	7074 7075	464 464			7020-11 7020-13
212	231-595-7	7647-01-0	7111	482	201-107-4	73-01-0	7076-01	464			7020-13
480	201 000 7	7017 01 0	7111-00	424	200-929-3	76-05-1	7079	464			7020-22
373	215-185-5	1310-73-2	7112	144	202-046-9	91-17-8	7080	464			7020-23
323			7113	426	203-613-3	108-75-8	7081	464			7020-26
489			7113-00	488			7084-05	464			020-27
121	237-029-5	10294-42-5	7114	139	202-423-8	95-48-7	7085	464			7020-33
483 483			7114-00 7114-06	462 462			7086-00 7086-01	464 464			020-40 020-41
114	233-140-8	10043-52-4	7116	462			7086-02	464			020-42
483	200-140-0	10045-02-4	7116-00	462			7086-02	464			020-42
483			7116-06	462			7086-06	464			7021-00
483			7116-39	462			7086-07	464			7021-02
111	200-540-9		7117	462			7086-08	464			7021-03
488			7118-01	462			7086-22	464			7021-07
488			7118-04	462			7086-23	488			7022-01
488	004 744 0	7007 07 0	7118-06	462			7086-26	488			7022-02
281 469	231-714-2	7697-37-2	7119 7119-01	462 462			7086-28 7086-40	488 488			7022-04 7022-06
469			7120-03	462			7086-41	147	201-557-4	84-74-2	7024
324			7120-03	462			7086-41	481	201-007-4	O4-/4-Z	7024-00
469			7121-01	465			7087-00	481			7024-01
469			7121-03	465			7087-01	481			024-02
469			7121-04	465			7087-02	481			7024-05
469			7121-05	465			7087-03	481			7024-09
469			7121-06	465			7087-06	147	202-425-9	95-50-1	7025
469	047 475 0	1700 05 4	7121-08	465	004 505 7	7047.04.0	7087-26	482			7025-00
71 469	217-175-6	1762-95-4	7122 7122-00	212 463	231-595-7	7647-01-0	7088 7088-00	482 151	203-868-0	111-42-2	7025-01 7026
95			7123	463			7088-01	482			7026-00
163	205-358-3		7123	463			7088-02	481			7020-00
164	205-358-3	139-33-3	7125	463			7088-03	155	203-305-9	105-53-3	7028
164			7126	463			7088-09	481			7028-00
164			7127	463			7088-13	481			7028-01
226	231-441-9		7128	462			7089-03	156	204-826-4	127-19-5	7030
483			7128-02	463			7090-00	157	200-679-5	68-12-2	7032
483 483			7128-04 7128-08	463 463			7090-01 7090-03	158 372	200-664-3 215-185-5	67-68-5 1310-73-2	7033 7036
322	215-181-3	1310-58-3	7120-06	463			7090-03	175	203-473-3	107-21-1	7030
486			7129-01	463			7090-29	481			7037-00
322	215-181-3	1310-58-3	7129-01	462			7091-00	213	231-595-7	7647-01-0	7037-00
324	5 .51 0		7131	462			7091-01	178	200-001-8	50-00-0	7040
324			7132	462			7091-03	482			7040-00
486			7132-01	163	231-791-2		7092	178	200-001-8	50-00-0	7041
									200-842-0	75-12-7	7042
308	231-829-8		7133	488			7092-00	180	200 0 12 0		7042-00
486	231-829-8		7133 7133-01	488			7092-04	481	200 012 0		
486 486		7770 50 0	7133 7133-01 7133-06	488 488	200 024 0	07.00.5	7092-04 7092-06	481 481		FC 04 F	7043-00
486	231-829-8	7778-50-9	7133 7133-01	488	200-664-3	67-68-5	7092-04	481	200-289-5	56-81-5	7043-00 7044
486 486 314 486	231-906-6		7133 7133-01 7133-06 7134 7134-01	488 488 158 464	200-664-3	67-68-5	7092-04 7092-06 7093 7094-00	481 481 186 481		56-81-5	7043-00 7044 7044-00
486 486 314		7778-50-9 7783-35-9	7133 7133-01 7133-06 7134	488 488 158	200-664-3	67-68-5	7092-04 7092-06 7093	481 481 186		56-81-5	7043-00 7044 7044-00 7045
486 486 314 486 254 267 485	231-906-6		7133 7133-01 7133-06 7134 7134-01 7135 7136 7136-01	488 488 158 464 464 464 464			7092-04 7092-06 7093 7094-00 7094-01 7094-03 7094-06	481 481 186 481 325 481 481	200-289-5		7043-00 7044 7044-00 7045 7045-00 7046-00
486 486 314 486 254 267 485 269	231-906-6		7133 7133-01 7133-06 7134-01 7134-01 7135 7136 7136-01 7137	488 488 158 464 464 464 464 349	200-664-3	67-68-5 7761-88-8	7092-04 7092-06 7093 7094-00 7094-01 7094-03 7094-06 7095	481 481 186 481 325 481 481 220		56-81-5 7722-84-1	7043-00 7044 7044-00 7045 7045-00 7046-00 7047
486 486 314 486 254 267 485 269 486	231-906-6 231-992-5		7133 7133-01 7133-06 7134-01 7134-01 7135 7136-01 7137-01	488 488 158 464 464 464 464 349 465			7092-04 7092-06 7093 7094-00 7094-01 7094-03 7094-06 7095 7095-00	481 481 186 481 325 481 481 220 481	200-289-5		7043-00 7044 7044-00 7045 7045-00 7046-00 7047 7047-00
486 486 314 486 254 267 485 269 486 356	231-906-6 231-992-5 207-838-8		7133 7133-01 7133-06 7134 7134-01 7135 7136 7136-01 7137-01 7137	488 488 158 464 464 464 464 464 349 465 465			7092-04 7092-06 7093 7094-00 7094-01 7094-03 7094-06 7095 7095-00 7095-01	481 481 186 481 325 481 481 220 481 481	200-289-5	7722-84-1	7043-00 7044 7044-00 7045 7045-00 7046-00 7047 7047-00 7047-01
486 486 314 486 254 267 485 269 486 356 357	231-906-6 231-992-5 207-838-8 207-838-8		7133 7133-01 7133-06 7134 7134-01 7135 7136 7136-01 7137 7137-01 7138 7139	488 488 158 464 464 464 464 464 349 465 465			7092-04 7092-06 7093 7094-00 7094-01 7094-03 7094-06 7095-00 7095-01 7095-03	481 481 186 481 325 481 481 220 481 481 218	200-289-5		7043-00 7044 7044-00 7045 7045-00 7046-00 7047 7047-00 7047-01 7048
486 486 314 486 254 267 485 269 486 356 357 361	231-906-6 231-992-5 207-838-8 207-838-8 231-598-3		7133 7133-01 7133-06 7134 7134-01 7135 7136 7136-01 7137 7137-01 7138 7139 7140	488 488 158 464 464 464 464 349 465 465 465			7092-04 7092-06 7093 7094-00 7094-01 7094-03 7094-06 7095-00 7095-01 7095-03 7095-05	481 481 186 481 325 481 481 220 481 481 218	200-289-5	7722-84-1	7043-00 7044 7044-00 7045 7045-00 7046-00 7047 7047-01 7048 7048-01
486 486 314 486 254 267 485 269 486 356 357	231-906-6 231-992-5 207-838-8 207-838-8		7133 7133-01 7133-06 7134 7134-01 7135 7136 7136-01 7137 7137-01 7138 7139	488 488 158 464 464 464 464 464 349 465 465			7092-04 7092-06 7093 7094-00 7094-01 7094-03 7094-06 7095-00 7095-01 7095-03	481 481 186 481 325 481 481 220 481 481 218	200-289-5	7722-84-1	7043-00 7044 7044-00 7045-00 7045-00 7046-00 7047-7047-01 7048 7048-01 7049-01
486 486 314 486 254 267 485 269 486 356 357 361 285 294	231-906-6 231-992-5 207-838-8 207-838-8 231-598-3		7133 7133-01 7133-06 7134 7134-01 7135 7136 7136-01 7137 7137-01 7138 7139 7140 7141 7142	488 488 158 464 464 464 464 465 465 465 465 465 349	231-853-9	7761-88-8	7092-04 7092-06 7093 7094-00 7094-01 7094-03 7094-06 7095-00 7095-01 7095-01 7095-03 7095-09 7096-09	481 481 186 481 325 481 481 220 481 481 481 481 482 482	200-289-5	7722-84-1	7043-00 7044 7044-00 7045-00 7045-00 7046-00 7047-01 7047-01 7048-01 7048-01 7049-01
486 486 314 486 254 267 485 269 486 356 351 285 294	231-906-6 231-992-5 207-838-8 207-838-8 231-598-3		7133 7133-01 7133-06 7134 7134-01 7135 7136 7136-01 7137 7137-01 7138 7139 7140 7141 7142	488 488 158 464 464 464 464 465 465 465 465 465	231-853-9	7761-88-8	7092-04 7092-06 7093 7094-00 7094-01 7094-03 7094-06 7095-00 7095-01 7095-03 7096-06 7095-09	481 481 186 481 325 481 481 220 481 481 218 481 481	200-289-5	7722-84-1	043-00 7044 0044-00 7045 045-00 046-00 7047 047-00 047-01 7048 048-01 050-01
486 486 314 486 254 267 485 269 486 356 357 361 285 294 300 100 103	231-906-6 231-992-5 207-838-8 207-838-8 231-598-3		7133 7133-01 7133-06 7134 7134-01 7135 7136 7136-01 7137 7137-01 7138 7139 7140 7141 7142 7143 7144 7145	488 488 488 158 464 464 464 464 465 465 465 465	231-853-9 231-853-9 215-185-5	7761-88-8 7761-88-8 1310-73-2	7092-04 7092-06 7093 7094-00 7094-01 7094-03 7094-06 7095-00 7095-01 7096-03 7096-09 7096-00 7096-00 7097-01	481 481 186 481 325 481 481 220 481 481 481 482 482 482	200-289-5 231-765-0 231-765-0 200-838-9	7722-84-1 7722-84-1 75-09-2	7043-00 7044-00 7044-00 7045-00 7046-00 7047-00 7047-01 7048-01 7049-01 7050-01 7051-01 7052-01 7053
486 486 314 486 254 267 485 269 486 357 361 285 294 300 100 103 351	231-906-6 231-992-5 207-838-8 207-838-8 231-598-3		7133 7133-01 7133-06 7134 7134-01 7135 7136 7136 7136-01 7137 7137-01 7138 7139 7140 7141 7142 7143 7144 7145 7147	488 488 488 158 464 464 464 464 464 465 465 465	231-853-9	7761-88-8 7761-88-8	7092-04 7092-06 7093 7094-01 7094-01 7094-03 7094-06 7095-00 7095-01 7095-03 7096-00 7096-00 7097-01 7097-01 7098	481 481 186 481 325 481 481 220 481 481 481 482 482 482 481 481 481 50 60 60 60 60 60 60 60 60 60 6	200-289-5 231-765-0 231-765-0	7722-84-1 7722-84-1	7043-00 7044 7044-00 7045-00 7045-00 7047-00 7047-01 7048-01 7049-01 7050-01 7051-01 7053 7055
486 486 314 486 254 485 269 486 356 357 361 285 294 300 100 103 351 351	231-906-6 231-992-5 207-838-8 207-838-8 231-598-3		7133 7133-01 7133-06 7134 7134-01 7135 7136 7136-01 7137 7137-01 7138 7139 7140 7141 7142 7143 7144 7145 7147 7148	488 488 488 158 464 464 464 464 465 465 465 465	231-853-9 231-853-9 215-185-5	7761-88-8 7761-88-8 1310-73-2	7092-04 7092-06 7093 7094-01 7094-01 7094-03 7094-06 7095-00 7095-01 7096-03 7096-00 7096-00 7097 7097 7097-01 7098 7098-00	481 481 186 481 325 481 220 481 218 481 481 482 482 482 481 150 171 481	200-289-5 231-765-0 231-765-0 200-838-9 205-483-3	7722-84-1 7722-84-1 75-09-2 141-43-5	7043-00 7044 7044-00 7045-7045-00 7046-00 7047-01 7047-01 7048-01 7048-01 7050-01 7051-01 7052-01 7053 7055-01
486 486 3144 486 254 485 269 486 356 357 361 285 294 300 103 351 351 351	231-906-6 231-992-5 207-838-8 207-838-8 231-598-3 205-634-3	7783-35-9	7133 7133-01 7133-06 7134 7134-01 7135 7136 7136-01 7137 7137-01 7138 7139 7140 7141 7142 7143 7144 7145 7147 7148 7149	488 488 488 464 464 464 464 465 465 465 465	231-853-9 231-853-9 215-185-5 215-185-5	7761-88-8 7761-88-8 1310-73-2 1310-73-2	7092-04 7092-06 7093 7094-01 7094-01 7094-03 7094-06 7095-00 7095-01 7096-03 7096-09 7096-00 7097-01 7098-01 7098-01	481 481 186 481 325 481 481 220 481 481 482 482 482 481 150 171 481 329	200-289-5 231-765-0 231-765-0 200-838-9	7722-84-1 7722-84-1 75-09-2	7043-00 7044 7044-00 7045 7045-00 7046-00 7046-00 7047-01 7047-01 7049-01 7050-01 7051-01 7052-01 7053 7065 7055-01 7057
486 486 314 486 254 267 485 269 486 356 357 294 300 100 351 351 351 351 351	231-906-6 231-992-5 207-838-8 207-838-8 231-598-3 205-634-3	7783-35-9 67-66-3	7133 7133-01 7133-06 7134 7134-01 7135 7136 7136 7136 7137 7137-01 7138 7139 7140 7141 7142 7143 7144 7145 7147 7148 7149 7152	488 488 488 464 464 464 464 464	231-853-9 231-853-9 215-185-5 215-185-5	7761-88-8 7761-88-8 1310-73-2 1310-73-2	7092-04 7092-06 7093 7094-01 7094-01 7094-03 7094-06 7095-00 7095-01 7095-03 7096-00 7096-00 7097-01 7098 7098-00 7098-01 7098-01 7098-01 7099-01	481 481 186 481 325 481 481 220 481 218 481 218 481 482 482 482 481 481 150 171 481 329 481	200-289-5 231-765-0 231-765-0 200-838-9 205-483-3	7722-84-1 7722-84-1 75-09-2 141-43-5	7043-00 7044 7044-00 7045 7045-00 7046-00 7047-01 7047-01 7048-01 7049-01 7050-01 7051-01 7055-01 7057-01
486 486 3144 486 254 485 269 486 356 357 361 285 294 300 103 351 351 351	231-906-6 231-992-5 207-838-8 207-838-8 231-598-3 205-634-3	7783-35-9	7133 7133-01 7133-06 7134 7134-01 7135 7136 7136-01 7137 7137-01 7138 7139 7140 7141 7142 7143 7144 7145 7147 7148 7149	488 488 488 464 464 464 464 465 465 465 465	231-853-9 231-853-9 215-185-5 215-185-5	7761-88-8 7761-88-8 1310-73-2 1310-73-2	7092-04 7092-06 7093 7094-01 7094-01 7094-03 7094-06 7095-00 7095-01 7096-03 7096-09 7096-00 7097-01 7098-01 7098-01	481 481 186 481 325 481 481 220 481 481 482 482 482 481 150 171 481 329	200-289-5 231-765-0 231-765-0 200-838-9 205-483-3	7722-84-1 7722-84-1 75-09-2 141-43-5	7043-00 7044 7044-00 7045-00 7046-00 7047-01 7047-01 7048-01 7049-01 7050-01 7052-01 7053 7055 7050-07



PROD. NO.	CAS NO.	EINECS	PAGE	PROD. NO.	CAS NO.	EINECS	PAGE	PROD. NO.	CAS NO.	EINECS	PAGE
7157-00			485	7216	1310-73-2	215-185-5	374	7279-00			483
7157-00			485	7216-06	1310-73-2	2 10- 100-0	463	7280-00			483
7157-05			485	7217-06			463	7280-43			483
7158	56-81-5	200-289-5	186	7218	7733-02-0	231-793-3	444	7283-00			486
7158-05			483	7218-02			463	7283-02			486
7158-06			483	7219	7647-01-0	231-595-7	212	7283-05			486
7 159			100	7219-07			463	7285-00			485
7159-00			483	7221			294	7285-02			485
7159-05	0040.05.4	200 201 2	483	7221-03			466	7285-05			485
7160	8012-95-1	232-384-2	287	7222-06			466	7291-01			485
7160-05			480	7225-04			466	7296-39			484
7161-05			480	7225-04			466	7297-00			484
7165	127-18-4	204-825-9	407	7225-06			466	7297-43			484
7165-00			489	7225-24			466	7297-47			484
7165-01			489	7227			148	7298-39			483
7165-02			489	7230	7647-01-0	231-595-7	212	7305	75-09-2	200-838-9	150
7165-03			489	7232	10043-52-4	233-140-8	114	7306	56-23-5	200-262-8	408
7165-04			489	7233			115	7309	71-55-6	200-756-3	422
7165-26			489	7238-06			463	7309-00			489
7165-30			489	7239-09			463	7314-01			487
7169	57-55-6	200-338-0	339	7241-00			468	7314-02			487
7176	25322-68-3	203-473-3	306	7242	7664-93-9	231-639-5	399	7315-01			487
7178	25322-68-3	203-473-3	306	7243-00	7004-33-3	231-000-0	482	7315-02			487
7179	25322-69-4	200-338-0	307	7244			145	7323-00			468
7179-00	LOOLL GO !	200 000 0	486	7245	7761-88-8	231-853-9	350	7327-00			468
7179-02			486	7246	7761-88-8	231-853-9	349	7328-01			467
7179-05			486	7247	7722-64-7	231-760-3	329	7328-03			467
7180-00			486	7248	10326-27-9		86	7328-06			467
7180-02			486	7248-00			486	7329-01			467
7180-05			486	7248-02			486	7329-03			467
7181-00			485	7248-05			486	7329-06			467
7181-00			485	7246-05			314	7329-00	14634-91-4	238-676-6	176
7181-02			485	7250			313	7334-01	14004-31-4	230-070-0	464
7182-00			485	7251			139	7334-03			464
7182-00			485	7251-00			486	7334-04			464
7182-05			485	7251-00			486	7334-06			464
7183-02			486	7251-05			486	7334-07			464
7184-02			485	7252			253	7334-08			464
7189-01			465	7252-00			485	7337-01			462
7189-02			465	7252-02			485	7337-07			462
7400.00			105	7050.05			405	7007.00			400
7189-03			465	7252-05			485	7337-08			462
7189-04			465 465	7253 7254	7440 00 0	221 702 2	354 444	7340 7341			81
7189-06 7191-00			486	7254-00	7446-20-0	231-793-3	485	7348			81 342
7191-00			486	7254-02			485	7354	75-25-2	200-854-6	95
7191-05			486	7254-05			485	7356	75-25-2	200-004-0	74
7192-00			486	7255	121-54-0	204-479-9	201	7368	106-93-4	203-444-5	146
7192-02			486	7261			99	7374	9005-64-5		429
7192-05			486	7262			100	7385		200-001-8	179
7194	544-76-3	208-878-9	194	7263			101	7386	67-66-3	200-663-8	126
7199-00			481	7263-00			485	7390	50-00-0	200-001-8	179
7199-00			481	7263-00			485	7391	72-17-3	200-001-8	377
7201	1310-73-2	215-185-5	373	7263-05			485	7394	9005-65-6		429
7202	1310-73-2	215-185-5	371	7264			101	7397			119
7203	1310-58-3	215-181-3	322	7264-00			486	7400	68-12-2	200-679-5	157
7206			324	7264-01			486	7400-00			489
7207-00			486	7264-02			486	7400-01			489
7207-02			486	7264-05			486	7401-00			489
7207-05			486	7265			102	7401-01			489
7208-00			468	7267			104	7401-29			489
7210		231-441-9	226	7268			104	7401-30			489
7210	7647-01-0	231-595-7	213	7269			104	7402-00			489
7211-00			463	7269-00			485	7404-07			489
7211-01			463	7269-02			485	7404-08			489
7211-03			463	7269-05			485	7404-37			489
7212	7647-01-0	231-595-7	212	7270			105	7404-38			489
7212-00			464	7272			130	7413	865-49-6	212-742-4	127
7212-01			464	7272-00			483	7418-06			465
7212-03			464	7272-43			483	7418-07			465
7213-03			462	7273			103	7419-06			464
7213-07			462	7273-00			464	7419-07			464
7213-08			462	7273-01			464	7420	56-23-5	200-262-8	409
7213-23			462	7273-02			464	7420-06			462
7213-27	1015	045 :	462	7273-03			464	7420-07			462
7214	1310-73-2	215-185-5	373	7274-00			483	7421-00			468
7214-03			462	7275-00			483	7423-00	E0 00 0	200 001 0	468
7214-04			462	7276-00			483	7424	50-00-0	200-001-8	179
7214-07			462	7276-06			483	7424-00			468
7214-27 7215	7664-93-9	231-639-5	462 402	7277-00 7278-00			483 483	7425-00 7425-01			468 468
/210	7004-33-3	201-000-0	402	1210-00			400	/420-UT			400



PROD. NO.	CAS NO.	EINECS	PAGE	PROD. NO.	CAS NO.	EINECS	PAGE	PROD. NO.	CAS NO.	EINECS	PAGE
7426-00			468	7519-02			465	8049	108-10-1	203-550-1	267
7427-00			468	7519-02			465	8050	108-87-2	203-624-3	264
7429-00			468	7519-03			465	8052	78-93-3	201-159-0	265
										201-159-0	
7430-00			468	7519-25			465	8053	78-93-3	201-159-0	265
7432-00			468	7523-00			465	8055-06			477
7433-00			468	7523-02			465	8055-07			477
7435-00			468	7523-05			465	8056	78-83-1	201-148-0	232
7436-00			468	7523-22			465	8056-06			477
7437-00			468	7524-04			466	8057-06			477
7439	108-32-7	203-572-1	338	7525	877-24-7	212-889-4	319	8058-06			477
7441-00			480	7526			304	8060-06			477
7441-01			480	7530-00			482	8061-06			477
7441-02			480	7534			342	8066	71-23-8	200-746-9	333
7441-03			480	7535			342	8067	67-63-0	200-661-7	336
7441-08			480	7541			97	8068	67-63-0	200-661-7	336
7442-05			480	7557	1310-73-2	215-185-5	372	8068-06	07 00 0	200 001 7	477
7442-03			480	7567-00	1310-73-2	2 10- 100-0	485				477
								8068-07			
7446-01			480	7567-02			485	8070-01	400.00.0	000 500 0	477
7448-01			480	7567-05			485	8072	108-20-3	203-560-6	155
7461			101	7568-01			485	8072-06			477
7462			102	7568-02			485	8072-07			477
7463			226	7575-06			462	8073	110-86-1	203-809-9	340
7464		231-992-5	254	7575-07			462	8074	110-86-1	203-809-9	341
7465-00			482	7577	7553-56-2	231-441-9	226	8075	109-99-9	203-726-8	411
7465-01			482	7598	1310-73-2	215-185-5	371	8077	108-88-3	203-625-9	419
7465-02			482	7606-01	1010 70 2	210 100 0	473	8078	108-88-3	203-625-9	419
7466-00			465	7606-02			473	8079	121-44-8	204-469-4	424
7466-01			465	7606-03			473	8080	1330-20-7	215-535-7	436
7466-03 7466-04			465 465	7606-04 7606-06			473 473	8083 8084	110-80-5 109-86-4	203-804-1 203-713-7	171 262
7466-06			465	7606-07			473	8093-01			477
7466-07			465	7606-08			473	8094-06			477
7466-08			465	7606-09			473	8095-06			477
7466-10			465	7606-11			473	8096-02			478
7466-12			465	7606-31			475	8097-06			478
7466-22			465	7609-01			465	8098	64-17-5	200-578-6	167
7466-27			465	7609-02			465	8098-01	04-17-5	200-370-0	478
7400-27		231-639-5	399		7664-93-9	231-639-5	399	8099-06			477
		231-039-0		7610							
7472-02 7473-00			485 485	7613 7620	139-33-3 1310-73-2	205-358-3 215-185-5	164 371	8100-06 8101-01			477 478
7473-02			485	7621	1310-73-2	215-185-5	372	8102-01			478
7473-05			485	7622	7697-37-2	231-714-2	282	8102-04			478
7474-06			483	7624	7664-93-9	231-639-5	399	8103	78-92-2	201-158-5	107
7475	1310-73-2	215-185-5	374	7631			314	8104	110-82-7	203-806-2	141
7475-39			483	7638	7664-93-9	231-639-5	402	8108-00			472
7476	7779-88-6	231-943-8	442	7639	1310-58-3	215-181-3	322	8108-01			472
7477	7647-01-0	231-595-7	213	7641			351	8108-02			472
7478	7758-98-7	231-847-6	176	7642	7783-35-9		254	8108-03			472
7479	6381-59-5	205-698-2	176	7645	10326-27-9		86	8108-04			472
7480	1310-73-2	215-185-5	374	7650-07	10020 27 0		466	8108-06			472
7481		231-791-2	164	7651	7647-01-0	231-595-7	213	8108-07			472
		231-791-2									
7486	7004.00.0	004 400 0	310	7655	7647-01-0	231-595-7	211	8108-08			472
7488	7664-93-9	231-163-9	342	7704-06			466	8108-09			472
7490-07			466	8001	67-64-1	200-662-2	45	8108-11			472
7490-08			466	8002	67-64-1	200-662-2	46	8108-31			475
7491			301	8004	75-05-8	200-835-2	50	8109-00			471
7493	7647-01-0	231-595-7	213	8006	64-17-5	200-578-6	167	8109-01			471
7495-04			466	8007	64-17-5	200-578-6	169	8109-02			471
7495-18			466	8009	71-41-0	200-752-1	72	8109-03			471
7496-00			489	8010	123-51-3	204-633-5	231	8109-04			471
7496-01			489	8012	75-85-4	200-908-9	72	8109-06			471
7496-25			489	8013	100-66-3	202-876-1	73	8109-00			471
7496-28					71-43-2	200-753-7	88				
			489	8014				8109-08			471
7496-29			489	8017	71-36-3	200-751-6	106	8109-09			471
7496-30			489	8018	123-86-4	204-658-1	106	8109-11			471
7504-00			480	8019	75-65-0	200-889-7	107	8109-31			475
7505-00			480	8023	108-90-7	203-628-5	124	8110	142-82-5	205-563-8	192
7506-00			480	8026	110-82-7	203-806-2	140	8111	142-82-5	205-563-8	193
7511-04			466	8028	108-94-1	203-806-2	141	8111-00			472
7511-06			466	8030	109-89-7	203-716-3	152	8111-01			472
7512-00			468	8031	123-91-1	204-661-8	160	8111-02			472
7512-00			468	8033	60-29-7	200-467-2	152	8111-03			472
7514-00			468	8037	141-78-6	205-500-4	173	8111-04			472
7515-00			468	8038	141-78-6	205-500-4	174	8111-06			472
7516-00			468	8042	107-06-2	203-458-1	148	8111-07			472
						2U3-4U8- I					
7516-01			468	8043	64742-49-0	202 777 0	200	8111-08			472
7517-00			468	8044	110-54-3	203-777-6	199	8111-09			472
7518-08			466	8045	67-56-1	200-659-6	259	8111-11			472
7519-00			465	8046	67-56-1	200-659-6	257	8111-31			475
7519-01			465	8047	67-56-1	200-659-6	260	8113	78-78-4	201-142-8	263



	PROD. NO.	CAS NO.	EINECS	PAGE	PROD. NO.	CAS NO.	EINECS	PAGE	PROD. NO.	CAS NO.	EINECS	PAGE
_	8113-00			472	8154-04			473	8170-11			470
	8113-01			472	8154-06			473	8170-31			475
	8113-02			472	8154-07			473	8172	112-40-3	203-967-9	162
	8113-03			472	8154-08			473	8172-00			471
	8113-04 8113-06			472 472	8154-09 8154-11			473 473	8172-01 8172-02			471 471
	8113-07			472	8154-31			475	8172-03			471
	8113-08			472	8156-00			474	8172-06			471
	8113-09			472	8156-01			474	8172-07			471
	8113-11			472	8156-02			474	8172-08			471
	8113-31	109-66-0	202 002 4	475	8156-03			474	8172-09			471
	8114 8115	8032-32-4	203-692-4 232-453-7	290 296	8156-04 8156-07			474 474	8172-11 8172-31			471 475
	8115-00	0002 02 4	202 400 7	472	8156-11			474	8174-00			474
	8115-01			472	8156-31			476	8174-04			474
	8115-02			472	8157-02			473	8174-06			474
	8115-03			472	8157-03			473	8174-07			474
	8115-04			472	8157-07			473	8174-08	07.00.0	200 001 7	474
	8115-06 8115-07			472 472	8157-08 8157-09			473 473	8175 8175-00	67-63-0	200-661-7	334 474
	8115-08			472	8157-31	111 65 0	202 002 1	475	8175-01			474
	8115-09 8115-11			472 472	8159 8159-00	111-65-9	203-892-1	283 471	8175-02 8175-03			474 474
	8115-31			472	8159-03			471	8175-04			474
	8116	8032-32-4	232-453-7	297	8159-07			471	8175-06			474
	8117	109-99-9	203-726-8	411	8159-09			471	8175-07			474
	8118	1330-20-7	215-535-7	437	8160-00			474	8175-08			474
	8118-00			479	8160-02			474	8175-09			474
	8119 8119-01	67-63-0	200-661-7	336 479	8160-04 8160-06			474 474	8175-11 8175-31			474 475
	8120-01			479	8160-07			474	8176-31			475
	8121-01 8122-01			479 479	8160-09 8160-31			474 476	8178-31 8180-31			476 475
	8123	93-58-3	202-259-7	262	8163-00			470	8182-31			475
	8123-01	00 00 0	202 200 7	479	8163-01			470	8184-00			473
	8124-01			479	8163-02			470	8184-01			473
	8125-01			479	8163-03			470	8184-02			473
	8126-01			479	8163-04			470	8184-03			473
	8127-01 8128-01			479 479	8163-06 8163-07			470 470	8184-04 8184-06			473 473
	0.20 0.			.,,,	0.000 07				0.0.0			
	8129-00			479	8163-08			470	8184-07			473
	8130-01			479	8163-09			470	8184-08			473
	8131-96	100-42-5	202 0E1 E	479	8163-11			470 476	8184-09			473
	8133 8134	75-05-8	202-851-5 200-835-2	391 49	8163-31 8165-00			470	8184-11 8184-31			473 475
	8137-01	73-03-0	200-033-2	469	8165-01			470	8188-96			479
	8137-10			469	8165-02			470	8192	629-59-4	211-096-0	409
	8138-01			469	8165-03			470	8196-00			471
	8138-10			469	8165-06			470	8196-01			471
	8142	67-64-1	200-662-2	45	8165-07			470	8196-02			471
	8144	75-05-8	200-835-2	50	8165-08			470	8196-03			471
	8146	108-38-3	203-576-3	437	8165-09			470	8196-04			471
	8147	95-47-6	202-422-2	438	8165-11			470	8196-06			471
	8150-00			479	8165-31			475 471	8196-07			471 471
	8151-00 8151-01			473 473	8167-00 8167-02			471	8196-08 8196-09			471
	8151-02			473	8167-03			471	8196-11			471
	8151-03			473	8167-07			471	8196-31			475
	8151-04			473	8167-11			471	8197-24			479
	8151-06			473	8168	124-18-5	204-686-4	144	8204	141-78-6	205-500-4	173
	8151-07			473	8168-00			470	8205	110-54-3	203-777-6	197
	8151-08			473	8168-01			470	8210	78-83-1	201-148-0	232
	8151-09			473	8168-02			470	8212	75-65-0	200-889-7	107
	8151-11			473	8168-03			470	8217	540-84-1	208-759-1	426
	8151-31			475	8168-04			470	8224	108-10-1	203-550-1	267
	8153-00 8153-01			474 474	8168-06 8168-07			470 470	8228 8229	64-17-5 64-17-5	200-578-6 200-578-6	168 169
	8153-02			474	8168-08			470	8239	109-66-0	203-692-4	288
	8153-03			474	8168-09			470	8240	8032-32-4	232-453-7	296
	8153-04			474	8168-11			470	8241	8032-32-4	232-453-7	297
	8153-06			474	8168-31			476	8242	8032-32-4	232-453-7	297
	8153-07			474	8170-00			470	8243	8032-32-4	232-453-7	298
	8153-08			474	8170-01			470	8254	60-29-7	200-467-2	154
	8153-09			474 474	8170-02			470 470	8257 8271	75-05-8	200-835-2 214-300-6	49 430
	8153-11 8153-31			474	8170-03 8170-04			470	8271	1120-21-4 111-84-2	203-913-4	283
	8154-00			478	8170-04			470	8287	109-66-0	203-692-4	289
	8154-01			473	8170-07			470	8387	64742-96-7	265-200-4	234
	8154-02			473	8170-08			470	8401	60-29-7	200-467-2	154
	8154-03			473	8170-09			470	8402	67-56-1	200-659-6	256



PROD. NO.	CAS NO.	EINECS	PAGE	PROD. NO.	CAS NO.	EINECS	PAGE	PROD. NO.	CAS NO.	EINECS	PAGE
8404	67-56-1	200-659-6	257	9262	110-54-3	203-777-6	197				
8409	107-15-3	203-468-6	174	9263	67-56-1	200-659-6	255				
8414	111-76-2	203-905-0	105	9264	75-09-2	200-838-9	149				
8417	112-34-5	203-961-6	105	9265	8032-32-4	232-453-7	295				
8430 8434	120-82-1 123-91-1	204-428-0 204-661-8	422 159	9266 9267	8032-32-4 64742-49-0	232-453-7	295 200				
8462	64-17-5	200-578-6	166	9270	8032-32-4	232-453-7	296				
8564	111-46-6	203-872-2	286	9272	8032-32-4	232-453-7	295				
8578	111-92-2	203-921-8	146	9276	141-78-6	205-500-4	173				
8599	77-73-6	201-052-9	151	9277	110-54-3	203-777-6	197				
0000	140.00 5	005 500 0	100	0000	144 70 0	005 500 4	170				
8662 8668	142-82-5 110-54-3	205-563-8 203-777-6	193 196	9282 9292	141-78-6 110-82-7	205-500-4 203-806-2	172 140				
8669	110-54-3	203-777-6	199	9295	75-09-2	200-838-9	149				
8685	109-66-0	203-692-4	290	9304	110-54-3	203-777-6	198				
8689	112-27-6	203-953-2	424	9305	64742-49-0		199				
8691	78-30-8	201-103-5	423	9315	75-09-2	200-838-9	149				
8699	64-18-6	200-579-1	181	9316	75-09-2	200-838-9	150				
8706 8707	110-82-7	203-806-2	141	9331 9333	109-66-0	203-692-4	289				
8710	107-06-2 110-54-3	203-458-1 205-563-8	148 198	9334	109-66-0 67-63-0	203-692-4 200-661-7	289 334				
0/10	110-54-5	200-300-0	130	3334	07-03-0	200-001-7	334				
8714	108-88-3	203-625-9	418	9335	540-84-1	208-759-1	425				
8715	540-84-1	208-759-1	426	9336	108-88-3	203-625-9	418				
8782	766-09-6	212-161-6	175	9338	142-82-5	205-563-8	191				
8792	100.00.0	204 622 0	233	9344	68-12-2	200-679-5	157				
8819 8821	123-38-6 75-31-0	204-623-0 200-860-9	337 233	9345 9351	872-50-4 108-88-3	212-828-1 203-625-9	268 418				
8822	75-56-9	200-879-2	339	9360	127-18-4	204-825-9	407				
8834			233	9364	108-88-3	203-625-9	419				
8844			204	9365	142-82-5	205-563-8	192				
8845			204	9393	110-86-1	203-809-9	340				
8855			203	0420	100.00.0	202 726 0	409				
8860			203	9439 9440	109-99-9 109-99-9	203-726-8 203-726-8	410				
8861			203	9441	109-99-9	203-726-8	410				
8862			202	9444	120-82-1	204-428-0	421				
8863			202	9446	109-99-9	203-726-8	410				
8890			202	9448			143				
8891			201	9449			143				
8892 8898	67-56-1	200-659-6	202 203	9451 9470	76-05-1	200-929-3	143				
8899	07-50-1	200-059-0	203	9478	75-05-8	200-929-3	424 52				
9003	7726-95-6	231-778-1	94	9479			52				
9006	7439-97-6	231-106-7	250	9480	540-84-1	208-759-1	425				
9008	7723-14-0 75-05-8	231-768-7 200-835-2	303 48	9485 9488			118 285				
9015	1313-60-6	215-209-4	381	9497			117				
9017	75-05-8	200-835-2	48	9502	7429-90-5	231-072-3	53				
9019	7550-45-0	231-441-9	417	9503	7440-50-8	231-159-6	134				
9020	75-07-0	200-836-8	42	9504	7440-50-8	231-159-6	134				
9035	75-05-8	200-835-2	49	9509	7440-31-5	231-141-8	414				
9042	1634-04-4	216-653-1	263	9510			118				
9043	1634-04-4	216-653-1	263	9512			118				
9048-03	78-83-1	201-148-0	231	9513	8002-74-2	232-315-6	288				
9049	67-56-1	200-659-6	258	9516	1330-20-7	215-535-7	436				
9051	7705-08-0	231-729-4	229	9518			143				
9077	67-56-1	200-659-6	255 258	9519	64 10 7	200-580-7	285 43				
9091	67-56-1 7726-95-6	200-659-6 231-778-1	94	9524 9530	64-19-7 7647-01-0	231-595-7	207				
9097	67-56-1	200-659-6	257	9563	7664-39-3	231-634-8	216				
9098	16940-66-2	241-004-4	356	9598	7697-37-2	231-714-2	279				
9099	26628-22-8	247-852-1	355	9600	7647-01-0	231-595-7	210				
9111	121-44-8	204-469-4	423	9653	7601-90-3	231-512-4	293				
9132	76-13-1	200-936-1	423	9820	64-18-6	200-579-1	180				
9161	16940-66-2	241-004-4	355	9821	75-05-8	200-835-2	47				
9174	67-66-3	200-663-8	124	9822	67-56-1	200-659-6	256				
9175	67-66-3	200-663-8	125	9823	7732-18-5	231-791-2	434				
9177	142-82-5	205-563-8	191	9824			183				
9185	142-82-5	205-563-8	192	9826			182				
9213 9214	68-12-2 78-93-3	200-679-5 201-159-0	156 265	9829 9830	67-63-0	200-661-7	57 334				
9233	95-50-1	202-425-9	147	9831	141-78-6	205-500-4	172				
9234	67-68-5	200-664-3	158	9832	64-18-6	000 : -	182				
9237	60-29-7	200-467-2	152	A426	9012-36-6	232-731-8	52				
9250 9254	60-29-7 67-64-1	200-467-2 200-662-2	153 45	E343 F780	7440-44-0 27565-41-9	231-153-3	51 130				
9255	75-05-8	200-835-2	45	G045-01	110-83-8	203-807-8	142				
9257	67-66-3	200-663-8	124	L034-07	112-53-8	203-982-0	162				
9258	110-82-7	203-806-2	140	M834-09	107-22-2	203-474-9	187				
9259	60-29-7	200-467-2	153	CON-2			489				
9260	141-78-6	205-500-4	172 268								
9261	872-50-4	212-828-1	200				ı				